

Introductory Offer : 2011 ONLY . . . or until withdrawn

Shown Prices are actually USD List Prices

AUSTRALIAN LIST PRICE - SAME \$\$\$ Amount - NO Import Freight



THE NE-1000

THE PERFECT PUMP AT THE PERFECT PRICE

Feature-Loaded & Computer Controllable

The NE-1000 programmable syringe pump was developed by long-time experts in the industry. Our full line of pumps are among the most feature-loaded, stand-alone or (RS-232) computer controllable pumps on the market and cost significantly less.

The Perfect Price

Most syringe pumps on the market have retail prices of over \$2000. Most labs and companies, however need to obtain authorization to purchase expensive items (over \$1000). The NE-1000 is strategically priced at \$995 to make it easier for customers to purchase the pump through their credit cards.

2-Year Warranty & Unlimited Technical Support

After the sale of a pump, never worry about it again. If there are any problems, we offer an extensive two-year warranty on all parts. We also provide excellent technical support to tackle any questions for the life of the pump.

NEW ERA PUMP SYSTEMS



NE-9000 PROGRAMMABLE PERISTALTIC PUMP

New Era's newest edition. Built-in programmability. The *Learn and Repeat™* mode automates any small dispensing need from .5ml/hr to 1250ml/min

Model: NE-9000
Price: \$500



NE-300 "JUST INFUSION"™ SYRINGE PUMP

A basic single syringe pump for customers that just want infusion at an affordable price. Has adjustable rates and dispensed volume display.

Model: NE-300
Retail Price: \$275



NE-1000 SINGLE SYRINGE PUMP

Full featured single syringe programmable pump. Built-in easy to use programming features and an optional wide range of accessories, this syringe pump accepts one syringe up to 60cc. Operates stand-alone or can be computer controlled via RS-232.

Model: NE-1000
Retail Price: \$750

Other models include: NE-1010, NE-1002X, NE-1050



NE-4000 DOUBLE SYRINGE PUMP

Holds 2 syringes up to 60cc. Has all of the same full featured programmable capabilities as the NE-1000.


Model: NE-4000
Retail Price: \$975



NE-1600/NE-1800/1200 MULTI-SYRINGE PUMP

NE-1600 holds 2 or 4 syringes up to 140cc or 6 up to 60cc each
NE-1800 holds 2, 4, 6, or 8 syringes up to 10 cc
NE-1200 holds 2, 4, 6, 8, 10, 12 syringes up to 3cc each
Has all of the same full featured capabilities as the NE-1000.

Model: NE-1600/1800/1200
Retail Price: \$1500
Model: NE-1650/1850/1250 High Speed
Retail Price: \$1500



PUMP UPGRADES

X-Upgrade: Smooth linear increasing or decreasing gradient available on all NE-1000 family of pumps
X2-Upgrade: 340 programmable phases & X-upgrade

Retail Prices:
X Upgrade NE-1000 Series: \$105
X Upgrade NE-500 Series: \$55
X2 Upgrade NE-1000 Series: \$205
X2 Upgrade NE-500 Series: \$110

NEW ERA PUMP SYSTEMS



NE-500 OEM SYRINGE PUMP

The OEM design of the NE-500 is designed for easy system integration. With the addition of a starter kit, the NE-500 can function similar to the NE-1000.

Model: NE-500
Retail Price: \$495 (Starter kit: \$25)

Other Models include NE-501, NE-510, NE-511, NE-4500, NE-4501



DUAL NE-1000

Two NE-1000 pumps connected into one system where one pump controls the other. For continuous and simultaneous flow and emulsification applications.

Model: DUAL-NE-1000
Retail Price: \$1510

***Limited Inventory Price is enacted when 3 or more dual systems are purchased!!!**
****High Profit prices are enacted when 5 or more dual systems are purchased!!!**

Other Models include: Dualing any two same or different pumps in the NE-1000 family and the NE-500 series.



SYRINGE HEATER KIT

The syringe Heater Kit can heat liquids to over 100° C. Temperature sensor in heating pad. Multiple syringe heater pads can be attached to the control unit.

Model: Syringe-Heater-Kit
Retail Price: \$399



SYRINGE PUMP ACCESSORIES, SYRINGES & PLUMBING SUPPLIES

Accessories:

Syringes:

Plumbing Supplies:

Cables:

Starter Kits:

Syringe Kits:

Syringes : NormJect - always in stock Melb

Tubing, plastic fittings, Double Check Valves from ValuePlastics also in stock