

SRI Model 310MM Gas Chromatograph with Built-in Hydrogen generator and air compressor

May 2020

The SRI Model 310MM Gas Chromatograph (part# 0310-0091) is the perfect GC for many labs and individuals.

Measure Cannabinoids, Terpenes and Residual Solvents

Low Cost. A fraction of the cost of most other GCs.

Portable. You can fly with it as airline baggage.

Rugged. Use it in the lab or the field.

Fast. 10 minutes for a cannabinoid analysis. 16 minutes for terpenes and cannabinoids in a single test.

No gas cylinders required. The built-in Hydrogen generator and air compressor provide all required gases. Runs on just distilled water.

Powerful. The two columns let you separate gases or liquids without changing any hardware. The sensitive FID detector measures down to .01% THC.

Easy. Anybody can learn to operate this GC in a few hours. **You don't have to be a scientist.**

Software Included. And you can install the software on as many computers as you like. Updates are free.

Two Year Warranty.. One year on H2 cell.


Free Training and US based Tech Support.

Videos See the many you-tube videos. Click the icon at the bottom of srigc.com


To Order:

SRI part# 0310-0091

Model 310MM GC




12 x 14 x 15 inches
30 pounds



Built-in H2 generator

Terpenes and Cannabinoids in
a single analysis



SRI Tech Support: 310-214-5092

www.srigc.com

310MMFlyerMay2020

Page 1