



# Multiple Gas Analyzer #2 GC System

These chromatograms show the separation of a 1000ppm C<sub>1</sub>-C<sub>6</sub> standard plus room air on a MG#2 GC equipped with a TCD and FID detectors.



This customized MG#2 GC is equipped with:

- TCD Detector
- HID Detector
- FID-Methanizer
- Two Packed Columns
- One Capillary Column
- 10-port Gas Sampling Valve and Dual Loops
- 4 channel PeakSimple Data System
- ...on the compact 8610C chassis

**8610-0270 Multiple Gas Analyzer #2 GC with TCD detector and 1 channel PeakSimple data system**

**8610-0271 Multiple Gas Analyzer #2 GC with TCD, Methanizer, FID, 4 channel PeakSimple data system and built-in air compressor**

**8610-0272 Multiple Gas Analyzer #2 GC with TCD and HID detectors, and 4 channel PeakSimple data system**



OPTIONS & UPGRADES: 6 channel USB PeakSimple data system, capillary column, H<sub>2</sub>-50XR hydrogen generator, Thermal Desorber  
(VOLTAGE: for 110VAC, use "part number-1" [ex: 8610-0270-1] for 220VAC, use "part number-2")