

Model 8610C Gas Chromatograph

- *Mounts up to Six Detectors and Five Injectors*
- *Ambient to 400°C Temperature Programmable Column Ovens*
- *Dimensions: 19" wide x 13.5" high x 14.5" deep*
- *Implement virtually any EPA or ASTM method*

The Model 8610C Gas Chromatograph is our most versatile and popular model. While it is very compact next to comparable laboratory GCs from other manufacturers, it is large and flexible enough to perform an amazing variety of applications. See our Preconfigured GC section (starting on page 3) for examples of the 8610C adaptability.



Up to six detectors can be installed on the same 8610C GC.



ECD

PID

FID / DELCD

FPD