



# SunShell RP Guard Filter



<Cartridge Type, Bonded with C18 and End-Capped with TMS>


Available as a guard column for reversed phase



- ✓ The filter is made of porous glass sized 4 mm i.d. and 4 mm thickness.
- ✓ Pore diameter is 2  $\mu$ m.
- ✓ Low dead volume structure
- ✓ Back pressure on glass filter is ca. 0.1 MPa at 1.0 mL/min of flow rate.
- ✓ Upper pressure limit is more than 60 MPa
- ✓ Available for 2.1 mm i.d to 4.6 mm i.d. column



## Evaluation of SunShell RP Guard Filter



## Price of SunShell RP Guard Filter

Name	quantity	Part number	Photo
SunShell RP Guard Filter For exchange	5 pieces	CBGAAC	
SunShell RP Guard Filter Holder	1 piece	CBGAAH	