



# Syringe Filters (now with luer lock inlet)



[Simplify Sample Prep and Protect Analytical Equipment Using Syringe Filters](#) (103k PDF)

Looking for Syringes?

Restek offers the inert, non-reactive NORM-JECT® Polypropylene Syringes.

[NORM-JECT® Polypropylene Syringes](#)

- Variety of filter types, porosities, and diameters.
- Color coded for easy identification.
- Polypropylene housing.
- Reusable storage container.
- Quantity break pricing, for greater savings.

Price breaks now available!

Price per 100-pack. Price breaks are available at 5 and 10 packs.

Your correct price will be reflected on your invoice.

FREE sample packs available. Use these handy packs for method development or to compare with your current brand. Request yours today by adding -248 to the part number. Sample pack orders cannot be placed online—please call. Limit one sample pack per customer.

Description	Porosity	Diameter	Color	qty.	Cat.#
Syringe Filter, Cellulose Acetate	0.22µm	13mm	green	100-pk.	26156
Syringe Filter, Cellulose Acetate	0.45µm	13mm	blue	100-pk.	26155
Syringe Filter, Cellulose Acetate	0.22µm	25mm	green	100-pk.	26158
Syringe Filter, Cellulose Acetate	0.45µm	25mm	blue	100-pk.	26157
Syringe Filter, Nylon	0.22µm	13mm	yellow	100-pk.	26146
Syringe Filter, Nylon	0.45µm	13mm	pink	100-pk.	26147
Syringe Filter, Nylon	0.22µm	25mm	yellow	100-pk.	26148
Syringe Filter, Nylon	0.45µm	25mm	pink	100-pk.	26149
Syringe Filter, PTFE (polytetrafluoroethylene)	0.22µm	13mm	purple	100-pk.	26142
Syringe Filter, PTFE (polytetrafluoroethylene)	0.45µm	13mm	orange	100-pk.	26143
Syringe Filter, PTFE (polytetrafluoroethylene)	0.22µm	25mm	purple	100-pk.	26144
Syringe Filter, PTFE (polytetrafluoroethylene)	0.45µm	25mm	orange	100-pk.	26145
Syringe Filter, PVDF (polyvinylidifluoride)	0.22µm	13mm	brown	100-pk.	26150
Syringe Filter, PVDF (polyvinylidifluoride)	0.45µm	13mm	red	100-pk.	26151
Syringe Filter, PVDF (polyvinylidifluoride)	0.22µm	25mm	brown	100-pk.	26152
Syringe Filter, PVDF (polyvinylidifluoride)	0.45µm	25mm	red	100-pk.	26153

Cellulose Acetate, Nylon, PVDF - hydrophilic applications

PTFE - hydrophobic applications, resistant to acid and base, used mostly with nonaqueous samples

Click thumbnails to see enlarged view.



26146  
13mm Syringe Filter,  
0.22µm, Nylon, 100-pk.



26144  
25mm Syringe Filter,  
0.22µm, PTFE, 100-pk.



26153  
25mm Syringe Filter,  
0.45µm, PVDF, 100-pk.

Australian Distributors  
**HROMalytic**  
**ECH nology**  
 Chromatography Products '08  
 www.chromtech.net.au E-mail : info@chromtech.net.au Tel : +61 3 9762 2034 Fax : +61 3 9761 1169