














Treated Swagelok® and Valco® Fittings

Siltek®/Sulfinert® Treated and Silcosteel®-CR Treated Swagelok® Fittings

- Full line of treated 1/16", 1/8", and 1/4" fittings.
- Siltek®/Sulfinert® treatment ensures ultimate inertness; Silcosteel®-CR treatment enhances corrosion resistance by 10X, or more.
- Custom treatment available for any Swagelok® fitting, or other system parts.

Fitting Type	Size	Swagelok® #	Siltek®/Sulfinert® Treated*		Silcosteel®-CR Treated	
			qty.	cat.#	qty.	cat.#
 Union	1/16"	SS-100-6	ea.	22540	ea.	22575
	1/8"	SS-200-6	ea.	22541	ea.	22576
	1/4"	SS-400-6	ea.	22542	ea.	22577
	3/8"	SS-600-6	ea.	22909 new!	ea.	22904 new!
 Tee	1/16"	SS-100-3	ea.	22543	ea.	22578
	1/8"	SS-200-3	ea.	22544	ea.	22579
	1/4"	SS-400-3	ea.	22545	ea.	22580
	3/8"	SS-600-3	ea.	22910 new!	ea.	22905 new!
 Reducing Union	1/8" to 1/16"	SS-200-6-1	ea.	22546	ea.	22581
	1/4" to 1/16"	SS-400-6-1	ea.	22547	ea.	22582
	1/4" to 1/8"	SS-400-6-2	ea.	22548	ea.	22583
	3/8" to 1/4"	SS-600-6-4	ea.	22911 new!	ea.	22906 new!
 Elbow	1/8"	SS-200-9	ea.	22549	ea.	22584
	1/4"	SS-400-9	ea.	22550	ea.	22585
 Plug	1/16"	SS-100-P	ea.	22572	ea.	22619
	1/8"	SS-200-P	ea.	22573	ea.	22620
	1/4"	SS-400-P	ea.	22574	ea.	22597
 Cross	1/8"	SS-200-4	ea.	22551	ea.	22586
	1/4"	SS-400-4	ea.	22552	ea.	22587
 Tube End Reducer	1/8" tube to 1/16"	SS-100-R-2	ea.	22553	ea.	22588
	1/4" tube to 1/16"	SS-100-R-4	ea.	22554	ea.	22589
	1/8" tube to 1/4"	SS-400-R-2	ea.	22555	ea.	22590
	1/4" tube to 1/8"	SS-200-R-4	ea.	22556	ea.	22591
 Port Connector	1/8"	SS-201-PC	ea.	22557	ea.	22592
	1/4"	SS-401-PC	ea.	22558	ea.	22593
	1/8" tube to 1/4"	SS-401-PC-2	ea.	22559	ea.	22594
 Male Connector	1/8" to 1/8" NPT	SS-200-1-2	ea.	22561	ea.	22595
	1/4" to 1/4" NPT	SS-400-1-4	ea.	22562	ea.	22596
	1/16" to 1/8" NPT	SS-100-1-2	ea.	22563	ea.	22610
	1/4" to 1/4" NPT	SS-200-1-4	ea.	22564	ea.	22611
	1/4" to 1/8" NPT	SS-400-1-2	ea.	22565	ea.	22612
	3/8" to 3/8" NPT	SS-600-1-6	ea.	22912 > new!	ea.	22907 > new!
	3/8" to 1/4" NPT	SS-600-1-4	ea.	22913	ea.	22908
 Female Connector	1/8" to 1/8" NPT	SS-200-7-2	ea.	22566	ea.	22613
	1/4" to 1/4" NPT	SS-400-7-4	ea.	22567	ea.	22614
	1/4" to 1/8" NPT	SS-400-7-2	ea.	22568	ea.	22615
	1/8" to 1/4" NPT	SS-200-7-4	ea.	22569	ea.	22616
 Bulkhead Union	1/8"	SS-200-61	ea.	22570	ea.	22617
	1/4"	SS-400-61	ea.	22571	ea.	22618

Siltek®/Sulfinert® Treated and Silcosteel® Treated Valco® Fittings

Fitting Type	Size	Siltek®/Sulfinert® Treated*		Silcosteel® Treated	
		qty.	cat.#	qty.	cat.#
 Valco® Zero Dead Volume Tee	1/16"	ea.	22534	ea.	20581
	1/8"	ea.	22535	ea.	20583
 Valco® Zero Dead Volume Union	1/16"	ea.	22532	ea.	20580
	1/8"	ea.	22533	ea.	20582

*See Frequently Asked Questions on page 384

Siltek®/Sulfinert® Treated, Silcosteel® Treated, and Silcosteel®-CR Treated Parker Fittings

A full line of 1/16", 1/8" and 1/4" fittings is available with Siltek®/Sulfinert®, Silcosteel®, or Silcosteel®-CR treatment. Because of expanding applications for these coatings, we have received many requests for a broader product offering. If you still do not see everything you need, contact us for information on custom coating services.

Fitting Type	Size	Parker #	Siltek®/Sulfinert® Treated*		Silcosteel® Treated		Silcosteel®-CR Treated	
			qty.	cat.#	qty.	cat.#	qty.	cat.#
	1/16"	1 SC 1	ea.	22520	ea.	20510	ea.	22863
	1/8"	2 SC 2	ea.	22521	ea.	20511	ea.	22864
	1/4"	4 SC 4	ea.	22522	ea.	20512	ea.	22865
	1/16"	1 ET 1	ea.	22526	ea.	20513	ea.	22866
	1/8"	2 ET 2	ea.	22527	ea.	20514	ea.	22867
	1/4"	4 ET 4	ea.	22528	ea.	20515	ea.	22868
	1/8" to 1/16"	2 RU 1	ea.	22523	ea.	20519	ea.	22869
	1/4" to 1/16"	4 RU 1	ea.	22524	ea.	20520	ea.	22870
	1/4" to 1/8"	4 RU 2	ea.	22525	ea.	20521	ea.	22871
	1/16"	1 EE 1	—	—	—	—	ea.	22874
	1/8"	1 EE 2	ea.	22530	ea.	20517	ea.	22875
	1/4"	1 EE 4	ea.	22531	ea.	20518	ea.	22876
	1/16"	1 BLP 1	ea.	21539	ea.	21518	ea.	22877
	1/8"	2 BLP 2	ea.	21540	ea.	21519	ea.	22878
	1/4"	4 BLP 4	ea.	21541	ea.	21520	ea.	22879
	1/8"	2 ECR 2	ea.	21542	ea.	21521	ea.	22872
	1/4"	4 ECR 4	ea.	21543	ea.	21522	ea.	22873
	1/8" tube to 1/16"	2 TUR 1	ea.	21544	ea.	21523	ea.	22880
	1/4" tube to 1/16"	4 TUR 1	ea.	21545	ea.	21524	ea.	22881
	1/8" tube to 1/4"	2 TUR 4	ea.	21546	ea.	21525	ea.	22882
	1/4" tube to 1/8"	4 TUR 2	ea.	21547	ea.	21526	ea.	22883
	1/8"	2 PC 2	ea.	21548	ea.	21527	ea.	22884
	1/4"	4 PC 4	ea.	21549	ea.	21528	ea.	22885
	1/8" tube to 1/4"	2 PC 4	ea.	21550	ea.	21529	ea.	22886
	1/8" to 1/8" NPT	2 MSC 2N	ea.	21551	ea.	21530	ea.	22887
	1/4" to 1/4" NPT	4 MSC 4N	ea.	21552	ea.	21531	ea.	22888
	1/16" to 1/8" NPT	1 MSC 2N	ea.	21553	ea.	21532	ea.	22889
	1/8" to 1/4" NPT	2 MSC 4N	ea.	21554	ea.	21533	ea.	22890
	1/4" to 1/8" NPT	4 MSC 2N	ea.	21555	ea.	21534	ea.	22891
	1/8" to 1/8" NPT	2 FSC 2N	ea.	21556	ea.	21535	ea.	22892
	1/4" to 1/4" NPT	4 FSC 4N	ea.	21557	ea.	21536	ea.	22893
	1/4" to 1/8" NPT	4 FSC 2N	ea.	21558	ea.	21537	ea.	22894
	1/8" to 1/4" NPT	2 FSC 4N	ea.	21559	ea.	21538	ea.	22895
	1/8"	2A PR4 VT SS	ea.	21586	ea.	21576	—	—
	1/4"	4A PR4 VT SS	ea.	21587	ea.	21577	—	—
	1/8"	2A B2LJ2 SSP	ea.	21588	ea.	21578	—	—
	1/4"	4A B2LJ2 SSP	ea.	21589	ea.	21579	—	—

new!

For the newest developments, visit our Performance Coatings Division website: www.restekcoatings.com

*See Frequently Asked Questions on page 384.

Please note: Nuts and ferrules are not treated unless requested (custom parts). Nuts and ferrules normally are not in contact with samples, and thus do not require coating.