

VICI JOUR - Low Pressure Fittings

VICI JOUR PRODUCTS

Tubing for chromatography

High pressure fittings

- Nanovolume
- Stainless steel
- Polymeric

Low pressure fittings

Filters & mobile phase filters

Mobile phase handling systems

Frits

PEEK columns & guard columns

Sample loops

Accessories/miscellaneous

Low Pressure Bulkhead Unions

- Panel mounting
- Easy connection of low pressure solvent lines
- Low dead volume
- Completely inert and biocompatible



Low Pressure Bulkhead Unions

Pressure rating: up to 105 bar (1500 psi) for 1/16", 35 bar (500 psi) for 1/8"

Description	Product No.
PEEK, 0.75 mm bore, with ETFE ferrules and PPS flangeless nuts for 1/16" tubing	JR- CBUMPK
PEEK, 1.50 mm bore, with ETFE ferrules and PPS flangeless nuts for 1/8" tubing	JR-CBULPK

Spare Parts

Description	Qty/pkg	Product No.	
Nuts	PPS, flangeless, 1/16", black	10	JR-55070-10
	PPS, flangeless, 1/8", black	10	JR-55071-10
Ferrules	ETFE, 1/16"	10	JR-041-10
	ETFE, 1/8"	10	JR-051-10
Adapters	PP fingertight sleeve, yellow	10	JR-55086-10
	PP fingertight sleeve, blue	10	JR-55090-10