



VICI JOUR - Mobile Phase Handling Systems

VICI JOUR PRODUCTS

- Tubing for chromatography
- High pressure fittings
 - Nanovolume
 - Stainless steel
 - Polymeric
- Low pressure fittings
- Filters & mobile phase filters
- Mobile phase handling systems
- Frits
- PEEK columns & guard columns
- Sample loops
- Accessories/miscellaneous

CUSTOMER SUPPORT

Safety-coated Glass Bottles

- Special design for helium sparging of LC mobile phases
- Safety coating minimizes risk of injury



VICI Jour Safety-Coated Glass Bottles are available in 1 and 2 liter sizes. Schott Duran® glass bottles with GL45 screw threads are coated on the outside with a sturdy polymer layer. The chances of breaking the bottle by dropping it are decreased, and if a bottle should break, the glass fragments are kept inside the polymer coating. The coating will contain the solvent for a short while – long enough for safe disposal – and also protects solvents against UV light up to 385 nm.

Note:

The safety coating is not intended to replace safe working conditions and precautions in laboratories. The bottles are not for use with pressures above 0.8 bar (11.6 psi). Use only for continuous helium sparging.

Safety-coated Glass Bottles

Pressure rating: 0.8 bar (11.6 psi)

Description	Product No.
1 liter bottle	JR-9000-0005
2 liter bottle	JR-9000-0010

TOP ▲