

Striped color-coded PEEK tubing


SPECS

- Maximum recommended working temp: 100°C (continuous) for 1/16" OD tubing with ID up to 0.75 mm
- Max working pressures: see page 94
- Tubing tolerances: see page 95

We recommend striped PEEK tubing for easy color-coded identification of tubing ID, with maximum chemical resistance and biocompatibility. The stripe is added to the outside, so dye never contacts the fluid stream.

All colorants used in the manufacture of this tubing are RoHS-compliant. (Restriction of Hazardous Substances)

Qty/pkg:		1 m		3 m	
Tubing OD	Color	ID (µm)	Product No.	Product No.	Product No.
360 µm	Orange	50	JR-T-3650-M1	JR-T-3650-M3	
	Green	75	JR-T-3675-M1	JR-T-3675-M3	
	Grey	100	JR-T-36100-M1	JR-T-36100-M3	
	Red	125	JR-T-36125-M1	JR-T-36125-M3	
	Black	150	JR-T-36150-M1	JR-T-36150-M3	



1/32"	Color	ID (mm)	1 m Product No.	3 m Product No.
	Orange	0.050	JR-T-5702-M1	JR-T-5702-M3
Pink	0.064	JR-T-5764-M1	JR-T-5764-M3	
Green	0.075	JR-T-5703-M1	JR-T-5703-M3	
Black	0.100	JR-T-5704-M1	JR-T-5704-M3	
Red	0.13	JR-T-5705-M1	JR-T-5705-M3	
Yellow	0.18	JR-T-5707-M1	JR-T-5707-M3	
Orange	0.20	JR-T-5608-M1	JR-T-5608-M3	
Blue	0.25	JR-T-5710-M1	JR-T-5710-M3	
Orange	0.50	JR-T-5720-M1	JR-T-5720-M3	

1/16"	Color	ID (mm)	3 m Product No.	5m Product No.
	Orange	0.050	JR-T-5802-M3	JR-T-5802-M5
Pink	0.064	JR-T-5998C-M3	JR-T-5998C-M5	
Green	0.075	JR-T-5803-M3	JR-T-5803-M5	
Black	0.10	JR-T-5804-M3	JR-T-5804-M5	
Red	0.13	JR-T-5999-M3	JR-T-5999-M5	
Yellow	0.18	JR-T-6000-M3	JR-T-6000-M5	
Blue	0.25	JR-T-6001-M3	JR-T-6001-M5	
Orange	0.50	JR-T-6002-M3	JR-T-6002-M5	
Green	0.75	JR-T-6003-M3	JR-T-6003-M5	
Grey	1.00	JR-T-60031-M3	JR-T-60031-M5	
Black	1.40	JR-T-60032-M3	JR-T-60032-M5	

	Color	ID (mm)	10 m Product No.	25m Product No.
	Orange	0.050	JR-T-5802-M10	JR-T-5802-M25
Pink	0.064	JR-T-5998C-M10	JR-T-5998C-M25	
Green	0.075	JR-T-5803-M10	JR-T-5803-M25	
Black	0.100	JR-T-5804-M10	JR-T-5804-M25	
Red	0.13	JR-T-5999-M10	JR-T-5999-M25	
Yellow	0.18	JR-T-6000-M10	JR-T-6000-M25	
Blue	0.25	JR-T-6001-M10	JR-T-6001-M25	
Orange	0.50	JR-T-6002-M10	JR-T-6002-M25	
Green	0.75	JR-T-6003-M10	JR-T-6003-M25	
Grey	1.00	JR-T-60031-M10	JR-T-60031-M25	
Black	1.40	JR-T-60032-M10	JR-T-60032-M25	

TECH TIP

VICI Jour sample loops for analytical applications are made of PEEK striped color-coded tubing.
PEEK loops 81

OPTIONS

Contact your distributor or VICI about:
 ■ Longer bulk tubing lengths
 ■ Customized ODs and IDs (30 m minimum quantity)

MORE INFORMATION

About PEEK® 99
 Chemical resistance
 PEEK 96-98
 PEEK at higher temperatures 96
 Maximum working pressure 94
 Tubing ID chart: volume per length ... 93
 Tubing tolerances 95

RELATED PRODUCTS

Clean-cut tubing cutter 86
 Easy-Flange kit 85
 PEEK fittings
 High pressure ... 19-29
 Low pressure . 42,44-47
 PEEK tubing elbows ... 86