

## VICI caps

NEW

S40 threads

## INCLUDES

1 cap, S40 thread  
 1 set: 1/16" nut/ ferrule/  
 color-sleeve adapter  
 2 sets:  
 1/8" nut/ferrule/  
 color-sleeve adapter



- Fit all bottles with S40 thread
- 3 ports with 1/4-28 threads for 1/8" or 1/16" tubing

The latest addition to the Jour Safety product line, the S40 VICI cap seals to the bottle with an O-ring (EPDM) which is resistant to commonly used HPLC solvents.

	No. of ports	Product No.
VICI cap, S40	3	JR-S-14002

## SPECS

- Connections: 1/4-28
- Materials  
 Closure: PP\* (Polypropylene)  
 Insert: PTFE  
 O-ring: EPDM  
 Nuts: PPS\*  
 Ferrules: ETFE  
 Sleeves: PP\*  
 \*not a wetted part

## VICI caps

DIN50 threads

## INCLUDES

1 cap, DIN50 thread  
 1 O-ring,  
 1 set: 1/16" nut/ ferrule/  
 color-sleeve adapter  
 1, 2, or 3 sets:  
 1/8" nut/ferrule/ color-  
 sleeve adapter



- Fit all bottles or canisters with DIN50 thread
- 2, 3, or 4 ports with 1/4-28 threads for 1/8" or 1/16" tubing

DIN50 safety caps seal to the bottle/canister with an O-ring behind a PTFE lip, so the O-ring is therefore not a wetted part.

	No. of ports	Product No.
VICI cap, DIN50	2	JR-S-15001
	3	JR-S-15002
	4	JR-S-15003



## SPECS

- Connections: 1/4-28
- Materials  
 Closure: PP (Polypropylene)  
 Insert: PTFE  
 Stopcock: CTFE  
 O-ring: EPDM/  
 NBR  
 Nuts: PPS  
 Ferrules: ETFE  
 Sleeves: PP

## VICI caps

GL38 threads

## INCLUDES

1 cap, GL37/430 thread  
 1 O-ring  
 1 set: 1/16" nut/ ferrule/  
 color-sleeve adapter  
 1, 2, or 3 sets:  
 1/8" nut/ferrule/  
 color-sleeve adapter  
 1 adapter ring



- Fit all bottles or canisters with GL38/430 thread (See NOTES, below)
- 2, 3, or 4 ports with 1/4-28 threads for 1/8" or 1/16" tubing

Color-coded PP sleeves fit on all ports on 2 and 3 port caps. With 4 port caps, the use of PP sleeves is limited to two sleeves per cap. (Four colors are included.)

	No. of ports	Product No.
VICI cap, GL38	2	JR-S-13001
	3	JR-S-13002
	4	JR-S-13003
O-ring, FEP-coated, GL38		JR-S-40009

## SPECS

- Connections: 1/4-28
- Materials  
 Closure: PP (Polypropylene)  
 Insert: PTFE  
 O-ring: EPDM (alternative: FEP-coated)  
 Adapter ring: PA/PP  
 Nuts: PPS  
 Ferrules: ETFE  
 Sleeves: PP

## TECH TIP

PFA tubing provides less re-gassing of the mobile phase.

## MORE INFORMATION

Chemical resistance 96-98

## RELATED PRODUCTS

Color-Sleeve adapters . 17  
 No-Met mobile phase filters . . . . . 52  
 PFA tubing . . . . . 9  
 Safety-coated bottles (for additional sizes) . 70  
 Safety air inlet valves . . 59

## NOTES

GL38 thread is available in different variations. Please check the inside diameter of the neck:

- For 21-24 mm ID (e.g. Wheaton® bottles), use the GL38 cap as is, without an adapter.
- For 24-28 mm ID (e.g. Nalgene® bottles), use the adapter ring (JR-S-40006) which is included with this VICI cap.