

**Eco cap kit****INCLUDES**

- 1 cap, GL45 thread, with three ports
- 1 safety-coated glass bottle (option: 1 L or 2 L)
- 1 male Luer plug, polypropylene

In addition to an Eco cap with polypropylene male Luer plug, this kit offers a safety bottle (1 or 2L).

The Eco cap consists of a three-hole ETFE insert with a polyethylene collar, compatible with all standard GL45 thread bottles. The Eco cap has no threading for fittings on the inside of the tubing ports – the tubing is just inserted.

*Product No.*

Eco cap kit	
1 Liter bottle	JR-9000-0013
2 Liter bottle	JR-9000-0014

*GL45 threads***SPECS**

- Materials
- Screw collar: PP (polypropylene)
- Pouring ring: PP
- Insert: ETFE
- O-rings: EPDM
- Luer plug: PP

**Eco cap****INCLUDES**

- 1 cap, GL45 thread, with three ports
- 1 male Luer plug, polypropylene

- Competitively priced alternative to sparge and deliver mobile phases
- Also available as bottle kit

The Eco cap consists of a three-hole ETFE insert with a polyethylene collar, compatible with all standard GL45 thread bottles. The Eco cap has no threading for fittings on the inside of the tubing ports – the tubing is just inserted.

*Product No.*

Eco cap	JR-9000-0007
---------	--------------

*GL45 threads***SPECS**

- Materials
- Screw collar: PP (polypropylene)
- Pouring ring: PP
- Insert: ETFE
- O-rings: EPDM
- Luer plug: PP

**Spare parts***for Eco cap**Product No.*

Male Luer plug, polypropylene 10/pkg	JR-070350-10
---	--------------

**NOTE**

These caps are not usable for building up a helium atmosphere within in the solvent bottle, but only for continuous helium sparging.

**TECH TIP**

Once the mobile phase is degassed, we recommend our impermeable tubing between bottle cap and pump to prevent regassing of helium degassed solvents.

**RELATED PRODUCTS**

Impermeable tubing . . . 13  
Polymeric tubing . . . 9-13  
Waste caps . . . . . 66-68