

Columns

PEEK COLUMNS ■ GUARD COLUMNS

Columns index

PEEK columns

2.1 mm ID	page 77
4.0 mm ID	77

In-line biocompatible

guard columns	78
----------------------------	----

This PEEK guard column system extends the life of valuable analytical columns by trapping contaminants from the sample and mobile phase. Only PEEK and titanium contact the mobile phase and sample.

Column couplers and plugs

Column coupler, one-piece fingertight	79
Column end plugs	79
Column end plugs, logotype	79

Packing kits	78
---------------------------	----

VICI Jour column packing kits are ready-to-use tools to pack our PEEK columns and PEEK guard columns. Includes needed slurry reservoir, packing adapter, and packing seals.

