



VICI JOUR - Low Pressure Fittings

VICI JOUR PRODUCTS

Tubing for chromatography

High pressure fittings

- Nanovolume
- Stainless steel
- Polymeric

Low pressure fittings

Filters & mobile phase filters

Mobile phase handling systems

Frits

PEEK columns & guard columns

Sample loops

Accessories/miscellaneous

CUSTOMER SUPPORT

Peristaltic Tubing Adapter Set

- Connects soft tubing to PTFE and PTFE-like tubing
- Works with peristaltic tubing up to 4.50 mm OD
- Biocompatible
- Low dead volume
- For 1/16" and 1/8" OD PTFE and PTFE-like tubing

This "one-piece-fits-all" set consists of everything needed to connect 1/16" and 1/8" OD PTFE and PTFE-like tubing to soft tubing like silicone up to 4.50 mm OD. The adapter comes with three different bushings (for the different soft tubing OD's) and two PPS fittings with inverted ETFE ferrules (one each for 1/16" and 1/8" OD tubing). The body and bushings are PEEK.

Dimensions

Length (without fitting): 28 mm; diameter: 12 mm; bore: 0.50 mm



Peristaltic Tubing Adapter Set

Description	Product No.
PEEK adapter for peristaltic to PTFE tubing connection	JR-PTAS

TOP ▲