



NEW!

Micro Valves for GC and LC

- 200 psi, .060" bore
- Compact 1" design
- Convenient panel mount
- Variety of configurations

Simplify your liquid or gas handling application with a VICI Precision Sampling Micro valve. The unique design of the fitting detail allows a leak-free seal with no potential for rotor damage from overtightening. Internal parts are PEEK® and PTFE.



(Fittings not included.
Order on pages 42 and 44.)

SPECS

200 psi
.060" bore
1/4-28 fitting detail
All polymer-based materials

Micro valves for GC and LC

1/4-28

	Prod No	3 PORT	4 PORT		
"T" flow path	3 ports			"T"	
	4 ports				660110
180° flow path	2 ports			180°	
	4 ports				660210
90° flow path	2 ports			4 PORT	
	3 ports				660310
	4 ports				660320