

**Analytical valves,  
1/16" Valco fittings, 0.40 mm ports (.016")**

*Model C2*

**SPECS**  
**5000 psi liq**  
**75°C max**  
 Metal stator  
 Valcon H rotor  
  
**5000 psi liq**  
**50°C max**  
 PAEK stator  
 Valcon E rotor

Includes stainless steel nuts and ferrules of the stator material. Valves with PAEK stators have PEEK nuts and ferrules.

Standard electric actuator:  
 110 VAC for USA  
 110/230 VAC to 24 VDC power supply for international.  
 Microelectric actuator:  
 24 VDC, with 110/230 VAC to 24 VDC power supply.

\* The 6 port valve includes a 20 µl loop of the stator material.

**5,000 psi**

**Analytical**

**1/16"**

**0.40 mm**

**OPTIONS**

- Continuous flow version is available as Model C6. See page 163.
- Hastelloy C stators
- Semi-prep version with 0.75 mm ports (.030") available
- Loop fill port assembly for injection from front of the valve. See page 41.



Order loops from page 159.

**OPTIONAL FLOWPATH**

Model C2 6 port valves can also be ordered with a dual 3-way rotor, as described in EPA Method 555.

To specify this flowpath, substitute "6X" for "6" in the valve or rotor product number.



**AUTOSAMPLER REPLACEMENT VALVES**

The Cheminert Model C2 6 port valve is an excellent replacement for the valve originally supplied in many autosamplers, including autosamplers manufactured by Beckman, Gilson, Spark-Holland, CTC, Thermo Fisher, and Varian.

Call technical support to determine which replacement is best for your application.



**4 Port**

Prod No Price



**6 Port\***

Prod No Price



**8 Port**

Prod No Price



**10 Port**

Prod No Price

**N60 stainless stator**

Manual  
 With pneumatic actuator  
  
 With standard electric actuator  
 With microelectric actuator  
  
 Replacement valve  
 Replacement rotor  
 Replacement stator

C2-2004  
 C2-2004A  
  
 C2-2004E  
 C2-2004EH  
  
 C2-2004D  
 C2-20R4  
 C-2C04

C2-2006  
 C2-2006A  
  
 C2-2006E  
 C2-2006EH  
  
 C2-2006D  
 C2-20R6  
 C-2C06

C2H-2008  
 C2H-2008A  
  
 C2H-2008E  
 C2H-2008EH  
  
 C2H-2008D  
 C2-20R8H  
 C-2C08H

C2H-2000  
 C2H-2000A  
  
 C2H-2000E  
 C2H-2000EH  
  
 C2H-2000D  
 C2-20R0H  
 C-2C00H

**PAEK stator**

Manual  
 With pneumatic actuator  
  
 With standard electric actuator  
 With microelectric actuator  
  
 Replacement valve  
 Replacement rotor  
 Replacement stator

C2-2344  
 C2-2344A  
  
 C2-2344E  
 C2-2344EH  
  
 C2-2344D  
 C2-23R4  
 C-2C44

C2-2346  
 C2-2346A  
  
 C2-2346E  
 C2-2346EH  
  
 C2-2346D  
 C2-23R6  
 C-2C46

C2H-2348  
 C2H-2348A  
  
 C2H-2348E  
 C2H-2348EH  
  
 C2H-2348D  
 C2-23R8H  
 C-2C48H

C2H-2340  
 C2H-2340A  
  
 C2H-2340E  
 C2H-2340EH  
  
 C2H-2340D  
 C2-23R0H  
 C-2C40H

**Titanium stator**

Manual  
 With pneumatic actuator  
  
 With standard electric actuator  
 With microelectric actuator  
  
 Replacement valve  
 Replacement rotor  
 Replacement stator

C2-2034  
 C2-2034A  
  
 C2-2034E  
 C2-2034EH  
  
 C2-2034D  
 C2-20R4  
 C-2C34

C2-2036  
 C2-2036A  
  
 C2-2036E  
 C2-2036EH  
  
 C2-2036D  
 C2-20R6  
 C-2C36

C2H-2038  
 C2H-2038A  
  
 C2H-2038E  
 C2H-2038EH  
  
 C2H-2038D  
 C2-20R8H  
 C-2C38H

C2H-2030  
 C2H-2030A  
  
 C2H-2030E  
 C2H-2030EH  
  
 C2H-2030D  
 C2-20R0H  
 C-2C30H



**Model C2  
1/16" ZDV fittings**