




### NEW Integrated motor/stream selector, 1/16" Valco ZDV fittings, 0.40 mm ports (.016")

Model C55

- 5,000 psi**
- Integrated**
- Stream selector**
- 10-32 ZDV**
- 1/16"**
**0.40 mm**

Includes stainless steel nuts and ferrules of the stator material.

Valves with PAEK stators have PEEK nuts and ferrules.



#### SPECS

**5000 psi liq**  
**50°C max**  
Metal stator  
Valcon H rotor



**5000 psi liq**  
**50°C max**  
PAEK stator  
Valcon E rotor



	<b>4 Position</b>	<b>6 Position</b>	<b>8 Position</b>	<b>10 Position</b>
	<i>Prod No</i>	<i>Prod No</i>	<i>Prod No</i>	<i>Prod No</i>
<b>N60 stainless stator</b>				
With integrated actuator	C55-2004I	C55-2006I	C55-2008I	C55-2000I
With motor/sensor only	C55-2004I-S	C55-2006I-S	C55-2008I-S	C55-2000I-S
With motor only	C55-2004IX	C55-2006IX	C55-2008IX	C55-2000IX
Replacement rotor	C5-20R4	C5-20R6	C5-20R8H	C5-20R0H
Replacement stator	C55-2C04	C55-2C06	C55-2C08	C55-2C00
<b>PAEK stator</b>				
With integrated actuator	C55-2344I	C55-2346I	C55-2348I	C55-2340I
With motor/sensor only	C55-2344I-S	C55-2346I-S	C55-2348I-S	C55-2340I-S
With motor only	C55-2344IX	C55-2346IX	C55-2348IX	C55-2340IX
Replacement rotor	C5-23R4	C5-23R6	C5-23R8H	C5-23R0H
Replacement stator	C55-2C44	C55-2C46	C55-2C48	C55-2C40

#### OPTIONS

- Optional bore:  
0.25 mm (.010")  
0.75 mm (.030")
- 4 and 8 positions available



### NEW Integrated motor/stream selector, 1/16" Valco ZDV fittings, 0.75 mm ports (.030")

Model C65Z

<b>SPECS</b> 100 psi gas/ 250 psi liq 50°C max PPS stator Valcon E2 rotor	Includes Valco ZDV PEEK nuts and ferrules.			<b>Low pressure</b>
				<b>Integrated</b>
			<b>Stream selector</b>	
			<b>10-32 ZDV</b>	
			<b>1/16" 0.75 mm</b>	
	<b>4 Position</b> Prod No	<b>6 Position</b> Prod No	<b>8 Position</b> Prod No	<b>10 Position</b> Prod No
With integrated actuator	C65Z-3184I	C65Z-3186I	C65Z-3188I	C65Z-3180I
With motor and sensor only	C65Z-3184I-S	C65Z-3186I-S	C65Z-3188I-S	C65Z-3180I-S

### NEW Integrated motor/stream stream selector, 1/4-28 fittings for 1/16" tubing, 0.75 mm ports (.030")

Model C65

<b>SPECS</b> 100 psi gas/ 250 psi liq 50°C max PPS stator Valcon E2 rotor	Includes multicolored Cheminert flangeless fittings for 1/16" tubing.	See photo on facing page.		<b>Low pressure</b>
				<b>Integrated</b>
			<b>Stream selector</b>	
			<b>1/4-28 Internal</b>	
			<b>1/16" 0.75 mm</b>	
	<b>4 Position</b> Prod No	<b>6 Position</b> Prod No	<b>8 Position</b> Prod No	<b>10 Position</b> Prod No
With integrated actuator	C65-3184I	C65-3186I	C65-3188I	C65-3180I
With motor and sensor only	C65-3184I-S	C65-3186I-S	C65-3188I-S	C65-3180I-S

### NEW Integrated motor/stream stream selector, 1/4-28 fittings for 1/8" tubing, 1.50 mm ports (.060")

Model C65

<b>SPECS</b> 100 psi gas/ 250 psi liq 50°C max PPS stator Valcon E2 rotor	Includes multicolored Cheminert flangeless fittings for 1/8" tubing.	See photo on facing page.		<b>Low pressure</b>
				<b>Integrated</b>
			<b>Stream selector</b>	
			<b>1/4-28 Internal</b>	
			<b>1/8" 1.50 mm</b>	
	<b>4 Position</b> Prod No	<b>6 Position</b> Prod No	<b>8 Position</b> Prod No	<b>10 Position</b> Prod No
With integrated actuator	C65-6184I	C65-6186I	C65-6188I	C65-6180I
With motor and sensor only	C65-6184I-S	C65-6186I-S	C65-6188I-S	C65-6180I-S