

Post-Column Reaction Tee Fittings


2009 #60

Post-column reaction tee fitting

The tee column end fitting (TCEF) has a third connection perpendicular to the normal flowpath. The TCEF permits post-column derivation, or may be used as a curtain flow column inlet fitting. Standard material is Type 316 stainless.



Column OD	Cone OD	Inlet/outlet OD	Bore	Side OD	Bore	Prod No	Price
1/16"	1.0 mm	1/32"	0.25 mm	1/32"	0.25 mm	TCEF1.5.5C	
1/16"	1.0 mm	1/32"	0.90 mm	1/32"	0.25 mm	TCEF1.5.5T	
1/16"	1.0 mm	1/16"	0.25 mm	1/16"	0.25 mm	TCEF111	
1/8"	1.0 mm	1/16"	0.50 mm	1/16"	0.50 mm	TCEF211	
1/8"	1.0 mm	1/16"	1.65 mm	1/16"	0.40 mm	TCEF211T	
1/4"	4.6 mm	1/16"	0.25 mm	1/16"	0.25 mm	TCEF411C	
1/4"	4.6 mm	1/16"	0.75 mm	1/16"	0.75 mm	TCEF411	
1/4"	4.6 mm	1/16"	1.65 mm	1/16"	0.75 mm	TCEF411T	
1/4"	4.6 mm	1/8"	0.75 mm	1/16"	0.75 mm	TCEF421	
3/8"	6.0 mm	1/16"	0.75 mm	1/16"	0.75 mm	TCEF611	
3/8"	6.0 mm	1/16"	1.65 mm	1/16"	0.75 mm	TCEF611T	
1/2"	9.0 mm	1/16"	0.75 mm	1/16"	0.75 mm	TCEF811	
1/2"	9.0 mm	1/16"	1.65 mm	1/16"	0.75 mm	TCEF811T	



Post-column reaction fitting (TCEF411)

0.25 mm	=	.010"
0.50 mm	=	.020"
0.75 mm	=	.030"
1.0 mm	=	.040"
1.5 mm	=	.060"
2.0 mm	=	.080"
4.6 mm	=	.180"
6.0 mm	=	.236"
6.4 mm	=	.253"
7.0 mm	=	.275"
10.0 mm	=	.400"
27.0 mm	=	1.08"
1/32"	=	0.8 mm
1/16"	=	1.6 mm
1/8"	=	3.2 mm
1/4"	=	6.4 mm
3/8"	=	9.5 mm
1/2"	=	12.7 mm

CHROMalytic TECHNOlogy Pty Ltd AUSTRALIAN Distributors e-mail: sales@chromtech.net.au Tel: 03 9762 2034