


Low Pressure Bulkhead Unions

Bulkhead unions *Cheminert to Cheminert* *1/4-28 to 1/4-28*

Includes flangeless 1/4-28 fittings for tubing OD indicated.

Tubing OD	Bore	PEEK		CTFE		316 Stainless	
		Prod No	Price	Prod No	Price	Prod No	Price
1/16"	0.25 mm	CBUCPK		CBUCKF		CBUCS6	
1/16"	0.50 mm	CBUPK		CBUKF		CBUS6	
1/16"	0.75 mm	CBUMPK		CBUMKF		CBUMS6	
1/8"	1.50 mm	CBULPK		CBULKF		CBULS6	



Bulkhead unions *Cheminert to 1/16" ZDV* *1/4-28 to 10-32*

Includes flangeless 1/4-28 and ZDV 10-32 fittings for 1/16" OD tubing.

Tubing OD	Bore	PEEK		CTFE		316 Stainless	
		Prod No	Price	Prod No	Price	Prod No	Price
1/16"	0.25 mm	CZBUCPK		CZBUCKF		CZBUCS6	
1/16"	0.50 mm	CZBUPK		CZBUKF		CZBUS6	
1/16"	0.75 mm	CZBUMPK		CZBUMKF		CZBUMS6	

