

Adapters

Tube adapters

1/4-28

Tube adapters have male 1/4-28 threads going to 1/4" or 1/8" OD tubing.

Tubing		PEEK		CTFE		316 Stainless	
OD	Bore	Prod No	Price	Prod No	Price	Prod No	Price
1/8"	1.5 mm	CTA2PK		CTA2KF		CTA2S6	
1/4"	1.5 mm	CTA4PK		CTA4KF		CTA4S6	



Luer adapters

Luer to 1/4-28 or 10-32

Luer adapters make a leak-tight connection from male or female luer to 1/4-28 threads.

Description	Bore	PEEK		CTFE		PFA	
		Prod No	Price	Prod No	Price	Prod No	Price
Female luer to 1/4-28	1.50 mm	CFLAPK		CFLAKF		CFLAPFA	
to 10-32	0.75 mm	ZUFLPK		ZUFLKF		-	
Male luer to 1/4-28	1.50 mm	CMLAPK		CMLAKF		CMLAPFA	



Luer adapter bulkhead unions

Luer to 1/4-28 or 10-32

Our luer adapter bulkhead union connects a male or female luer to 1/4-28 or 10-32 fittings. These are the ideal fittings for through-the-panel syringe injections. The 1/4-28 versions include flangeless fittings for 1/16" OD tubing. Versions with 10-32 connections (for 1/16" OD tubing) include a fingertight PEEK nut and a ferrule of the same material as the union.



Description	Bore	PEEK		CTFE	
		Prod No	Price	Prod No	Price
Female luer to 1/4-28	1.50 mm	CBUFLPK		CBUFLKF	
to 10-32	1.00 mm	ZBUFLPK		ZBUFLKF	
Male luer to 10-32	1.00 mm	ZBUMLPK		ZBUMLKF	