



2009 #60



### Mobile phase or solvent reservoirs 1/4-28

These high density polyethylene reservoirs for in-line solvent use come with polypropylene caps, 1/4-28 flangeless fittings, and 1/8" PTFE tubes for one or two lines plus vent. Plugs are included for conversion to solvent storage when the reservoir is removed from the system. Optional PTFE filters with titanium frits are available on the facing page.

Capacity	Cap	Prod No	Price
0.5 liter	2-hole	C-MPR2	
0.5 liter	3-hole	C-MPR3	
0.5 liter	plain	C-BOT16	
1.0 liter	plain	C-BOT32	

### VICI-cap



The VICI-cap is the most economical way to helium sparge and deliver HPLC mobile phases. The insert is manufactured from PTFE with an EPDM\* O-ring and a polypropylene screw cap.

The VICI-cap is available for either GL45 or S40 threaded bottles. It has a 1/4-28 female port and three ports for tubing insertion: two 1/8" tubing ports and a 1/16" tubing port. The tubing ports are made so that you push the tubing through the hole, while 1/4-28 fittings provide the best connection. Unused ports can be plugged as required.

	Prod No	Price
VICI-cap GL-45	JR-9000-0001	
VICI-cap S40	JR-9000-0006	

### Valves for vials

The screw-cap Mininert is available in a variety of sizes. The crimp-top valve for 13 mm ID glassware slides into the neck of the vial and features a threaded flange, which is turned to provide a leak-tight fit.



Pkg/12:	Cap/thread size	Prod No	Price
	13 mm-425	PS-614158	
	15 mm-425	PS-614160	
	18 mm-400	PS-614161	
	20 mm-400	PS-614170	
	24 mm-400	PS-614163	
	Crimp top	PS-614250	

#### MORE INFORMATION

Bulkhead connectors . . . page 73  
 Flangeless fittings . . . . 68  
 Plugs . . . . . 71  
 Polymeric tubing . . . . . 90

#### TECH TIP

The VICI-cap is not usable for building up a helium atmosphere within the solvent bottle. It is only designed for continuous helium sparging.

\*Ethylene Propylene Diene Monomer