

## Polymeric Fittings – High Pressure

### PEEK Color Coded & natural One-Piece Fingertight Fittings

#### PEEK Color Coded & natural One-Piece Fingertight Fittings

- Fingertight up to 350 bar (5000 psi)
- Excellent chemical resistance
- Biocompatible

To identify various tubing lines the PEEK fingertight fittings are available in six colors. PEEK fingertights are rated to 350 bar (5000 psi) and can be used in virtually any type of HPLC fitting detail on the market with 10-32 thread.



#### SPECS

**Material**  
PEEK

**Threads**  
10-32

**Dimensions**  
For OD 1/16" Tubing

**Tolerances**  
ID +/- 0.05 mm (.002")

**Pressure rating**  
< 350 bar (< 5000 psi)

**Special Info**  
Fingertight Fittings – works with all 1/16" OD tubing and female 10-32 ports

Color Coded – six colors to identify tubing lines

Biocompatible

The coloration pigments are carefully chosen not to contain any heavy metals and hazardous chemicals

#### PEEK Color Coded One-Piece Fingertight Fittings-molded version

Part No.	Description	Qty/pkg
JR-55020-5	Fitting, PEEK, one-piece natural, 10-32	5
JR-55021-5	Fitting, PEEK, one-piece black, 10-32	5
JR-55022-5	Fitting, PEEK, one-piece red, 10-32	5
JR-55023-5	Fitting, PEEK, one-piece yellow, 10-32	5
JR-55024-5	Fitting, PEEK, one-piece blue, 10-32	5
JR-55025-5	Fitting, PEEK, one-piece green, 10-32	5



Part No. JR-55025-5