


## VICI JOUR - Low Pressure Fittings

### VICI JOUR PRODUCTS

- Tubing for chromatography
- High pressure fittings
  - Nanovolume
  - Stainless steel
  - Polymeric
- Low pressure fittings
- Filters & mobile phase filters
- Mobile phase handling systems
- Frits
- PEEK columns & guard columns
- Sample loops
- Accessories/miscellaneous

### Low Pressure Bulkhead Unions

- Panel mounting
- Easy connection of low pressure solvent lines
- Low dead volume
- Completely inert and biocompatible



### Low Pressure Bulkhead Unions

Pressure rating: up to 105 bar (1500 psi) for 1/16", 35 bar (500 psi) for 1/8"

Description	Product No.
PEEK, 0.75 mm bore, with ETFE ferrules and PPS flangeless nuts for 1/16" tubing	JR-CBUMPK
PEEK, 1.50 mm bore, with ETFE ferrules and PPS flangeless nuts for 1/8" tubing	JR-CBULPK

### Spare Parts

Description		Qty/pkg	Product No.
Nuts	PPS, flangeless, 1/16", black	10	JR-55070-10
	PPS, flangeless, 1/8", black	10	JR-55071-10
Ferrules	ETFE, 1/16"	10	JR-041-10
	ETFE, 1/8"	10	JR-051-10
Adapters	PP fingertight sleeve, yellow	10	JR-55086-10
	PP fingertight sleeve, blue	10	JR-55090-10