



## VICI JOUR - High Pressure Fittings

### VICI JOUR PRODUCTS

Tubing for chromatography

High pressure fittings

- Nanovolume
- Stainless steel
- Polymeric

Low pressure fittings

Filters & mobile phase filters

Mobile phase handling systems

Frits

PEEK columns & guard columns

Sample loops

Accessories/miscellaneous

### ValvTool

- For making connections in those hard-to-reach areas
- Ideal for nuts with 1/4" and 5/16" heads



The new ValvTool is a great timesaver, providing easy access to many hard-to-reach areas. The unique design with its slotted wrench allows you to tighten nuts where an interfering loop or capillary might otherwise make it difficult. The ValvTool is the ideal for Rheodyne® and Valco valves and fittings as well as for most of the HPLC and PEEK fittings on the market with 1/4" and 5/16" heads.

### CUSTOMER SUPPORT

Description	Product No.
ValvTool	JR-800

### Further Reference

- [File](#)
- [Open-end wrenches](#)
- [Tweezer](#)
- [Ferrule removal kit](#)
- [Hex key set](#)
- [Drill set](#)

[Home](#) | [Products](#) | [Search](#) | [Contact VICI Jour](#) | [About VICI Jour](#)

VICI® is a registered trademark of Valco Instruments Co. Inc. and VICI AG  
© 2005 VICI