

BIOPROCESS

FITTINGS

PharmaLok™ Clamps
PG-19



Bag Ports
PG-20



Mini Flange Fittings 3/4"
PG-22



Maxi Flange Fittings 1-1/2"
PG-23



Mini Flange Accessories 3/4"
PG-24



Maxi Flange Accessories 1-1/2"
PG-24



PHARMALOK™ Series

Value Plastics' new sanitary fitting clamp was created with the end user in mind. Easily closed with one hand, the clamp's unique and innovative design is one-of-a-kind!

The biopharmaceutical industry has made a strong move towards single-use technology so Value Plastics has developed a ground-breaking design that simplifies the assembly process by making it effortless to securely attach the fittings to bag ports and filters. Complies with ASME BPE 2009 DT-10.4 Hygienic Clamp Ferrule: Design Criteria.

Features:

- ASME BPE 2009 DT-10.4 compliant
- Multiple colors for color coding applications
- Connects both plastic and stainless steel fittings
- Zip tie compatible
- Easy to use, one-handed closure (no tools needed)
- 360° compression
- Tamper evident
- SIP/CIP compatible

Specifications:

Material: USP Class VI & ISO 11137 Certified glass-filled nylon

Sterilization: Gamma Radiation up to 50 kGy,
Autoclavable up to 273°F (134° C) for a max of 30 min.

Size: .75" & 1.5", for use with "Mini" & "Maxi" flange sanitary fittings

Operating Temperature: -94° (-70° C) to 140° (60° C)

Maximum Pressure: 80 psi