

MALE LUER: Slip

MTLS210 (-1 -6005 -9 -J1A)

Male Luer Slip to 200 Series Barb, 1/16" (1.6 mm) ID Tubing



MSL004 (-1 -6005 -9 -9002)

Male Luer Slip to 500 Series Barb, 1/16" (1.6 mm) ID Tubing



MTLS220 (-1 -3 -5 -6005 -81 -9 -J1A)

Male Luer Slip to 200 Series Barb, 3/32" (2.4 mm) ID Tubing



MSL007 (-1 -6005 -9 -9002)

Male Luer Slip to 500 Series Barb, 3/32" (2.4 mm) ID Tubing



MTLS230 (-1 -3 -5 -6005 -81 -9 -9002 -J1A)

Male Luer Slip to 200 Series Barb, 1/8" (3.2 mm) ID Tubing



MSL013 (-1 -6005 -9 -9002)

Male Luer Slip to 500 Series Barb, 1/8" (3.2 mm) ID Tubing



MTLS240 (-1 -2 -5 -6005 -81 -9 -J1A)

Male Luer Slip to 200 Series Barb, 5/32" (4.0 mm) ID Tubing



MSL025 (-1 -6005 -9 -9002)

Male Luer Slip to 500 Series Barb, 5/32" (4.0 mm) ID Tubing



MSL035 (-1 -6005 -9 -9002 -J1A)

Male Luer Slip to 500 Series Barb, 3/16" (4.8 mm) ID Tubing



MSL055 (-1 -6005 -9 -9002)

Male Luer Slip to 500 Series Barb, 1/4" (6.4 mm) ID Tubing



FEMALE LUER: Thread Style

FTLL004 (-1 -2 -3 -4 -5 -06 -6005 -07 -9 -9002 -J1A)

Female Luer Thread Style to 500 Series Barb, 1/16" (1.6 mm) ID Tubing



FTLL210 (-1 -2 -3 -4 -5 -6005 -9 -9002 -J1A)

Female Luer Thread Style to 200 Series Barb, 1/16" (1.6 mm) ID Tubing



FTLL410 (-1 -6005 -9 -9002 -J1A)

Female Luer Thread Style to 400 Series Barb, 1/16" (1.6 mm) ID Tubing



FTLL007 (-1 -2 -3 -4 -5 -06 -6005 -07 -9 -9002)

Female Luer Thread Style to 500 Series Barb, 3/32" (2.4 mm) ID Tubing



FTLL220 (-1 -2 -3 -5 -6005 -9 -9002 -J1A)

Female Luer Thread Style to 200 Series Barb, 3/32" (2.4 mm) ID Tubing



FTLL420 (-1 -6005 -9 -9002 -J1A)

Female Luer Thread Style to 400 Series Barb, 3/32" (2.4 mm) ID Tubing



FTLL013 (-1 -2 -3 -4 -5 -06 -6005 -07 -9 -9002)

Female Luer Thread Style to 500 Series Barb, 1/8" (3.2 mm) ID Tubing



FTLL230 (-1 -2 -3 -5 -6005 -9 -9002 -J1A)

Female Luer Thread Style to 200 Series Barb, 1/8" (3.2 mm) ID Tubing



FTLL430 (-1 -6005 -9 -9002 -J1A)

Female Luer Thread Style to 400 Series Barb, 1/8" (3.2 mm) ID Tubing



29

Material Reference

-001 (various, see part description)

-06 Orange Nylon

-07 Yellow Nylon

-0 Natural Nylon

-1 White Nylon

-2 Black Nylon

-3 Red Nylon

-4 Green Nylon

-5 Blue Nylon

-7 Glass Filled Natural Nylon

-8 Gray Nylon

-8B Gray Nylon

-9 Clear Polycarbonate

-10 Acetal

-40 Polysulfone

-81 Lace White ABS

-91 White Polycarbonate

-1006 Natural Acetal

-1007 Black Acetal

-6005 Animal-Free Polypropylene

-6006 White Polypropylene

-8003 Natural ABS

-8012 White ABS

-9002 Rad. Stable Polycarbonate

-9010 Rad. Stable Polycarbonate

-B Buna-N O-ring

-CM030 Dark Blue Nylon

-E PTEE/EPDM

-EP EPDM O-ring

-J1A Natural Kynar PVDF

-J1A2 Black Kynar PVDF

-PX Platinum-Cured Silicone

-U71 White Class VI Glass Filled Nylon

-U75 Blue Class VI Glass Filled Nylon

-VP1 Animal-free Polyethylene

-X2 Buna-N Rubber

-XO Stainless Steel

LUER
FITTINGS

Resin Reference

Pg 82 Chemical Resistance

Pg 83 Regulatory Status

Pg 84 Sterilization Stability

Pg 85 Typical Properties

025a