

**LE7230 (-1 -2 -6005 -9 -9010)**  
Male Luer to 200 Series Barb Elbow 1/8" (3.2 mm) ID Tubing (May be used with separate rotating lock ring; FSLLR)



**CA091 (-001) Gamma Stable Polycarbonate**  
Male Luer Elbow to 3/16" Classic Series Barb, (CF361-9002) with Rotating Lock Ring (FSLLR-9002)



LUER  
FITTINGS

- 001 (various, see part description)
- 06 Orange Nylon
- 07 Yellow Nylon
- 0 Natural Nylon
- 1 White Nylon
- 2 Black Nylon
- 3 Red Nylon
- 4 Green Nylon
- 5 Blue Nylon
- 7 Glass Filled Natural Nylon
- 8 Gray Nylon
- 8B Gray Nylon
- 9 Clear Polycarbonate
- 10 Acetal
- 40 Polysulfone
- 81 Lace White ABS
- 91 White Polycarbonate
- 1006 Natural Acetal
- 1007 Black Acetal
- 6005 Animal-Free Polypropylene
- 6006 White Polypropylene
- 8003 Natural ABS
- 8012 White ABS
- 9002 Rad. Stable Polycarbonate
- 9010 Rad. Stable Polycarbonate
- B Buna-N O-ring
- CM030 Dark Blue Nylon
- E PTEE/EPDM
- EP EPDM O-ring
- J1A Natural Kynar PVDF
- J1A2 Black Kynar PVDF
- PX Platinum-Cured Silicone
- U71 White Class VI Glass Filled Nylon
- U75 Blue Class VI Glass Filled Nylon
- VP1 Animal-free Polyethylene
- X2 Buna-N Rubber
- XO Stainless Steel

### Need samples for evaluation?

You can call to request most of our standard product in quantities of 5, or less, for free application evaluation.  
**+1 (970) 267-5200 Toll Free: (888) 404-5837**

### Need it fast?

You can order any of our 4,000+ parts online now with no order \$ minimums, 24 hours a day, 7 days a week!

[www.valueplastics.com](http://www.valueplastics.com)

**MTLLP (-1 -2 -3 -6005 -07 -9 -J1A)**  
Male Luer Integral Lock Ring Plug, Closed at Luer Tip



**LP4 (-1 -2 -3 -4 -5 -6005 -9 -9002 -J1A)**  
Male Luer Integral Lock Ring Plug, Closed at Grip



**VLP4 (-1 -3 -6005 -81 -J1A)**  
Male Luer Integral Lock Ring Vented Plug, Closed at Grip



**MLRLP (-1 -6005 -9)**  
Male Luer Plug (May be used with separate rotating lock ring; FSLLR)



**MLSPL (-1 -6005 -9)**  
Male Luer Plug (May be used with separate stationary lock ring; FSLLR)



**MTLP (-1 -6005 -9 -J1A)**  
Male Luer Plug, Closed at Tether Loop (May be used with separate stationary lock ring; MLLR)



**MTLSP (-1 -6005 -81 -9 -J1A)**  
Male Luer Slip Plug, Closed at Grip



**FTLP (-1 -6005 -9 -J1A)**  
Female Luer Lug Style Cap with Tether Loop



**FTLLP (-1 -2 -3 -5 -6005 -07 -81 -9 -9002 -J1A)**  
Female Luer Thread Style Cap



**LP34 (-1 -3 -5 -6005 -9 -9002 -J1A)**  
Male Luer Integral Lock Ring to Female Luer Thread Style Plug



**Pg 82** Chemical Resistance

**Pg 83** Regulatory Status

**Pg 84** Sterilization Stability

**Pg 85** Typical Properties