

**BDMT210 (-40 -81)**

Bondable Tee 200 Series Barb, 1/16" (1.6 mm) ID Tubing

**BDMT410 (-40 -81)**

Bondable Tee 400 Series Barb, 1/16" (1.6 mm) ID Tubing

**BDMT220 (-40 -81)**

Bondable Tee 200 Series Barb, 3/32" (2.4 mm) ID Tubing

**BDMT420 (-40 -81)**

Bondable Tee 400 Series Barb, 3/32" (2.4 mm) ID Tubing

**BDMT230 (-40 -81)**

Bondable Tee 200 Series Barb, 1/8" (3.2 mm) ID Tubing

**BDMT430 (-40 -81)**

Bondable Tee 400 Series Barb, 1/8" (3.2 mm) ID Tubing

**BDMX (-40 -81)**

Bondable 10-32 Special Tapered Thread with 1/4" Hex

**BDMK (-40 -81)**

Bondable 10-32 UNF Thread with 1/4" Hex

**BDSM (-81)**

Bondable 1/4-28 UNF Thread with 5/16" Hex

**BDM6 (-40 -81)**

Bondable M6x1 Thread with 5/16" Hex

**BDM16 (-40 -81)**

Bondable 1/16-27 NPT Thread with 5/16" Hex

**BDM18 (-40 -81)**

Bondable 1/8-27 NPT Thread with 7/16" Hex

**BDF16 (-40 -81)**

Bondable 1/16-27 NPT Thread with 5/16" Hex

**BDF18 (-40 -81)**

Bondable 1/8-27 NPT Thread with 7/16" Hex

**BDFM6 (-40 -81)**

Bondable M6x1 Thread with 5/16" Hex

**BDFS (-40 -81)**

Bondable 1/4-28 UNF Thread with 5/16" Hex

**-001** (various, see part description)**-06** Orange Nylon**-07** Yellow Nylon**-0** Natural Nylon**-1** White Nylon**-2** Black Nylon**-3** Red Nylon**-4** Green Nylon**-5** Blue Nylon**-7** Glass Filled Natural Nylon**-8** Gray Nylon**-8B** Gray Nylon**-9** Clear Polycarbonate**-10** Acetal**-40** Polysulfone**-81** Lace White ABS**-91** White Polycarbonate**-1006** Natural Acetal**-1007** Black Acetal**-6005** Animal-Free Polypropylene**-6006** White Polypropylene**-8003** Natural ABS**-8012** White ABS**-9002** Rad. Stable Polycarbonate**-9010** Rad. Stable Polycarbonate**-B** Buna-N O-ring**-CM030** Dark Blue Nylon**-E** PTEE/EPDM**-EP** EPDM O-ring**-J1A** Natural Kynar PVDF**-J1A2** Black Kynar PVDF**-PX** Platinum-Cured Silicone**-U71** White Class VI Glass Filled Nylon**-U75** Blue Class VI Glass Filled Nylon**-VP1** Animal-free Polyethylene**-X2** Buna-N Rubber**-XO** Stainless SteelBUILD-A  
PART**Need it fast?**

You can order any of our 4,000+ parts online now with no order \$ minimums, 24 hours a day, 7 days a week!

[www.valueplastics.com](http://www.valueplastics.com)**Resin Reference****Pg 82** Chemical Resistance**Pg 83** Regulatory Status**Pg 84** Sterilization Stability**Pg 85** Typical Properties