

OTHER COMPONENTS



Check Valves
PG-72



Clamps
PG-72



Injection Sites
PG-73



Stopcocks
PG-73



Spikes
PG-73



Tubing
PG-74



Designer Kits
PG-75



Tubesetter Assembly Tool
PG-77



OTHER COMPONENTS: Components

VPS5401068N (Polycarbonate with Silicone Disk)

Check Valve, Male/Female Luer, 420 mL/min flow rate at 1 M head 30 psi back pressure, 1/4 psi inlet cracking pressure



VPS5401020N (Polycarbonate with Silicone Disk)

Check Valve, Male/Female Luer, 129 mL/min flow rate at 100mmHg, 45 psi back pressure, 1-5 psi inlet cracking pressure



VPS5401036SN (Polycarbonate with Silicone Disk)

Reflux Check Valve, Male/Female Luer, 350mL/min flow Rate at 1 M head, 45 psi back pressure, opens at 284 mmHg



VPS5401001N (ABS Housing with Silicone Disk)

Check Valve, Male Slip Luer, 420 mL/min flow rate at 1M head, 30 psi back pressure, 1/4 psi inlet cracking pressure



VPM2930601N (White ABS)

Tube Clamp with split, 1/8" (3.2 mm) OD Tubing

