

**VPS3901074N (Polyisoprene with ABS Acrylic Body)**

IN-4000 Intermittent Injection Site to Male Luer Lock

**VPM1060200N (Clear Low Density Polyethylene)**

Long Spike Guard Closed at Tip

**VPM0480201N (White Polyethylene)**

Luer Tip Syringe Cap

**VPB0850400N (Clear ABS)**

Female Luer Thread Style Coupler

**VPM1080303N (Blue Polyethylene)**

Male Luer Integral Lock Ring Plug

**VPB1000050N (Polycarbonate Body with Polyethylene Handle)**

One Way Stopcock, Capped Male Luer Lock to Capped Female Luer Thread Style

**VPB1000055N (Polycarbonate Body with Polyethylene Handle)**

One way Stopcock, Male Slip Luer to Capped Female Luer Thread Style

**VPB1000079N (Polycarbonate Body with Polyethylene Handle)**

Three Way Stopcock, (2) Capped Female Luer Thread Style and Capped Male Luer Lock

**VPB1000062N (Polycarbonate Body with Polyethylene Handle)**

Three Way Stopcock, (2) Female Luer Thread Style and Male Slip Luer

**VPB1000099N (Polycarbonate Body with Polyethylene Handle)**

Four Way Stopcock, (2) Capped Female Luer Thread Style and Male Slip Luer

**VPB1000110N (Polycarbonate Body with Polyethylene Handle)**

Four Way Stopcock, (2) Capped Female Luer Thread Style &amp; Capped Male Luer Lock

**VPS5069200N (White ABS Spike with Clear LDPE Guard)**

Bag Spike with vented guard (replaces VPS5069700N)

**-001** (various, see part description)**-06** Orange Nylon**-07** Yellow Nylon**-0** Natural Nylon**-1** White Nylon**-2** Black Nylon**-3** Red Nylon**-4** Green Nylon**-5** Blue Nylon**-7** Glass Filled Natural Nylon**-8** Gray Nylon**-8B** Gray Nylon**-9** Clear Polycarbonate**-10** Acetal**-40** Polysulfone**-81** Lace White ABS**-91** White Polycarbonate**-1006** Natural Acetal**-1007** Black Acetal**-6005** Animal-Free Polypropylene**-6006** White Polypropylene**-8003** Natural ABS**-8012** White ABS**-9002** Rad. Stable Polycarbonate**-9010** Rad. Stable Polycarbonate**-B** Buna-N O-ring**-CM030** Dark Blue Nylon**-E** PTEE/EPDM**-EP** EPDM O-ring**-J1A** Natural Kynar PVDF**-J1A2** Black Kynar PVDF**-PX** Platinum-Cured Silicone**-U71** White Class VI Glass Filled Nylon**-U75** Blue Class VI Glass Filled Nylon**-VP1** Animal-free Polyethylene**-X2** Buna-N Rubber**-XO** Stainless SteelOTHER  
COMPONENTS

## Resin Reference

**Pg 82** Chemical Resistance**Pg 83** Regulatory Status**Pg 84** Sterilization Stability**Pg 85** Typical Properties