

**SCV SERIES:**  
One-Way Check Valve

**SCV SERIES:**  
Anti-Siphon Valve

**SCV21053 (Clear and White ABS w/Silicone Diaphragm)**

Check Valve, Socket for .118" (3.0 mm) ID x .161" (4.1 mm) OD Tubing x 2



Cracking Pressure:  $\leq .174$  psig  
Max. Back Pressure: 87.023 psig  
Flow Rate:  $\geq 90$  ml/min  
Flow Direction: Left to right

**SCV22050 (Clear SAN and Clear ABS w/Silicone Diaphragm)**

Check Valve, Socket for .118" (3.0 mm) ID x .161" (4.1 mm) OD to Male Locking Luer



Cracking Pressure:  $\leq .174$  psig  
Max. Back Pressure: 87.023 psig  
Flow Rate:  $\geq 90$  ml/min  
Flow Direction: Left to right

**SCV23050 (Clear SAN w/Silicone Diaphragm)**

Check Valve, Female Locking Luer to Male Locking Luer



Cracking Pressure:  $\leq .174$  psig  
Max. Back Pressure: 87.023 psig  
Flow Rate:  $\geq 90$  ml/min  
Flow Direction: Left to right

**SCV25052 (Clear ABS and Blue SAN w/Silicone Diaphragm)**

Check Valve, Stem for .157" (4.0 mm) ID Tubing to Stem for .118" (3.0 mm) ID Tubing



Cracking Pressure:  $\leq .174$  psig  
Max. Back Pressure: 87.023 psig  
Flow Rate:  $\geq 90$  ml/min  
Flow Direction: Left to right

**SCV SERIES:**  
Anti-Siphon Valve

**SCV03052 (Clear and Blue ABS w/Silicone Diaphragm)**

Check Valve, Pocket for .165" (4.2 mm) OD Tubing to Socket Fitting .122" (3.1 mm) ID x .165" (4.2 mm) OD Tubing



Cracking Pressure: 1.450-4.351 psig  
Max. Back Pressure: 304.579 psig, max 30 sec.  
Flow Rate: max 200 ml/min  
Flow Direction: Left to right

**SCV04052 (Clear and Blue SAN w/Silicone Diaphragm)**

Check Valve, Female Locking Luer to Male Locking Luer



Cracking Pressure: 1.450-4.351 psig  
Max. Back Pressure: 304.579 psig, max 30 sec.  
Flow Rate: max 200 ml/min  
Flow Direction: Left to right

**SCV05052 (Clear SAN and Blue Acrylic w/Silicone Diaphragm)**

Check Valve, Pocket for .165" (4.2 mm) OD Tubing to Male Locking Luer



Cracking Pressure: 1.450-4.351 psig  
Max. Back Pressure: 304.579 psig, max 30 sec.  
Flow Rate: max 200 ml/min  
Flow Direction: Left to right

**SCV06052 (Clear and Blue ABS w/Silicone Diaphragm)**

Check Valve, Female Locking Luer to Socket for .122" (3.1 mm) ID to .165" (4.2 mm) OD Tubing



Cracking Pressure: 1.450-4.351 psig  
Max. Back Pressure: 304.579 psig, max 30 sec.  
Flow Rate: max 200 ml/min  
Flow Direction: Left to right

**SCV07052 (Clear ABS and Blue SAN w/Silicone Diaphragm)**

Check Valve, Pocket for .122" (3.1 mm) OD Tubing to Male Locking Luer



Cracking Pressure: 1.450-4.351 psig  
Max. Back Pressure: 87.023 psig  
Flow Rate: max 200 ml/min  
Flow Direction: Left to right

**SCV09093 (Clear and White Polycarbonate w/Silicone Diaphragm)**

Check Valve, Female Locking Luer to Pocket for .118" (3.0 mm) OD Tubing



Cracking Pressure: 1.450-4.351 psig  
Max. Back Pressure: 87.023 psig  
Flow Rate: max 200 ml/min  
Flow Direction: Left to right

**SCV11095 (Clear and Green Polycarbonate w/Silicone Diaphragm)**

Check Valve, Female Locking Luer to Pocket for .098" (2.5 mm) OD Tubing



Cracking Pressure: 1.450-4.351 psig  
Max. Back Pressure: 87.023 psig  
Flow Rate: max 200 ml/min  
Flow Direction: Left to right

**SCV12052 (White and Blue SAN w/Silicone Diaphragm)**

Check Valve, Female Locking Luer to Male Locking Luer



Cracking Pressure: 1.450-4.351 psig  
Max. Back Pressure: 87.023 psig  
Flow Rate: max 200 ml/min  
Flow Direction: Left to right

**SCV15050 (Clear ABS and Clear SAN w/Silicone Diaphragm)**

Check Valve, Pocket for .165" (4.2 mm) OD Tubing to Male Locking Luer



Cracking Pressure: 1.450-4.351 psig  
Max. Back Pressure: 87.023 psig  
Flow Rate: max 200 ml/min  
Flow Direction: Left to right

**SCV16090 (Clear Polycarbonate w/Silicone Diaphragm)**

Check Valve, Female Locking Luer to Pocket for .118" (3.0 mm) OD Tubing



Cracking Pressure: 1.450-4.351 psig  
Max. Back Pressure: 87.023 psig  
Flow Rate: max 200 ml/min  
Flow Direction: Left to right