

Material Suffix Key

-1 White Nylon	-10 Acetal	-2 Black Nylon	-3 Red Nylon
-4 Green Nylon	-5 Blue Nylon	-06 Orange Nylon	-07 Yellow Nylon
-6 Natural Polypropylene	-6005 Animal-Free Polypropylene	-8003 Natural ABS	-81 Lace White ABS
-9 Clear Polycarbonate	-9002 Radiation Stable Polycarbonate	-J1A Natural Kynar PVDF	

MTLLP (-1 -2 -3 -07 -6 -6005 -9 -J1A)

Male Luer Integral Lock Ring Plug, Closed at Luer Tip



LP4 (-1 -2 -3 -4 -5 -6 -9 -9002 -J1A)

Male Luer Integral Lock Ring Plug, Closed at Grip



VLP4 (-1 -3 -6 -81 -J1A)

Male Luer Integral Lock Ring Vented Plug, Closed at Grip



MLRLP (-1 -6 -9)

Male Luer Plug (May be used with separate rotating lock ring; FSLLR)



MLSLP (-1 -6 -9)

Male Luer Plug (May be used with separate stationary lock ring; FSLLR)



MTLP (-1 -6 -9 -J1A)

Male Luer Plug, Closed at Tether Loop (May be used with separate stationary lock ring; MLRL)



MTLSP (-1 -6 -81 -9 -J1A)

Male Luer Slip Plug, Closed at Grip



FTLP (-1 -6 -9 -J1A)

Female Luer Lug Style Cap with Tether Loop



FTLLP (-1 -2 -3 -5 -07 -6 -6005 -81 -9 -9002 -J1A)

Female Luer Thread Style Cap



LP34 (-1 -3 -5 -6 -9 -9002 -J1A)

Male Luer Integral Lock Ring to Female Luer Thread Style Plug



LUER TO THREAD:

Male

ABSML-6 UNDER (-1 -2 -3 -4 -5 -06 -07)

Male Luer with Rotating Lock Ring and 5/16" Hex to 1/4-28 UNF Bottom Sealing Thread



ABSML-J1A UNDER (-1 -2 -3 -4 -5 -06 -07)

Male Luer with Rotating Lock Ring and 5/16" Hex to 1/4-28 UNF Bottom Sealing Thread



XMTLL (-1 -6 -9 -J1A)

Male Luer Integral Lock Ring to 10-32 Special Tapered Thread



SMTLL (-1 -6 -9 -J1A)

Male Luer Integral Lock Ring to 1/4-28 UNF Thread



18MTLL (-1 -2 -6 -J1A)

Male Luer Integral Lock Ring with 7/16" Hex to 1/8-27 NPT Thread



VALUE PLASTICS, INC.

HRM *plastic* +61(0)3 9762 2034
EChnology Pty Ltd

Australian Distributors:
Importers & Manufacturers