

# BLOOD PRESSURE



**BPF SERIES**  
PG-30



**BAYONET MALE CUFF CONNECTORS & DINAMAP® MONITOR ADAPTERS**  
PG-31



**BPLF 1/4 TURN CUFF CONNECTORS**  
PG-32



**DINAMAP® CUFF CONNECTORS**  
PG-33



## THE **BPF** SERIES

From Value Plastics, the reference standard for healthcare fluid management products: The BPF Series bayonet-style female quick connector. BPF series quick connect fittings are designed to withstand abuse, provide economy and ease of use over other alternatives, and comply with ISO Standard BS EN 1060-3.

### **Ergonomic Shape**

- Easy to hold, even with gloves
- Won't irritate patients in case of incidental contact

### **Large Guarded Button**

- Easy to use, even with gloves
- Easy to understand, convenient operation
- Protected against accidental disconnect

### **Push To Connect**

- Easiest connector to use, even with gloves

### **O-Ring Seal**

- Never gets more difficult to put together
- Consistent performance throughout product life

### **Bayonet Rotates In Connector**

- Prevents kinked tubes
- Facilitates easy connections to patient cuff

### **Unique Spring-Loaded Bayonet Retainer**

- Clicks when assembled to signal secure connection
- Resists accidental disconnects

### **Rugged Polymer Construction**

- Withstands tens of thousands of cycles
- Withstands hospital sanitizing procedures

### **Unimpeded Flow Path**

- Better acoustics for accurate readings



VALUE PLASTICS, INC.

**HRM** *aiytic* +61(0)3 9762 2034  
**ECH** *nology* Pty Ltd

Australian Distributors:  
Importers & Manufacturers