

Material Suffix Key

- | | | | |
|--------------------------|--------------------------------------|--------------------------------------|-------------------------|
| -1 White Nylon | -10 Acetal | -2 Black Nylon | -3 Red Nylon |
| -4 Green Nylon | -5 Blue Nylon | -06 Orange Nylon | -07 Yellow Nylon |
| -6 Natural Polypropylene | -6005 Animal-Free Polypropylene | -8003 Natural ABS | -81 Lace White ABS |
| -9 Clear Polycarbonate | -9002 Radiation Stable Polycarbonate | -9010 Radiation Stable Polycarbonate | -J1A Natural Kynar PVDF |

N007/004 (-1 -2 -6 -9 -J1A)
Straight Through Reduction Connector with 500 Series Barbs, 3/32" (2.4 mm) and 1/16" (1.6 mm) ID Tubing



N220/210 (-1 -2 -6 -9 -J1A)
Straight Through Reduction Connector with 200 Series Barbs, 3/32" (2.4 mm) and 1/16" (1.6 mm) ID Tubing



N420/410 (-1 -2 -6 -10 -J1A)
Straight Through Reduction Connector with 400 Series Barbs, 3/32" (2.4 mm) and 1/16" (1.6 mm) ID Tubing



AD (-1 -2 -6 -9 -J1A)
Straight Through Reduction Connector with Classic Series Barbs, 3/32" (2.4 mm) and 1/16" (1.6 mm) ID Tubing



N013/004 (-1 -2 -6 -9 -J1A)
Straight Through Reduction Connector with 500 Series Barbs, 1/8" (3.2 mm) and 1/16" (1.6 mm) ID Tubing



N230/210 (-1 -2 -6 -9 -J1A)
Straight Through Reduction Connector with 200 Series Barbs, 1/8" (3.2 mm) and 1/16" (1.6 mm) ID Tubing



N430/410 (-1 -2 -6 -J1A)
Straight Through Reduction Connector with 400 Series Barbs, 1/8" (3.2 mm) and 1/16" (1.6 mm) ID Tubing



AC (-0 -1 -2 -6 -9 -J1A)
Straight Through Reduction Connector with Classic Series Barbs, 1/8" (3.2 mm) and 1/16" (1.6 mm) ID Tubing



N025/004 (-1 -2 -6 -6005 -9 -J1A)
Straight Through Reduction Connector with 500 Series Barbs, 5/32" (4.0 mm) and 1/16" (1.6 mm) ID Tubing




T040 (-1 -2 -6 -9 -J1A)
Straight Through Reduction Connector with Classic Series Barbs, 5/32" (4.0 mm) and 1/16" (1.6 mm) ID Tubing



I050 (-1 -2 -6 -9 -J1A)
Straight Through Reduction Connector with Classic Series Barbs, 3/16" (4.8 mm) and 1/16" (1.6 mm) ID Tubing



N013/007 (-1 -2 -6 -9 -9002 -J1A)
Straight Through Reduction Connector with 500 Series Barbs, 1/8" (3.2 mm) and 3/32" (2.4 mm) ID Tubing



N230/220 (-1 -2 -6 -9 -J1A)
Straight Through Reduction Connector with 200 Series Barbs, 1/8" (3.2 mm) and 3/32" (2.4 mm) ID Tubing



N430/420 (-1 -2 -6 -J1A)
Straight Through Reduction Connector with 400 Series Barbs, 1/8" (3.2 mm) and 3/32" (2.4 mm) ID Tubing



DC (-1 -2 -6 -9)
Straight Through Reduction Connector with Classic Series Barbs, 1/8" (3.2 mm) and 3/32" (2.4 mm) ID Tubing



N025/007 (-1 -2 -6 -6005 -9 -J1A)
Straight Through Reduction Connector with 500 Series Barbs, 5/32" (4.0 mm) and 3/32" (2.4 mm) ID Tubing



T040 (-1 -2 -6 -9 -J1A)
Straight Through Reduction Connector with Classic Series Barbs, 5/32" (4.0 mm) and 3/32" (2.4 mm) ID Tubing



N035/007 (-1 -2 -6 -9 -J1A)
Straight Through Reduction Connector with 500 Series Barbs, 3/16" (4.8 mm) and 3/32" (2.4 mm) ID Tubing

