

Material Suffix Key

- | | | | |
|--------------------------|--------------------------------------|--------------------------------------|-------------------------|
| -1 White Nylon | -10 Acetal | -2 Black Nylon | -3 Red Nylon |
| -4 Green Nylon | -5 Blue Nylon | -06 Orange Nylon | -07 Yellow Nylon |
| -6 Natural Polypropylene | -6005 Animal-Free Polypropylene | -8003 Natural ABS | -81 Lace White ABS |
| -9 Clear Polycarbonate | -9002 Radiation Stable Polycarbonate | -9010 Radiation Stable Polycarbonate | -J1A Natural Kynar PVDF |

T013 (-6005)

Tee Connector with 500 Series Barbs, 1/8" (3.2 mm) ID Tubing




T30 (-1 -2 -6 -9 -J1A)

Tee Connector with Classic Series Barbs, 1/8" (3.2 mm) ID Tubing




T40 (-1 -2 -6 -J1A)

Tee Connector with Classic Series Barbs, 5/32" (4 mm) ID Tubing




T50 (-1 -2 -6 -J1A)

Tee Connector with Classic Series Barbs, 3/16" (4.8 mm) ID Tubing




T055 (-6005)

Tee Connector with 500 Series Barbs, 1/4" (6.4 mm) ID Tubing



T60 (-1 -2 -6 -J1A)

Tee Connector with Classic Series Barbs, 1/4" (6.4 mm) ID Tubing



T210/220 (-1 -2 -6 -9 -J1A)

Tee Reduction Connector with 200 Series Barbs, 1/16" (1.6 mm) and 3/32" (2.4 mm) ID Tubing




T410/420 (-1 -2 -6 -J1A)

Tee Reduction Connector with 400 Series Barbs, 1/16" (1.6 mm) and 3/32" (2.4 mm) ID Tubing



T420/410 (-1 -2 -6 -J1A)

Tee Reduction Connector with 400 Series Barbs, 3/32" (2.4 mm) and 1/16" (1.6 mm) ID Tubing



T220/210 (-1 -2 -6 -9 -J1A)

Tee Reduction Connector with 200 Series Barbs, 3/32" (2.4 mm) and 1/16" (1.6 mm) ID Tubing



T220R210 (-1 -2 -6 -9 -J1A)

Tee Reduction Connector with 200 Series Barbs, 3/32" (2.4 mm) and 1/16" (1.6 mm) ID Tubing



T420R410 (-1 -2 -6 -J1A)

Tee Reduction Connector with 400 Series Barbs, 3/32" (2.4 mm) and 1/16" (1.6 mm) ID Tubing



T430/410 (-1 -2 -6 -J1A)

Tee Reduction Connector with 400 Series Barbs, 1/8" (3.2 mm) and 1/16" (1.6 mm) ID Tubing



T230/210 (-1 -2 -6 -9 -J1A)

Tee Reduction Connector with 200 Series Barbs, 1/8" (3.2 mm) and 1/16" (1.6 mm) ID Tubing

