

LARGE TUBE-TO-TUBE



STRAIGHT THROUGH
PG-53



**STRAIGHT THROUGH
REDUCER**
PG-54



Y
PG-55



PRESS-IN PLUGS
PG-55



ELBOW
PG-56



ELBOW REDUCER
PG-56



TEE
PG-57



TEE REDUCER
PG-57



LARGE TUBE-TO-TUBE

Large tube-to-tube fittings encompass sizes fitting tubing IDs between 3/8" (9.5 mm) and 1.0" (25.4 mm). They are available from stock in 3 barb styles, most commonly in animal-free and natural polypropylene and polycarbonate resins, with PVDF available as an option and other materials available on special order.

AVAILABLE BARB STYLES

Tube ID	500 Series	600 Series	900 Series
3/8" (9.5 mm)	X	X	X
1/2" (12.7 mm)	X	X	X
5/8" (16.0 mm)	X	X	
3/4" (19.0 mm)	X	X	
1" (25.4 mm)	X	X	



VALUE PLASTICS, INC.

HRM *alytic* +61(0)3 9762 2034
ECH *nology* Pty Ltd

Australian Distributors:
Importers & Manufacturers