

## LARGE TUBE-TO-TUBE: Straight Through

### Material Suffix Key

-0 Natural Nylon  
-J1A Natural Kynnar PVDF


-6 Natural Polypropylene

-6005 Animal-Free Polypropylene

-9002 Radiation Stable Polycarbonate



**SHOWN: N670**



**SHOWN: N110**

Minimum package quantities vary by tube IDs, please verify minimum package sizes when placing orders

**N070 (-0 -6 -9)**

Straight Through Connector with 500 Series Barbs, 3/8" (9.5 mm) ID Tubing

**N090 (-0 -6 -9)**

Straight Through Connector with 500 Series Barbs, 5/8" (16.0 mm) ID Tubing

**N670 (-6 -6005 -9 -J1A)**

Straight Through Connector with 600 Series Barbs, 3/8" (9.5 mm) ID Tubing

**N690 (-0 -6 -9)**

Straight Through Connector with 600 Series Barbs, 5/8" (16.0 mm) ID Tubing

**N970 (-9)**

Straight Through Connector with 900 Series Barbs, 3/8" (9.5 mm) ID Tubing

**N100 (-0 -6)**

Straight Through Connector with 500 Series Barbs, 3/4" (19.0 mm) ID Tubing

**N080 (-0 -6 -9)**

Straight Through Connector with 500 Series Barbs, 1/2" (12.7 mm) ID Tubing

**N110 (-0 -6)**

Straight Through Connector with 500 Series Barbs, 1" (25.4 mm) ID Tubing

**N680 (-6 -9 -J1A)**

Straight Through Connector with 600 Series Barbs, 1/2" (12.7 mm) ID Tubing

**N975 (-9)**

Straight Through Connector with 900 Series Barbs, 1/2" (12.7 mm) ID Tubing



our online ordering for component parts, now available 24 hours a day, 7 days a week, with no \$ minimums! You can find every catalog part at:

**www.valueplastics.com**

(Online ordering is not available in countries served by our regional partners)

Website : [www.chromtech.net.au](http://www.chromtech.net.au) E-mail : [info@chromtech.net.au](mailto:info@chromtech.net.au) TelNo : 03 9762 2034 . . . in AUSTRALIA

100 piece minimums, unless noted in red text • ■ Red Text = Special order items • Volume pricing may be available for quantities 5,000 and over • Call for details!



VALUE PLASTICS, INC.

HRM alytic +61(0)3 9762 2034  
ECHnology Pty Ltd

Australian Distributors:  
Importers & Manufacturers