

THREAD:**10-32 UNF, Special Tapered Thread****Material Suffix Key**

-1 White Nylon
-4 Green Nylon
-6 Natural Polypropylene
-9 Clear Polycarbonate

-10 Acetal
-5 Blue Nylon
-6005 Animal-Free Polypropylene
-9002 Radiation Stable Polycarbonate

-2 Black Nylon
-06 Orange Nylon
-8003 Natural ABS
-J1A Natural Kynar PVDF

-3 Red Nylon
-07 Yellow Nylon
-81 Lace White ABS

X210 (-1 -2 -3 -5 -6 -J1A)

10-32 Special Tapered Thread with 1/4" Hex to 200 Series Barb, 1/16" (1.6 mm) ID Tubing

**K210 (-1 -2 -3 -5 -6 -J1A)**

10-32 UNF Thread with 1/4" Hex to 200 Series Barb, 1/16" (1.6 mm) ID Tubing

**X410 (-1 -2 -6 -J1A)**

10-32 Special Tapered Thread with 1/4" Hex to 400 Series Barb, 1/16" (1.6 mm) ID Tubing

**K410 (-1 -2 -6 -J1A)**

10-32 UNF Thread with 1/4" Hex to 400 Series Barb, 1/16" (1.6 mm) ID Tubing

**X10 (-1 -2 -6 -J1A)**

10-32 Special Tapered Thread with 1/4" Hex to Classic Series Barb, 1/16" (1.6 mm) ID Tubing

**K10 (-1 -2 -6 -J1A)**

10-32 UNF Thread with 1/4" Hex to Classic Series Barb, 1/16" (1.6 mm) ID Tubing

**X220 (-1 -2 -3 -6 -J1A)**

10-32 Special Tapered Thread with 1/4" Hex to 200 Series Barb, 3/32" (2.4 mm) ID Tubing

**K220 (-1 -2 -6 -J1A)**

10-32 UNF Thread with 1/4" Hex to 200 Series Barb, 3/32" (2.4 mm) ID Tubing

**X420 (-1 -2 -6 -J1A)**

10-32 Special Tapered Thread with 1/4" Hex to 400 Series Barb, 3/32" (2.4 mm) ID Tubing

**K420 (-1 -2 -6 -J1A)**

10-32 UNF Thread with 1/4" Hex to 400 Series Barb, 3/32" (2.4 mm) ID Tubing

**D (-1 -2 -3 -4 -5 -6 -J1A)**

10-32 UNF Thread with 1/4" Hex to Classic Series Barb, 3/32" (2.4 mm) ID Tubing

**X230 (-1 -2 -3 -5 -6 -J1A)**

10-32 Special Tapered Thread with 1/4" Hex to 200 Series Barb, 1/8" (3.2 mm) ID Tubing

**K230 (-1 -2 -6 -J1A)**

10-32 UNF Thread with 1/4" Hex to 200 Series Barb, 1/8" (3.2 mm) ID Tubing

**X430 (-1 -2 -3 -6 -J1A)**

10-32 Special Tapered Thread with 1/4" Hex to 400 Series Barb, 1/8" (3.2 mm) ID Tubing

**K430 (-1 -2 -6 -J1A)**

10-32 UNF Thread with 1/4" Hex to 400 Series Barb, 1/8" (3.2 mm) ID Tubing

**B (-1 -2 -3 -4 -5 -6 -J1A)**

10-32 UNF Thread with 1/4" Hex to Classic Series Barb, 1/8" (3.2 mm) ID Tubing

**C (-1 -2 -3 -4 -5 -6 -J1A)**

10-32 UNF Thread with 1/4" Hex to Classic Series Barb, Short Neck, 1/8" (3.2 mm) ID Tubing

**K240 (-1 -2 -6 -J1A)**

10-32 UNF Thread with 1/4" Hex to 200 Series Barb, 5/32" (4.0 mm) ID Tubing

**X240 (-1 -2 -6 -J1A)**

10-32 Special Tapered Thread with 1/4" Hex to 200 Series Barb, 5/32" (4.0 mm) ID Tubing

**X445 (-1 -2 -6 -J1A)**

10-32 Special Tapered Thread with 1/4" Hex to 400 Series Barb, .170" (4.3 mm) ID Tubing

**K445 (-1 -2 -6 -J1A)**

10-32 UNF Thread with 1/4" Hex to 400 Series Barb, .170" (4.3 mm) ID Tubing



Website : www.chromtech.net.au E-mail : info@chromtech.net.au TelNo : 03 9762 2034 . . . in AUSTRALIA

100 piece minimums, unless noted in red text • ■ Red Text = Special order items • Volume pricing may be available for quantities 5,000 and over • Call for details!



VALUE PLASTICS, INC.

CHROM *metric* +61(0)3 9762 2034
EChnology Pty Ltd

Australian Distributors:
Importers & Manufacturers