



BOTTLE THREAD ADAPTER

- Adapt GL45 caps to fit GL38 solvent bottles
- Polypropylene construction

In lab applications where various bottle types are used, there may be occasions when the use of separate bottle caps is inconvenient. Save time and hassle by using the new BC-151 solvent bottle adapter, which will convert Vaplock 'Economy' GL45 bottle caps to standard 38-430 threaded solvent bottles.

BOTTLE THREAD ADAPTER

BC-151	ADAPTER, GL45 to GL38, Polypropylene	\$42.00
---------------	--------------------------------------	---------

Nut and Ferrule Assembly, 1/4-28

GL45 Vaplock Economy Cap BC-322N

BC-151 Bottle Cap Adapter

GL38 (or 38-430) glass reagent bottle



Vaplock™ is a trademark of Western Fluids Engineering & Mfg, LLC
 © 2007 Western Fluids Engineering and Manufacturing, LLC. All rights reserved.