

VAPLOCKTM
CLOSED SYSTEMS

SOLVENT RESERVOIRS

FOR STORAGE & CONTAINMENT

BOROSILICATE GLASS BOTTLES FOR LABORATORY SOLVENTS






SOLVENT RESERVOIRS

- GL45 threaded lab bottles
- Borosilicate glass, available in clear or amber
- Select from clear uncoated glass, clear with protective mesh, or plastic safety coated

Although many container types and sizes are used in the laboratory, one of the more commonly used standards is the "GL45" threaded bottle. The nomenclature is simple: "GL" denotes a glass bottle, and "45" is the thread dimension in millimeters when measured across the outermost diameter. This laboratory standard is available for purchase in clear or amber, in a variety of sizes and protective coverings. Pressure-tested 1 liter bottles are also available as an added safety measure.



GL45 BOTTLE
Measures 45 mm from outside thread to outside thread (drawing not to scale)

1 BAR PRESSURE-TESTED BOTTLES, GL45, 1 LITER

BT-4506	CLEAR, 1 Liter, pressure-tested	\$68.77
BT-4512	CLEAR 1 Liter, pressure-tested, plastic coated	\$120.19

GL45 STANDARD LAB BOTTLES

100 ML		
BT45-01	CLEAR 100mL	\$7.75
BT45-01M	CLEAR 100mL, with protective mesh	\$11.95
BT45-01P	CLEAR 100mL, plastic safety coated	\$24.25
BT45-01A	AMBER 100mL, non-coated	\$34.00
BT45-01AM	AMBER 100mL, with protective mesh	\$38.00
250 ML		
BT45-025	CLEAR 250mL	\$9.50
BT45-025M	CLEAR 250mL, with protective mesh	\$13.50
BT45-025P	CLEAR 250mL, plastic safety coated	\$30.50
BT45-025A	AMBER 250mL, non-coated	\$50.00
BT45-025AM	AMBER 250mL, with protective mesh	\$54.00
500 ML		
BT45-05	CLEAR 500mL	\$13.00
BT45-05M	CLEAR 500mL, with protective mesh	\$17.00
BT45-05P	CLEAR 500mL, plastic safety coated	\$37.00
BT45-05A	AMBER 500mL, non-coated	\$64.00
BT45-05AM	AMBER 500mL, with protective mesh	\$68.00
1 LITER		
BT45-1	CLEAR 1 Liter	\$20.00
BT45-1M	CLEAR 1 Liter, with protective mesh	\$25.00
BT45-1P	CLEAR 1 Liter, plastic safety coated	\$58.00
BT45-1A	AMBER 1 Liter, non-coated	\$93.00
BT45-1AM	AMBER 1 Liter, with protective mesh	\$98.00
2 LITER		
BT45-2	CLEAR 2 Liter	\$52.00
BT45-2M	CLEAR 2 Liter, with protective mesh	\$57.00
BT45-2P	CLEAR 2 Liter, plastic safety coated	\$96.00
BT45-2A	AMBER 2 Liter, non-coated	\$127.00
BT45-2AM	AMBER 2 Liter, with protective mesh	\$132.00
5 LITER		
BT45-5	CLEAR 5 Liter	\$220.00
BT45-5M	CLEAR 5 Liter, with protective mesh	\$225.00
BT45-5P	CLEAR 5 Liter, plastic safety coated	\$310.00
BT45-5A	AMBER 5 Liter, non-coated	\$370.00
BT45-5AM	AMBER 5 Liter, with protective mesh	\$375.00
10 LITER		
BT45-10	CLEAR 10 Liter	\$380.00
BT45-10M	CLEAR 10 Liter, with protective mesh	\$385.00
BT45-10P *	CLEAR 10 Liter, plastic safety coated	\$470.00
BT45-10A	AMBER 10 Liter, non-coated	\$695.00
BT45-10AM	AMBER 10 Liter, with protective mesh	\$700.00
20 LITER		
BT45-20	CLEAR 20 Liter	\$910.00
BT45-20M	CLEAR 20 Liter, with protective mesh	\$915.00
BT45-20A	AMBER 20 Liter, non-coated	\$1,595.00
BT45-20AM	AMBER 20 Liter, with protective mesh	\$1,600.00

* SPECIAL ORDER ONLY

USD Prices for AUD Enquire

PRICING SUBJECT TO CHANGE WITHOUT NOTICE

