

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

<u>Catalog No</u>	<u>Description</u>	<u>Category</u>	<u>AUD 1.00</u>	<u>Web Link</u>
T10C10D	1/16inOD x0.010inID x10cm 2-pk	TU - Stainless Steel	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T30C5-50	1/16inOD x0.005inID x30cm 50-pk	TU - Stainless Steel	\$332.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T100C40-50	1/16inOD x0.040inID x100cm 50-pk	TU - Stainless Steel	\$624.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T100C30-50	1/16inOD x0.030inID x100cm 50-pk	TU - Stainless Steel	\$624.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T100C20-50	1/16inOD x0.020inID x100cm 50-pk	TU - Stainless Steel	\$624.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T100C10-50	1/16inOD x0.010inID x100cm 50-pk	TU - Stainless Steel	\$624.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T50C40-50	1/16inOD x0.040inID x50cm 50-pk	TU - Stainless Steel	\$382.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T50C30-50	1/16inOD x0.030inID x50cm 50-pk	TU - Stainless Steel	\$382.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T50C20-50	1/16inOD x0.020inID x50cm 50-pk	TU - Stainless Steel	\$382.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T50C10-50	1/16inOD x0.010inID x50cm 50-pk	TU - Stainless Steel	\$382.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T50C5-50	1/16inOD x0.005inID x50cm50-pk	TU - Stainless Steel	\$382.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T30C40-50	1/16inOD x0.040inID x30cm Straight SS Tubing 50-pk	TU - Stainless Steel	\$332.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T30C30-50	1/16inOD x0.030inID x30cm Straight SS Tubing 50-pk	TU - Stainless Steel	\$332.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T5C40D	1/16inOD x0.040inID x5cm 2-pk	TU - Stainless Steel	\$10.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T30C10-50	1/16inOD x0.010inID x30cm	TU - Stainless Steel	\$332.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T5C20-100	1/16inOD x0.020inID x5cm 100-pk	TU - Stainless Steel	\$402.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T20C40-50	1/16inOD x0.030inID x20cm 50-pk	TU - Stainless Steel	\$306.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T20C30-50	1/16inOD x0.030inID x20cm 50-pk	TU - Stainless Steel	\$306.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T20C20-50	1/16inOD x0.020inID x20cm	TU - Stainless Steel	\$306.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T20C10-50	1/16inOD x0.010inID x20cm	TU - Stainless Steel	\$306.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T20C5-50	1/16inOD x0.005inID x20cm 50-pk	TU - Stainless Steel	\$306.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T10C40-50	1/16inOD x0.040inID x10cm 50-pk	TU - Stainless Steel	\$243.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T10C30-50	1/16inOD x0.030inID x10cm 50-pk	TU - Stainless Steel	\$243.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T10C20-50	1/16inOD x0.020inID x10cm 50-pk	TU - Stainless Steel	\$243.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T10C10-50	1/16inOD x0.010inID x10cm 50-pk	TU - Stainless Steel	\$243.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T10C5-50	1/16inOD x0.005inID x10cm	TU - Stainless Steel	\$243.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T5C40-50	1/16inOD x0.040inID x5cm 50-pk	TU - Stainless Steel	\$216.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T5C30-50	1/16inOD x0.030inID x5cm 50-pk	TU - Stainless Steel	\$216.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T30C20-50	1/16inOD x0.020inID x30cm	TU - Stainless Steel	\$332.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T20C40-100	1/16inOD x0.030inID x20cm 100-pk	TU - Stainless Steel	\$561.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T100C40-100	1/16inOD x0.040inID x100cm 100-pk	TU - Stainless Steel	\$1,198.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T100C30-100	1/16inOD x0.030inID x100cm 100-pk	TU - Stainless Steel	\$1,198.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T100C20-100	1/16inOD x0.020inID x100cm 100-pk	TU - Stainless Steel	\$1,198.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T100C10-100	1/16inOD x0.010inID x100cm 100-pk	TU - Stainless Steel	\$1,198.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

T50C40-100	1/16inOD x0.040inID x50cm 100-pk	TU - Stainless Steel	\$720.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T50C30-100	1/16inOD x0.030inID x50cm 100-pk	TU - Stainless Steel	\$720.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T50C20-100	1/16inOD x0.020inID x50cm 100-pk	TU - Stainless Steel	\$720.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T50C10-100	1/16inOD x0.010inID x50cm 100-pk	TU - Stainless Steel	\$720.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T50C5-100	1/16inOD x0.005inID x50cm 100-pk	TU - Stainless Steel	\$720.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T30C40-100	1/16inOD x0.040inID x30cm Straight SS Tubing 100-pk	TU - Stainless Steel	\$606.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T30C30-100	1/16inOD x0.030inID x30cm Straight SS Tubing 100-pk	TU - Stainless Steel	\$606.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T30C20-100	1/16inOD x0.020inID x30cm	TU - Stainless Steel	\$606.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T5C5-100	1/16inOD x0.005inID x5cm 100-pk	TU - Stainless Steel	\$402.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T30C5-100	1/16inOD x0.005inID x30cm 100-pk	TU - Stainless Steel	\$606.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T5C10-100	1/16inOD x0.010inID x5cm 100-pk	TU - Stainless Steel	\$402.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T20C30-100	1/16inOD x0.030inID x20cm 100-pk	TU - Stainless Steel	\$561.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T20C20-100	1/16inOD x0.020inID x20cm	TU - Stainless Steel	\$561.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T20C10-100	1/16inOD x0.010inID x20cm	TU - Stainless Steel	\$561.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T20C5-100	1/16inOD x0.005inID x20cm 100-pk	TU - Stainless Steel	\$561.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T10C40-100	1/16inOD x0.040inID x10cm 100-pk	TU - Stainless Steel	\$453.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T10C30-100	1/16inOD x0.030inID x10cm 100-pk	TU - Stainless Steel	\$453.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T10C20-100	1/16inOD x0.020inID x10cm 100-pk	TU - Stainless Steel	\$453.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T10C10-100	1/16inOD x0.010inID x10cm 100-pk	TU - Stainless Steel	\$453.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T10C5-100	1/16inOD x0.005inID x10cm	TU - Stainless Steel	\$453.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T5C40-100	1/16inOD x0.040inID x5cm 100-pk	TU - Stainless Steel	\$402.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T5C30-100	1/16inOD x0.030inID x5cm 100-pk	TU - Stainless Steel	\$402.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T5C5-50	1/16inOD x0.040inID x5cm 10-pk	TU - Stainless Steel	\$216.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T30C10-100	1/16inOD x0.010inID x30cm	TU - Stainless Steel	\$606.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T30C5D	1/16inOD x0.005inID x30cm Straight SS Tubing 2-pk	TU - Stainless Steel	\$15.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T5C20-50	1/16inOD x0.020inID x5cm50-pk	TU - Stainless Steel	\$216.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T100C40D	1/16inOD x0.040inID x100cm 2-pk	TU - Stainless Steel	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T100C30D	1/16inOD x0.030inID x100cm 2-pk	TU - Stainless Steel	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T100C20D	1/16inOD x0.020inID x100cm Straight SS Tubing 2-pk	TU - Stainless Steel	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T100C10D	1/16inOD x0.010inID x100cm Straight SS Tubing	TU - Stainless Steel	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T50C40D	1/16inOD x0.040inID x50cm 2-pk	TU - Stainless Steel	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T50C30D	1/16inOD x0.030inID x50cm 2-pk	TU - Stainless Steel	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T50C20D	1/16inOD x0.020inID x50cm Straight SS Tubing 2-pk	TU - Stainless Steel	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

T50C10D	1/16inOD x0.010inID x50cm Straight SS Tubing 2-pk	TU - Stainless Steel	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T50C5D	1/16inOD x0.005inID x50cm Straight SS Tubing 2-pk	TU - Stainless Steel	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T30C40D	1/16inOD x0.040inID x30cm Straight SS Tubing 2-pk	TU - Stainless Steel	\$15.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T30C30D	1/16inOD x0.030inID x30cm Straight SS Tubing 2-pk	TU - Stainless Steel	\$15.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T5C10-10	1/16inOD x0.010inID x5cm 10-pk	TU - Stainless Steel	\$48.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T30C10D	1/16inOD x0.010inID x20cm Straight SS Tubing	TU - Stainless Steel	\$15.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T5C20-10	1/16inOD x0.020inID x5cm 10-pk	TU - Stainless Steel	\$48.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T20C40D	1/16inOD x0.030inID x20cm 2-pk	TU - Stainless Steel	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T20C30D	1/16inOD x0.030inID x20cm 2-pk	TU - Stainless Steel	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T20C20D	1/16inOD x0.020inID x20cm	TU - Stainless Steel	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T20C10D	1/16inOD x0.010inID x20cm	TU - Stainless Steel	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T20C5D	1/16inOD x0.005inID x20cm 2-pk	TU - Stainless Steel	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T10C40D	1/16inOD x0.040inID x10cm 2-pk	TU - Stainless Steel	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T10C30D	1/16inOD x0.030inID x10cm 2-pk	TU - Stainless Steel	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T10C20D	1/16inOD x0.020inID x10cm 2-pk	TU - Stainless Steel	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T10C5D	1/16inOD x0.005inID x10cm	TU - Stainless Steel	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T5C30D	1/16inOD x0.030inID x5cm 2-pk	TU - Stainless Steel	\$10.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T5C20D	1/16inOD x0.020inID x5cm 2-pk	TU - Stainless Steel	\$10.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T5C10D	1/16inOD x0.010inID x5cm 2-pk	TU - Stainless Steel	\$10.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T30C20D	1/16inOD x0.020inID x20cm Straight SS Tubing 2-pk	TU - Stainless Steel	\$15.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T30C5-10	1/16inOD x0.005inID x30cm Straight SS Tubing	TU - Stainless Steel	\$70.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T5C5D	1/16inOD x0.005inID x5cm	TU - Stainless Steel	\$10.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T100C40-10	1/16inOD x0.040inID x100cm 10-pk	TU - Stainless Steel	\$128.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T100C30-10	1/16inOD x0.030inID x100cm 10-pk	TU - Stainless Steel	\$128.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T100C20-10	1/16inOD x0.020inID x100cm 2-pk	TU - Stainless Steel	\$128.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T100C10-10	1/16inOD x0.010inID x100cm 10-pk	TU - Stainless Steel	\$128.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T50C40-10	1/16inOD x0.040inID x50cm 10-pk	TU - Stainless Steel	\$84.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T50C30-10	1/16inOD x0.030inID x50cm 10-pk	TU - Stainless Steel	\$84.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T50C20-10	1/16inOD x0.020inID x50cm10-pk	TU - Stainless Steel	\$84.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T50C10-10	1/16inOD x0.010inID x50cm 10-pk	TU - Stainless Steel	\$84.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T50C5-10	1/16inOD x0.005inID x50cm Straight SS Tubing 2-pk	TU - Stainless Steel	\$84.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T30C40-10	1/16inOD x0.040inID x30cm Straight SS Tubing 10-pk	TU - Stainless Steel	\$70.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T30C30-10	1/16inOD x0.030inID x30cm Straight SS Tubing 10-pk	TU - Stainless Steel	\$70.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

T5C5-10	1/16inOD x0.005inID x5cm 10-pk	TU - Stainless Steel	\$48.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T30C10-10	1/16inOD x0.010inID x30cm	TU - Stainless Steel	\$70.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T5C10-50	1/16inOD x0.010inID x5cm 50-pk	TU - Stainless Steel	\$216.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T20C40-10	1/16inOD x0.040inID x20cm 10-pk	TU - Stainless Steel	\$62.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T20C30-10	1/16inOD x0.030inID x20cm 10-pk	TU - Stainless Steel	\$62.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T20C20-10	1/16inOD x0.020inID x20cm	TU - Stainless Steel	\$62.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T20C10-10	1/16inOD x0.010inID x20cm	TU - Stainless Steel	\$62.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T20C5-10	1/16inOD x0.005inID x20cm 10-pk	TU - Stainless Steel	\$62.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T10C40-10	1/16inOD x0.040inID x10cm 10-pk	TU - Stainless Steel	\$51.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T10C30-10	1/16inOD x0.030inID x10cm 10-pk	TU - Stainless Steel	\$51.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T10C20-10	1/16inOD x0.020inID x10cm 10-pk	TU - Stainless Steel	\$51.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T10C10-10	1/16inOD x0.010inID x10cm 10-pk	TU - Stainless Steel	\$51.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T10C5-10	1/16inOD x0.005inID x10cm	TU - Stainless Steel	\$51.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T5C40-10	1/16inOD x0.040inID x5cm 10-pk	TU - Stainless Steel	\$48.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T5C30-10	1/16inOD x0.030inID x5cm 10-pk	TU - Stainless Steel	\$48.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T30C20-10	1/16inOD x0.020inID x30cm	TU - Stainless Steel	\$70.00	<a href="http://www.chromtech.net.au/pdf2/vici09_93.pdf">http://www.chromtech.net.au/pdf2/vici09_93.pdf</a>
T10N20D	1/32 SS Tubing 0.020"ID (2 x 20cm)	Not Categorised	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
THC270	Hastelloy C 1/8inOD x.070inID	TU - Hastelloy C	\$26.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
TINCO261	Inconel 600 1/8inOD x 0.061in	TU - Inconel	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
TINCO282	Inconel 600 1/8inOD x 0.082in	TU - Inconel	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T50N5-100	1/32inOD x0.005inID x50cm	TU - Stainless Steel	\$828.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T100N20-50	1/32inOD x0.020inID x100cm	TU - Stainless Steel	\$714.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T20N5-100	1/32inOD x0.005inID x20cm	TU - Stainless Steel	\$638.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T5N5-100	1/32inOD x0.005inID x5cm 100-pk	TU - Stainless Steel	\$459.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T5N10-100	1/32inOD x0.010inID x5cm	TU - Stainless Steel	\$459.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T100N10-50	1/32inOD x0.010inID x100cm	TU - Stainless Steel	\$714.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T5N20-100	1/32inOD x0.020inID x5cm	TU - Stainless Steel	\$459.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T10N5-100	1/32inOD x0.005inID x10cm	TU - Stainless Steel	\$510.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T10N10-100	1/32inOD x0.010inID x10cm	TU - Stainless Steel	\$510.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T10N20-100	1/32inOD x0.020inID x10cm	TU - Stainless Steel	\$510.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T20N10-100	1/32inOD x0.010inID x20cm	TU - Stainless Steel	\$638.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T20N20-100	1/32inOD x0.020inID x20cm	TU - Stainless Steel	\$638.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T30N5-100	1/32inOD x0.005inID x30cm	TU - Stainless Steel	\$688.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T100N20-100	1/32inOD x0.020inID x100cm	TU - Stainless Steel	\$1,365.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T30N20-100	1/32inOD x0.020inID x30cm	TU - Stainless Steel	\$688.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

T50N20-100	1/32inOD x0.020inID x50cm	TU - Stainless Steel	\$828.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T100N10-100	1/32inOD x0.010inID x100cm	TU - Stainless Steel	\$1,365.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T50N20-50	1/32inOD x0.020inID x50cm	TU - Stainless Steel	\$453.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T30N20D	1/32inOD x0.020inID x20cm	TU - Stainless Steel	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T30N10-100	1/32inOD x0.010inID x30cm	TU - Stainless Steel	\$688.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T100N10D	1/32inOD x0.010inID x100cm	TU - Stainless Steel	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T50N10D	1/32inOD x0.010inID x50cm	TU - Stainless Steel	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T20N10-10	1/32inOD x0.010inID x20cm	TU - Stainless Steel	\$70.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T20N5-10	1/32inOD x0.005inID x20cm	TU - Stainless Steel	\$70.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T10N20-10	1/32inOD x0.020inID x10cm	TU - Stainless Steel	\$58.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T10N10-10	1/32inOD x0.010inID x10cm	TU - Stainless Steel	\$58.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T10N5-10	1/32inOD x0.005inID x10cm	TU - Stainless Steel	\$58.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T5N20-10	1/32inOD x0.020inID x5cm	TU - Stainless Steel	\$54.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T5N10-10	1/32inOD x0.010inID x5cm	TU - Stainless Steel	\$54.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T30N5-10	1/32inOD x0.005inID x30cm	TU - Stainless Steel	\$81.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T100N20D	1/32inOD x0.020inID x100cm	TU - Stainless Steel	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T30N10-10	1/32inOD x0.010inID x30cm	TU - Stainless Steel	\$81.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T50N20D	1/32inOD x0.020inID x50cm	TU - Stainless Steel	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T50N5D	1/32inOD x0.005inID x50cm	TU - Stainless Steel	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T30N10D	1/32inOD x0.010inID x20cm	TU - Stainless Steel	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T30N5D	1/32inOD x0.005inID x20cm	TU - Stainless Steel	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T20N20D	1/32inOD x0.020inID x20cm	TU - Stainless Steel	\$14.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T20N10D	1/32inOD x0.010inID x20cm	TU - Stainless Steel	\$14.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T20N5D	1/32inOD x0.005inID x20cm	TU - Stainless Steel	\$14.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T10N10D	1/32inOD x0.010inID x10cm	TU - Stainless Steel	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T10N5D	1/32inOD x0.005inID x10cm	TU - Stainless Steel	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T5N5-10	1/32inOD x0.005inID x5cm	TU - Stainless Steel	\$54.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T5N20-50	1/32inOD x0.020inID x5cm	TU - Stainless Steel	\$243.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T50N5-50	1/32inOD x0.005inID x50cm	TU - Stainless Steel	\$453.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T30N20-50	1/32inOD x0.020inID x30cm	TU - Stainless Steel	\$376.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T30N10-50	1/32inOD x0.010inID x30cm	TU - Stainless Steel	\$376.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T30N5-50	1/32inOD x0.005inID x30cm	TU - Stainless Steel	\$376.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T20N20-50	1/32inOD x0.020inID x20cm	TU - Stainless Steel	\$345.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T20N10-50	1/32inOD x0.010inID x20cm	TU - Stainless Steel	\$345.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T20N5-50	1/32inOD x0.005inID x20cm	TU - Stainless Steel	\$345.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>

VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

T10N20-50	1/32inOD x0.020inID x10cm	TU - Stainless Steel	\$267.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T20N20-10	1/32inOD x0.020inID x20cm	TU - Stainless Steel	\$70.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T10N5-50	1/32inOD x0.005inID x10cm	TU - Stainless Steel	\$117.00	
T50N10-50	1/32inOD x0.010inID x50cm	TU - Stainless Steel	\$453.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T5N10-50	1/32inOD x0.010inID x5cm	TU - Stainless Steel	\$243.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T5N5-50	1/32inOD x0.005inID x5cm 50-pk	TU - Stainless Steel	\$243.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T100N20-10	1/32inOD x0.020inID x100cm	TU - Stainless Steel	\$147.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T100N10-10	1/32inOD x0.010inID x100cm	TU - Stainless Steel	\$147.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T50N20-10	1/32inOD x0.020inID x50cm	TU - Stainless Steel	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T50N10-10	1/32inOD x0.010inID x50cm	TU - Stainless Steel	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T50N5-10	1/32inOD x0.005inID x50cm	TU - Stainless Steel	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T30N20-10	1/32inOD x0.020inID x30cm	TU - Stainless Steel	\$81.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
T10N10-50	1/32inOD x0.010inID x10cm	TU - Stainless Steel	\$267.00	<a href="http://www.chromtech.net.au/pdf2/vici09_92.pdf">http://www.chromtech.net.au/pdf2/vici09_92.pdf</a>
THC130	Hastelloy C 1/16inOD x.030inID	TU - Hastelloy C	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_91.pdf">http://www.chromtech.net.au/pdf2/vici09_91.pdf</a>
TNI120	Nickel 200 1/16inOD x.020inID	TU - Nickel	\$6.00	<a href="http://www.chromtech.net.au/pdf2/vici09_91.pdf">http://www.chromtech.net.au/pdf2/vici09_91.pdf</a>
TNI285	Nickel 200 1/8inOD x.085inID	TU - Nickel	\$9.00	<a href="http://www.chromtech.net.au/pdf2/vici09_91.pdf">http://www.chromtech.net.au/pdf2/vici09_91.pdf</a>
TNI130	Nickel 200 1/16inOD x.030inID	TU - Nickel	\$6.00	<a href="http://www.chromtech.net.au/pdf2/vici09_91.pdf">http://www.chromtech.net.au/pdf2/vici09_91.pdf</a>
TNI.520	Nickel 200 1/32inOD x.020inID	TU - Nickel	\$6.00	<a href="http://www.chromtech.net.au/pdf2/vici09_91.pdf">http://www.chromtech.net.au/pdf2/vici09_91.pdf</a>
TNI.510	Nickel 200 1/32inOD x.010inID	TU - Nickel	\$6.00	<a href="http://www.chromtech.net.au/pdf2/vici09_91.pdf">http://www.chromtech.net.au/pdf2/vici09_91.pdf</a>
TNI105	Nickel 200 1/16inOD x.005inID	TU - Nickel	\$8.00	<a href="http://www.chromtech.net.au/pdf2/vici09_91.pdf">http://www.chromtech.net.au/pdf2/vici09_91.pdf</a>
TNI140	Nickel 200 1/16inOD x.040inID	TU - Nickel	\$6.00	<a href="http://www.chromtech.net.au/pdf2/vici09_91.pdf">http://www.chromtech.net.au/pdf2/vici09_91.pdf</a>
TSS260	316 Stainless 1/8inOD x.060inID	TU - Stainless Steel	\$6.00	<a href="http://www.chromtech.net.au/pdf2/vici09_91.pdf">http://www.chromtech.net.au/pdf2/vici09_91.pdf</a>
TSS105	316 Stainless 1/16inOD x.005inID	TU - Stainless Steel	\$4.00	<a href="http://www.chromtech.net.au/pdf2/vici09_91.pdf">http://www.chromtech.net.au/pdf2/vici09_91.pdf</a>
TSS285	316 Stainless 1/8inOD x.085inID	TU - Stainless Steel	\$4.00	<a href="http://www.chromtech.net.au/pdf2/vici09_91.pdf">http://www.chromtech.net.au/pdf2/vici09_91.pdf</a>
TSS.505	316 Stainless 1/32inOD x.005inID	TU - Stainless Steel	\$4.00	<a href="http://www.chromtech.net.au/pdf2/vici09_91.pdf">http://www.chromtech.net.au/pdf2/vici09_91.pdf</a>
TSS.507	316 Stainless 1/32inOD x.007inID	TU - Stainless Steel	\$4.00	<a href="http://www.chromtech.net.au/pdf2/vici09_91.pdf">http://www.chromtech.net.au/pdf2/vici09_91.pdf</a>
TSS.510	316 Stainless 1/32inOD x.010inID	TU - Stainless Steel	\$4.00	<a href="http://www.chromtech.net.au/pdf2/vici09_91.pdf">http://www.chromtech.net.au/pdf2/vici09_91.pdf</a>
TSS110	316 Stainless 1/16inOD x.010inID	TU - Stainless Steel	\$4.00	<a href="http://www.chromtech.net.au/pdf2/vici09_91.pdf">http://www.chromtech.net.au/pdf2/vici09_91.pdf</a>
TSS115	316 Stainless 1/16inOD x.015inID	TU - Stainless Steel	\$4.00	<a href="http://www.chromtech.net.au/pdf2/vici09_91.pdf">http://www.chromtech.net.au/pdf2/vici09_91.pdf</a>
TSS.520	316 Stainless 1/32inOD x.020inID	TU - Stainless Steel	\$4.00	<a href="http://www.chromtech.net.au/pdf2/vici09_91.pdf">http://www.chromtech.net.au/pdf2/vici09_91.pdf</a>
TSS120	316 Stainless 1/16inOD x.020inID	TU - Stainless Steel	\$4.00	<a href="http://www.chromtech.net.au/pdf2/vici09_91.pdf">http://www.chromtech.net.au/pdf2/vici09_91.pdf</a>
TSS126	316 Stainless 1/16inOD x.026inID	TU - Stainless Steel	\$4.00	<a href="http://www.chromtech.net.au/pdf2/vici09_91.pdf">http://www.chromtech.net.au/pdf2/vici09_91.pdf</a>
TSS230	316 Stainless 1/8inOD x.030inID	TU - Stainless Steel	\$6.00	<a href="http://www.chromtech.net.au/pdf2/vici09_91.pdf">http://www.chromtech.net.au/pdf2/vici09_91.pdf</a>
TSS140	316 Stainless 1/16inOD x.040inID	TU - Stainless Steel	\$4.00	<a href="http://www.chromtech.net.au/pdf2/vici09_91.pdf">http://www.chromtech.net.au/pdf2/vici09_91.pdf</a>
TSS240	316 Stainless 1/8inOD x.040inID	TU - Stainless Steel	\$6.00	<a href="http://www.chromtech.net.au/pdf2/vici09_91.pdf">http://www.chromtech.net.au/pdf2/vici09_91.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

TSS130	316 Stainless 1/16inOD x.030inID	TU - Stainless Steel	\$4.00	<a href="http://www.chromtech.net.au/pdf2/vici09_91.pdf">http://www.chromtech.net.au/pdf2/vici09_91.pdf</a>
TSS267	316 Stainless 1/8inOD x.070inID	TU - Stainless Steel	\$6.00	<a href="http://www.chromtech.net.au/pdf2/vici09_91.pdf">http://www.chromtech.net.au/pdf2/vici09_91.pdf</a>
JR-797	Tool Clean-Cut Tubing Cutter	HPLC - Tools	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_90.pdf">http://www.chromtech.net.au/pdf2/vici09_90.pdf</a>
JR-798	Tool replacement blade for JR-797	HPLC - Tools	\$3.00	<a href="http://www.chromtech.net.au/pdf2/vici09_90.pdf">http://www.chromtech.net.au/pdf2/vici09_90.pdf</a>
TFEP260-100FT	1/8in FEP Tubing 0.060in 100ft	Polymeric Tubing	\$224.00	<a href="http://www.chromtech.net.au/pdf2/vici09_90.pdf">http://www.chromtech.net.au/pdf2/vici09_90.pdf</a>
TTF110-10FT	1/16in PTFE Tubing 0.010in 10ft	Polymeric Tubing	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_90.pdf">http://www.chromtech.net.au/pdf2/vici09_90.pdf</a>
JR-9001-5		Polymeric Tubing	\$15.00	<a href="http://www.chromtech.net.au/pdf2/vici09_90.pdf">http://www.chromtech.net.au/pdf2/vici09_90.pdf</a>
TFEP260-25FT	1/8in FEP Tubing 0.060in 25ft	Polymeric Tubing	\$57.00	<a href="http://www.chromtech.net.au/pdf2/vici09_90.pdf">http://www.chromtech.net.au/pdf2/vici09_90.pdf</a>
TTZ260-25FT	1/8in ETFE Tubing 0.060in 25ft	Polymeric Tubing	\$102.00	<a href="http://www.chromtech.net.au/pdf2/vici09_90.pdf">http://www.chromtech.net.au/pdf2/vici09_90.pdf</a>
TTF106-25FT	1/16in PTFE Tubing 0.006in 25ft	Polymeric Tubing	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_90.pdf">http://www.chromtech.net.au/pdf2/vici09_90.pdf</a>
TTF130-10FT	1/16in PTFE Tubing 0.030in 10ft	Polymeric Tubing	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_90.pdf">http://www.chromtech.net.au/pdf2/vici09_90.pdf</a>
TTF115-10FT	1/16in PTFE Tubing 0.015in 10ft	Polymeric Tubing	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_90.pdf">http://www.chromtech.net.au/pdf2/vici09_90.pdf</a>
TTF106-10FT	1/16in PTFE Tubing 0.006in 10ft	Polymeric Tubing	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_90.pdf">http://www.chromtech.net.au/pdf2/vici09_90.pdf</a>
TFEP260-10FT	1/8in FEP Tubing 0.060in 10ft	Polymeric Tubing	\$26.00	<a href="http://www.chromtech.net.au/pdf2/vici09_90.pdf">http://www.chromtech.net.au/pdf2/vici09_90.pdf</a>
TTZ260-10FT	1/8in ETFE Tubing 0.060in 10ft	Polymeric Tubing	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_90.pdf">http://www.chromtech.net.au/pdf2/vici09_90.pdf</a>
TTF110-25FT	1/16in PTFE Tubing 0.010in 25ft	Polymeric Tubing	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_90.pdf">http://www.chromtech.net.au/pdf2/vici09_90.pdf</a>
TTF120-10FT	1/16in PTFE Tubing 0.020in 10ft	Polymeric Tubing	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_90.pdf">http://www.chromtech.net.au/pdf2/vici09_90.pdf</a>
TTF130-100FT	1/16in PTFE Tubing 0.030in 100ft	Polymeric Tubing	\$165.00	<a href="http://www.chromtech.net.au/pdf2/vici09_90.pdf">http://www.chromtech.net.au/pdf2/vici09_90.pdf</a>
TFEP120-100FT	1/16in FEP Tubing 0.020in 100ft	Polymeric Tubing	\$200.00	<a href="http://www.chromtech.net.au/pdf2/vici09_90.pdf">http://www.chromtech.net.au/pdf2/vici09_90.pdf</a>
TTF285-100FT	1/8in PTFE Tubing 0.085in 100ft	Polymeric Tubing	\$216.00	<a href="http://www.chromtech.net.au/pdf2/vici09_90.pdf">http://www.chromtech.net.au/pdf2/vici09_90.pdf</a>
TFEP130-25FT	1/16in FEP Tubing 0.030in 25ft	Polymeric Tubing	\$50.00	<a href="http://www.chromtech.net.au/pdf2/vici09_90.pdf">http://www.chromtech.net.au/pdf2/vici09_90.pdf</a>
TTZ260-100FT	1/8in ETFE Tubing 0.060in 100ft	Polymeric Tubing	\$396.00	<a href="http://www.chromtech.net.au/pdf2/vici09_90.pdf">http://www.chromtech.net.au/pdf2/vici09_90.pdf</a>
TFEP110-25FT	1/16in FEP Tubing 0.010in 25ft	Polymeric Tubing	\$50.00	<a href="http://www.chromtech.net.au/pdf2/vici09_90.pdf">http://www.chromtech.net.au/pdf2/vici09_90.pdf</a>
TFEP130-10FT	1/16in FEP Tubing 0.030in 10ft	Polymeric Tubing	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_90.pdf">http://www.chromtech.net.au/pdf2/vici09_90.pdf</a>
TFEP130-100FT	1/16in FEP Tubing 0.030in 100ft	Polymeric Tubing	\$200.00	<a href="http://www.chromtech.net.au/pdf2/vici09_90.pdf">http://www.chromtech.net.au/pdf2/vici09_90.pdf</a>
TFEP110-10FT	1/16in FEP Tubing 0.010in 10ft	Polymeric Tubing	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_90.pdf">http://www.chromtech.net.au/pdf2/vici09_90.pdf</a>
TFEP110-100FT	1/16in FEP Tubing 0.010in 100ft	Polymeric Tubing	\$200.00	<a href="http://www.chromtech.net.au/pdf2/vici09_90.pdf">http://www.chromtech.net.au/pdf2/vici09_90.pdf</a>
TTF120-100FT	1/16in PTFE Tubing 0.020in 100ft	Polymeric Tubing	\$165.00	<a href="http://www.chromtech.net.au/pdf2/vici09_90.pdf">http://www.chromtech.net.au/pdf2/vici09_90.pdf</a>
TTF115-100FT	1/16in PTFE Tubing 0.015in 100ft	Polymeric Tubing	\$165.00	<a href="http://www.chromtech.net.au/pdf2/vici09_90.pdf">http://www.chromtech.net.au/pdf2/vici09_90.pdf</a>
TTF110-100FT	1/16in PTFE Tubing 0.010in 100ft	Polymeric Tubing	\$165.00	<a href="http://www.chromtech.net.au/pdf2/vici09_90.pdf">http://www.chromtech.net.au/pdf2/vici09_90.pdf</a>
TTF106-100FT	1/16in PTFE Tubing 0.006in 100ft	Polymeric Tubing	\$165.00	<a href="http://www.chromtech.net.au/pdf2/vici09_90.pdf">http://www.chromtech.net.au/pdf2/vici09_90.pdf</a>
TTF130-25FT	1/16in PTFE Tubing 0.030in 25ft	Polymeric Tubing	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_90.pdf">http://www.chromtech.net.au/pdf2/vici09_90.pdf</a>
TTF115-25FT	1/16in PTFE Tubing 0.015in 100ft	Polymeric Tubing	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_90.pdf">http://www.chromtech.net.au/pdf2/vici09_90.pdf</a>
TFEP120-10FT	1/16in FEP Tubing 0.020in 10ft	Polymeric Tubing	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_90.pdf">http://www.chromtech.net.au/pdf2/vici09_90.pdf</a>
TTZ130-100FT	1/16in ETFE Tubing 0.030in 100ft	Polymeric Tubing	\$267.00	<a href="http://www.chromtech.net.au/pdf2/vici09_90.pdf">http://www.chromtech.net.au/pdf2/vici09_90.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

TTF230-10FT	1/8in PTFE Tubing 0.030in 10ft	Polymeric Tubing	\$26.00	<a href="http://www.chromtech.net.au/pdf2/vici09_90.pdf">http://www.chromtech.net.au/pdf2/vici09_90.pdf</a>
TTF230-25FT	1/8in PTFE Tubing 0.030in 25ft	Polymeric Tubing	\$57.00	<a href="http://www.chromtech.net.au/pdf2/vici09_90.pdf">http://www.chromtech.net.au/pdf2/vici09_90.pdf</a>
TFEP120-25FT	1/16in FEP Tubing 0.020in 25ft	Polymeric Tubing	\$50.00	<a href="http://www.chromtech.net.au/pdf2/vici09_90.pdf">http://www.chromtech.net.au/pdf2/vici09_90.pdf</a>
TTF260-25FT	1/8in PTFE Tubing 0.060in 25ft	Polymeric Tubing	\$57.00	<a href="http://www.chromtech.net.au/pdf2/vici09_90.pdf">http://www.chromtech.net.au/pdf2/vici09_90.pdf</a>
TTZ110-10FT	1/16in ETFE Tubing 0.010in 10ft	Polymeric Tubing	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_90.pdf">http://www.chromtech.net.au/pdf2/vici09_90.pdf</a>
TTF285-25FT	1/8in PTFE Tubing 0.085in 25ft	Polymeric Tubing	\$57.00	<a href="http://www.chromtech.net.au/pdf2/vici09_90.pdf">http://www.chromtech.net.au/pdf2/vici09_90.pdf</a>
TTF260-10FT	1/8in PTFE Tubing 0.060in 10ft	Polymeric Tubing	\$26.00	<a href="http://www.chromtech.net.au/pdf2/vici09_90.pdf">http://www.chromtech.net.au/pdf2/vici09_90.pdf</a>
TTF260-100FT	1/8in PTFE Tubing 0.060in 100ft	Polymeric Tubing	\$216.00	<a href="http://www.chromtech.net.au/pdf2/vici09_90.pdf">http://www.chromtech.net.au/pdf2/vici09_90.pdf</a>
TTZ120-100FT	1/16in ETFE Tubing 0.020in 100ft	Polymeric Tubing	\$267.00	<a href="http://www.chromtech.net.au/pdf2/vici09_90.pdf">http://www.chromtech.net.au/pdf2/vici09_90.pdf</a>
TTZ110-100FT	1/16in ETFE Tubing 0.010in 100ft	Polymeric Tubing	\$267.00	<a href="http://www.chromtech.net.au/pdf2/vici09_90.pdf">http://www.chromtech.net.au/pdf2/vici09_90.pdf</a>
TTZ130-25FT	1/16in ETFE Tubing 0.030in 25ft	Polymeric Tubing	\$72.00	<a href="http://www.chromtech.net.au/pdf2/vici09_90.pdf">http://www.chromtech.net.au/pdf2/vici09_90.pdf</a>
TTZ120-25FT	1/16in ETFE Tubing 0.020in 25ft	Polymeric Tubing	\$72.00	<a href="http://www.chromtech.net.au/pdf2/vici09_90.pdf">http://www.chromtech.net.au/pdf2/vici09_90.pdf</a>
TTZ110-25FT	1/16in ETFE Tubing 0.010in 25ft	Polymeric Tubing	\$72.00	<a href="http://www.chromtech.net.au/pdf2/vici09_90.pdf">http://www.chromtech.net.au/pdf2/vici09_90.pdf</a>
TTZ130-10FT	1/16in ETFE Tubing 0.030in 10ft	Polymeric Tubing	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_90.pdf">http://www.chromtech.net.au/pdf2/vici09_90.pdf</a>
TTZ120-10FT	1/16in ETFE Tubing 0.020in 10ft	Polymeric Tubing	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_90.pdf">http://www.chromtech.net.au/pdf2/vici09_90.pdf</a>
TTF230-100FT	1/8in PTFE Tubing 0.030in 100ft	Polymeric Tubing	\$216.00	<a href="http://www.chromtech.net.au/pdf2/vici09_90.pdf">http://www.chromtech.net.au/pdf2/vici09_90.pdf</a>
JR-T-98075	1/16in Striped PEEK 0.030in Green Kit 5ea:5.10,20cm	PEEK Tubing Color-coded	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_89.pdf">http://www.chromtech.net.au/pdf2/vici09_89.pdf</a>
JR-T-98013	1/16in Striped PEEK 0.005in Red Kit 5ea:5.10,20cm	PEEK Tubing Color-coded	\$48.00	<a href="http://www.chromtech.net.au/pdf2/vici09_89.pdf">http://www.chromtech.net.au/pdf2/vici09_89.pdf</a>
JR-TD-6010	1/16in Dual Layer PEEK Tubing 0.020in OrangeOuter per ft	PEEK Tubing Color-coded	\$3.00	<a href="http://www.chromtech.net.au/pdf2/vici09_89.pdf">http://www.chromtech.net.au/pdf2/vici09_89.pdf</a>
JR-T-98200	1/16in Striped PEEK 0.005in Super Kit 75pieces allsizes	PEEK Tubing Color-coded	\$198.00	<a href="http://www.chromtech.net.au/pdf2/vici09_89.pdf">http://www.chromtech.net.au/pdf2/vici09_89.pdf</a>
JR-T-98050	1/16in Striped PEEK 0.020in Orange Kit 5ea:5.10,20cm	PEEK Tubing Color-coded	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_89.pdf">http://www.chromtech.net.au/pdf2/vici09_89.pdf</a>
JR-T-98017	1/16in Striped PEEK 0.007in Yellow Kit 5ea:5.10,20cm	PEEK Tubing Color-coded	\$48.00	<a href="http://www.chromtech.net.au/pdf2/vici09_89.pdf">http://www.chromtech.net.au/pdf2/vici09_89.pdf</a>
JR-T-60031	1/16in Striped PEEK 0.04in Grey; per metre	PEEK Tubing Color-coded	\$4.00	<a href="http://www.chromtech.net.au/pdf2/vici09_89.pdf">http://www.chromtech.net.au/pdf2/vici09_89.pdf</a>
JR-T-6003	Tubing, PEEK, 1/16" x 0.75 mm ID, green stripe; per metre	PEEK Tubing Color-coded	\$4.00	<a href="http://www.chromtech.net.au/pdf2/vici09_89.pdf">http://www.chromtech.net.au/pdf2/vici09_89.pdf</a>
JR-T-6001	Tubing PEEK 1/16" x 0.25 mm ID blue striped; per metre	PEEK Tubing Color-coded	\$4.00	<a href="http://www.chromtech.net.au/pdf2/vici09_89.pdf">http://www.chromtech.net.au/pdf2/vici09_89.pdf</a>
JR-T-6000	Tubing, PEEK, 1/16" x 018 mm ID, yellow striped; per metre	PEEK Tubing Color-coded	\$4.00	<a href="http://www.chromtech.net.au/pdf2/vici09_89.pdf">http://www.chromtech.net.au/pdf2/vici09_89.pdf</a>



VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

JR-TD-6007	1/16in Dual Layer PEEK Tubing 0.005in RedOuter per ft	PEEK Tubing Color-coded	\$4.00	<a href="http://www.chromtech.net.au/pdf2/vici09_89.pdf">http://www.chromtech.net.au/pdf2/vici09_89.pdf</a>
JR-T-98025	1/16in Striped PEEK 0.010in Blue Kit Sea:5.10,20cm	PEEK Tubing Color-coded	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_89.pdf">http://www.chromtech.net.au/pdf2/vici09_89.pdf</a>
JR-TD-6009	1/16in Dual Layer PEEK Tubing 0.010in BlueOuter per ft	PEEK Tubing Color-coded	\$3.00	<a href="http://www.chromtech.net.au/pdf2/vici09_89.pdf">http://www.chromtech.net.au/pdf2/vici09_89.pdf</a>
JR-TD-6008	1/16in Dual Layer PEEK Tubing 0.007in YellowOuter per ft	PEEK Tubing Color-coded	\$4.00	<a href="http://www.chromtech.net.au/pdf2/vici09_89.pdf">http://www.chromtech.net.au/pdf2/vici09_89.pdf</a>
JR-TD-6011	1/16in Dual Layer PEEK Tubing 0.030inGreenOuter per ft	PEEK Tubing Color-coded	\$3.00	<a href="http://www.chromtech.net.au/pdf2/vici09_89.pdf">http://www.chromtech.net.au/pdf2/vici09_89.pdf</a>
JR-TD-5804	1/16in Dual Layer PEEK Tubing 0.004in BlackOuter per ft	PEEK Tubing Color-coded	\$4.00	<a href="http://www.chromtech.net.au/pdf2/vici09_89.pdf">http://www.chromtech.net.au/pdf2/vici09_89.pdf</a>
JR-T-5804	1/16 PEEK Striped Black 0.004"	PEEK Tubing Natural	\$4.00	<a href="http://www.chromtech.net.au/pdf2/vici09_89.pdf">http://www.chromtech.net.au/pdf2/vici09_89.pdf</a>
JR-T-6002	Tubing PEEK 1/16" x 0.50 mm ID orange striped; per metre	PEEK Tubing Natural	\$4.00	<a href="http://www.chromtech.net.au/pdf2/vici09_89.pdf">http://www.chromtech.net.au/pdf2/vici09_89.pdf</a>
JR-T-5999	Tubing, PEEK, 1/16" x 0.005"Red striped; per metre	PEEK Tubing Natural	\$4.00	<a href="http://www.chromtech.net.au/pdf2/vici09_89.pdf">http://www.chromtech.net.au/pdf2/vici09_89.pdf</a>
TPK.102	360um PEEK Tubing .002in per ft	EF Nickel Tubing	\$6.00	<a href="http://www.chromtech.net.au/pdf2/vici09_88.pdf">http://www.chromtech.net.au/pdf2/vici09_88.pdf</a>
TPK130-25FT	1/16in PEEK Tubing .030in per 25ft	PEEK Tubing Natural	\$76.00	<a href="http://www.chromtech.net.au/pdf2/vici09_88.pdf">http://www.chromtech.net.au/pdf2/vici09_88.pdf</a>
TPK260-100FT	1/8in PEEK Tubing .060in per 100ft	PEEK Tubing Natural	\$471.00	<a href="http://www.chromtech.net.au/pdf2/vici09_88.pdf">http://www.chromtech.net.au/pdf2/vici09_88.pdf</a>
TPK.515-100FT	1/32in PEEK Tubing .015in per 100ft	PEEK Tubing Natural	\$306.00	<a href="http://www.chromtech.net.au/pdf2/vici09_88.pdf">http://www.chromtech.net.au/pdf2/vici09_88.pdf</a>
TPK.104	360um PEEK Tubing .004in per ft	PEEK Tubing Natural	\$6.00	<a href="http://www.chromtech.net.au/pdf2/vici09_88.pdf">http://www.chromtech.net.au/pdf2/vici09_88.pdf</a>
TPK.105	360um PEEK Tubing .005in per ft	PEEK Tubing Natural	\$6.00	<a href="http://www.chromtech.net.au/pdf2/vici09_88.pdf">http://www.chromtech.net.au/pdf2/vici09_88.pdf</a>
TPK.106	360um PEEK Tubing .006in per ft	PEEK Tubing Natural	\$6.00	<a href="http://www.chromtech.net.au/pdf2/vici09_88.pdf">http://www.chromtech.net.au/pdf2/vici09_88.pdf</a>
TPK.502-10FT	1/32in PEEK Tubing .002in per 10ft	PEEK Tubing Natural	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_88.pdf">http://www.chromtech.net.au/pdf2/vici09_88.pdf</a>
TPK.505-10FT	1/32in PEEK Tubing .005in per 10ft	PEEK Tubing Natural	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_88.pdf">http://www.chromtech.net.au/pdf2/vici09_88.pdf</a>
TPK.510-10FT	1/32in PEEK Tubing .010in per 10ft	PEEK Tubing Natural	\$33.00	<a href="http://www.chromtech.net.au/pdf2/vici09_88.pdf">http://www.chromtech.net.au/pdf2/vici09_88.pdf</a>
TPK.515-10FT	1/32in PEEK Tubing .015in per 10ft	PEEK Tubing Natural	\$33.00	<a href="http://www.chromtech.net.au/pdf2/vici09_88.pdf">http://www.chromtech.net.au/pdf2/vici09_88.pdf</a>
TPK.502-25FT	1/32in PEEK Tubing .002in per 25 ft	PEEK Tubing Natural	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_88.pdf">http://www.chromtech.net.au/pdf2/vici09_88.pdf</a>
TPK.505-25FT	1/32in PEEK Tubing .005in per 25ft	PEEK Tubing Natural	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_88.pdf">http://www.chromtech.net.au/pdf2/vici09_88.pdf</a>
TPK110-100FT	1/16in PEEK Tubing .010in per 100ft	PEEK Tubing Natural	\$294.00	<a href="http://www.chromtech.net.au/pdf2/vici09_88.pdf">http://www.chromtech.net.au/pdf2/vici09_88.pdf</a>
TPK.515-25FT	1/32in PEEK Tubing .015in per 25ft	PEEK Tubing Natural	\$80.00	<a href="http://www.chromtech.net.au/pdf2/vici09_88.pdf">http://www.chromtech.net.au/pdf2/vici09_88.pdf</a>
TPK260-25FT	1/8in PEEK Tubing .060in per 25ft	PEEK Tubing Natural	\$122.00	<a href="http://www.chromtech.net.au/pdf2/vici09_88.pdf">http://www.chromtech.net.au/pdf2/vici09_88.pdf</a>
TPK106-10FT	1/16in PEEK Tubing .006in per 10ft	PEEK Tubing Natural	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_88.pdf">http://www.chromtech.net.au/pdf2/vici09_88.pdf</a>
TPK110-10FT	1/16in PEEK Tubing .010in per 10ft	PEEK Tubing Natural	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_88.pdf">http://www.chromtech.net.au/pdf2/vici09_88.pdf</a>
TPK120-10FT	1/16in PEEK Tubing .020in per 10ft	PEEK Tubing Natural	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_88.pdf">http://www.chromtech.net.au/pdf2/vici09_88.pdf</a>

VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

TPK130-10FT	1/16in PEEK Tubing .030in per 10ft	PEEK Tubing Natural	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_88.pdf">http://www.chromtech.net.au/pdf2/vici09_88.pdf</a>
TPK106-25FT	1/16in PEEK Tubing .006in per 25ft	PEEK Tubing Natural	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_88.pdf">http://www.chromtech.net.au/pdf2/vici09_88.pdf</a>
TPK110-25FT	1/16in PEEK Tubing .010in per 25ft	PEEK Tubing Natural	\$76.00	<a href="http://www.chromtech.net.au/pdf2/vici09_88.pdf">http://www.chromtech.net.au/pdf2/vici09_88.pdf</a>
TPK120-25FT	1/16in PEEK Tubing .020in per 25ft	PEEK Tubing Natural	\$76.00	<a href="http://www.chromtech.net.au/pdf2/vici09_88.pdf">http://www.chromtech.net.au/pdf2/vici09_88.pdf</a>
TPK106-100FT	1/16in PEEK Tubing .006in per 100ft	PEEK Tubing Natural	\$357.00	<a href="http://www.chromtech.net.au/pdf2/vici09_88.pdf">http://www.chromtech.net.au/pdf2/vici09_88.pdf</a>
TPK120-100FT	1/16in PEEK Tubing .020in per 100ft	PEEK Tubing Natural	\$294.00	<a href="http://www.chromtech.net.au/pdf2/vici09_88.pdf">http://www.chromtech.net.au/pdf2/vici09_88.pdf</a>
TPK130-100FT	1/16in PEEK Tubing .030in per 100ft	PEEK Tubing Natural	\$294.00	<a href="http://www.chromtech.net.au/pdf2/vici09_88.pdf">http://www.chromtech.net.au/pdf2/vici09_88.pdf</a>
TPK260-10FT	1/8in PEEK Tubing .060in per 10ft	PEEK Tubing Natural	\$51.00	<a href="http://www.chromtech.net.au/pdf2/vici09_88.pdf">http://www.chromtech.net.au/pdf2/vici09_88.pdf</a>
TPK.510-25FT	1/32in PEEK Tubing .010in per 25t	PEEK Tubing Natural	\$80.00	<a href="http://www.chromtech.net.au/pdf2/vici09_88.pdf">http://www.chromtech.net.au/pdf2/vici09_88.pdf</a>
TEFNI130	1/16in EFNi Tubing .030in per ft,50ft max	EF Nickel Tubing	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_87.pdf">http://www.chromtech.net.au/pdf2/vici09_87.pdf</a>
TEFNI.502	1/32in EFNi Tubing .002in per ft,2ft max	EF Nickel Tubing	\$26.00	<a href="http://www.chromtech.net.au/pdf2/vici09_87.pdf">http://www.chromtech.net.au/pdf2/vici09_87.pdf</a>
TEFNI.110	360um EFNi Tubing .007in per ft,50ft max	EF Nickel Tubing	\$15.00	<a href="http://www.chromtech.net.au/pdf2/vici09_87.pdf">http://www.chromtech.net.au/pdf2/vici09_87.pdf</a>
TEFNI.105	360um EFNi Tubing .005in per ft,30ft max	EF Nickel Tubing	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_87.pdf">http://www.chromtech.net.au/pdf2/vici09_87.pdf</a>
TEFNI.505	1/32in EFNi Tubing .005in per ft,30ft max	EF Nickel Tubing	\$22.00	<a href="http://www.chromtech.net.au/pdf2/vici09_87.pdf">http://www.chromtech.net.au/pdf2/vici09_87.pdf</a>
TEFNI.104	360um EFNi Tubing .004in per ft,20ft max	EF Nickel Tubing	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_87.pdf">http://www.chromtech.net.au/pdf2/vici09_87.pdf</a>
TEFNI.504	1/32in EFNi Tubing .004in per ft,20ft max	EF Nickel Tubing	\$22.00	<a href="http://www.chromtech.net.au/pdf2/vici09_87.pdf">http://www.chromtech.net.au/pdf2/vici09_87.pdf</a>
TEFNI.102	360um EFNi Tubing .002in per ft,2ft max	EF Nickel Tubing	\$22.00	<a href="http://www.chromtech.net.au/pdf2/vici09_87.pdf">http://www.chromtech.net.au/pdf2/vici09_87.pdf</a>
TEFNI.101	360um EFNi Tubing .001in per ft, 1ft max	EF Nickel Tubing	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_87.pdf">http://www.chromtech.net.au/pdf2/vici09_87.pdf</a>
TEFNI120	1/16in EFNi Tubing .020in per ft,50ft max	EF Nickel Tubing	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_87.pdf">http://www.chromtech.net.au/pdf2/vici09_87.pdf</a>
TEFNI.520	1/32in EFNi Tubing .020in per ft,2ft max	EF Nickel Tubing	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_87.pdf">http://www.chromtech.net.au/pdf2/vici09_87.pdf</a>
TEFNI.515	1/32in EFNi Tubing .015in per ft,50ft max	EF Nickel Tubing	\$16.00	<a href="http://www.chromtech.net.au/pdf2/vici09_87.pdf">http://www.chromtech.net.au/pdf2/vici09_87.pdf</a>
TEFNI140	1/16in EFNi Tubing .040in per ft,50ft max	EF Nickel Tubing	\$15.00	<a href="http://www.chromtech.net.au/pdf2/vici09_87.pdf">http://www.chromtech.net.au/pdf2/vici09_87.pdf</a>
TEFNI.512	1/32in EFNi Tubing .012in per ft,50ft max	EF Nickel Tubing	\$16.00	<a href="http://www.chromtech.net.au/pdf2/vici09_87.pdf">http://www.chromtech.net.au/pdf2/vici09_87.pdf</a>
TEFNI.510	1/32in EFNi Tubing .010in per ft,50ft max	EF Nickel Tubing	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_87.pdf">http://www.chromtech.net.au/pdf2/vici09_87.pdf</a>
TEFNI.507	1/32in EFNi Tubing .007in per ft,50ft max	EF Nickel Tubing	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_87.pdf">http://www.chromtech.net.au/pdf2/vici09_87.pdf</a>
ASFV40K360	UHP Valve 360um 0.15mm	UHP Valves	\$1,122.00	<a href="http://www.chromtech.net.au/pdf2/vici09_85.pdf">http://www.chromtech.net.au/pdf2/vici09_85.pdf</a>
ASFVOD40K1	Dual On/off Valve 1/16in 0.15mm	UHP Valves	\$1,606.00	<a href="http://www.chromtech.net.au/pdf2/vici09_85.pdf">http://www.chromtech.net.au/pdf2/vici09_85.pdf</a>
ASFVOD40K.5	Dual On/off Valve 1/32in 0.15mm	UHP Valves	\$1,708.00	<a href="http://www.chromtech.net.au/pdf2/vici09_85.pdf">http://www.chromtech.net.au/pdf2/vici09_85.pdf</a>
ASFVOD40K360	Dual On/off Valve 360um 0.15mm	UHP Valves	\$1,989.00	<a href="http://www.chromtech.net.au/pdf2/vici09_85.pdf">http://www.chromtech.net.au/pdf2/vici09_85.pdf</a>
ASFVO40K360	Prime/Purge Valve 360um 0.15mm	UHP Valves	\$1,071.00	<a href="http://www.chromtech.net.au/pdf2/vici09_85.pdf">http://www.chromtech.net.au/pdf2/vici09_85.pdf</a>
ASFV40K.5	UHP Valve 132in 0.15mm	UHP Valves	\$1,008.00	<a href="http://www.chromtech.net.au/pdf2/vici09_85.pdf">http://www.chromtech.net.au/pdf2/vici09_85.pdf</a>
ASFVO40K1	Prime/Purge Valve 1/16in 0.15mm	UHP Valves	\$912.00	<a href="http://www.chromtech.net.au/pdf2/vici09_85.pdf">http://www.chromtech.net.au/pdf2/vici09_85.pdf</a>
ASFVO40K.5	Prime/Purge Valve 1/32in0.15mm	UHP Valves	\$969.00	<a href="http://www.chromtech.net.au/pdf2/vici09_85.pdf">http://www.chromtech.net.au/pdf2/vici09_85.pdf</a>
ASFV40K1	UHP Valve 1/16in 0.15mm	UHP Valves	\$950.00	<a href="http://www.chromtech.net.au/pdf2/vici09_85.pdf">http://www.chromtech.net.au/pdf2/vici09_85.pdf</a>
SPSS40	Ultra High Pressure Injector System	UHP Injector System	\$5,100.00	<a href="http://www.chromtech.net.au/pdf2/vici09_84.pdf">http://www.chromtech.net.au/pdf2/vici09_84.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

CP-CM1-P	Controller Microlynx-4	Liquid Handling Pumps	\$765.00	<a href="http://www.chromtech.net.au/pdf2/vici09_83.pdf">http://www.chromtech.net.au/pdf2/vici09_83.pdf</a>
CP-DSM2	M50 Pump Motor	Liquid Handling Pumps	\$702.00	<a href="http://www.chromtech.net.au/pdf2/vici09_83.pdf">http://www.chromtech.net.au/pdf2/vici09_83.pdf</a>
CP-DSM	M6 Pump Motor	Liquid Handling Pumps	\$669.00	<a href="http://www.chromtech.net.au/pdf2/vici09_83.pdf">http://www.chromtech.net.au/pdf2/vici09_83.pdf</a>
CP2-4841-100D	M6 Pump Only 50ul-50ml	Liquid Handling Pumps	\$1,964.00	<a href="http://www.chromtech.net.au/pdf2/vici09_83.pdf">http://www.chromtech.net.au/pdf2/vici09_83.pdf</a>
CP2-4841-100SM	M6 Pump w/Stepper Motor Only 10nl-10ml	Liquid Handling Pumps	\$2,626.00	<a href="http://www.chromtech.net.au/pdf2/vici09_83.pdf">http://www.chromtech.net.au/pdf2/vici09_83.pdf</a>
CP2-4841-100M1	M6 Pump w/ Controller Stepper Motor 10nl-10ml	Liquid Handling Pumps	\$3,494.00	<a href="http://www.chromtech.net.au/pdf2/vici09_83.pdf">http://www.chromtech.net.au/pdf2/vici09_83.pdf</a>
CD50-8182-M2A	M50 Dilutor/Dispenser	Dilutor/Dispensers,M Series	\$4,590.00	<a href="http://www.chromtech.net.au/pdf2/vici09_82.pdf">http://www.chromtech.net.au/pdf2/vici09_82.pdf</a>
CD10-4841-M1A	M10 Dilutor/Dispenser	Dilutor/Dispensers,M Series	\$4,462.00	<a href="http://www.chromtech.net.au/pdf2/vici09_82.pdf">http://www.chromtech.net.au/pdf2/vici09_82.pdf</a>
C-MPR2	Mobile phase reservoir -0.5 liter 2-hole 1/4-28	FIA Filters,Access	\$58.00	<a href="http://www.chromtech.net.au/pdf2/vici09_81.pdf">http://www.chromtech.net.au/pdf2/vici09_81.pdf</a>
C-BOT16	Mobile phase reservoir -0.5 liter plain	FIA Filters,Access	\$9.00	<a href="http://www.chromtech.net.au/pdf2/vici09_81.pdf">http://www.chromtech.net.au/pdf2/vici09_81.pdf</a>
JR-9000-0006	PTFE VICI-Cap S40 2 1/8-Port and vent	FIA Filters,Access	\$34.00	<a href="http://www.chromtech.net.au/pdf2/vici09_81.pdf">http://www.chromtech.net.au/pdf2/vici09_81.pdf</a>
C-BOT32	Mobile phase reservoir -1.0 liter plain	FIA Filters,Access	\$14.00	<a href="http://www.chromtech.net.au/pdf2/vici09_81.pdf">http://www.chromtech.net.au/pdf2/vici09_81.pdf</a>
JR-9000-0001	PTFE VICI-Cap GL-45 2 1/8-Port and vent	FIA Filters,Access	\$34.00	<a href="http://www.chromtech.net.au/pdf2/vici09_81.pdf">http://www.chromtech.net.au/pdf2/vici09_81.pdf</a>
C-MPR3	Mobile phase reservoir -0.5 liter 3-hole 1/4-28	FIA Filters,Access	\$76.00	<a href="http://www.chromtech.net.au/pdf2/vici09_81.pdf">http://www.chromtech.net.au/pdf2/vici09_81.pdf</a>
JR-9000-0602	Filter/Sparger PTFE Last Drop 2.5 um filter 10 um sparger	Filters,Perifits	\$62.00	<a href="http://www.chromtech.net.au/pdf2/vici09_79.pdf">http://www.chromtech.net.au/pdf2/vici09_79.pdf</a>
JR-9000-0530	Filter SS Last Drop 2 um stepped tubing connector	Filters,Perifits	\$27.00	<a href="http://www.chromtech.net.au/pdf2/vici09_79.pdf">http://www.chromtech.net.au/pdf2/vici09_79.pdf</a>
JR-9000-0640	Filter/Sparger SS Last Drop 2 um filter 10 um sparger	Filters,Perifits	\$62.00	<a href="http://www.chromtech.net.au/pdf2/vici09_79.pdf">http://www.chromtech.net.au/pdf2/vici09_79.pdf</a>
JR-9000-0520	Filter PTFE Last Drop 2.5 um stepped tubing connector	Filters,Perifits	\$27.00	<a href="http://www.chromtech.net.au/pdf2/vici09_79.pdf">http://www.chromtech.net.au/pdf2/vici09_79.pdf</a>
PS-614158	Valve for screw-cap vial with 13-425 mm bottle neck	MobilePhase Filters,Spargers	\$93.00	<a href="http://www.chromtech.net.au/pdf2/vici09_79.pdf">http://www.chromtech.net.au/pdf2/vici09_79.pdf</a>
PS-614163	Valve for screw-cap vial with 24-400 mm bottle neck	MobilePhase Filters,Spargers	\$112.00	<a href="http://www.chromtech.net.au/pdf2/vici09_79.pdf">http://www.chromtech.net.au/pdf2/vici09_79.pdf</a>
PS-614170	Valve for screw-cap vial with 20-400 mm bottle neck 12/	MobilePhase Filters,Spargers	\$106.00	<a href="http://www.chromtech.net.au/pdf2/vici09_79.pdf">http://www.chromtech.net.au/pdf2/vici09_79.pdf</a>
PS-614161	Valve for screw-cap vial with 18-400 mm bottle neck 12/	MobilePhase Filters,Spargers	\$100.00	<a href="http://www.chromtech.net.au/pdf2/vici09_79.pdf">http://www.chromtech.net.au/pdf2/vici09_79.pdf</a>
PS-614160	Valve for screw-cap vial with 15-425 mm bottle neck	MobilePhase Filters,Spargers	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_79.pdf">http://www.chromtech.net.au/pdf2/vici09_79.pdf</a>
PS-614250	Mininert Valve for 20mm Crimp Top Vial	MobilePhase Filters,Spargers	\$126.00	<a href="http://www.chromtech.net.au/pdf2/vici09_79.pdf">http://www.chromtech.net.au/pdf2/vici09_79.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

JR-367016-20	Filter SS mobile phase 20 um fittings for 1/16" tubing	Filters,Perifits	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_78.pdf">http://www.chromtech.net.au/pdf2/vici09_78.pdf</a>
JR-367008-2	Filter SS mobile phase 2 um fittings for 1/8" tubing	Filters,Perifits	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_78.pdf">http://www.chromtech.net.au/pdf2/vici09_78.pdf</a>
JR-367008-10	Filter SS mobile phase 10 um fittings for 1/8" tubing	Filters,Perifits	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_78.pdf">http://www.chromtech.net.au/pdf2/vici09_78.pdf</a>
C-MPFTI2		Filters,Perifits	\$26.00	<a href="http://www.chromtech.net.au/pdf2/vici09_78.pdf">http://www.chromtech.net.au/pdf2/vici09_78.pdf</a>
JR-367008-20	Filter SS mobile phase 20 um fittings for 1/8" tubing	Filters,Perifits	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_78.pdf">http://www.chromtech.net.au/pdf2/vici09_78.pdf</a>
JR-367016-10	Filter SS mobile phase 10 um fittings for 1/16" tubing	Filters,Perifits	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_78.pdf">http://www.chromtech.net.au/pdf2/vici09_78.pdf</a>
JR-367016-2	Filter SS mobile phase 2 um fittings for 1/16" tubing	Filters,Perifits	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_78.pdf">http://www.chromtech.net.au/pdf2/vici09_78.pdf</a>
C-MPFTI10		Filters,Perifits	\$26.00	<a href="http://www.chromtech.net.au/pdf2/vici09_78.pdf">http://www.chromtech.net.au/pdf2/vici09_78.pdf</a>
JR-32178	Filter PE mobile phase No-Met 20 um 1/8"	Filters,Perifits	\$16.00	<a href="http://www.chromtech.net.au/pdf2/vici09_78.pdf">http://www.chromtech.net.au/pdf2/vici09_78.pdf</a>
C-PFL	Perifit fittings -for tube 1.85 to 2.29mmID	Filters,Perifits	\$10.00	<a href="http://www.chromtech.net.au/pdf2/vici09_78.pdf">http://www.chromtech.net.au/pdf2/vici09_78.pdf</a>
C-PFM	Perifit fittings -for tube 1.12 to 1.65mmID	Filters,Perifits	\$9.00	<a href="http://www.chromtech.net.au/pdf2/vici09_78.pdf">http://www.chromtech.net.au/pdf2/vici09_78.pdf</a>
C-PFS	Perifit fittings -for tube 0.50 to 1.02mmID	Filters,Perifits	\$9.00	<a href="http://www.chromtech.net.au/pdf2/vici09_78.pdf">http://www.chromtech.net.au/pdf2/vici09_78.pdf</a>
JR-32179	Filter PE mobile phase replacement No-Met 20 um	Filters,Perifits	\$3.00	<a href="http://www.chromtech.net.au/pdf2/vici09_78.pdf">http://www.chromtech.net.au/pdf2/vici09_78.pdf</a>
CFE-S2	In-line filter -2u 1/4-28	Adaptors	\$8.00	<a href="http://www.chromtech.net.au/pdf2/vici09_77.pdf">http://www.chromtech.net.au/pdf2/vici09_77.pdf</a>
CFPA4KF	CTFE Pipe adapters -Female 1/4-28 to female NPT 1/4inx1.5mm bore	Adaptors	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_77.pdf">http://www.chromtech.net.au/pdf2/vici09_77.pdf</a>
C-F1.5TI		Adaptors	\$6.00	<a href="http://www.chromtech.net.au/pdf2/vici09_77.pdf">http://www.chromtech.net.au/pdf2/vici09_77.pdf</a>
ZU1FPK.5	Biocompatible Filte 1/16 0.5mm PEK	Adaptors	\$56.00	<a href="http://www.chromtech.net.au/pdf2/vici09_77.pdf">http://www.chromtech.net.au/pdf2/vici09_77.pdf</a>
JR-26504	Coupler PEEK universal 0.50 mm ID orange	Adaptors	\$26.00	<a href="http://www.chromtech.net.au/pdf2/vici09_77.pdf">http://www.chromtech.net.au/pdf2/vici09_77.pdf</a>
JR-26503	Coupler PEEK universal 0.25 mm ID blue	Adaptors	\$26.00	<a href="http://www.chromtech.net.au/pdf2/vici09_77.pdf">http://www.chromtech.net.au/pdf2/vici09_77.pdf</a>
JR-26502	Coupler PEEK universal 0.17 mm ID yellow	Adaptors	\$26.00	<a href="http://www.chromtech.net.au/pdf2/vici09_77.pdf">http://www.chromtech.net.au/pdf2/vici09_77.pdf</a>
JR-26501	Coupler PEEK universal 0.13 mm ID red	Adaptors	\$26.00	<a href="http://www.chromtech.net.au/pdf2/vici09_77.pdf">http://www.chromtech.net.au/pdf2/vici09_77.pdf</a>
ZLCA1PK	1/4-28 to 10-32 ZVD Port Adaptor	Adaptors	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_77.pdf">http://www.chromtech.net.au/pdf2/vici09_77.pdf</a>
CFE-S10	In-line filter -10u 1/4-28	Adaptors	\$8.00	<a href="http://www.chromtech.net.au/pdf2/vici09_77.pdf">http://www.chromtech.net.au/pdf2/vici09_77.pdf</a>
CEPA2PK	PEEK Pipe adapters -Male 1/4-28 to male NPT 1/8inODx1.5mm bore	Adaptors	\$22.00	<a href="http://www.chromtech.net.au/pdf2/vici09_77.pdf">http://www.chromtech.net.au/pdf2/vici09_77.pdf</a>
CFPA2KF	CTFE Pipe adapters -Female 1/4-28 to female NPT 1/8inODx1.5mm bore	Adaptors	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_77.pdf">http://www.chromtech.net.au/pdf2/vici09_77.pdf</a>
CFPA2PK	PEEK Pipe adapters -Female 1/4-28 to female NPT 1/8inODx1.5mm bore	Adaptors	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_77.pdf">http://www.chromtech.net.au/pdf2/vici09_77.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

CEPA4KF	CTFE Pipe adapters -Male 1/4-28 to male NPT 1/4inODx1.5mm bore	Adaptors	\$22.00	<a href="http://www.chromtech.net.au/pdf2/vici09_77.pdf">http://www.chromtech.net.au/pdf2/vici09_77.pdf</a>
CEPA4PK	PEEK Pipe adapters -Male 1/4-28 to male NPT 1/4inODx1.5mm bore	Adaptors	\$22.00	<a href="http://www.chromtech.net.au/pdf2/vici09_77.pdf">http://www.chromtech.net.au/pdf2/vici09_77.pdf</a>
CEPA2KF	CTFE Pipe adapters -Male 1/4-28 to male NPT 1/8inODx1.5mm bore	Adaptors	\$22.00	<a href="http://www.chromtech.net.au/pdf2/vici09_77.pdf">http://www.chromtech.net.au/pdf2/vici09_77.pdf</a>
CPA2PK	PEEK Pipe adapters -Female 1/4-28 to male NPT 1/8inODx1.5mm bore	Adaptors	\$22.00	<a href="http://www.chromtech.net.au/pdf2/vici09_77.pdf">http://www.chromtech.net.au/pdf2/vici09_77.pdf</a>
CPA4PK	PEEK Pipe adapters -Female 1/4-28 to male NPT 1/4inODx1.5mm bore	Adaptors	\$22.00	<a href="http://www.chromtech.net.au/pdf2/vici09_77.pdf">http://www.chromtech.net.au/pdf2/vici09_77.pdf</a>
CPA2KF	CTFE Pipe adapters -Female 1/4-28 to male NPT 1/8inODx1.5mm bore	Adaptors	\$22.00	<a href="http://www.chromtech.net.au/pdf2/vici09_77.pdf">http://www.chromtech.net.au/pdf2/vici09_77.pdf</a>
CPA4KF	CTFE Pipe adapters -Female 1/4-28 to male NPT 1/4inODx1.5mm bore	Adaptors	\$22.00	<a href="http://www.chromtech.net.au/pdf2/vici09_77.pdf">http://www.chromtech.net.au/pdf2/vici09_77.pdf</a>
CFE-S75	In-line filter -5u 1/4-28	Adaptors	\$8.00	<a href="http://www.chromtech.net.au/pdf2/vici09_77.pdf">http://www.chromtech.net.au/pdf2/vici09_77.pdf</a>
CFPA4PK	PEEK Pipe adapters -Female 1/4-28 to female NPT 1/4inODx1.5mm bore	Adaptors	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_77.pdf">http://www.chromtech.net.au/pdf2/vici09_77.pdf</a>
CMLAKF	CTFE Luer adapter-Male Luer to 1/4-28 1.50mm bore	Adaptors	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_76.pdf">http://www.chromtech.net.au/pdf2/vici09_76.pdf</a>
CMLAPFA	PFA Luer adapter-Male Luer to 1/4-28 1.50mm bore	Adaptors	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_76.pdf">http://www.chromtech.net.au/pdf2/vici09_76.pdf</a>
CBUFLPK	PEEK Luer adapter-bulkhead unions -PEEK Female Luer to 1/4-28 x 0.75mm bore	Adaptors	\$33.00	<a href="http://www.chromtech.net.au/pdf2/vici09_76.pdf">http://www.chromtech.net.au/pdf2/vici09_76.pdf</a>
CBUFLKF	Female Luer adapter bulkhead union CTFE 1/16" OD tube	Adaptors	\$33.00	<a href="http://www.chromtech.net.au/pdf2/vici09_76.pdf">http://www.chromtech.net.au/pdf2/vici09_76.pdf</a>
ZBUFLKF	CTFE Luer adapter-bulkhead unions -CTFE Female Luer to 10-32 x 0.75mm bore	Adaptors	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_76.pdf">http://www.chromtech.net.au/pdf2/vici09_76.pdf</a>
ZBUFLPK	PEEK Luer adapter-bulkhead unions -PEEK Female Luer to 10-32 x 0.75mm bore	Adaptors	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_76.pdf">http://www.chromtech.net.au/pdf2/vici09_76.pdf</a>
CMLAPK	PEEK Luer Adapter-Male Luer to 1/4-28 1.50mm bore	Adaptors	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_76.pdf">http://www.chromtech.net.au/pdf2/vici09_76.pdf</a>
ZBUMLKF	CTFE Luer adapter-bulkhead unions -CTFE Male Luer to 10-32 x 1.00mm bore	Adaptors	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_76.pdf">http://www.chromtech.net.au/pdf2/vici09_76.pdf</a>
CTA4PK	PEEK Tube adapters -Male 1/4-28 1/4inx1.5mm bore	Adaptors	\$10.00	<a href="http://www.chromtech.net.au/pdf2/vici09_76.pdf">http://www.chromtech.net.au/pdf2/vici09_76.pdf</a>
ZUFLKF	CTFE Luer adapter-Female Luer to 10/32 X 0.75mm bore	Adaptors	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_76.pdf">http://www.chromtech.net.au/pdf2/vici09_76.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

ZBUMLPK	PEEK Luer adapter-bulkhead unions -PEEK Male Luer to 10-32 x 1.00mm bore	Adaptors	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_76.pdf">http://www.chromtech.net.au/pdf2/vici09_76.pdf</a>
CTA2KF	CTFE Tube adapters -Male 1/4-28 1/8inx1.5mm bore	Adaptors	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_76.pdf">http://www.chromtech.net.au/pdf2/vici09_76.pdf</a>
CTA4KF	Tube adapters -Male 1/4-28 1/4inx1.5mm boreCTFE	Adaptors	\$10.00	<a href="http://www.chromtech.net.au/pdf2/vici09_76.pdf">http://www.chromtech.net.au/pdf2/vici09_76.pdf</a>
CFLAPFA	PFA Luer adapter-Female Luer to 1/4-28 1.50mm bore	Adaptors	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_76.pdf">http://www.chromtech.net.au/pdf2/vici09_76.pdf</a>
CTA2S6	316 Stainless Tube adapters -Male 1/4-28 1/8inx1.5mm bore	Adaptors	\$10.00	<a href="http://www.chromtech.net.au/pdf2/vici09_76.pdf">http://www.chromtech.net.au/pdf2/vici09_76.pdf</a>
ZUFLPK	PEEK Female Luer to 10-32 Fitting	Adaptors	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_76.pdf">http://www.chromtech.net.au/pdf2/vici09_76.pdf</a>
CTA4S6	316 Stainless Tube adapters -Male 1/4-28 1/4inx1.5mm bore	Adaptors	\$9.00	<a href="http://www.chromtech.net.au/pdf2/vici09_76.pdf">http://www.chromtech.net.au/pdf2/vici09_76.pdf</a>
CFLAKF	CTFE Luer adapter-Female Luer to 1/4-28 1.50mm bore	Adaptors	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_76.pdf">http://www.chromtech.net.au/pdf2/vici09_76.pdf</a>
CFLAPK	PEEK Luer adapter-Female Luer to 1/4-28 1.50mm bore	Adaptors	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_76.pdf">http://www.chromtech.net.au/pdf2/vici09_76.pdf</a>
CGC41	Acetal Glass Connectors -Complete assembly 1/4-24 female to 1/4inOD glass 1mmID glass tubing	Low Pressure T, X, Manifolds	\$50.00	<a href="http://www.chromtech.net.au/pdf2/vici09_75.pdf">http://www.chromtech.net.au/pdf2/vici09_75.pdf</a>
CGCN	Acetal Glass Connectors -1/4-24 female to 1/4inOD glass -Components Nut	Mixing Tee, Glass Connectors	\$14.00	<a href="http://www.chromtech.net.au/pdf2/vici09_75.pdf">http://www.chromtech.net.au/pdf2/vici09_75.pdf</a>
CM1XKF	Mixing Tee -1/16inx0.75mm bore inc 1/4-28 flangeless fits	Mixing Tee, Glass Connectors	\$82.00	<a href="http://www.chromtech.net.au/pdf2/vici09_75.pdf">http://www.chromtech.net.au/pdf2/vici09_75.pdf</a>
CTA2PK	PEEK Tube adapters -Male 1/4-28 1/8inx1.5mm bore	Mixing Tee, Glass Connectors	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_75.pdf">http://www.chromtech.net.au/pdf2/vici09_75.pdf</a>
CGCRKF	CTFE Glass Connectors -1/4-24 female to 1/4inOD glass -Split rings	Mixing Tee, Glass Connectors	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_75.pdf">http://www.chromtech.net.au/pdf2/vici09_75.pdf</a>
CGCR	Acetal Glass Connectors -1/4-24 female to 1/4inOD glass -Split rings	Mixing Tee, Glass Connectors	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_75.pdf">http://www.chromtech.net.au/pdf2/vici09_75.pdf</a>
CGCG42	Acetal Glass Connectors - 1/4-24 female to 1/4inOD glass -Glass rod 3-1/4inlong 2mmID	Mixing Tee, Glass Connectors	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_75.pdf">http://www.chromtech.net.au/pdf2/vici09_75.pdf</a>
CGCG41	Acetal Glass Connectors -1/4-24 female to 1/4inOD glass -Glass rod 3-1/4inlong 1mmID	Mixing Tee, Glass Connectors	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_75.pdf">http://www.chromtech.net.au/pdf2/vici09_75.pdf</a>
CGCNKF	CTFE Glass Connectors -1/4-24 female to 1/4inOD glass -Components Nut	Mixing Tee, Glass Connectors	\$14.00	<a href="http://www.chromtech.net.au/pdf2/vici09_75.pdf">http://www.chromtech.net.au/pdf2/vici09_75.pdf</a>
CM1XPK	Mixing Tees1/4-28 0.75mm PEEK	Mixing Tee, Glass Connectors	\$82.00	<a href="http://www.chromtech.net.au/pdf2/vici09_75.pdf">http://www.chromtech.net.au/pdf2/vici09_75.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

CGCB	Acetal Glass Connectors -1/4-24 female to 1/4inOD glass -Components Bushing	Mixing Tee, Glass Connectors	\$14.00	<a href="http://www.chromtech.net.au/pdf2/vici09_75.pdf">http://www.chromtech.net.au/pdf2/vici09_75.pdf</a>
CGC42KF	CTFE Glass Connectors -Complete assembly 1/4-24 female to 1/4inOD glass 2mmID glass tubing	Mixing Tee, Glass Connectors	\$56.00	<a href="http://www.chromtech.net.au/pdf2/vici09_75.pdf">http://www.chromtech.net.au/pdf2/vici09_75.pdf</a>
CGC42	Acetal Glass Connectors -Complete assembly 1/4-24 female to 1/4inOD glass 2mmID glass tubing	Mixing Tee, Glass Connectors	\$50.00	<a href="http://www.chromtech.net.au/pdf2/vici09_75.pdf">http://www.chromtech.net.au/pdf2/vici09_75.pdf</a>
CGC41KF	CTFE Glass Connectors -Complete assembly 1/4-24 female to 1/4inOD glass 1mmID glass tubing	Mixing Tee, Glass Connectors	\$56.00	<a href="http://www.chromtech.net.au/pdf2/vici09_75.pdf">http://www.chromtech.net.au/pdf2/vici09_75.pdf</a>
CM2XPK	Mixing Tee -1/8inx1.50mm bore inc 1/4-28 flangeless fits	Mixing Tee, Glass Connectors	\$82.00	<a href="http://www.chromtech.net.au/pdf2/vici09_75.pdf">http://www.chromtech.net.au/pdf2/vici09_75.pdf</a>
CGCBKF	CTFE Glass Connectors -1/4-24 female to 1/4inOD glass -Components Bushing	Mixing Tee, Glass Connectors	\$14.00	<a href="http://www.chromtech.net.au/pdf2/vici09_75.pdf">http://www.chromtech.net.au/pdf2/vici09_75.pdf</a>
CM2XKF	Mixing Tee -1/8inx1.50mm bore inc 1/4-28 flangeless fits	Mixing Tee, Glass Connectors	\$82.00	<a href="http://www.chromtech.net.au/pdf2/vici09_75.pdf">http://www.chromtech.net.au/pdf2/vici09_75.pdf</a>
CTMPK	TEES -1/16inx0.75mm bore inc 1/4-28 flangeless fits PEEK	Low Pressure T, X, Manifolds	\$42.00	<a href="http://www.chromtech.net.au/pdf2/vici09_74.pdf">http://www.chromtech.net.au/pdf2/vici09_74.pdf</a>
C5M1KF	Manifolds 5 ports -1/16inx0.75mm bore inc 1/4-28 flangeless fits	Low Pressure T, X, Manifolds	\$76.00	<a href="http://www.chromtech.net.au/pdf2/vici09_74.pdf">http://www.chromtech.net.au/pdf2/vici09_74.pdf</a>
CTCPK	TEES -1/16inx0.25mm bore inc 1/4-28 flangeless fits PEEK	Low Pressure T, X, Manifolds	\$50.00	<a href="http://www.chromtech.net.au/pdf2/vici09_74.pdf">http://www.chromtech.net.au/pdf2/vici09_74.pdf</a>
CTPK	TEES -1/16inx0.50mm bore inc 1/4-28 flangeless fits PEEK	Low Pressure T, X, Manifolds	\$42.00	<a href="http://www.chromtech.net.au/pdf2/vici09_74.pdf">http://www.chromtech.net.au/pdf2/vici09_74.pdf</a>
C9M2KF	Manifolds 9 ports -1/8inx1.50mm bore inc 1/4-28 flangeless fits	Low Pressure T, X, Manifolds	\$112.00	<a href="http://www.chromtech.net.au/pdf2/vici09_74.pdf">http://www.chromtech.net.au/pdf2/vici09_74.pdf</a>
C9M2PK	Manifolds 9 ports -1/8inx1.50mm bore inc 1/4-28 flangeless fits	Low Pressure T, X, Manifolds	\$112.00	<a href="http://www.chromtech.net.au/pdf2/vici09_74.pdf">http://www.chromtech.net.au/pdf2/vici09_74.pdf</a>
C9M1KF	Manifolds 9 ports -1/16inx0.75mm bore inc 1/4-28 flangeless fits	Low Pressure T, X, Manifolds	\$112.00	<a href="http://www.chromtech.net.au/pdf2/vici09_74.pdf">http://www.chromtech.net.au/pdf2/vici09_74.pdf</a>
C9M1PK	Manifolds 9 ports -1/16inx0.75mm bore inc 1/4-28 flangeless fits	Low Pressure T, X, Manifolds	\$112.00	<a href="http://www.chromtech.net.au/pdf2/vici09_74.pdf">http://www.chromtech.net.au/pdf2/vici09_74.pdf</a>
C5M2PK	Manifolds 5 ports -1/8inx1.50mm bore inc 1/4-28 flangeless fits	Low Pressure T, X, Manifolds	\$76.00	<a href="http://www.chromtech.net.au/pdf2/vici09_74.pdf">http://www.chromtech.net.au/pdf2/vici09_74.pdf</a>
C5M1PK	Manifolds 5 ports -1/16inx0.75mm bore inc 1/4-28 flangeless fits	Low Pressure T, X, Manifolds	\$76.00	<a href="http://www.chromtech.net.au/pdf2/vici09_74.pdf">http://www.chromtech.net.au/pdf2/vici09_74.pdf</a>
CXLKF	Crosses -1/8inx1.50mm bore inc 1/4-28 flangeless fits CTFE	Low Pressure T, X, Manifolds	\$50.00	<a href="http://www.chromtech.net.au/pdf2/vici09_74.pdf">http://www.chromtech.net.au/pdf2/vici09_74.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

CXLPK	Crosses -1/8inx1.50mm bore inc 1/4-28 flangeless fits PEEK	Low Pressure T, X, Manifolds	\$50.00	<a href="http://www.chromtech.net.au/pdf2/vici09_74.pdf">http://www.chromtech.net.au/pdf2/vici09_74.pdf</a>
CTLPK	TEES -1/8inx1.50mm bore inc 1/4-28 flangeless fits PEEK	Low Pressure T, X, Manifolds	\$42.00	<a href="http://www.chromtech.net.au/pdf2/vici09_74.pdf">http://www.chromtech.net.au/pdf2/vici09_74.pdf</a>
CTCKF	TEES -1/16inx0.25mm bore inc 1/4-28 flangeless fits CTFE	Low Pressure T, X, Manifolds	\$50.00	<a href="http://www.chromtech.net.au/pdf2/vici09_74.pdf">http://www.chromtech.net.au/pdf2/vici09_74.pdf</a>
C5M2KF	Manifolds 5 ports -1/8inx1.50mm bore inc 1/4-28 flangeless fits	Low Pressure T, X, Manifolds	\$76.00	<a href="http://www.chromtech.net.au/pdf2/vici09_74.pdf">http://www.chromtech.net.au/pdf2/vici09_74.pdf</a>
CTMKF	TEES -1/16inx0.75mm bore inc 1/4-28 flangeless fits CTFE	Low Pressure T, X, Manifolds	\$42.00	<a href="http://www.chromtech.net.au/pdf2/vici09_74.pdf">http://www.chromtech.net.au/pdf2/vici09_74.pdf</a>
CXMKF	Crosses -1/16inx0.75mm bore inc 1/4-28 flangeless fits CTFE	Low Pressure T, X, Manifolds	\$50.00	<a href="http://www.chromtech.net.au/pdf2/vici09_74.pdf">http://www.chromtech.net.au/pdf2/vici09_74.pdf</a>
CTLKF	TEES -1/8inx1.50mm bore inc 1/4-28 flangeless fits CTFE	Low Pressure T, X, Manifolds	\$42.00	<a href="http://www.chromtech.net.au/pdf2/vici09_74.pdf">http://www.chromtech.net.au/pdf2/vici09_74.pdf</a>
CXCPK	Crosses -1/16inx0.25mm bore inc 1/4-28 flangeless fits PEEK	Low Pressure T, X, Manifolds	\$56.00	<a href="http://www.chromtech.net.au/pdf2/vici09_74.pdf">http://www.chromtech.net.au/pdf2/vici09_74.pdf</a>
CXCKF	Crosses -1/16inx0.25mm bore inc 1/4-28 flangeless fits CTFE	Low Pressure T, X, Manifolds	\$56.00	<a href="http://www.chromtech.net.au/pdf2/vici09_74.pdf">http://www.chromtech.net.au/pdf2/vici09_74.pdf</a>
CXPK	Crosses -1/16inx0.50mm bore inc 1/4-28 flangeless fits PEEK	Low Pressure T, X, Manifolds	\$50.00	<a href="http://www.chromtech.net.au/pdf2/vici09_74.pdf">http://www.chromtech.net.au/pdf2/vici09_74.pdf</a>
CXKF	Crosses -1/16inx0.50mm bore inc 1/4-28 flangeless fits CTFE	Low Pressure T, X, Manifolds	\$50.00	<a href="http://www.chromtech.net.au/pdf2/vici09_74.pdf">http://www.chromtech.net.au/pdf2/vici09_74.pdf</a>
CXMPK	Crosses -1/16inx0.75mm bore inc 1/4-28 flangeless fits PEEK	Low Pressure T, X, Manifolds	\$50.00	<a href="http://www.chromtech.net.au/pdf2/vici09_74.pdf">http://www.chromtech.net.au/pdf2/vici09_74.pdf</a>
CTKF	TEES -1/16inx0.50mm bore inc 1/4-28 flangeless fits CTFE	Low Pressure T, X, Manifolds	\$42.00	<a href="http://www.chromtech.net.au/pdf2/vici09_74.pdf">http://www.chromtech.net.au/pdf2/vici09_74.pdf</a>
CZBUS6	316 Stainless Bulkhead Unions -1/4-28 to 10-32 1/16inx0.50mm bore	Low Press Bulkhead Union	\$26.00	<a href="http://www.chromtech.net.au/pdf2/vici09_73.pdf">http://www.chromtech.net.au/pdf2/vici09_73.pdf</a>
CBULPK	PEEK Bulkhead Unions -1/4-28 1/8inx1.50mm bore	Low Press Bulkhead Union	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_73.pdf">http://www.chromtech.net.au/pdf2/vici09_73.pdf</a>
CZBUMS6	316 Stainless Bulkhead Unions -1/4-28 to 10-32 1/16inx0.75mm bore	Low Press Bulkhead Union	\$26.00	<a href="http://www.chromtech.net.au/pdf2/vici09_73.pdf">http://www.chromtech.net.au/pdf2/vici09_73.pdf</a>
CZBUMKF	CTFE Bulkhead Unions -1/4-28 to 10-32 1/16inx0.75mm bore	Low Press Bulkhead Union	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_73.pdf">http://www.chromtech.net.au/pdf2/vici09_73.pdf</a>
CZBUMPK	PEEK Bulkhead Unions -1/4-28 to 10-32 1/16inx0.75mm bore	Low Press Bulkhead Union	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_73.pdf">http://www.chromtech.net.au/pdf2/vici09_73.pdf</a>



## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

CZBUKF	CTFE Bulkhead Unions -1/4-28 to 10-32 1/16inx0.50mm bore	Low Press Bulkhead Union	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_73.pdf">http://www.chromtech.net.au/pdf2/vici09_73.pdf</a>
CZBUCS6	316 Stainless Bulkhead Unions -1/4-28 to 10-32 1/16inx0.25mm bore	Low Press Bulkhead Union	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_73.pdf">http://www.chromtech.net.au/pdf2/vici09_73.pdf</a>
CZBUCKF	CTFE Bulkhead Unions -1/4-28 to 10-32 1/16inx0.25mm bore	Low Press Bulkhead Union	\$36.00	<a href="http://www.chromtech.net.au/pdf2/vici09_73.pdf">http://www.chromtech.net.au/pdf2/vici09_73.pdf</a>
CZBUCPK	PEEK Bulkhead Unions -1/4-28 to 10-32 1/16inx0.25mm bore	Low Press Bulkhead Union	\$36.00	<a href="http://www.chromtech.net.au/pdf2/vici09_73.pdf">http://www.chromtech.net.au/pdf2/vici09_73.pdf</a>
CZBUPK	PEEK Bulkhead Unions -1/4-28 to 10-32 1/16inx0.50mm bore	Low Press Bulkhead Union	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_73.pdf">http://www.chromtech.net.au/pdf2/vici09_73.pdf</a>
CBULKF	CTFE Bulkhead Unions -1/4-28 1/8inx1.50mm bore	Low Press Bulkhead Union	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_73.pdf">http://www.chromtech.net.au/pdf2/vici09_73.pdf</a>
CBUCPK	PEEK Bulkhead Unions -1/4-28 1/16inx0.25mm bore	Low Press Bulkhead Union	\$36.00	<a href="http://www.chromtech.net.au/pdf2/vici09_73.pdf">http://www.chromtech.net.au/pdf2/vici09_73.pdf</a>
CBUMS6	316 Stainless Bulkhead Unions -1/4-28 1/16inx0.75mm bore	Low Press Bulkhead Union	\$26.00	<a href="http://www.chromtech.net.au/pdf2/vici09_73.pdf">http://www.chromtech.net.au/pdf2/vici09_73.pdf</a>
CBUMKF	CTFE Bulkhead Unions -1/4-28 1/16inx0.75mm bore	Low Press Bulkhead Union	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_73.pdf">http://www.chromtech.net.au/pdf2/vici09_73.pdf</a>
CBUMPK	PEEK Bulkhead Unions -1/4-28 1/16inx0.75mm bore	Low Press Bulkhead Union	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_73.pdf">http://www.chromtech.net.au/pdf2/vici09_73.pdf</a>
CBUS6	316 Stainless Bulkhead Unions -1/4-28 1/16inx0.50mm bore	Low Press Bulkhead Union	\$26.00	<a href="http://www.chromtech.net.au/pdf2/vici09_73.pdf">http://www.chromtech.net.au/pdf2/vici09_73.pdf</a>
CBUKF	CTFE Bulkhead Unions -1/4-28 1/16inx0.50mm bore	Low Press Bulkhead Union	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_73.pdf">http://www.chromtech.net.au/pdf2/vici09_73.pdf</a>
CBUPK	PEEK Bulkhead Unions -1/4-28 1/16inx0.50mm bore	Low Press Bulkhead Union	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_73.pdf">http://www.chromtech.net.au/pdf2/vici09_73.pdf</a>
CBUCS6	316 Stainless Bulkhead Unions -1/4-28 1/16inx0.25mm bore	Low Press Bulkhead Union	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_73.pdf">http://www.chromtech.net.au/pdf2/vici09_73.pdf</a>
CBUCKF	CTFE Bulkhead Unions -1/4-28 1/16inx0.25mm bore	Low Press Bulkhead Union	\$36.00	<a href="http://www.chromtech.net.au/pdf2/vici09_73.pdf">http://www.chromtech.net.au/pdf2/vici09_73.pdf</a>
CBULS6	316 Stainless Bulkhead Unions -1/4-28 1/8inx1.50mm bore	Low Press Bulkhead Union	\$26.00	<a href="http://www.chromtech.net.au/pdf2/vici09_73.pdf">http://www.chromtech.net.au/pdf2/vici09_73.pdf</a>
CUMKF	Union 1/4-28-1/4-28,CTFE,0.75mm bore/1/16in Tubing	Low Press Plugs/Caps	\$26.00	<a href="http://www.chromtech.net.au/pdf2/vici09_72.pdf">http://www.chromtech.net.au/pdf2/vici09_72.pdf</a>
CZUMKF	Union 1/4-28-10-32 1/16 0.75mm CTFE	Low Press Unions	\$26.00	<a href="http://www.chromtech.net.au/pdf2/vici09_72.pdf">http://www.chromtech.net.au/pdf2/vici09_72.pdf</a>
CUMPK	Union 1/4-28-1/4-28,PEEK,0.75mm bore/1/16in Tubing	Low Press Unions	\$26.00	<a href="http://www.chromtech.net.au/pdf2/vici09_72.pdf">http://www.chromtech.net.au/pdf2/vici09_72.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

CULPK	Union 1/4-28-1/4-28,PEEK,1.50mm bore/1/8in Tubing	Low Press Unions	\$26.00	<a href="http://www.chromtech.net.au/pdf2/vici09_72.pdf">http://www.chromtech.net.au/pdf2/vici09_72.pdf</a>
CUTPK	1/8" Flanged Butt Union PEEK	Low Press Unions	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_72.pdf">http://www.chromtech.net.au/pdf2/vici09_72.pdf</a>
CUTKF	1/8" Flanged Butt Union CTFE	Low Press Unions	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_72.pdf">http://www.chromtech.net.au/pdf2/vici09_72.pdf</a>
CUTPP*	1/8" Flanged Butt Union PP	Low Press Unions	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_72.pdf">http://www.chromtech.net.au/pdf2/vici09_72.pdf</a>
CZUCPK	Union 1/4-28-10-32 1/16 0.25mm PEEK	Low Press Unions	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_72.pdf">http://www.chromtech.net.au/pdf2/vici09_72.pdf</a>
CZUCKF	Union 1/4-28-10-32 1/16 0.25mm CTFE	Low Press Unions	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_72.pdf">http://www.chromtech.net.au/pdf2/vici09_72.pdf</a>
CZUCS6	Union 1/4-28-10-32 1/16 0.25mm 316SS	Low Press Unions	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_72.pdf">http://www.chromtech.net.au/pdf2/vici09_72.pdf</a>
CZUPK	Union 1/4-28-10-32 1/16 0.50mm PEEK	Low Press Unions	\$26.00	<a href="http://www.chromtech.net.au/pdf2/vici09_72.pdf">http://www.chromtech.net.au/pdf2/vici09_72.pdf</a>
CZUKF	Union 1/4-28-10-32 1/16 0.50mm CTFE	Low Press Unions	\$26.00	<a href="http://www.chromtech.net.au/pdf2/vici09_72.pdf">http://www.chromtech.net.au/pdf2/vici09_72.pdf</a>
CULKF	Union 1/4-28-1/4-28,CTFE,1.50mm bore/1/8in Tubing	Low Press Unions	\$26.00	<a href="http://www.chromtech.net.au/pdf2/vici09_72.pdf">http://www.chromtech.net.au/pdf2/vici09_72.pdf</a>
CZUMPK	Union 1/4-28-10-32 1/16 0.75mm PEEK	Low Press Unions	\$26.00	<a href="http://www.chromtech.net.au/pdf2/vici09_72.pdf">http://www.chromtech.net.au/pdf2/vici09_72.pdf</a>
CUCPK	Union 1/4-28-1/4-28,PEEK, 0.25mm bore,1/16in Tubing	Low Press Unions	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_72.pdf">http://www.chromtech.net.au/pdf2/vici09_72.pdf</a>
CZUMS6	316 SS Union Fitting 1/4-28 to 10-32 for 1/16" OD Tubing 0.75mm bore	Low Press Unions	\$26.00	<a href="http://www.chromtech.net.au/pdf2/vici09_72.pdf">http://www.chromtech.net.au/pdf2/vici09_72.pdf</a>
CU4LPK	Union 1/4-28-1/2-20,PEEK, 1.50mm bore	Low Press Unions	\$57.00	<a href="http://www.chromtech.net.au/pdf2/vici09_72.pdf">http://www.chromtech.net.au/pdf2/vici09_72.pdf</a>
CU4LKF	Union 1/4-28-1/2-20, CTFE, 1.50mm bore	Low Press Unions	\$57.00	<a href="http://www.chromtech.net.au/pdf2/vici09_72.pdf">http://www.chromtech.net.au/pdf2/vici09_72.pdf</a>
CFL-4KF	Nut 1/2-20 CTFE 1/4-28-1/2-20 Union	Low Press Unions	\$10.00	<a href="http://www.chromtech.net.au/pdf2/vici09_72.pdf">http://www.chromtech.net.au/pdf2/vici09_72.pdf</a>
CFL-4D	Nut 1/2-20 Delrin 1/4-28-1/2-20 Union	Low Press Unions	\$8.00	<a href="http://www.chromtech.net.au/pdf2/vici09_72.pdf">http://www.chromtech.net.au/pdf2/vici09_72.pdf</a>
CFL-CB4KF-S	CTFE Ferrule 1/4-28-1/2-20 Union	Low Press Unions	\$6.00	<a href="http://www.chromtech.net.au/pdf2/vici09_72.pdf">http://www.chromtech.net.au/pdf2/vici09_72.pdf</a>
CUKF	Union 1/4-28-1/4-28,CTFE, 0.50mm bore,1/16in Tubing	Low Press Unions	\$26.00	<a href="http://www.chromtech.net.au/pdf2/vici09_72.pdf">http://www.chromtech.net.au/pdf2/vici09_72.pdf</a>
CUPK	Union 1/4-28-1/4-28,CTFE,0.50mm bore/1/16in Tubing	Low Press Unions	\$26.00	<a href="http://www.chromtech.net.au/pdf2/vici09_72.pdf">http://www.chromtech.net.au/pdf2/vici09_72.pdf</a>
CUCKF	Union 1/4-28-1/4-28,CTFE, 0.25mm bore,1/16in Tubing	Low Press Unions	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_72.pdf">http://www.chromtech.net.au/pdf2/vici09_72.pdf</a>
CZUS6	Union 1/4-8-10-32 1/16 0.75mm 316SS	Low Press Unions	\$26.00	<a href="http://www.chromtech.net.au/pdf2/vici09_72.pdf">http://www.chromtech.net.au/pdf2/vici09_72.pdf</a>
JR-201540	Easy-Flange kit in plastic case	HPLC - Easy-Flange	\$216.00	<a href="http://www.chromtech.net.au/pdf2/vici09_71.pdf">http://www.chromtech.net.au/pdf2/vici09_71.pdf</a>
CCKF-5	CTFE Caps 1/4-28	Low Press Plugs/Caps	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_71.pdf">http://www.chromtech.net.au/pdf2/vici09_71.pdf</a>
CPPK	PEEK Plugs 1/4-28	Low Press Plugs/Caps	\$16.00	<a href="http://www.chromtech.net.au/pdf2/vici09_71.pdf">http://www.chromtech.net.au/pdf2/vici09_71.pdf</a>
CPKF	CTFE Plugs 1/4-28	Low Press Plugs/Caps	\$16.00	<a href="http://www.chromtech.net.au/pdf2/vici09_71.pdf">http://www.chromtech.net.au/pdf2/vici09_71.pdf</a>
MZP1PK	PEEK 1/16in x.500in long hex Plugs	Low Press Plugs/Caps	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_71.pdf">http://www.chromtech.net.au/pdf2/vici09_71.pdf</a>
LZP1PK	PEEK 1/16inx.880in long hex	Low Press Plugs/Caps	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_71.pdf">http://www.chromtech.net.au/pdf2/vici09_71.pdf</a>
ZP1FPK	PEEK 1/16in fingertight Plugs	Low Press Plugs/Caps	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_71.pdf">http://www.chromtech.net.au/pdf2/vici09_71.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

CCPK-5	PEEK Caps 1/4-28	Low Press Plugs/Caps	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_71.pdf">http://www.chromtech.net.au/pdf2/vici09_71.pdf</a>
CFT-H	Accessories -Tubing holder	Flanging Tool	\$26.00	<a href="http://www.chromtech.net.au/pdf2/vici09_70.pdf">http://www.chromtech.net.au/pdf2/vici09_70.pdf</a>
CFT1K-110	Starter Kit -1/16in tubing	Flanging Tool	\$420.00	<a href="http://www.chromtech.net.au/pdf2/vici09_70.pdf">http://www.chromtech.net.au/pdf2/vici09_70.pdf</a>
CFT1K-220	Starter Kit -1/16in tubing	Flanging Tool	\$420.00	<a href="http://www.chromtech.net.au/pdf2/vici09_70.pdf">http://www.chromtech.net.au/pdf2/vici09_70.pdf</a>
CFT2K-220	Starter Kit -1/8in tubing	Flanging Tool	\$420.00	<a href="http://www.chromtech.net.au/pdf2/vici09_70.pdf">http://www.chromtech.net.au/pdf2/vici09_70.pdf</a>
CFT-R	Accessories -Razor blades	Flanging Tool	\$8.00	<a href="http://www.chromtech.net.au/pdf2/vici09_70.pdf">http://www.chromtech.net.au/pdf2/vici09_70.pdf</a>
CFT-220	Flanging tools -230 VAC	Flanging Tool	\$274.00	<a href="http://www.chromtech.net.au/pdf2/vici09_70.pdf">http://www.chromtech.net.au/pdf2/vici09_70.pdf</a>
CFT2K-110	Starter Kit -1/8in tubing	Flanging Tool	\$420.00	<a href="http://www.chromtech.net.au/pdf2/vici09_70.pdf">http://www.chromtech.net.au/pdf2/vici09_70.pdf</a>
CFT-TXL	Accessories -Flanging tips -tubing ID <2.00mm	Flanging Tool	\$36.00	<a href="http://www.chromtech.net.au/pdf2/vici09_70.pdf">http://www.chromtech.net.au/pdf2/vici09_70.pdf</a>
CFT-TL	Accessories -Flanging tips -tubing ID <1.50mm	Flanging Tool	\$36.00	<a href="http://www.chromtech.net.au/pdf2/vici09_70.pdf">http://www.chromtech.net.au/pdf2/vici09_70.pdf</a>
CFT-TM	Accessories -Flanging tips -tubing ID <1.00mm	Flanging Tool	\$36.00	<a href="http://www.chromtech.net.au/pdf2/vici09_70.pdf">http://www.chromtech.net.au/pdf2/vici09_70.pdf</a>
CFT-TXC	Accessories -Flanging tips -tubing ID <0.25mm	Flanging Tool	\$36.00	<a href="http://www.chromtech.net.au/pdf2/vici09_70.pdf">http://www.chromtech.net.au/pdf2/vici09_70.pdf</a>
CFT-110	Flanging tools -110 VAC	Flanging Tool	\$274.00	<a href="http://www.chromtech.net.au/pdf2/vici09_70.pdf">http://www.chromtech.net.au/pdf2/vici09_70.pdf</a>
CFT-TC	Accessories -Flanging tips -tubing ID <0.75mm	Flanging Tool	\$36.00	<a href="http://www.chromtech.net.au/pdf2/vici09_70.pdf">http://www.chromtech.net.au/pdf2/vici09_70.pdf</a>
CF-2R	Red-pp End fittings -1/8inOD	LowPress Std TubeEnd Fits/Ext Nuts	\$16.00	<a href="http://www.chromtech.net.au/pdf2/vici09_69.pdf">http://www.chromtech.net.au/pdf2/vici09_69.pdf</a>
CF-2A	Assorted-pp End fittings -1/8inOD	LowPress Std TubeEnd Fits/Ext Nuts	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_69.pdf">http://www.chromtech.net.au/pdf2/vici09_69.pdf</a>
CF-1Y	Yellow-pp End fittings -1/16inOD	LowPress Std TubeEnd Fits/Ext Nuts	\$16.00	<a href="http://www.chromtech.net.au/pdf2/vici09_69.pdf">http://www.chromtech.net.au/pdf2/vici09_69.pdf</a>
CF-2Y	Yellow-pp End fittings -1/8inOD	LowPress Std TubeEnd Fits/Ext Nuts	\$16.00	<a href="http://www.chromtech.net.au/pdf2/vici09_69.pdf">http://www.chromtech.net.au/pdf2/vici09_69.pdf</a>
CF-1A	Assorted-pp End fittings -1/16inOD	LowPress Std TubeEnd Fits/Ext Nuts	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_69.pdf">http://www.chromtech.net.au/pdf2/vici09_69.pdf</a>
CF-2W	White-pp End fittings -1/8inOD	LowPress Std TubeEnd Fits/Ext Nuts	\$16.00	<a href="http://www.chromtech.net.au/pdf2/vici09_69.pdf">http://www.chromtech.net.au/pdf2/vici09_69.pdf</a>
CF-W1	End fittings -1/16inOD -Washers	LowPress Std TubeEnd Fits/Ext Nuts	\$4.00	<a href="http://www.chromtech.net.au/pdf2/vici09_69.pdf">http://www.chromtech.net.au/pdf2/vici09_69.pdf</a>
CF-W2	End fittings -1/8inOD Washers	LowPress Std TubeEnd Fits/Ext Nuts	\$4.00	<a href="http://www.chromtech.net.au/pdf2/vici09_69.pdf">http://www.chromtech.net.au/pdf2/vici09_69.pdf</a>
CEN1PK	PEEK External 1/16in Nuts for flanged tube ends	LowPress Std TubeEnd Fits/Ext Nuts	\$16.00	<a href="http://www.chromtech.net.au/pdf2/vici09_69.pdf">http://www.chromtech.net.au/pdf2/vici09_69.pdf</a>
CEN2PK	PEEK External 1/8in Nuts for flanged tube ends	LowPress Std TubeEnd Fits/Ext Nuts	\$16.00	<a href="http://www.chromtech.net.au/pdf2/vici09_69.pdf">http://www.chromtech.net.au/pdf2/vici09_69.pdf</a>
CEN2KF	CTFE External 1/8in Nuts for flanged tube ends	LowPress Std TubeEnd Fits/Ext Nuts	\$16.00	<a href="http://www.chromtech.net.au/pdf2/vici09_69.pdf">http://www.chromtech.net.au/pdf2/vici09_69.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

CF-2P	Purple-pp End fittings -1/8inOD	LowPress Std TubeEnd Fits/Ext Nuts	\$16.00	<a href="http://www.chromtech.net.au/pdf2/vici09_69.pdf">http://www.chromtech.net.au/pdf2/vici09_69.pdf</a>
CF-1R	Red-pp End fittings -1/16inOD	LowPress Std TubeEnd Fits/Ext Nuts	\$16.00	<a href="http://www.chromtech.net.au/pdf2/vici09_69.pdf">http://www.chromtech.net.au/pdf2/vici09_69.pdf</a>
CEN1KF	CTFE External 1/16in Nuts for flanged tube ends	LowPress Std TubeEnd Fits/Ext Nuts	\$16.00	<a href="http://www.chromtech.net.au/pdf2/vici09_69.pdf">http://www.chromtech.net.au/pdf2/vici09_69.pdf</a>
CF-2DG	Dark gray-pp End fittings -1/8inOD	LowPress Std TubeEnd Fits/Ext Nuts	\$16.00	<a href="http://www.chromtech.net.au/pdf2/vici09_69.pdf">http://www.chromtech.net.au/pdf2/vici09_69.pdf</a>
CF-1W	White-pp End fittings -1/16inOD	LowPress Std TubeEnd Fits/Ext Nuts	\$16.00	<a href="http://www.chromtech.net.au/pdf2/vici09_69.pdf">http://www.chromtech.net.au/pdf2/vici09_69.pdf</a>
CF-1P	Purple-pp End fittings -1/16inOD	LowPress Std TubeEnd Fits/Ext Nuts	\$16.00	<a href="http://www.chromtech.net.au/pdf2/vici09_69.pdf">http://www.chromtech.net.au/pdf2/vici09_69.pdf</a>
CF-2BK	Black-pp End fittings -1/8inOD	LowPress Std TubeEnd Fits/Ext Nuts	\$16.00	<a href="http://www.chromtech.net.au/pdf2/vici09_69.pdf">http://www.chromtech.net.au/pdf2/vici09_69.pdf</a>
CF-1BE	Blue-pp End fittings -1/16inOD	LowPress Std TubeEnd Fits/Ext Nuts	\$16.00	<a href="http://www.chromtech.net.au/pdf2/vici09_69.pdf">http://www.chromtech.net.au/pdf2/vici09_69.pdf</a>
CF-2BE	Blue-pp End fittings -1/8inOD	LowPress Std TubeEnd Fits/Ext Nuts	\$16.00	<a href="http://www.chromtech.net.au/pdf2/vici09_69.pdf">http://www.chromtech.net.au/pdf2/vici09_69.pdf</a>
CF-1BR	Brown-pp End fittings -1/16inOD	LowPress Std TubeEnd Fits/Ext Nuts	\$16.00	<a href="http://www.chromtech.net.au/pdf2/vici09_69.pdf">http://www.chromtech.net.au/pdf2/vici09_69.pdf</a>
CF-1DG	Dark gray-pp End fittings -1/16inOD	LowPress Std TubeEnd Fits/Ext Nuts	\$16.00	<a href="http://www.chromtech.net.au/pdf2/vici09_69.pdf">http://www.chromtech.net.au/pdf2/vici09_69.pdf</a>
CF-1G	Green -pp End fittings -1/16inOD	LowPress Std TubeEnd Fits/Ext Nuts	\$16.00	<a href="http://www.chromtech.net.au/pdf2/vici09_69.pdf">http://www.chromtech.net.au/pdf2/vici09_69.pdf</a>
CF-1E	Orange-pp End fittings -1/16inOD	LowPress Std TubeEnd Fits/Ext Nuts	\$16.00	<a href="http://www.chromtech.net.au/pdf2/vici09_69.pdf">http://www.chromtech.net.au/pdf2/vici09_69.pdf</a>
CF-1L	Lavendar-pp End fittings -1/16inOD	LowPress Std TubeEnd Fits/Ext Nuts	\$16.00	<a href="http://www.chromtech.net.au/pdf2/vici09_69.pdf">http://www.chromtech.net.au/pdf2/vici09_69.pdf</a>
CF-2L	Lavendar-pp End fittings -1/8inOD	LowPress Std TubeEnd Fits/Ext Nuts	\$16.00	<a href="http://www.chromtech.net.au/pdf2/vici09_69.pdf">http://www.chromtech.net.au/pdf2/vici09_69.pdf</a>
CF-1N	Natural-pp End fittings -1/16inOD	LowPress Std TubeEnd Fits/Ext Nuts	\$16.00	<a href="http://www.chromtech.net.au/pdf2/vici09_69.pdf">http://www.chromtech.net.au/pdf2/vici09_69.pdf</a>
CF-2N	Natural-pp End fittings -1/8inOD	LowPress Std TubeEnd Fits/Ext Nuts	\$16.00	<a href="http://www.chromtech.net.au/pdf2/vici09_69.pdf">http://www.chromtech.net.au/pdf2/vici09_69.pdf</a>
CF-2G	Green -pp End fittings -1/8inOD	LowPress Std TubeEnd Fits/Ext Nuts	\$16.00	<a href="http://www.chromtech.net.au/pdf2/vici09_69.pdf">http://www.chromtech.net.au/pdf2/vici09_69.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

CF-2E	Orange-pp End fittings -1/8inOD	LowPress Std TubeEnd Fits/Ext Nuts	\$16.00	<a href="http://www.chromtech.net.au/pdf2/vici09_69.pdf">http://www.chromtech.net.au/pdf2/vici09_69.pdf</a>
CF-2BR	Brown-pp End fittings -1/8inOD	LowPress Std TubeEnd Fits/Ext Nuts	\$16.00	<a href="http://www.chromtech.net.au/pdf2/vici09_69.pdf">http://www.chromtech.net.au/pdf2/vici09_69.pdf</a>
CF-1BK	Black-pp End fittings -1/16inOD	LowPress Std TubeEnd Fits/Ext Nuts	\$16.00	<a href="http://www.chromtech.net.au/pdf2/vici09_69.pdf">http://www.chromtech.net.au/pdf2/vici09_69.pdf</a>
CFL-2A	Tube end fittings -1/16in Assorted(1 of each color)	LowPress Flangeless TubeEnd Fits	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_68.pdf">http://www.chromtech.net.au/pdf2/vici09_68.pdf</a>
CFL-1W	Tube end fittings -1/16in OD White	LowPress Flangeless TubeEnd Fits	\$14.00	<a href="http://www.chromtech.net.au/pdf2/vici09_68.pdf">http://www.chromtech.net.au/pdf2/vici09_68.pdf</a>
CFL-2W	Tube end fittings -1/8in OD White	LowPress Flangeless TubeEnd Fits	\$14.00	<a href="http://www.chromtech.net.au/pdf2/vici09_68.pdf">http://www.chromtech.net.au/pdf2/vici09_68.pdf</a>
CFL-1PK	Tube end fittings -PEEK nuts 1/16in	LowPress Flangeless TubeEnd Fits	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_68.pdf">http://www.chromtech.net.au/pdf2/vici09_68.pdf</a>
CFL-1Y	Tube end fittings -1/16in OD Yellow	LowPress Flangeless TubeEnd Fits	\$14.00	<a href="http://www.chromtech.net.au/pdf2/vici09_68.pdf">http://www.chromtech.net.au/pdf2/vici09_68.pdf</a>
CFL-2Y	Tube end fittings -1/8in OD Yellow	LowPress Flangeless TubeEnd Fits	\$14.00	<a href="http://www.chromtech.net.au/pdf2/vici09_68.pdf">http://www.chromtech.net.au/pdf2/vici09_68.pdf</a>
CFL-1A	Tube end fittings -1/16in Assorted(1 of each color)	LowPress Flangeless TubeEnd Fits	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_68.pdf">http://www.chromtech.net.au/pdf2/vici09_68.pdf</a>
CFL-2R	Tube end fittings -1/8in OD Red	LowPress Flangeless TubeEnd Fits	\$14.00	<a href="http://www.chromtech.net.au/pdf2/vici09_68.pdf">http://www.chromtech.net.au/pdf2/vici09_68.pdf</a>
CFL-1A-PK	Flangeless Tube End Fittings12-Colors with PEEK Collapsible Ferrules 1/16,	LowPress Flangeless TubeEnd Fits	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_68.pdf">http://www.chromtech.net.au/pdf2/vici09_68.pdf</a>
CFL-1R	Tube end fittings -1/16in OD Red	LowPress Flangeless TubeEnd Fits	\$14.00	<a href="http://www.chromtech.net.au/pdf2/vici09_68.pdf">http://www.chromtech.net.au/pdf2/vici09_68.pdf</a>
CFL-CB1PK	1/16" Collapsible PEEK Ferrules	LowPress Flangeless TubeEnd Fits	\$15.00	<a href="http://www.chromtech.net.au/pdf2/vici09_68.pdf">http://www.chromtech.net.au/pdf2/vici09_68.pdf</a>
CFL-CB2PK	1/8" Collapsible PEEK Ferrules	LowPress Flangeless TubeEnd Fits	\$15.00	<a href="http://www.chromtech.net.au/pdf2/vici09_68.pdf">http://www.chromtech.net.au/pdf2/vici09_68.pdf</a>
CFL-CB2KF	1/8in Collapsible CTFE Ferrule	LowPress Flangeless TubeEnd Fits	\$15.00	<a href="http://www.chromtech.net.au/pdf2/vici09_68.pdf">http://www.chromtech.net.au/pdf2/vici09_68.pdf</a>
CFL-2PK	Tube end fittings -PEEK nuts 1/16in	LowPress Flangeless TubeEnd Fits	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_68.pdf">http://www.chromtech.net.au/pdf2/vici09_68.pdf</a>
CST	Tube end fittings -Setting tool 1/16in	LowPress Flangeless TubeEnd Fits	\$8.00	<a href="http://www.chromtech.net.au/pdf2/vici09_68.pdf">http://www.chromtech.net.au/pdf2/vici09_68.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

CFL-CB1KF	1/16" Collapsible CTFE Ferrules	LowPress Flangeless TubeEnd Fits	\$15.00	<a href="http://www.chromtech.net.au/pdf2/vici09_68.pdf">http://www.chromtech.net.au/pdf2/vici09_68.pdf</a>
CFL-1BE	Tube end fittings -1/16in OD Blue	LowPress Flangeless TubeEnd Fits	\$14.00	<a href="http://www.chromtech.net.au/pdf2/vici09_68.pdf">http://www.chromtech.net.au/pdf2/vici09_68.pdf</a>
CFL-2P	Tube end fittings -1/8in OD Purple	LowPress Flangeless TubeEnd Fits	\$14.00	<a href="http://www.chromtech.net.au/pdf2/vici09_68.pdf">http://www.chromtech.net.au/pdf2/vici09_68.pdf</a>
CFL-2A-PK	Tube end fittings -1/8in OD PEEK(1 of each color)	LowPress Flangeless TubeEnd Fits	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_68.pdf">http://www.chromtech.net.au/pdf2/vici09_68.pdf</a>
CFL-2BK	Tube end fittings -1/8in OD Black	LowPress Flangeless TubeEnd Fits	\$14.00	<a href="http://www.chromtech.net.au/pdf2/vici09_68.pdf">http://www.chromtech.net.au/pdf2/vici09_68.pdf</a>
CFL-2BE	Tube end fittings -1/8in OD Blue	LowPress Flangeless TubeEnd Fits	\$14.00	<a href="http://www.chromtech.net.au/pdf2/vici09_68.pdf">http://www.chromtech.net.au/pdf2/vici09_68.pdf</a>
CFL-1BR	Tube end fittings -1/16in OD Brown	LowPress Flangeless TubeEnd Fits	\$14.00	<a href="http://www.chromtech.net.au/pdf2/vici09_68.pdf">http://www.chromtech.net.au/pdf2/vici09_68.pdf</a>
CFL-2BR	Tube end fittings -1/8in OD Brown	LowPress Flangeless TubeEnd Fits	\$14.00	<a href="http://www.chromtech.net.au/pdf2/vici09_68.pdf">http://www.chromtech.net.au/pdf2/vici09_68.pdf</a>
CFL-1DG	Tube end fittings -1/16in OD Dark gray	LowPress Flangeless TubeEnd Fits	\$14.00	<a href="http://www.chromtech.net.au/pdf2/vici09_68.pdf">http://www.chromtech.net.au/pdf2/vici09_68.pdf</a>
CFL-2DG	Tube end fittings -1/8in OD	LowPress Flangeless TubeEnd Fits	\$14.00	<a href="http://www.chromtech.net.au/pdf2/vici09_68.pdf">http://www.chromtech.net.au/pdf2/vici09_68.pdf</a>
CFL-1E	Tube end fittings -1/16in OD Orange	LowPress Flangeless TubeEnd Fits	\$14.00	<a href="http://www.chromtech.net.au/pdf2/vici09_68.pdf">http://www.chromtech.net.au/pdf2/vici09_68.pdf</a>
CFL-1BK	Tube end fittings -1/16in OD Black	LowPress Flangeless TubeEnd Fits	\$14.00	<a href="http://www.chromtech.net.au/pdf2/vici09_68.pdf">http://www.chromtech.net.au/pdf2/vici09_68.pdf</a>
CFL-1G	Tube end fittings -1/16in OD Green	LowPress Flangeless TubeEnd Fits	\$14.00	<a href="http://www.chromtech.net.au/pdf2/vici09_68.pdf">http://www.chromtech.net.au/pdf2/vici09_68.pdf</a>
CFL-2E	Tube end fittings -1/8in OD Orange	LowPress Flangeless TubeEnd Fits	\$14.00	<a href="http://www.chromtech.net.au/pdf2/vici09_68.pdf">http://www.chromtech.net.au/pdf2/vici09_68.pdf</a>
CFL-1P	Tube end fittings -1/16in OD Purple	LowPress Flangeless TubeEnd Fits	\$14.00	<a href="http://www.chromtech.net.au/pdf2/vici09_68.pdf">http://www.chromtech.net.au/pdf2/vici09_68.pdf</a>
CFL-2N	Tube end fittings -1/8in OD Natural	LowPress Flangeless TubeEnd Fits	\$14.00	<a href="http://www.chromtech.net.au/pdf2/vici09_68.pdf">http://www.chromtech.net.au/pdf2/vici09_68.pdf</a>
CFL-1N	Tube end fittings -1/16in OD Natural	LowPress Flangeless TubeEnd Fits	\$14.00	<a href="http://www.chromtech.net.au/pdf2/vici09_68.pdf">http://www.chromtech.net.au/pdf2/vici09_68.pdf</a>
CFL-2L	Tube end fittings -1/8in OD Lavendar	LowPress Flangeless TubeEnd Fits	\$14.00	<a href="http://www.chromtech.net.au/pdf2/vici09_68.pdf">http://www.chromtech.net.au/pdf2/vici09_68.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

CFL-1L	Tube end fittings -1/16in OD Lavendar	LowPress Flangeless TubeEnd Fits	\$14.00	<a href="http://www.chromtech.net.au/pdf2/vici09_68.pdf">http://www.chromtech.net.au/pdf2/vici09_68.pdf</a>
CFL-2G	Tube end fittings -1/8in OD Green	LowPress Flangeless TubeEnd Fits	\$14.00	<a href="http://www.chromtech.net.au/pdf2/vici09_68.pdf">http://www.chromtech.net.au/pdf2/vici09_68.pdf</a>
CNFT	Tightening Tool for 1/32in PEEK Fittings Green	CI Fittings Starter Kit	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_67.pdf">http://www.chromtech.net.au/pdf2/vici09_67.pdf</a>
JR-35P	PEEK Starter Kit	CI Fittings Starter Kit	\$309.00	<a href="http://www.chromtech.net.au/pdf2/vici09_67.pdf">http://www.chromtech.net.au/pdf2/vici09_67.pdf</a>
C360ET	Tightening Tool for 360um PEEK Fittings Red	CI Fittings Starter Kit	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_67.pdf">http://www.chromtech.net.au/pdf2/vici09_67.pdf</a>
ZNFT	Tightening Tool for 1/16in PEEK Fittings Blue	Microbore UHPLC 1/16	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_67.pdf">http://www.chromtech.net.au/pdf2/vici09_67.pdf</a>
JR-55025-5	Fitting PEEK one-piece green 10-32 5/pkg	High Pressure PEEK Fittings	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_66.pdf">http://www.chromtech.net.au/pdf2/vici09_66.pdf</a>
JR-55016-12	Color-It fingertight adapter multi-color assorted 2 pcs. each color 12/pkg	High Pressure PEEK Fittings	\$10.00	<a href="http://www.chromtech.net.au/pdf2/vici09_66.pdf">http://www.chromtech.net.au/pdf2/vici09_66.pdf</a>
JR-55015-5	Color-It fingertight adapter red 5/pkg	High Pressure PEEK Fittings	\$4.00	<a href="http://www.chromtech.net.au/pdf2/vici09_66.pdf">http://www.chromtech.net.au/pdf2/vici09_66.pdf</a>
JR-55014-5	Color-It fingertight adapter white 5/pkg	High Pressure PEEK Fittings	\$4.00	<a href="http://www.chromtech.net.au/pdf2/vici09_66.pdf">http://www.chromtech.net.au/pdf2/vici09_66.pdf</a>
JR-55013-5	Color-It fingertight adapter black 5/pkg	High Pressure PEEK Fittings	\$4.00	<a href="http://www.chromtech.net.au/pdf2/vici09_66.pdf">http://www.chromtech.net.au/pdf2/vici09_66.pdf</a>
JR-55012-5	Color-It fingertight adapter green 5/pkg	High Pressure PEEK Fittings	\$4.00	<a href="http://www.chromtech.net.au/pdf2/vici09_66.pdf">http://www.chromtech.net.au/pdf2/vici09_66.pdf</a>
JR-55011-5	Color-It fingertight adapter yellow 5/pkg	High Pressure PEEK Fittings	\$4.00	<a href="http://www.chromtech.net.au/pdf2/vici09_66.pdf">http://www.chromtech.net.au/pdf2/vici09_66.pdf</a>
JR-55024-5	Fitting PEEK one-piece blue 10-32 5/pkg	High Pressure PEEK Fittings	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_66.pdf">http://www.chromtech.net.au/pdf2/vici09_66.pdf</a>
JR-55023-5	Fitting PEEK one-piece yellow 10-32 5/pkg	High Pressure PEEK Fittings	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_66.pdf">http://www.chromtech.net.au/pdf2/vici09_66.pdf</a>
JR-55022-5	Fitting PEEK one-piece red 10-32 5/pkg	High Pressure PEEK Fittings	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_66.pdf">http://www.chromtech.net.au/pdf2/vici09_66.pdf</a>
JR-55021-5	Fitting PEEK one-piece black 10-32 5/pkg	High Pressure PEEK Fittings	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_66.pdf">http://www.chromtech.net.au/pdf2/vici09_66.pdf</a>
JR-55020-5	Fitting PEEK one-piece natural 10-32 5/pkg	High Pressure PEEK Fittings	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_66.pdf">http://www.chromtech.net.au/pdf2/vici09_66.pdf</a>
JR-55010-5	Color-It fingertight adapter blue 5/pkg	High Pressure PEEK Fittings	\$4.00	<a href="http://www.chromtech.net.au/pdf2/vici09_66.pdf">http://www.chromtech.net.au/pdf2/vici09_66.pdf</a>
JR-5508-5	Fitting PEEK 1/16"OD tubing one-piece fingertight hex- head 10-32 5/pkg	High Pressure PEEK Fittings	\$38.00	<a href="http://www.chromtech.net.au/pdf2/vici09_66.pdf">http://www.chromtech.net.au/pdf2/vici09_66.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

ZRU1.5TFPK	PEEK 1/16in to 1/32in 1/32in bore Standard	High Pressure PEEK Fittings	\$36.00	<a href="http://www.chromtech.net.au/pdf2/vici09_65.pdf">http://www.chromtech.net.au/pdf2/vici09_65.pdf</a>
ZRU2.5TFPK	1/8into 1/32in x1/32in -Standard	High Pressure PEEK Fittings	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_65.pdf">http://www.chromtech.net.au/pdf2/vici09_65.pdf</a>
ZBU2LPK	PEEK 1/8in fingertight x2.0mm bore -Bulkhead	High Pressure PEEK Fittings	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_65.pdf">http://www.chromtech.net.au/pdf2/vici09_65.pdf</a>
ZU2TPK	PEEK 1/8inhex x1/8in bore -standard	High Pressure PEEK Fittings	\$26.00	<a href="http://www.chromtech.net.au/pdf2/vici09_65.pdf">http://www.chromtech.net.au/pdf2/vici09_65.pdf</a>
ZBU2TPK	PEEK Internal Unions -1/8in hex 1/8in bulkhead	High Pressure PEEK Fittings	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_65.pdf">http://www.chromtech.net.au/pdf2/vici09_65.pdf</a>
ZRU1.5FPK	1/16in to 1/32in x0.25mm bore -Standard	High Pressure PEEK Fittings	\$36.00	<a href="http://www.chromtech.net.au/pdf2/vici09_65.pdf">http://www.chromtech.net.au/pdf2/vici09_65.pdf</a>
ZBRU1.5FPK	1/16in to 1/32in x0.25mm bore -Bulkhead	High Pressure PEEK Fittings	\$42.00	<a href="http://www.chromtech.net.au/pdf2/vici09_65.pdf">http://www.chromtech.net.au/pdf2/vici09_65.pdf</a>
ZBRU1.5LFPK	1/16in to 1/32in x0.50mm bore -Bulkhead	High Pressure PEEK Fittings	\$42.00	<a href="http://www.chromtech.net.au/pdf2/vici09_65.pdf">http://www.chromtech.net.au/pdf2/vici09_65.pdf</a>
ZBRU1.5TFPK	1/16in to 1/32in x1/32in -Bulkhead	High Pressure PEEK Fittings	\$42.00	<a href="http://www.chromtech.net.au/pdf2/vici09_65.pdf">http://www.chromtech.net.au/pdf2/vici09_65.pdf</a>
ZRU2.5FPK	1/8into 1/32in x0.25mm bore -Standard	High Pressure PEEK Fittings	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_65.pdf">http://www.chromtech.net.au/pdf2/vici09_65.pdf</a>
ZBRU2.5FPK	PEEK 1/8in to 1/32in 0.25mm bore Bulkhead	High Pressure PEEK Fittings	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_65.pdf">http://www.chromtech.net.au/pdf2/vici09_65.pdf</a>
ZU2LPK	PEEK Internal Unions -1/8in hex 2.0 mm standard	High Pressure PEEK Fittings	\$26.00	<a href="http://www.chromtech.net.au/pdf2/vici09_65.pdf">http://www.chromtech.net.au/pdf2/vici09_65.pdf</a>
ZBRU2.5LFPK	1/8into 1/32in x0.50mm bore -Bulkhead	High Pressure PEEK Fittings	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_65.pdf">http://www.chromtech.net.au/pdf2/vici09_65.pdf</a>
ZBRU21CFPK	1/8into 1/16in x0.25mm bore -Bulkhead	High Pressure PEEK Fittings	\$40.00	<a href="http://www.chromtech.net.au/pdf2/vici09_65.pdf">http://www.chromtech.net.au/pdf2/vici09_65.pdf</a>
ZBRU2.5TFPK	1/8into 1/32in x1/32in -Bulkhead	High Pressure PEEK Fittings	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_65.pdf">http://www.chromtech.net.au/pdf2/vici09_65.pdf</a>
ZRU21CFPK	1/8into 1/16in x0.25mm bore -Standard	High Pressure PEEK Fittings	\$33.00	<a href="http://www.chromtech.net.au/pdf2/vici09_65.pdf">http://www.chromtech.net.au/pdf2/vici09_65.pdf</a>
ZRU21FPK	PEEK 1/8in to 1/16in 0.75mm bore Bore Standard	High Pressure PEEK Fittings	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_65.pdf">http://www.chromtech.net.au/pdf2/vici09_65.pdf</a>
ZRU21LFPK	1/8into 1/16in x1.00mm bore -Standard	High Pressure PEEK Fittings	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_65.pdf">http://www.chromtech.net.au/pdf2/vici09_65.pdf</a>



## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

ZBRU21LFPK	PEEK 1/8in to 1/16in 1.00mm bore Bulkhead	High Pressure PEEK Fittings	\$36.00	<a href="http://www.chromtech.net.au/pdf2/vici09_65.pdf">http://www.chromtech.net.au/pdf2/vici09_65.pdf</a>
ZRU21TFPK	1/8into 1/16in x1/16in -Standard	High Pressure PEEK Fittings	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_65.pdf">http://www.chromtech.net.au/pdf2/vici09_65.pdf</a>
ZBRU21TFPK	1/8into 1/16in x1/16in -Bulkhead	High Pressure PEEK Fittings	\$36.00	<a href="http://www.chromtech.net.au/pdf2/vici09_65.pdf">http://www.chromtech.net.au/pdf2/vici09_65.pdf</a>
ZRU1.5LFPK	PEEK 1/16in to 1/32in 0.50mm bore Bore Standard	High Pressure PEEK Fittings	\$36.00	<a href="http://www.chromtech.net.au/pdf2/vici09_65.pdf">http://www.chromtech.net.au/pdf2/vici09_65.pdf</a>
ZU.5FPK	PEEK 1/32in fingertight x0.25mm bore -standard	High Pressure PEEK Fittings	\$33.00	<a href="http://www.chromtech.net.au/pdf2/vici09_65.pdf">http://www.chromtech.net.au/pdf2/vici09_65.pdf</a>
ZRU2.5LFPK	1/8into 1/32in x0.50mm bore -Standard	High Pressure PEEK Fittings	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_65.pdf">http://www.chromtech.net.au/pdf2/vici09_65.pdf</a>
ZBU1MFPK	PEEK 1/16in fingertight x0.50mm bore -Bulkhead	High Pressure PEEK Fittings	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_65.pdf">http://www.chromtech.net.au/pdf2/vici09_65.pdf</a>
ZBU2PK	PEEK 1/16in fingertight x0.75mm bore -Bulkhead	High Pressure PEEK Fittings	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_65.pdf">http://www.chromtech.net.au/pdf2/vici09_65.pdf</a>
ZBU.5LFPK	1/32in to 0.50mm-Std	High Pressure PEEK Fittings	\$46.00	<a href="http://www.chromtech.net.au/pdf2/vici09_65.pdf">http://www.chromtech.net.au/pdf2/vici09_65.pdf</a>
ZU.5TFPK	PEEK Internal Unions -1/32in fingertight 0.50mm standard	High Pressure PEEK Fittings	\$33.00	<a href="http://www.chromtech.net.au/pdf2/vici09_65.pdf">http://www.chromtech.net.au/pdf2/vici09_65.pdf</a>
ZBU.5TFPK	1/32into1/32in Std	High Pressure PEEK Fittings	\$46.00	<a href="http://www.chromtech.net.au/pdf2/vici09_65.pdf">http://www.chromtech.net.au/pdf2/vici09_65.pdf</a>
ZU1CFPK	PEEK Internal Unions -1/16in fingertight 0.25mm standard	High Pressure PEEK Fittings	\$27.00	<a href="http://www.chromtech.net.au/pdf2/vici09_65.pdf">http://www.chromtech.net.au/pdf2/vici09_65.pdf</a>
ZU.5LFPK	PEEK 1/32in fingertight x0.25mm bore -standard	High Pressure PEEK Fittings	\$33.00	<a href="http://www.chromtech.net.au/pdf2/vici09_65.pdf">http://www.chromtech.net.au/pdf2/vici09_65.pdf</a>
ZBU1CFPK	PEEK Internal Unions -1/16in fingertight 0.25mm bulkhead	High Pressure PEEK Fittings	\$33.00	<a href="http://www.chromtech.net.au/pdf2/vici09_65.pdf">http://www.chromtech.net.au/pdf2/vici09_65.pdf</a>
ZBRU21FPK	PEEK 1/8in to 1/16in 0.75mm bore Bulkhead	High Pressure PEEK Fittings	\$36.00	<a href="http://www.chromtech.net.au/pdf2/vici09_65.pdf">http://www.chromtech.net.au/pdf2/vici09_65.pdf</a>
ZU1MFPK	PEEK Internal Unions -1/16in fingertight 0.50mm standard	High Pressure PEEK Fittings	\$24.00	<a href="http://www.chromtech.net.au/pdf2/vici09_65.pdf">http://www.chromtech.net.au/pdf2/vici09_65.pdf</a>
ZU1FPK	PEEK Internal Unions -1/16in fingertight 0.75mm standard	High Pressure PEEK Fittings	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_65.pdf">http://www.chromtech.net.au/pdf2/vici09_65.pdf</a>
ZBU1FPK	PEEK Internal Unions -1/16in fingertight 0.75mm bulkhead	High Pressure PEEK Fittings	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_65.pdf">http://www.chromtech.net.au/pdf2/vici09_65.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

ZBU1PK	PEEK 1/16in fingertight x0.75mm bore -Bulkhead	High Pressure PEEK Fittings	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_65.pdf">http://www.chromtech.net.au/pdf2/vici09_65.pdf</a>
ZU2PK	PEEK Internal Unions -1/8in hex 0.75mm standard	High Pressure PEEK Fittings	\$26.00	<a href="http://www.chromtech.net.au/pdf2/vici09_65.pdf">http://www.chromtech.net.au/pdf2/vici09_65.pdf</a>
ZBU.5FPK	1/32in to 0.25mm-Std	High Pressure PEEK Fittings	\$46.00	<a href="http://www.chromtech.net.au/pdf2/vici09_65.pdf">http://www.chromtech.net.au/pdf2/vici09_65.pdf</a>
ZU1TPK	PEEK Internal Unions -1/16in hex 1/16in standard	High Pressure PEEK Fittings	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_65.pdf">http://www.chromtech.net.au/pdf2/vici09_65.pdf</a>
ZU1TFPK	PEEK Internal Unions -1/16in fingertight 1/16in standard	High Pressure PEEK Fittings	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_65.pdf">http://www.chromtech.net.au/pdf2/vici09_65.pdf</a>
ZU1PK	PEEK 1/16in fingertight x0.75mm bore -standard	High Pressure PEEK Fittings	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_65.pdf">http://www.chromtech.net.au/pdf2/vici09_65.pdf</a>
ZBU1MPK	PEEK 1/16in fingertight x0.75mm bore -Bulkhead	High Pressure PEEK Fittings	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_65.pdf">http://www.chromtech.net.au/pdf2/vici09_65.pdf</a>
ZU1MPK	PEEK Internal Unions -1/16in hex 0.50mm standard	High Pressure PEEK Fittings	\$24.00	<a href="http://www.chromtech.net.au/pdf2/vici09_65.pdf">http://www.chromtech.net.au/pdf2/vici09_65.pdf</a>
ZBU1CPK	PEEK Internal Unions -1/16in hex 0.25mm bulkhead	High Pressure PEEK Fittings	\$33.00	<a href="http://www.chromtech.net.au/pdf2/vici09_65.pdf">http://www.chromtech.net.au/pdf2/vici09_65.pdf</a>
ZU1CPK	PEEK Internal Unions -1/16in hex 0.25mm standard	High Pressure PEEK Fittings	\$27.00	<a href="http://www.chromtech.net.au/pdf2/vici09_65.pdf">http://www.chromtech.net.au/pdf2/vici09_65.pdf</a>
ZBU1TFPK	PEEK Internal Unions -1/16in fingertight 1/16in bulkhead	High Pressure PEEK Fittings	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_65.pdf">http://www.chromtech.net.au/pdf2/vici09_65.pdf</a>
ZBU1TPK	PEEK 1/16in fingertight x1/16in -bulkhead	High Pressure PEEK Fittings	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_65.pdf">http://www.chromtech.net.au/pdf2/vici09_65.pdf</a>
ZT2LPK	PEEK 1/8in 2.00mm bore Bore	High Pressure PEEK Fittings	\$46.00	<a href="http://www.chromtech.net.au/pdf2/vici09_64.pdf">http://www.chromtech.net.au/pdf2/vici09_64.pdf</a>
ZT.5LFPK	PEEK 1/32in 0.50mm bore Bore	High Pressure PEEK Fittings	\$54.00	<a href="http://www.chromtech.net.au/pdf2/vici09_64.pdf">http://www.chromtech.net.au/pdf2/vici09_64.pdf</a>
ZX.5LFPK	PEEK 1/32in 0.50mm bore Bore	High Pressure PEEK Fittings	\$66.00	<a href="http://www.chromtech.net.au/pdf2/vici09_64.pdf">http://www.chromtech.net.au/pdf2/vici09_64.pdf</a>
ZT1CFPK	PEEK 1/16in 0.25mm bore Bore	High Pressure PEEK Fittings	\$50.00	<a href="http://www.chromtech.net.au/pdf2/vici09_64.pdf">http://www.chromtech.net.au/pdf2/vici09_64.pdf</a>
ZX1CFPK	PEEK 1/16in 0.25mm bore Bore	High Pressure PEEK Fittings	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_64.pdf">http://www.chromtech.net.au/pdf2/vici09_64.pdf</a>
ZT1MFPK	PEEK 1/16in 0.50mm bore Bore	High Pressure PEEK Fittings	\$46.00	<a href="http://www.chromtech.net.au/pdf2/vici09_64.pdf">http://www.chromtech.net.au/pdf2/vici09_64.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

ZX1MFPK	PEEK 1/16in 0.50mm bore Bore	High Pressure PEEK Fittings	\$58.00	<a href="http://www.chromtech.net.au/pdf2/vici09_64.pdf">http://www.chromtech.net.au/pdf2/vici09_64.pdf</a>
ZT1FPK	PEEK 1/16in 0.725mm bore Bore	High Pressure PEEK Fittings	\$42.00	<a href="http://www.chromtech.net.au/pdf2/vici09_64.pdf">http://www.chromtech.net.au/pdf2/vici09_64.pdf</a>
ZX1FPK	PEEK 1/16in 0.75mm bore Bore	High Pressure PEEK Fittings	\$56.00	<a href="http://www.chromtech.net.au/pdf2/vici09_64.pdf">http://www.chromtech.net.au/pdf2/vici09_64.pdf</a>
ZT1LFPK	PEEK 1/16in 1.00mm bore Bore	High Pressure PEEK Fittings	\$42.00	<a href="http://www.chromtech.net.au/pdf2/vici09_64.pdf">http://www.chromtech.net.au/pdf2/vici09_64.pdf</a>
ZX1LFPK	PEEK 1/16in 1.00mm bore Bore	High Pressure PEEK Fittings	\$56.00	<a href="http://www.chromtech.net.au/pdf2/vici09_64.pdf">http://www.chromtech.net.au/pdf2/vici09_64.pdf</a>
ZX.5FPK	PEEK 1/32in 0.25mm bore Bore	High Pressure PEEK Fittings	\$66.00	<a href="http://www.chromtech.net.au/pdf2/vici09_64.pdf">http://www.chromtech.net.au/pdf2/vici09_64.pdf</a>
ZX2LPK	PEEK 1/8in 2.00mm bore Bore	High Pressure PEEK Fittings	\$58.00	<a href="http://www.chromtech.net.au/pdf2/vici09_64.pdf">http://www.chromtech.net.au/pdf2/vici09_64.pdf</a>
ZC1PK	PEEK 1/16inx.500in long hex	High Pressure PEEK Fittings	\$15.00	<a href="http://www.chromtech.net.au/pdf2/vici09_64.pdf">http://www.chromtech.net.au/pdf2/vici09_64.pdf</a>
ZT2PK	PEEK 1/8in 0.75mm bore Bore	High Pressure PEEK Fittings	\$46.00	<a href="http://www.chromtech.net.au/pdf2/vici09_64.pdf">http://www.chromtech.net.au/pdf2/vici09_64.pdf</a>
ZT.5FPK	PEEK 1/32in 0.25mm bore Bore	High Pressure PEEK Fittings	\$54.00	<a href="http://www.chromtech.net.au/pdf2/vici09_64.pdf">http://www.chromtech.net.au/pdf2/vici09_64.pdf</a>
C-ZP1FPK	1/16in FT PEEK Plugs	High Pressure PEEK Fittings	\$15.00	<a href="http://www.chromtech.net.au/pdf2/vici09_64.pdf">http://www.chromtech.net.au/pdf2/vici09_64.pdf</a>
C-LZP1PK	1/16in Long 0.880in Hex PEEK Plugs	High Pressure PEEK Fittings	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_64.pdf">http://www.chromtech.net.au/pdf2/vici09_64.pdf</a>
C-MZP1PK	1/16in Hex PEEK Plugs	High Pressure PEEK Fittings	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_64.pdf">http://www.chromtech.net.au/pdf2/vici09_64.pdf</a>
ZC2PK	PEEK 1/8in hex	High Pressure PEEK Fittings	\$16.00	<a href="http://www.chromtech.net.au/pdf2/vici09_64.pdf">http://www.chromtech.net.au/pdf2/vici09_64.pdf</a>
LZP1PK	PEEK 1/16inx.880in long hex	High Pressure PEEK Fittings	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_64.pdf">http://www.chromtech.net.au/pdf2/vici09_64.pdf</a>
MZP1PK	PEEK 1/16in x.500in long hex Plugs	High Pressure PEEK Fittings	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_64.pdf">http://www.chromtech.net.au/pdf2/vici09_64.pdf</a>
ZC1FPK	PEEK 1/16in fingertight	High Pressure PEEK Fittings	\$16.00	<a href="http://www.chromtech.net.au/pdf2/vici09_64.pdf">http://www.chromtech.net.au/pdf2/vici09_64.pdf</a>
ZP1FPK	PEEK 1/16in fingertight Plugs	High Pressure PEEK Fittings	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_64.pdf">http://www.chromtech.net.au/pdf2/vici09_64.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

LZP.5FPK	1/32in FT PEEK Plugs	High Pressure PEEK Fittings	\$14.00	<a href="http://www.chromtech.net.au/pdf2/vici09_64.pdf">http://www.chromtech.net.au/pdf2/vici09_64.pdf</a>
ZC.5FPK	PEEK 1/32in fingertight	High Pressure PEEK Fittings	\$16.00	<a href="http://www.chromtech.net.au/pdf2/vici09_64.pdf">http://www.chromtech.net.au/pdf2/vici09_64.pdf</a>
ZP.5FPK	PEEK 1/32in fingertight	High Pressure PEEK Fittings	\$14.00	<a href="http://www.chromtech.net.au/pdf2/vici09_64.pdf">http://www.chromtech.net.au/pdf2/vici09_64.pdf</a>
ZX2PK	PEEK 1/8in 0.75mm bore Bore	High Pressure PEEK Fittings	\$58.00	<a href="http://www.chromtech.net.au/pdf2/vici09_64.pdf">http://www.chromtech.net.au/pdf2/vici09_64.pdf</a>
ZP2PK	PEEK 1/8in hex Plugs	High Pressure PEEK Fittings	\$16.00	<a href="http://www.chromtech.net.au/pdf2/vici09_64.pdf">http://www.chromtech.net.au/pdf2/vici09_64.pdf</a>
ZN1FPK-10	PEEK 1/16in fingertight nut	High Pressure PEEK Fittings	\$42.00	<a href="http://www.chromtech.net.au/pdf2/vici09_63.pdf">http://www.chromtech.net.au/pdf2/vici09_63.pdf</a>
ZGF.5PK-10	1/32in Ferrules High Pressure PEEK-Grooved 10-pk	High Pressure PEEK Fittings	\$56.00	<a href="http://www.chromtech.net.au/pdf2/vici09_63.pdf">http://www.chromtech.net.au/pdf2/vici09_63.pdf</a>
LZN1PK-10	PEEK 1/16in x.880in long hex nut	High Pressure PEEK Fittings	\$42.00	<a href="http://www.chromtech.net.au/pdf2/vici09_63.pdf">http://www.chromtech.net.au/pdf2/vici09_63.pdf</a>
ZGF1PK-10	1/16in Ferrules High Pressure PEEK-Grooved 0-pk	High Pressure PEEK Fittings	\$50.00	<a href="http://www.chromtech.net.au/pdf2/vici09_63.pdf">http://www.chromtech.net.au/pdf2/vici09_63.pdf</a>
MZN1PK-10	PEEK 1/16in x.500in long hex nut	High Pressure PEEK Fittings	\$36.00	<a href="http://www.chromtech.net.au/pdf2/vici09_63.pdf">http://www.chromtech.net.au/pdf2/vici09_63.pdf</a>
ZN1PK-10	PEEK Internal Nuts-High Pressure -1/16in hex	High Pressure PEEK Fittings	\$36.00	<a href="http://www.chromtech.net.au/pdf2/vici09_63.pdf">http://www.chromtech.net.au/pdf2/vici09_63.pdf</a>
MZNF1PK-5	10-32 1/16 PEEK,Hex,Med PEEK Ferrule 5-pk	High Pressure PEEK Fittings	\$54.00	<a href="http://www.chromtech.net.au/pdf2/vici09_63.pdf">http://www.chromtech.net.au/pdf2/vici09_63.pdf</a>
ZN2PK-10	PEEK Internal Nuts-High Pressure -1/8in hex	High Pressure PEEK Fittings	\$42.00	<a href="http://www.chromtech.net.au/pdf2/vici09_63.pdf">http://www.chromtech.net.au/pdf2/vici09_63.pdf</a>
ZNF1PK-5	10-32 FT PEEK, PEEK Ferrule 5-pk	High Pressure PEEK Fittings	\$54.00	<a href="http://www.chromtech.net.au/pdf2/vici09_63.pdf">http://www.chromtech.net.au/pdf2/vici09_63.pdf</a>
LZN.5FPK-10	1/32in FT PEEK Internal HighP PEEK Ferrule	High Pressure PEEK Fittings	\$50.00	<a href="http://www.chromtech.net.au/pdf2/vici09_63.pdf">http://www.chromtech.net.au/pdf2/vici09_63.pdf</a>
MZNF1PKG-5	10-32 1/16 PEEK,Hex,Med PEEK/GF Ferrule 5-pk	High Pressure PEEK Fittings	\$56.00	<a href="http://www.chromtech.net.au/pdf2/vici09_63.pdf">http://www.chromtech.net.au/pdf2/vici09_63.pdf</a>
MZNF1KF-5	10-32 1/16 PEEK,Hex,Med CTFE/GF Ferrule 5-pk	High Pressure PEEK Fittings	\$54.00	<a href="http://www.chromtech.net.au/pdf2/vici09_63.pdf">http://www.chromtech.net.au/pdf2/vici09_63.pdf</a>
ZNF1PKG-5	10-32 FT PEEK PEEK/GF Ferrule 5-pk	High Pressure PEEK Fittings	\$57.00	<a href="http://www.chromtech.net.au/pdf2/vici09_63.pdf">http://www.chromtech.net.au/pdf2/vici09_63.pdf</a>

VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

LZNF1PK-5	10-32 1/16 PEEK,Hex,Long PEEK Ferrule 5-pk	High Pressure PEEK Fittings	\$56.00	<a href="http://www.chromtech.net.au/pdf2/vici09_63.pdf">http://www.chromtech.net.au/pdf2/vici09_63.pdf</a>
LZNF1KF-5	10-32 1/16 PEEK,Hex,Long CTFE Ferrule 5-pk	High Pressure PEEK Fittings	\$56.00	<a href="http://www.chromtech.net.au/pdf2/vici09_63.pdf">http://www.chromtech.net.au/pdf2/vici09_63.pdf</a>
ZNF1FPKG-5	10-32 1/16 PEEK,Hex,Short PEEK/GF Ferrule 5-pk	High Pressure PEEK Fittings	\$57.00	<a href="http://www.chromtech.net.au/pdf2/vici09_63.pdf">http://www.chromtech.net.au/pdf2/vici09_63.pdf</a>
ZNF1FPK-5	10-32 1/16 PEEK,Hex,Short PEEK Ferrule 5-pk	High Pressure PEEK Fittings	\$56.00	<a href="http://www.chromtech.net.au/pdf2/vici09_63.pdf">http://www.chromtech.net.au/pdf2/vici09_63.pdf</a>
ZNF1FKF-5	10-32 1/16 PEEK,Hex,Short CTFE Ferrule 5-pk	High Pressure PEEK Fittings	\$56.00	<a href="http://www.chromtech.net.au/pdf2/vici09_63.pdf">http://www.chromtech.net.au/pdf2/vici09_63.pdf</a>
ZN.5FPK-10	PEEK Internal Nuts-High Pressure -1/32in fingertight	High Pressure PEEK Fittings	\$42.00	<a href="http://www.chromtech.net.au/pdf2/vici09_63.pdf">http://www.chromtech.net.au/pdf2/vici09_63.pdf</a>
ZNF1KF-5	10-32 FT PEEK, CTFE Ferrule 5-pk	High Pressure PEEK Fittings	\$54.00	<a href="http://www.chromtech.net.au/pdf2/vici09_63.pdf">http://www.chromtech.net.au/pdf2/vici09_63.pdf</a>
C-NEF.5FPK.35S1	1/32in 100um 325-375um for FS inv liners	Nanovolume Col EndFittings for FS	\$81.00	<a href="http://www.chromtech.net.au/pdf2/vici09_62.pdf">http://www.chromtech.net.au/pdf2/vici09_62.pdf</a>
C-NLS1.25	FS Cap Replac Liners 27mm 225-275um	Nanovolume Col EndFittings for FS	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_62.pdf">http://www.chromtech.net.au/pdf2/vici09_62.pdf</a>
C-NLS1.35	FS Cap Replac Liners 27mm 325-375um	Nanovolume Col EndFittings for FS	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_62.pdf">http://www.chromtech.net.au/pdf2/vici09_62.pdf</a>
C-NLS1.30	FS Cap Replac Liners 27mm 275-325um	Nanovolume Col EndFittings for FS	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_62.pdf">http://www.chromtech.net.au/pdf2/vici09_62.pdf</a>
C-NEF.5FPK.20S1	1/32in 100um 175-225um for FS inv liners	Nanovolume Col EndFittings for FS	\$81.00	<a href="http://www.chromtech.net.au/pdf2/vici09_62.pdf">http://www.chromtech.net.au/pdf2/vici09_62.pdf</a>
C-NEF.5XFPK.15S1	1/32in 100um 125-175um for FS inv liners	Nanovolume Col EndFittings for FS	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_62.pdf">http://www.chromtech.net.au/pdf2/vici09_62.pdf</a>
C-NLS1.20	FS Cap Replac Liners 27mm 175-225um	Nanovolume Col EndFittings for FS	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_62.pdf">http://www.chromtech.net.au/pdf2/vici09_62.pdf</a>
C-NEF.5XFPK.20S1	1/32in 100um 175-225um for FS inv liners	Nanovolume Col EndFittings for FS	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_62.pdf">http://www.chromtech.net.au/pdf2/vici09_62.pdf</a>
C-NEF.5XFPK.25S1	1/32in 100um 225-275um for FS inv liners	Nanovolume Col EndFittings for FS	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_62.pdf">http://www.chromtech.net.au/pdf2/vici09_62.pdf</a>
C-NEF.5FPK.25S1	1/32in 100um 225-275um for FS inv liners	Nanovolume Col EndFittings for FS	\$81.00	<a href="http://www.chromtech.net.au/pdf2/vici09_62.pdf">http://www.chromtech.net.au/pdf2/vici09_62.pdf</a>
C-NEF.5XFPK.30S1	1/32in 100um 275-325um for FS inv liners	Nanovolume Col EndFittings for FS	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_62.pdf">http://www.chromtech.net.au/pdf2/vici09_62.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C-NEF.5FPK.30S1	1/32in 100um 275-325um for FS inv liners	Nanovolume Col EndFittings for FS	\$81.00	<a href="http://www.chromtech.net.au/pdf2/vici09_62.pdf">http://www.chromtech.net.au/pdf2/vici09_62.pdf</a>
C-NEF.5XFPK.35S1	1/32in 100um 325-375um for FS inv liners	Nanovolume Col EndFittings for FS	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_62.pdf">http://www.chromtech.net.au/pdf2/vici09_62.pdf</a>
C-NLS1.15	FS Cap Replac Liners 27mm 125-175um	Nanovolume Col EndFittings for FS	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_62.pdf">http://www.chromtech.net.au/pdf2/vici09_62.pdf</a>
C-NEF.5FPK.15S1	1/32in 100um 125-175um for FS inv liners	Nanovolume Col EndFittings for FS	\$81.00	<a href="http://www.chromtech.net.au/pdf2/vici09_62.pdf">http://www.chromtech.net.au/pdf2/vici09_62.pdf</a>
C-NL.30S-5	Liners for 1/32in Conn FS 10mm 325-375um	Nanovolume Unions/Liners for FS	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_61.pdf">http://www.chromtech.net.au/pdf2/vici09_61.pdf</a>
C-NL.35S-5	Liners for 1/32in Conn FS 10mm 325-375um	Nanovolume Unions/Liners for FS	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_61.pdf">http://www.chromtech.net.au/pdf2/vici09_61.pdf</a>
C-NL.35L-5	Liners for 1/32in Conn FS 27mm 325-375um	Nanovolume Unions/Liners for FS	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_61.pdf">http://www.chromtech.net.au/pdf2/vici09_61.pdf</a>
C-NL.20S-5	Liners for 1/32in Conn FS 10mm 175-225um	Nanovolume Unions/Liners for FS	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_61.pdf">http://www.chromtech.net.au/pdf2/vici09_61.pdf</a>
C-NL.15S-5	Liners for 1/32in Conn FS 10mm 125-175um	Nanovolume Unions/Liners for FS	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_61.pdf">http://www.chromtech.net.au/pdf2/vici09_61.pdf</a>
C-NEU.5FPK.25	Unions 100um FS 225-275um	Nanovolume Unions/Liners for FS	\$57.00	<a href="http://www.chromtech.net.au/pdf2/vici09_61.pdf">http://www.chromtech.net.au/pdf2/vici09_61.pdf</a>
C-NEU.5FPK.20	Unions 100um FS 175-225um	Nanovolume Unions/Liners for FS	\$57.00	<a href="http://www.chromtech.net.au/pdf2/vici09_61.pdf">http://www.chromtech.net.au/pdf2/vici09_61.pdf</a>
C-NL.15L-5	Liners for 1/32in Conn FS 27mm 125-175um	Nanovolume Unions/Liners for FS	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_61.pdf">http://www.chromtech.net.au/pdf2/vici09_61.pdf</a>
C-NEU.5FPK.35	Unions 100um FS 325-375um	Nanovolume Unions/Liners for FS	\$57.00	<a href="http://www.chromtech.net.au/pdf2/vici09_61.pdf">http://www.chromtech.net.au/pdf2/vici09_61.pdf</a>
C-NEU.5XFPK.35	Unions 150um FS 325-375um	Nanovolume Unions/Liners for FS	\$69.00	<a href="http://www.chromtech.net.au/pdf2/vici09_61.pdf">http://www.chromtech.net.au/pdf2/vici09_61.pdf</a>
C-NEU.5XFPK.20	Unions 150um FS 175-225um	Nanovolume Unions/Liners for FS	\$69.00	<a href="http://www.chromtech.net.au/pdf2/vici09_61.pdf">http://www.chromtech.net.au/pdf2/vici09_61.pdf</a>
C-NEU.5XFPK.30	Unions 150um FS 275-325um	Nanovolume Unions/Liners for FS	\$69.00	<a href="http://www.chromtech.net.au/pdf2/vici09_61.pdf">http://www.chromtech.net.au/pdf2/vici09_61.pdf</a>
C-NEU.5XFPK.25	Unions 150um FS 225-275um	Nanovolume Unions/Liners for FS	\$69.00	<a href="http://www.chromtech.net.au/pdf2/vici09_61.pdf">http://www.chromtech.net.au/pdf2/vici09_61.pdf</a>
C-NEU.5FPK.15	Unions 100um FS 125-175um	Nanovolume Unions/Liners for FS	\$57.00	<a href="http://www.chromtech.net.au/pdf2/vici09_61.pdf">http://www.chromtech.net.au/pdf2/vici09_61.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C-NEU.5XFPK.15	Unions 150um FS 125-175um	Nanovolume Unions/Liners for FS	\$69.00	<a href="http://www.chromtech.net.au/pdf2/vici09_61.pdf">http://www.chromtech.net.au/pdf2/vici09_61.pdf</a>
C-NL.25S-5	Liners for 1/32in Conn FS 10mm 225-275um	Nanovolume Unions/Liners for FS	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_61.pdf">http://www.chromtech.net.au/pdf2/vici09_61.pdf</a>
C-NL.30L-5	Liners for 1/32in Conn FS 27mm 275-325um	Nanovolume Unions/Liners for FS	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_61.pdf">http://www.chromtech.net.au/pdf2/vici09_61.pdf</a>
C-NL.20L-5	Liners for 1/32in Conn FS 27mm 175-225um	Nanovolume Unions/Liners for FS	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_61.pdf">http://www.chromtech.net.au/pdf2/vici09_61.pdf</a>
C-NL.25L-5	Liners for 1/32in Conn FS 27mm 225-275um	Nanovolume Unions/Liners for FS	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_61.pdf">http://www.chromtech.net.au/pdf2/vici09_61.pdf</a>
C-NEU.5FPK.30	Unions 100um FS 275-325um	Nanovolume Unions/Liners for FS	\$57.00	<a href="http://www.chromtech.net.au/pdf2/vici09_61.pdf">http://www.chromtech.net.au/pdf2/vici09_61.pdf</a>
.2FR.5-5	Nanovolumeume SS Frits 1/32in	1/32" Nanovolume Frits	\$22.00	<a href="http://www.chromtech.net.au/pdf2/vici09_60.pdf">http://www.chromtech.net.au/pdf2/vici09_60.pdf</a>
.5FR.5-5	Nanovolumeume SS Frits 1/32in	1/32" Nanovolume Frits	\$22.00	<a href="http://www.chromtech.net.au/pdf2/vici09_60.pdf">http://www.chromtech.net.au/pdf2/vici09_60.pdf</a>
C-NEU.5FPK	1/32in Union 100um for PEEK/Tit	1/32" Nanovolume Red Unions	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_60.pdf">http://www.chromtech.net.au/pdf2/vici09_60.pdf</a>
C-NYXFPK	1/32in Y for SS /FS tubing 100um	1/32" Nanovolume T,Y X's	\$134.00	<a href="http://www.chromtech.net.au/pdf2/vici09_60.pdf">http://www.chromtech.net.au/pdf2/vici09_60.pdf</a>
C-NVISF	1/32 Fill ports for CN2 Injectors	1/32" Nanovolume T,Y X's	\$16.00	<a href="http://www.chromtech.net.au/pdf2/vici09_60.pdf">http://www.chromtech.net.au/pdf2/vici09_60.pdf</a>
C-NXFPK	1/32in Cross SS /FS tubing 150um	1/32" Nanovolume T,Y X's	\$128.00	<a href="http://www.chromtech.net.au/pdf2/vici09_60.pdf">http://www.chromtech.net.au/pdf2/vici09_60.pdf</a>
C-NYFPK	1/32in Y for SS /FS tubing 150um	1/32" Nanovolume T,Y X's	\$106.00	<a href="http://www.chromtech.net.au/pdf2/vici09_60.pdf">http://www.chromtech.net.au/pdf2/vici09_60.pdf</a>
C-NXXFPK	1/32in Cross SS /FS tubing 100um	1/32" Nanovolume T,Y X's	\$134.00	<a href="http://www.chromtech.net.au/pdf2/vici09_60.pdf">http://www.chromtech.net.au/pdf2/vici09_60.pdf</a>
C-NTXFPK	1/32in Tee SS /FS tubing 100um	1/32" Nanovolume T,Y X's	\$134.00	<a href="http://www.chromtech.net.au/pdf2/vici09_60.pdf">http://www.chromtech.net.au/pdf2/vici09_60.pdf</a>

VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

ZERU1.5FPK	PEEK Reducing Union for 1/32"	1/32" Nanovolume T,Y X's	\$33.00	<a href="http://www.chromtech.net.au/pdf2/vici09_60.pdf">http://www.chromtech.net.au/pdf2/vici09_60.pdf</a>
C-NTFPK	1/32in Union 150um for PEEK/Tit	1/32" Nanovolume T,Y X's	\$106.00	<a href="http://www.chromtech.net.au/pdf2/vici09_60.pdf">http://www.chromtech.net.au/pdf2/vici09_60.pdf</a>
C-EN.5FPKB	360um External Nut	1/32" Nanovolume Fittings	\$6.00	<a href="http://www.chromtech.net.au/pdf2/vici09_59.pdf">http://www.chromtech.net.au/pdf2/vici09_59.pdf</a>
ZGF.5PK	1/32 Collaps Ferrule	1/32" Nanovolume Fittings	\$6.00	<a href="http://www.chromtech.net.au/pdf2/vici09_59.pdf">http://www.chromtech.net.au/pdf2/vici09_59.pdf</a>
C-NNFLFPK	1/32in Int Nut w/Collaps Ferrule,10p valve	1/32" Nanovolume Fittings	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_59.pdf">http://www.chromtech.net.au/pdf2/vici09_59.pdf</a>
C-NNFFPK	1/32in Int Nut w/Collaps Ferrule,4,6p valve	1/32" Nanovolume Fittings	\$10.00	<a href="http://www.chromtech.net.au/pdf2/vici09_59.pdf">http://www.chromtech.net.au/pdf2/vici09_59.pdf</a>
C-NPFPK	1/32in Internal Plug	1/32" Nanovolume Fittings	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_59.pdf">http://www.chromtech.net.au/pdf2/vici09_59.pdf</a>
C360RU.5S62	360um Red Union 1/32 to 360um 50um	CI 360um Nanovolume Fittings	\$62.00	<a href="http://www.chromtech.net.au/pdf2/vici09_58.pdf">http://www.chromtech.net.au/pdf2/vici09_58.pdf</a>
C360RU.5FS4	360um Red Union 1/32 to 360um 100um FS	CI 360um Nanovolume Fittings	\$87.00	<a href="http://www.chromtech.net.au/pdf2/vici09_58.pdf">http://www.chromtech.net.au/pdf2/vici09_58.pdf</a>
C360US62	360um Union 50um	CI 360um Nanovolume Fittings	\$76.00	<a href="http://www.chromtech.net.au/pdf2/vici09_58.pdf">http://www.chromtech.net.au/pdf2/vici09_58.pdf</a>
C360NFS6	360um Nut/Ferrule	CI 360um Nanovolume Fittings	\$15.00	<a href="http://www.chromtech.net.au/pdf2/vici09_58.pdf">http://www.chromtech.net.au/pdf2/vici09_58.pdf</a>
C360C	SS 360um Cap	CI 360um Nanovolume Fittings	\$26.00	<a href="http://www.chromtech.net.au/pdf2/vici09_58.pdf">http://www.chromtech.net.au/pdf2/vici09_58.pdf</a>
C360UFS4	360um Union 100um FS	CI 360um Nanovolume Fittings	\$102.00	<a href="http://www.chromtech.net.au/pdf2/vici09_58.pdf">http://www.chromtech.net.au/pdf2/vici09_58.pdf</a>
C360UFS2	360um Union 50um FS	CI 360um Nanovolume Fittings	\$106.00	<a href="http://www.chromtech.net.au/pdf2/vici09_58.pdf">http://www.chromtech.net.au/pdf2/vici09_58.pdf</a>
C360RU.5FS6	360um Red Union 1/32 to 360um 1 50um FS	CI 360um Nanovolume Fittings	\$81.00	<a href="http://www.chromtech.net.au/pdf2/vici09_58.pdf">http://www.chromtech.net.au/pdf2/vici09_58.pdf</a>
C360CFS	SS 360um Cap FS Tubing	CI 360um Nanovolume Fittings	\$40.00	<a href="http://www.chromtech.net.au/pdf2/vici09_58.pdf">http://www.chromtech.net.au/pdf2/vici09_58.pdf</a>
C360US64	360um Union 100um	CI 360um Nanovolume Fittings	\$72.00	<a href="http://www.chromtech.net.au/pdf2/vici09_58.pdf">http://www.chromtech.net.au/pdf2/vici09_58.pdf</a>



## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C360NFFS	SS 360um Nut/Ferrule FS Tubing	CI 360um Nanovolume Fittings	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_58.pdf">http://www.chromtech.net.au/pdf2/vici09_58.pdf</a>
C360RU1FS6	360um Red Union 1/32 to 360um 150um FS	CI 360um Nanovolume Fittings	\$76.00	<a href="http://www.chromtech.net.au/pdf2/vici09_58.pdf">http://www.chromtech.net.au/pdf2/vici09_58.pdf</a>
C360IZR.5S6PK	Internal Red for360um Tubing 1/32 to 360um PEEK	CI 360um Nanovolume Fittings	\$44.00	<a href="http://www.chromtech.net.au/pdf2/vici09_58.pdf">http://www.chromtech.net.au/pdf2/vici09_58.pdf</a>
C360IZR.5S6	Internal Red for360um Tubing 1/32 to 360um SS	CI 360um Nanovolume Fittings	\$44.00	<a href="http://www.chromtech.net.au/pdf2/vici09_58.pdf">http://www.chromtech.net.au/pdf2/vici09_58.pdf</a>
C360RU1S66	360um Red Union 1/32 to 360um 150um	CI 360um Nanovolume Fittings	\$46.00	<a href="http://www.chromtech.net.au/pdf2/vici09_58.pdf">http://www.chromtech.net.au/pdf2/vici09_58.pdf</a>
C360RU.5S66	360um Red Union 1/32 to 360um 1 50um	CI 360um Nanovolume Fittings	\$51.00	<a href="http://www.chromtech.net.au/pdf2/vici09_58.pdf">http://www.chromtech.net.au/pdf2/vici09_58.pdf</a>
C360RU.5S64	360um Red Union 1/32 to 360um 100um	CI 360um Nanovolume Fittings	\$56.00	<a href="http://www.chromtech.net.au/pdf2/vici09_58.pdf">http://www.chromtech.net.au/pdf2/vici09_58.pdf</a>
C360UFS6	360um Union 150um FS	CI 360um Nanovolume Fittings	\$98.00	<a href="http://www.chromtech.net.au/pdf2/vici09_58.pdf">http://www.chromtech.net.au/pdf2/vici09_58.pdf</a>
C360US66	360um Union 150um	CI 360um Nanovolume Fittings	\$66.00	<a href="http://www.chromtech.net.au/pdf2/vici09_58.pdf">http://www.chromtech.net.au/pdf2/vici09_58.pdf</a>
C360RU.5FS2	360um Red Union 1/32 to 360um 50um FS	CI 360um Nanovolume Fittings	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_58.pdf">http://www.chromtech.net.au/pdf2/vici09_58.pdf</a>
C360QUPKG2	360um Union, Quick Mount 50um	CI 360um Nanovolume Fittings	\$90.00	<a href="http://www.chromtech.net.au/pdf2/vici09_57.pdf">http://www.chromtech.net.au/pdf2/vici09_57.pdf</a>
C360UPKG4	360um Union 100um	CI 360um Nanovolume Fittings	\$62.00	<a href="http://www.chromtech.net.au/pdf2/vici09_57.pdf">http://www.chromtech.net.au/pdf2/vici09_57.pdf</a>
C360NFPKG	360um Nut/Ferrule	CI 360um Nanovolume Fittings	\$15.00	<a href="http://www.chromtech.net.au/pdf2/vici09_57.pdf">http://www.chromtech.net.au/pdf2/vici09_57.pdf</a>
C360CPKG	360um Cap	CI 360um Nanovolume Fittings	\$15.00	<a href="http://www.chromtech.net.au/pdf2/vici09_57.pdf">http://www.chromtech.net.au/pdf2/vici09_57.pdf</a>
C360PPK	360um Plug	CI 360um Nanovolume Fittings	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_57.pdf">http://www.chromtech.net.au/pdf2/vici09_57.pdf</a>
C360ET	360um Tightening Tool	CI 360um Nanovolume Fittings	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_57.pdf">http://www.chromtech.net.au/pdf2/vici09_57.pdf</a>
C360QTPKG4	360um Tee, Quick Mount 100um	CI 360um Nanovolume Fittings	\$105.00	<a href="http://www.chromtech.net.au/pdf2/vici09_57.pdf">http://www.chromtech.net.au/pdf2/vici09_57.pdf</a>
C360QXPKG4	360um Cross, Quick Mount 100um	CI 360um Nanovolume Fittings	\$128.00	<a href="http://www.chromtech.net.au/pdf2/vici09_57.pdf">http://www.chromtech.net.au/pdf2/vici09_57.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C360QXPKG6	360um Cross, Quick Mount 150um	CI 360um Nanovolume Fittings	\$120.00	<a href="http://www.chromtech.net.au/pdf2/vici09_57.pdf">http://www.chromtech.net.au/pdf2/vici09_57.pdf</a>
C360QUPKG6	360um Union, Quick Mount 150um	CI 360um Nanovolume Fittings	\$76.00	<a href="http://www.chromtech.net.au/pdf2/vici09_57.pdf">http://www.chromtech.net.au/pdf2/vici09_57.pdf</a>
C360QTPKG6	360um Tee, Quick Mount 150um	CI 360um Nanovolume Fittings	\$98.00	<a href="http://www.chromtech.net.au/pdf2/vici09_57.pdf">http://www.chromtech.net.au/pdf2/vici09_57.pdf</a>
C360QUPKG4	360um Union, Quick Mount 100um	CI 360um Nanovolume Fittings	\$82.00	<a href="http://www.chromtech.net.au/pdf2/vici09_57.pdf">http://www.chromtech.net.au/pdf2/vici09_57.pdf</a>
C360QTPKG2	360um Tee, Quick Mount 50um	CI 360um Nanovolume Fittings	\$112.00	<a href="http://www.chromtech.net.au/pdf2/vici09_57.pdf">http://www.chromtech.net.au/pdf2/vici09_57.pdf</a>
C360RU1PK6	360um Reducing Union 1/16 to 360um	CI 360um Nanovolume Fittings	\$40.00	<a href="http://www.chromtech.net.au/pdf2/vici09_57.pdf">http://www.chromtech.net.au/pdf2/vici09_57.pdf</a>
C360UPKG6	360um Union 150um	CI 360um Nanovolume Fittings	\$56.00	<a href="http://www.chromtech.net.au/pdf2/vici09_57.pdf">http://www.chromtech.net.au/pdf2/vici09_57.pdf</a>
C360UPKG2	360um Union 50um	CI 360um Nanovolume Fittings	\$66.00	<a href="http://www.chromtech.net.au/pdf2/vici09_57.pdf">http://www.chromtech.net.au/pdf2/vici09_57.pdf</a>
C360QXPKG2	360um Cross, Quick Mount 50um	CI 360um Nanovolume Fittings	\$135.00	<a href="http://www.chromtech.net.au/pdf2/vici09_57.pdf">http://www.chromtech.net.au/pdf2/vici09_57.pdf</a>
TEMPLATE1	Template	Tools	\$4.00	<a href="http://www.chromtech.net.au/pdf2/vici09_55.pdf">http://www.chromtech.net.au/pdf2/vici09_55.pdf</a>
PV	PinVise Drill Index -.0135in to .039in	Tools	\$38.00	<a href="http://www.chromtech.net.au/pdf2/vici09_55.pdf">http://www.chromtech.net.au/pdf2/vici09_55.pdf</a>
OEW-2	Open end Wrench -3/8inx7/16in for 1/8in nuts	Tools	\$15.00	<a href="http://www.chromtech.net.au/pdf2/vici09_55.pdf">http://www.chromtech.net.au/pdf2/vici09_55.pdf</a>
OEW	Open end Wrench -3/16inx1/4in for 1/32in and 1/16in nuts	Tools	\$8.00	<a href="http://www.chromtech.net.au/pdf2/vici09_55.pdf">http://www.chromtech.net.au/pdf2/vici09_55.pdf</a>
OEW-3	Open end Wrench -1/2inx9/16in for 1/4in nuts	Tools	\$15.00	<a href="http://www.chromtech.net.au/pdf2/vici09_55.pdf">http://www.chromtech.net.au/pdf2/vici09_55.pdf</a>
FRK1	Ferrule removal kit	Frits/Screens for Filters	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_54.pdf">http://www.chromtech.net.au/pdf2/vici09_54.pdf</a>
SWH4	Custom Socket Wrench -1/4in slotted	Tools	\$9.00	<a href="http://www.chromtech.net.au/pdf2/vici09_54.pdf">http://www.chromtech.net.au/pdf2/vici09_54.pdf</a>
HKS	Hex Key Set -0.05in to 5/32in x8	Tools	\$15.00	<a href="http://www.chromtech.net.au/pdf2/vici09_54.pdf">http://www.chromtech.net.au/pdf2/vici09_54.pdf</a>
2FR4-10	Stainless steel 1/4in frits 2u x1.00mm	Frits/Screens for Filters	\$15.00	<a href="http://www.chromtech.net.au/pdf2/vici09_53.pdf">http://www.chromtech.net.au/pdf2/vici09_53.pdf</a>
.5FR1HC-10	Hastelloy C 1/16in frits 0.5u x0.75mmID	Frits/Screens for Filters	\$27.00	<a href="http://www.chromtech.net.au/pdf2/vici09_53.pdf">http://www.chromtech.net.au/pdf2/vici09_53.pdf</a>
2FR1-10	Stainless steel 1/16in frits 2u x0.75mm	Frits/Screens for Filters	\$15.00	<a href="http://www.chromtech.net.au/pdf2/vici09_53.pdf">http://www.chromtech.net.au/pdf2/vici09_53.pdf</a>
2FR1HC-10	Hastelloy C 1/16in frits 2u x0.75mm	Frits/Screens for Filters	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_53.pdf">http://www.chromtech.net.au/pdf2/vici09_53.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

2FR1TI-10	Titanium 1/16in frits 2u x0.75mm	Frits/Screens for Filters	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_53.pdf">http://www.chromtech.net.au/pdf2/vici09_53.pdf</a>
10FR1-10	Stainless steel 1/16in frits 10u x0.75mmID	Frits/Screens for Filters	\$15.00	<a href="http://www.chromtech.net.au/pdf2/vici09_53.pdf">http://www.chromtech.net.au/pdf2/vici09_53.pdf</a>
.5FR2-10	Stainless steel 1/8in frits 0.5u x1.00mm	Frits/Screens for Filters	\$15.00	<a href="http://www.chromtech.net.au/pdf2/vici09_53.pdf">http://www.chromtech.net.au/pdf2/vici09_53.pdf</a>
.5FR2HC-10	Replac Frits Hastelloy C	Frits/Screens for Filters	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_53.pdf">http://www.chromtech.net.au/pdf2/vici09_53.pdf</a>
1FR2-10	Stainless steel 1/16in frits 0.5u x1.00mm	Frits/Screens for Filters	\$15.00	<a href="http://www.chromtech.net.au/pdf2/vici09_53.pdf">http://www.chromtech.net.au/pdf2/vici09_53.pdf</a>
1FR2HC-10	Has Cl 1/16in frits 0.5u x1.00mm	Frits/Screens for Filters	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_53.pdf">http://www.chromtech.net.au/pdf2/vici09_53.pdf</a>
2FR2-10	Stainless steel 1/8in frits 2u x1.00mm	Frits/Screens for Filters	\$15.00	<a href="http://www.chromtech.net.au/pdf2/vici09_53.pdf">http://www.chromtech.net.au/pdf2/vici09_53.pdf</a>
2FR2HC-10	Hastelloy C 1/8in frits 2u x1.00mm	Frits/Screens for Filters	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_53.pdf">http://www.chromtech.net.au/pdf2/vici09_53.pdf</a>
2FR2TI-10	Titanium 1/8in frits 2u x1.00mm	Frits/Screens for Filters	\$57.00	<a href="http://www.chromtech.net.au/pdf2/vici09_53.pdf">http://www.chromtech.net.au/pdf2/vici09_53.pdf</a>
.5FR1-10	Stainless steel 1/16in frits 0.5u x0.75mmID	Frits/Screens for Filters	\$15.00	<a href="http://www.chromtech.net.au/pdf2/vici09_53.pdf">http://www.chromtech.net.au/pdf2/vici09_53.pdf</a>
.5FR4-10	Stainless steel 1/4in frits 0.5u x1.00mmID	Frits/Screens for Filters	\$15.00	<a href="http://www.chromtech.net.au/pdf2/vici09_53.pdf">http://www.chromtech.net.au/pdf2/vici09_53.pdf</a>
10FR4HC-10	Hastelloy C 1/4in frits 10u x1.00mmID	Frits/Screens for Filters	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_53.pdf">http://www.chromtech.net.au/pdf2/vici09_53.pdf</a>
2FR4HC-10	Hastelloy C 1/4in frits 2u x1.00mmID	Frits/Screens for Filters	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_53.pdf">http://www.chromtech.net.au/pdf2/vici09_53.pdf</a>
2FR4TI-10	Titanium 1/4in frits 2u x1.00mm	Frits/Screens for Filters	\$76.00	<a href="http://www.chromtech.net.au/pdf2/vici09_53.pdf">http://www.chromtech.net.au/pdf2/vici09_53.pdf</a>
10FR4-10	Stainless steel 1/4in frits 10u x1.00mmID	Frits/Screens for Filters	\$15.00	<a href="http://www.chromtech.net.au/pdf2/vici09_53.pdf">http://www.chromtech.net.au/pdf2/vici09_53.pdf</a>
.5SR.5-10	Stainless steel 1/32in screens 0.5u x1.00mm	Frits/Screens for Filters	\$26.00	<a href="http://www.chromtech.net.au/pdf2/vici09_53.pdf">http://www.chromtech.net.au/pdf2/vici09_53.pdf</a>
2SR.5-10	Stainless steel 1/32in screens 2u x0.075mm	Frits/Screens for Filters	\$15.00	<a href="http://www.chromtech.net.au/pdf2/vici09_53.pdf">http://www.chromtech.net.au/pdf2/vici09_53.pdf</a>
10SR.5-10	Stainless steel 1/32in screens 10u x0.125mm	Frits/Screens for Filters	\$15.00	<a href="http://www.chromtech.net.au/pdf2/vici09_53.pdf">http://www.chromtech.net.au/pdf2/vici09_53.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

.5SR1-10	Stainless steel 1/16in screens 0.5u x1.00mm	Frits/Screens for Filters	\$26.00	<a href="http://www.chromtech.net.au/pdf2/vici09_53.pdf">http://www.chromtech.net.au/pdf2/vici09_53.pdf</a>
1SR1-10	Stainless steel 1/32in screens 1u x0.05mm	Frits/Screens for Filters	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_53.pdf">http://www.chromtech.net.au/pdf2/vici09_53.pdf</a>
2SR1-10	Stainless steel 1/16in screens 2u x0.075mm	Frits/Screens for Filters	\$15.00	<a href="http://www.chromtech.net.au/pdf2/vici09_53.pdf">http://www.chromtech.net.au/pdf2/vici09_53.pdf</a>
10SR1-10	Stainless steel 1/16in screens 10u x0.125mm	Frits/Screens for Filters	\$15.00	<a href="http://www.chromtech.net.au/pdf2/vici09_53.pdf">http://www.chromtech.net.au/pdf2/vici09_53.pdf</a>
.5SR2-10	Stainless steel 1/8in screens 0.5um x1.00mm	Frits/Screens for Filters	\$26.00	<a href="http://www.chromtech.net.au/pdf2/vici09_53.pdf">http://www.chromtech.net.au/pdf2/vici09_53.pdf</a>
1SR2-10	Stainless steel 1/8in screens 1u x0.05mm	Frits/Screens for Filters	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_53.pdf">http://www.chromtech.net.au/pdf2/vici09_53.pdf</a>
10FR2-10	Stainless steel 1/8in frits 10u x1.00mmID	Frits/Screens for Filters	\$15.00	<a href="http://www.chromtech.net.au/pdf2/vici09_53.pdf">http://www.chromtech.net.au/pdf2/vici09_53.pdf</a>
1SR.5-10	Stainless steel 1/32in screens 1u x0.05mm	Frits/Screens for Filters	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_53.pdf">http://www.chromtech.net.au/pdf2/vici09_53.pdf</a>
2SR2-10	Stainless steel 1/8in screens 2u x0.075mm	Frits/Screens for Filters	\$15.00	<a href="http://www.chromtech.net.au/pdf2/vici09_53.pdf">http://www.chromtech.net.au/pdf2/vici09_53.pdf</a>
10SR2-10	Stainless steel 1/8in screens 10u x0.125mm	Frits/Screens for Filters	\$15.00	<a href="http://www.chromtech.net.au/pdf2/vici09_53.pdf">http://www.chromtech.net.au/pdf2/vici09_53.pdf</a>
.5SR4-10	Stainless steel 1/8in screens 0.5um x1.00mm	Frits/Screens for Filters	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_53.pdf">http://www.chromtech.net.au/pdf2/vici09_53.pdf</a>
10SR4-10	Stainless steel 1/4in screens 10u x0.125mm	Frits/Screens for Filters	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_53.pdf">http://www.chromtech.net.au/pdf2/vici09_53.pdf</a>
2SR4-10	Stainless steel 1/4in screens 2u x0.075mm	Frits/Screens for Filters	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_53.pdf">http://www.chromtech.net.au/pdf2/vici09_53.pdf</a>
1SR4-10	Stainless steel 1/4in screens 1u x0.05mm	Frits/Screens for Filters	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_53.pdf">http://www.chromtech.net.au/pdf2/vici09_53.pdf</a>
ZBRUFR21	1/8in to 1/16in -0.75mm bore	Filters	\$46.00	<a href="http://www.chromtech.net.au/pdf2/vici09_52.pdf">http://www.chromtech.net.au/pdf2/vici09_52.pdf</a>
ZRUFR42	1/4in to 1/8inx2.00mm bore	Filters	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_52.pdf">http://www.chromtech.net.au/pdf2/vici09_52.pdf</a>
ZBRUFR42	1/4in to 1/8inx2.00mm bore	Filters	\$51.00	<a href="http://www.chromtech.net.au/pdf2/vici09_52.pdf">http://www.chromtech.net.au/pdf2/vici09_52.pdf</a>
ZBRUFR41	1/4in to 1/16inx1.00mm bore	Filters	\$51.00	<a href="http://www.chromtech.net.au/pdf2/vici09_52.pdf">http://www.chromtech.net.au/pdf2/vici09_52.pdf</a>
ZRUFR41	1/4in to 1/16inx1.00mm bore	Filters	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_52.pdf">http://www.chromtech.net.au/pdf2/vici09_52.pdf</a>
ZBUFR2	1/8in to 1/8inx2.00mm bore	Filters	\$46.00	<a href="http://www.chromtech.net.au/pdf2/vici09_52.pdf">http://www.chromtech.net.au/pdf2/vici09_52.pdf</a>
ZUFR2	1/8in to 1/8inx2.00mm bore	Filters	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_52.pdf">http://www.chromtech.net.au/pdf2/vici09_52.pdf</a>
ZRUFR21	1/8in to 1/16in -0.75mm bore	Filters	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_52.pdf">http://www.chromtech.net.au/pdf2/vici09_52.pdf</a>
ZBUFR1	1/16in to 1/16in -0.50mm bore	Filters	\$46.00	<a href="http://www.chromtech.net.au/pdf2/vici09_52.pdf">http://www.chromtech.net.au/pdf2/vici09_52.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

ZUFR1	1/16in to 1/16in -0.50mm bore	Filters	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_52.pdf">http://www.chromtech.net.au/pdf2/vici09_52.pdf</a>
ZUFR1C	1/16in to 1/16in -0.25mm bore	Filters	\$51.00	<a href="http://www.chromtech.net.au/pdf2/vici09_52.pdf">http://www.chromtech.net.au/pdf2/vici09_52.pdf</a>
ZRUFR1.5	1/16in to 1/32in -0.25mm bore	Filters	\$46.00	<a href="http://www.chromtech.net.au/pdf2/vici09_52.pdf">http://www.chromtech.net.au/pdf2/vici09_52.pdf</a>
ZUFR.5	1/32in to 1/32in -0.25mm bore	Filters	\$66.00	<a href="http://www.chromtech.net.au/pdf2/vici09_52.pdf">http://www.chromtech.net.au/pdf2/vici09_52.pdf</a>
ZBRUFR1.5	1/16in to 1/32in -0.25mm bore	Filters	\$51.00	<a href="http://www.chromtech.net.au/pdf2/vici09_52.pdf">http://www.chromtech.net.au/pdf2/vici09_52.pdf</a>
ZBUFR.5	1/32in to 1/32in -0.25mm bore	Filters	\$72.00	<a href="http://www.chromtech.net.au/pdf2/vici09_52.pdf">http://www.chromtech.net.au/pdf2/vici09_52.pdf</a>
ZBUFR1C	1/16in to 1/16in -0.25mm bore	Filters	\$57.00	<a href="http://www.chromtech.net.au/pdf2/vici09_52.pdf">http://www.chromtech.net.au/pdf2/vici09_52.pdf</a>
ZRUFR1.5F	1/16in to 1/32in -0.25mm bore	Filters	\$46.00	<a href="http://www.chromtech.net.au/pdf2/vici09_51.pdf">http://www.chromtech.net.au/pdf2/vici09_51.pdf</a>
ZBUFR2F	1/8in to 1/8in -2.00mm bore B/H	Filters	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_51.pdf">http://www.chromtech.net.au/pdf2/vici09_51.pdf</a>
ZUFR.5F	1/32in to 1/32in -0.25mm bore	Filters	\$66.00	<a href="http://www.chromtech.net.au/pdf2/vici09_51.pdf">http://www.chromtech.net.au/pdf2/vici09_51.pdf</a>
ZBUFR.5F	1/32in to 1/32in -0.25mm bore B/H	Filters	\$72.00	<a href="http://www.chromtech.net.au/pdf2/vici09_51.pdf">http://www.chromtech.net.au/pdf2/vici09_51.pdf</a>
ZRUFR42F	1/4in to 1/8in -2.00mm bore	Filters	\$42.00	<a href="http://www.chromtech.net.au/pdf2/vici09_51.pdf">http://www.chromtech.net.au/pdf2/vici09_51.pdf</a>
ZRUFR41F	1/4in to 1/16in -1.00mm bore	Filters	\$42.00	<a href="http://www.chromtech.net.au/pdf2/vici09_51.pdf">http://www.chromtech.net.au/pdf2/vici09_51.pdf</a>
ZBRUFR42F	1/4in to 1/8in -2.00mm bore B/H	Filters	\$48.00	<a href="http://www.chromtech.net.au/pdf2/vici09_51.pdf">http://www.chromtech.net.au/pdf2/vici09_51.pdf</a>
ZUFR2F	1/8in to 1/8in -2.00mm bore	Filters	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_51.pdf">http://www.chromtech.net.au/pdf2/vici09_51.pdf</a>
ZBRUFR21F	1/8in to 1/16in -0.75mm bore B/H	Filters	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_51.pdf">http://www.chromtech.net.au/pdf2/vici09_51.pdf</a>
ZRUFR21F	1/8in to 1/16in -0.75mm bore	Filters	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_51.pdf">http://www.chromtech.net.au/pdf2/vici09_51.pdf</a>
ZBUFR1F	1/16in to 1/16in -0.50mm bore B/H	Filters	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_51.pdf">http://www.chromtech.net.au/pdf2/vici09_51.pdf</a>
ZUFR1F	1/16in to 1/16in -0.50mm bore	Filters	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_51.pdf">http://www.chromtech.net.au/pdf2/vici09_51.pdf</a>
ZBUFR1CF	1/16in to 1/16in -0.25mm bore B/H	Filters	\$57.00	<a href="http://www.chromtech.net.au/pdf2/vici09_51.pdf">http://www.chromtech.net.au/pdf2/vici09_51.pdf</a>
ZUFR1CF	1/16in to 1/16in -0.25mm bore	Filters	\$51.00	<a href="http://www.chromtech.net.au/pdf2/vici09_51.pdf">http://www.chromtech.net.au/pdf2/vici09_51.pdf</a>
ZBRUFR1.5F	1/16in to 1/32in -0.25mm bore B/H	Filters	\$54.00	<a href="http://www.chromtech.net.au/pdf2/vici09_51.pdf">http://www.chromtech.net.au/pdf2/vici09_51.pdf</a>
ZBRUFR41F	1/4in to 1/16in -1.00mm bore B/H	Filters	\$48.00	<a href="http://www.chromtech.net.au/pdf2/vici09_51.pdf">http://www.chromtech.net.au/pdf2/vici09_51.pdf</a>
ZRUF21	1/8in to 1/16in -0.75mm bore	Filters	\$27.00	<a href="http://www.chromtech.net.au/pdf2/vici09_50.pdf">http://www.chromtech.net.au/pdf2/vici09_50.pdf</a>
ZUF4	1/4in to 1/4in -4.60mm bore	Filters	\$27.00	<a href="http://www.chromtech.net.au/pdf2/vici09_50.pdf">http://www.chromtech.net.au/pdf2/vici09_50.pdf</a>
ZBRUF42	1/4in to 1/8in -2.00mm bore	Filters	\$33.00	<a href="http://www.chromtech.net.au/pdf2/vici09_50.pdf">http://www.chromtech.net.au/pdf2/vici09_50.pdf</a>
ZRUF42	1/4in to 1/8in -2.00mm bore	Filters	\$27.00	<a href="http://www.chromtech.net.au/pdf2/vici09_50.pdf">http://www.chromtech.net.au/pdf2/vici09_50.pdf</a>
ZBUF2	1/8in to 1/8in -0.75mm bore	Filters	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_50.pdf">http://www.chromtech.net.au/pdf2/vici09_50.pdf</a>
ZBRUF21	1/8in to 1/16in -0.75mm bore	Filters	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_50.pdf">http://www.chromtech.net.au/pdf2/vici09_50.pdf</a>
ZBUF4	1/4in to 1/4in -4.60mm bore	Filters	\$33.00	<a href="http://www.chromtech.net.au/pdf2/vici09_50.pdf">http://www.chromtech.net.au/pdf2/vici09_50.pdf</a>
ZBUF1	1/16in to 1/16in -0.75mm bore	Filters	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_50.pdf">http://www.chromtech.net.au/pdf2/vici09_50.pdf</a>
ZUF1	1/16in to 1/16in -0.75mm bore	Filters	\$22.00	<a href="http://www.chromtech.net.au/pdf2/vici09_50.pdf">http://www.chromtech.net.au/pdf2/vici09_50.pdf</a>
ZBRUF1.5	1/16in to 1/32in -0.25mm bore	Filters	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_50.pdf">http://www.chromtech.net.au/pdf2/vici09_50.pdf</a>
ZRUF1.5	1/16in to 1/32in -0.25mm bore	Filters	\$38.00	<a href="http://www.chromtech.net.au/pdf2/vici09_50.pdf">http://www.chromtech.net.au/pdf2/vici09_50.pdf</a>
ZUF2	1/8in to 1/8in -0.75mm bore	Filters	\$27.00	<a href="http://www.chromtech.net.au/pdf2/vici09_50.pdf">http://www.chromtech.net.au/pdf2/vici09_50.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

TCEF811T	1/2in 9.0mm IO-ODx1/16in -bore x1.65mm Side-ODx1/16in -bore x0.75mm	Pre-Column Guard	\$87.00	<a href="http://www.chromtech.net.au/pdf2/vici09_47.pdf">http://www.chromtech.net.au/pdf2/vici09_47.pdf</a>
PCS415F	1/4in x 5cm precolumn system	Pre-Column Guard	\$70.00	<a href="http://www.chromtech.net.au/pdf2/vici09_47.pdf">http://www.chromtech.net.au/pdf2/vici09_47.pdf</a>
SFEC42	0.25ml(4.0mmIDx2cm) Cartridge	Pre-Column Guard	\$26.00	<a href="http://www.chromtech.net.au/pdf2/vici09_47.pdf">http://www.chromtech.net.au/pdf2/vici09_47.pdf</a>
SFECH412	0.25ml(4.0mmIDx2cm) Cartridge Assembly	Pre-Column Guard	\$74.00	<a href="http://www.chromtech.net.au/pdf2/vici09_47.pdf">http://www.chromtech.net.au/pdf2/vici09_47.pdf</a>
PCL412	1/4in x 2cm precolumn insert	Pre-Column Guard	\$38.00	<a href="http://www.chromtech.net.au/pdf2/vici09_47.pdf">http://www.chromtech.net.au/pdf2/vici09_47.pdf</a>
PCS412F	1/4in x 2cm precolumn system	Pre-Column Guard	\$56.00	<a href="http://www.chromtech.net.au/pdf2/vici09_47.pdf">http://www.chromtech.net.au/pdf2/vici09_47.pdf</a>
PCL415	1/4in x 5cm precolumn insert	Pre-Column Guard	\$51.00	<a href="http://www.chromtech.net.au/pdf2/vici09_47.pdf">http://www.chromtech.net.au/pdf2/vici09_47.pdf</a>
TCEF1.5.5C	1/16in 1.0mm IO-ODx1/32in -bore x0.25mm Side-ODx1/32in bore0.25mm	HPLC Column EndFittings	\$87.00	<a href="http://www.chromtech.net.au/pdf2/vici09_46.pdf">http://www.chromtech.net.au/pdf2/vici09_46.pdf</a>
TCEF1.5.5T	1/16in 1.0mm IO-ODx1/32in -bore x0.90mm Side-ODx1/32in -bore x0.25mm	Post Column Teaction Tee	\$87.00	<a href="http://www.chromtech.net.au/pdf2/vici09_46.pdf">http://www.chromtech.net.au/pdf2/vici09_46.pdf</a>
TCEF111	1/16in 1.0mm IO-ODx1/16in -bore x0.25mm Side-ODx1/16in -bore x0.25mm	Post Column Teaction Tee	\$62.00	<a href="http://www.chromtech.net.au/pdf2/vici09_46.pdf">http://www.chromtech.net.au/pdf2/vici09_46.pdf</a>
TCEF211	1/8in 1.0mm IO-ODx1/16in -bore x0.50mm Side-ODx1/16in -bore x0.50mm	Post Column Teaction Tee	\$62.00	<a href="http://www.chromtech.net.au/pdf2/vici09_46.pdf">http://www.chromtech.net.au/pdf2/vici09_46.pdf</a>
TCEF211T	1/8in 1.0mm IO-ODx1/16in -bore x1.65mm Side-ODx1/16in -bore x0.40mm	Post Column Teaction Tee	\$62.00	<a href="http://www.chromtech.net.au/pdf2/vici09_46.pdf">http://www.chromtech.net.au/pdf2/vici09_46.pdf</a>
TCEF411	1/4in 4.6mm IO-ODx1/16in -bore x0.75mm Side-ODx1/16in -bore x0.75mm	Post Column Teaction Tee	\$62.00	<a href="http://www.chromtech.net.au/pdf2/vici09_46.pdf">http://www.chromtech.net.au/pdf2/vici09_46.pdf</a>
TCEF411T	1/4in 4.6mm IO-ODx1/16in -bore x1.65mm Side-ODx1/16in -bore x0.75mm	Post Column Teaction Tee	\$62.00	<a href="http://www.chromtech.net.au/pdf2/vici09_46.pdf">http://www.chromtech.net.au/pdf2/vici09_46.pdf</a>
TCEF421	1/4in 4.6mm IO-ODx1/8in -bore x0.75mm Side-ODx1/16in -bore x0.75mm	Post Column Teaction Tee	\$62.00	<a href="http://www.chromtech.net.au/pdf2/vici09_46.pdf">http://www.chromtech.net.au/pdf2/vici09_46.pdf</a>
TCEF611	3/8in 6.0mm IO-ODx1/16in -bore x0.75mm Side-ODx1/16in -bore x0.75mm	Post Column Teaction Tee	\$76.00	<a href="http://www.chromtech.net.au/pdf2/vici09_46.pdf">http://www.chromtech.net.au/pdf2/vici09_46.pdf</a>
TCEF611T	3/8in 6.0mm IO-ODx1/16in -bore x1.65mm Side-ODx1/16in -bore x0.75mm	Post Column Teaction Tee	\$76.00	<a href="http://www.chromtech.net.au/pdf2/vici09_46.pdf">http://www.chromtech.net.au/pdf2/vici09_46.pdf</a>
TCEF811	1/2in 9.0mm IO-ODx1/16in -bore x0.75mm Side-ODx1/16in -bore x0.75mm	Post Column Teaction Tee	\$87.00	<a href="http://www.chromtech.net.au/pdf2/vici09_46.pdf">http://www.chromtech.net.au/pdf2/vici09_46.pdf</a>
TCEF411C	1/4in 4.6mm IO-ODx1/16in -bore x0.25mm Side-ODx1/16in -bore x0.25mm	Post Column Teaction Tee	\$70.00	<a href="http://www.chromtech.net.au/pdf2/vici09_46.pdf">http://www.chromtech.net.au/pdf2/vici09_46.pdf</a>
2FR8	Stainless 1/2in frits 2u x1.00mm	HPLC Column EndFittings	\$2.00	<a href="http://www.chromtech.net.au/pdf2/vici09_45.pdf">http://www.chromtech.net.au/pdf2/vici09_45.pdf</a>
2FR1KTI	Titanium 1/16in frits 2u x0.75mm	HPLC Column EndFittings	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_45.pdf">http://www.chromtech.net.au/pdf2/vici09_45.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

2FR4-10	Stainless steel 1/4in frits 2u x1.00mm	HPLC Column EndFittings	\$15.00	<a href="http://www.chromtech.net.au/pdf2/vici09_45.pdf">http://www.chromtech.net.au/pdf2/vici09_45.pdf</a>
2FR4HC-10	Hastelloy C 1/4in frits 2u x1.00mmID	HPLC Column EndFittings	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_45.pdf">http://www.chromtech.net.au/pdf2/vici09_45.pdf</a>
2FR4TI-10	Titanium 1/4in frits 2u x1.00mm	HPLC Column EndFittings	\$76.00	<a href="http://www.chromtech.net.au/pdf2/vici09_45.pdf">http://www.chromtech.net.au/pdf2/vici09_45.pdf</a>
10FR4-10	Stainless steel 1/4in frits 10u x1.00mmID	HPLC Column EndFittings	\$15.00	<a href="http://www.chromtech.net.au/pdf2/vici09_45.pdf">http://www.chromtech.net.au/pdf2/vici09_45.pdf</a>
.5FR4-10	Stainless steel 1/4in frits 0.5u x1.00mmID	HPLC Column EndFittings	\$15.00	<a href="http://www.chromtech.net.au/pdf2/vici09_45.pdf">http://www.chromtech.net.au/pdf2/vici09_45.pdf</a>
2FR1KHC	Hastelloy C 1/16in frits 2u x0.75mm	HPLC Column EndFittings	\$6.00	<a href="http://www.chromtech.net.au/pdf2/vici09_45.pdf">http://www.chromtech.net.au/pdf2/vici09_45.pdf</a>
2FR1K	Stainless steel 1.0in frits 2u x1.5mm	HPLC Column EndFittings	\$3.00	<a href="http://www.chromtech.net.au/pdf2/vici09_45.pdf">http://www.chromtech.net.au/pdf2/vici09_45.pdf</a>
10FR2-10	Stainless steel 1/8in frits 10u x1.00mmID	HPLC Column EndFittings	\$15.00	<a href="http://www.chromtech.net.au/pdf2/vici09_45.pdf">http://www.chromtech.net.au/pdf2/vici09_45.pdf</a>
2FR8HC	Hastelloy C 1/2in frits 2u x1.00mmID	HPLC Column EndFittings	\$4.00	<a href="http://www.chromtech.net.au/pdf2/vici09_45.pdf">http://www.chromtech.net.au/pdf2/vici09_45.pdf</a>
2FR1TI-10	Titanium 1/16in frits 2u x0.75mm	HPLC Column EndFittings	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_45.pdf">http://www.chromtech.net.au/pdf2/vici09_45.pdf</a>
2FR6TI	Titanium 3/8in frits 2u x1.00mmID	HPLC Column EndFittings	\$8.00	<a href="http://www.chromtech.net.au/pdf2/vici09_45.pdf">http://www.chromtech.net.au/pdf2/vici09_45.pdf</a>
2FR6HC	Hastelloy C 3/8in frits 2u x1.00mmID	HPLC Column EndFittings	\$3.00	<a href="http://www.chromtech.net.au/pdf2/vici09_45.pdf">http://www.chromtech.net.au/pdf2/vici09_45.pdf</a>
2FR6	Stainless 3/8in frits 2u x1.00mm	HPLC Column EndFittings	\$2.00	<a href="http://www.chromtech.net.au/pdf2/vici09_45.pdf">http://www.chromtech.net.au/pdf2/vici09_45.pdf</a>
2FR8TI	Titanium1/2in frits 2u x1.00mmID	HPLC Column EndFittings	\$9.00	<a href="http://www.chromtech.net.au/pdf2/vici09_45.pdf">http://www.chromtech.net.au/pdf2/vici09_45.pdf</a>
.5FR1-10	Stainless steel 1/16in frits 0.5u x0.75mmID	HPLC Column EndFittings	\$15.00	<a href="http://www.chromtech.net.au/pdf2/vici09_45.pdf">http://www.chromtech.net.au/pdf2/vici09_45.pdf</a>
ECEF616.0	External 3/8in to 1/16in 0.40mm bore -6.0mmID	HPLC Column EndFittings	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_45.pdf">http://www.chromtech.net.au/pdf2/vici09_45.pdf</a>
10FR4HC-10	Hastelloy C 1/4in frits 10u x1.00mmID	HPLC Column EndFittings	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_45.pdf">http://www.chromtech.net.au/pdf2/vici09_45.pdf</a>
ECEF617.0	External 3/8in to 1/16in 0.40mm bore 7.0mmID	HPLC Column EndFittings	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_45.pdf">http://www.chromtech.net.au/pdf2/vici09_45.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

ECEF617.0F	External 3/8in to 1/16in 0.40mm bore -7.0mmID	HPLC Column EndFittings	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_45.pdf">http://www.chromtech.net.au/pdf2/vici09_45.pdf</a>
ECEF819.0	External 1/2in to 1/16in 0.75mm bore 9.0mmID	HPLC Column EndFittings	\$33.00	<a href="http://www.chromtech.net.au/pdf2/vici09_45.pdf">http://www.chromtech.net.au/pdf2/vici09_45.pdf</a>
ECEF8110.0F	External 1/2in to 1/16in 0.75mm bore 10.0mmID	HPLC Column EndFittings	\$36.00	<a href="http://www.chromtech.net.au/pdf2/vici09_45.pdf">http://www.chromtech.net.au/pdf2/vici09_45.pdf</a>
.5FR2-10	Stainless steel 1/8in frits 0.5u x1.00mm	HPLC Column EndFittings	\$15.00	<a href="http://www.chromtech.net.au/pdf2/vici09_45.pdf">http://www.chromtech.net.au/pdf2/vici09_45.pdf</a>
ECEF1K1F	External 1in to 1/16in 0.75mm bore 20.0mmID	HPLC Column EndFittings	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_45.pdf">http://www.chromtech.net.au/pdf2/vici09_45.pdf</a>
2FR2TI-10	Titanium 1/8in frits 2u x1.00mm	HPLC Column EndFittings	\$57.00	<a href="http://www.chromtech.net.au/pdf2/vici09_45.pdf">http://www.chromtech.net.au/pdf2/vici09_45.pdf</a>
.5FR1HC-10	Hastelloy C 1/16in frits 0.5u x0.75mmID	HPLC Column EndFittings	\$27.00	<a href="http://www.chromtech.net.au/pdf2/vici09_45.pdf">http://www.chromtech.net.au/pdf2/vici09_45.pdf</a>
2FR1-10	Stainless steel 1/16in frits 2u x0.75mm	HPLC Column EndFittings	\$15.00	<a href="http://www.chromtech.net.au/pdf2/vici09_45.pdf">http://www.chromtech.net.au/pdf2/vici09_45.pdf</a>
2FR1HC-10	Hastelloy C 1/16in frits 2u x0.75mm	HPLC Column EndFittings	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_45.pdf">http://www.chromtech.net.au/pdf2/vici09_45.pdf</a>
10FR1-10	Stainless steel 1/16in frits 10u x0.75mmID	HPLC Column EndFittings	\$15.00	<a href="http://www.chromtech.net.au/pdf2/vici09_45.pdf">http://www.chromtech.net.au/pdf2/vici09_45.pdf</a>
2FR2-10	Stainless steel 1/8in frits 2u x1.00mm	HPLC Column EndFittings	\$15.00	<a href="http://www.chromtech.net.au/pdf2/vici09_45.pdf">http://www.chromtech.net.au/pdf2/vici09_45.pdf</a>
2FR2HC-10	Hastelloy C 1/8in frits 2u x1.00mm	HPLC Column EndFittings	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_45.pdf">http://www.chromtech.net.au/pdf2/vici09_45.pdf</a>
ECEF1K1	External 1in to 1/16in 0.75mm bore 20.0mmID	HPLC Column EndFittings	\$90.00	<a href="http://www.chromtech.net.au/pdf2/vici09_45.pdf">http://www.chromtech.net.au/pdf2/vici09_45.pdf</a>
ECEF819.0F	External 1/2in to 1/16in 0.75mm bore 9.0mmID	HPLC Column EndFittings	\$36.00	<a href="http://www.chromtech.net.au/pdf2/vici09_45.pdf">http://www.chromtech.net.au/pdf2/vici09_45.pdf</a>
ECEF616.0F	External 3/8in to 1/16in 0.40mm bore 6.0mmID	HPLC Column EndFittings	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_45.pdf">http://www.chromtech.net.au/pdf2/vici09_45.pdf</a>
ECEF8110.0	External 1/2in to 1/16in 0.75mm bore 10.0mmID	HPLC Column EndFittings	\$33.00	<a href="http://www.chromtech.net.au/pdf2/vici09_45.pdf">http://www.chromtech.net.au/pdf2/vici09_45.pdf</a>
CEF414.0	Internal w/o frit 1/4in to 1/16in 0.4mm bore -4.0mmID	HPLC Column EndFittings	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_44.pdf">http://www.chromtech.net.au/pdf2/vici09_44.pdf</a>
CEF414.6	Internal w/o frit 1/4in to 1/16in 0.4mm bore -4.6mmID	HPLC Column EndFittings	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_44.pdf">http://www.chromtech.net.au/pdf2/vici09_44.pdf</a>



## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

CEF212.0	Internal w/o frit 1/8in to 1/16in 0.25mm bore -2.0mmID	HPLC Column EndFittings	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_44.pdf">http://www.chromtech.net.au/pdf2/vici09_44.pdf</a>
CEF413.0	Internal w/o frit 1/4in to 1/16in 0.4mm bore -3.0 mmID	HPLC Column EndFittings	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_44.pdf">http://www.chromtech.net.au/pdf2/vici09_44.pdf</a>
CEF414.0F	Internal w/2u frit 1/4in to 1/16in 0.4mm bore -4.0mmID	HPLC Column EndFittings	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_44.pdf">http://www.chromtech.net.au/pdf2/vici09_44.pdf</a>
CEF414.6F	Internal w/2u frit 1/4in to 1/16in 0.4mm bore -4.6mmID	HPLC Column EndFittings	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_44.pdf">http://www.chromtech.net.au/pdf2/vici09_44.pdf</a>
CEF413.0F	Internal w/2u frit 1/4in to 1/16in 0.4mm bore -3.0mmID	HPLC Column EndFittings	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_44.pdf">http://www.chromtech.net.au/pdf2/vici09_44.pdf</a>
CEF412.1F	Internal Rem 2u frit 1/4in to 1/16in 0.4mm bore -2.1mmID	HPLC Column EndFittings	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_44.pdf">http://www.chromtech.net.au/pdf2/vici09_44.pdf</a>
CEF412.1	Internalw/o frit 1/4in to 1/16in 0.4mm bore -2.1mmID	HPLC Column EndFittings	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_44.pdf">http://www.chromtech.net.au/pdf2/vici09_44.pdf</a>
CEF212.0F	Internal Rem 2u frit 1/8in to 1/16in 0.25mm bore -2.0mmID	HPLC Column EndFittings	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_44.pdf">http://www.chromtech.net.au/pdf2/vici09_44.pdf</a>
CEF211.0F	Internal Rem 2u frit 1/8in to 1/16in 0.25mm bore -1.0mmID	HPLC Column EndFittings	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_43.pdf">http://www.chromtech.net.au/pdf2/vici09_43.pdf</a>
ECEF211.0	External w/o frit 1/8in to 1/16in 0.25mm bore -1.0mmID	HPLC Column EndFittings	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_43.pdf">http://www.chromtech.net.au/pdf2/vici09_43.pdf</a>
ECEF111.0F	External Rem 2u frit 1/16intol/16in 0.25mm bore -1.0mmID	HPLC Column EndFittings	\$22.00	<a href="http://www.chromtech.net.au/pdf2/vici09_43.pdf">http://www.chromtech.net.au/pdf2/vici09_43.pdf</a>
ECEF111.0	External w/o frit 1/16intol/16in 0.25mm bore -1.0mmID	HPLC Column EndFittings	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_43.pdf">http://www.chromtech.net.au/pdf2/vici09_43.pdf</a>
ECEF211.0F	External Rem 2u frit 1/8in to 1/16in 0.25mm bore -1.0mm ID	HPLC Column EndFittings	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_43.pdf">http://www.chromtech.net.au/pdf2/vici09_43.pdf</a>
CEF211.0	Internal w/o frit 1/8in to 1/16in 0.25mm bore -1.0mmID	HPLC Column EndFittings	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_43.pdf">http://www.chromtech.net.au/pdf2/vici09_43.pdf</a>
CEF2.51.0F	Internal Rem 2u frit 1/8in to 1/16in 0.25mm bore -1.0mmID	HPLC Column EndFittings	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_43.pdf">http://www.chromtech.net.au/pdf2/vici09_43.pdf</a>
CEF2.51.0	Internal w/o frit 1/8in to 1/32in 0.25mm bore -1.0mmID	HPLC Column EndFittings	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_43.pdf">http://www.chromtech.net.au/pdf2/vici09_43.pdf</a>
CEF1F	Internal Rem 2u frit 1/16in to 1/16in 0.25mm bore -1.0mmID	HPLC Column EndFittings	\$22.00	<a href="http://www.chromtech.net.au/pdf2/vici09_43.pdf">http://www.chromtech.net.au/pdf2/vici09_43.pdf</a>
CEF1	Internal w/o frit 1/16in to 1/16in 0.25mmID -1.0mmID	HPLC Column EndFittings	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_43.pdf">http://www.chromtech.net.au/pdf2/vici09_43.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

CEF1.5F	Internal 1/16 to 1/32 Rem 2u frit 0.25mm bore -1.0mmID	HPLC Column EndFittings	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_43.pdf">http://www.chromtech.net.au/pdf2/vici09_43.pdf</a>
CEF1.5	Internal 1/16 to 1/32in w/o frit	HPLC Column EndFittings	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_43.pdf">http://www.chromtech.net.au/pdf2/vici09_43.pdf</a>
SI4G	Septum, Low Bleed.10-pk	Valco Fittings- SyringeAdaptors	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_41.pdf">http://www.chromtech.net.au/pdf2/vici09_41.pdf</a>
ZF2SI	Septum injector nut -Replacement support	Valco Fittings- SyringeAdaptors	\$8.00	<a href="http://www.chromtech.net.au/pdf2/vici09_41.pdf">http://www.chromtech.net.au/pdf2/vici09_41.pdf</a>
EN2SI	Septum injector nut -w/support	Valco Fittings- SyringeAdaptors	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_41.pdf">http://www.chromtech.net.au/pdf2/vici09_41.pdf</a>
C-LFP	Loop Fill Port Assembly	Valco Fittings- SyringeAdaptors	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_41.pdf">http://www.chromtech.net.au/pdf2/vici09_41.pdf</a>
VISL-1	Fill port Liner for VISF-1	Valco Fittings- PipeAdaptors	\$3.00	<a href="http://www.chromtech.net.au/pdf2/vici09_40.pdf">http://www.chromtech.net.au/pdf2/vici09_40.pdf</a>
VISL-2	Fill port Liner for VISF-2	Valco Fittings- PipeAdaptors	\$3.00	<a href="http://www.chromtech.net.au/pdf2/vici09_40.pdf">http://www.chromtech.net.au/pdf2/vici09_40.pdf</a>
VISF-A	Fill port:use w/> 2in long 22g needle Syringe 1/8in fitting	Valco Fittings- PipeAdaptors	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_40.pdf">http://www.chromtech.net.au/pdf2/vici09_40.pdf</a>
VISF-2	for 2in 22# needle 1/16in fitting metal valves	Valco Fittings- PipeAdaptors	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_40.pdf">http://www.chromtech.net.au/pdf2/vici09_40.pdf</a>
VISF.5FPK	Fill Port 1/32 Injectors, blunt 26G	Valco Fittings- PipeAdaptors	\$15.00	<a href="http://www.chromtech.net.au/pdf2/vici09_40.pdf">http://www.chromtech.net.au/pdf2/vici09_40.pdf</a>
VISF-1	for .75in 22# needle 1/16in fitting metal valves	Valco Fittings- PipeAdaptors	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_40.pdf">http://www.chromtech.net.au/pdf2/vici09_40.pdf</a>
ZVISF-1PEH	Zero Dead Volume Fill Port	Valco Fittings- PipeAdaptors	\$51.00	<a href="http://www.chromtech.net.au/pdf2/vici09_40.pdf">http://www.chromtech.net.au/pdf2/vici09_40.pdf</a>
VISL-1H	Fill port Liner for VISF-1H	Valco Fittings- PipeAdaptors	\$3.00	<a href="http://www.chromtech.net.au/pdf2/vici09_40.pdf">http://www.chromtech.net.au/pdf2/vici09_40.pdf</a>
ZVISF-1PFAH	Zero Dead Volume Fill Port	Valco Fittings- PipeAdaptors	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_40.pdf">http://www.chromtech.net.au/pdf2/vici09_40.pdf</a>
C-VISF-1	Fill Port C22Z,C62Z Injectors	Valco Fittings- PipeAdaptors	\$14.00	<a href="http://www.chromtech.net.au/pdf2/vici09_40.pdf">http://www.chromtech.net.au/pdf2/vici09_40.pdf</a>
C-VISF-1H	Fill Port C22Z,C62Z Injectors	Valco Fittings- PipeAdaptors	\$14.00	<a href="http://www.chromtech.net.au/pdf2/vici09_40.pdf">http://www.chromtech.net.au/pdf2/vici09_40.pdf</a>
ZVISF-1PE	Ferrule for VISF-1 or -2	Valco Fittings- PipeAdaptors	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_40.pdf">http://www.chromtech.net.au/pdf2/vici09_40.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

ZF1VISF	Fill port Ferrule for VISF-1 or VISF-2	Valco Fittings- PipeAdaptors	\$6.00	<a href="http://www.chromtech.net.au/pdf2/vici09_40.pdf">http://www.chromtech.net.au/pdf2/vici09_40.pdf</a>
PEA44	1/4in NPT male to 1/4in external fitting -4.6mm bore	Valco Fittings- PipeAdaptors	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_39.pdf">http://www.chromtech.net.au/pdf2/vici09_39.pdf</a>
PEA22	1/8in NPT male to 1/8in external fitting -2.0mm bore	Valco Fittings- PipeAdaptors	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_39.pdf">http://www.chromtech.net.au/pdf2/vici09_39.pdf</a>
PEA42	1/4in NPT male to 1/8in external fitting -2.0mm bore	Valco Fittings- PipeAdaptors	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_39.pdf">http://www.chromtech.net.au/pdf2/vici09_39.pdf</a>
FPEA84	1/8in NPT female to 1/4in external fitting -4.6mm bore	Valco Fittings- PipeAdaptors	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_39.pdf">http://www.chromtech.net.au/pdf2/vici09_39.pdf</a>
PEA82	1/2in NPT male to 1/8in external fitting -2.0mm bore	Valco Fittings- PipeAdaptors	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_39.pdf">http://www.chromtech.net.au/pdf2/vici09_39.pdf</a>
PEA84	1/2in NPT male to 1/4in external fitting -4.6mm bore	Valco Fittings- PipeAdaptors	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_39.pdf">http://www.chromtech.net.au/pdf2/vici09_39.pdf</a>
FPEA22	1/8in NPT female to 1/8in external fitting -2.0mm bore	Valco Fittings- PipeAdaptors	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_39.pdf">http://www.chromtech.net.au/pdf2/vici09_39.pdf</a>
FPEA24	1/8in NPT female to 1/4in external fitting -4.6mm bore	Valco Fittings- PipeAdaptors	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_39.pdf">http://www.chromtech.net.au/pdf2/vici09_39.pdf</a>
FPEA42	1/4in NPT female to 1/8in external fitting -2.0mm bore	Valco Fittings- PipeAdaptors	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_39.pdf">http://www.chromtech.net.au/pdf2/vici09_39.pdf</a>
FPEA44	1/8in NPT female to 1/4in external fitting -4.6mm bore	Valco Fittings- PipeAdaptors	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_39.pdf">http://www.chromtech.net.au/pdf2/vici09_39.pdf</a>
FPEA82	1/2in NPT female to:1/8in external fitting -2.0mm bore	Valco Fittings- PipeAdaptors	\$38.00	<a href="http://www.chromtech.net.au/pdf2/vici09_39.pdf">http://www.chromtech.net.au/pdf2/vici09_39.pdf</a>
PEA24	1/8in NPT male to 1/4in external fitting -2.0mm bore	Valco Fittings- PipeAdaptors	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_39.pdf">http://www.chromtech.net.au/pdf2/vici09_39.pdf</a>
PZA82	1/2in NPT male to 1/8in ZDV fitting -1.0mm bore	Valco Fittings- PipeAdaptors	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_38.pdf">http://www.chromtech.net.au/pdf2/vici09_38.pdf</a>
PZA21	1/8in NPT male to 1/16in ZDV fitting -1.0mm bore	Valco Fittings- PipeAdaptors	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_38.pdf">http://www.chromtech.net.au/pdf2/vici09_38.pdf</a>
PZA21T	1/8in NPT male to 1/16in ZDV fitting -1/16in bore	Valco Fittings- PipeAdaptors	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_38.pdf">http://www.chromtech.net.au/pdf2/vici09_38.pdf</a>
FPZA42	1/4in NPT female to 1/8in ZDV fitting -1.0mm bore	Valco Fittings- PipeAdaptors	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_38.pdf">http://www.chromtech.net.au/pdf2/vici09_38.pdf</a>
FPZA41	1/4in NPT female to 1/16in ZDV fitting -1.0mm bore	Valco Fittings- PipeAdaptors	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_38.pdf">http://www.chromtech.net.au/pdf2/vici09_38.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

FPZA22L	1/8in NPT female to 1/8in ZDV fitting -2.0mm bore	Valco Fittings-PipeAdaptors	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_38.pdf">http://www.chromtech.net.au/pdf2/vici09_38.pdf</a>
FPZA22	1/8in NPT female to 1/8in ZDV fitting -1.0mm bore	Valco Fittings-PipeAdaptors	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_38.pdf">http://www.chromtech.net.au/pdf2/vici09_38.pdf</a>
FPZA21	1/8in NPT female to 1/16in ZDV fitting -1.0mm bore	Valco Fittings-PipeAdaptors	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_38.pdf">http://www.chromtech.net.au/pdf2/vici09_38.pdf</a>
PZA82L	1/2in NPT male to 1/8in ZDV fitting -2.0mm bore	Valco Fittings-PipeAdaptors	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_38.pdf">http://www.chromtech.net.au/pdf2/vici09_38.pdf</a>
FPZA82L	1/2in NPT female to 1/8in ZDV fitting -2.0mm bore	Valco Fittings-PipeAdaptors	\$38.00	<a href="http://www.chromtech.net.au/pdf2/vici09_38.pdf">http://www.chromtech.net.au/pdf2/vici09_38.pdf</a>
PZA81	1/2in NPT male to 1/16in ZDV fitting -1.0mm bore	Valco Fittings-PipeAdaptors	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_38.pdf">http://www.chromtech.net.au/pdf2/vici09_38.pdf</a>
PZA44L	1/4in NPT male to 1/4in ZDV fitting -4.6mm bore	Valco Fittings-PipeAdaptors	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_38.pdf">http://www.chromtech.net.au/pdf2/vici09_38.pdf</a>
PZA42L	1/4in NPT male to 1/8in ZDV fitting -2.0mm bore	Valco Fittings-PipeAdaptors	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_38.pdf">http://www.chromtech.net.au/pdf2/vici09_38.pdf</a>
PZA42	1/4in NPT male to 1/8in ZDV fitting -1.0mm bore	Valco Fittings-PipeAdaptors	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_38.pdf">http://www.chromtech.net.au/pdf2/vici09_38.pdf</a>
PZA41	1/4in NPT male to 1/16in ZDV fitting -1.0mm bore	Valco Fittings-PipeAdaptors	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_38.pdf">http://www.chromtech.net.au/pdf2/vici09_38.pdf</a>
PZA84L	1/2in NPT male to 1/4in ZDV fitting -4.6mm bore	Valco Fittings-PipeAdaptors	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_38.pdf">http://www.chromtech.net.au/pdf2/vici09_38.pdf</a>
FPZA44L	1/4in NPT female to 1/4in ZDV fitting -4.6mm bore	Valco Fittings-PipeAdaptors	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_38.pdf">http://www.chromtech.net.au/pdf2/vici09_38.pdf</a>
FPZA81	1/2in NPT female to 1/16in ZDV fitting -1.0mm bore	Valco Fittings-PipeAdaptors	\$38.00	<a href="http://www.chromtech.net.au/pdf2/vici09_38.pdf">http://www.chromtech.net.au/pdf2/vici09_38.pdf</a>
FPZA82	1/2in NPT female to 1/8in ZDV fitting -1.0mm bore	Valco Fittings-PipeAdaptors	\$38.00	<a href="http://www.chromtech.net.au/pdf2/vici09_38.pdf">http://www.chromtech.net.au/pdf2/vici09_38.pdf</a>
FPZA42L	1/4in NPT female to 1/8in ZDV fitting -2.0mm bore	Valco Fittings-PipeAdaptors	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_38.pdf">http://www.chromtech.net.au/pdf2/vici09_38.pdf</a>
FPZA84L	1/2in NPT female to 1/4in ZDV fitting -4.6mm bore	Valco Fittings-PipeAdaptors	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_38.pdf">http://www.chromtech.net.au/pdf2/vici09_38.pdf</a>
PZA22	1/8in NPT male to 1/8in ZDV fitting -1.0mm bore	Valco Fittings-PipeAdaptors	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_38.pdf">http://www.chromtech.net.au/pdf2/vici09_38.pdf</a>
FP2Z3M21	Two 1/8in female pipe to three 1/16inZDV fittings -1.0mm bore	Valco Fittings-Manifold Pipe	\$98.00	<a href="http://www.chromtech.net.au/pdf2/vici09_37.pdf">http://www.chromtech.net.au/pdf2/vici09_37.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

FP1Z3M44	One 1/4in female pipe to three 1/4inZDV fittings -4.6mm bore	Valco Fittings-Manifold Pipe	\$80.00	<a href="http://www.chromtech.net.au/pdf2/vici09_37.pdf">http://www.chromtech.net.au/pdf2/vici09_37.pdf</a>
FP2Z3M42	Two 1/4in female pipe to three 1/8inZDV fittings -2.0mm bore	Valco Fittings-Manifold Pipe	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_37.pdf">http://www.chromtech.net.au/pdf2/vici09_37.pdf</a>
FP2Z3M41	Two 1/4in female pipe to three 1/16inZDV fittings -1.0mm bore	Valco Fittings-Manifold Pipe	\$98.00	<a href="http://www.chromtech.net.au/pdf2/vici09_37.pdf">http://www.chromtech.net.au/pdf2/vici09_37.pdf</a>
FP2Z3M24	Two 1/8in female pipe to three 1/4inZDV fittings -4.6mm bore	Valco Fittings-Manifold Pipe	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_37.pdf">http://www.chromtech.net.au/pdf2/vici09_37.pdf</a>
FP2Z3M44	Two 1/4in female pipe to three 1/4inZDV fittings -4.6mm bore	Valco Fittings-Manifold Pipe	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_37.pdf">http://www.chromtech.net.au/pdf2/vici09_37.pdf</a>
FP1Z3M42	One 1/4in female pipe to three 1/8inZDV fittings -2.0mm bore	Valco Fittings-Manifold Pipe	\$80.00	<a href="http://www.chromtech.net.au/pdf2/vici09_37.pdf">http://www.chromtech.net.au/pdf2/vici09_37.pdf</a>
FP1Z3M41	One 1/4in female pipe to three 1/8inZDV fittings -1.0mm bore	Valco Fittings-Manifold Pipe	\$86.00	<a href="http://www.chromtech.net.au/pdf2/vici09_37.pdf">http://www.chromtech.net.au/pdf2/vici09_37.pdf</a>
FP1Z3M24	One 1/8in female pipe to three 1/4inZDV fittings -4.6mm bore	Valco Fittings-Manifold Pipe	\$80.00	<a href="http://www.chromtech.net.au/pdf2/vici09_37.pdf">http://www.chromtech.net.au/pdf2/vici09_37.pdf</a>
FP1Z3M22	One 1/8in female pipe to three 1/8inZDV fittings -2.0mm bore	Valco Fittings-Manifold Pipe	\$80.00	<a href="http://www.chromtech.net.au/pdf2/vici09_37.pdf">http://www.chromtech.net.au/pdf2/vici09_37.pdf</a>
FP1Z3M21	One 1/8in female pipe to three 1/8inZDV fittings -1.0mm bore	Valco Fittings-Manifold Pipe	\$86.00	<a href="http://www.chromtech.net.au/pdf2/vici09_37.pdf">http://www.chromtech.net.au/pdf2/vici09_37.pdf</a>
FP2Z3M22	Two 1/8in female pipe to three 1/8inZDV fittings -2.0mm bore	Valco Fittings-Manifold Pipe	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_37.pdf">http://www.chromtech.net.au/pdf2/vici09_37.pdf</a>
ZLTA41	1/4in to 1/16in x2.075in long -1/16in bore	Valco Fittings-Special	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_36.pdf">http://www.chromtech.net.au/pdf2/vici09_36.pdf</a>
ZTA41	1/4in to 1/16in x0.975in long -1/16in bore	Valco Fittings-Special	\$15.00	<a href="http://www.chromtech.net.au/pdf2/vici09_36.pdf">http://www.chromtech.net.au/pdf2/vici09_36.pdf</a>
ZXLTA41	1/4in to 1/16in x2.800in long -1/16in bore	Valco Fittings-Special	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_36.pdf">http://www.chromtech.net.au/pdf2/vici09_36.pdf</a>
ZbAA1	Bulkhead union -Standard 1/8 Aerosol to 1/16in ZDV	Valco Fittings-Special	\$38.00	<a href="http://www.chromtech.net.au/pdf2/vici09_36.pdf">http://www.chromtech.net.au/pdf2/vici09_36.pdf</a>
EZR42	1/4in ext to 1/8in int x1.00mm bore	Valco Fittings-Int Red	\$33.00	<a href="http://www.chromtech.net.au/pdf2/vici09_35.pdf">http://www.chromtech.net.au/pdf2/vici09_35.pdf</a>
EZR42T	1/4in ext to 1/8in int x1/8in bore	Valco Fittings-Int Red	\$33.00	<a href="http://www.chromtech.net.au/pdf2/vici09_35.pdf">http://www.chromtech.net.au/pdf2/vici09_35.pdf</a>
EZR41T	1/4in ext to 1/16in int x1/16in bore	Valco Fittings-Int Red	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_35.pdf">http://www.chromtech.net.au/pdf2/vici09_35.pdf</a>
EZR41	1/4in ext to 1/16in int x1.00mm bore	Valco Fittings-Int Red	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_35.pdf">http://www.chromtech.net.au/pdf2/vici09_35.pdf</a>
EZR21T	1/8in ext to 1/16in int x1/16in bore	Valco Fittings-Int Red	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_35.pdf">http://www.chromtech.net.au/pdf2/vici09_35.pdf</a>
EZR2.5	1/8in ext to 1/32in int x0.25mm bore	Valco Fittings-Int Red	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_35.pdf">http://www.chromtech.net.au/pdf2/vici09_35.pdf</a>
EZR1.5T	1/16in ext to 1/32in int x1/32in bore	Valco Fittings-Int Red	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_35.pdf">http://www.chromtech.net.au/pdf2/vici09_35.pdf</a>
EZR1.5	1/16in ext to 1/32in int x0.25mm bore	Valco Fittings-Int Red	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_35.pdf">http://www.chromtech.net.au/pdf2/vici09_35.pdf</a>
EZA1	Injector/Detector 1/16in ext to 1/16in int	Valco Fittings-Int Red	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_35.pdf">http://www.chromtech.net.au/pdf2/vici09_35.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

EZA2	Injector/Detector 1/8in ext to 1/8in int	Valco Fittings-Int Red	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_35.pdf">http://www.chromtech.net.au/pdf2/vici09_35.pdf</a>
EZR21	1/8in ext to 1/16in int x0.50mm bore	Valco Fittings-Int Red	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_35.pdf">http://www.chromtech.net.au/pdf2/vici09_35.pdf</a>
IZR42	W/o frit 1/4in to 1/8in -1.00mm bore	Valco Fittings-Int Red	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_34.pdf">http://www.chromtech.net.au/pdf2/vici09_34.pdf</a>
IZR42F	w/2u frit 1/4in to 1/8in -1.00mm bore	Valco Fittings-Int Red	\$26.00	<a href="http://www.chromtech.net.au/pdf2/vici09_34.pdf">http://www.chromtech.net.au/pdf2/vici09_34.pdf</a>
C360IZR.5S6PKG	1/32 to 360um PEEK	Valco Fittings-Int Red	\$44.00	<a href="http://www.chromtech.net.au/pdf2/vici09_34.pdf">http://www.chromtech.net.au/pdf2/vici09_34.pdf</a>
IZR21CF	w/2u frit 1/8in to 1/16in -0.25mm bore	Valco Fittings-Int Red	\$27.00	<a href="http://www.chromtech.net.au/pdf2/vici09_34.pdf">http://www.chromtech.net.au/pdf2/vici09_34.pdf</a>
IZR1.5LF	w/2u frit 1/16in to 1/32in -0.25mm bore	Valco Fittings-Int Red	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_34.pdf">http://www.chromtech.net.au/pdf2/vici09_34.pdf</a>
IZR1.5F	w/2u frit 1/16in to 1/32in -0.25mm bore	Valco Fittings-Int Red	\$36.00	<a href="http://www.chromtech.net.au/pdf2/vici09_34.pdf">http://www.chromtech.net.au/pdf2/vici09_34.pdf</a>
IZR42L	W/o frit 1/4in to 1/8in -2.00mm bore	Valco Fittings-Int Red	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_34.pdf">http://www.chromtech.net.au/pdf2/vici09_34.pdf</a>
IZR21LF	w/2u frit 1/8in to 1/16in -1.00mm bore	Valco Fittings-Int Red	\$24.00	<a href="http://www.chromtech.net.au/pdf2/vici09_34.pdf">http://www.chromtech.net.au/pdf2/vici09_34.pdf</a>
IZR41	W/o frit 1/4in to 1/16in -1.00mm bore	Valco Fittings-Int Red	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_34.pdf">http://www.chromtech.net.au/pdf2/vici09_34.pdf</a>
IZR21T	W/o frit 1/8in to 1/16in -1/16in bore	Valco Fittings-Int Red	\$16.00	<a href="http://www.chromtech.net.au/pdf2/vici09_34.pdf">http://www.chromtech.net.au/pdf2/vici09_34.pdf</a>
IZR21L	W/o frit 1/8in to 1/16in -1.00mm bore	Valco Fittings-Int Red	\$16.00	<a href="http://www.chromtech.net.au/pdf2/vici09_34.pdf">http://www.chromtech.net.au/pdf2/vici09_34.pdf</a>
IZR21	Internal Reducing Fitting W/o frit 1/8in to 1/16in -0.50mm bore	Valco Fittings-Int Red	\$16.00	<a href="http://www.chromtech.net.au/pdf2/vici09_34.pdf">http://www.chromtech.net.au/pdf2/vici09_34.pdf</a>
IZR21C	W/o frit 1/8in to 1/16in -0.25mm bore	Valco Fittings-Int Red	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_34.pdf">http://www.chromtech.net.au/pdf2/vici09_34.pdf</a>
IZR1.5T	W/o frit 1/16in to 1/32in -1/32in bore	Valco Fittings-Int Red	\$26.00	<a href="http://www.chromtech.net.au/pdf2/vici09_34.pdf">http://www.chromtech.net.au/pdf2/vici09_34.pdf</a>
IZR1.5	W/o frit 1/16in to 1/32in -0.25mm bore	Valco Fittings-Int Red	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_34.pdf">http://www.chromtech.net.au/pdf2/vici09_34.pdf</a>
C360IZR.5S6	1/32 to 360um Stainless	Valco Fittings-Int Red	\$44.00	<a href="http://www.chromtech.net.au/pdf2/vici09_34.pdf">http://www.chromtech.net.au/pdf2/vici09_34.pdf</a>
IZR42LF	w/2u frit 1/4in to 1/8in -2.00mm bore	Valco Fittings-Int Red	\$26.00	<a href="http://www.chromtech.net.au/pdf2/vici09_34.pdf">http://www.chromtech.net.au/pdf2/vici09_34.pdf</a>
IZR41F	w/2u frit 1/4in to 1/16in -1.00mm bore	Valco Fittings-Int Red	\$26.00	<a href="http://www.chromtech.net.au/pdf2/vici09_34.pdf">http://www.chromtech.net.au/pdf2/vici09_34.pdf</a>
IZR1.5L	w/2u frit 1/16in to 1/32in -0.50mm bore	Valco Fittings-Int Red	\$26.00	<a href="http://www.chromtech.net.au/pdf2/vici09_34.pdf">http://www.chromtech.net.au/pdf2/vici09_34.pdf</a>
IZR21F	w/2u frit 1/8in to 1/16in -0.50mm bore	Valco Fittings-Int Red	\$24.00	<a href="http://www.chromtech.net.au/pdf2/vici09_34.pdf">http://www.chromtech.net.au/pdf2/vici09_34.pdf</a>
Z4M1PK	PEEK 4 inlets 0.25mm inlet x0.75mm outlet bores	Valco Fittings-Manifolds	\$147.00	<a href="http://www.chromtech.net.au/pdf2/vici09_33.pdf">http://www.chromtech.net.au/pdf2/vici09_33.pdf</a>
Z12M1PK	PEEK 12 inlets 0.25mm inlet x0.75mm outlet bores	Valco Fittings-Manifolds	\$294.00	<a href="http://www.chromtech.net.au/pdf2/vici09_33.pdf">http://www.chromtech.net.au/pdf2/vici09_33.pdf</a>
Z8M1	Stainless steel 8 inlets 0.25mm inlet x0.75mm outlet bores	Valco Fittings-Manifolds	\$198.00	<a href="http://www.chromtech.net.au/pdf2/vici09_33.pdf">http://www.chromtech.net.au/pdf2/vici09_33.pdf</a>
Z14M1PK	PEEK 14 inlets 0.25mm inlet x0.75mm outlet bores	Valco Fittings-Manifolds	\$332.00	<a href="http://www.chromtech.net.au/pdf2/vici09_33.pdf">http://www.chromtech.net.au/pdf2/vici09_33.pdf</a>
Z6M1PK	PEEK 6 inlets 0.25mm inlet x0.75mm outlet bores	Valco Fittings-Manifolds	\$198.00	<a href="http://www.chromtech.net.au/pdf2/vici09_33.pdf">http://www.chromtech.net.au/pdf2/vici09_33.pdf</a>
Z8M1PK	PEEK 8 inlets 0.25mm inlet x0.75mm outlet bores	Valco Fittings-Manifolds	\$224.00	<a href="http://www.chromtech.net.au/pdf2/vici09_33.pdf">http://www.chromtech.net.au/pdf2/vici09_33.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

Z6M1	Stainless steel 6 inlets 0.25mm inlet x0.75mm outlet bores	Valco Fittings-Manifolds	\$147.00	<a href="http://www.chromtech.net.au/pdf2/vici09_33.pdf">http://www.chromtech.net.au/pdf2/vici09_33.pdf</a>
Z12M2	300 SS 12 inlets 2.00mm inlet x2.00mm outlet bores	Valco Fittings-Manifolds	\$294.00	<a href="http://www.chromtech.net.au/pdf2/vici09_33.pdf">http://www.chromtech.net.au/pdf2/vici09_33.pdf</a>
Z4M1	Stainless steel 4 inlets 0.25mm inlet x0.75mm outlet bores	Valco Fittings-Manifolds	\$111.00	<a href="http://www.chromtech.net.au/pdf2/vici09_33.pdf">http://www.chromtech.net.au/pdf2/vici09_33.pdf</a>
Z8M2	300 SS 8 inlets 2.00mm inlet x2.00mm outlet bores	Valco Fittings-Manifolds	\$224.00	<a href="http://www.chromtech.net.au/pdf2/vici09_33.pdf">http://www.chromtech.net.au/pdf2/vici09_33.pdf</a>
Z16M1PK	PEEK 16 inlets 0.25mm inlet x0.75mm outlet bores	Valco Fittings-Manifolds	\$369.00	<a href="http://www.chromtech.net.au/pdf2/vici09_33.pdf">http://www.chromtech.net.au/pdf2/vici09_33.pdf</a>
Z4M2	300 SS 4 inlets 2.00mm inlet x2.00mm outlet bores	Valco Fittings-Manifolds	\$147.00	<a href="http://www.chromtech.net.au/pdf2/vici09_33.pdf">http://www.chromtech.net.au/pdf2/vici09_33.pdf</a>
Z6M2	300 SS 6 inlets 2.00mm inlet x2.00mm outlet bores	Valco Fittings-Manifolds	\$198.00	<a href="http://www.chromtech.net.au/pdf2/vici09_33.pdf">http://www.chromtech.net.au/pdf2/vici09_33.pdf</a>
Z10M1	Stainless steel 10 inlets 0.25mm inlet x0.75mm outlet bores	Valco Fittings-Manifolds	\$224.00	<a href="http://www.chromtech.net.au/pdf2/vici09_33.pdf">http://www.chromtech.net.au/pdf2/vici09_33.pdf</a>
Z10M2	300 SS 10 inlets 2.00mm inlet x2.00mm outlet bores	Valco Fittings-Manifolds	\$255.00	<a href="http://www.chromtech.net.au/pdf2/vici09_33.pdf">http://www.chromtech.net.au/pdf2/vici09_33.pdf</a>
Z14M1	Stainless steel 14 inlets 0.25mm inlet x0.75mm outlet bores	Valco Fittings-Manifolds	\$294.00	<a href="http://www.chromtech.net.au/pdf2/vici09_33.pdf">http://www.chromtech.net.au/pdf2/vici09_33.pdf</a>
Z12M1	Stainless steel 12 inlets 0.25mm inlet x0.75mm outlet bores	Valco Fittings-Manifolds	\$255.00	<a href="http://www.chromtech.net.au/pdf2/vici09_33.pdf">http://www.chromtech.net.au/pdf2/vici09_33.pdf</a>
Z10M1PK	PEEK 10 inlets 0.25mm inlet x0.75mm outlet bores	Valco Fittings-Manifolds	\$255.00	<a href="http://www.chromtech.net.au/pdf2/vici09_33.pdf">http://www.chromtech.net.au/pdf2/vici09_33.pdf</a>
ZX1M	300 Series SS 1/16inx0.50mm bore, Cross	Valco Fittings-Cross	\$51.00	<a href="http://www.chromtech.net.au/pdf2/vici09_32.pdf">http://www.chromtech.net.au/pdf2/vici09_32.pdf</a>
ZX4L	300 Series SS 1/4inx4.00mm bore, Cross	Valco Fittings-Cross	\$51.00	<a href="http://www.chromtech.net.au/pdf2/vici09_32.pdf">http://www.chromtech.net.au/pdf2/vici09_32.pdf</a>
ZX1C	300 Series SS 1/32inx0.25mm bore, Cross	Valco Fittings-Cross	\$58.00	<a href="http://www.chromtech.net.au/pdf2/vici09_32.pdf">http://www.chromtech.net.au/pdf2/vici09_32.pdf</a>
ZX1	300 Series SS 1/16inx0.75mm bore, Cross	Valco Fittings-Cross	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_32.pdf">http://www.chromtech.net.au/pdf2/vici09_32.pdf</a>
ZX1L	300 Series SS 1/16inx1.00mm bore, Cross	Valco Fittings-Cross	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_32.pdf">http://www.chromtech.net.au/pdf2/vici09_32.pdf</a>
ZX2	300 Series SS 1/8inx0.75mm bore, Cross	Valco Fittings-Cross	\$46.00	<a href="http://www.chromtech.net.au/pdf2/vici09_32.pdf">http://www.chromtech.net.au/pdf2/vici09_32.pdf</a>
ZX2L	300 Series SS 1/8inx2.00mm bore, Cross	Valco Fittings-Cross	\$46.00	<a href="http://www.chromtech.net.au/pdf2/vici09_32.pdf">http://www.chromtech.net.au/pdf2/vici09_32.pdf</a>
ZX4	300 Series SS 1/4inx1.00mm bore, Cross	Valco Fittings-Cross	\$51.00	<a href="http://www.chromtech.net.au/pdf2/vici09_32.pdf">http://www.chromtech.net.au/pdf2/vici09_32.pdf</a>
ZX.5	300 Series SS 1/32inx0.25mm bore, Cross	Valco Fittings-Cross	\$58.00	<a href="http://www.chromtech.net.au/pdf2/vici09_32.pdf">http://www.chromtech.net.au/pdf2/vici09_32.pdf</a>
ZX.5L	300 Series SS 1/16inx0.50mm bore, Cross	Valco Fittings-Cross	\$58.00	<a href="http://www.chromtech.net.au/pdf2/vici09_32.pdf">http://www.chromtech.net.au/pdf2/vici09_32.pdf</a>
ZT1M	300 Series SS 1/16inx0.50mm bore, Tee	Valco Fittings-Tee	\$36.00	<a href="http://www.chromtech.net.au/pdf2/vici09_32.pdf">http://www.chromtech.net.au/pdf2/vici09_32.pdf</a>
ZT.5	300 Series SS 1/32inx.25mm bore, Tee	Valco Fittings-Tee	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_32.pdf">http://www.chromtech.net.au/pdf2/vici09_32.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

ZT1C	300 Series SS 1/16x0.25mm bore, Tee	Valco Fittings-Tee	\$40.00	<a href="http://www.chromtech.net.au/pdf2/vici09_32.pdf">http://www.chromtech.net.au/pdf2/vici09_32.pdf</a>
ZT1	300 Series SS 1/16inx0.75mm bore, Tee	Valco Fittings-Tee	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_32.pdf">http://www.chromtech.net.au/pdf2/vici09_32.pdf</a>
ZT1L	300 Series SS 1/16inx1.00mm bore, Tee	Valco Fittings-Tee	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_32.pdf">http://www.chromtech.net.au/pdf2/vici09_32.pdf</a>
ZT2	300 Series SS 1/8inx0.75mm bore, Tee	Valco Fittings-Tee	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_32.pdf">http://www.chromtech.net.au/pdf2/vici09_32.pdf</a>
ZT2L	300 Series SS 1/8inx2.00mm bore, Tee	Valco Fittings-Tee	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_32.pdf">http://www.chromtech.net.au/pdf2/vici09_32.pdf</a>
ZT4	300 Series SS 1/4inx1.00mm bore, Tee	Valco Fittings-Tee	\$36.00	<a href="http://www.chromtech.net.au/pdf2/vici09_32.pdf">http://www.chromtech.net.au/pdf2/vici09_32.pdf</a>
ZT4L	300 Series SS 1/4inx4.60mm bore, Tee	Valco Fittings-Tee	\$36.00	<a href="http://www.chromtech.net.au/pdf2/vici09_32.pdf">http://www.chromtech.net.au/pdf2/vici09_32.pdf</a>
ZT.5L	300 Series SS 1/32inx.50mm bore, Tee	Valco Fittings-Tee	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_32.pdf">http://www.chromtech.net.au/pdf2/vici09_32.pdf</a>
ZbrU42	1/4in to 1/8in x1.0mm bore -Bulkhead	Valco Fittings-Red Unions	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_31.pdf">http://www.chromtech.net.au/pdf2/vici09_31.pdf</a>
EZrU41t	1/4in to 1/16in x1/16in bore -standard	Valco Fittings-Red Unions	\$27.00	<a href="http://www.chromtech.net.au/pdf2/vici09_31.pdf">http://www.chromtech.net.au/pdf2/vici09_31.pdf</a>
EZbrU1.5l	1/16in to 1/32in x0.50mm bore -bulkhead	Valco Fittings-Red Unions	\$51.00	<a href="http://www.chromtech.net.au/pdf2/vici09_31.pdf">http://www.chromtech.net.au/pdf2/vici09_31.pdf</a>
EZbrU42l	1/4in to 1/8in x2.0mm bore -bulkhead	Valco Fittings-Red Unions	\$36.00	<a href="http://www.chromtech.net.au/pdf2/vici09_31.pdf">http://www.chromtech.net.au/pdf2/vici09_31.pdf</a>
EZbrU42t	1/4in to 1/8in x1/8in bore -bulkhead	Valco Fittings-Red Unions	\$36.00	<a href="http://www.chromtech.net.au/pdf2/vici09_31.pdf">http://www.chromtech.net.au/pdf2/vici09_31.pdf</a>
ZbrU21t	1/8in to 1/16in x1/16in -Bulkhead	Valco Fittings-Red Unions	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_31.pdf">http://www.chromtech.net.au/pdf2/vici09_31.pdf</a>
ZbrU41C	1/4in to 1/16in x0.25mm bore -Bulkhead	Valco Fittings-Red Unions	\$36.00	<a href="http://www.chromtech.net.au/pdf2/vici09_31.pdf">http://www.chromtech.net.au/pdf2/vici09_31.pdf</a>
EZbrU41	1/4in to 1/16in x0.75mm bore -bulkhead	Valco Fittings-Red Unions	\$33.00	<a href="http://www.chromtech.net.au/pdf2/vici09_31.pdf">http://www.chromtech.net.au/pdf2/vici09_31.pdf</a>
ZbrU41t	1/4in to 1/16in x1/16in -Bulkhead	Valco Fittings-Red Unions	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_31.pdf">http://www.chromtech.net.au/pdf2/vici09_31.pdf</a>
EZbrU41t	1/4in to 1/16in x1/16in bore -bulkhead	Valco Fittings-Red Unions	\$33.00	<a href="http://www.chromtech.net.au/pdf2/vici09_31.pdf">http://www.chromtech.net.au/pdf2/vici09_31.pdf</a>
EZrU.5l	1/16in to 1/32in x0.50mm bore -Standard	Valco Fittings-Red Unions	\$38.00	<a href="http://www.chromtech.net.au/pdf2/vici09_31.pdf">http://www.chromtech.net.au/pdf2/vici09_31.pdf</a>
EZbrU.5l	1/16in to 1/32in x0.25mm bore -Bulkhead	Valco Fittings-Red Unions	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_31.pdf">http://www.chromtech.net.au/pdf2/vici09_31.pdf</a>
EZrU.5l	1/16in to 1/32in x0.50mm bore -Standard	Valco Fittings-Red Unions	\$38.00	<a href="http://www.chromtech.net.au/pdf2/vici09_31.pdf">http://www.chromtech.net.au/pdf2/vici09_31.pdf</a>
EZbrU.5l	1/16in to 1/32in x0.50mm bore -Bulkhead	Valco Fittings-Red Unions	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_31.pdf">http://www.chromtech.net.au/pdf2/vici09_31.pdf</a>



## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

EZrU.51t	1/16in to 1/32in x1/32in bore -Standard	Valco Fittings-Red Unions	\$38.00	<a href="http://www.chromtech.net.au/pdf2/vici09_31.pdf">http://www.chromtech.net.au/pdf2/vici09_31.pdf</a>
EZbrU.51t	1/16in to 1/32in x1/32in bore -Bulkhead	Valco Fittings-Red Unions	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_31.pdf">http://www.chromtech.net.au/pdf2/vici09_31.pdf</a>
EZbrU42	1/4in to 1/8in x1.0mm bore -bulkhead	Valco Fittings-Red Unions	\$36.00	<a href="http://www.chromtech.net.au/pdf2/vici09_31.pdf">http://www.chromtech.net.au/pdf2/vici09_31.pdf</a>
EZrU42t	1/4in to 1/8in x1/8in bore -standard	Valco Fittings-Red Unions	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_31.pdf">http://www.chromtech.net.au/pdf2/vici09_31.pdf</a>
ZbrU41	1/4in to 1/16in x0.75mm bore -Bulkhead	Valco Fittings-Red Unions	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_31.pdf">http://www.chromtech.net.au/pdf2/vici09_31.pdf</a>
EZbrU2.5l	1/8in to 1/32in x0.50mm bore -bulkhead	Valco Fittings-Red Unions	\$51.00	<a href="http://www.chromtech.net.au/pdf2/vici09_31.pdf">http://www.chromtech.net.au/pdf2/vici09_31.pdf</a>
EZrU1.5l	1/16in to 1/32in x0.50mm bore -standard	Valco Fittings-Red Unions	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_31.pdf">http://www.chromtech.net.au/pdf2/vici09_31.pdf</a>
EZrU1.5	1/16in to 1/32 x0.25mm bore -standard	Valco Fittings-Red Unions	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_31.pdf">http://www.chromtech.net.au/pdf2/vici09_31.pdf</a>
EZrU1.5t	1/16in to 1/32in x1/32in bore -standard	Valco Fittings-Red Unions	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_31.pdf">http://www.chromtech.net.au/pdf2/vici09_31.pdf</a>
EZbrU1.5t	1/16in to 1/32in x1/32in bore -bulkhead	Valco Fittings-Red Unions	\$51.00	<a href="http://www.chromtech.net.au/pdf2/vici09_31.pdf">http://www.chromtech.net.au/pdf2/vici09_31.pdf</a>
EZrU42l	1/4in to 1/8in x2.0mm bore -standard	Valco Fittings-Red Unions	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_31.pdf">http://www.chromtech.net.au/pdf2/vici09_31.pdf</a>
EZrU2.5l	1/8in to 1/32in x0.50mm bore -standard	Valco Fittings-Red Unions	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_31.pdf">http://www.chromtech.net.au/pdf2/vici09_31.pdf</a>
EZrU41	1/4in to 1/16in x0.75mm bore -standard	Valco Fittings-Red Unions	\$27.00	<a href="http://www.chromtech.net.au/pdf2/vici09_31.pdf">http://www.chromtech.net.au/pdf2/vici09_31.pdf</a>
EZrU2.5t	1/8in to 1/32in x1/32in bore -standard	Valco Fittings-Red Unions	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_31.pdf">http://www.chromtech.net.au/pdf2/vici09_31.pdf</a>
EZbrU2.5t	1/8in to 1/32in x1/32in bore -bulkhead	Valco Fittings-Red Unions	\$51.00	<a href="http://www.chromtech.net.au/pdf2/vici09_31.pdf">http://www.chromtech.net.au/pdf2/vici09_31.pdf</a>
EZbrU21t	1/8in to 1/16in x1/16in bore -bulkhead	Valco Fittings-Red Unions	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_31.pdf">http://www.chromtech.net.au/pdf2/vici09_31.pdf</a>
EZrU2.5	1/8in to 1/32in x0.25mm bore -bulkhead	Valco Fittings-Red Unions	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_31.pdf">http://www.chromtech.net.au/pdf2/vici09_31.pdf</a>
EZrU41C	1/4in to 1/16in x0.25mm bore -standard	Valco Fittings-Red Unions	\$36.00	<a href="http://www.chromtech.net.au/pdf2/vici09_31.pdf">http://www.chromtech.net.au/pdf2/vici09_31.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

EZrU21C	1/8in to 1/16in x0.25mm bore -standard	Valco Fittings-Red Unions	\$33.00	<a href="http://www.chromtech.net.au/pdf2/vici09_31.pdf">http://www.chromtech.net.au/pdf2/vici09_31.pdf</a>
EZrU21t	1/8in to 1/16in x1/16in bore -standard	Valco Fittings-Red Unions	\$26.00	<a href="http://www.chromtech.net.au/pdf2/vici09_31.pdf">http://www.chromtech.net.au/pdf2/vici09_31.pdf</a>
EZbrU21	1/8in to 1/16in x0.75mm bore -bulkhead	Valco Fittings-Red Unions	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_31.pdf">http://www.chromtech.net.au/pdf2/vici09_31.pdf</a>
EZrU21	1/8in to 1/16in x0.75mm bore -standard	Valco Fittings-Red Unions	\$26.00	<a href="http://www.chromtech.net.au/pdf2/vici09_31.pdf">http://www.chromtech.net.au/pdf2/vici09_31.pdf</a>
EZrU42	n1/4in to 1/8in x1.0mm bore -standard	Valco Fittings-Red Unions	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_31.pdf">http://www.chromtech.net.au/pdf2/vici09_31.pdf</a>
ErU21	1/8in to 1/16in x0.75mm bore -standard	Valco Fittings-Red Unions	\$27.00	<a href="http://www.chromtech.net.au/pdf2/vici09_30.pdf">http://www.chromtech.net.au/pdf2/vici09_30.pdf</a>
EbrU12l	1/8in to 1/16in x1.00mm bore -bulkhead	Valco Fittings-Red Unions	\$33.00	<a href="http://www.chromtech.net.au/pdf2/vici09_30.pdf">http://www.chromtech.net.au/pdf2/vici09_30.pdf</a>
EbrU24l	1/4in to 1/8in x2.0mm bore -bulkhead	Valco Fittings-Red Unions	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_30.pdf">http://www.chromtech.net.au/pdf2/vici09_30.pdf</a>
EbrU14t	1/4in to 1/16in x1/16in bore -bulkhead	Valco Fittings-Red Unions	\$38.00	<a href="http://www.chromtech.net.au/pdf2/vici09_30.pdf">http://www.chromtech.net.au/pdf2/vici09_30.pdf</a>
EbrU14l	1/4in to 1/16in x1.00mm bore -bulkhead	Valco Fittings-Red Unions	\$38.00	<a href="http://www.chromtech.net.au/pdf2/vici09_30.pdf">http://www.chromtech.net.au/pdf2/vici09_30.pdf</a>
EbrU12t	1/8in to 1/16in x1/16in bore-bulkhead	Valco Fittings-Red Unions	\$33.00	<a href="http://www.chromtech.net.au/pdf2/vici09_30.pdf">http://www.chromtech.net.au/pdf2/vici09_30.pdf</a>
ErU42t	1/4in to 1/8in x1/8in bore -standard	Valco Fittings-Red Unions	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_30.pdf">http://www.chromtech.net.au/pdf2/vici09_30.pdf</a>
ErU42l	1/4in to 1/8in x2.0mm bore -standard	Valco Fittings-Red Unions	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_30.pdf">http://www.chromtech.net.au/pdf2/vici09_30.pdf</a>
ErU41t	1/4in to 1/16in x1/16in bore -standard	Valco Fittings-Red Unions	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_30.pdf">http://www.chromtech.net.au/pdf2/vici09_30.pdf</a>
ZbrU21	1/8in to 1/16in x0.75mm bore -Bulkhead	Valco Fittings-Red Unions	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_30.pdf">http://www.chromtech.net.au/pdf2/vici09_30.pdf</a>
ErU21l	1/8in to 1/16in x1.00mm bore -standard	Valco Fittings-Red Unions	\$27.00	<a href="http://www.chromtech.net.au/pdf2/vici09_30.pdf">http://www.chromtech.net.au/pdf2/vici09_30.pdf</a>
ErU42	1/4in to 1/8in x1.0mm bore -standard	Valco Fittings-Red Unions	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_30.pdf">http://www.chromtech.net.au/pdf2/vici09_30.pdf</a>
ErU21t	1/8in to 1/16in x1/16in bore-standard	Valco Fittings-Red Unions	\$27.00	<a href="http://www.chromtech.net.au/pdf2/vici09_30.pdf">http://www.chromtech.net.au/pdf2/vici09_30.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

ErU41	1/4in to 1/16in x0.75mm bore -standard	Valco Fittings-Red Unions	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_30.pdf">http://www.chromtech.net.au/pdf2/vici09_30.pdf</a>
ZrU1.5	1/16in to 1/32in x0.25mm bore -Standard	Valco Fittings-Red Unions	\$33.00	<a href="http://www.chromtech.net.au/pdf2/vici09_29.pdf">http://www.chromtech.net.au/pdf2/vici09_29.pdf</a>
ZrU41t	1/4in to 1/16in x1/16in -Standard	Valco Fittings-Red Unions	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_29.pdf">http://www.chromtech.net.au/pdf2/vici09_29.pdf</a>
ZrU41	1/4in to 1/16in x0.75mm bore -Standard	Valco Fittings-Red Unions	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_29.pdf">http://www.chromtech.net.au/pdf2/vici09_29.pdf</a>
ZrU41C	1/4in to 1/16in x0.25mm bore -Standard	Valco Fittings-Red Unions	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_29.pdf">http://www.chromtech.net.au/pdf2/vici09_29.pdf</a>
ZrU21t	1/8in to 1/16in x1/16in -Standard	Valco Fittings-Red Unions	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_29.pdf">http://www.chromtech.net.au/pdf2/vici09_29.pdf</a>
ZrU21	1/8in to 1/16in x0.75mm bore -Standard	Valco Fittings-Red Unions	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_29.pdf">http://www.chromtech.net.au/pdf2/vici09_29.pdf</a>
ZrU21C	1/8in to 1/16in x0.25mm bore -Standard	Valco Fittings-Red Unions	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_29.pdf">http://www.chromtech.net.au/pdf2/vici09_29.pdf</a>
ZrU2.5t	1/8in to 1/32in x1/32in -Standard	Valco Fittings-Red Unions	\$38.00	<a href="http://www.chromtech.net.au/pdf2/vici09_29.pdf">http://www.chromtech.net.au/pdf2/vici09_29.pdf</a>
ZrU1.5I	Internal Reducing Union -SS 1/16 to 1/32 0.50mm Bore	Valco Fittings-Red Unions	\$33.00	<a href="http://www.chromtech.net.au/pdf2/vici09_29.pdf">http://www.chromtech.net.au/pdf2/vici09_29.pdf</a>
ZrU1.5XC	1/16in to 1/32in x0.15mm bore	Valco Fittings-Red Unions	\$48.00	<a href="http://www.chromtech.net.au/pdf2/vici09_29.pdf">http://www.chromtech.net.au/pdf2/vici09_29.pdf</a>
ZrU2.5I	1/8in to 1/32in x0.50mm bore -Standard	Valco Fittings-Red Unions	\$38.00	<a href="http://www.chromtech.net.au/pdf2/vici09_29.pdf">http://www.chromtech.net.au/pdf2/vici09_29.pdf</a>
ZrU42	1/4in to 1/8in x1.0mm bore -Standard	Valco Fittings-Red Unions	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_29.pdf">http://www.chromtech.net.au/pdf2/vici09_29.pdf</a>
ZrU2.5	1/8in to 1/32in x.25mm bore -Standard	Valco Fittings-Red Unions	\$38.00	<a href="http://www.chromtech.net.au/pdf2/vici09_29.pdf">http://www.chromtech.net.au/pdf2/vici09_29.pdf</a>
ZbrU2.5I	1/8in to 1/32in x0.50mm bore -Bulkhead	Valco Fittings-Red Unions	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_29.pdf">http://www.chromtech.net.au/pdf2/vici09_29.pdf</a>
ZrU1.5t	1/16in to 1/32in x1/32in -Standard	Valco Fittings-Red Unions	\$33.00	<a href="http://www.chromtech.net.au/pdf2/vici09_29.pdf">http://www.chromtech.net.au/pdf2/vici09_29.pdf</a>
ZbrU2.5t	1/8in to 1/32in x1/32in -Bulkhead	Valco Fittings-Red Unions	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_29.pdf">http://www.chromtech.net.au/pdf2/vici09_29.pdf</a>
ZbrU2.5	Bulkhead Int Red 1/8 to 1/32in 0.25mm bore	Valco Fittings-Red Unions	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_29.pdf">http://www.chromtech.net.au/pdf2/vici09_29.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

ZbrU1.5t	1/16in to 1/32in x1/32in -Bulkhead	Valco Fittings-Red Unions	\$40.00	<a href="http://www.chromtech.net.au/pdf2/vici09_29.pdf">http://www.chromtech.net.au/pdf2/vici09_29.pdf</a>
ZbrU1.5l	1/16in to 1/32in x0.50mm bore -Bulkhead	Valco Fittings-Red Unions	\$40.00	<a href="http://www.chromtech.net.au/pdf2/vici09_29.pdf">http://www.chromtech.net.au/pdf2/vici09_29.pdf</a>
ZbrU1.5	1/16in to 1/32in x0.25mm bore -Bulkhead	Valco Fittings-Red Unions	\$40.00	<a href="http://www.chromtech.net.au/pdf2/vici09_29.pdf">http://www.chromtech.net.au/pdf2/vici09_29.pdf</a>
ZrU42t	1/4in to 1/8in x1/8in -Standard	Valco Fittings-Red Unions	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_29.pdf">http://www.chromtech.net.au/pdf2/vici09_29.pdf</a>
ZrU42l	1/4in to 1/8in x2.0mm bore -Standard	Valco Fittings-Red Unions	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_29.pdf">http://www.chromtech.net.au/pdf2/vici09_29.pdf</a>
ZbrU21C	1/8in to 1/16in x0.25mm bore -Bulkhead	Valco Fittings-Red Unions	\$36.00	<a href="http://www.chromtech.net.au/pdf2/vici09_29.pdf">http://www.chromtech.net.au/pdf2/vici09_29.pdf</a>
EZU2	1/8in x1.0mm bore -Standard	Valco Fittings-Unions	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_27.pdf">http://www.chromtech.net.au/pdf2/vici09_27.pdf</a>
EZU.5l	1/32in x0.50mm bore	Valco Fittings-Unions	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_27.pdf">http://www.chromtech.net.au/pdf2/vici09_27.pdf</a>
EZU1C	1/16in x0.25mm bore	Valco Fittings-Unions	\$36.00	<a href="http://www.chromtech.net.au/pdf2/vici09_27.pdf">http://www.chromtech.net.au/pdf2/vici09_27.pdf</a>
EZbU1C	1/16in x0.25mm bore	Valco Fittings-Unions	\$42.00	<a href="http://www.chromtech.net.au/pdf2/vici09_27.pdf">http://www.chromtech.net.au/pdf2/vici09_27.pdf</a>
EZU1M	1/16in x0.50mm bore	Valco Fittings-Unions	\$36.00	<a href="http://www.chromtech.net.au/pdf2/vici09_27.pdf">http://www.chromtech.net.au/pdf2/vici09_27.pdf</a>
EZU1	1/32in x0.75mm bore	Valco Fittings-Unions	\$27.00	<a href="http://www.chromtech.net.au/pdf2/vici09_27.pdf">http://www.chromtech.net.au/pdf2/vici09_27.pdf</a>
EZU.5	1/32in x0.25mm bore	Valco Fittings-Unions	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_27.pdf">http://www.chromtech.net.au/pdf2/vici09_27.pdf</a>
EZbU1t	1/16in x1/16in bore	Valco Fittings-Unions	\$33.00	<a href="http://www.chromtech.net.au/pdf2/vici09_27.pdf">http://www.chromtech.net.au/pdf2/vici09_27.pdf</a>
EZbU1M	1/16in x0.50mm bore	Valco Fittings-Unions	\$42.00	<a href="http://www.chromtech.net.au/pdf2/vici09_27.pdf">http://www.chromtech.net.au/pdf2/vici09_27.pdf</a>
EZbU2	1/8in x1.0mm bore -Bulkhead	Valco Fittings-Unions	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_27.pdf">http://www.chromtech.net.au/pdf2/vici09_27.pdf</a>
EZU2l	1/8in x2.0mm bore -Standard	Valco Fittings-Unions	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_27.pdf">http://www.chromtech.net.au/pdf2/vici09_27.pdf</a>
EZbU2l	1/8in x2.0mm bore -Bulkhead	Valco Fittings-Unions	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_27.pdf">http://www.chromtech.net.au/pdf2/vici09_27.pdf</a>
EZU2t	Through-type Union 1/8in external to internal union	Valco Fittings-Unions	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_27.pdf">http://www.chromtech.net.au/pdf2/vici09_27.pdf</a>
EZbU2t	1/8in x1/8in bore -Bulkhead	Valco Fittings-Unions	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_27.pdf">http://www.chromtech.net.au/pdf2/vici09_27.pdf</a>
EZU1t	1/16in x1/16in bore	Valco Fittings-Unions	\$27.00	<a href="http://www.chromtech.net.au/pdf2/vici09_27.pdf">http://www.chromtech.net.au/pdf2/vici09_27.pdf</a>
EbU2l	1/8in x2.0mm bore -Bulkhead	Valco Fittings-Unions	\$27.00	<a href="http://www.chromtech.net.au/pdf2/vici09_27.pdf">http://www.chromtech.net.au/pdf2/vici09_27.pdf</a>
EZbU1	1/16in x0.75mm bore	Valco Fittings-Unions	\$33.00	<a href="http://www.chromtech.net.au/pdf2/vici09_27.pdf">http://www.chromtech.net.au/pdf2/vici09_27.pdf</a>
EU2l	1/8in x2.0mm bore -Standard	Valco Fittings-Unions	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_27.pdf">http://www.chromtech.net.au/pdf2/vici09_27.pdf</a>
EU.5t	1/32in Ultra low mass external union -Bore 1/32in	Valco Fittings-Unions	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_27.pdf">http://www.chromtech.net.au/pdf2/vici09_27.pdf</a>
EU2t	Through-type Union 1/8in external union	Valco Fittings-Unions	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_27.pdf">http://www.chromtech.net.au/pdf2/vici09_27.pdf</a>
EbU2t	1/8in x1/8in bore -Bulkhead	Valco Fittings-Unions	\$27.00	<a href="http://www.chromtech.net.au/pdf2/vici09_27.pdf">http://www.chromtech.net.au/pdf2/vici09_27.pdf</a>
EU4	1/4in x2.0mm bore -Standard	Valco Fittings-Unions	\$22.00	<a href="http://www.chromtech.net.au/pdf2/vici09_27.pdf">http://www.chromtech.net.au/pdf2/vici09_27.pdf</a>
EbU4	1/4in x2.0mm bore -Bulkhead	Valco Fittings-Unions	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_27.pdf">http://www.chromtech.net.au/pdf2/vici09_27.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

EbU4I	1/4in x4.6mm bore -Bulkhead	Valco Fittings-Unions	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_27.pdf">http://www.chromtech.net.au/pdf2/vici09_27.pdf</a>
EU4t	Through-type Union 1/4in external union	Valco Fittings-Unions	\$22.00	<a href="http://www.chromtech.net.au/pdf2/vici09_27.pdf">http://www.chromtech.net.au/pdf2/vici09_27.pdf</a>
EbU4t	1/4in x1/4in bore -Bulkhead	Valco Fittings-Unions	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_27.pdf">http://www.chromtech.net.au/pdf2/vici09_27.pdf</a>
EU.5	1/32in Ultra low mass external union -Bore 0.25mm	Valco Fittings-Unions	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_27.pdf">http://www.chromtech.net.au/pdf2/vici09_27.pdf</a>
EU.5I	1/32in Ultra low mass external union -Bore 0.50mm	Valco Fittings-Unions	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_27.pdf">http://www.chromtech.net.au/pdf2/vici09_27.pdf</a>
EU4I	1/4in x4.6mm bore -Standard	Valco Fittings-Unions	\$22.00	<a href="http://www.chromtech.net.au/pdf2/vici09_27.pdf">http://www.chromtech.net.au/pdf2/vici09_27.pdf</a>
EU2	1/8in x1.0mm bore -Standard	Valco Fittings-Unions	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_27.pdf">http://www.chromtech.net.au/pdf2/vici09_27.pdf</a>
ZU.5	SS-303 1/32in x0.25mm bore -standard	Valco Fittings-Unions	\$33.00	<a href="http://www.chromtech.net.au/pdf2/vici09_26.pdf">http://www.chromtech.net.au/pdf2/vici09_26.pdf</a>
ZU.5XC	SS-303 1/32in x0.15mm bore -standard	Valco Fittings-Unions	\$40.00	<a href="http://www.chromtech.net.au/pdf2/vici09_26.pdf">http://www.chromtech.net.au/pdf2/vici09_26.pdf</a>
ZU.5I	SS-303 1/32in x0.50mm bore -standard	Valco Fittings-Unions	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_26.pdf">http://www.chromtech.net.au/pdf2/vici09_26.pdf</a>
ZU.5t	SS-303 1/32in x0.50mm bore -standard	Valco Fittings-Unions	\$33.00	<a href="http://www.chromtech.net.au/pdf2/vici09_26.pdf">http://www.chromtech.net.au/pdf2/vici09_26.pdf</a>
ZU1XC	SS-303 1/16in x0.15mm bore -standard	Valco Fittings-Unions	\$40.00	<a href="http://www.chromtech.net.au/pdf2/vici09_26.pdf">http://www.chromtech.net.au/pdf2/vici09_26.pdf</a>
ZU1C	SS-303 1/16in x0.25mm bore -standard	Valco Fittings-Unions	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_26.pdf">http://www.chromtech.net.au/pdf2/vici09_26.pdf</a>
ZU1M	Union Fitting SS-303 1/16" x 0.50mm bore - Standard	Valco Fittings-Unions	\$26.00	<a href="http://www.chromtech.net.au/pdf2/vici09_26.pdf">http://www.chromtech.net.au/pdf2/vici09_26.pdf</a>
ZU1	SS-303 1/16in x0.75mm bore -standard	Valco Fittings-Unions	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_26.pdf">http://www.chromtech.net.au/pdf2/vici09_26.pdf</a>
ZU1I	SS-303 1/16in x1.0mm bore -standard	Valco Fittings-Unions	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_26.pdf">http://www.chromtech.net.au/pdf2/vici09_26.pdf</a>
ZU2	SS-303 1/8in x0.75mm bore -standard	Valco Fittings-Unions	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_26.pdf">http://www.chromtech.net.au/pdf2/vici09_26.pdf</a>
ZU2I	SS-303 1/8in x2.0mm bore -standard	Valco Fittings-Unions	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_26.pdf">http://www.chromtech.net.au/pdf2/vici09_26.pdf</a>
ZU2t	Through-type Union 1/8in internal union	Valco Fittings-Unions	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_26.pdf">http://www.chromtech.net.au/pdf2/vici09_26.pdf</a>
ZbU.5XC	SS-303 1/32in x0.15mm bore -Bulkhead	Valco Fittings-Unions	\$48.00	<a href="http://www.chromtech.net.au/pdf2/vici09_26.pdf">http://www.chromtech.net.au/pdf2/vici09_26.pdf</a>
ZbU1XC	SS-303 1/16in x0.15in -bulkhead	Valco Fittings-Unions	\$48.00	<a href="http://www.chromtech.net.au/pdf2/vici09_26.pdf">http://www.chromtech.net.au/pdf2/vici09_26.pdf</a>
ZbU1I	SS-303 1/16in x1.0mm bore -Bulkhead	Valco Fittings-Unions	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_26.pdf">http://www.chromtech.net.au/pdf2/vici09_26.pdf</a>
ZU4I	SS-303 1/4in x4.6mm bore -Bulkhead	Valco Fittings-Unions	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_26.pdf">http://www.chromtech.net.au/pdf2/vici09_26.pdf</a>
ZU4t	Through-type Union 1/4in internal union	Valco Fittings-Unions	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_26.pdf">http://www.chromtech.net.au/pdf2/vici09_26.pdf</a>
ZU4	SS-303 1/4in x0.75mm bore -standard	Valco Fittings-Unions	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_26.pdf">http://www.chromtech.net.au/pdf2/vici09_26.pdf</a>
ZbU.5	SS-303 1/32in x0.25mm bore -Bulkhead	Valco Fittings-Unions	\$40.00	<a href="http://www.chromtech.net.au/pdf2/vici09_26.pdf">http://www.chromtech.net.au/pdf2/vici09_26.pdf</a>
ZbU.5t	SS-303 1/32in x0.50mm bore -Bulkhead	Valco Fittings-Unions	\$40.00	<a href="http://www.chromtech.net.au/pdf2/vici09_26.pdf">http://www.chromtech.net.au/pdf2/vici09_26.pdf</a>
ZbU.5I	SS-303 1/32in x0.50mm bore -Bulkhead	Valco Fittings-Unions	\$40.00	<a href="http://www.chromtech.net.au/pdf2/vici09_26.pdf">http://www.chromtech.net.au/pdf2/vici09_26.pdf</a>
ZbU1	SS-303 1/16in x0.75mm bore -Bulkhead	Valco Fittings-Unions	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_26.pdf">http://www.chromtech.net.au/pdf2/vici09_26.pdf</a>
ZbU4t	SS-303 1/4in x1/4in bore -Bulkhead	Valco Fittings-Unions	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_26.pdf">http://www.chromtech.net.au/pdf2/vici09_26.pdf</a>
ZbU1t	SS-303 1/16in x1/16in -bulkhead	Valco Fittings-Unions	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_26.pdf">http://www.chromtech.net.au/pdf2/vici09_26.pdf</a>
ZbU2	SS-303 1/8in x0.75mm bore -Bulkhead	Valco Fittings-Unions	\$27.00	<a href="http://www.chromtech.net.au/pdf2/vici09_26.pdf">http://www.chromtech.net.au/pdf2/vici09_26.pdf</a>
ZbU2I	SS-303 1/8in x2.0mm bore -Bulkhead	Valco Fittings-Unions	\$27.00	<a href="http://www.chromtech.net.au/pdf2/vici09_26.pdf">http://www.chromtech.net.au/pdf2/vici09_26.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

ZbU2t	SS-303 1/8in x1/8in bore -Bulkhead	Valco Fittings-Unions	\$27.00	<a href="http://www.chromtech.net.au/pdf2/vici09_26.pdf">http://www.chromtech.net.au/pdf2/vici09_26.pdf</a>
ZbU4	SS-303 1/4in x0.75mm bore -Bulkhead	Valco Fittings-Unions	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_26.pdf">http://www.chromtech.net.au/pdf2/vici09_26.pdf</a>
ZbU4I	SS-303 1/4in x4.6mm bore -standard	Valco Fittings-Unions	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_26.pdf">http://www.chromtech.net.au/pdf2/vici09_26.pdf</a>
ZbU1C	SS-303 1/16in x0.25mm bore -Bulkhead	Valco Fittings-Unions	\$36.00	<a href="http://www.chromtech.net.au/pdf2/vici09_26.pdf">http://www.chromtech.net.au/pdf2/vici09_26.pdf</a>
ZbU1M	SS-303 1/16in x0.50mm bore -Bulkhead	Valco Fittings-Unions	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_26.pdf">http://www.chromtech.net.au/pdf2/vici09_26.pdf</a>
MU.5CS6	Stainless steel Union -1/32inx0.25mm	Valco Fittings-MicroVolume	\$56.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
Mx1CPK	PEEK Cross -1/16inx0.25mm	Valco Fittings-MicroVolume	\$76.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
Mx1CKF	CTFE Cross -1/16inx0.25mm	Valco Fittings-MicroVolume	\$76.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MRx.5S6	1/32in connectors -SS ring for union tee or cross	Valco Fittings-MicroVolume	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MRx1S6	1/16in connectors -PEEK ring for union tee or cross	Valco Fittings-MicroVolume	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MRY.5S6	1/32in connectors -SS ring for Y-1piece	Valco Fittings-MicroVolume	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MRY1S6	1/16in connectors -SS ring for Y-1piece	Valco Fittings-MicroVolume	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MRx.5PK	1/32in connectors -PEEK ring for union tee or cross	Valco Fittings-MicroVolume	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MRx1PK	1/16in connectors -PEEK ring for union tee or cross	Valco Fittings-MicroVolume	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MRY.5PK	1/32in connectors -PEEK ring for Y-1piece	Valco Fittings-MicroVolume	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MRY1PK	1/16in connectors -PEEK ring for Y-1piece	Valco Fittings-MicroVolume	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MU.5CHC	Hastelloy C Union -1/32inx0.25mm	Valco Fittings-MicroVolume	\$86.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
ZN1	1/16in connectors -Nuts for SS ring	Valco Fittings-MicroVolume	\$2.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
Mx.5xCKF	CTFE Cross -1/32inx0.15mm	Valco Fittings-MicroVolume	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
Mx.5xCPK	PEEK Cross -1/32inx0.15mm	Valco Fittings-MicroVolume	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
Mx.5xCTI	Titanium Cross -1/32inx0.15mm	Valco Fittings-MicroVolume	\$141.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

Mx.5xCHC	Hastelloy C Cross -1/32inx0.15mm	Valco Fittings-MicroVolume	\$134.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MU.5CTI	MU.5CTI Union -1/32inx0.25mm	Valco Fittings-MicroVolume	\$93.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
Mx1CTI	Titanium Cross -1/16inx0.25mm	Valco Fittings-MicroVolume	\$108.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
Mx.5xCS6	Stainless steel Cross -1/32inx0.15mm	Valco Fittings-MicroVolume	\$99.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MT1CPK	PEEK Tee -1/16inx0.25mm	Valco Fittings-MicroVolume	\$56.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MY.5xCKF	CTFE Y-1piece -1/32inx0.15mm bore	Valco Fittings-MicroVolume	\$76.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MU.5CPK	PEEK Union -1/32inx0.25mm	Valco Fittings-MicroVolume	\$56.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MY.5xCPK	PEEK Y-1piece -1/32inx0.15mm bore	Valco Fittings-MicroVolume	\$76.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
ZN.5	1/32in connectors -Nuts for SS ring	Valco Fittings-MicroVolume	\$4.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MY.5xCHC	Hastelloy C Y-1piece -1/32inx0.15mm bore	Valco Fittings-MicroVolume	\$114.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MU.5xCS6	Stainless steel Tee -1/32inx0.15mm bore	Valco Fittings-MicroVolume	\$70.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MT1CS6	Stainless steel Tee -1/16inx0.25mm	Valco Fittings-MicroVolume	\$56.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MU.5xCHC	Hastelloy C Union -1/32inx0.15mm bore	Valco Fittings-MicroVolume	\$99.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MU.5xCTI	Titanium Union 1/32inx0.15mm bore	Valco Fittings-MicroVolume	\$108.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MU.5xCPK	PEEK Union -1/32inx0.15mm bore	Valco Fittings-MicroVolume	\$62.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MU.5xCKF	CTFE Union -1/32inx0.15mm bore	Valco Fittings-MicroVolume	\$62.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MT.5xCS6	Stainless steel Tee -1/32inx0.15mm bore	Valco Fittings-MicroVolume	\$87.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MT.5xCHC	Hastelloy C Tee -1/32inx0.15mm bore	Valco Fittings-MicroVolume	\$114.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

MT.5xCTI	Titanium Tee -1/32inx0.15mm bore	Valco Fittings-MicroVolume	\$122.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MT.5xCPK	PEEK Tee -1/32inx0.15mm bore	Valco Fittings-MicroVolume	\$76.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MY1CS6	Stainless steel Y-1piece -1/16inx0.25mm	Valco Fittings-MicroVolume	\$56.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MY.5xCS6	Stainless steel Y-1piece -1/32inx0.15mm bore	Valco Fittings-MicroVolume	\$87.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
Mx1CHC	Hastelloy C Cross -1/16inx0.25mm	Valco Fittings-MicroVolume	\$102.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MY.5xCTI	Titanium Y-1piece -1/32inx0.15mm bore	Valco Fittings-MicroVolume	\$122.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MT1CHC	Hastelloy C Tee -1/16inx0.25mm	Valco Fittings-MicroVolume	\$86.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MT1CTI	Titanium Tee -1/16inx0.25mm	Valco Fittings-MicroVolume	\$93.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MT.5CHC	Hastelloy C Tee -1/32inx0.25mm	Valco Fittings-MicroVolume	\$93.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MT1CKF	CTFE Tee -1/16inx0.25mm	Valco Fittings-MicroVolume	\$56.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MT.5CS6	Stainless steel Tee -1/32inx0.25mm	Valco Fittings-MicroVolume	\$62.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MY1CHC	Hastelloy C Y-1piece -1/16inx0.25mm	Valco Fittings-MicroVolume	\$87.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MY1CTI	Titanium Y-1piece -1/16inx0.25mm	Valco Fittings-MicroVolume	\$93.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MY1CPK	PEEK Y-1piece -1/16inx0.25mm	Valco Fittings-MicroVolume	\$56.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MY1CKF	CTFE Y-1piece -1/16inx0.25mm	Valco Fittings-MicroVolume	\$56.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
Mx1CS6	Stainless steel Cross -1/16inx0.25mm	Valco Fittings-MicroVolume	\$76.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MT.5xCKF	CTFE Tee -1/32inx0.15mm bore	Valco Fittings-MicroVolume	\$76.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MU.5CKF	CTFE Union -1/32inx0.25mm	Valco Fittings-MicroVolume	\$56.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>



## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

MT1xCPK	PEEK Tee -1/16inx0.15mm bore	Valco Fittings-MicroVolume	\$70.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MT1xCKF	CTFE Tee -1/16inx0.15mm bore	Valco Fittings-MicroVolume	\$70.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MY1xCS6	Stainless steel Y-1piece -1/16inx0.15mm bore	Valco Fittings-MicroVolume	\$76.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MY1xCHC	Hastelloy C Y-1piece -1/16inx0.15mm bore	Valco Fittings-MicroVolume	\$108.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MY1xCPK	PEEK Y-1piece -1/16inx0.15mm bore	Valco Fittings-MicroVolume	\$70.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MU1CKF	CTFE Union -1/16inx0.25mm	Valco Fittings-MicroVolume	\$46.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
Mx1xCS6	Stainless steel Cross -1/16inx0.15mm	Valco Fittings-MicroVolume	\$93.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MT1xCTI	Titanium Tee -1/16x0.15mm bore	Valco Fittings-MicroVolume	\$114.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
Mx1xCTI	Titanium Cross -1/16inx0.15mm	Valco Fittings-MicroVolume	\$134.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MY1xCTI	Titanium Y-1piece -1/16inx0.15mm bore	Valco Fittings-MicroVolume	\$114.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
ZN1FPK	1/16in connectors -Nuts for PEEK ring	Valco Fittings-MicroVolume	\$4.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
ZN.5FPK	1/32in connectors -Nuts for PEEK ring	Valco Fittings-MicroVolume	\$4.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
Mx1xCPK	PEEK Cross -1/16inx0.15mm	Valco Fittings-MicroVolume	\$87.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
Mx1xCKF	CTFE Y-1piece -1/16inx0.25mm	Valco Fittings-MicroVolume	\$87.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MU1CS6	Stainless steel Union -1/16inx0.25mm	Valco Fittings-MicroVolume	\$46.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MU1CHC	Hastelloy C Union -1/16inx0.25mm	Valco Fittings-MicroVolume	\$78.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MU1CTI	Titanium Union -1/16inx0.25mm	Valco Fittings-MicroVolume	\$87.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MU1CPK	PEEK Union -1/16inx0.25mm	Valco Fittings-MicroVolume	\$46.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

Mx1xCHC	Hastelloy C Cross -1/16inx0.15mm	Valco Fittings-MicroVolume	\$128.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
Mx.5CTI	Titanium Cross -1/32inx0.25mm	Valco Fittings-MicroVolume	\$114.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MY.5CKF	CTFE Y-1piece -1/32inx0.25mm	Valco Fittings-MicroVolume	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MY.5CPK	PEEK Y-1piece -1/32inx0.25mm	Valco Fittings-MicroVolume	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MY1xCKF	CTFE Y-1piece -1/16inx0.15mm bore	Valco Fittings-MicroVolume	\$70.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MY.5CTI	Titanium Y-1piece -1/32inx0.25mm	Valco Fittings-MicroVolume	\$108.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MT1xCHC	Hastelloy C Tee -1/16x0.15mm bore	Valco Fittings-MicroVolume	\$108.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MY.5CHC	Hastelloy C Y-1piece -1/32inx0.25mm	Valco Fittings-MicroVolume	\$93.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MY.5CS6	Stainless steel Y-1piece -1/32inx0.25mm	Valco Fittings-MicroVolume	\$62.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MT.5CKF	CTFE Tee -1/32inx0.25mm	Valco Fittings-MicroVolume	\$62.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
Mx.5CHC	Hastelloy C Cross -1/32inx0.25mm	Valco Fittings-MicroVolume	\$108.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
Mx.5CS6	Stainless steel Cross -1/32inx0.25mm	Valco Fittings-MicroVolume	\$87.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
Mx.5CPK	PEEK Cross -1/32inx0.25mm	Valco Fittings-MicroVolume	\$82.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
Mx.5CKF	CTFE Cross -1/32inx0.25mm	Valco Fittings-MicroVolume	\$82.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MU1xCS6	Stainless steel Tee -1/16inx0.15mm bore	Valco Fittings-MicroVolume	\$62.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MT.5CTI	Titanium Tee -1/32inx0.25mm	Valco Fittings-MicroVolume	\$108.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MU1xCHC	Hastelloy C Union -1/16inx0.15mm bore	Valco Fittings-MicroVolume	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MU1xCTI	Titanium Union 1/16inx0.15mm bore	Valco Fittings-MicroVolume	\$99.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

MU1xCPK	PEEK Union -1/16inx0.15mm bore	Valco Fittings-MicroVolume	\$56.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MU1xCKF	CTFE Union -1/16inx0.15mm bore	Valco Fittings-MicroVolume	\$56.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MT1xCS6	Stainless steel Tee -1/16inx0.15mm bore	Valco Fittings-MicroVolume	\$76.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
MT.5CPK	PEEK Tee -1/32inx0.25mm	Valco Fittings-MicroVolume	\$62.00	<a href="http://www.chromtech.net.au/pdf2/vici09_23.pdf">http://www.chromtech.net.au/pdf2/vici09_23.pdf</a>
HPM	Miniature Helium Purifier	Instrumentation	\$702.00	<a href="http://www.chromtech.net.au/pdf2/vici09_227.pdf">http://www.chromtech.net.au/pdf2/vici09_227.pdf</a>
TCD2-NIFED	TCD Thermal conductivity detector, Entire unit (cell and electronics)	Instrumentation	\$1,479.00	<a href="http://www.chromtech.net.au/pdf2/vici09_227.pdf">http://www.chromtech.net.au/pdf2/vici09_227.pdf</a>
TCD2-NIFE-220	TCD Thermal conductivity detector 220VAC, Entire unit (cell and electronics)	Instrumentation	\$4,437.00	<a href="http://www.chromtech.net.au/pdf2/vici09_227.pdf">http://www.chromtech.net.au/pdf2/vici09_227.pdf</a>
TCD2-NIFE	TCD Thermal conductivity detector 220VAC, Entire unit (cell and electronics)	Instrumentation	\$4,437.00	<a href="http://www.chromtech.net.au/pdf2/vici09_227.pdf">http://www.chromtech.net.au/pdf2/vici09_227.pdf</a>
NPM-220	Mini Nitrogen Gas Purifier 220VAC	Instrumentation	\$702.00	<a href="http://www.chromtech.net.au/pdf2/vici09_227.pdf">http://www.chromtech.net.au/pdf2/vici09_227.pdf</a>
NPM	Mini Nitrogen Gas Purifier	Instrumentation	\$702.00	<a href="http://www.chromtech.net.au/pdf2/vici09_227.pdf">http://www.chromtech.net.au/pdf2/vici09_227.pdf</a>
TCD2-C-220	TCD Thermal conductivity detector 220VAC, Controller Only	Instrumentation	\$2,868.00	<a href="http://www.chromtech.net.au/pdf2/vici09_227.pdf">http://www.chromtech.net.au/pdf2/vici09_227.pdf</a>
TCD2-NIFED-220	TCD Thermal conductivity detector, 220VAC, Entire unit (cell and electronics)	Instrumentation	\$1,479.00	<a href="http://www.chromtech.net.au/pdf2/vici09_227.pdf">http://www.chromtech.net.au/pdf2/vici09_227.pdf</a>
TCD2-C	TCD Thermal conductivity detector, Controller Only	Instrumentation	\$2,868.00	<a href="http://www.chromtech.net.au/pdf2/vici09_227.pdf">http://www.chromtech.net.au/pdf2/vici09_227.pdf</a>
HPM-220	Miniature Helium Purifier 220VAC	Instrumentation	\$702.00	<a href="http://www.chromtech.net.au/pdf2/vici09_227.pdf">http://www.chromtech.net.au/pdf2/vici09_227.pdf</a>
NP2-220	Nitrogen purifier 230 VAC	Instrumentation	\$1,167.00	<a href="http://www.chromtech.net.au/pdf2/vici09_226.pdf">http://www.chromtech.net.au/pdf2/vici09_226.pdf</a>
I-23572HP2	Helium purifier, Replacement getter assembly, He	Instrumentation	\$586.00	<a href="http://www.chromtech.net.au/pdf2/vici09_226.pdf">http://www.chromtech.net.au/pdf2/vici09_226.pdf</a>
HP2-220	Helium purifier 220 VAC	Instrumentation	\$1,167.00	<a href="http://www.chromtech.net.au/pdf2/vici09_226.pdf">http://www.chromtech.net.au/pdf2/vici09_226.pdf</a>
NP2	Nitrogen purifier 110 VAC	Instrumentation	\$1,167.00	<a href="http://www.chromtech.net.au/pdf2/vici09_226.pdf">http://www.chromtech.net.au/pdf2/vici09_226.pdf</a>
HP2	Helium purifier 110 VAC	Instrumentation	\$1,167.00	<a href="http://www.chromtech.net.au/pdf2/vici09_226.pdf">http://www.chromtech.net.au/pdf2/vici09_226.pdf</a>
I-23572NP2	Helium purifier, Replacement getter assembly, N2	Instrumentation	\$586.00	<a href="http://www.chromtech.net.au/pdf2/vici09_226.pdf">http://www.chromtech.net.au/pdf2/vici09_226.pdf</a>
D-4-I-SH17-R	PDD Model D-5 for He PID - Shimadzu GC17 uses existing FID Electrometer	Instrumentation	\$4,972.00	<a href="http://www.chromtech.net.au/pdf2/vici09_223.pdf">http://www.chromtech.net.au/pdf2/vici09_223.pdf</a>
D-4-I-SH17-R-220	PDD Model D-5 for He PID - Shimadzu GC17 uses existing FID Electrometer 220VAC	Instrumentation	\$4,972.00	<a href="http://www.chromtech.net.au/pdf2/vici09_223.pdf">http://www.chromtech.net.au/pdf2/vici09_223.pdf</a>
D-4-I-TQ-R-220	PDD Model D-5 for He PID -Therm Trace GC uses existing FID Electrometer 220VAC	Instrumentation	\$4,972.00	<a href="http://www.chromtech.net.au/pdf2/vici09_223.pdf">http://www.chromtech.net.au/pdf2/vici09_223.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

D-4-I-VA38-R-220	PDD Model D-5 for He PID -Varian 3800 GC uses existing FID Electrometer 220VAC	Instrumentation	\$4,972.00	<a href="http://www.chromtech.net.au/pdf2/vici09_223.pdf">http://www.chromtech.net.au/pdf2/vici09_223.pdf</a>
D-4-I-SH14-R	PDD Model D-5 for He PID - Shimadzu GC14 uses existing FID Electrometer	Instrumentation	\$4,972.00	<a href="http://www.chromtech.net.au/pdf2/vici09_223.pdf">http://www.chromtech.net.au/pdf2/vici09_223.pdf</a>
D-4-I	PDD Model D-5 for He PID - all other GCs	Instrumentation	\$6,502.00	<a href="http://www.chromtech.net.au/pdf2/vici09_223.pdf">http://www.chromtech.net.au/pdf2/vici09_223.pdf</a>
D-4-I-SH14-R-220	PDD Model D-5 for He PID - Shimadzu GC14 uses existing FID Electrometer 220VAC	Instrumentation	\$4,972.00	<a href="http://www.chromtech.net.au/pdf2/vici09_223.pdf">http://www.chromtech.net.au/pdf2/vici09_223.pdf</a>
D-4-I-220	PDD Model D-5 for He PID - all other GCs 220VAC	Instrumentation	\$6,502.00	<a href="http://www.chromtech.net.au/pdf2/vici09_223.pdf">http://www.chromtech.net.au/pdf2/vici09_223.pdf</a>
D-4-I-TQ-R	PDD Model D-5 for He PID -Therm Trace GC uses existing FID Electrometer	Instrumentation	\$4,972.00	<a href="http://www.chromtech.net.au/pdf2/vici09_223.pdf">http://www.chromtech.net.au/pdf2/vici09_223.pdf</a>
D-4-I-VA38-R	PDD Model D-5 for He PID -Varian 3800 GC uses existing FID Electrometer	Instrumentation	\$4,972.00	<a href="http://www.chromtech.net.au/pdf2/vici09_223.pdf">http://www.chromtech.net.au/pdf2/vici09_223.pdf</a>
D-3-I-HP	PDD Model D-3 Plug-In system for Agilent 6890	Instrumentation	\$7,917.00	<a href="http://www.chromtech.net.au/pdf2/vici09_222.pdf">http://www.chromtech.net.au/pdf2/vici09_222.pdf</a>
D-4-I-HP58-220	PDD Model D-4 for He PID HP5890 220VAC	Instrumentation	\$6,630.00	<a href="http://www.chromtech.net.au/pdf2/vici09_222.pdf">http://www.chromtech.net.au/pdf2/vici09_222.pdf</a>
D-4-I-HP58	PDD Model D-4 for He PID HP5890	Instrumentation	\$6,630.00	<a href="http://www.chromtech.net.au/pdf2/vici09_222.pdf">http://www.chromtech.net.au/pdf2/vici09_222.pdf</a>
D-5-6890-220	PDD Model D-5 for He PID Plug-in system for Agilent 6890 220-VAC	Instrumentation	\$8,925.00	<a href="http://www.chromtech.net.au/pdf2/vici09_222.pdf">http://www.chromtech.net.au/pdf2/vici09_222.pdf</a>
D-5-6890	PDD Model D-5 for He PID Plug-in system for Agilent 6890	Instrumentation	\$8,925.00	<a href="http://www.chromtech.net.au/pdf2/vici09_222.pdf">http://www.chromtech.net.au/pdf2/vici09_222.pdf</a>
D-3-I-7890-220	PDD Model D-3 Plug-In system for Agilent 7890 220VAC	Instrumentation	\$7,395.00	<a href="http://www.chromtech.net.au/pdf2/vici09_222.pdf">http://www.chromtech.net.au/pdf2/vici09_222.pdf</a>
D-3-I-7890	PDD Model D-3 Plug-In system for Agilent 7890	Instrumentation	\$7,395.00	<a href="http://www.chromtech.net.au/pdf2/vici09_222.pdf">http://www.chromtech.net.au/pdf2/vici09_222.pdf</a>
D-3-I-HP-220	PDD Model D-3 Plug-In system for Agilent 6890 220VAC	Instrumentation	\$7,917.00	<a href="http://www.chromtech.net.au/pdf2/vici09_222.pdf">http://www.chromtech.net.au/pdf2/vici09_222.pdf</a>
D-2-I-220	PDD model D-2 helium photoionization mode, optimum trace packed column used 220VAC	Instrumentation	\$6,630.00	<a href="http://www.chromtech.net.au/pdf2/vici09_221.pdf">http://www.chromtech.net.au/pdf2/vici09_221.pdf</a>
D-2-I	PDD model D-2 helium photoionization mode, optimum trace packed column used 110V	Instrumentation	\$6,630.00	<a href="http://www.chromtech.net.au/pdf2/vici09_221.pdf">http://www.chromtech.net.au/pdf2/vici09_221.pdf</a>
D-2-220	PDD model D-2 Mode-selectable universal electron capture / photoionization 220VAC electron capture / photoionization detector system	Instrumentation	\$7,012.00	<a href="http://www.chromtech.net.au/pdf2/vici09_221.pdf">http://www.chromtech.net.au/pdf2/vici09_221.pdf</a>
PR50A60Z1	Valco internal fittings 1/16in 0-60 psi	Flow,Pressure,On/Off Devices	\$274.00	<a href="http://www.chromtech.net.au/pdf2/vici09_219.pdf">http://www.chromtech.net.au/pdf2/vici09_219.pdf</a>
PR50A15Z1	Valco internal fittings 1/16in 0-15 psi	Flow,Pressure,On/Off Devices	\$274.00	<a href="http://www.chromtech.net.au/pdf2/vici09_219.pdf">http://www.chromtech.net.au/pdf2/vici09_219.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

PR50A15Z2	Valco internal 1/8in 0-15 psi	Flow,Pressure,On/Off Devices	\$274.00	<a href="http://www.chromtech.net.au/pdf2/vici09_219.pdf">http://www.chromtech.net.au/pdf2/vici09_219.pdf</a>
PR50A15E2	External fittings 1/8in 0-15 psi	Flow,Pressure,On/Off Devices	\$274.00	<a href="http://www.chromtech.net.au/pdf2/vici09_219.pdf">http://www.chromtech.net.au/pdf2/vici09_219.pdf</a>
PR50A30Z1	Valco internal fittings 1/16in 0-30 psi	Flow,Pressure,On/Off Devices	\$274.00	<a href="http://www.chromtech.net.au/pdf2/vici09_219.pdf">http://www.chromtech.net.au/pdf2/vici09_219.pdf</a>
PR50A30E2	External fittings 1/8in 0-30 psi	Flow,Pressure,On/Off Devices	\$274.00	<a href="http://www.chromtech.net.au/pdf2/vici09_219.pdf">http://www.chromtech.net.au/pdf2/vici09_219.pdf</a>
PR50A60Z2	Valco internal 1/8in 0-60 psi	Flow,Pressure,On/Off Devices	\$274.00	<a href="http://www.chromtech.net.au/pdf2/vici09_219.pdf">http://www.chromtech.net.au/pdf2/vici09_219.pdf</a>
PR50A30Z2	Valco internal 1/8in 0-30 psi	Flow,Pressure,On/Off Devices	\$274.00	<a href="http://www.chromtech.net.au/pdf2/vici09_219.pdf">http://www.chromtech.net.au/pdf2/vici09_219.pdf</a>
PR50A60E2	External fittings 1/8in 0-60 psi	Flow,Pressure,On/Off Devices	\$274.00	<a href="http://www.chromtech.net.au/pdf2/vici09_219.pdf">http://www.chromtech.net.au/pdf2/vici09_219.pdf</a>
ZBNV1F-D	1/16in micrometering w/Valco Valves ZDV fittings Viton non-lubricated-low flow 2-175ml/min	Flow,Pressure,On/Off Devices	\$224.00	<a href="http://www.chromtech.net.au/pdf2/vici09_218.pdf">http://www.chromtech.net.au/pdf2/vici09_218.pdf</a>
ZBNV1-KZ	1/16in micrometering w/Valco Valves ZDV fittings Kalrez non-lubricated-std 2-225ml/min	Flow,Pressure,On/Off Devices	\$230.00	<a href="http://www.chromtech.net.au/pdf2/vici09_218.pdf">http://www.chromtech.net.au/pdf2/vici09_218.pdf</a>
ZBNV1-D	1/16in micrometering w/Valco Valves ZDV fittings Viton non-lubricated-std 2-225ml/min	Flow,Pressure,On/Off Devices	\$198.00	<a href="http://www.chromtech.net.au/pdf2/vici09_218.pdf">http://www.chromtech.net.au/pdf2/vici09_218.pdf</a>
ZBNV1	1/16in micrometering w/Valco Valves ZDV fittings Viton lubricated-low flow 2-175ml/min	Flow,Pressure,On/Off Devices	\$165.00	<a href="http://www.chromtech.net.au/pdf2/vici09_218.pdf">http://www.chromtech.net.au/pdf2/vici09_218.pdf</a>
BNV1-KZ	1/16in micrometering w/18in tubes Kalrez non-lubricated-std 2-225ml/min	Flow,Pressure,On/Off Devices	\$230.00	<a href="http://www.chromtech.net.au/pdf2/vici09_218.pdf">http://www.chromtech.net.au/pdf2/vici09_218.pdf</a>
ZBNV1F-KZ	1/16in micrometering w/Valco Valves ZDV fittings Viton non-lubricated-low flow 2-175ml/min	Flow,Pressure,On/Off Devices	\$249.00	<a href="http://www.chromtech.net.au/pdf2/vici09_218.pdf">http://www.chromtech.net.au/pdf2/vici09_218.pdf</a>
ZBNV1LF	1/16in micrometering w/Valco Valves ZDV fittings Viton lubricated -low flow 2-175ml/min	Flow,Pressure,On/Off Devices	\$172.00	<a href="http://www.chromtech.net.au/pdf2/vici09_218.pdf">http://www.chromtech.net.au/pdf2/vici09_218.pdf</a>
ZBNV1LF-D	1/16in micrometering w/Valco Valves ZDV fittings Kalrez non-lubricated -low flow 2-175ml/min	Flow,Pressure,On/Off Devices	\$204.00	<a href="http://www.chromtech.net.au/pdf2/vici09_218.pdf">http://www.chromtech.net.au/pdf2/vici09_218.pdf</a>
ZBNV1LF-KZ	1/16in micrometering w/Valco Valves ZDV fittings Kalrez non-lubricated -low flow 2-175ml/min	Flow,Pressure,On/Off Devices	\$236.00	<a href="http://www.chromtech.net.au/pdf2/vici09_218.pdf">http://www.chromtech.net.au/pdf2/vici09_218.pdf</a>
BNV1	1/16in micrometering w/18in tubes Viton lubricated-std 2-225ml/min	Flow,Pressure,On/Off Devices	\$165.00	<a href="http://www.chromtech.net.au/pdf2/vici09_218.pdf">http://www.chromtech.net.au/pdf2/vici09_218.pdf</a>
BNV1-D	1/16in micrometering w/18in tubes Viton non-lubricated-std 2-225ml/min	Flow,Pressure,On/Off Devices	\$198.00	<a href="http://www.chromtech.net.au/pdf2/vici09_218.pdf">http://www.chromtech.net.au/pdf2/vici09_218.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

BNV1LF	1/16in micrometering w/18in tubes Viton lubricated -low flow 2-175ml/min	Flow,Pressure,On/Off Devices	\$172.00	<a href="http://www.chromtech.net.au/pdf2/vici09_218.pdf">http://www.chromtech.net.au/pdf2/vici09_218.pdf</a>
BNV1LF-KZ	1/16 Micrometering Valve Kelrez non-lubricated Flow with 18" Tubes-low flow 2-175ml/min	Flow,Pressure,On/Off Devices	\$236.00	<a href="http://www.chromtech.net.au/pdf2/vici09_218.pdf">http://www.chromtech.net.au/pdf2/vici09_218.pdf</a>
BNV1LF-D	1/16 Micrometering Valve Viton Non-lubricated Low Flow with 18" Tubes-low flow 2-175ml/min	Flow,Pressure,On/Off Devices	\$204.00	<a href="http://www.chromtech.net.au/pdf2/vici09_218.pdf">http://www.chromtech.net.au/pdf2/vici09_218.pdf</a>
FC30SS2S3	Fixed span upstream ref SS body SS diaphragm 0-300 mL w/Spectrol 2-digit dial	Flow,Pressure,On/Off Devices	\$573.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC30SS4S4	Fixed span upstream ref SS body SS diaphragm 0-1.6 L w/Spectrol 4-digit dial	Flow,Pressure,On/Off Devices	\$586.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC31AV1	Fixed span upstream ref SS body Viton diaphragm 0-750 L w/Scrwdriver adjust	Flow,Pressure,On/Off Devices	\$280.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC30AV4S3	Fixed span upstream ref Aluminum body Viton diaphragm 0-1.6 L w/Spectrol 2-digit dial	Flow,Pressure,On/Off Devices	\$318.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC30SS3S3	Fixed span upstream ref SS body SS diaphragm 0-800 mL w/Spectrol 2-digit dial	Flow,Pressure,On/Off Devices	\$573.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC30SV3S3	Fixed span upstream ref SS body Viton diaphragm 0-800 mL w/Spectrol 2-digit dial	Flow,Pressure,On/Off Devices	\$555.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC30SV4S4	Fixed span upstream ref SS body Viton diaphragm 0-1.6 L w/Spectrol 4-digit dial	Flow,Pressure,On/Off Devices	\$573.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC30AV3S3	Fixed span upstream ref Aluminum body Viton diaphragm 0-800 mL w/Spectrol 2-digit dial	Flow,Pressure,On/Off Devices	\$318.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC30SV3S4	Fixed span upstream ref SS body Viton diaphragm 0-800 mL w/Spectrol 4-digit dial	Flow,Pressure,On/Off Devices	\$573.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC30SV2S3	Fixed span upstream ref SS body Viton diaphragm 0-300 mL w/Spectrol 2-digit dial	Flow,Pressure,On/Off Devices	\$555.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC30AS2S3	Fixed span upstream ref Aluminum body SS diaphragm 0-300 mL w/Spectrol 2-digit dial	Flow,Pressure,On/Off Devices	\$332.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC30AV2S3	Fixed span upstream ref Aluminum body Viton diaphragm 0-300 mL w/Spectrol 2-digit dial	Flow,Pressure,On/Off Devices	\$318.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC30AS3S3	Fixed span upstream ref Aluminum body SS diaphragm 0-800 mL w/Spectrol 2-digit dial	Flow,Pressure,On/Off Devices	\$332.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC30SV2S4	Fixed span upstream ref SS body Viton diaphragm 0-300 mL w/Spectrol 4-digit dial	Flow,Pressure,On/Off Devices	\$573.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC30SV4S3	Fixed span upstream ref SS body Viton diaphragm 0-1.6 L w/Spectrol 2-digit dial	Flow,Pressure,On/Off Devices	\$555.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

FC30SS4S3	Fixed span upstream ref SS body SS diaphragm 0-1.6 L w/Spectrol 2-digit dial	Flow,Pressure,On/Off Devices	\$573.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC30AV1S4	Fixed span upstream ref Aluminum body Viton diaphragm 0-200 mL w/Spectrol 4-digit dial	Flow,Pressure,On/Off Devices	\$332.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC30AS1S4	Fixed span upstream ref Aluminum body SS diaphragm 0-200 mL w/Spectrol 4-digit dial	Flow,Pressure,On/Off Devices	\$351.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC30SV1S4	Fixed span upstream ref SS body Viton diaphragm 0-200 mL w/Spectrol 4-digit dial	Flow,Pressure,On/Off Devices	\$573.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC30SS1S4	Fixed span upstream ref SS body SS diaphragm 0-200 mL w/Spectrol 4-digit dial	Flow,Pressure,On/Off Devices	\$586.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC30AV4S4	Fixed span upstream ref Aluminum body Viton diaphragm 0-1.6 L w/Spectrol 4-digit dial	Flow,Pressure,On/Off Devices	\$332.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC30AS2S4	Fixed span upstream ref Aluminum body SS diaphragm 0-300 mL w/Spectrol 4-digit dial	Flow,Pressure,On/Off Devices	\$351.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC30AS4S4	Fixed span upstream ref Aluminum body SS diaphragm 0-1.6 L w/Spectrol 4-digit dial	Flow,Pressure,On/Off Devices	\$351.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC30SS2S4	Fixed span upstream ref SS body SS diaphragm 0-300 mL w/Spectrol 4-digit dial	Flow,Pressure,On/Off Devices	\$586.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC30AV3S4	Fixed span upstream ref Aluminum body Viton diaphragm 0-800 mL w/Spectrol 4-digit dial	Flow,Pressure,On/Off Devices	\$332.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC30AS3S4	Fixed span upstream ref Aluminum body SS diaphragm 0-800 mL w/Spectrol 4-digit dial	Flow,Pressure,On/Off Devices	\$351.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC30SS4K	Fixed span upstream ref SS body SS diaphragm 0-1.6 L w/standard control knob	Flow,Pressure,On/Off Devices	\$530.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC30SS3S4	Fixed span upstream ref SS body SS diaphragm 0-800 mL w/Spectrol 4-digit dial	Flow,Pressure,On/Off Devices	\$586.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC30SS1S3	Fixed span upstream ref SS body SS diaphragm 0-200 mL w/Spectrol 2-digit dial	Flow,Pressure,On/Off Devices	\$573.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC30AV2S4	Fixed span upstream ref Aluminum body Viton diaphragm 0-300 mL w/Spectrol 4-digit dial	Flow,Pressure,On/Off Devices	\$332.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC30SS1K	Fixed span upstream ref SS body SS diaphragm 0-200mL w/standard control knob	Flow,Pressure,On/Off Devices	\$530.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC22SS1S3	Adjustable span-upstream ref SS body Viton diaphragm 1/8inF-NPT w/Spectrol 3-digit dial	Flow,Pressure,On/Off Devices	\$618.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC22AV1S4	Adjustable span-upstream ref Aluminum body Viton diaphragm 1/8inF-NPT w/Spectrol 4-digit dial	Flow,Pressure,On/Off Devices	\$363.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

FC22AS1S4	Adjustable span-upstream ref Aluminum body SS diaphragm 1/8inF-NPT w/Spectrol 4-digit dial	Flow,Pressure,On/Off Devices	\$382.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC22SV1S4	Adjustable span-upstream ref SS body Viton diaphragm 1/8inF-NPT w/Spectrol 4-digit dial	Flow,Pressure,On/Off Devices	\$618.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC22SS1S4	Adjustable span-upstream ref SS body Viton diaphragm 1/8inF-NPT w/Spectrol 4-digit dial	Flow,Pressure,On/Off Devices	\$638.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC30AV1K	Fixed span upstream ref Aluminum body Viton diaphragm 0-200mL w/standard control knob	Flow,Pressure,On/Off Devices	\$280.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC30AS1S3	Fixed span upstream ref Aluminum body SS diaphragm 0-200 mL w/Spectrol 2-digit dial	Flow,Pressure,On/Off Devices	\$332.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC30SV1K	Fixed span upstream ref SS body Viton diaphragm SS body 0-200mL w/standard control knob	Flow,Pressure,On/Off Devices	\$516.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC22AV1S3	Adjustable span-upstream ref Aluminum body 1/8inF-NPT Viton diaphragm w/Spectrol 3-digit dial	Flow,Pressure,On/Off Devices	\$351.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC30AV2K	Fixed span upstream ref Aluminum body Viton diaphragm 0-300 mL w/standard control knob	Flow,Pressure,On/Off Devices	\$280.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC30AS2K	Fixed span upstream ref Aluminum body SS diaphragm 0-200 mL w/standard control knob	Flow,Pressure,On/Off Devices	\$294.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC30SV2K	Fixed span upstream ref Aluminum body SS diaphragm 0-300 mL w/standard control knob	Flow,Pressure,On/Off Devices	\$516.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC30SS2K	Fixed span upstream ref SS body SS diaphragm 0-300 mL w/standard control knob	Flow,Pressure,On/Off Devices	\$530.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC30AV3K	Fixed span upstream ref Aluminum body Viton diaphragm 0-800 mL w/standard control knob	Flow,Pressure,On/Off Devices	\$280.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC30AS3K	Fixed span upstream ref Aluminum body SS diaphragm 0-800 mL w/standard control knob	Flow,Pressure,On/Off Devices	\$294.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC30SV3K	Fixed span upstream ref SS body Viton diaphragm 0-800 mL w/standard control knob	Flow,Pressure,On/Off Devices	\$516.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC30AS1K	Fixed span upstream ref Aluminum body SS diaphragm 0-200mL w/standard control knob	Flow,Pressure,On/Off Devices	\$294.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC22AS1K	Adjustable span-upstream ref Aluminum body 1/8inF-NPT SS diaphragm w/standard control knob	Flow,Pressure,On/Off Devices	\$318.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC30AS4S3	Fixed span upstream ref Aluminum body SS diaphragm 0-1.6 L w/Spectrol 2-digit dial	Flow,Pressure,On/Off Devices	\$332.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC30AV1S3	Fixed span upstream ref Aluminum body Viton diaphragm 0-200 mL w/Spectrol 2-digit dial	Flow,Pressure,On/Off Devices	\$318.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>



## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

FC30SV4K	Fixed span upstream ref SS body Viton diaphragm 0-1.6 L w/standard control knob	Flow,Pressure,On/Off Devices	\$516.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC30AS4K	Fixed span upstream ref Aluminum body SS diaphragm 0-1.6 L w/standard control knob	Flow,Pressure,On/Off Devices	\$294.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC30AV4K	Fixed span upstream ref Aluminum body Viton diaphragm 0-1.6 L w/standard control knob	Flow,Pressure,On/Off Devices	\$280.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC30SS3K	Fixed span upstream ref SS body SS diaphragm 0-800 mL w/standard control knob	Flow,Pressure,On/Off Devices	\$530.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC22SV1S3	Adjustable span-upstream ref SS body Viton 1/8inF-NPT diaphragm w/Spectrol 3-digit dial	Flow,Pressure,On/Off Devices	\$600.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC22AV1K	Adjustable span-upstream ref Aluminum body 1/8inF-NPT Viton diaphragm w/standard control knob	Flow,Pressure,On/Off Devices	\$306.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC22AS1S3	Adjustable span-upstream ref Aluminum body 1/8inF-NPT SS diaphragm w/Spectrol 3-digit dial	Flow,Pressure,On/Off Devices	\$363.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC22SV1K	Adjustable span-upstream ref SS body Viton 1/8inF-NPT diaphragm w/standard control knob	Flow,Pressure,On/Off Devices	\$561.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC22SS1K	Adjustable span-upstream ref SS body Viton 1/8inF-NPT diaphragm w/standard control knob	Flow,Pressure,On/Off Devices	\$580.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC30SV1S3	Fixed span upstream ref SS body Viton diaphragm 0-200 mL w/Spectrol 2-digit dial	Flow,Pressure,On/Off Devices	\$555.00	<a href="http://www.chromtech.net.au/pdf2/vici09_217.pdf">http://www.chromtech.net.au/pdf2/vici09_217.pdf</a>
FC10SS1S4	Fixed span upstream ref SS body SS diaphragm 0-150ml w/Spectrol 4-digit dial	Flow,Pressure,On/Off Devices	\$586.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10SS3S4	Fixed span upstream ref SS body SS diaphragm 0-850ml w/Spectrol 4-digit dial	Flow,Pressure,On/Off Devices	\$586.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10SV2S4	Fixed span upstream ref SS body Viton diaphragm 0-250ml w/Spectrol 4-digit dial	Flow,Pressure,On/Off Devices	\$573.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10AV4S4	Fixed span upstream ref Aluminum body Viton diaphragm 0-1.2 L w/Spectrol 4-digit dial	Flow,Pressure,On/Off Devices	\$332.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10AS4S4	Fixed span upstream ref Aluminum body SS diaphragm 0-1.2 L w/Spectrol 4-digit dial	Flow,Pressure,On/Off Devices	\$351.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10SV3S4	Fixed span upstream ref SS body Viton diaphragm 0-850ml w/Spectrol 4-digit dial	Flow,Pressure,On/Off Devices	\$573.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10AS3S4	Fixed span upstream ref Aluminum body SS diaphragm 0-850ml w/Spectrol 4-digit dial	Flow,Pressure,On/Off Devices	\$351.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10SS2S4	Fixed span upstream ref SS body SS diaphragm 0-250ml w/Spectrol 4-digit dial	Flow,Pressure,On/Off Devices	\$586.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

FC10AS2S4	Fixed span upstream ref Aluminum body SS diaphragm 0-250ml w/Spectrol 4-digit dial	Flow,Pressure,On/Off Devices	\$351.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10AV2S4	Fixed span upstream ref Aluminum body Viton diaphragm 0-250ml w/Spectrol 4-digit dial	Flow,Pressure,On/Off Devices	\$332.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10SV4S4	Fixed span upstream ref SS body Viton diaphragm 0-1.2 L w/Spectrol 4-digit dial	Flow,Pressure,On/Off Devices	\$573.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10SV6S4	Fixed span upstream ref SS body Viton diaphragm 0-10.0L w/Spectrol 4-digit dial	Flow,Pressure,On/Off Devices	\$573.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10AV3S4	Fixed span upstream ref Aluminum body Viton diaphragm 0-850ml w/Spectrol 4-digit dial	Flow,Pressure,On/Off Devices	\$332.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10SS4S4	Fixed span upstream ref SS body SS diaphragm0-1.2 L w/Spectrol 4-digit dial	Flow,Pressure,On/Off Devices	\$586.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10AV5S4	Fixed span upstream ref Aluminum body Viton diaphragm 0-4.5 L w/Spectrol 4-digit dial	Flow,Pressure,On/Off Devices	\$332.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10AS5S4	Fixed span upstream ref Aluminum body SS diaphragm 0-4.5 L w/Spectrol 4-digit dial	Flow,Pressure,On/Off Devices	\$351.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10SV5S4	Fixed span upstream ref SS body Viton diaphragm 0-4.5 L w/Spectrol 4-digit dial	Flow,Pressure,On/Off Devices	\$573.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10SS5S4	Fixed span upstream ref SS body SS diaphragm0-4.5 L w/Spectrol 4-digit dial	Flow,Pressure,On/Off Devices	\$586.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10SV1S4	Fixed span upstream ref SS body Viton diaphragm 0-150ml w/Spectrol 4-digit dial	Flow,Pressure,On/Off Devices	\$573.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10AS6S4	Fixed span upstream ref Aluminum body SS diaphragm 0-10.0L w/Spectrol 4-digit dial	Flow,Pressure,On/Off Devices	\$351.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10AS4S3	Fixed span upstream ref Aluminum body SS diaphragm 0-1.2 L w/Spectrol 3-digit dial	Flow,Pressure,On/Off Devices	\$338.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10SS6S4	Fixed span upstream ref SS body SS diaphragm0-10.0L w/Spectrol 4-digit dial	Flow,Pressure,On/Off Devices	\$586.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10AV4K	Fixed span upstream ref Aluminum body Viton diaphragm 0-1.2 L w/standard control knob	Flow,Pressure,On/Off Devices	\$280.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10SS3K	Fixed span upstream ref SS body SS diaphragm 0-850ml w/standard control knob	Flow,Pressure,On/Off Devices	\$530.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10SV3K	Fixed span upstream ref SS body Viton diaphragm 0-850ml w/standard control knob	Flow,Pressure,On/Off Devices	\$516.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10AV6S4	Fixed span upstream ref Aluminum body Viton diaphragm 0-10.0L w/Spectrol 4-digit dial	Flow,Pressure,On/Off Devices	\$332.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

FC10SV6K	Fixed span upstream ref SS body Viton diaphragm 0-10.0 L w/standard control knob	Flow,Pressure,On/Off Devices	\$516.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10SS4S3	Fixed span upstream ref SS body SS diaphragm 0-1.2 L w/Spectrol 3-digit dial	Flow,Pressure,On/Off Devices	\$573.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10AS2S3	Fixed span upstream ref Aluminum body SS diaphragm 0-250ml w/Spectrol 3-digit dial	Flow,Pressure,On/Off Devices	\$338.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10AV2S3	Fixed span upstream ref Aluminum body Viton diaphragm 0-250ml w/Spectrol 3-digit dial	Flow,Pressure,On/Off Devices	\$318.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10SS1S3	Fixed span upstream ref SS body SS diaphragm 0-150ml w/Spectrol 3-digit dial	Flow,Pressure,On/Off Devices	\$573.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10SV1S3	Fixed span upstream ref SS body Viton diaphragm 0-150ml w/Spectrol 3-digit dial	Flow,Pressure,On/Off Devices	\$555.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10AS1S3	Fixed span upstream ref Aluminum body SS diaphragm 0-150ml w/Spectrol 3-digit dial	Flow,Pressure,On/Off Devices	\$338.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10SS2S3	Fixed span upstream ref SS body SS diaphragm 0-250ml w/Spectrol 3-digit dial	Flow,Pressure,On/Off Devices	\$573.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10SS6K	Fixed span upstream ref SS body SS diaphragm 0-10.0 L w/standard control knob	Flow,Pressure,On/Off Devices	\$530.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10AV3S3	Fixed span upstream ref Aluminum body Viton diaphragm 0-850ml w/Spectrol 3-digit dial	Flow,Pressure,On/Off Devices	\$318.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10SS5K	Fixed span upstream ref SS body SS diaphragm 0-4.5 L w/standard control knob	Flow,Pressure,On/Off Devices	\$530.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10SV5K	Fixed span upstream ref SS body Viton diaphragm 0-4.5 L w/standard control knob	Flow,Pressure,On/Off Devices	\$516.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10AS5K	Fixed span upstream ref Aluminum body SS diaphragm 0-4.5 L w/standard control knob	Flow,Pressure,On/Off Devices	\$294.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10AV5K	Fixed span upstream ref Aluminum body Viton diaphragm 0-4.5 L w/standard control knob	Flow,Pressure,On/Off Devices	\$280.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10SS4K	Fixed span upstream ref SS body SS diaphragm 0-1.2 L w/standard control knob	Flow,Pressure,On/Off Devices	\$530.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10SV4K	Fixed span upstream ref SS body Viton diaphragm 0-1.2 L w/standard control knob	Flow,Pressure,On/Off Devices	\$516.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10AS4K	Fixed span upstream ref Aluminum body SS diaphragm 0-1.2 L w/standard control knob	Flow,Pressure,On/Off Devices	\$294.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10AV1S3	Fixed span upstream ref Aluminum body Viton diaphragm 0-150ml w/Spectrol 3-digit dial	Flow,Pressure,On/Off Devices	\$318.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

FC10AV5S3	Fixed span upstream ref Aluminum body Viton diaphragm 0-4.5 L w/Spectrol 3-digit dial	Flow,Pressure,On/Off Devices	\$318.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10AV1S4	Fixed span upstream ref Aluminum body Viton diaphragm 0-150ml w/Spectrol 4-digit dial	Flow,Pressure,On/Off Devices	\$332.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10SS6S3	Fixed span upstream ref SS body SS diaphragm0-10.0 L w/Spectrol 3-digit dial	Flow,Pressure,On/Off Devices	\$573.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10SV6S3	Fixed span upstream ref SS body Viton diaphragm 0-10.0 L w/Spectrol 3-digit dial	Flow,Pressure,On/Off Devices	\$555.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10AS6S3	Fixed span upstream ref Aluminum body SS diaphragm 0-10.0 L w/Spectrol 3-digit dial	Flow,Pressure,On/Off Devices	\$338.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10AV6S3	Fixed span upstream ref Aluminum body Viton diaphragm 0-10.0 L w/Spectrol 3-digit dial	Flow,Pressure,On/Off Devices	\$318.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10SS5S3	Fixed span upstream ref SS body SS diaphragm0-4.5 L w/Spectrol 3-digit dial	Flow,Pressure,On/Off Devices	\$573.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10SV2S3	Fixed span upstream ref SS body Viton diaphragm 0-250ml w/Spectrol 3-digit dial	Flow,Pressure,On/Off Devices	\$555.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10AS5S3	Fixed span upstream ref Aluminum body SS diaphragm 0-4.5 L w/Spectrol 3-digit dial	Flow,Pressure,On/Off Devices	\$338.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10AS1S4	Fixed span upstream ref Aluminum body SS diaphragm 0-150ml w/Spectrol 4-digit dial	Flow,Pressure,On/Off Devices	\$351.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10AS3K	Fixed span upstream ref Aluminum body SS diaphragm 0-850ml w/standard control knob	Flow,Pressure,On/Off Devices	\$294.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10SV4S3	Fixed span upstream ref SS body Viton diaphragm 0-1.2 L w/Spectrol 3-digit dial	Flow,Pressure,On/Off Devices	\$555.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10AV6K	Fixed span upstream ref Aluminum body Viton diaphragm 0-10.0 L w/standard control knob	Flow,Pressure,On/Off Devices	\$280.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10AV4S3	Fixed span upstream ref Aluminum body Viton diaphragm 0-1.2 L w/Spectrol 3-digit dial	Flow,Pressure,On/Off Devices	\$318.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10SS3S3	Fixed span upstream ref SS body SS diaphragm0-850ml w/Spectrol 3-digit dial	Flow,Pressure,On/Off Devices	\$573.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10SV3S3	Fixed span upstream ref SS body Viton diaphragm 0-850ml w/Spectrol 3-digit dial	Flow,Pressure,On/Off Devices	\$555.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10AS3S3	Fixed span upstream ref Aluminum body SS diaphragm 0-850ml w/Spectrol32-digit dial	Flow,Pressure,On/Off Devices	\$338.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10SV5S3	Fixed span upstream ref SS body Viton diaphragm 0-4.5 L w/Spectrol 3-digit dial	Flow,Pressure,On/Off Devices	\$555.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

FC10SS2K	Fixed span upstream ref SS body SS diaphragm 0-250ml w/standard control knob	Flow,Pressure,On/Off Devices	\$530.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10SV2K	Fixed span upstream ref Aluminum body SS diaphragm 0-250ml w/standard control knob	Flow,Pressure,On/Off Devices	\$516.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10AS2K	Fixed span upstream ref Aluminum body SS diaphragm 0-250ml w/standard control knob	Flow,Pressure,On/Off Devices	\$294.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10AV2K	Fixed span upstream ref Aluminum body Viton diaphragm 0-250ml w/standard control knob	Flow,Pressure,On/Off Devices	\$280.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10AV3K	Fixed span upstream ref Aluminum body Viton diaphragm 0-850ml w/standard control knob	Flow,Pressure,On/Off Devices	\$280.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10SS1K	Fixed span upstream ref SS body SS diaphragm 0-150mL w/standard control knob	Flow,Pressure,On/Off Devices	\$530.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10SV1K	Fixed span upstream ref SS body Viton diaphragm SS body 0-150mL w/standard control knob	Flow,Pressure,On/Off Devices	\$516.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10AS1K	Fixed span upstream ref Aluminum body SS diaphragm 0-150mL w/standard control knob	Flow,Pressure,On/Off Devices	\$294.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10AV1K	Fixed span upstream ref Aluminum body Viton diaphragm 0-150mL w/standard control knob	Flow,Pressure,On/Off Devices	\$280.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
FC10AS6K	Fixed span upstream ref Aluminum body SS diaphragm 0-10.0 L w/standard control knob	Valco Valves	\$294.00	<a href="http://www.chromtech.net.au/pdf2/vici09_216.pdf">http://www.chromtech.net.au/pdf2/vici09_216.pdf</a>
CVA50RS1	Condyne Combo valves 1/16in Valco 500ml/minRed	Flow,Pressure,On/Off Devices	\$128.00	<a href="http://www.chromtech.net.au/pdf2/vici09_215.pdf">http://www.chromtech.net.au/pdf2/vici09_215.pdf</a>
CVA10GS1	Condyne Combo valves 1/16in Valco 10ml/min Green	Flow,Pressure,On/Off Devices	\$128.00	<a href="http://www.chromtech.net.au/pdf2/vici09_215.pdf">http://www.chromtech.net.au/pdf2/vici09_215.pdf</a>
CVA1KYS1	Condyne Combo valves 1/16in Valco 1l/min Yellow	Flow,Pressure,On/Off Devices	\$128.00	<a href="http://www.chromtech.net.au/pdf2/vici09_215.pdf">http://www.chromtech.net.au/pdf2/vici09_215.pdf</a>
CVA50RS2	Condyne Combo valves 1/8in Valco 50ml/min Red	Flow,Pressure,On/Off Devices	\$128.00	<a href="http://www.chromtech.net.au/pdf2/vici09_215.pdf">http://www.chromtech.net.au/pdf2/vici09_215.pdf</a>
CVA500BS2	Condyne Combo valves 1/8in Valco 500ml/min Black	Flow,Pressure,On/Off Devices	\$128.00	<a href="http://www.chromtech.net.au/pdf2/vici09_215.pdf">http://www.chromtech.net.au/pdf2/vici09_215.pdf</a>
CVA500BS1	Condyne Combo valves 1/16in Valco 500ml/minBlack	Flow,Pressure,On/Off Devices	\$128.00	<a href="http://www.chromtech.net.au/pdf2/vici09_215.pdf">http://www.chromtech.net.au/pdf2/vici09_215.pdf</a>
CVA1KYS2	Condyne Combo valves 1/8in Valco 1l/min Yellow	Flow,Pressure,On/Off Devices	\$128.00	<a href="http://www.chromtech.net.au/pdf2/vici09_215.pdf">http://www.chromtech.net.au/pdf2/vici09_215.pdf</a>
CVA150US1	Condyne Combo valves 1/16in Valco 150ml/min Blue	Flow,Pressure,On/Off Devices	\$128.00	<a href="http://www.chromtech.net.au/pdf2/vici09_215.pdf">http://www.chromtech.net.au/pdf2/vici09_215.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

CVA10GS2	Condyne Combo valves 1/8in Valco 10ml/min Green	Flow,Pressure,On/Off Devices	\$128.00	<a href="http://www.chromtech.net.au/pdf2/vici09_215.pdf">http://www.chromtech.net.au/pdf2/vici09_215.pdf</a>
CVA150US2	Condyne Combo valves 1/8in Valco 150ml/min Blue	Flow,Pressure,On/Off Devices	\$128.00	<a href="http://www.chromtech.net.au/pdf2/vici09_215.pdf">http://www.chromtech.net.au/pdf2/vici09_215.pdf</a>
CNVEKB	Optional colored knobs-Black Standard (.62in)	Flow,Pressure,On/Off Devices	\$8.00	<a href="http://www.chromtech.net.au/pdf2/vici09_214.pdf">http://www.chromtech.net.au/pdf2/vici09_214.pdf</a>
CNVEKLB	Optional colored knobs-Black Long(1.25in)	Flow,Pressure,On/Off Devices	\$10.00	<a href="http://www.chromtech.net.au/pdf2/vici09_214.pdf">http://www.chromtech.net.au/pdf2/vici09_214.pdf</a>
CNVEKLS	Optional colored knobs-Silver Long(1.25in)	Flow,Pressure,On/Off Devices	\$10.00	<a href="http://www.chromtech.net.au/pdf2/vici09_214.pdf">http://www.chromtech.net.au/pdf2/vici09_214.pdf</a>
CNVEKS	Optional colored knobs-Silver Standard (.62in)	Flow,Pressure,On/Off Devices	\$8.00	<a href="http://www.chromtech.net.au/pdf2/vici09_214.pdf">http://www.chromtech.net.au/pdf2/vici09_214.pdf</a>
CNVEKLU	Optional colored knobs-Blue Long(1.25in)	Flow,Pressure,On/Off Devices	\$10.00	<a href="http://www.chromtech.net.au/pdf2/vici09_214.pdf">http://www.chromtech.net.au/pdf2/vici09_214.pdf</a>
CNVEKU	Optional colored knobs-Blue Standard(.62in)	Flow,Pressure,On/Off Devices	\$8.00	<a href="http://www.chromtech.net.au/pdf2/vici09_214.pdf">http://www.chromtech.net.au/pdf2/vici09_214.pdf</a>
CNVEKLR	Optional colored knobs-Red Long(1.25in)	Flow,Pressure,On/Off Devices	\$10.00	<a href="http://www.chromtech.net.au/pdf2/vici09_214.pdf">http://www.chromtech.net.au/pdf2/vici09_214.pdf</a>
CNVEKR	Optional colored knobs-Red Standard (.62in)	Flow,Pressure,On/Off Devices	\$8.00	<a href="http://www.chromtech.net.au/pdf2/vici09_214.pdf">http://www.chromtech.net.au/pdf2/vici09_214.pdf</a>
CNVEKLG	Optional colored knobs-Green Long (1.25in)	Flow,Pressure,On/Off Devices	\$10.00	<a href="http://www.chromtech.net.au/pdf2/vici09_214.pdf">http://www.chromtech.net.au/pdf2/vici09_214.pdf</a>
CNVEKG	Optional colored knobs-Green Standard (.62in)	Flow,Pressure,On/Off Devices	\$8.00	<a href="http://www.chromtech.net.au/pdf2/vici09_214.pdf">http://www.chromtech.net.au/pdf2/vici09_214.pdf</a>
CNV1A10S1	10ml/min max @ 40 psi He or N2 Aluminum body	Flow,Pressure,On/Off Devices	\$108.00	<a href="http://www.chromtech.net.au/pdf2/vici09_214.pdf">http://www.chromtech.net.au/pdf2/vici09_214.pdf</a>
CNV1A500S1	500ml/min @ max 40 psi He or N2 Aluminum body	Flow,Pressure,On/Off Devices	\$108.00	<a href="http://www.chromtech.net.au/pdf2/vici09_214.pdf">http://www.chromtech.net.au/pdf2/vici09_214.pdf</a>
CNV1S250S1	250ml/min @ max 40 psi He or N2 SS body	Flow,Pressure,On/Off Devices	\$172.00	<a href="http://www.chromtech.net.au/pdf2/vici09_214.pdf">http://www.chromtech.net.au/pdf2/vici09_214.pdf</a>
CNV1A250S1	250ml/min @ max 40 psi He or N2 Aluminum body	Flow,Pressure,On/Off Devices	\$108.00	<a href="http://www.chromtech.net.au/pdf2/vici09_214.pdf">http://www.chromtech.net.au/pdf2/vici09_214.pdf</a>
CNV1S150S1	150 ml/min max @ 40 psi He or N2 SS body	Flow,Pressure,On/Off Devices	\$16.00	<a href="http://www.chromtech.net.au/pdf2/vici09_214.pdf">http://www.chromtech.net.au/pdf2/vici09_214.pdf</a>
CNV1A150S1	150ml/min max @ 40 psi He or N2 Aluminum body	Flow,Pressure,On/Off Devices	\$108.00	<a href="http://www.chromtech.net.au/pdf2/vici09_214.pdf">http://www.chromtech.net.au/pdf2/vici09_214.pdf</a>

VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

CNV1S50S1	50ml/min @ max 40 psi He or N2 SS body	Flow,Pressure,On/Off Devices	\$172.00	<a href="http://www.chromtech.net.au/pdf2/vici09_214.pdf">http://www.chromtech.net.au/pdf2/vici09_214.pdf</a>
CNV1A50S1	50 ml/min max @ 40 psi He or N2 Aluminum body	Flow,Pressure,On/Off Devices	\$108.00	<a href="http://www.chromtech.net.au/pdf2/vici09_214.pdf">http://www.chromtech.net.au/pdf2/vici09_214.pdf</a>
CNV1S10S1	10 ml/min max @ 40 psi He or N2 SS body	Flow,Pressure,On/Off Devices	\$172.00	<a href="http://www.chromtech.net.au/pdf2/vici09_214.pdf">http://www.chromtech.net.au/pdf2/vici09_214.pdf</a>
CNV1S500S1	500ml/min @ max 40 psi He or N2 SS body	Flow,Pressure,On/Off Devices	\$172.00	<a href="http://www.chromtech.net.au/pdf2/vici09_214.pdf">http://www.chromtech.net.au/pdf2/vici09_214.pdf</a>
ASFV2HT	w/2in standoff High temperature/high pressure 1/8in 1.50mm	Flow,Pressure,On/Off Devices	\$484.00	<a href="http://www.chromtech.net.au/pdf2/vici09_213.pdf">http://www.chromtech.net.au/pdf2/vici09_213.pdf</a>
ASFVL	w/2in standoff Standard temperature 1/16in 0.75mm	Flow,Pressure,On/Off Devices	\$382.00	<a href="http://www.chromtech.net.au/pdf2/vici09_213.pdf">http://www.chromtech.net.au/pdf2/vici09_213.pdf</a>
ASFVHT4	w/2in standoff High temperature/high pressure 1/16in 0.75mm	Flow,Pressure,On/Off Devices	\$484.00	<a href="http://www.chromtech.net.au/pdf2/vici09_213.pdf">http://www.chromtech.net.au/pdf2/vici09_213.pdf</a>
ASFVHT	w/2in standoff High temperature/high pressure 1/16in 0.50mm	Flow,Pressure,On/Off Devices	\$408.00	<a href="http://www.chromtech.net.au/pdf2/vici09_213.pdf">http://www.chromtech.net.au/pdf2/vici09_213.pdf</a>
SFVL	Prime/Purge Valve Manual Standard temperature 1/16in 0.75mm	Flow,Pressure,On/Off Devices	\$255.00	<a href="http://www.chromtech.net.au/pdf2/vici09_213.pdf">http://www.chromtech.net.au/pdf2/vici09_213.pdf</a>
ASFVLHT4	w/4in standoff High temperature/high pressure 1/16in	Flow,Pressure,On/Off Devices	\$510.00	<a href="http://www.chromtech.net.au/pdf2/vici09_213.pdf">http://www.chromtech.net.au/pdf2/vici09_213.pdf</a>
ASFV2HT4	w/4in standoff High temperature/high pressure 1/8in 1.50mm	Flow,Pressure,On/Off Devices	\$561.00	<a href="http://www.chromtech.net.au/pdf2/vici09_213.pdf">http://www.chromtech.net.au/pdf2/vici09_213.pdf</a>
ASFVOLHT4	w/4in standoff High temperature/high pressure 1/16in 0.75mm	Flow,Pressure,On/Off Devices	\$510.00	<a href="http://www.chromtech.net.au/pdf2/vici09_213.pdf">http://www.chromtech.net.au/pdf2/vici09_213.pdf</a>
ASFVLHT	w/2in standoff High temperature/high pressure 1/16in 0.75mm	Flow,Pressure,On/Off Devices	\$465.00	<a href="http://www.chromtech.net.au/pdf2/vici09_213.pdf">http://www.chromtech.net.au/pdf2/vici09_213.pdf</a>
ASFV	w/2in standoff Standard temperature 1/16in 0.50mm	Flow,Pressure,On/Off Devices	\$345.00	<a href="http://www.chromtech.net.au/pdf2/vici09_213.pdf">http://www.chromtech.net.au/pdf2/vici09_213.pdf</a>
SFV	Prime/Purge Valve, Manual Standard temperature 1/16in 0.50mm	Flow,Pressure,On/Off Devices	\$216.00	<a href="http://www.chromtech.net.au/pdf2/vici09_213.pdf">http://www.chromtech.net.au/pdf2/vici09_213.pdf</a>
ASFVO2HT4	w/4in standoff High temperature/high pressure 1/8in 1.50mm	Flow,Pressure,On/Off Devices	\$561.00	<a href="http://www.chromtech.net.au/pdf2/vici09_213.pdf">http://www.chromtech.net.au/pdf2/vici09_213.pdf</a>
ASFVOLHT	w/2in standoff High temperature/high pressure 1/16in 0.75mm	Flow,Pressure,On/Off Devices	\$465.00	<a href="http://www.chromtech.net.au/pdf2/vici09_213.pdf">http://www.chromtech.net.au/pdf2/vici09_213.pdf</a>
ASFVO	w/2in standoff Standard temperature 1/16in 0.50mm	Flow,Pressure,On/Off Devices	\$345.00	<a href="http://www.chromtech.net.au/pdf2/vici09_213.pdf">http://www.chromtech.net.au/pdf2/vici09_213.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

SFVOL	On/Off Valve Standard temperature 1/16in 0.75mm	Flow,Pressure,On/Off Devices	\$255.00	<a href="http://www.chromtech.net.au/pdf2/vici09_213.pdf">http://www.chromtech.net.au/pdf2/vici09_213.pdf</a>
ASFVOL	w/2in standoff Standard temperature 1/16in 0.75mm	Flow,Pressure,On/Off Devices	\$382.00	<a href="http://www.chromtech.net.au/pdf2/vici09_213.pdf">http://www.chromtech.net.au/pdf2/vici09_213.pdf</a>
SFVOHT	On/Off Valve High Temp,High Pressure,1/16inx0.5mm	Flow,Pressure,On/Off Devices	\$286.00	<a href="http://www.chromtech.net.au/pdf2/vici09_213.pdf">http://www.chromtech.net.au/pdf2/vici09_213.pdf</a>
SFVOHT4	On/Off Valve High Temp,High Pressure,1/16inx0.5mm,Manual	Flow,Pressure,On/Off Devices	\$332.00	<a href="http://www.chromtech.net.au/pdf2/vici09_213.pdf">http://www.chromtech.net.au/pdf2/vici09_213.pdf</a>
ASFVOHT	w/2in standoff High temperature/high pressure 1/16in 0.50mm	Flow,Pressure,On/Off Devices	\$408.00	<a href="http://www.chromtech.net.au/pdf2/vici09_213.pdf">http://www.chromtech.net.au/pdf2/vici09_213.pdf</a>
ASFVO2HT	w/2in standoff High temperature/high pressure 1/8in 1.50mm	Flow,Pressure,On/Off Devices	\$484.00	<a href="http://www.chromtech.net.au/pdf2/vici09_213.pdf">http://www.chromtech.net.au/pdf2/vici09_213.pdf</a>
SFVO	On/Off Valve Standard temperature 1/16in 0.50mm	Flow,Pressure,On/Off Devices	\$210.00	<a href="http://www.chromtech.net.au/pdf2/vici09_213.pdf">http://www.chromtech.net.au/pdf2/vici09_213.pdf</a>
ASFVOHT4	w/4in standoff High temperature/high pressure 1/16in 0.50mm	Flow,Pressure,On/Off Devices	\$484.00	<a href="http://www.chromtech.net.au/pdf2/vici09_213.pdf">http://www.chromtech.net.au/pdf2/vici09_213.pdf</a>
VSH	Valve Spanner Handle	Actuators,Accessories	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_211.pdf">http://www.chromtech.net.au/pdf2/vici09_211.pdf</a>
TEMPLATE1	Template	Actuators,Accessories	\$4.00	<a href="http://www.chromtech.net.au/pdf2/vici09_211.pdf">http://www.chromtech.net.au/pdf2/vici09_211.pdf</a>
PV	PinVise Drill Index -.0135in to .039in	Actuators,Accessories	\$38.00	<a href="http://www.chromtech.net.au/pdf2/vici09_211.pdf">http://www.chromtech.net.au/pdf2/vici09_211.pdf</a>
SFVO	Standard temperature 1/16in 0.50mm	Flow,Pressure,On/Off Devices	\$210.00	<a href="http://www.chromtech.net.au/pdf2/vici09_211.pdf">http://www.chromtech.net.au/pdf2/vici09_211.pdf</a>
OEW-3	Open end Wrench -1/2inx9/16in for 1/4in nuts	Actuators,Accessories	\$15.00	<a href="http://www.chromtech.net.au/pdf2/vici09_210.pdf">http://www.chromtech.net.au/pdf2/vici09_210.pdf</a>
PM	Pencil Magnet	Actuators,Accessories	\$6.00	<a href="http://www.chromtech.net.au/pdf2/vici09_210.pdf">http://www.chromtech.net.au/pdf2/vici09_210.pdf</a>
SWH4	Custom Socket Wrench -1/4in slotted	Actuators,Accessories	\$9.00	<a href="http://www.chromtech.net.au/pdf2/vici09_210.pdf">http://www.chromtech.net.au/pdf2/vici09_210.pdf</a>
HKS	Hex Key Set -0.05in to 5/32in x8	Actuators,Accessories	\$15.00	<a href="http://www.chromtech.net.au/pdf2/vici09_210.pdf">http://www.chromtech.net.au/pdf2/vici09_210.pdf</a>
OEW	Open end Wrench -3/16inx1/4in for 1/32in and 1/16in nuts	Actuators,Accessories	\$8.00	<a href="http://www.chromtech.net.au/pdf2/vici09_210.pdf">http://www.chromtech.net.au/pdf2/vici09_210.pdf</a>
OEW-2	Open end Wrench -3/8inx7/16in for 1/8in nuts	Actuators,Accessories	\$15.00	<a href="http://www.chromtech.net.au/pdf2/vici09_210.pdf">http://www.chromtech.net.au/pdf2/vici09_210.pdf</a>
IZERA1.5	Int/Ext Reducer 1/16in to 1/32in 1.0mm bore	Valco Fittings-FS	\$56.00	<a href="http://www.chromtech.net.au/pdf2/vici09_21.pdf">http://www.chromtech.net.au/pdf2/vici09_21.pdf</a>
IZERA1.5M	Int/Ext Reducer 1/16in to 1/32in 0.50mm bore	Valco Fittings-FS	\$58.00	<a href="http://www.chromtech.net.au/pdf2/vici09_21.pdf">http://www.chromtech.net.au/pdf2/vici09_21.pdf</a>
IZERA1.5C	Int/Ext Reducer 1/16in to 1/32in 0.25mm bore	Valco Fittings-FS	\$68.00	<a href="http://www.chromtech.net.au/pdf2/vici09_21.pdf">http://www.chromtech.net.au/pdf2/vici09_21.pdf</a>
CMH11L	Closemount hardware -for Air or standard electric actuators -for -for Cheminert valves low pressure design	Actuators,Accessories	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_208.pdf">http://www.chromtech.net.au/pdf2/vici09_208.pdf</a>
CMHMP	Closemount hardware -for Air or standard electric actuators -for Valco two position valves w/no mounting holes	Actuators,Accessories	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_208.pdf">http://www.chromtech.net.au/pdf2/vici09_208.pdf</a>



## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

CMHMP	Closemount hardware -for Air or standard electric actuators -for Valco two position valves w/no mounting holes	Actuators,Accessories	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_208.pdf">http://www.chromtech.net.au/pdf2/vici09_208.pdf</a>
CMH13L	Closemount hardware -for Microelectric actuators For Cheminert multiposition valves low pressure design	Actuators,Accessories	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_208.pdf">http://www.chromtech.net.au/pdf2/vici09_208.pdf</a>
CMH13H	Closemount hardware -for Microelectric actuators For Cheminert multiposition valves high pressure design	Actuators,Accessories	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_208.pdf">http://www.chromtech.net.au/pdf2/vici09_208.pdf</a>
CMH12L	Closemount hardware -for Microelectric actuators For Cheminert two position valves low pressure design	Actuators,Accessories	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_208.pdf">http://www.chromtech.net.au/pdf2/vici09_208.pdf</a>
CMH12H	Closemount hardware -for Microelectric actuators For Cheminert two position valves high pressure design	Actuators,Accessories	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_208.pdf">http://www.chromtech.net.au/pdf2/vici09_208.pdf</a>
CMH13	Closemount hardware -for Microelectric actuators For Valco multiposition valves (UW Type only)	Actuators,Accessories	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_208.pdf">http://www.chromtech.net.au/pdf2/vici09_208.pdf</a>
CMH12H	Closemount hardware -for Microelectric actuators For Cheminert two position valves high pressure design	Actuators,Accessories	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_208.pdf">http://www.chromtech.net.au/pdf2/vici09_208.pdf</a>
WMMA	Closemount hardware -for manual valves w/threaded mounting holes	Actuators,Accessories	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_208.pdf">http://www.chromtech.net.au/pdf2/vici09_208.pdf</a>
CMH	Closemount hardware -for Air or standard electric actuators -for Valco two position valves w/1 or 2 mounting holes	Actuators,Accessories	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_208.pdf">http://www.chromtech.net.au/pdf2/vici09_208.pdf</a>
WMMA10	Closemount hardware -for manual valves w/unthreaded mounting holes	Actuators,Accessories	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_208.pdf">http://www.chromtech.net.au/pdf2/vici09_208.pdf</a>
CMH11H	Closemount hardware -for Air or standard electric actuators -for Cheminert valves high pressure design	Actuators,Accessories	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_208.pdf">http://www.chromtech.net.au/pdf2/vici09_208.pdf</a>
CMH12H	Closemount hardware -for Microelectric actuators For Cheminert two position valves high pressure design	Actuators,Accessories	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_208.pdf">http://www.chromtech.net.au/pdf2/vici09_208.pdf</a>
HWSC-SC8-8B	Screws -for Standard electric actuators Valco multiposition valves	Actuators,Accessories	\$0.00	<a href="http://www.chromtech.net.au/pdf2/vici09_205.pdf">http://www.chromtech.net.au/pdf2/vici09_205.pdf</a>
HWSC-SC8-8B	Screws -for Standard electric actuators Valco multiposition valves	Actuators,Accessories	\$0.00	<a href="http://www.chromtech.net.au/pdf2/vici09_205.pdf">http://www.chromtech.net.au/pdf2/vici09_205.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

2SOAMMP	Standoff assembly -for Microelectric actuators Cheminert multiposition valves	Actuators,Accessories	\$90.00	<a href="http://www.chromtech.net.au/pdf2/vici09_205.pdf">http://www.chromtech.net.au/pdf2/vici09_205.pdf</a>
CR10	Clamp ring -for Microelectric actuators Cheminert multiposition valves	Actuators,Accessories	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_205.pdf">http://www.chromtech.net.au/pdf2/vici09_205.pdf</a>
HWSC-SC8-8TDH	Screws -for Microelectric actuators Cheminert multiposition valves	Actuators,Accessories	\$0.00	<a href="http://www.chromtech.net.au/pdf2/vici09_205.pdf">http://www.chromtech.net.au/pdf2/vici09_205.pdf</a>
2SOA	Standoff assembly -for Standard electric actuators Valco 2-position valves w/1 or 2 mounting holes	Actuators,Accessories	\$57.00	<a href="http://www.chromtech.net.au/pdf2/vici09_205.pdf">http://www.chromtech.net.au/pdf2/vici09_205.pdf</a>
HWSC-SC8-8B	Screws -for Standard electric actuators Valco multiposition valves	Actuators,Accessories	\$0.00	<a href="http://www.chromtech.net.au/pdf2/vici09_205.pdf">http://www.chromtech.net.au/pdf2/vici09_205.pdf</a>
2SOAMP	Standoff assembly -for Microelectric actuators Cheminert 2-position valves	Actuators,Accessories	\$57.00	<a href="http://www.chromtech.net.au/pdf2/vici09_205.pdf">http://www.chromtech.net.au/pdf2/vici09_205.pdf</a>
CR3	Clamp ring -for Air actuators Cheminert valves	Actuators,Accessories	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_205.pdf">http://www.chromtech.net.au/pdf2/vici09_205.pdf</a>
HWSC-SC8-8B	Screws -for Standard electric actuators Valco multiposition valves	Actuators,Accessories	\$0.00	<a href="http://www.chromtech.net.au/pdf2/vici09_205.pdf">http://www.chromtech.net.au/pdf2/vici09_205.pdf</a>
CR8	Clamp ring -for Microelectric actuators Cheminert 2-position valves	Actuators,Accessories	\$14.00	<a href="http://www.chromtech.net.au/pdf2/vici09_205.pdf">http://www.chromtech.net.au/pdf2/vici09_205.pdf</a>
CR3	Clamp ring -for Air actuators Cheminert valves	Actuators,Accessories	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_205.pdf">http://www.chromtech.net.au/pdf2/vici09_205.pdf</a>
CR3	Clamp ring -for Air actuators Cheminert valves	Actuators,Accessories	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_205.pdf">http://www.chromtech.net.au/pdf2/vici09_205.pdf</a>
4SOWK	Standoff assembly -for manual Valco W Type two position valves w/1 or 2 mounting holes	Actuators,Accessories	\$90.00	<a href="http://www.chromtech.net.au/pdf2/vici09_205.pdf">http://www.chromtech.net.au/pdf2/vici09_205.pdf</a>
4SOWKMP	Standoff assembly -for manual Cheminert valves	Actuators,Accessories	\$90.00	<a href="http://www.chromtech.net.au/pdf2/vici09_205.pdf">http://www.chromtech.net.au/pdf2/vici09_205.pdf</a>
4SOUTH	Standoff assembly -for manual Valco UW Type two position valves w/1 or 2 mounting holes	Actuators,Accessories	\$82.00	<a href="http://www.chromtech.net.au/pdf2/vici09_205.pdf">http://www.chromtech.net.au/pdf2/vici09_205.pdf</a>
4SOUTHMP	Standoff assembly -for manual Valco UW Type two position valves w/no mounting holes	Actuators,Accessories	\$82.00	<a href="http://www.chromtech.net.au/pdf2/vici09_205.pdf">http://www.chromtech.net.au/pdf2/vici09_205.pdf</a>
4SOWKMP	Standoff assembly -for manual Cheminert valves	Actuators,Accessories	\$90.00	<a href="http://www.chromtech.net.au/pdf2/vici09_205.pdf">http://www.chromtech.net.au/pdf2/vici09_205.pdf</a>
CR3	Clamp ring -for Air actuators Cheminert valves	Actuators,Accessories	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_205.pdf">http://www.chromtech.net.au/pdf2/vici09_205.pdf</a>
HWSC-SC8-8B	Screws -for Standard electric actuators Valco multiposition valves	Actuators,Accessories	\$0.00	<a href="http://www.chromtech.net.au/pdf2/vici09_205.pdf">http://www.chromtech.net.au/pdf2/vici09_205.pdf</a>
2SOAMP	Standoff assembly -for Microelectric actuators Cheminert 2-position valves	Actuators,Accessories	\$57.00	<a href="http://www.chromtech.net.au/pdf2/vici09_205.pdf">http://www.chromtech.net.au/pdf2/vici09_205.pdf</a>
HWSC-SC8-6	Screws -for Air actuators Cheminert valves	Actuators,Accessories	\$0.00	<a href="http://www.chromtech.net.au/pdf2/vici09_205.pdf">http://www.chromtech.net.au/pdf2/vici09_205.pdf</a>
4SOA	Standoff assembly -for Air actuators Valco 2-position valves w/1 or 2 mounting holes	Actuators,Accessories	\$57.00	<a href="http://www.chromtech.net.au/pdf2/vici09_205.pdf">http://www.chromtech.net.au/pdf2/vici09_205.pdf</a>
CR3	Clamp ring -for Air actuators Cheminert valves	Actuators,Accessories	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_205.pdf">http://www.chromtech.net.au/pdf2/vici09_205.pdf</a>
HWSC-SC8-6	Screws -for Air actuators Cheminert valves	Actuators,Accessories	\$0.00	<a href="http://www.chromtech.net.au/pdf2/vici09_205.pdf">http://www.chromtech.net.au/pdf2/vici09_205.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

4SOAMP	Standoff assembly -for Air actuators Cheminert valves	Actuators,Accessories	\$57.00	<a href="http://www.chromtech.net.au/pdf2/vici09_205.pdf">http://www.chromtech.net.au/pdf2/vici09_205.pdf</a>
2SOAMP	Standoff assembly -for Microelectric actuators Cheminert 2-position valves	Actuators,Accessories	\$57.00	<a href="http://www.chromtech.net.au/pdf2/vici09_205.pdf">http://www.chromtech.net.au/pdf2/vici09_205.pdf</a>
CR3	Clamp ring -for Air actuators Cheminert valves	Actuators,Accessories	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_205.pdf">http://www.chromtech.net.au/pdf2/vici09_205.pdf</a>
2SOAMP	Standoff assembly -for Microelectric actuators Cheminert 2-position valves	Actuators,Accessories	\$57.00	<a href="http://www.chromtech.net.au/pdf2/vici09_205.pdf">http://www.chromtech.net.au/pdf2/vici09_205.pdf</a>
4SOAMP	Standoff assembly -for Air actuators Cheminert valves	Actuators,Accessories	\$57.00	<a href="http://www.chromtech.net.au/pdf2/vici09_205.pdf">http://www.chromtech.net.au/pdf2/vici09_205.pdf</a>
CR3	Clamp ring -for Air actuators Cheminert valves	Actuators,Accessories	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_205.pdf">http://www.chromtech.net.au/pdf2/vici09_205.pdf</a>
HWSC-SC8-6	Screws -for Air actuators Cheminert valves	Actuators,Accessories	\$0.00	<a href="http://www.chromtech.net.au/pdf2/vici09_205.pdf">http://www.chromtech.net.au/pdf2/vici09_205.pdf</a>
4SOAMP	Standoff assembly -for Air actuators Cheminert valves	Actuators,Accessories	\$57.00	<a href="http://www.chromtech.net.au/pdf2/vici09_205.pdf">http://www.chromtech.net.au/pdf2/vici09_205.pdf</a>
CR3	Clamp ring -for Air actuators Cheminert valves	Actuators,Accessories	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_205.pdf">http://www.chromtech.net.au/pdf2/vici09_205.pdf</a>
2SOAMMP	Standoff assembly -for Microelectric actuators Cheminert multiposition valves	Actuators,Accessories	\$90.00	<a href="http://www.chromtech.net.au/pdf2/vici09_205.pdf">http://www.chromtech.net.au/pdf2/vici09_205.pdf</a>
CR10	Clamp ring -for Microelectric actuators Cheminert multiposition valves	Actuators,Accessories	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_205.pdf">http://www.chromtech.net.au/pdf2/vici09_205.pdf</a>
HWSC-SC8-6	Screws -for Air actuators Cheminert valves	Actuators,Accessories	\$0.00	<a href="http://www.chromtech.net.au/pdf2/vici09_205.pdf">http://www.chromtech.net.au/pdf2/vici09_205.pdf</a>
HWSC-SC8-8B	Screws -for Standard electric actuators Valco multiposition valves	Actuators,Accessories	\$0.00	<a href="http://www.chromtech.net.au/pdf2/vici09_205.pdf">http://www.chromtech.net.au/pdf2/vici09_205.pdf</a>
CR8	Clamp ring -for Microelectric actuators Cheminert 2-position valves	Actuators,Accessories	\$14.00	<a href="http://www.chromtech.net.au/pdf2/vici09_205.pdf">http://www.chromtech.net.au/pdf2/vici09_205.pdf</a>
HWSC-SC8-8B	Screws -for Standard electric actuators Valco multiposition valves	Actuators,Accessories	\$0.00	<a href="http://www.chromtech.net.au/pdf2/vici09_205.pdf">http://www.chromtech.net.au/pdf2/vici09_205.pdf</a>
CR8	Clamp ring -for Microelectric actuators Cheminert 2-position valves	Actuators,Accessories	\$14.00	<a href="http://www.chromtech.net.au/pdf2/vici09_205.pdf">http://www.chromtech.net.au/pdf2/vici09_205.pdf</a>
2SOA	Standoff assembly -for Standard electric actuators Valco 2-position valves w/1 or 2 mounting holes	Actuators,Accessories	\$57.00	<a href="http://www.chromtech.net.au/pdf2/vici09_205.pdf">http://www.chromtech.net.au/pdf2/vici09_205.pdf</a>
2SOAMP	Standoff assembly -for Microelectric actuators Cheminert 2-position valves	Actuators,Accessories	\$57.00	<a href="http://www.chromtech.net.au/pdf2/vici09_205.pdf">http://www.chromtech.net.au/pdf2/vici09_205.pdf</a>
HWSC-SC8-8TDH	Screws -for Microelectric actuators Cheminert multiposition valves	Actuators,Accessories	\$0.00	<a href="http://www.chromtech.net.au/pdf2/vici09_205.pdf">http://www.chromtech.net.au/pdf2/vici09_205.pdf</a>
RAD	Right Angle Drive w/2in standoff assemblyClosemount h'ware	Actuators,Accessories	\$306.00	<a href="http://www.chromtech.net.au/pdf2/vici09_204.pdf">http://www.chromtech.net.au/pdf2/vici09_204.pdf</a>
2RAD	Right Angle Drive w/Closemount h'ware	Actuators,Accessories	\$363.00	<a href="http://www.chromtech.net.au/pdf2/vici09_204.pdf">http://www.chromtech.net.au/pdf2/vici09_204.pdf</a>
WMMASO	Knob for a W type valve	Actuators,Accessories	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_204.pdf">http://www.chromtech.net.au/pdf2/vici09_204.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

SOUTH	T-handle for a UW type valve	Actuators,Accessories	\$26.00	<a href="http://www.chromtech.net.au/pdf2/vici09_204.pdf">http://www.chromtech.net.au/pdf2/vici09_204.pdf</a>
I-21014-01	Replace Thermocouple	Actuators,Accessories	\$27.00	<a href="http://www.chromtech.net.au/pdf2/vici09_203.pdf">http://www.chromtech.net.au/pdf2/vici09_203.pdf</a>
ITC10399-220	Heated Enclosures and Controllers -ITC Instrumentation temperature controller 0C to 399C 230 VAC	Actuators,Accessories	\$516.00	<a href="http://www.chromtech.net.au/pdf2/vici09_203.pdf">http://www.chromtech.net.au/pdf2/vici09_203.pdf</a>
ITC10399	Heated Enclosures and Controllers -ITC Instrumentation temperature controller 0C to 399C 110 VAC	Actuators,Accessories	\$516.00	<a href="http://www.chromtech.net.au/pdf2/vici09_203.pdf">http://www.chromtech.net.au/pdf2/vici09_203.pdf</a>
CM	Heated Enclosures and Controllers -Column mandrel w/o heater	Actuators,Accessories	\$82.00	<a href="http://www.chromtech.net.au/pdf2/vici09_203.pdf">http://www.chromtech.net.au/pdf2/vici09_203.pdf</a>
HCE2	Heated Enclosures and Controllers -Heated column enclosure 2-columns 8inx5-1/4inx4ind 65W/350C	Actuators,Accessories	\$382.00	<a href="http://www.chromtech.net.au/pdf2/vici09_203.pdf">http://www.chromtech.net.au/pdf2/vici09_203.pdf</a>
HCEC	Heated Enclosures and Controllers -Heated column enclosure 1-column 8inx8inx6ind 65W/350C	Actuators,Accessories	\$465.00	<a href="http://www.chromtech.net.au/pdf2/vici09_203.pdf">http://www.chromtech.net.au/pdf2/vici09_203.pdf</a>
HCE1	Heated Enclosures and Controllers -Heated column enclosure 1-column 4inx4-1/4inx3-5/8ind 65W/350C	Actuators,Accessories	\$312.00	<a href="http://www.chromtech.net.au/pdf2/vici09_203.pdf">http://www.chromtech.net.au/pdf2/vici09_203.pdf</a>
HCEB	Heated Enclosures and Controllers -Heated column enclosure 1-column 4-1/4inx5-1/8inx3-5/8ind 65W/350C	Actuators,Accessories	\$312.00	<a href="http://www.chromtech.net.au/pdf2/vici09_203.pdf">http://www.chromtech.net.au/pdf2/vici09_203.pdf</a>
HA3	Heater assembly,use with HVE3	Actuators,Accessories	\$224.00	<a href="http://www.chromtech.net.au/pdf2/vici09_202.pdf">http://www.chromtech.net.au/pdf2/vici09_202.pdf</a>
I-21208-06	Heater Cartridges,65W,220VAC	Actuators,Accessories	\$66.00	<a href="http://www.chromtech.net.au/pdf2/vici09_202.pdf">http://www.chromtech.net.au/pdf2/vici09_202.pdf</a>
HA2	Heater assembly,use with HVE2	Actuators,Accessories	\$134.00	<a href="http://www.chromtech.net.au/pdf2/vici09_202.pdf">http://www.chromtech.net.au/pdf2/vici09_202.pdf</a>
HA6	Heater assembly,use with HVEC6	Actuators,Accessories	\$3.00	<a href="http://www.chromtech.net.au/pdf2/vici09_202.pdf">http://www.chromtech.net.au/pdf2/vici09_202.pdf</a>
HBS	Heater Block-Low Mass	Actuators,Accessories	\$72.00	<a href="http://www.chromtech.net.au/pdf2/vici09_202.pdf">http://www.chromtech.net.au/pdf2/vici09_202.pdf</a>
HB	Heater Block-Std, 1 valve	Actuators,Accessories	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_202.pdf">http://www.chromtech.net.au/pdf2/vici09_202.pdf</a>
HB1N	Heater Block-Std, 1 NV type (1/32) valve	Actuators,Accessories	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_202.pdf">http://www.chromtech.net.au/pdf2/vici09_202.pdf</a>
I-21208-32	Heater Cartridges,100W,110VAC	Actuators,Accessories	\$52.00	<a href="http://www.chromtech.net.au/pdf2/vici09_202.pdf">http://www.chromtech.net.au/pdf2/vici09_202.pdf</a>
I-21208-33	Heater Cartridges,100W,220VAC	Actuators,Accessories	\$52.00	<a href="http://www.chromtech.net.au/pdf2/vici09_202.pdf">http://www.chromtech.net.au/pdf2/vici09_202.pdf</a>
I-21208-05	Heater Cartridges,65W,110VAC	Actuators,Accessories	\$58.00	<a href="http://www.chromtech.net.au/pdf2/vici09_202.pdf">http://www.chromtech.net.au/pdf2/vici09_202.pdf</a>
6HA1	Heater assembly,use with HVEA or HVEB	Actuators,Accessories	\$134.00	<a href="http://www.chromtech.net.au/pdf2/vici09_202.pdf">http://www.chromtech.net.au/pdf2/vici09_202.pdf</a>
HA1T	Heater assembly,use with HVEC	Actuators,Accessories	\$134.00	<a href="http://www.chromtech.net.au/pdf2/vici09_202.pdf">http://www.chromtech.net.au/pdf2/vici09_202.pdf</a>
HVEA	Heated valve enclosures for two position and multiposition valves 1-valve 4in x 4-1/4in x 3-5/8ind 65W/350C	Actuators,Accessories	\$312.00	<a href="http://www.chromtech.net.au/pdf2/vici09_201.pdf">http://www.chromtech.net.au/pdf2/vici09_201.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

HVEC	Heated valve enclosures for two position and multiposition valves 1-valve 8in x 8in x 6ind 100W/350C	Actuators,Accessories	\$465.00	<a href="http://www.chromtech.net.au/pdf2/vici09_201.pdf">http://www.chromtech.net.au/pdf2/vici09_201.pdf</a>
HVEBN	Heated valve enclosures for two position and multiposition valves 1x1/32in valve 4-1/4in x 5-1/8in x 3-5/8ind 65W/350 C	Actuators,Accessories	\$312.00	<a href="http://www.chromtech.net.au/pdf2/vici09_201.pdf">http://www.chromtech.net.au/pdf2/vici09_201.pdf</a>
HVEAN	Heated valve enclosures for two position and multiposition valves 1x1/32in valve 4in x 4-1/4in x 3-5/8ind 65W/350C	Actuators,Accessories	\$312.00	<a href="http://www.chromtech.net.au/pdf2/vici09_201.pdf">http://www.chromtech.net.au/pdf2/vici09_201.pdf</a>
HVE2	Heated valve enclosures for two position and multiposition valves 2-valves 8in x 5-1/4in x 4ind 125W/350	Actuators,Accessories	\$382.00	<a href="http://www.chromtech.net.au/pdf2/vici09_201.pdf">http://www.chromtech.net.au/pdf2/vici09_201.pdf</a>
HVE6	Heated valve enclosures for two position and multiposition valves 6-valves 13-3/4in x 8in x 6ind 300W/C	Actuators,Accessories	\$702.00	<a href="http://www.chromtech.net.au/pdf2/vici09_201.pdf">http://www.chromtech.net.au/pdf2/vici09_201.pdf</a>
HVE3	Heated valve enclosures for two position and multiposition valves 3-valves 13-1/2in x 5-3/4in x 4ind 150W/350	Actuators,Accessories	\$465.00	<a href="http://www.chromtech.net.au/pdf2/vici09_201.pdf">http://www.chromtech.net.au/pdf2/vici09_201.pdf</a>
HVEB	Heated valve enclosures for two position and multiposition valves 1-valve 4-1/4in x 5-1/8in x 3-5/8ind 65W/350C	Actuators,Accessories	\$312.00	<a href="http://www.chromtech.net.au/pdf2/vici09_201.pdf">http://www.chromtech.net.au/pdf2/vici09_201.pdf</a>
FSMUAS1.5	Make-up adaptor 1/16in to 1/32in x1.5"L x 0.75mm bore	Valco Fittings-FS	\$118.00	<a href="http://www.chromtech.net.au/pdf2/vici09_20.pdf">http://www.chromtech.net.au/pdf2/vici09_20.pdf</a>
FSMUA1.5	Make-up adaptor 1/16in to 1/32in x 3.5"L x 0.75mm bore	Valco Fittings-FS	\$118.00	<a href="http://www.chromtech.net.au/pdf2/vici09_20.pdf">http://www.chromtech.net.au/pdf2/vici09_20.pdf</a>
FSMUAS1.5M	Make-up adaptor 1/16in to 1/32in x 1.5"L x 0.50mm bore	Valco Fittings-FS	\$134.00	<a href="http://www.chromtech.net.au/pdf2/vici09_20.pdf">http://www.chromtech.net.au/pdf2/vici09_20.pdf</a>
FSMUAS1.5L	Make-up adaptor 1/16in to 1/32in x 1.5"L x 1.0mm bore	Valco Fittings-FS	\$118.00	<a href="http://www.chromtech.net.au/pdf2/vici09_20.pdf">http://www.chromtech.net.au/pdf2/vici09_20.pdf</a>
PFW60	Position Feedback for manual valves Cheminert W type 6 port valves	Actuators,Accessories	\$10.00	<a href="http://www.chromtech.net.au/pdf2/vici09_199.pdf">http://www.chromtech.net.au/pdf2/vici09_199.pdf</a>
DVI	Digital valve interface DVI for 2 position air actuators 110 VAC	Actuators,Accessories	\$312.00	<a href="http://www.chromtech.net.au/pdf2/vici09_199.pdf">http://www.chromtech.net.au/pdf2/vici09_199.pdf</a>
DVI-220	Digital valve interface DVI for 2 position air actuators 230 VAC	Actuators,Accessories	\$312.00	<a href="http://www.chromtech.net.au/pdf2/vici09_199.pdf">http://www.chromtech.net.au/pdf2/vici09_199.pdf</a>
PFW90	Position Feedback for manual valves Cheminert W type 4 8 and 10 port valves	Actuators,Accessories	\$108.00	<a href="http://www.chromtech.net.au/pdf2/vici09_199.pdf">http://www.chromtech.net.au/pdf2/vici09_199.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

PFW36	8,10 Port Position Feedback for Valco W type valves	Actuators,Accessories	\$108.00	<a href="http://www.chromtech.net.au/pdf2/vici09_199.pdf">http://www.chromtech.net.au/pdf2/vici09_199.pdf</a>
PFC2	Position Feedback for manual valves for Cheminert C2 and C4 series valves	Actuators,Accessories	\$108.00	<a href="http://www.chromtech.net.au/pdf2/vici09_199.pdf">http://www.chromtech.net.au/pdf2/vici09_199.pdf</a>
PFAF	PFAF Position feedback for two position air actuators	Actuators,Accessories	\$192.00	<a href="http://www.chromtech.net.au/pdf2/vici09_199.pdf">http://www.chromtech.net.au/pdf2/vici09_199.pdf</a>
41E1-120VAC	4-Way Solenoid air valve 41E1 for multiposition air actuators 110 VAC	Actuators,Accessories	\$159.00	<a href="http://www.chromtech.net.au/pdf2/vici09_198.pdf">http://www.chromtech.net.au/pdf2/vici09_198.pdf</a>
41E1-24VAC	4-Way Solenoid air valve 41E1 for multiposition air actuators 24 VAC	Actuators,Accessories	\$159.00	<a href="http://www.chromtech.net.au/pdf2/vici09_198.pdf">http://www.chromtech.net.au/pdf2/vici09_198.pdf</a>
41E1-12VDC	4-Way Solenoid air valve 41E1 for multiposition air actuators 12VDC	Actuators,Accessories	\$159.00	<a href="http://www.chromtech.net.au/pdf2/vici09_198.pdf">http://www.chromtech.net.au/pdf2/vici09_198.pdf</a>
41E1-24VDC	4-Way Solenoid air valve 41E1 for multiposition air actuators 24 VDC	Actuators,Accessories	\$159.00	<a href="http://www.chromtech.net.au/pdf2/vici09_198.pdf">http://www.chromtech.net.au/pdf2/vici09_198.pdf</a>
MSVA-110VAC	Manifold 3-way solenoid assembly MSVA for two position air actuators 110 VAC	Actuators,Accessories	\$300.00	<a href="http://www.chromtech.net.au/pdf2/vici09_198.pdf">http://www.chromtech.net.au/pdf2/vici09_198.pdf</a>
MSVA-220VAC	Manifold 3-way solenoid assembly MSVA for two position air actuators 230 VAC	Actuators,Accessories	\$300.00	<a href="http://www.chromtech.net.au/pdf2/vici09_198.pdf">http://www.chromtech.net.au/pdf2/vici09_198.pdf</a>
MSVA-24VAC	Manifold 3-way solenoid assembly MSVA for two position air actuators 24 VAC	Actuators,Accessories	\$300.00	<a href="http://www.chromtech.net.au/pdf2/vici09_198.pdf">http://www.chromtech.net.au/pdf2/vici09_198.pdf</a>
MSVA-12VDC	Manifold 3-way solenoid assembly MSVA for two position air actuators 12 VDC	Actuators,Accessories	\$300.00	<a href="http://www.chromtech.net.au/pdf2/vici09_198.pdf">http://www.chromtech.net.au/pdf2/vici09_198.pdf</a>
MSVA-24VDC	Manifold 3-way solenoid assembly MSVA for two position air actuators 24 VDC	Actuators,Accessories	\$300.00	<a href="http://www.chromtech.net.au/pdf2/vici09_198.pdf">http://www.chromtech.net.au/pdf2/vici09_198.pdf</a>
HSSA	High speed switching accessory HSSA for all Air Actuators	Actuators,Accessories	\$153.00	<a href="http://www.chromtech.net.au/pdf2/vici09_198.pdf">http://www.chromtech.net.au/pdf2/vici09_198.pdf</a>
41E1-220VAC	4-Way Solenoid air valve 41E1 for multiposition air actuators 230 VAC	Actuators,Accessories	\$159.00	<a href="http://www.chromtech.net.au/pdf2/vici09_198.pdf">http://www.chromtech.net.au/pdf2/vici09_198.pdf</a>
SVI	SVI Serial valve interface for all air and standard electric actuators 110 VAC	Actuators,Accessories	\$930.00	<a href="http://www.chromtech.net.au/pdf2/vici09_197.pdf">http://www.chromtech.net.au/pdf2/vici09_197.pdf</a>
SVI-220	SVI Serial valve interface for all air and standard electric actuators 230 VAC	Actuators,Accessories	\$930.00	<a href="http://www.chromtech.net.au/pdf2/vici09_197.pdf">http://www.chromtech.net.au/pdf2/vici09_197.pdf</a>
I-22239	DVI/SVI Interface Cable	Actuators,Accessories	\$70.00	<a href="http://www.chromtech.net.au/pdf2/vici09_197.pdf">http://www.chromtech.net.au/pdf2/vici09_197.pdf</a>
DVSP4	Digital sequence programmer for all air and electric actuators 4 intervals 110VAC	Actuators,Accessories	\$1,148.00	<a href="http://www.chromtech.net.au/pdf2/vici09_196.pdf">http://www.chromtech.net.au/pdf2/vici09_196.pdf</a>
DVSP4-220	Digital sequence programmer for all air and electric actuators 4 intervals 230VAC	Actuators,Accessories	\$1,148.00	<a href="http://www.chromtech.net.au/pdf2/vici09_196.pdf">http://www.chromtech.net.au/pdf2/vici09_196.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

A45S	Standard Valve Two position Air actuators8 ports 45 deg rotation use w/existing standoff	Actuators,Accessories	\$27.00	<a href="http://www.chromtech.net.au/pdf2/vici09_195.pdf">http://www.chromtech.net.au/pdf2/vici09_195.pdf</a>
AT36S	High temperature Two position Air actuators10 ports 36 deg rotation use w/existing standoff	Actuators,Accessories	\$267.00	<a href="http://www.chromtech.net.au/pdf2/vici09_195.pdf">http://www.chromtech.net.au/pdf2/vici09_195.pdf</a>
AT30	High temperature Two position Air actuators12 ports 30 deg rotation w/+losemount assembly	Actuators,Accessories	\$318.00	<a href="http://www.chromtech.net.au/pdf2/vici09_195.pdf">http://www.chromtech.net.au/pdf2/vici09_195.pdf</a>
AT90	High temperature Two position Air Actuators 3-4 ports 90 deg rotation w/+losemount assembly	Actuators,Accessories	\$318.00	<a href="http://www.chromtech.net.au/pdf2/vici09_195.pdf">http://www.chromtech.net.au/pdf2/vici09_195.pdf</a>
AT902	High temperature Two position Air Actuators 3-4 ports 90 deg rotation w/2in standoff assembly	Actuators,Accessories	\$318.00	<a href="http://www.chromtech.net.au/pdf2/vici09_195.pdf">http://www.chromtech.net.au/pdf2/vici09_195.pdf</a>
AT90S	High temperature Two position Air Actuators 3-4 ports 90 deg rotation use w/existing standoff	Actuators,Accessories	\$267.00	<a href="http://www.chromtech.net.au/pdf2/vici09_195.pdf">http://www.chromtech.net.au/pdf2/vici09_195.pdf</a>
AT60	High temperature Two position Air actuators6 ports 60 deg rotation w/+losemount assembly	Actuators,Accessories	\$318.00	<a href="http://www.chromtech.net.au/pdf2/vici09_195.pdf">http://www.chromtech.net.au/pdf2/vici09_195.pdf</a>
AT602	High temperature Two position Air actuators6 ports 60 deg rotation w/2in standoff assembly	Actuators,Accessories	\$318.00	<a href="http://www.chromtech.net.au/pdf2/vici09_195.pdf">http://www.chromtech.net.au/pdf2/vici09_195.pdf</a>
AT60S	High temperature Two position Air actuators6 ports 60 deg rotation use w/existing standoff	Actuators,Accessories	\$267.00	<a href="http://www.chromtech.net.au/pdf2/vici09_195.pdf">http://www.chromtech.net.au/pdf2/vici09_195.pdf</a>
A304	Standard Valve Two position Air actuators12 ports 30 deg rotation w/4in standoff assembly	Actuators,Accessories	\$318.00	<a href="http://www.chromtech.net.au/pdf2/vici09_195.pdf">http://www.chromtech.net.au/pdf2/vici09_195.pdf</a>
AT452	High temperature Two position Air actuators8 ports 45 deg rotation w/2in standoff assembly	Actuators,Accessories	\$318.00	<a href="http://www.chromtech.net.au/pdf2/vici09_195.pdf">http://www.chromtech.net.au/pdf2/vici09_195.pdf</a>
A30	Standard Valve Two position Air actuators12 ports 30 deg rotation closemount assembly	Actuators,Accessories	\$318.00	<a href="http://www.chromtech.net.au/pdf2/vici09_195.pdf">http://www.chromtech.net.au/pdf2/vici09_195.pdf</a>
A90	Standard Valve Two position Air actuators 3-4 ports 90 deg rotation w/closemount assembly	Actuators,Accessories	\$318.00	<a href="http://www.chromtech.net.au/pdf2/vici09_195.pdf">http://www.chromtech.net.au/pdf2/vici09_195.pdf</a>
A904	Standard Valve Two position Air actuators 3-4 ports 90 deg rotation w/4in standoff assembly	Actuators,Accessories	\$318.00	<a href="http://www.chromtech.net.au/pdf2/vici09_195.pdf">http://www.chromtech.net.au/pdf2/vici09_195.pdf</a>
A90S	Standard Valve Two position Air actuators 3-4 ports 90 deg rotation use w/existing standoff	Actuators,Accessories	\$267.00	<a href="http://www.chromtech.net.au/pdf2/vici09_195.pdf">http://www.chromtech.net.au/pdf2/vici09_195.pdf</a>
A60	Standard Valve Two position Air actuators6 ports 60 deg rotation w/closemount assembly	Actuators,Accessories	\$318.00	<a href="http://www.chromtech.net.au/pdf2/vici09_195.pdf">http://www.chromtech.net.au/pdf2/vici09_195.pdf</a>
A604	Standard Valve Two position Air actuators6 ports 60 deg rotation w/4in standoff assembly	Actuators,Accessories	\$318.00	<a href="http://www.chromtech.net.au/pdf2/vici09_195.pdf">http://www.chromtech.net.au/pdf2/vici09_195.pdf</a>
A60S	Standard Valve Two position Air actuators6 ports 60 deg rotation use w/existing standoff	Actuators,Accessories	\$27.00	<a href="http://www.chromtech.net.au/pdf2/vici09_195.pdf">http://www.chromtech.net.au/pdf2/vici09_195.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

OR	Two position Air actuators Replacement Standard O-rings	Actuators,Accessories	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_195.pdf">http://www.chromtech.net.au/pdf2/vici09_195.pdf</a>
A45	Standard Valve Two position Air actuators8 ports 45 deg rotation w/closemount assembly	Actuators,Accessories	\$318.00	<a href="http://www.chromtech.net.au/pdf2/vici09_195.pdf">http://www.chromtech.net.au/pdf2/vici09_195.pdf</a>
AT302	High temperature Two position Air actuators12 ports 30 deg rotation w/2in standoff assembly	Actuators,Accessories	\$318.00	<a href="http://www.chromtech.net.au/pdf2/vici09_195.pdf">http://www.chromtech.net.au/pdf2/vici09_195.pdf</a>
A454	Standard Valve Two position Air actuators8 ports 45 deg rotation w/4in standoff assembly	Actuators,Accessories	\$318.00	<a href="http://www.chromtech.net.au/pdf2/vici09_195.pdf">http://www.chromtech.net.au/pdf2/vici09_195.pdf</a>
AT45	High temperature Two position Air actuators8 ports 45 deg rotation w/+losemount assembly	Actuators,Accessories	\$318.00	<a href="http://www.chromtech.net.au/pdf2/vici09_195.pdf">http://www.chromtech.net.au/pdf2/vici09_195.pdf</a>
A36S	Standard Valve Two position Air actuators10 ports 36 deg rotation use w/existing standoff	Actuators,Accessories	\$267.00	<a href="http://www.chromtech.net.au/pdf2/vici09_195.pdf">http://www.chromtech.net.au/pdf2/vici09_195.pdf</a>
A364	Standard Valve Two position Air actuators10 ports 36 deg rotation w/4in standoff assembly	Actuators,Accessories	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_195.pdf">http://www.chromtech.net.au/pdf2/vici09_195.pdf</a>
A36	Standard Valve Two position Air actuators10 ports 36 deg rotation w/closemount assembly	Actuators,Accessories	\$318.00	<a href="http://www.chromtech.net.au/pdf2/vici09_195.pdf">http://www.chromtech.net.au/pdf2/vici09_195.pdf</a>
AT30S	High temperature Two position Air actuators12 ports 30 deg rotation use w/existing standoff	Actuators,Accessories	\$267.00	<a href="http://www.chromtech.net.au/pdf2/vici09_195.pdf">http://www.chromtech.net.au/pdf2/vici09_195.pdf</a>
ORT	Two position Air actuators Replacement High temp O-rings	Actuators,Accessories	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_195.pdf">http://www.chromtech.net.au/pdf2/vici09_195.pdf</a>
F-TCF	Actuator Compression Fitting	Actuators,Accessories	\$4.00	<a href="http://www.chromtech.net.au/pdf2/vici09_195.pdf">http://www.chromtech.net.au/pdf2/vici09_195.pdf</a>
A30S	Standard Valve Two position Air actuators12 ports 30 deg rotation use w/existing standoff	Actuators,Accessories	\$267.00	<a href="http://www.chromtech.net.au/pdf2/vici09_195.pdf">http://www.chromtech.net.au/pdf2/vici09_195.pdf</a>
AT45S	High temperature Two position Air actuators8 ports 45 deg rotation use w/existing standoff	Valco Valves	\$267.00	<a href="http://www.chromtech.net.au/pdf2/vici09_195.pdf">http://www.chromtech.net.au/pdf2/vici09_195.pdf</a>
AT362	High temperature Two position Air actuators10 ports 36 deg rotation w/2in standoff assembly	Valco Valves	\$318.00	<a href="http://www.chromtech.net.au/pdf2/vici09_195.pdf">http://www.chromtech.net.au/pdf2/vici09_195.pdf</a>
AT36	High temperature Two position Air actuators10 ports 36 deg rotation w/+losemount assembly	Valco Valves	\$318.00	<a href="http://www.chromtech.net.au/pdf2/vici09_195.pdf">http://www.chromtech.net.au/pdf2/vici09_195.pdf</a>
A8	Standard Multiposition Valve Air actuators 8 position w/closemount assembly	Actuators,Accessories	\$516.00	<a href="http://www.chromtech.net.au/pdf2/vici09_194.pdf">http://www.chromtech.net.au/pdf2/vici09_194.pdf</a>
A162	Standard Multiposition Valve Air actuators 16 position w/2in standoff assembly	Actuators,Accessories	\$516.00	<a href="http://www.chromtech.net.au/pdf2/vici09_194.pdf">http://www.chromtech.net.au/pdf2/vici09_194.pdf</a>
A16	Standard Multiposition Valve Air actuators 16 position closemount assembly	Actuators,Accessories	\$516.00	<a href="http://www.chromtech.net.au/pdf2/vici09_194.pdf">http://www.chromtech.net.au/pdf2/vici09_194.pdf</a>
A12S	Standard Multiposition Valve Air actuators 12 position use w/Existing standoff	Actuators,Accessories	\$465.00	<a href="http://www.chromtech.net.au/pdf2/vici09_194.pdf">http://www.chromtech.net.au/pdf2/vici09_194.pdf</a>



## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

A122	Standard Multiposition Valve Air actuators 12 position w/2in standoff assembly	Actuators,Accessories	\$516.00	<a href="http://www.chromtech.net.au/pdf2/vici09_194.pdf">http://www.chromtech.net.au/pdf2/vici09_194.pdf</a>
A12	Standard Multiposition Valve Air actuators 12 position closemount assembly	Actuators,Accessories	\$516.00	<a href="http://www.chromtech.net.au/pdf2/vici09_194.pdf">http://www.chromtech.net.au/pdf2/vici09_194.pdf</a>
A10S	Standard Multiposition Valve Air actuators 10 position use w/Existing standoff	Actuators,Accessories	\$465.00	<a href="http://www.chromtech.net.au/pdf2/vici09_194.pdf">http://www.chromtech.net.au/pdf2/vici09_194.pdf</a>
A102	Standard Multiposition Valve Air actuators 10 position w/2in standoff assembly	Actuators,Accessories	\$516.00	<a href="http://www.chromtech.net.au/pdf2/vici09_194.pdf">http://www.chromtech.net.au/pdf2/vici09_194.pdf</a>
A10	Standard Multiposition Valve Air actuators 10 position w/closemount assembly	Actuators,Accessories	\$516.00	<a href="http://www.chromtech.net.au/pdf2/vici09_194.pdf">http://www.chromtech.net.au/pdf2/vici09_194.pdf</a>
A16S	Standard Multiposition Valve Air actuators 16 position use w/Existing standoff	Actuators,Accessories	\$465.00	<a href="http://www.chromtech.net.au/pdf2/vici09_194.pdf">http://www.chromtech.net.au/pdf2/vici09_194.pdf</a>
A82	Standard Multiposition Valve Air actuators 8 position w/2in standoff assembly	Actuators,Accessories	\$516.00	<a href="http://www.chromtech.net.au/pdf2/vici09_194.pdf">http://www.chromtech.net.au/pdf2/vici09_194.pdf</a>
A6S	Standard Multiposition Valve Air actuators 6 position use w/Existing standoff	Actuators,Accessories	\$465.00	<a href="http://www.chromtech.net.au/pdf2/vici09_194.pdf">http://www.chromtech.net.au/pdf2/vici09_194.pdf</a>
A62	Standard Multiposition Valve Air actuators 6 position w/2in standoff assembly	Actuators,Accessories	\$516.00	<a href="http://www.chromtech.net.au/pdf2/vici09_194.pdf">http://www.chromtech.net.au/pdf2/vici09_194.pdf</a>
A6	Standard Multiposition Valve Air actuators 6 position w/closemount assembly	Actuators,Accessories	\$516.00	<a href="http://www.chromtech.net.au/pdf2/vici09_194.pdf">http://www.chromtech.net.au/pdf2/vici09_194.pdf</a>
A4S	Standard Multiposition Valve Air actuators 4 position use w/Existing standoff	Actuators,Accessories	\$465.00	<a href="http://www.chromtech.net.au/pdf2/vici09_194.pdf">http://www.chromtech.net.au/pdf2/vici09_194.pdf</a>
A42	Standard Multiposition Valve Air actuators 4 position w/2in standoff assembly	Actuators,Accessories	\$516.00	<a href="http://www.chromtech.net.au/pdf2/vici09_194.pdf">http://www.chromtech.net.au/pdf2/vici09_194.pdf</a>
A8S	Standard Multiposition Valve Air actuators 8 position use w/Existing standoff	Actuators,Accessories	\$465.00	<a href="http://www.chromtech.net.au/pdf2/vici09_194.pdf">http://www.chromtech.net.au/pdf2/vici09_194.pdf</a>
AT16	High temperature Multiposition Air actuators 16 position w/+losemount assembly	Actuators,Accessories	\$516.00	<a href="http://www.chromtech.net.au/pdf2/vici09_194.pdf">http://www.chromtech.net.au/pdf2/vici09_194.pdf</a>
AT84	High temperature Multiposition Air actuators 8 position w/4in standoff assembly	Actuators,Accessories	\$516.00	<a href="http://www.chromtech.net.au/pdf2/vici09_194.pdf">http://www.chromtech.net.au/pdf2/vici09_194.pdf</a>
AT8S	High temperature Multiposition Air actuators 8 position use w/Existing standoff	Actuators,Accessories	\$465.00	<a href="http://www.chromtech.net.au/pdf2/vici09_194.pdf">http://www.chromtech.net.au/pdf2/vici09_194.pdf</a>
AT10	High temperature Multiposition Air actuators 10 position w/+losemount assembly	Actuators,Accessories	\$516.00	<a href="http://www.chromtech.net.au/pdf2/vici09_194.pdf">http://www.chromtech.net.au/pdf2/vici09_194.pdf</a>
AT104	High temperature Multiposition Air actuators 10 position w/4in standoff assembly	Actuators,Accessories	\$516.00	<a href="http://www.chromtech.net.au/pdf2/vici09_194.pdf">http://www.chromtech.net.au/pdf2/vici09_194.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

AT10S	High temperature Multiposition Air actuators 10 position use w/Existing standoff	Actuators,Accessories	\$465.00	<a href="http://www.chromtech.net.au/pdf2/vici09_194.pdf">http://www.chromtech.net.au/pdf2/vici09_194.pdf</a>
AT12	High temperature Multiposition Air actuators 12 position w/+losemount assembly	Actuators,Accessories	\$516.00	<a href="http://www.chromtech.net.au/pdf2/vici09_194.pdf">http://www.chromtech.net.au/pdf2/vici09_194.pdf</a>
AT44	High temperature Multiposition Air actuators 4 position w/4in standoff assembly	Actuators,Accessories	\$516.00	<a href="http://www.chromtech.net.au/pdf2/vici09_194.pdf">http://www.chromtech.net.au/pdf2/vici09_194.pdf</a>
AT12S	High temperature Multiposition Air actuators 12 position use w/Existing standoff	Actuators,Accessories	\$465.00	<a href="http://www.chromtech.net.au/pdf2/vici09_194.pdf">http://www.chromtech.net.au/pdf2/vici09_194.pdf</a>
AT4	High temperature Multiposition Air actuators 4 position w/+losemount assembly	Actuators,Accessories	\$516.00	<a href="http://www.chromtech.net.au/pdf2/vici09_194.pdf">http://www.chromtech.net.au/pdf2/vici09_194.pdf</a>
AT164	High temperature Multiposition Air actuators 16 position w/4in standoff assembly	Actuators,Accessories	\$516.00	<a href="http://www.chromtech.net.au/pdf2/vici09_194.pdf">http://www.chromtech.net.au/pdf2/vici09_194.pdf</a>
AT16S	High temperature Multiposition Air actuators 16 position use w/Existing standoff	Actuators,Accessories	\$465.00	<a href="http://www.chromtech.net.au/pdf2/vici09_194.pdf">http://www.chromtech.net.au/pdf2/vici09_194.pdf</a>
AT6S	High temperature Multiposition Air actuators 6 position use w/Existing standoff	Actuators,Accessories	\$465.00	<a href="http://www.chromtech.net.au/pdf2/vici09_194.pdf">http://www.chromtech.net.au/pdf2/vici09_194.pdf</a>
AT64	High temperature Multiposition Air actuators 6 position w/4in standoff assembly	Actuators,Accessories	\$516.00	<a href="http://www.chromtech.net.au/pdf2/vici09_194.pdf">http://www.chromtech.net.au/pdf2/vici09_194.pdf</a>
AT6	High temperature Multiposition Air actuators 6 position w/+losemount assembly	Actuators,Accessories	\$516.00	<a href="http://www.chromtech.net.au/pdf2/vici09_194.pdf">http://www.chromtech.net.au/pdf2/vici09_194.pdf</a>
AT4S	High temperature Multiposition Air actuators 4 position use w/Existing standoff	Actuators,Accessories	\$465.00	<a href="http://www.chromtech.net.au/pdf2/vici09_194.pdf">http://www.chromtech.net.au/pdf2/vici09_194.pdf</a>
A4	Standard Multiposition Valve Air actuators 4 position w/closemount assembly	Actuators,Accessories	\$516.00	<a href="http://www.chromtech.net.au/pdf2/vici09_194.pdf">http://www.chromtech.net.au/pdf2/vici09_194.pdf</a>
AT124	High temperature Multiposition Air actuators 12 position w/4in standoff assembly	Actuators,Accessories	\$516.00	<a href="http://www.chromtech.net.au/pdf2/vici09_194.pdf">http://www.chromtech.net.au/pdf2/vici09_194.pdf</a>
AT8	High temperature Multiposition Air actuators 8 position w/+losemount assembly	Actuators,Accessories	\$516.00	<a href="http://www.chromtech.net.au/pdf2/vici09_194.pdf">http://www.chromtech.net.au/pdf2/vici09_194.pdf</a>
E362	Standard Valve electric actuator 2-position valves 10 ports 36 deg rotation w/ 2in standoff assembly 110v	Actuators,Accessories	\$726.00	<a href="http://www.chromtech.net.au/pdf2/vici09_193.pdf">http://www.chromtech.net.au/pdf2/vici09_193.pdf</a>
E36	Standard Valve electric actuator 2-position valves 10 ports 36 deg rotation w/ closemount assembly 110v	Actuators,Accessories	\$726.00	<a href="http://www.chromtech.net.au/pdf2/vici09_193.pdf">http://www.chromtech.net.au/pdf2/vici09_193.pdf</a>
E45S	Standard Valve electric actuator 2-position valves 8 ports 45 deg rotation use w/existing standoff 110v	Actuators,Accessories	\$675.00	<a href="http://www.chromtech.net.au/pdf2/vici09_193.pdf">http://www.chromtech.net.au/pdf2/vici09_193.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

E452	Standard Valve electric actuator 2-position valves 8 ports 45 deg rotation w/ 2in standoff assembly 110v	Actuators,Accessories	\$726.00	<a href="http://www.chromtech.net.au/pdf2/vici09_193.pdf">http://www.chromtech.net.au/pdf2/vici09_193.pdf</a>
E4X22	Standard Electric actuator Multiposition valves 4x2 position valves w/ 2in standoff assembly 110v	Actuators,Accessories	\$1,032.00	<a href="http://www.chromtech.net.au/pdf2/vici09_193.pdf">http://www.chromtech.net.au/pdf2/vici09_193.pdf</a>
E6X22	Standard Electric actuator Multiposition valves 6x2 position valves w/ 2in standoff assembly 110v	Actuators,Accessories	\$1,032.00	<a href="http://www.chromtech.net.au/pdf2/vici09_193.pdf">http://www.chromtech.net.au/pdf2/vici09_193.pdf</a>
E16S	Standard Electric actuator Multiposition valves 16 position valves use w/existing standoff	Actuators,Accessories	\$981.00	<a href="http://www.chromtech.net.au/pdf2/vici09_193.pdf">http://www.chromtech.net.au/pdf2/vici09_193.pdf</a>
E162	Standard Electric actuator Multiposition valves 16 position valves w/ 2in standoff assembly 110v	Actuators,Accessories	\$1,032.00	<a href="http://www.chromtech.net.au/pdf2/vici09_193.pdf">http://www.chromtech.net.au/pdf2/vici09_193.pdf</a>
E16	Standard Electric actuator Multiposition valves 16 position valves w/closemount assembly 110v	Actuators,Accessories	\$1,032.00	<a href="http://www.chromtech.net.au/pdf2/vici09_193.pdf">http://www.chromtech.net.au/pdf2/vici09_193.pdf</a>
E12S	Standard Electric actuator Multiposition valves 12 position valves use w/existing standoff	Actuators,Accessories	\$981.00	<a href="http://www.chromtech.net.au/pdf2/vici09_193.pdf">http://www.chromtech.net.au/pdf2/vici09_193.pdf</a>
E122	Standard Electric actuator Multiposition valves 12 position valves w/ 2in standoff assembly 110v	Actuators,Accessories	\$1,032.00	<a href="http://www.chromtech.net.au/pdf2/vici09_193.pdf">http://www.chromtech.net.au/pdf2/vici09_193.pdf</a>
E12	Standard Electric actuator Multiposition valves 12 position valves w/closemount assembly 110v	Actuators,Accessories	\$1,032.00	<a href="http://www.chromtech.net.au/pdf2/vici09_193.pdf">http://www.chromtech.net.au/pdf2/vici09_193.pdf</a>
E10S	Standard Electric actuator Multiposition valves 10 position valves use w/existing standoff	Actuators,Accessories	\$981.00	<a href="http://www.chromtech.net.au/pdf2/vici09_193.pdf">http://www.chromtech.net.au/pdf2/vici09_193.pdf</a>
E102	Standard Electric actuator Multiposition valves 10 position valves w/ 2in standoff assembly 110v	Actuators,Accessories	\$1,032.00	<a href="http://www.chromtech.net.au/pdf2/vici09_193.pdf">http://www.chromtech.net.au/pdf2/vici09_193.pdf</a>
E10	Standard Electric actuator Multiposition valves 10 position valves w/ closemount assembly 110v	Actuators,Accessories	\$1,032.00	<a href="http://www.chromtech.net.au/pdf2/vici09_193.pdf">http://www.chromtech.net.au/pdf2/vici09_193.pdf</a>
E8S	Standard Electric actuator Multiposition valves 8 position valves use w/existing standoff	Actuators,Accessories	\$981.00	<a href="http://www.chromtech.net.au/pdf2/vici09_193.pdf">http://www.chromtech.net.au/pdf2/vici09_193.pdf</a>
E82	Standard Electric actuator Multiposition valves 8 position valves w/ 2in standoff assembly 110v	Actuators,Accessories	\$1,032.00	<a href="http://www.chromtech.net.au/pdf2/vici09_193.pdf">http://www.chromtech.net.au/pdf2/vici09_193.pdf</a>
E42	Standard Electric actuator Multiposition valves 4 position valves w/ 2in standoff assembly 110v	Actuators,Accessories	\$1,032.00	<a href="http://www.chromtech.net.au/pdf2/vici09_193.pdf">http://www.chromtech.net.au/pdf2/vici09_193.pdf</a>
E6X2S	Standard Electric actuator Multiposition valves 6x2 position valves use w/existing standoff	Actuators,Accessories	\$981.00	<a href="http://www.chromtech.net.au/pdf2/vici09_193.pdf">http://www.chromtech.net.au/pdf2/vici09_193.pdf</a>
E36S	Standard Valve electric actuator 2-position valves 10 ports 36 deg rotation use w/existing standoff	Actuators,Accessories	\$675.00	<a href="http://www.chromtech.net.au/pdf2/vici09_193.pdf">http://www.chromtech.net.au/pdf2/vici09_193.pdf</a>
E6X2	Standard Electric actuator Multiposition valves 6x2 position valves w/ +losemount assembly 110v	Actuators,Accessories	\$1,032.00	<a href="http://www.chromtech.net.au/pdf2/vici09_193.pdf">http://www.chromtech.net.au/pdf2/vici09_193.pdf</a>

VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

E6S	Standard Electric actuator Multiposition valves 6 position valves use w/existing standoff	Actuators,Accessories	\$981.00	<a href="http://www.chromtech.net.au/pdf2/vici09_193.pdf">http://www.chromtech.net.au/pdf2/vici09_193.pdf</a>
E62	Standard Electric actuator Multiposition valves 6 position valves w/ 2in standoff assembly 110v	Actuators,Accessories	\$1,032.00	<a href="http://www.chromtech.net.au/pdf2/vici09_193.pdf">http://www.chromtech.net.au/pdf2/vici09_193.pdf</a>
E6	Standard Electric actuator Multiposition valves 6 position valves w/ +losemount assembly 110v	Actuators,Accessories	\$1,032.00	<a href="http://www.chromtech.net.au/pdf2/vici09_193.pdf">http://www.chromtech.net.au/pdf2/vici09_193.pdf</a>
E4X2S	Standard Electric actuator Multiposition valves 4x2 position valves use w/existing standoff	Actuators,Accessories	\$981.00	<a href="http://www.chromtech.net.au/pdf2/vici09_193.pdf">http://www.chromtech.net.au/pdf2/vici09_193.pdf</a>
E4	Standard Electric actuator Multiposition valves 4 position valves w/ +losemount assembly 110v	Actuators,Accessories	\$1,032.00	<a href="http://www.chromtech.net.au/pdf2/vici09_193.pdf">http://www.chromtech.net.au/pdf2/vici09_193.pdf</a>
E4X2	Standard Electric actuator Multiposition valves 4x2 position valves w/ +losemount assembly 110v	Actuators,Accessories	\$1,032.00	<a href="http://www.chromtech.net.au/pdf2/vici09_193.pdf">http://www.chromtech.net.au/pdf2/vici09_193.pdf</a>
E45	Standard Valve electric actuator 2-position valves 8 ports 45 deg rotation w/ closemount assembly 110v	Actuators,Accessories	\$726.00	<a href="http://www.chromtech.net.au/pdf2/vici09_193.pdf">http://www.chromtech.net.au/pdf2/vici09_193.pdf</a>
E30S	Standard Valve electric actuator 2-position valves 12 ports 30 deg rotation use w/existing standoff	Actuators,Accessories	\$675.00	<a href="http://www.chromtech.net.au/pdf2/vici09_193.pdf">http://www.chromtech.net.au/pdf2/vici09_193.pdf</a>
E302	Standard Valve electric actuator 2-position valves 12 ports 30 deg rotation w/ 2in standoff assembly 110v	Actuators,Accessories	\$726.00	<a href="http://www.chromtech.net.au/pdf2/vici09_193.pdf">http://www.chromtech.net.au/pdf2/vici09_193.pdf</a>
E30	Standard Valve electric actuator 2-position valves 12 ports 30 deg rotation closemount assembly 110v	Actuators,Accessories	\$726.00	<a href="http://www.chromtech.net.au/pdf2/vici09_193.pdf">http://www.chromtech.net.au/pdf2/vici09_193.pdf</a>
E8	Standard Electric actuator Multiposition valves 8 position valves w/ +losemount assembly 110v	Actuators,Accessories	\$1,032.00	<a href="http://www.chromtech.net.au/pdf2/vici09_193.pdf">http://www.chromtech.net.au/pdf2/vici09_193.pdf</a>
E602	Standard Valve electric actuator 2-position valves 6 ports 60 deg rotation w/ 2in standoff assembly 110v	Actuators,Accessories	\$726.00	<a href="http://www.chromtech.net.au/pdf2/vici09_193.pdf">http://www.chromtech.net.au/pdf2/vici09_193.pdf</a>
E60	Standard Valve electric actuator 2-position valves 6 ports 60 deg rotation w/ closemount assembly 110v	Actuators,Accessories	\$726.00	<a href="http://www.chromtech.net.au/pdf2/vici09_193.pdf">http://www.chromtech.net.au/pdf2/vici09_193.pdf</a>
E90S	Standard Valve electric actuator 2-position valves 3-4 ports 90 deg rotation use w/existing standoff 110v	Actuators,Accessories	\$675.00	<a href="http://www.chromtech.net.au/pdf2/vici09_193.pdf">http://www.chromtech.net.au/pdf2/vici09_193.pdf</a>
E902	Standard Valve electric actuator 2-position valves 3-4 ports 90 deg rotation w/ 2in standoff assembly 110v	Actuators,Accessories	\$726.00	<a href="http://www.chromtech.net.au/pdf2/vici09_193.pdf">http://www.chromtech.net.au/pdf2/vici09_193.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

E90	Standard Valve electric actuator 2-position valves 3-4 ports 90 deg rotation w/ closemount assembly 110v	Actuators,Accessories	\$726.00	<a href="http://www.chromtech.net.au/pdf2/vici09_193.pdf">http://www.chromtech.net.au/pdf2/vici09_193.pdf</a>
E4S	Standard Electric actuator Multiposition valves 4 position valves use w/existing standoff	Actuators,Accessories	\$981.00	<a href="http://www.chromtech.net.au/pdf2/vici09_193.pdf">http://www.chromtech.net.au/pdf2/vici09_193.pdf</a>
E60S	Standard Valve electric actuator 2-position valves 6 ports 60 deg rotation use w/existing standoff	Actuators,Accessories	\$675.00	<a href="http://www.chromtech.net.au/pdf2/vici09_193.pdf">http://www.chromtech.net.au/pdf2/vici09_193.pdf</a>
EUDC	OEM Universal Actuators, Medium torque, w/BCD	Actuators,Accessories	\$1,032.00	<a href="http://www.chromtech.net.au/pdf2/vici09_192.pdf">http://www.chromtech.net.au/pdf2/vici09_192.pdf</a>
EUDB	OEM Universal Actuators, Medium torque, w/USB	Actuators,Accessories	\$1,032.00	<a href="http://www.chromtech.net.au/pdf2/vici09_192.pdf">http://www.chromtech.net.au/pdf2/vici09_192.pdf</a>
EUTC	OEM Universal Actuators, High torque,w/BCD	Actuators,Accessories	\$1,059.00	<a href="http://www.chromtech.net.au/pdf2/vici09_192.pdf">http://www.chromtech.net.au/pdf2/vici09_192.pdf</a>
EUTB	OEM Universal Actuators, High torque,w/USB	Actuators,Accessories	\$1,059.00	<a href="http://www.chromtech.net.au/pdf2/vici09_192.pdf">http://www.chromtech.net.au/pdf2/vici09_192.pdf</a>
EUHC	OEM Universal Actuators, High Speed, w/BCD	Actuators,Accessories	\$1,008.00	<a href="http://www.chromtech.net.au/pdf2/vici09_192.pdf">http://www.chromtech.net.au/pdf2/vici09_192.pdf</a>
EUH	OEM Universal Actuators, High Speed, w/out Interface	Actuators,Accessories	\$969.00	<a href="http://www.chromtech.net.au/pdf2/vici09_192.pdf">http://www.chromtech.net.au/pdf2/vici09_192.pdf</a>
EUHB	OEM Universal Actuators, High Speed, w/USB	Actuators,Accessories	\$1,008.00	<a href="http://www.chromtech.net.au/pdf2/vici09_192.pdf">http://www.chromtech.net.au/pdf2/vici09_192.pdf</a>
EUTA	OEM Universal Actuators, Medium torque, w/out Interface	Actuators,Accessories	\$1,059.00	<a href="http://www.chromtech.net.au/pdf2/vici09_192.pdf">http://www.chromtech.net.au/pdf2/vici09_192.pdf</a>
EUDA	OEM Universal Actuators, High Speed, Medium torque	Actuators,Accessories	\$1,032.00	<a href="http://www.chromtech.net.au/pdf2/vici09_192.pdf">http://www.chromtech.net.au/pdf2/vici09_192.pdf</a>
EUHA	OEM Universal Actuators, High Speed, w/RS-232/485	Actuators,Accessories	\$1,008.00	<a href="http://www.chromtech.net.au/pdf2/vici09_192.pdf">http://www.chromtech.net.au/pdf2/vici09_192.pdf</a>
EUT	OEM Universal Actuators, High torque, w/out Interface	Actuators,Accessories	\$1,020.00	<a href="http://www.chromtech.net.au/pdf2/vici09_192.pdf">http://www.chromtech.net.au/pdf2/vici09_192.pdf</a>
EUD	OEM Universal Actuators, Medium torque, w/out Interface	Actuators,Accessories	\$994.00	<a href="http://www.chromtech.net.au/pdf2/vici09_192.pdf">http://www.chromtech.net.au/pdf2/vici09_192.pdf</a>
I-22897-08	Multi-position Microelectric actuator Multi-drop cable for 8 Microelectric actuators	Actuators,Accessories	\$198.00	<a href="http://www.chromtech.net.au/pdf2/vici09_191.pdf">http://www.chromtech.net.au/pdf2/vici09_191.pdf</a>
I-22897-05	Multi-position Microelectric actuator Multi-drop cable for 5 Microelectric actuators	Actuators,Accessories	\$141.00	<a href="http://www.chromtech.net.au/pdf2/vici09_191.pdf">http://www.chromtech.net.au/pdf2/vici09_191.pdf</a>
I-22897-03	Multi-position Microelectric actuator Multi-drop cable for3 Microelectric actuators	Actuators,Accessories	\$111.00	<a href="http://www.chromtech.net.au/pdf2/vici09_191.pdf">http://www.chromtech.net.au/pdf2/vici09_191.pdf</a>
EMT	Multi-position Microelectric actuator High torque actuator closemount assembly	Actuators,Accessories	\$1,314.00	<a href="http://www.chromtech.net.au/pdf2/vici09_191.pdf">http://www.chromtech.net.au/pdf2/vici09_191.pdf</a>
I-22897-04	Multi-position Microelectric actuator Multi-drop cable for 4 Microelectric actuators	Actuators,Accessories	\$123.00	<a href="http://www.chromtech.net.au/pdf2/vici09_191.pdf">http://www.chromtech.net.au/pdf2/vici09_191.pdf</a>
I-22897-06	Multi-position Microelectric actuator Multi-drop cable for 6 Microelectric actuators	Actuators,Accessories	\$153.00	<a href="http://www.chromtech.net.au/pdf2/vici09_191.pdf">http://www.chromtech.net.au/pdf2/vici09_191.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

EMTS	Multi-position Microelectric actuator High torque actuator use w/existing standoff	Actuators,Accessories	\$1,263.00	<a href="http://www.chromtech.net.au/pdf2/vici09_191.pdf">http://www.chromtech.net.au/pdf2/vici09_191.pdf</a>
EMH	Multi-position Microelectric actuator High speed actuator closamount assembly	Actuators,Accessories	\$1,236.00	<a href="http://www.chromtech.net.au/pdf2/vici09_191.pdf">http://www.chromtech.net.au/pdf2/vici09_191.pdf</a>
EMH2	Multi-position Microelectric actuator High speed actuator 2in standoff assembly	Actuators,Accessories	\$1,236.00	<a href="http://www.chromtech.net.au/pdf2/vici09_191.pdf">http://www.chromtech.net.au/pdf2/vici09_191.pdf</a>
EMHS	Multi-position Microelectric actuator High speed actuator use w/existing standoff	Actuators,Accessories	\$1,185.00	<a href="http://www.chromtech.net.au/pdf2/vici09_191.pdf">http://www.chromtech.net.au/pdf2/vici09_191.pdf</a>
I-22897-02	Multi-position Microelectric actuator Multi-drop cable for 2 Microelectric actuators	Actuators,Accessories	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_191.pdf">http://www.chromtech.net.au/pdf2/vici09_191.pdf</a>
EMT2	Multi-position Microelectric actuator High torque actuator 2in standoff assembly	Actuators,Accessories	\$1,314.00	<a href="http://www.chromtech.net.au/pdf2/vici09_191.pdf">http://www.chromtech.net.au/pdf2/vici09_191.pdf</a>
EU.5L	1/32in Ultra low mass external union -Bore 0.50mm	Valco Fittings	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_19.pdf">http://www.chromtech.net.au/pdf2/vici09_19.pdf</a>
EU.5	1/32in Ultra low mass external union -Bore 0.25mm	Valco Fittings	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_19.pdf">http://www.chromtech.net.au/pdf2/vici09_19.pdf</a>
EU.5T	1/32in Ultra low mass external union -Bore 1/32in	Valco Fittings	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_19.pdf">http://www.chromtech.net.au/pdf2/vici09_19.pdf</a>
SEAL1-HP	Injector nut system HP5890- 6890 Series 1 and 11 -Replacement Seal Disk	Valco Fittings	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_19.pdf">http://www.chromtech.net.au/pdf2/vici09_19.pdf</a>
FSZNA-HP	Injector nut system HP5890- 6890 Series 1 and 11 -Replacement Nut	Valco Fittings	\$51.00	<a href="http://www.chromtech.net.au/pdf2/vici09_19.pdf">http://www.chromtech.net.au/pdf2/vici09_19.pdf</a>
FSZA-HP	Injector nut system HP5890 and 6890 Series 1 and 11 -Nut and seal disk	Valco Fittings	\$80.00	<a href="http://www.chromtech.net.au/pdf2/vici09_19.pdf">http://www.chromtech.net.au/pdf2/vici09_19.pdf</a>
I-22897-04	Multi-position Microelectric actuator Multi-drop cable for 4 Microelectric actuators	Actuators,Accessories	\$123.00	<a href="http://www.chromtech.net.au/pdf2/vici09_189.pdf">http://www.chromtech.net.au/pdf2/vici09_189.pdf</a>
EDS	2-position valves Microelectric actuator High torque actuator w/existing standoff	Actuators,Accessories	\$969.00	<a href="http://www.chromtech.net.au/pdf2/vici09_189.pdf">http://www.chromtech.net.au/pdf2/vici09_189.pdf</a>
ET2	2-position valves Microelectric actuator High torque actuator 2in standoff assembly	Actuators,Accessories	\$1,020.00	<a href="http://www.chromtech.net.au/pdf2/vici09_189.pdf">http://www.chromtech.net.au/pdf2/vici09_189.pdf</a>
ET2	2-position valves Microelectric actuator High torque actuator 2in standoff assembly	Actuators,Accessories	\$1,020.00	<a href="http://www.chromtech.net.au/pdf2/vici09_189.pdf">http://www.chromtech.net.au/pdf2/vici09_189.pdf</a>
I-22897-03	Multi-position Microelectric actuator Multi-drop cable for 2 Microelectric actuators	Actuators,Accessories	\$111.00	<a href="http://www.chromtech.net.au/pdf2/vici09_189.pdf">http://www.chromtech.net.au/pdf2/vici09_189.pdf</a>
I-22897-08	Multi-position Microelectric actuator Multi-drop cable for 8 Microelectric actuators	Actuators,Accessories	\$198.00	<a href="http://www.chromtech.net.au/pdf2/vici09_189.pdf">http://www.chromtech.net.au/pdf2/vici09_189.pdf</a>
ED2	2-position valves Microelectric actuator High torque actuator w/ 2in standoff assembly	Actuators,Accessories	\$1,020.00	<a href="http://www.chromtech.net.au/pdf2/vici09_189.pdf">http://www.chromtech.net.au/pdf2/vici09_189.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

I-22897-05	Multi-position Microelectric actuator Multi-drop cable for 5 Microelectric actuators	Actuators,Accessories	\$141.00	<a href="http://www.chromtech.net.au/pdf2/vici09_189.pdf">http://www.chromtech.net.au/pdf2/vici09_189.pdf</a>
ETS	2-position valves Microelectric actuator High torque actuator use w/existing standoff	Actuators,Accessories	\$969.00	<a href="http://www.chromtech.net.au/pdf2/vici09_189.pdf">http://www.chromtech.net.au/pdf2/vici09_189.pdf</a>
I-22897-06	Multi-position Microelectric actuator Multi-drop cable for 6 Microelectric actuators	Actuators,Accessories	\$153.00	<a href="http://www.chromtech.net.au/pdf2/vici09_189.pdf">http://www.chromtech.net.au/pdf2/vici09_189.pdf</a>
EQS	2-position valves Microelectric actuator Highest speed actuator use w/existing standoff	Actuators,Accessories	\$892.00	<a href="http://www.chromtech.net.au/pdf2/vici09_189.pdf">http://www.chromtech.net.au/pdf2/vici09_189.pdf</a>
ED	2-position valves Microelectric actuator High torque actuator closamount assembly	Actuators,Accessories	\$1,020.00	<a href="http://www.chromtech.net.au/pdf2/vici09_189.pdf">http://www.chromtech.net.au/pdf2/vici09_189.pdf</a>
EQ2	2-position valves Microelectric actuator Highest speed actuator 2in standoff assembly	Actuators,Accessories	\$944.00	<a href="http://www.chromtech.net.au/pdf2/vici09_189.pdf">http://www.chromtech.net.au/pdf2/vici09_189.pdf</a>
EH	2-position valves Microelectric actuator High speed actuator closamount assembly	Actuators,Accessories	\$944.00	<a href="http://www.chromtech.net.au/pdf2/vici09_189.pdf">http://www.chromtech.net.au/pdf2/vici09_189.pdf</a>
EH2	2-position valves Microelectric actuator High speed actuator 2in standoff assembly	Actuators,Accessories	\$944.00	<a href="http://www.chromtech.net.au/pdf2/vici09_189.pdf">http://www.chromtech.net.au/pdf2/vici09_189.pdf</a>
EHS	2-position valves Microelectric actuator High speed actuator use w/existing standoff	Actuators,Accessories	\$892.00	<a href="http://www.chromtech.net.au/pdf2/vici09_189.pdf">http://www.chromtech.net.au/pdf2/vici09_189.pdf</a>
I-22897-02	Multi-position Microelectric actuator Multi-drop cable for 5 Microelectric actuators	Actuators,Accessories	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_189.pdf">http://www.chromtech.net.au/pdf2/vici09_189.pdf</a>
EP	2-position valves Microelectric actuator Medium torque actuator closamount assembly	Actuators,Accessories	\$944.00	<a href="http://www.chromtech.net.au/pdf2/vici09_189.pdf">http://www.chromtech.net.au/pdf2/vici09_189.pdf</a>
EP2	2-position valves Microelectric actuator Medium torque actuator 2in standoff assembly	Actuators,Accessories	\$944.00	<a href="http://www.chromtech.net.au/pdf2/vici09_189.pdf">http://www.chromtech.net.au/pdf2/vici09_189.pdf</a>
EPS	2-position valves Microelectric actuator Medium torque actuator use w/existing standoff	Actuators,Accessories	\$892.00	<a href="http://www.chromtech.net.au/pdf2/vici09_189.pdf">http://www.chromtech.net.au/pdf2/vici09_189.pdf</a>
EQ	2-position valves Microelectric actuator Highest speed actuator closamount assembly	Actuators,Accessories	\$944.00	<a href="http://www.chromtech.net.au/pdf2/vici09_189.pdf">http://www.chromtech.net.au/pdf2/vici09_189.pdf</a>
C65-3186I	1/16inx0.75mm ports,PPS Stator, Integrated Motor/Stream Selector,6 Position,w/Integrated actuator	Cheminert Valves	\$855.00	<a href="http://www.chromtech.net.au/pdf2/vici09_185.pdf">http://www.chromtech.net.au/pdf2/vici09_185.pdf</a>
C65-3188I	1/16inx0.75mm ports,PPS Stator, Integrated Motor/Stream Selector,6 Position,w/Integrated actuator	Cheminert Valves	\$975.00	<a href="http://www.chromtech.net.au/pdf2/vici09_185.pdf">http://www.chromtech.net.au/pdf2/vici09_185.pdf</a>
C65-3180I	1/16inx0.75mm ports,PPS Stator, Integrated Motor/Stream Selector,10 Position,w/Integrated actuator	Cheminert Valves	\$1,000.00	<a href="http://www.chromtech.net.au/pdf2/vici09_185.pdf">http://www.chromtech.net.au/pdf2/vici09_185.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C65-3186I-S	1/16inx0.75mm ports,PPS Stator, Integrated Motor/Stream Selector,6 Position,w/Motor,sensor Only	Cheminert Valves	\$663.00	<a href="http://www.chromtech.net.au/pdf2/vici09_185.pdf">http://www.chromtech.net.au/pdf2/vici09_185.pdf</a>
C65-3188I-S	1/16inx0.75mm ports,PPS Stator, Integrated Motor/Stream Selector,8 Position,w/Motor,sensor Only	Cheminert Valves	\$784.00	<a href="http://www.chromtech.net.au/pdf2/vici09_185.pdf">http://www.chromtech.net.au/pdf2/vici09_185.pdf</a>
C65-3184I-S	1/16inx0.75mm ports,PPS Stator, Integrated Motor/Stream Selector,10 Position,w/Motor,sensor Only	Cheminert Valves	\$663.00	<a href="http://www.chromtech.net.au/pdf2/vici09_185.pdf">http://www.chromtech.net.au/pdf2/vici09_185.pdf</a>
C65-3184I	1/16inx0.75mm ports,PPS Stator, Integrated Motor/Stream Selector,10 Position,w/Integrated actuator	Cheminert Valves	\$855.00	<a href="http://www.chromtech.net.au/pdf2/vici09_185.pdf">http://www.chromtech.net.au/pdf2/vici09_185.pdf</a>
C65Z-3180I-S	1/16inx0.75mm ports,PPS Stator, Integrated Motor/Stream Selector,10 Position,w/Motor,sensor Only	Cheminert Valves	\$740.00	<a href="http://www.chromtech.net.au/pdf2/vici09_185.pdf">http://www.chromtech.net.au/pdf2/vici09_185.pdf</a>
C65Z-3186I-S	1/16inx0.75mm ports,PPS Stator, Integrated Motor/Stream Selector,6 Position,w/Motor,sensor Only	Cheminert Valves	\$663.00	<a href="http://www.chromtech.net.au/pdf2/vici09_185.pdf">http://www.chromtech.net.au/pdf2/vici09_185.pdf</a>
C65-6184I	1/8inx1.50mm ports,PPS Stator, Integrated Motor/Stream Selector,4 Position,w/Integrated actuator	Cheminert Valves	\$855.00	<a href="http://www.chromtech.net.au/pdf2/vici09_185.pdf">http://www.chromtech.net.au/pdf2/vici09_185.pdf</a>
C65Z-3184I-S	1/16inx0.75mm ports,PPS Stator, Integrated Motor/Stream Selector,4 Position,w/Motor,sensor Only	Cheminert Valves	\$663.00	<a href="http://www.chromtech.net.au/pdf2/vici09_185.pdf">http://www.chromtech.net.au/pdf2/vici09_185.pdf</a>
C65Z-3180I	1/16inx0.75mm ports,PPS Stator, Integrated Motor/Stream Selector,10 Position,w/Integrated actuator	Cheminert Valves	\$930.00	<a href="http://www.chromtech.net.au/pdf2/vici09_185.pdf">http://www.chromtech.net.au/pdf2/vici09_185.pdf</a>
C65Z-3188I	1/16inx0.75mm ports,PPS Stator, Integrated Motor/Stream Selector,8 Position,w/Integrated actuator	Cheminert Valves	\$906.00	<a href="http://www.chromtech.net.au/pdf2/vici09_185.pdf">http://www.chromtech.net.au/pdf2/vici09_185.pdf</a>
C65Z-3186I	1/16inx0.75mm ports,PPS Stator, Integrated Motor/Stream Selector,6 Position,w/Integrated actuator	Cheminert Valves	\$855.00	<a href="http://www.chromtech.net.au/pdf2/vici09_185.pdf">http://www.chromtech.net.au/pdf2/vici09_185.pdf</a>
C65Z-3184I	1/16inx0.75mm ports,PPS Stator, Integrated Motor/Stream Selector,4 Position,w/Integrated actuator	Cheminert Valves	\$855.00	<a href="http://www.chromtech.net.au/pdf2/vici09_185.pdf">http://www.chromtech.net.au/pdf2/vici09_185.pdf</a>



## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C65-6186I	1/8inx1.50mm ports,PPS Stator, Integrated Motor/Stream Selector,6 Position,w/Integrated actuator	Cheminert Valves	\$855.00	<a href="http://www.chromtech.net.au/pdf2/vici09_185.pdf">http://www.chromtech.net.au/pdf2/vici09_185.pdf</a>
C65Z-3188I-S	1/16inx0.75mm ports,PPS Stator, Integrated Motor/Stream Selector,8 Position,w/Motor,sensor Only	Cheminert Valves	\$714.00	<a href="http://www.chromtech.net.au/pdf2/vici09_185.pdf">http://www.chromtech.net.au/pdf2/vici09_185.pdf</a>
C65-6188I	1/8inx1.50mm ports,PPS Stator, Integrated Motor/Stream Selector,8 Position,w/Integrated actuator	Cheminert Valves	\$975.00	<a href="http://www.chromtech.net.au/pdf2/vici09_185.pdf">http://www.chromtech.net.au/pdf2/vici09_185.pdf</a>
C65-6180I	1/8inx1.50mm ports,PPS Stator, Integrated Motor/Stream Selector,10 Position,w/Integrated actuator	Cheminert Valves	\$1,000.00	<a href="http://www.chromtech.net.au/pdf2/vici09_185.pdf">http://www.chromtech.net.au/pdf2/vici09_185.pdf</a>
C65-6184I-S	1/8inx1.50mm ports,PPS Stator, Integrated Motor/Stream Selector,4 Position,w/Motor,sensor Only	Cheminert Valves	\$663.00	<a href="http://www.chromtech.net.au/pdf2/vici09_185.pdf">http://www.chromtech.net.au/pdf2/vici09_185.pdf</a>
C65-6186I-S	1/8inx1.50mm ports,PPS Stator, Integrated Motor/Stream Selector,6 Position,w/Motor,sensor Only	Cheminert Valves	\$663.00	<a href="http://www.chromtech.net.au/pdf2/vici09_185.pdf">http://www.chromtech.net.au/pdf2/vici09_185.pdf</a>
C65-6188I-S	1/8inx1.50mm ports,PPS Stator, Integrated Motor/Stream Selector,8 Position,w/Motor,sensor Only	Cheminert Valves	\$784.00	<a href="http://www.chromtech.net.au/pdf2/vici09_185.pdf">http://www.chromtech.net.au/pdf2/vici09_185.pdf</a>
C65-6180I-S	1/8inx1.50mm ports,PPS Stator, Integrated Motor/Stream Selector,10 Position,w/Motor,sensor Only	Cheminert Valves	\$810.00	<a href="http://www.chromtech.net.au/pdf2/vici09_185.pdf">http://www.chromtech.net.au/pdf2/vici09_185.pdf</a>
C65-3180I-S	1/16inx0.75mm ports,PPS Stator, Integrated Motor/Stream Selector,10 Position,w/Motor,sensor Only	Cheminert Valves	\$810.00	<a href="http://www.chromtech.net.au/pdf2/vici09_185.pdf">http://www.chromtech.net.au/pdf2/vici09_185.pdf</a>
C55-2000I	1/16inx0.40mm ports,N60 SS Stator, Integrated Stream Selector,10 Position,w/Integrated actuator	Cheminert Valves	\$1,052.00	<a href="http://www.chromtech.net.au/pdf2/vici09_184.pdf">http://www.chromtech.net.au/pdf2/vici09_184.pdf</a>
C55-2004I	1/16inx0.40mm ports,N60 SS Stator, Integrated Stream Selector,4 Position,w/Integrated actuator	Cheminert Valves	\$1,052.00	<a href="http://www.chromtech.net.au/pdf2/vici09_184.pdf">http://www.chromtech.net.au/pdf2/vici09_184.pdf</a>
C55-2008I	1/16inx0.40mm ports,N60 SS Stator, Integrated Stream Selector,8 Position,w/Integrated actuator	Cheminert Valves	\$1,116.00	<a href="http://www.chromtech.net.au/pdf2/vici09_184.pdf">http://www.chromtech.net.au/pdf2/vici09_184.pdf</a>
C55-2000IX	1/16inx0.40mm ports,PAEK Stator, Integrated Stream Selector,10 Position,w/Motor Only	Cheminert Valves	\$924.00	<a href="http://www.chromtech.net.au/pdf2/vici09_184.pdf">http://www.chromtech.net.au/pdf2/vici09_184.pdf</a>
C55-2C00	1/16inx0.40mm ports,N60 SS Stator, Integrated Stream Selector,10 Position,Replace Stator	Cheminert Valves	\$478.00	<a href="http://www.chromtech.net.au/pdf2/vici09_184.pdf">http://www.chromtech.net.au/pdf2/vici09_184.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C55-2C08	1/16inx0.40mm ports,N60 SS Stator, Integrated Stream Selector,8 Position,Replace Stator	Cheminert Valves	\$408.00	<a href="http://www.chromtech.net.au/pdf2/vici09_184.pdf">http://www.chromtech.net.au/pdf2/vici09_184.pdf</a>
C55-2C06	1/16inx0.40mm ports,N60 SS Stator, Integrated Stream Selector,6 Position,Replace Stator	Cheminert Valves	\$338.00	<a href="http://www.chromtech.net.au/pdf2/vici09_184.pdf">http://www.chromtech.net.au/pdf2/vici09_184.pdf</a>
C55-2C04	1/16inx0.40mm ports,N60 SS Stator, Integrated Stream Selector,4 Position,Replace Stator	Cheminert Valves	\$338.00	<a href="http://www.chromtech.net.au/pdf2/vici09_184.pdf">http://www.chromtech.net.au/pdf2/vici09_184.pdf</a>
C5-20R0H	1/16inx0.40mm ports,N60 SS Stator, Integrated Stream Selector,10 Position, Replace Rotor	Cheminert Valves	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_184.pdf">http://www.chromtech.net.au/pdf2/vici09_184.pdf</a>
C5-20R8H	1/16inx0.40mm ports,PAEK Stator, Integrated Stream Selector,8 Position, Replace Rotor	Cheminert Valves	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_184.pdf">http://www.chromtech.net.au/pdf2/vici09_184.pdf</a>
C55-2C40	1/16inx0.40mm ports,PAEK Stator, Integrated Stream Selector,10 Position,Replace Stator	Cheminert Valves	\$618.00	<a href="http://www.chromtech.net.au/pdf2/vici09_184.pdf">http://www.chromtech.net.au/pdf2/vici09_184.pdf</a>
C5-20R4	1/16inx0.40mm ports,N60 SS Stator, Integrated Stream Selector,4 Position, Replace Rotor	Cheminert Valves	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_184.pdf">http://www.chromtech.net.au/pdf2/vici09_184.pdf</a>
C55-2348I	1/16inx0.40mm ports,PAEK Stator, Integrated Stream Selector,8 Position,w/Integrated actuator	Cheminert Valves	\$1,244.00	<a href="http://www.chromtech.net.au/pdf2/vici09_184.pdf">http://www.chromtech.net.au/pdf2/vici09_184.pdf</a>
C55-2008IX	1/16inx0.40mm ports,N60 SS Stator, Integrated Stream Selector,8 Position,w/Motor Only	Cheminert Valves	\$861.00	<a href="http://www.chromtech.net.au/pdf2/vici09_184.pdf">http://www.chromtech.net.au/pdf2/vici09_184.pdf</a>
C55-2000IX	1/16inx0.40mm ports,N60 SS Stator, Integrated Stream Selector,10 Position,w/Motor Only	Cheminert Valves	\$796.00	<a href="http://www.chromtech.net.au/pdf2/vici09_184.pdf">http://www.chromtech.net.au/pdf2/vici09_184.pdf</a>
C55-2004IX	1/16inx0.40mm ports,N60 SS Stator, Integrated Stream Selector,4 Position,w/Motor Only	Cheminert Valves	\$796.00	<a href="http://www.chromtech.net.au/pdf2/vici09_184.pdf">http://www.chromtech.net.au/pdf2/vici09_184.pdf</a>
C55-2004I-S	1/16inx0.40mm ports,N60 SS Stator, Integrated Stream Selector,4 Position,w/Motor,sensor Only	Cheminert Valves	\$861.00	<a href="http://www.chromtech.net.au/pdf2/vici09_184.pdf">http://www.chromtech.net.au/pdf2/vici09_184.pdf</a>
C55-2000I-S	1/16inx0.40mm ports,PAEK Stator, Integrated Stream Selector,10 Position,w/Motor,sensor Only	Cheminert Valves	\$988.00	<a href="http://www.chromtech.net.au/pdf2/vici09_184.pdf">http://www.chromtech.net.au/pdf2/vici09_184.pdf</a>
C5-23R0H	1/16inx0.40mm ports,PAEK Stator, Integrated Stream Selector,10 Position, Replace Rotor	Cheminert Valves	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_184.pdf">http://www.chromtech.net.au/pdf2/vici09_184.pdf</a>
C55-2008I-S	1/16inx0.40mm ports,N60 SS Stator, Integrated Stream Selector,8 Position,w/Motor,sensor Only	Cheminert Valves	\$924.00	<a href="http://www.chromtech.net.au/pdf2/vici09_184.pdf">http://www.chromtech.net.au/pdf2/vici09_184.pdf</a>
C5-20R6	1/16inx0.40mm ports,N60 SS Stator, Integrated Stream Selector,6 Position, Replace Rotor	Cheminert Valves	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_184.pdf">http://www.chromtech.net.au/pdf2/vici09_184.pdf</a>
C5-23R8H	1/16inx0.40mm ports,PAEK Stator, Integrated Stream Selector,8 Position, Replace Rotor	Cheminert Valves	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_184.pdf">http://www.chromtech.net.au/pdf2/vici09_184.pdf</a>
C5-23R6	1/16inx0.40mm ports,PAEK Stator, Integrated Stream Selector,6 Position, Replace Rotor	Cheminert Valves	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_184.pdf">http://www.chromtech.net.au/pdf2/vici09_184.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C55-2C46	1/16inx0.40mm ports,PAEK Stator, Integrated Stream Selector,6 Position,Replace Stator	Cheminert Valves	\$478.00	<a href="http://www.chromtech.net.au/pdf2/vici09_184.pdf">http://www.chromtech.net.au/pdf2/vici09_184.pdf</a>
C5-23R4	1/16inx0.40mm ports,PAEK Stator, Integrated Stream Selector,4 Position, Replace Rotor	Cheminert Valves	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_184.pdf">http://www.chromtech.net.au/pdf2/vici09_184.pdf</a>
C55-2C48	1/16inx0.40mm ports,PAEK Stator, Integrated Stream Selector,8 Position,Replace Stator	Cheminert Valves	\$549.00	<a href="http://www.chromtech.net.au/pdf2/vici09_184.pdf">http://www.chromtech.net.au/pdf2/vici09_184.pdf</a>
C55-2340IX	1/16inx0.40mm ports,PAEK Stator, Integrated Stream Selector,8 Position,w/Motor Only	Cheminert Valves	\$1,052.00	<a href="http://www.chromtech.net.au/pdf2/vici09_184.pdf">http://www.chromtech.net.au/pdf2/vici09_184.pdf</a>
C55-2348IX	1/16inx0.40mm ports,PAEK Stator, Integrated Stream Selector,8 Position,w/Motor Only	Cheminert Valves	\$988.00	<a href="http://www.chromtech.net.au/pdf2/vici09_184.pdf">http://www.chromtech.net.au/pdf2/vici09_184.pdf</a>
C55-2346IX	1/16inx0.40mm ports,PAEK Stator, Integrated Stream Selector,6 Position,w/Motor Only	Cheminert Valves	\$924.00	<a href="http://www.chromtech.net.au/pdf2/vici09_184.pdf">http://www.chromtech.net.au/pdf2/vici09_184.pdf</a>
C55-2344I	1/16inx0.40mm ports,PAEK Stator Integrated Stream Selector,4 Position,w/Integrated actuator	Cheminert Valves	\$1,179.00	<a href="http://www.chromtech.net.au/pdf2/vici09_184.pdf">http://www.chromtech.net.au/pdf2/vici09_184.pdf</a>
C55-2344IX	1/16inx0.40mm ports,PAEK Stator, Integrated Stream Selector,4 Position,w/Motor Only	Cheminert Valves	\$924.00	<a href="http://www.chromtech.net.au/pdf2/vici09_184.pdf">http://www.chromtech.net.au/pdf2/vici09_184.pdf</a>
C55-2346I	1/16inx0.40mm ports,PAEK Stator, Integrated Stream Selector,6 Position,w/Integrated actuator	Cheminert Valves	\$1,179.00	<a href="http://www.chromtech.net.au/pdf2/vici09_184.pdf">http://www.chromtech.net.au/pdf2/vici09_184.pdf</a>
C55-2340I-S	1/16inx0.40mm ports,PAEK Stator, Integrated Stream Selector,8 Position,w/Motor,sensor Only	Cheminert Valves	\$1,116.00	<a href="http://www.chromtech.net.au/pdf2/vici09_184.pdf">http://www.chromtech.net.au/pdf2/vici09_184.pdf</a>
C55-2348I-S	1/16inx0.40mm ports,PAEK Stator, Integrated Stream Selector,8 Position,w/Motor,sensor Only	Cheminert Valves	\$105.00	<a href="http://www.chromtech.net.au/pdf2/vici09_184.pdf">http://www.chromtech.net.au/pdf2/vici09_184.pdf</a>
C55-2000I	1/16inx0.40mm ports,PAEK Stator, Integrated Stream Selector,10 Position,w/Integrated actuator	Cheminert Valves	\$1,179.00	<a href="http://www.chromtech.net.au/pdf2/vici09_184.pdf">http://www.chromtech.net.au/pdf2/vici09_184.pdf</a>
C55-2346I-S	1/16inx0.40mm ports,PAEK Stator, Integrated Stream Selector,6 Position,w/Motor,sensor Only	Cheminert Valves	\$988.00	<a href="http://www.chromtech.net.au/pdf2/vici09_184.pdf">http://www.chromtech.net.au/pdf2/vici09_184.pdf</a>
C55-2344I-S	1/16inx0.40mm ports,PAEK Stator, Integrated Stream Selector,4 Position,w/Motor,sensor Only	Cheminert Valves	\$988.00	<a href="http://www.chromtech.net.au/pdf2/vici09_184.pdf">http://www.chromtech.net.au/pdf2/vici09_184.pdf</a>
C55-2340I	1/16inx0.40mm ports,PAEK Stator, Integrated Stream Selector,8 Position,w/Integrated actuator	Cheminert Valves	\$1,306.00	<a href="http://www.chromtech.net.au/pdf2/vici09_184.pdf">http://www.chromtech.net.au/pdf2/vici09_184.pdf</a>
C55-2C44	1/16inx0.40mm ports,PAEK Stator, Integrated Stream Selector,4 Position,Replace Stator	Cheminert Valves	\$478.00	<a href="http://www.chromtech.net.au/pdf2/vici09_184.pdf">http://www.chromtech.net.au/pdf2/vici09_184.pdf</a>
C55-2000I-S	1/16inx0.40mm ports,N60 SS Stator, Integrated Stream Selector,10 Position,w/Motor,sensor Only	Cheminert Valves	\$861.00	<a href="http://www.chromtech.net.au/pdf2/vici09_184.pdf">http://www.chromtech.net.au/pdf2/vici09_184.pdf</a>
C62-3186I-S	1/4-26 1/16in Tubing 0.75mm 6 Ports Integrated Motor/valve,w/Motor,sensor Only	Cheminert Valves	\$618.00	<a href="http://www.chromtech.net.au/pdf2/vici09_183.pdf">http://www.chromtech.net.au/pdf2/vici09_183.pdf</a>
CFSL500TF	PTFE Sample Loops-500ul	Cheminert Valves	\$26.00	<a href="http://www.chromtech.net.au/pdf2/vici09_183.pdf">http://www.chromtech.net.au/pdf2/vici09_183.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C62-6186I-S	1/4-26 1/8in Tubing1.50mm 6 Ports Integrated Motor/valve,w/Motor,sensor Only	Cheminert Valves	\$618.00	<a href="http://www.chromtech.net.au/pdf2/vici09_183.pdf">http://www.chromtech.net.au/pdf2/vici09_183.pdf</a>
CFSL1KFEP	FEP Sample Loops-1ml	Cheminert Valves	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_183.pdf">http://www.chromtech.net.au/pdf2/vici09_183.pdf</a>
CFSL1KTF	PTFE Sample Loops-1ml	Cheminert Valves	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_183.pdf">http://www.chromtech.net.au/pdf2/vici09_183.pdf</a>
CFSL1KPK	PEEK Sample Loops-1ml	Cheminert Valves	\$51.00	<a href="http://www.chromtech.net.au/pdf2/vici09_183.pdf">http://www.chromtech.net.au/pdf2/vici09_183.pdf</a>
CFSL2KFEP	FEP Sample Loops-2ml	Cheminert Valves	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_183.pdf">http://www.chromtech.net.au/pdf2/vici09_183.pdf</a>
CFSL2KTF	FEP Sample Loops-2ml	Cheminert Valves	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_183.pdf">http://www.chromtech.net.au/pdf2/vici09_183.pdf</a>
C62-6184I	1/4-26 1/8in Tubing1.50mm 4 Ports Integrated Motor/valve,w/Motor,sensor,controller	Cheminert Valves	\$810.00	<a href="http://www.chromtech.net.au/pdf2/vici09_183.pdf">http://www.chromtech.net.au/pdf2/vici09_183.pdf</a>
C62-3188I-S	1/4-26 1/16in Tubing 0.75mm 8 Ports Integrated Motor/valve,w/Motor,sensor Only	Cheminert Valves	\$726.00	<a href="http://www.chromtech.net.au/pdf2/vici09_183.pdf">http://www.chromtech.net.au/pdf2/vici09_183.pdf</a>
CFSL500FEP	FEP Sample Loops-250ul	Cheminert Valves	\$26.00	<a href="http://www.chromtech.net.au/pdf2/vici09_183.pdf">http://www.chromtech.net.au/pdf2/vici09_183.pdf</a>
C62-3184I-S	1/4-26 1/16in Tubing 0.75mm 4 Ports Integrated Motor/valve,w/Motor,sensor Only	Cheminert Valves	\$618.00	<a href="http://www.chromtech.net.au/pdf2/vici09_183.pdf">http://www.chromtech.net.au/pdf2/vici09_183.pdf</a>
C62-6184I-S	1/4-26 1/8in Tubing1.50mm 4 Ports Integrated Motor/valve,w/Motor,sensor Only	Cheminert Valves	\$618.00	<a href="http://www.chromtech.net.au/pdf2/vici09_183.pdf">http://www.chromtech.net.au/pdf2/vici09_183.pdf</a>
C62-6180I	1/4-26 1/8in Tubing1.50mm 10 Ports Integrated Motor/valve,w/Motor,sensor,controller	Cheminert Valves	\$944.00	<a href="http://www.chromtech.net.au/pdf2/vici09_183.pdf">http://www.chromtech.net.au/pdf2/vici09_183.pdf</a>
C62-6188I	1/4-26 1/8in Tubing1.50mm 8 Ports Integrated Motor/valve,w/Motor,sensor,controller	Cheminert Valves	\$918.00	<a href="http://www.chromtech.net.au/pdf2/vici09_183.pdf">http://www.chromtech.net.au/pdf2/vici09_183.pdf</a>
C62-3180I-S	1/4-26 1/16in Tubing 0.75mm 10 Ports Integrated Motor/valve,w/Motor,sensor Only	Cheminert Valves	\$753.00	<a href="http://www.chromtech.net.au/pdf2/vici09_183.pdf">http://www.chromtech.net.au/pdf2/vici09_183.pdf</a>
C62-6186I	1/4-26 1/8in Tubing1.50mm 6 Ports Integrated Motor/valve,w/Motor,sensor,controller	Cheminert Valves	\$810.00	<a href="http://www.chromtech.net.au/pdf2/vici09_183.pdf">http://www.chromtech.net.au/pdf2/vici09_183.pdf</a>
CFSL2KPK	PEEK Sample Loops-2ml	Cheminert Valves	\$80.00	<a href="http://www.chromtech.net.au/pdf2/vici09_183.pdf">http://www.chromtech.net.au/pdf2/vici09_183.pdf</a>
C62-3186I	1/4-26 1/16in Tubing 0.75mm 6 Ports Integrated Motor/valve,w/Motor,sensor,controller	Cheminert Valves	\$810.00	<a href="http://www.chromtech.net.au/pdf2/vici09_183.pdf">http://www.chromtech.net.au/pdf2/vici09_183.pdf</a>
CFSL20TF	PTFE Sample Loops-20ul	Cheminert Valves	\$22.00	<a href="http://www.chromtech.net.au/pdf2/vici09_183.pdf">http://www.chromtech.net.au/pdf2/vici09_183.pdf</a>
CFSL20PK	PEEK Sample Loops-20ul	Cheminert Valves	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_183.pdf">http://www.chromtech.net.au/pdf2/vici09_183.pdf</a>
CFSL50FEP	FEP Sample Loops-50ul	Cheminert Valves	\$22.00	<a href="http://www.chromtech.net.au/pdf2/vici09_183.pdf">http://www.chromtech.net.au/pdf2/vici09_183.pdf</a>
CFSL50TF	PTFE Sample Loops-50ul	Cheminert Valves	\$22.00	<a href="http://www.chromtech.net.au/pdf2/vici09_183.pdf">http://www.chromtech.net.au/pdf2/vici09_183.pdf</a>
C62-3180I	1/4-26 1/16in Tubing 0.75mm 10 Ports Integrated Motor/valve,w/Motor,sensor,controller	Cheminert Valves	\$944.00	<a href="http://www.chromtech.net.au/pdf2/vici09_183.pdf">http://www.chromtech.net.au/pdf2/vici09_183.pdf</a>
CFSL500PK	PEEK Sample Loops-500ul	Cheminert Valves	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_183.pdf">http://www.chromtech.net.au/pdf2/vici09_183.pdf</a>
C62-6180I-S	1/4-26 1/8in Tubing1.50mm 10 Ports Integrated Motor/valve,w/Motor,sensor Only	Cheminert Valves	\$753.00	<a href="http://www.chromtech.net.au/pdf2/vici09_183.pdf">http://www.chromtech.net.au/pdf2/vici09_183.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

CFSL20FEP	FEP Sample Loops-20ul	Cheminert Valves	\$22.00	<a href="http://www.chromtech.net.au/pdf2/vici09_183.pdf">http://www.chromtech.net.au/pdf2/vici09_183.pdf</a>
C62-6188I-S	1/4-26 1/8in Tubing1.50mm 8 Ports Integrated Motor/valve,w/Motor,sensor Only	Cheminert Valves	\$726.00	<a href="http://www.chromtech.net.au/pdf2/vici09_183.pdf">http://www.chromtech.net.au/pdf2/vici09_183.pdf</a>
C62-3188I	1/4-26 1/16in Tubing 0.75mm 8 Ports Integrated Motor/valve,w/Motor,sensor,controller	Cheminert Valves	\$918.00	<a href="http://www.chromtech.net.au/pdf2/vici09_183.pdf">http://www.chromtech.net.au/pdf2/vici09_183.pdf</a>
C62-3184I	1/4-26 1/16in Tubing 0.75mm 4 Ports Integrated Motor/valve,w/Motor,sensor,controller	Cheminert Valves	\$810.00	<a href="http://www.chromtech.net.au/pdf2/vici09_183.pdf">http://www.chromtech.net.au/pdf2/vici09_183.pdf</a>
CFSL50PK	PEEK Sample Loops-50ul	Cheminert Valves	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_183.pdf">http://www.chromtech.net.au/pdf2/vici09_183.pdf</a>
CFSL100FEP	FEP Sample Loops-100ul	Cheminert Valves	\$22.00	<a href="http://www.chromtech.net.au/pdf2/vici09_183.pdf">http://www.chromtech.net.au/pdf2/vici09_183.pdf</a>
CFSL100TF	PTFE Sample Loops-100ul	Cheminert Valves	\$22.00	<a href="http://www.chromtech.net.au/pdf2/vici09_183.pdf">http://www.chromtech.net.au/pdf2/vici09_183.pdf</a>
CFSL100PK	PEEK Sample Loops-100ul	Cheminert Valves	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_183.pdf">http://www.chromtech.net.au/pdf2/vici09_183.pdf</a>
CFSL250FEP	FEP Sample Loops-250ul	Cheminert Valves	\$22.00	<a href="http://www.chromtech.net.au/pdf2/vici09_183.pdf">http://www.chromtech.net.au/pdf2/vici09_183.pdf</a>
CFSL250TF	PTFE Sample Loops-250ul	Cheminert Valves	\$22.00	<a href="http://www.chromtech.net.au/pdf2/vici09_183.pdf">http://www.chromtech.net.au/pdf2/vici09_183.pdf</a>
CFSL250PK	PEEK Sample Loops-250ul	Cheminert Valves	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_183.pdf">http://www.chromtech.net.au/pdf2/vici09_183.pdf</a>
C62Z-3184I-S	1/16inx0.75mm Integrated Motor Valve, 4 Port, w/Motor,sensor Only	Cheminert Valves	\$618.00	<a href="http://www.chromtech.net.au/pdf2/vici09_182.pdf">http://www.chromtech.net.au/pdf2/vici09_182.pdf</a>
C62Z-3188I-S	1/16inx0.75mm Integrated Motor Valve, 10 Port, w/Motor,sensor Only	Cheminert Valves	\$657.00	<a href="http://www.chromtech.net.au/pdf2/vici09_182.pdf">http://www.chromtech.net.au/pdf2/vici09_182.pdf</a>
CZSL2KFEP	2ml FEP Sample loop	Cheminert Valves	\$48.00	<a href="http://www.chromtech.net.au/pdf2/vici09_182.pdf">http://www.chromtech.net.au/pdf2/vici09_182.pdf</a>
C62Z-3186I-S	1/16inx0.75mm Integrated Motor Valve, 6 Port, w/Motor,sensor Only	Cheminert Valves	\$618.00	<a href="http://www.chromtech.net.au/pdf2/vici09_182.pdf">http://www.chromtech.net.au/pdf2/vici09_182.pdf</a>
CZSL500FEP	500ul FEP Sample loop	Cheminert Valves	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_182.pdf">http://www.chromtech.net.au/pdf2/vici09_182.pdf</a>
CZSL500TF	500ul PTFE Sample loop	Cheminert Valves	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_182.pdf">http://www.chromtech.net.au/pdf2/vici09_182.pdf</a>
CZSL500PK	500ul PEEK Sample loop	Cheminert Valves	\$48.00	<a href="http://www.chromtech.net.au/pdf2/vici09_182.pdf">http://www.chromtech.net.au/pdf2/vici09_182.pdf</a>
CZSL1KFEP	1ml FEP Sample loop	Cheminert Valves	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_182.pdf">http://www.chromtech.net.au/pdf2/vici09_182.pdf</a>
CZSL250TF	250ul PTFE Sample loop	Cheminert Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_182.pdf">http://www.chromtech.net.au/pdf2/vici09_182.pdf</a>
CZSL1KPK	1ml PEEK Sample loop	Cheminert Valves	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_182.pdf">http://www.chromtech.net.au/pdf2/vici09_182.pdf</a>
CZSL250FEP	250ul FEP Sample loop	Cheminert Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_182.pdf">http://www.chromtech.net.au/pdf2/vici09_182.pdf</a>
CZSL2KTF	2ml PTFE Sample loop	Cheminert Valves	\$48.00	<a href="http://www.chromtech.net.au/pdf2/vici09_182.pdf">http://www.chromtech.net.au/pdf2/vici09_182.pdf</a>
CZSL2KPK	2ml PEEK Sample loop	Cheminert Valves	\$87.00	<a href="http://www.chromtech.net.au/pdf2/vici09_182.pdf">http://www.chromtech.net.au/pdf2/vici09_182.pdf</a>
C62Z-3184I	1/16inx0.75mm Integrated Motor Valve, 4 Port, w/Motor,sensor,controller	Cheminert Valves	\$810.00	<a href="http://www.chromtech.net.au/pdf2/vici09_182.pdf">http://www.chromtech.net.au/pdf2/vici09_182.pdf</a>
C62Z-3186I	1/16inx0.75mm Integrated Motor Valve, 6 Port, w/Motor,sensor,controller	Cheminert Valves	\$810.00	<a href="http://www.chromtech.net.au/pdf2/vici09_182.pdf">http://www.chromtech.net.au/pdf2/vici09_182.pdf</a>
C62Z-3188I	1/16inx0.75mm Integrated Motor Valve, 10 Port, w/Motor,sensor,controller	Cheminert Valves	\$848.00	<a href="http://www.chromtech.net.au/pdf2/vici09_182.pdf">http://www.chromtech.net.au/pdf2/vici09_182.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C62Z-3180I	1/16inx0.75mm Integrated Motor Valve, 10 Port, w/Motor,sensor,controller	Cheminert Valves	\$873.00	<a href="http://www.chromtech.net.au/pdf2/vici09_182.pdf">http://www.chromtech.net.au/pdf2/vici09_182.pdf</a>
CZSL1KTF	1ml PTFE Sample loop	Cheminert Valves	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_182.pdf">http://www.chromtech.net.au/pdf2/vici09_182.pdf</a>
CZSL20FEP	20ul FEP Sample loop	Cheminert Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_182.pdf">http://www.chromtech.net.au/pdf2/vici09_182.pdf</a>
C62Z-3180I-S	1/16inx0.75mm Integrated Motor Valve, 10 Port, w/Motor,sensor Only	Cheminert Valves	\$682.00	<a href="http://www.chromtech.net.au/pdf2/vici09_182.pdf">http://www.chromtech.net.au/pdf2/vici09_182.pdf</a>
CZSL5FEP	5ul FEP Sample loop	Cheminert Valves	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_182.pdf">http://www.chromtech.net.au/pdf2/vici09_182.pdf</a>
CZSL5TF	5ul PTFE Sample loop	Cheminert Valves	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_182.pdf">http://www.chromtech.net.au/pdf2/vici09_182.pdf</a>
CZSL5PK	5ul PEEK Sample loops for C1 C2 and C3 inc fits	Cheminert Valves	\$40.00	<a href="http://www.chromtech.net.au/pdf2/vici09_182.pdf">http://www.chromtech.net.au/pdf2/vici09_182.pdf</a>
CZSL10FEP	100ul FEP Sample loop	Cheminert Valves	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_182.pdf">http://www.chromtech.net.au/pdf2/vici09_182.pdf</a>
CZSL250PK	250ul PEEK Sample loops for C1 C2 and C3 inc fits	Cheminert Valves	\$40.00	<a href="http://www.chromtech.net.au/pdf2/vici09_182.pdf">http://www.chromtech.net.au/pdf2/vici09_182.pdf</a>
CZSL10PK	10ul PEEK Sample loops for C1 C2 and C3 inc fits	Cheminert Valves	\$40.00	<a href="http://www.chromtech.net.au/pdf2/vici09_182.pdf">http://www.chromtech.net.au/pdf2/vici09_182.pdf</a>
CZSL50FEP	50ul FEP Sample loop	Cheminert Valves	\$24.00	<a href="http://www.chromtech.net.au/pdf2/vici09_182.pdf">http://www.chromtech.net.au/pdf2/vici09_182.pdf</a>
CZSL20PK	20ul PEEK Sample loops for C1 C2 and C3 inc fits	Cheminert Valves	\$34.00	<a href="http://www.chromtech.net.au/pdf2/vici09_182.pdf">http://www.chromtech.net.au/pdf2/vici09_182.pdf</a>
CZSL50TF	50ul PTFE Sample loop	Cheminert Valves	\$24.00	<a href="http://www.chromtech.net.au/pdf2/vici09_182.pdf">http://www.chromtech.net.au/pdf2/vici09_182.pdf</a>
CZSL50PK	50ul PEEK Sample loop	Cheminert Valves	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_182.pdf">http://www.chromtech.net.au/pdf2/vici09_182.pdf</a>
CZSL100FEP	100ul FEP Sample loop	Cheminert Valves	\$24.00	<a href="http://www.chromtech.net.au/pdf2/vici09_182.pdf">http://www.chromtech.net.au/pdf2/vici09_182.pdf</a>
CZSL100TF	100ul PTFE Sample loop	Cheminert Valves	\$24.00	<a href="http://www.chromtech.net.au/pdf2/vici09_182.pdf">http://www.chromtech.net.au/pdf2/vici09_182.pdf</a>
CZSL100PK	100ul PEEK Sample loop	Cheminert Valves	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_182.pdf">http://www.chromtech.net.au/pdf2/vici09_182.pdf</a>
CZSL10TF	10ul PTFE Sample loop	Cheminert Valves	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_182.pdf">http://www.chromtech.net.au/pdf2/vici09_182.pdf</a>
CZSL20TF	20ul PEEK Sample loop	Cheminert Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_182.pdf">http://www.chromtech.net.au/pdf2/vici09_182.pdf</a>
C3-2006A	Analytical 20ul 0.40mm ports N60 stainless stator Pneumatic actuator 1/16in fittings	Cheminert Valves	\$879.00	<a href="http://www.chromtech.net.au/pdf2/vici09_181.pdf">http://www.chromtech.net.au/pdf2/vici09_181.pdf</a>
C3-2346	Analytical 20ul 0.40mm ports PAEK stator Manual 1/16in fittings	Cheminert Valves	\$702.00	<a href="http://www.chromtech.net.au/pdf2/vici09_181.pdf">http://www.chromtech.net.au/pdf2/vici09_181.pdf</a>
C3-2006	Analytical 20ul 0.40mm ports N60 stainless stator Manual 1/16in fittings	Cheminert Valves	\$561.00	<a href="http://www.chromtech.net.au/pdf2/vici09_181.pdf">http://www.chromtech.net.au/pdf2/vici09_181.pdf</a>
C2-20R6	Titanium stator Replacement rotor 6 Port 1/16in fittings 0.40mm ports	Cheminert Valves	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_181.pdf">http://www.chromtech.net.au/pdf2/vici09_181.pdf</a>
C3-2C46	Analytical 20ul 0.40mm ports PAEK stator Replacement stator 1/16in fittings	Cheminert Valves	\$490.00	<a href="http://www.chromtech.net.au/pdf2/vici09_181.pdf">http://www.chromtech.net.au/pdf2/vici09_181.pdf</a>
C3-2346A	Analytical 20ul 0.40mm ports PAEK stator Pneumatic actuator 1/16in fittings	Cheminert Valves	\$1,020.00	<a href="http://www.chromtech.net.au/pdf2/vici09_181.pdf">http://www.chromtech.net.au/pdf2/vici09_181.pdf</a>
C3-2006E	Analytical 20ul 0.40mm ports N60 stainless stator Electric actuator 1/16in fittings	Cheminert Valves	\$1,287.00	<a href="http://www.chromtech.net.au/pdf2/vici09_181.pdf">http://www.chromtech.net.au/pdf2/vici09_181.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C3-2346E	Analytical 20ul 0.40mm ports PAEK stator Electric actuator 1/16in fittings	Cheminert Valves	\$1,428.00	<a href="http://www.chromtech.net.au/pdf2/vici09_181.pdf">http://www.chromtech.net.au/pdf2/vici09_181.pdf</a>
C3-2006EH	Analytical 20ul 0.40mm ports N60 stainless stator Microelectric actuator 1/16in fittings	Cheminert Valves	\$1,504.00	<a href="http://www.chromtech.net.au/pdf2/vici09_181.pdf">http://www.chromtech.net.au/pdf2/vici09_181.pdf</a>
C3-2346EH	Analytical 20ul 0.40mm ports PAEK stator Microelectric actuator 1/16in fittings	Cheminert Valves	\$1,644.00	<a href="http://www.chromtech.net.au/pdf2/vici09_181.pdf">http://www.chromtech.net.au/pdf2/vici09_181.pdf</a>
C3-2006D	Analytical 20ul 0.40mm ports N60 stainless stator Replacement valve 1/16in fittings	Cheminert Valves	\$561.00	<a href="http://www.chromtech.net.au/pdf2/vici09_181.pdf">http://www.chromtech.net.au/pdf2/vici09_181.pdf</a>
C3-2346D	Analytical 20ul 0.40mm ports PAEK stator Replacement valve 1/16in fittings	Cheminert Valves	\$702.00	<a href="http://www.chromtech.net.au/pdf2/vici09_181.pdf">http://www.chromtech.net.au/pdf2/vici09_181.pdf</a>
C2-20R6	Titanium stator Replacement rotor 6 Port 1/16in fittings 0.40mm ports	Cheminert Valves	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_181.pdf">http://www.chromtech.net.au/pdf2/vici09_181.pdf</a>
C2V-2C06	1/16inx0.40mm 6 Port Vertical Injector N60 SS Stator,Replace Stator	Cheminert Valves	\$351.00	<a href="http://www.chromtech.net.au/pdf2/vici09_181.pdf">http://www.chromtech.net.au/pdf2/vici09_181.pdf</a>
C3-2C06	Analytical 20ul 0.40mm ports N60 stainless stator Replacement stator 1/16in fittings	Cheminert Valves	\$351.00	<a href="http://www.chromtech.net.au/pdf2/vici09_181.pdf">http://www.chromtech.net.au/pdf2/vici09_181.pdf</a>
C2V-2C46	1/16inx0.40mm 6 Port Vertical Injector PAEKStator,Replace Stator	Cheminert Valves	\$490.00	<a href="http://www.chromtech.net.au/pdf2/vici09_181.pdf">http://www.chromtech.net.au/pdf2/vici09_181.pdf</a>
C2V-2006	1/16inx0.40mm 6 Port Vertical Injector N60 SS Stator,Manual	Cheminert Valves	\$561.00	<a href="http://www.chromtech.net.au/pdf2/vici09_181.pdf">http://www.chromtech.net.au/pdf2/vici09_181.pdf</a>
C2V-2346	1/16inx0.40mm 6 Port Vertical Injector PAEK Stator,Manual	Cheminert Valves	\$702.00	<a href="http://www.chromtech.net.au/pdf2/vici09_181.pdf">http://www.chromtech.net.au/pdf2/vici09_181.pdf</a>
C2V-2006A	1/16inx0.40mm 6 Port Vertical Injector N60 SS Stator,w/Pneumatic actuator	Cheminert Valves	\$879.00	<a href="http://www.chromtech.net.au/pdf2/vici09_181.pdf">http://www.chromtech.net.au/pdf2/vici09_181.pdf</a>
C2V-2346A	1/16inx0.40mm 6 Port Vertical Injector PAEK Stator,w/Pneumatic actuator	Cheminert Valves	\$1,020.00	<a href="http://www.chromtech.net.au/pdf2/vici09_181.pdf">http://www.chromtech.net.au/pdf2/vici09_181.pdf</a>
C2V-2006E	1/16inx0.40mm 6 Port Vertical Injector N60 SS Stator,w/Std Electric actuator	Cheminert Valves	\$1,287.00	<a href="http://www.chromtech.net.au/pdf2/vici09_181.pdf">http://www.chromtech.net.au/pdf2/vici09_181.pdf</a>
C2V-2346E	1/16inx0.40mm 6 Port Vertical Injector PAEK Stator,w/Std Electric actuator	Cheminert Valves	\$1,428.00	<a href="http://www.chromtech.net.au/pdf2/vici09_181.pdf">http://www.chromtech.net.au/pdf2/vici09_181.pdf</a>
C2V-2006EH	1/16inx0.40mm 6 Port Vertical Injector N60 SS Stator,w/Microelectric actuator	Cheminert Valves	\$1,504.00	<a href="http://www.chromtech.net.au/pdf2/vici09_181.pdf">http://www.chromtech.net.au/pdf2/vici09_181.pdf</a>
C2V-2346EH	1/16inx0.40mm 6 Port Vertical Injector N60 SS Stator,w/Microelectric actuator	Cheminert Valves	\$1,644.00	<a href="http://www.chromtech.net.au/pdf2/vici09_181.pdf">http://www.chromtech.net.au/pdf2/vici09_181.pdf</a>
C2V-2006D	1/16inx0.40mm 6 Port Vertical Injector N60 SS Stator,Replace Valve	Cheminert Valves	\$561.00	<a href="http://www.chromtech.net.au/pdf2/vici09_181.pdf">http://www.chromtech.net.au/pdf2/vici09_181.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C2V-2346D	1/16inx0.40mm 6 Port Vertical Injector PAEK Stator,Replace Valve	Cheminert Valves	\$702.00	<a href="http://www.chromtech.net.au/pdf2/vici09_181.pdf">http://www.chromtech.net.au/pdf2/vici09_181.pdf</a>
C2-23R6	Analytical 20ul 0.40mm ports PAEK stator Replacement rotor 1/16in fittings	Cheminert Valves	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_181.pdf">http://www.chromtech.net.au/pdf2/vici09_181.pdf</a>
C2-23R6	Analytical 20ul 0.40mm ports PAEK stator Replacement rotor 1/16in fittings	Cheminert Valves	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_181.pdf">http://www.chromtech.net.au/pdf2/vici09_181.pdf</a>
CZSL500PK	500ul PEEK Sample loop	Cheminert Valves	\$48.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
C52-2340I-S	PAEK Stator 1/16inx0.4mm 8 Port Integrated Motor,w/Motor Sensor Only	Cheminert Valves	\$1,020.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
C52-2346IX	PAEK Stator 1/16inx0.4mm 6 Port Integrated Motor,w/Motor Only	Cheminert Valves	\$828.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
CZSL250PK	250ul PEEK Sample loops for C1 C2 and C3 inc fits	Cheminert Valves	\$40.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
CSL500	500ul 1/16in Stainless steel loops inc fits	Cheminert Valves	\$40.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
C52-2344IX	PAEK Stator 1/16inx0.4mm 4 Port Integrated Motor,w/Motor Only	Cheminert Valves	\$828.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
CSL1K	1ml Stainless Steel Sample loops for C1 C2 and C3 inc fits	Cheminert Valves	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
CZSL1KPK	1ml PEEK Sample loop	Cheminert Valves	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
CSL2K	2ml 1/16in Stainless steel loops inc fits	Cheminert Valves	\$57.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
CSL5K	5ml 1/16in Stainless steel loops inc fits	Cheminert Valves	\$68.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
C52-2348IX	PAEK Stator 1/16inx0.4mm 8 Port Integrated Motor,w/Motor Only	Cheminert Valves	\$892.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
CZSL5KPK	5ml PEEKSample loops for C1 C2 and C3 inc fits	Cheminert Valves	\$135.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
CSL10K	10ml 1/16in Stainless steel loops inc fits	Cheminert Valves	\$87.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
CSL50	50ul Stainless Steel Sample loops for C1 C2 and C3 inc fits	Cheminert Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
CZSL2KPK	2ml PEEK Sample loop	Cheminert Valves	\$87.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
CSL20	20ul 1/16in Stainless steel loops inc fits	Cheminert Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
C52-2008IX	N60 SS Stator 1/16inx0.4mm 8 Port Integrated Motor,w/Motor Only	Cheminert Valves	\$765.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
C52-2346I	PAEK Stator 1/16inx0.4mm 6 Port Integrated Motor,w/Integated actuator	Cheminert Valves	\$1,083.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
C2-20R4	N60 SS stator Replacement rotor 4 Port 1/16in fittings 0.40mm ports	Cheminert Valves	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
C2-20R6	N60 SS stator Replacement rotor 6 Port 1/16in fittings 0.40mm ports	Cheminert Valves	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>



VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C2-20R8H	N60 SS stator Replacement rotor 8 Port 1/16in fittings 0.40mm ports	Cheminert Valves	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
C2-20R0H	N60 SS stator Replacement rotor 4 Port 1/16in fittings 0.40mm ports	Cheminert Valves	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
C52-2C04	N60 SS Stator 1/16inx0.4mm 10 Port Integrated Motor,Replace Stator	Cheminert Valves	\$280.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
C52-2C06	N60 SS Stator 1/16inx0.4mm 10 Port Integrated Motor,Replace Stator	Cheminert Valves	\$280.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
C52-2C08	N60 SS Stator 1/16inx0.4mm 8 Port Integrated Motor,Replace Stator	Cheminert Valves	\$351.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
C52-2004IX	N60 SS Stator 1/16inx0.4mm 4 Port Integrated Motor,w/Motor Only	Cheminert Valves	\$702.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
C52-2C00	N60 SS Stator 1/16inx0.4mm 10 Port Integrated Motor,w/Replace Stator	Cheminert Valves	\$420.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
C52-2006IX	N60 SS Stator 1/16inx0.4mm 6 Port Integrated Motor,w/Motor Only	Cheminert Valves	\$702.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
C52-2344I	PAEK Stator 1/16inx0.4mm 4 Port Integrated Motor,w/Integated actuator	Cheminert Valves	\$1,083.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
CZSL50PK	50ul PEEK Sample loop	Cheminert Valves	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
CSL100	100ul 1/16in Stainless steel loops inc fits	Cheminert Valves	\$40.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
CZSL100PK	100ul PEEK Sample loop	Cheminert Valves	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
CSL250	250ul Stainless Steel Sample loops for C1 C2 and C3 inc fits	Cheminert Valves	\$40.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
C52-2346I-S	PAEK Stator 1/16inx0.4mm 6 Port Integrated Motor,w/Motor Sensor Only	Cheminert Valves	\$892.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
C52-2344I-S	PAEK Stator 1/16inx0.4mm 4 Port Integrated Motor,w/Motor Sensor Only	Cheminert Valves	\$892.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
C52-2340I	PAEK Stator 1/16inx0.4mm 8 Port Integrated Motor,w/Integated actuator	Cheminert Valves	\$1,212.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
C52-2348I-S	PAEK Stator 1/16inx0.4mm 8 Port Integrated Motor,w/Motor Sensor Only	Cheminert Valves	\$957.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
C52-2348I	PAEK Stator 1/16inx0.4mm 8 Port Integrated Motor,w/Integated actuator	Cheminert Valves	\$1,148.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
C52-2340IX	PAEK Stator 1/16inx0.4mm 8 Port Integrated Motor,w/Motor Only	Cheminert Valves	\$957.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
CZSL10PK	10ul PEEK Sample loops for C1 C2 and C3 inc fits	Cheminert Valves	\$40.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C2-23R4	PAEK Stator Replacement rotor 4 Port 1/16in fittings 0.40mm ports	Cheminert Valves	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
C2-23R6	PAEK Stator Replacement rotor 6 Port 1/16in fittings 0.40mm ports	Cheminert Valves	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
C2-23R8	PAEK Stator Replacement rotor 8 Port 1/16in fittings 0.40mm ports	Cheminert Valves	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
C2-23R0H	N60 SS stator Replacement rotor 10 Port 1/16in fittings 0.40mm ports	Cheminert Valves	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
C52-2C44	PAEK Stator 1/16inx0.4mm 4 Port Integrated Motor,Replace Stator	Cheminert Valves	\$420.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
C52-2C46	PAEK Stator 1/16inx0.4mm 6 Port Integrated Motor,Replace Stator	Cheminert Valves	\$420.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
C52-2C48	PAEK Stator 1/16inx0.4mm 8 Port Integrated Motor,Replace Stator	Cheminert Valves	\$490.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
C52-2C40	PAEK Stator 1/16inx0.4mm 10 Port Integrated Motor,Replace Stator	Cheminert Valves	\$561.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
CSL2	2ul Stainless Steel Sample loops for C1 C2 and C3 inc fits	Cheminert Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
CZSL2PK	2ul PEEK Sample loops for C1 C2 and C3 inc fits	Cheminert Valves	\$40.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
CSL5	5ul Stainless Steel Sample loops for C1 C2 and C3 inc fits	Cheminert Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
C52-2000IX	N60 SS Stator 1/16inx0.4mm 10 Port Integrated Motor,w/Motor Only	Cheminert Valves	\$828.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
CSL10	10ul Stainless Steel Sample loops for C1 C2 and C3 inc fits	Cheminert Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
C52-2000I-S	N60 SS Stator 1/16inx0.4mm 10 Port Integrated Motor,w/Motor Sensor Only	Cheminert Valves	\$892.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
CZSL20PK	20ul PEEK Sample loops for C1 C2 and C3 inc fits	Cheminert Valves	\$34.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
C52-2004I	N60 SS Stator 1/16inx0.4mm 4 Port Integrated Motor,w/Integated actuator	Cheminert Valves	\$957.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
C52-2006I	N60 SS Stator 1/16inx0.4mm 6 Port Integrated Motor,w/Integated actuator	Cheminert Valves	\$957.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
C52-2008I	N60 SS Stator 1/16inx0.4mm 8 Port Integrated Motor,w/Integated actuator	Cheminert Valves	\$1,020.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
C52-2000I	N60 SS Stator 1/16inx0.4mm 10 Port Integrated Motor,w/Integated actuator	Cheminert Valves	\$1,083.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C52-2004I-S	N60 SS Stator 1/16inx0.4mm 4 Port Integrated Motor,w/Motor Sensor Only	Cheminert Valves	\$765.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
C52-2006I-S	N60 SS Stator 1/16inx0.4mm 6 Port Integrated Motor,w/Motor Sensor Only	Cheminert Valves	\$765.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
C52-2008I-S	N60 SS Stator 1/16inx0.4mm 8 Port Integrated Motor,w/Motor Sensor Only	Cheminert Valves		<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
CZSL5PK	5ul PEEK Sample loops for C1 C2 and C3 inc fits	Cheminert Valves	\$40.00	<a href="http://www.chromtech.net.au/pdf2/vici09_180.pdf">http://www.chromtech.net.au/pdf2/vici09_180.pdf</a>
FS1.35-R10	Union -Replacement ferrules:Column ID:..45mm	Valco Fittings	\$74.00	<a href="http://www.chromtech.net.au/pdf2/vici09_18.pdf">http://www.chromtech.net.au/pdf2/vici09_18.pdf</a>
FSKZU1	Union	Valco Fittings	\$118.00	<a href="http://www.chromtech.net.au/pdf2/vici09_18.pdf">http://www.chromtech.net.au/pdf2/vici09_18.pdf</a>
FS1.45-R10	Union -Replacement ferrules:Column ID:.32mm	Valco Fittings	\$74.00	<a href="http://www.chromtech.net.au/pdf2/vici09_18.pdf">http://www.chromtech.net.au/pdf2/vici09_18.pdf</a>
FSZN1	Union -Replacement nut	Valco Fittings	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_18.pdf">http://www.chromtech.net.au/pdf2/vici09_18.pdf</a>
FS1.75-R10	Union -Replacement ferrules:Column ID:.53mm	Valco Fittings	\$74.00	<a href="http://www.chromtech.net.au/pdf2/vici09_18.pdf">http://www.chromtech.net.au/pdf2/vici09_18.pdf</a>
FSQ1	Union -Replacement liner	Valco Fittings	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_18.pdf">http://www.chromtech.net.au/pdf2/vici09_18.pdf</a>
C2V-1006EH	Microbore 1/16mmx 0.25mm Vert Port /Injectors N60 SS stator, w/Microelectric actuator	Cheminert Valves	\$1,714.00	<a href="http://www.chromtech.net.au/pdf2/vici09_179.pdf">http://www.chromtech.net.au/pdf2/vici09_179.pdf</a>
C3-1346EH	Microbore 5ul 0.25mm ports PAEK stator Microelectric actuator 1/16in fittings	Cheminert Valves	\$1,983.00	<a href="http://www.chromtech.net.au/pdf2/vici09_179.pdf">http://www.chromtech.net.au/pdf2/vici09_179.pdf</a>
C2V-1006E	Microbore 1/16mmx 0.25mm Vert Port /Injectors N60 SS stator, w/Std Electric actuator	Cheminert Valves	\$1,498.00	<a href="http://www.chromtech.net.au/pdf2/vici09_179.pdf">http://www.chromtech.net.au/pdf2/vici09_179.pdf</a>
C2V-1006D	Microbore 1/16mmx 0.25mm Vert Port /Injectors N60 SS stator, Replace Valve	Cheminert Valves	\$771.00	<a href="http://www.chromtech.net.au/pdf2/vici09_179.pdf">http://www.chromtech.net.au/pdf2/vici09_179.pdf</a>
C2V-1006A	Microbore 1/16mmx 0.25mm Vert Port /Injectors N60 SS stator, w/Pneumatic actuator	Cheminert Valves	\$1,090.00	<a href="http://www.chromtech.net.au/pdf2/vici09_179.pdf">http://www.chromtech.net.au/pdf2/vici09_179.pdf</a>
C2V-1346	Microbore 1/16mmx 0.25mm Vert Port /Injectors PAEK stator, Manual	Cheminert Valves	\$912.00	<a href="http://www.chromtech.net.au/pdf2/vici09_179.pdf">http://www.chromtech.net.au/pdf2/vici09_179.pdf</a>
C2V-1006	Microbore 1/16mmx 0.25mm Vert Port /Injectors N60 SS stator, Manual	Cheminert Valves	\$771.00	<a href="http://www.chromtech.net.au/pdf2/vici09_179.pdf">http://www.chromtech.net.au/pdf2/vici09_179.pdf</a>
C3-1C46	Microbore 5ul 0.25mm ports PAEK stator Replacement stator 1/16in fittings	Cheminert Valves	\$702.00	<a href="http://www.chromtech.net.au/pdf2/vici09_179.pdf">http://www.chromtech.net.au/pdf2/vici09_179.pdf</a>
C3-1C06	Microbore 5ul 0.25mm ports N60 stainless stator Replacement stator 1/16in fittings	Cheminert Valves	\$561.00	<a href="http://www.chromtech.net.au/pdf2/vici09_179.pdf">http://www.chromtech.net.au/pdf2/vici09_179.pdf</a>
C2V-1346EH	Microbore 1/16mmx 0.25mm Vert Port /Injectors PAEK stator, w/Microelectric actuator	Cheminert Valves	\$1,983.00	<a href="http://www.chromtech.net.au/pdf2/vici09_179.pdf">http://www.chromtech.net.au/pdf2/vici09_179.pdf</a>
C3-1006D	Microbore 5ul 0.25mm ports N60 stainless stator Replacement valve 1/16in fittings	Cheminert Valves	\$771.00	<a href="http://www.chromtech.net.au/pdf2/vici09_179.pdf">http://www.chromtech.net.au/pdf2/vici09_179.pdf</a>
C2V-1346E	Microbore 1/16mmx 0.25mm Vert Port /Injectors PAEK stator, w/Std Electric actuator	Cheminert Valves	\$1,638.00	<a href="http://www.chromtech.net.au/pdf2/vici09_179.pdf">http://www.chromtech.net.au/pdf2/vici09_179.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C3-1006EH	Microbore 5ul 0.25mm ports N60 stainless stator Microelectric actuator 1/16in fittings	Cheminert Valves	\$1,714.00	<a href="http://www.chromtech.net.au/pdf2/vici09_179.pdf">http://www.chromtech.net.au/pdf2/vici09_179.pdf</a>
C3-1346E	Microbore 5ul 0.25mm ports PAEK stator Electric actuator 1/16in fittings	Cheminert Valves	\$1,638.00	<a href="http://www.chromtech.net.au/pdf2/vici09_179.pdf">http://www.chromtech.net.au/pdf2/vici09_179.pdf</a>
C3-1006E	Microbore 5ul 0.25mm ports N60 stainless stator Electric actuator 1/16in fittings	Cheminert Valves	\$1,498.00	<a href="http://www.chromtech.net.au/pdf2/vici09_179.pdf">http://www.chromtech.net.au/pdf2/vici09_179.pdf</a>
C3-1346A	Microbore 5ul 0.25mm ports PAEK stator Pneumatic actuator 1/16in fittings	Cheminert Valves	\$1,230.00	<a href="http://www.chromtech.net.au/pdf2/vici09_179.pdf">http://www.chromtech.net.au/pdf2/vici09_179.pdf</a>
C3-1006A	Microbore 5ul 0.25mm ports N60 stainless stator Pneumatic actuator 1/16in fittings	Cheminert Valves	\$1,090.00	<a href="http://www.chromtech.net.au/pdf2/vici09_179.pdf">http://www.chromtech.net.au/pdf2/vici09_179.pdf</a>
C3-1346	Microbore 5ul 0.25mm ports PAEK stator Manual 1/16in fittings	Cheminert Valves	\$912.00	<a href="http://www.chromtech.net.au/pdf2/vici09_179.pdf">http://www.chromtech.net.au/pdf2/vici09_179.pdf</a>
C3-1006	Microbore 5ul 0.25mm ports N60 stainless stator Manual 1/16in fittings	Cheminert Valves	\$771.00	<a href="http://www.chromtech.net.au/pdf2/vici09_179.pdf">http://www.chromtech.net.au/pdf2/vici09_179.pdf</a>
C3-1346D	Microbore 5ul 0.25mm ports PAEK stator Replacement valve 1/16in fittings	Cheminert Valves	\$912.00	<a href="http://www.chromtech.net.au/pdf2/vici09_179.pdf">http://www.chromtech.net.au/pdf2/vici09_179.pdf</a>
C2V-1346D	Microbore 1/16mmx 0.25mm Vert Port /Injectors PAEK stator, Replace Valve	Cheminert Valves	\$912.00	<a href="http://www.chromtech.net.au/pdf2/vici09_179.pdf">http://www.chromtech.net.au/pdf2/vici09_179.pdf</a>
C2V-1346A	Microbore 1/16mmx 0.25mm Vert Port /Injectors PAEK stator, w/Pneumatic actuator	Cheminert Valves	\$1,230.00	<a href="http://www.chromtech.net.au/pdf2/vici09_179.pdf">http://www.chromtech.net.au/pdf2/vici09_179.pdf</a>
C2-10R6	Microbore 1/16mmx 0.25mm Vert Port /Injectors N60 SS stator, Replacement rotor	Cheminert Valves	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_179.pdf">http://www.chromtech.net.au/pdf2/vici09_179.pdf</a>
C2-13R6	Microbore 1/16mmx 0.25mm Vert Port /Injectors PAEK stator, Replacement rotor	Cheminert Valves	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_179.pdf">http://www.chromtech.net.au/pdf2/vici09_179.pdf</a>
C2V-1C06	Microbore 1/16mmx 0.25mm Vert Port /Injectors N60 SS stator, Replace stator	Cheminert Valves	\$561.00	<a href="http://www.chromtech.net.au/pdf2/vici09_179.pdf">http://www.chromtech.net.au/pdf2/vici09_179.pdf</a>
C2V-1C46	Microbore 1/16mmx 0.25mm Vert Port /Injectors PAEK stator, Replace stator	Cheminert Valves	\$702.00	<a href="http://www.chromtech.net.au/pdf2/vici09_179.pdf">http://www.chromtech.net.au/pdf2/vici09_179.pdf</a>
CZSL1KPK	1ml PEEK Sample loop	Cheminert Valves	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
CSL50	50ul Stainless Steel Sample loops for C1 C2 and C3 inc fits	Cheminert Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
CZSL50PK	50ul PEEK Sample loop	Cheminert Valves	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
CSL100	100ul 1/16in Stainless steel loops inc fits	Cheminert Valves	\$40.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
CZSL100PK	100ul PEEK Sample loop	Cheminert Valves	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
CSL500	500ul 1/16in Stainless steel loops inc fits	Cheminert Valves	\$40.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
CZSL10PK	10ul PEEK Sample loops for C1 C2 and C3 inc fits	Cheminert Valves	\$40.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

CSL1K	1ml Stainless Steel Sample loops for C1 C2 and C3 inc fits	Cheminert Valves	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
C52-1348IX	Microbore 1/16inx 0.25mm PAEK stator 8 Port w/Motor Only	Cheminert Valves	\$1,083.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
CSL2K	2ml 1/16in Stainless steel loops inc fits	Cheminert Valves	\$57.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
CZSL2KPK	2ml PEEK Sample loop	Cheminert Valves	\$87.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
CSL5K	5ml 1/16in Stainless steel loops inc fits	Cheminert Valves	\$68.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
CSL10K	10ml 1/16in Stainless steel loops inc fits	Cheminert Valves	\$87.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
C52-1340I-S	Microbore 1/16inx 0.25mm PAEK stator 10 Port w/Motor Sensor Only	Cheminert Valves	\$1,212.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
CZSL500PK	500ul PEEK Sample loop	Cheminert Valves	\$48.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
C52-1C44	Microbore 1/16inx 0.25mm PAEK stator 4 Port Replace stator	Cheminert Valves	\$632.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
CZSL5PK	5ul PEEK Sample loops for C1 C2 and C3 inc fits	Cheminert Valves	\$40.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
CSL5	5ul Stainless Steel Sample loops for C1 C2 and C3 inc fits	Cheminert Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
CZSL2PK	2ul PEEK Sample loops for C1 C2 and C3 inc fits	Cheminert Valves	\$40.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
CSL2	2ul Stainless Steel Sample loops for C1 C2 and C3 inc fits	Cheminert Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
C52-1C40	Microbore 1/16inx 0.25mm PAEK stator 10 Port Replace stator	Cheminert Valves	\$771.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
C52-1344IX	Microbore 1/16inx 0.25mm PAEK stator 4 Port w/Motor Only	Cheminert Valves	\$1,020.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
C52-1C46	Microbore 1/16inx 0.25mm PAEK stator 6 Port Replace stator	Cheminert Valves	\$632.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
C52-1346IX	Microbore 1/16inx 0.25mm PAEK stator 6 Port w/Motor Only	Cheminert Valves	\$1,020.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
C2-13R0	Microbore 1/16inx 0.25mm PAEK stator 10 Port Replacement rotor	Cheminert Valves	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
C2-13R8H	Microbore 1/16inx 0.25mm PAEK stator 8 Port Replacement rotor	Cheminert Valves	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
C2-13R6	Microbore 1/16inx 0.25mm PAEK stator 6 Port Replacement rotor	Cheminert Valves	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
C2-13R4	Microbore 1/16inx 0.25mm PAEK stator 4 Port Replacement rotor	Cheminert Valves	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
C52-1340IX	Microbore 1/16inx 0.25mm PAEK stator 10 Port w/Motor Only	Cheminert Valves	\$1,148.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C52-1344I-S	Microbore 1/16inx 0.25mm PEEK stator 4 Port w/Motor Sensor Only	Cheminert Valves	\$1,083.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
C52-1C48	Microbore 1/16inx 0.25mm PEEK stator 8 Port Replace stator	Cheminert Valves	\$702.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
C52-1000IX	Microbore 1/16inx 0.25mm N60 SS stator 10 Port w/Motor Only	Cheminert Valves	\$1,020.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
CSL10	10ul Stainless Steel Sample loops for C1 C2 and C3 inc fits	Cheminert Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
CZSL5KPK	5ml PEEK Sample loops for C1 C2 and C3 inc fits	Cheminert Valves	\$135.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
CSL20	20ul 1/16in Stainless steel loops inc fits	Cheminert Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
CZSL20PK	20ul PEEK Sample loops for C1 C2 and C3 inc fits	Cheminert Valves	\$34.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
C52-1004I	Microbore 1/16inx 0.25mm N60 SS stator 10 Port w/Integrated actuator	Cheminert Valves	\$1,148.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
C52-1006I	Microbore 1/16inx 0.25mm N60 SS stator 6 Port w/Integrated actuator	Cheminert Valves	\$1,148.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
C52-1008I	Microbore 1/16inx 0.25mm N60 SS stator 8 Port w/Integrated actuator	Cheminert Valves	\$1,212.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
C52-1000I	Microbore 1/16inx 0.25mm N60 SS stator 10 Port w/Integrated actuator	Cheminert Valves	\$1,275.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
C52-1004I-S	Microbore 1/16inx 0.25mm N60 SS stator 10 Port w/Motor Sensor Only	Cheminert Valves	\$957.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
C52-1006I-S	Microbore 1/16inx 0.25mm N60 SS stator 6 Port w/Motor Sensor Only	Cheminert Valves	\$957.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
C52-1008I-S	Microbore 1/16inx 0.25mm N60 SS stator 8 Port w/Motor Sensor Only	Cheminert Valves	\$1,020.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
C52-1000I-S	Microbore 1/16inx 0.25mm N60 SS stator 10 Port w/Motor Sensor Only	Cheminert Valves	\$1,083.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
C52-1004IX	Microbore 1/16inx 0.25mm N60 SS stator 10 Port w/Motor Only	Cheminert Valves	\$892.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
C52-1348I-S	Microbore 1/16inx 0.25mm PEEK stator 8 Port w/Motor Sensor Only	Cheminert Valves	\$1,148.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
C52-1C06	Microbore 1/16inx 0.25mm N60 SS stator 6 Port, Replace stator	Cheminert Valves	\$490.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
C52-1340I	Microbore 1/16inx 0.25mm PEEK stator 10 Port w/Integrated actuator	Cheminert Valves	\$1,402.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
C52-1348I	Microbore 1/16inx 0.25mm PEEK stator 8 Port w/Integrated actuator	Cheminert Valves	\$1,338.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C52-1346I	Microbore 1/16inx 0.25mm PAEK stator 6 Port w/Integrated actuator	Cheminert Valves	\$1,275.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
C52-1344I	Microbore 1/16inx 0.25mm PAEK stator 4 Port w/Integrated actuator	Cheminert Valves	\$1,275.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
C52-1006IX	Microbore 1/16inx 0.25mm N60 SS stator 6 Port w/Motor Only	Cheminert Valves	\$892.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
C52-1C08	Microbore 1/16inx 0.25mm N60 SS stator 8 Port, Replace stator	Cheminert Valves	\$561.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
C52-1008IX	Microbore 1/16inx 0.25mm N60 SS stator 8 Port w/Motor Only	Cheminert Valves	\$957.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
C52-1C04	Microbore 1/16inx 0.25mm N60 SS stator 4 Port, Replace stator	Cheminert Valves	\$490.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
C2-10R0	Microbore 1/16inx 0.25mm N60 SS stator 10 Port Replacement rotor	Cheminert Valves	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
C2-10R8H	Microbore 1/16inx 0.25mm N60 SS stator 8 Port Replacement rotor	Cheminert Valves	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
C2-10R6	Microbore 1/16inx 0.25mm N60 SS stator 6 Port Replacement rotor	Cheminert Valves	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
C2-10R4	Microbore 1/16inx 0.25mm N60 SS stator 4 Port Replacement rotor	Cheminert Valves	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
C52-1346I-S	Microbore 1/16inx 0.25mm PAEK stator 6 Port w/Motor Sensor Only	Cheminert Valves	\$1,083.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
C52-1C00	Microbore 1/16inx 0.25mm N60 SS stator 10 Port, Replace stator	Cheminert Valves	\$632.00	<a href="http://www.chromtech.net.au/pdf2/vici09_178.pdf">http://www.chromtech.net.au/pdf2/vici09_178.pdf</a>
C45-9784A	4 position Pneumatic act. 1/2-20 fittings for 1/4in tubing 4.6mm ports	Cheminert Valves	\$1,185.00	<a href="http://www.chromtech.net.au/pdf2/vici09_177.pdf">http://www.chromtech.net.au/pdf2/vici09_177.pdf</a>
CFL-4PPS	Fittings for C45 Valve, PPS	Cheminert Valves	\$9.00	<a href="http://www.chromtech.net.au/pdf2/vici09_177.pdf">http://www.chromtech.net.au/pdf2/vici09_177.pdf</a>
CFL-4KF	Fittings for C45 Valve, CTFE	Cheminert Valves	\$10.00	<a href="http://www.chromtech.net.au/pdf2/vici09_177.pdf">http://www.chromtech.net.au/pdf2/vici09_177.pdf</a>
CFL-4D	Fittings for C45 Valve, Delrin	Cheminert Valves	\$8.00	<a href="http://www.chromtech.net.au/pdf2/vici09_177.pdf">http://www.chromtech.net.au/pdf2/vici09_177.pdf</a>
C45-97R6	6 position Replace rotor 1/2-20 fittings for 1/4in tubing 4.6mm ports	Cheminert Valves	\$141.00	<a href="http://www.chromtech.net.au/pdf2/vici09_177.pdf">http://www.chromtech.net.au/pdf2/vici09_177.pdf</a>
C45-9786D	6 position Replacement valve 1/2-20 fittings for 1/4in tubing 4.6mm ports	Cheminert Valves	\$740.00	<a href="http://www.chromtech.net.au/pdf2/vici09_177.pdf">http://www.chromtech.net.au/pdf2/vici09_177.pdf</a>
C45-9786EMT	6 position Microelectric act. 1/2-20 fittings for 1/4in tubing 4.6mm ports	Cheminert Valves	\$2,052.00	<a href="http://www.chromtech.net.au/pdf2/vici09_177.pdf">http://www.chromtech.net.au/pdf2/vici09_177.pdf</a>
C45-9784EMT	4 position Microelectric act. 1/2-20 fittings for 1/4in tubing 4.6mm ports	Cheminert Valves	\$1,983.00	<a href="http://www.chromtech.net.au/pdf2/vici09_177.pdf">http://www.chromtech.net.au/pdf2/vici09_177.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C45-9786E	6 position Electric act. 1/2-20 fittings for 1/4in tubing 4.6mm ports	Cheminert Valves	\$1,773.00	<a href="http://www.chromtech.net.au/pdf2/vici09_177.pdf">http://www.chromtech.net.au/pdf2/vici09_177.pdf</a>
C45-9786A	6 position Pneumatic act. 1/2-20 fittings for 1/4in tubing 4.6mm ports	Cheminert Valves	\$1,256.00	<a href="http://www.chromtech.net.au/pdf2/vici09_177.pdf">http://www.chromtech.net.au/pdf2/vici09_177.pdf</a>
C45-97R4	4 position Replace rotor 1/2-20 fittings for 1/4in tubing 4.6mm ports	Cheminert Valves	\$141.00	<a href="http://www.chromtech.net.au/pdf2/vici09_177.pdf">http://www.chromtech.net.au/pdf2/vici09_177.pdf</a>
CFL-CB4KF-S	CTFE Ferrules for C45 Valve	Cheminert Valves	\$6.00	<a href="http://www.chromtech.net.au/pdf2/vici09_177.pdf">http://www.chromtech.net.au/pdf2/vici09_177.pdf</a>
C45-9784E	4 position Electric act. 1/2-20 fittings for 1/4in tubing 4.6mm ports	Cheminert Valves	\$1,702.00	<a href="http://www.chromtech.net.au/pdf2/vici09_177.pdf">http://www.chromtech.net.au/pdf2/vici09_177.pdf</a>
C45-9784D	4 position Replacement valve 1/2-20 fittings for 1/4in tubing 4.6mm ports	Cheminert Valves	\$669.00	<a href="http://www.chromtech.net.au/pdf2/vici09_177.pdf">http://www.chromtech.net.au/pdf2/vici09_177.pdf</a>
C35Z-31R20	20 Position Replacement rotor 1/16in ZDV fittings 0.75mm ports	Cheminert Valves	\$192.00	<a href="http://www.chromtech.net.au/pdf2/vici09_176.pdf">http://www.chromtech.net.au/pdf2/vici09_176.pdf</a>
C35Z-31R20	20 Position Replacement rotor 1/16in ZDV fittings 0.75mm ports	Cheminert Valves	\$192.00	<a href="http://www.chromtech.net.au/pdf2/vici09_176.pdf">http://www.chromtech.net.au/pdf2/vici09_176.pdf</a>
C35Z-31826D	26 Position Replacement valve 1/16in ZDV fittings 0.75mm ports	Cheminert Valves	\$873.00	<a href="http://www.chromtech.net.au/pdf2/vici09_176.pdf">http://www.chromtech.net.au/pdf2/vici09_176.pdf</a>
C35Z-31824D	20 Position Replacement valve 1/16in ZDV fittings 0.75mm ports	Cheminert Valves	\$873.00	<a href="http://www.chromtech.net.au/pdf2/vici09_176.pdf">http://www.chromtech.net.au/pdf2/vici09_176.pdf</a>
C35Z-31820D	24 Position Replacement valve 1/16in ZDV fittings 0.75mm ports	Cheminert Valves	\$804.00	<a href="http://www.chromtech.net.au/pdf2/vici09_176.pdf">http://www.chromtech.net.au/pdf2/vici09_176.pdf</a>
C35Z-31824EMT	24 Position Microelectric actuator 1/16in ZDV fittings 0.75mm ports	Cheminert Valves	\$2,187.00	<a href="http://www.chromtech.net.au/pdf2/vici09_176.pdf">http://www.chromtech.net.au/pdf2/vici09_176.pdf</a>
C35Z-31820EMT	20 Position Microelectric actuator 1/16in ZDV fittings 0.75mm ports	Cheminert Valves	\$2,116.00	<a href="http://www.chromtech.net.au/pdf2/vici09_176.pdf">http://www.chromtech.net.au/pdf2/vici09_176.pdf</a>
C35Z-31R20	20 Position Replacement rotor 1/16in ZDV fittings 0.75mm ports	Cheminert Valves	\$192.00	<a href="http://www.chromtech.net.au/pdf2/vici09_176.pdf">http://www.chromtech.net.au/pdf2/vici09_176.pdf</a>
C35Z-31826EMT	26 Position Microelectric actuator 1/16in ZDV fittings 0.75mm ports	Cheminert Valves	\$2,187.00	<a href="http://www.chromtech.net.au/pdf2/vici09_176.pdf">http://www.chromtech.net.au/pdf2/vici09_176.pdf</a>
C25-6186EMH	6 position Microelectric act. 1/4-28 fittings for 1/8in tubing 1.50mm ports	Cheminert Valves	\$1,620.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25-6180E	10 position Electric act. 1/4-28 fittings for 1/8in tubing 1.50mm ports	Cheminert Valves	\$1,575.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25-618	8 position Replacement stator 1/4-28 fittings for 1/8in tubing 1.50mm ports	Cheminert Valves	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25-610	10 position Replacement stator 1/4-28 fittings for 1/8in tubing 1.50mm ports	Cheminert Valves	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>



## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C25-684	4 position Replacement stator 1/4-28 fittings for 1/8in tubing 1.50mm ports	Cheminert Valves	\$230.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25-686	6 position Replacement stator 1/4-28 fittings for 1/8in tubing 1.50mm ports	Cheminert Valves	\$230.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25-688	8 position Replacement stator 1/4-28 fittings for 1/8in tubing 1.50mm ports	Cheminert Valves	\$357.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25-680	10 position Replacement stator 1/4-28 fittings for 1/8in tubing 1.50mm ports	Cheminert Valves	\$388.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25-6188D	8 position Replacement valve 1/4-28 fittings for 1/8in tubing 1.50mm ports	Cheminert Valves	\$516.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25-314	4 position Replacement rotor 1/4-28 1/16inTubing 0.75mm ports	Cheminert Valves	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25-6186D	6 position Replacement valve 1/4-28 fittings for 1/8in tubing 1.50mm ports	Cheminert Valves	\$382.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25-6184D	4 position Replacement valve 1/4-28 fittings for 1/8in tubing 1.50mm ports	Cheminert Valves	\$382.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25-614	4 position Replacement stator 1/4-28 fittings for 1/8in tubing 1.50mm ports	Cheminert Valves	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25-6188EMH	8 position Microelectric act. 1/4-28 fittings for 1/8in tubing 1.50mm ports	Cheminert Valves	\$1,754.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25-388	8 position Replacement stator 1/4-28 fittings for 1/16in tubing 0.75mm ports	Cheminert Valves	\$357.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25-6184EMH	4 position Microelectric act. 1/4-28 fittings for 1/8in tubing 1.50mm ports	Cheminert Valves	\$1,620.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25-6184	4 position Manual 1/4-28 fittings for 1/8in tubing 1.50mm ports	Cheminert Valves	\$382.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25-6186	6 position Manual 1/4-28 fittings for 1/8in tubing 1.50mm ports	Cheminert Valves	\$382.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25-6188	8 position Manual 1/4-28 fittings for 1/8in tubing 1.50mm ports	Cheminert Valves	\$516.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25-6180	10 position Manual 1/4-28 fittings for 1/8in tubing 1.50mm ports	Cheminert Valves	\$542.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25-6184A	4 position Pneumatic act. 1/4-28 fittings for 1/8in tubing 1.50mm ports	Cheminert Valves	\$898.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25-6186A	6 position Pneumatic act. 1/4-28 fittings for 1/8in tubing 1.50mm ports	Cheminert Valves	\$898.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C25-6188A	8 position Pneumatic act. 1/4-28 fittings for 1/8in tubing 1.50mm ports	Cheminert Valves	\$1,032.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25-6180A	10 position Pneumatic act. 1/4-28 fittings for 1/8in tubing 1.50mm ports	Cheminert Valves	\$1,059.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25-6184E	4 position Electric act. 1/4-28 fittings for 1/8in tubing 1.50mm ports	Cheminert Valves	\$1,416.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25-6186E	6 position Electric act. 1/4-28 fittings for 1/8in tubing 1.50mm ports	Cheminert Valves	\$1,416.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25-6188E	8 position Electric act. 1/4-28 fittings for 1/8in tubing 1.50mm ports	Cheminert Valves	\$1,550.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25-6180EMH	10 position Microelectric act. 1/4-28 fittings for 1/8in tubing 1.50mm ports	Cheminert Valves	\$1,779.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25-3186A	6 position Pneumatic act. 1/4-28 fittings for 1/16in tubing 0.75mm ports	Cheminert Valves	\$898.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25-3180D	10 position Replacement valve 1/4-28 fittings for 1/16in tubing 0.75mm ports	Cheminert Valves	\$542.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25-3188D	8 position Replacement valve 1/4-28 fittings for 1/16in tubing 0.75mm ports	Cheminert Valves	\$516.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25-3186D	6 position Replacement valve 1/4-28 fittings for 1/16in tubing 0.75mm ports	Cheminert Valves	\$382.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25-3184D	4 position Replacement valve 1/4-28 fittings for 1/16in tubing 0.75mm ports	Cheminert Valves	\$382.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25-3180EMH	10 position Microelectric act. 1/4-28 fittings for 1/16in tubing 0.75mm ports	Cheminert Valves	\$1,779.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25-3188EMH	8 position Microelectric act. 1/4-28 fittings for 1/16in tubing 0.75mm ports	Cheminert Valves	\$1,754.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25-3186EMH	6 position Microelectric act. 1/4-28 fittings for 1/16in tubing 0.75mm ports	Cheminert Valves	\$1,620.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25-3184EMH	4 position Microelectric act. 1/4-28 fittings for 1/16in tubing 0.75mm ports	Cheminert Valves	\$1,620.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25-3180E	10 position Electric act. 1/4-28 fittings for 1/16in tubing 0.75mm ports	Cheminert Valves	\$1,575.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25-3188E	8 position Electric act. 1/4-28 fittings for 1/16in tubing 0.75mm ports	Cheminert Valves	\$1,550.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25-3186E	6 position Electric act. 1/4-28 fittings for 1/16in tubing 0.75mm ports	Cheminert Valves	\$1,416.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C25-3184E	4 position Electric act. 1/4-28 fittings for 1/16in tubing 0.75mm ports	Cheminert Valves	\$1,416.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25-616	6 position Replacement stator 1/4-28 fittings for 1/8in tubing 1.50mm ports	Cheminert Valves	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25-3188A	8 position Pneumatic act. 1/4-28 fittings for 1/16in tubing 0.75mm ports	Cheminert Valves	\$1,032.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25-380	10 position Replacement stator 1/4-28 fittings for 1/16in tubing 0.75mm ports	Cheminert Valves	\$388.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25-3184A	4 position Pneumatic act. 1/4-28 fittings for 1/16in tubing 0.75mm ports	Cheminert Valves	\$898.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25-3180	10 position Manual 1/4-28 fittings for 1/16in tubing 0.75mm ports	Cheminert Valves	\$542.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25-3188	8 position Manual 1/4-28 fittings for 1/16in tubing 0.75mm ports	Cheminert Valves	\$516.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25-3186	6 position Manual 1/4-28 fittings for 1/16in tubing 0.75mm ports	Cheminert Valves	\$382.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25-3184	6 position Manual 1/4-28 fittings for 1/16in tubing 0.75mm ports	Cheminert Valves	\$382.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25-316	6 position Replacement rotor 1/4-28 1/16inTubing 0.75mm ports	Cheminert Valves	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25-318	8 position Replacement rotor 1/4-28 1/16inTubing 0.75mm ports	Cheminert Valves	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25-310	10 position Replacement rotor 1/4-28 1/16inTubing 0.75mm ports	Cheminert Valves	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25-384	4 position Replacement stator 1/4-28 fittings for 1/16in tubing 0.75mm ports	Cheminert Valves	\$230.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25-6180D	10 position Replacement valve 1/4-28 fittings for 1/8in tubing 1.50mm ports	Cheminert Valves	\$542.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25-386	6 position Replacement stator 1/4-28 fittings for 1/16in tubing 0.75mm ports	Cheminert Valves	\$230.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25-3180A	10 position Pneumatic act.1/4-28 fittings for 1/16in tubing 0.75mm ports	Cheminert Valves	\$1,059.00	<a href="http://www.chromtech.net.au/pdf2/vici09_175.pdf">http://www.chromtech.net.au/pdf2/vici09_175.pdf</a>
C25Z-3186D	6 position Replacement valve 1/16in ZDV fittings 0.75mm ports	Cheminert Valves	\$382.00	<a href="http://www.chromtech.net.au/pdf2/vici09_174.pdf">http://www.chromtech.net.au/pdf2/vici09_174.pdf</a>
C25Z-380	10 position Replacement stator 1/16in ZDV fittings 0.75mm ports	Cheminert Valves	\$306.00	<a href="http://www.chromtech.net.au/pdf2/vici09_174.pdf">http://www.chromtech.net.au/pdf2/vici09_174.pdf</a>

VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C25Z-388	8 position Replacement stator 1/16in ZDV fittings 0.75mm ports	Cheminert Valves	\$280.00	<a href="http://www.chromtech.net.au/pdf2/vici09_174.pdf">http://www.chromtech.net.au/pdf2/vici09_174.pdf</a>
C25Z-3186	6 position Manual 1/16in ZDV fittings 0.75mm ports	Cheminert Valves	\$382.00	<a href="http://www.chromtech.net.au/pdf2/vici09_174.pdf">http://www.chromtech.net.au/pdf2/vici09_174.pdf</a>
C25Z-3180	10 position Manual 1/16in ZDV fittings 0.75mm ports	Cheminert Valves	\$465.00	<a href="http://www.chromtech.net.au/pdf2/vici09_174.pdf">http://www.chromtech.net.au/pdf2/vici09_174.pdf</a>
C25Z-31814	14 position Manual 1/16in ZDV fittings 0.75mm ports	Cheminert Valves	\$606.00	<a href="http://www.chromtech.net.au/pdf2/vici09_174.pdf">http://www.chromtech.net.au/pdf2/vici09_174.pdf</a>
C25Z-325	14 position Replace rotor 1/16in ZDV fittings 0.75mm ports	Cheminert Valves	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_174.pdf">http://www.chromtech.net.au/pdf2/vici09_174.pdf</a>
C15-310	10 position Replacement rotor 1/4-28 fittings for 1/16in tubing 0.75mm ports	Cheminert Valves	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_174.pdf">http://www.chromtech.net.au/pdf2/vici09_174.pdf</a>
C15-310	10 position Replacement rotor 1/4-28 fittings for 1/16in tubing 0.75mm ports	Cheminert Valves	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_174.pdf">http://www.chromtech.net.au/pdf2/vici09_174.pdf</a>
C15-310	10 position Replacement rotor 1/4-28 fittings for 1/16in tubing 0.75mm ports	Cheminert Valves	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_174.pdf">http://www.chromtech.net.au/pdf2/vici09_174.pdf</a>
C25Z-31814D	14 position Replacement valve 1/16in ZDV fittings 0.75mm ports	Cheminert Valves	\$606.00	<a href="http://www.chromtech.net.au/pdf2/vici09_174.pdf">http://www.chromtech.net.au/pdf2/vici09_174.pdf</a>
C25Z-38-14	14 position Replacement stator 1/16in ZDV fittings 0.75mm ports	Cheminert Valves	\$447.00	<a href="http://www.chromtech.net.au/pdf2/vici09_174.pdf">http://www.chromtech.net.au/pdf2/vici09_174.pdf</a>
C25Z-3188D	8 position Replacement valve 1/16in ZDV fittings 0.75mm ports	Cheminert Valves	\$434.00	<a href="http://www.chromtech.net.au/pdf2/vici09_174.pdf">http://www.chromtech.net.au/pdf2/vici09_174.pdf</a>
C25Z-386	6 position Replacement stator 1/16in ZDV fittings 0.75mm ports	Cheminert Valves	\$230.00	<a href="http://www.chromtech.net.au/pdf2/vici09_174.pdf">http://www.chromtech.net.au/pdf2/vici09_174.pdf</a>
C25Z-31814EMH	14 position Microelectric act. 1/16in ZDV fittings 0.75mm ports	Cheminert Valves	\$1,970.00	<a href="http://www.chromtech.net.au/pdf2/vici09_174.pdf">http://www.chromtech.net.au/pdf2/vici09_174.pdf</a>
C25Z-3180EMH	10 position Microelectric act. 1/16in ZDV fittings 0.75mm ports	Cheminert Valves	\$1,702.00	<a href="http://www.chromtech.net.au/pdf2/vici09_174.pdf">http://www.chromtech.net.au/pdf2/vici09_174.pdf</a>
C25Z-3188EMH	8 position Microelectric act. 1/16in ZDV fittings 0.75mm ports	Cheminert Valves	\$1,671.00	<a href="http://www.chromtech.net.au/pdf2/vici09_174.pdf">http://www.chromtech.net.au/pdf2/vici09_174.pdf</a>
C25Z-3186EMH	6 position Microelectric act. 1/16in ZDV fittings 0.75mm ports	Cheminert Valves	\$1,620.00	<a href="http://www.chromtech.net.au/pdf2/vici09_174.pdf">http://www.chromtech.net.au/pdf2/vici09_174.pdf</a>
C25Z-31814E	14 position Electric act. 1/16in ZDV fittings 0.75mm ports	Cheminert Valves	\$1,638.00	<a href="http://www.chromtech.net.au/pdf2/vici09_174.pdf">http://www.chromtech.net.au/pdf2/vici09_174.pdf</a>
C25Z-3180E	10 position Electric act. 1/16in ZDV fittings 0.75mm ports	Cheminert Valves	\$1,498.00	<a href="http://www.chromtech.net.au/pdf2/vici09_174.pdf">http://www.chromtech.net.au/pdf2/vici09_174.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C25Z-3188E	8 position Electric act.1/16in ZDV fittings 0.75mm ports	Cheminert Valves	\$1,467.00	<a href="http://www.chromtech.net.au/pdf2/vici09_174.pdf">http://www.chromtech.net.au/pdf2/vici09_174.pdf</a>
C25Z-3186E	6 position Electric act. 1/16in ZDV fittings 0.75mm ports	Cheminert Valves	\$1,416.00	<a href="http://www.chromtech.net.au/pdf2/vici09_174.pdf">http://www.chromtech.net.au/pdf2/vici09_174.pdf</a>
C25Z-31814A	14 position Pneumatic act.1/16in ZDV fittings 0.75mm ports	Cheminert Valves	\$1,122.00	<a href="http://www.chromtech.net.au/pdf2/vici09_174.pdf">http://www.chromtech.net.au/pdf2/vici09_174.pdf</a>
C25Z-3180A	10 position Pneumatic act.1/16in ZDV fittings 0.75mm ports	Cheminert Valves	\$981.00	<a href="http://www.chromtech.net.au/pdf2/vici09_174.pdf">http://www.chromtech.net.au/pdf2/vici09_174.pdf</a>
C25Z-3188A	8 position Pneumatic act. 1/16in ZDV fittings 0.75mm ports	Cheminert Valves	\$950.00	<a href="http://www.chromtech.net.au/pdf2/vici09_174.pdf">http://www.chromtech.net.au/pdf2/vici09_174.pdf</a>
C25Z-3186A	6 position Pneumatic act. 1/16in ZDV fittings 0.75mm ports	Cheminert Valves	\$898.00	<a href="http://www.chromtech.net.au/pdf2/vici09_174.pdf">http://www.chromtech.net.au/pdf2/vici09_174.pdf</a>
C25Z-3180D	10 position Replacement valve 1/16in ZDV fittings 0.75mm ports	Cheminert Valves	\$465.00	<a href="http://www.chromtech.net.au/pdf2/vici09_174.pdf">http://www.chromtech.net.au/pdf2/vici09_174.pdf</a>
C25Z-3188	8 position Manual 1/16in ZDV fittings 0.75mm ports	Cheminert Valves	\$434.00	<a href="http://www.chromtech.net.au/pdf2/vici09_174.pdf">http://www.chromtech.net.au/pdf2/vici09_174.pdf</a>
C5-2346EMTD	PAEK Stator w/Microelectric actuator 6 Position 1/16in ZDV fittings 0.40mm ports	Cheminert Valves	\$3,321.00	<a href="http://www.chromtech.net.au/pdf2/vici09_173.pdf">http://www.chromtech.net.au/pdf2/vici09_173.pdf</a>
C5H-2008D	N60 stainless stator Replacement Valve 8 Column 1/16in fittings 0.40mm ports	Cheminert Valves	\$669.00	<a href="http://www.chromtech.net.au/pdf2/vici09_173.pdf">http://www.chromtech.net.au/pdf2/vici09_173.pdf</a>
C5H-2000D	N60 stainless stator Replacement Valve 6 Column 1/16in fittings 0.40mm ports	Cheminert Valves	\$740.00	<a href="http://www.chromtech.net.au/pdf2/vici09_173.pdf">http://www.chromtech.net.au/pdf2/vici09_173.pdf</a>
C5H-2008EMTD	N60 stainless stator w/Microelectric actuator 8 Column 1/16in fittings 0.40mm ports	Cheminert Valves	\$3,258.00	<a href="http://www.chromtech.net.au/pdf2/vici09_173.pdf">http://www.chromtech.net.au/pdf2/vici09_173.pdf</a>
C5-2C06	N60 stainless stator Replacement stator 10 Position 1/16in ZDV fittings 0.40mm ports	Cheminert Valves	\$338.00	<a href="http://www.chromtech.net.au/pdf2/vici09_173.pdf">http://www.chromtech.net.au/pdf2/vici09_173.pdf</a>
C5-2C00H	N60 stainless stator Replacement stator 6 Position 1/16in ZDV fittings 0.40mm ports	Cheminert Valves	\$478.00	<a href="http://www.chromtech.net.au/pdf2/vici09_173.pdf">http://www.chromtech.net.au/pdf2/vici09_173.pdf</a>
C5H-2348EMTD	PAEK Stator w/Microelectric actuator 8Column 1/16in fittings 0.40mm ports	Cheminert Valves	\$3,538.00	<a href="http://www.chromtech.net.au/pdf2/vici09_173.pdf">http://www.chromtech.net.au/pdf2/vici09_173.pdf</a>
C5H-2340EMTD	PAEK Stator w/Microelectric actuator 10 Column 1/16in fittings 0.40mm ports	Cheminert Valves	\$3,678.00	<a href="http://www.chromtech.net.au/pdf2/vici09_173.pdf">http://www.chromtech.net.au/pdf2/vici09_173.pdf</a>
C5-2346D	PAEK Stator Replacement valve 6 Position 1/16in ZDV fittings 0.40mm ports	Cheminert Valves	\$740.00	<a href="http://www.chromtech.net.au/pdf2/vici09_173.pdf">http://www.chromtech.net.au/pdf2/vici09_173.pdf</a>
C5H-2348D	PAEK Stator Replacement Valve 8 Column 1/16in fittings 0.40mm ports	Cheminert Valves	\$810.00	<a href="http://www.chromtech.net.au/pdf2/vici09_173.pdf">http://www.chromtech.net.au/pdf2/vici09_173.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C5H-2340D	PAEK Stator Replacement Valve 10 Column 1/16in fittings 0.40mm ports	Cheminert Valves	\$879.00	<a href="http://www.chromtech.net.au/pdf2/vici09_173.pdf">http://www.chromtech.net.au/pdf2/vici09_173.pdf</a>
C5-23R6	PAEK Stator Replacement rotor 6 Port 1/16in fittings 0.75mm ports	Cheminert Valves	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_173.pdf">http://www.chromtech.net.au/pdf2/vici09_173.pdf</a>
C5-23R8H	PAEK Stator Replacement rotor 8 Position 1/16in ZDV fittings 0.40mm ports	Cheminert Valves	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_173.pdf">http://www.chromtech.net.au/pdf2/vici09_173.pdf</a>
C5-23R0H	PAEK Stator Replacement rotor 10 Position 1/16in ZDV fittings 0.40mm ports	Cheminert Valves	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_173.pdf">http://www.chromtech.net.au/pdf2/vici09_173.pdf</a>
C5-2C46	PAEK Stator Replacement stator 6 Port 1/16in fittings 0.40mm ports	Cheminert Valves	\$478.00	<a href="http://www.chromtech.net.au/pdf2/vici09_173.pdf">http://www.chromtech.net.au/pdf2/vici09_173.pdf</a>
C5-2C48H	PAEK Stator Replacement stator 6 Port 1/16in fittings 0.40mm ports	Cheminert Valves	\$549.00	<a href="http://www.chromtech.net.au/pdf2/vici09_173.pdf">http://www.chromtech.net.au/pdf2/vici09_173.pdf</a>
C5-2C40H	PAEK Stator Replacement stator 10 Port 1/16in fittings 0.40mm ports	Cheminert Valves	\$618.00	<a href="http://www.chromtech.net.au/pdf2/vici09_173.pdf">http://www.chromtech.net.au/pdf2/vici09_173.pdf</a>
I-22697	2-position valves Microelectric actuator RS-232 interface cable	Cheminert Valves	\$38.00	<a href="http://www.chromtech.net.au/pdf2/vici09_173.pdf">http://www.chromtech.net.au/pdf2/vici09_173.pdf</a>
C5-2C08H	N60 stainless stator Replacement stator 8 Position 1/16in ZDV fittings 0.40mm ports	Cheminert Valves	\$408.00	<a href="http://www.chromtech.net.au/pdf2/vici09_173.pdf">http://www.chromtech.net.au/pdf2/vici09_173.pdf</a>
C5H-2000EMTD	N60 stainless stator w/Microelectric actuator 10 Column 1/16in fittings 0.40mm ports	Cheminert Valves	\$3,398.00	<a href="http://www.chromtech.net.au/pdf2/vici09_173.pdf">http://www.chromtech.net.au/pdf2/vici09_173.pdf</a>
C5H-2030E	Titanium stator w/Std Electric actuator 10 Position 1/16in fittings 0.4mm ports	Cheminert Valves	\$2,142.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-2036A	Titanium stator Pneumatic actuator 6 Port 1/16in fittings 0.75mm ports	Cheminert Valves	\$1,485.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5H-2030D	Titanium stator Replace Valve 10 Position 1/16in fittings 0.4mm ports	Cheminert Valves	\$1,110.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5H-2038A	Titanium stator w/Pneumatic actuator 8 Position 1/16in fittings 0.4mm ports	Cheminert Valves	\$1,556.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5H-2038D	Titanium stator Replace Valve 8 Position 1/16in fittings 0.4mm ports	Cheminert Valves	\$1,040.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5H-2030A	Titanium stator w/Pneumatic actuator 10 Position 1/16in fittings 0.4mm ports	Cheminert Valves	\$1,626.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-2034E	Titanium stator Electric actuator 4 Position 1/16in ZDV fittings 0.40mm ports	Cheminert Valves	\$2,001.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-2036EMH	Titanium stator Microelectric actuator 6 Port 1/16in fittings 0.75mm ports	Cheminert Valves	\$2,205.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C5H-2038E	Titanium stator w/Std Electric actuator 8 Position 1/16in fittings 0.4mm ports	Cheminert Valves	\$2,072.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-2034EMH	Titanium stator Microelectric actuator 4 Position 1/16in ZDV fittings 0.40mm ports	Cheminert Valves	\$2,205.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-2034A	Titanium stator Pneumatic actuator 4 Position 1/16in ZDV fittings 0.40mm ports	Cheminert Valves	\$1,485.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-2036D	Titanium stator Replacement Valve Stream Selector 6 Port 1/16in fittings 0.75mm ports	Cheminert Valves	\$969.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5H-2000A	N60 SS Stator w/Pneumatic actuator 10 Position 1/16in fittings 0.4mm ports	Cheminert Valves	\$1,256.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-2034D	Titanium stator Replacement Valve Stream Selector 4 Position 1/16in ZDV fittings 0.40mm ports	Cheminert Valves	\$969.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5H-2030EMT	Titanium stator w/Microelectric actuator 10 Position 1/16in fittings 0.4mm ports	Cheminert Valves	\$2,422.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5H-2038EMT	Titanium stator w/Microelectric actuator 8 Position 1/16in fittings 0.4mm ports	Cheminert Valves	\$2,352.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-2036E	Titanium stator Electric actuator 6 Port 1/16in fittings 0.75mm ports	Cheminert Valves	\$2,001.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5H-2008	N60 SS Stator Manual 8 Position 1/16in fittings 0.4mm ports	Cheminert Valves	\$669.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-2C46	PAEK Stator Replacement stator 6 Port 1/16in fittings 0.40mm ports	Cheminert Valves	\$478.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-2C48H	PAEK Stator Replacement stator 6 Port 1/16in fittings 0.40mm ports	Cheminert Valves	\$549.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-2C40H	PAEK stator Replacement stator 10 Position 1/16in ZDV fittings 0.40mm ports	Cheminert Valves	\$618.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-2034	Titanium stator Manual 4 Position 1/16in ZDV fittings 0.40mm ports	Cheminert Valves	\$969.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-20R0H	Replace Titanium rotor 6 Port 1/16in fittings 0.4mm ports	Cheminert Valves	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-2036	Titanium stator Manual 4 Position 1/16in ZDV fittings 0.40mm ports	Cheminert Valves	\$969.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-2006E	N60 stainless stator Electric actuator 6 Port 1/16in fittings 0.75mm ports	Cheminert Valves	\$1,632.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-2006	N60 stainless stator Manual 4 Position 1/16in ZDV fittings 0.40mm ports	Cheminert Valves	\$600.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C5H-2030	Titanium stator Manual 6 Position 1/16in fittings 0.4mm ports	Cheminert Valves	\$1,110.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5H-2000	N60 SS Stator Manual 6 Position 1/16in fittings 0.4mm ports	Cheminert Valves	\$740.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-2004A	N60 stainless stator Pneumatic actuator 4 Position 1/16in ZDV fittings 0.40mm ports	Cheminert Valves	\$1,116.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-2006A	N60 stainless stator Pneumatic actuator 6 Port 1/16in fittings 0.75mm ports	Cheminert Valves	\$1,116.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5H-2008A	N60 SS Stator w/Pneumatic actuator 8 Position 1/16in fittings 0.4mm ports	Cheminert Valves	\$1,185.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-2004E	N60 stainless stator Electric actuator 4 Position 1/16in ZDV fittings 0.40mm ports	Cheminert Valves	\$1,632.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5H-2008E	N60 SS Stator w/Std Electric actuator 8 Position 1/16in fittings 0.4mm ports	Cheminert Valves	\$1,702.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5H-2000E	N60 SS Stator w/Std Electric actuator 10 Position 1/16in fittings 0.4mm ports	Cheminert Valves	\$1,773.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-2004	N60 stainless stator Manual 4 Position 1/16in ZDV fittings 0.40mm ports	Cheminert Valves	\$600.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-2C06	N60 stainless stator Replacement stator 6 Position 1/16in ZDV fittings 0.40mm ports	Cheminert Valves	\$338.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-23R6	PAEK Stator Replacement rotor 6 Port 1/16in fittings 0.75mm ports	Cheminert Valves	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-23R8H	Replace PAEK rotor 8 Port 1/16in fittings 0.4mm ports	Cheminert Valves	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-23R0H	Replace PAEK rotor 10 Port 1/16in fittings 0.4mm ports	Cheminert Valves	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-2346	PAEK Stator Manual 4 Position 1/16in ZDV fittings 0.40mm ports	Cheminert Valves	\$740.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-2344	PAEK Stator Manual 4 Position 1/16in ZDV fittings 0.40mm ports	Cheminert Valves	\$740.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-2006EMT	N60 stainless stator w/Microelectric actuator Stream Selector 6 Port 1/16in fittings 0.40mm ports	Cheminert Valves	\$3,040.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-20R4	Titanium stator Replacement rotor 4 Position 1/16in ZDV fittings 0.40mm ports	Cheminert Valves	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-2C08H	N60 stainless stator Replacement stator 8 Position 1/16in ZDV fittings 0.40mm ports	Cheminert Valves	\$408.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>



## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C5H-2348D	PAEK Stator Replace Valve 8 Position 1/16in fittings 0.4mm ports	Cheminert Valves	\$810.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-2C04	N60 stainless stator Replacement stator 4 Position 1/16in ZDV fittings 0.40mm ports	Cheminert Valves	\$338.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-20R8H	Replace Titanium rotor 8 Port 1/16in fittings 0.4mm ports	Cheminert Valves	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-20R6	6 Column N60 stainless stator Replacement rotor 1/16in ZDV fittings 0.40mm ports	Cheminert Valves	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-20R4	Titanium stator Replacement rotor 4 Position 1/16in ZDV fittings 0.40mm ports	Cheminert Valves	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5H-2000D	N60 SS Stator Replace Valve 10 Position 1/16in fittings 0.4mm ports	Cheminert Valves	\$740.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5H-2008D	N60 SS Stator Replace Valve 8 Position 1/16in fittings 0.4mm ports	Cheminert Valves	\$669.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-2C00H	N60 stainless stator Replacement stator 10 Position 1/16in ZDV fittings 0.40mm ports	Cheminert Valves	\$478.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5H-2340E	Titanium stator w/Std Electric actuator 8 Position 1/16in fittings 0.4mm ports	Cheminert Valves	\$1,912.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5H-2340	Titanium stator Manual 8 Position 1/16in fittings 0.4mm ports	Cheminert Valves	\$879.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-2344A	PAEK Stator Pneumatic actuator 4 Position 1/16in ZDV fittings 0.40mm ports	Cheminert Valves	\$1,256.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-2346A	PAEK Stator Pneumatic actuator 6 Port 1/16in fittings 0.75mm ports	Cheminert Valves	\$1,256.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5H-2348A	PAEK Stator w/Pneumatic actuator 8 Position 1/16in fittings 0.4mm ports	Cheminert Valves	\$1,326.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5H-2340A	Titanium stator w/Pneumatic actuator 8 Position 1/16in fittings 0.4mm ports	Cheminert Valves	\$1,396.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-2344E	PAEK Stator Electric actuator 4 Position 1/16in ZDV fittings 0.40mm ports	Cheminert Valves	\$1,773.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-23R4	PAEK Stator Replacement rotor 4 Position 1/16in ZDV fittings 0.40mm ports	Cheminert Valves	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5H-2348E	PAEK Stator w/Std Electric actuator 8 Position 1/16in fittings 0.4mm ports	Cheminert Valves	\$1,970.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5H-2340D	Titanium stator Replace Valve 8 Position 1/16in fittings 0.4mm ports	Cheminert Valves	\$879.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C5-2344EMH	PAEK Stator Microelectric actuator 4 Position 1/16in ZDV fittings 0.40mm ports	Cheminert Valves	\$1,977.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-2346EMH	PAEK Stator Microelectric actuator 6 Port 1/16in fittings 0.75mm ports	Cheminert Valves	\$1,977.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5H-2348EMT	PAEK Stator w/Microelectric actuator 8 Position 1/16in fittings 0.4mm ports	Cheminert Valves	\$2,122.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5H-2340EMT	Titanium stator w/Microelectric actuator 8 Position 1/16in fittings 0.4mm ports	Cheminert Valves	\$2,193.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-2344D	PAEK Stator Replacement Valve Stream Selector 4 Position 1/16in ZDV fittings 0.40mm ports	Cheminert Valves	\$740.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-2346D	6 Column PAEK stator Replacement valve 1/16in ZDV fittings 0.40mm ports	Cheminert Valves	\$740.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5H-2000EMT	N60 SS Stator w/Microelectric actuator 10 Position 1/16in fittings 0.4mm ports	Cheminert Valves	\$2,052.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-2346E	PAEK Stator w/Std Electric actuator 6 Port 1/16in fittings 0.75mm ports	Cheminert Valves	\$1,773.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-2C30H	Titanium stator Replacement stator 10 Position 1/16in ZDV fittings 0.40mm ports	Cheminert Valves	\$822.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-2006D	N60 stainless stator Replacement Valve Stream Selector 6 Port 1/16in fittings 0.75mm ports	Cheminert Valves	\$600.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-20R8H	Titanium stator Replacement rotor 8 Position 1/16in ZDV fittings 0.40mm ports	Cheminert Valves	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-2004D	N60 stainless stator Replacement Valve Stream Selector 4 Position 1/16in ZDV fittings 0.40mm ports	Cheminert Valves	\$600.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-20R6	Titanium stator Replacement rotor 6 Position 1/16in ZDV fittings 0.40mm ports	Cheminert Valves	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-2006D	N60 stainless stator Replacement Valve Stream Selector 6 Port 1/16in fittings 0.40mm ports	Cheminert Valves	\$600.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-2C44	PAEK Stator Replacement stator 4 Position 1/16in ZDV fittings 0.40mm ports	Cheminert Valves	\$478.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-20R0H	Titanium stator Replacement rotor 10 Position 1/16in ZDV fittings 0.40mm ports	Cheminert Valves	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5H-2038	Titanium stator Manual 8 Position 1/16in fittings 0.4mm ports	Cheminert Valves	\$1,040.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-2C38H	Titanium stator Replacement stator 8 Position 1/16in ZDV fittings 0.40mm ports	Cheminert Valves	\$777.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C5-2C36	Titanium stator Replacement stator 6 Port 1/16in fittings 0.40mm ports	Cheminert Valves	\$708.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-2C34	Titanium stator Replacement stator 4 Position 1/16in ZDV fittings 0.40mm ports	Cheminert Valves	\$708.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5H-2348	PAEK Stator Manual 8 Position 1/16in fittings 0.4mm ports	Cheminert Valves	\$810.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-2004EMH	N60 stainless stator Microelectric actuator 4 Position 1/16in ZDV fittings 0.40mm ports	Cheminert Valves	\$1,964.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-2006EMH	N60 stainless stator Microelectric actuator 6 Port 1/16in fittings 0.75mm ports	Cheminert Valves	\$1,964.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5H-2008EMT	N60 SS Stator w/Microelectric actuator 8 Position 1/16in fittings 0.4mm ports	Cheminert Valves	\$1,983.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C5-20R6	6 Column N60 stainless stator Replacement rotor 1/16in ZDV fittings 0.40mm ports	Cheminert Valves	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_172.pdf">http://www.chromtech.net.au/pdf2/vici09_172.pdf</a>
C75H-1690EMT	1/16inx250um 10 Position,w/Microelectric actuator	Cheminert Valves	\$2,544.00	<a href="http://www.chromtech.net.au/pdf2/vici09_171.pdf">http://www.chromtech.net.au/pdf2/vici09_171.pdf</a>
C75H-1698EMT	1/16inx250um 8 Position,w/Microelectric actuator	Cheminert Valves	\$2,480.00	<a href="http://www.chromtech.net.au/pdf2/vici09_171.pdf">http://www.chromtech.net.au/pdf2/vici09_171.pdf</a>
C75H-1698D	1/16inx250um 8 Position,Replace Valve	Cheminert Valves	\$1,167.00	<a href="http://www.chromtech.net.au/pdf2/vici09_171.pdf">http://www.chromtech.net.au/pdf2/vici09_171.pdf</a>
C75H-1696D	1/16inx250um 6 Position,Replace Valve	Cheminert Valves	\$1,102.00	<a href="http://www.chromtech.net.au/pdf2/vici09_171.pdf">http://www.chromtech.net.au/pdf2/vici09_171.pdf</a>
C75H-1690D	1/16inx250um 10 Position,Replace Valve	Cheminert Valves	\$1,230.00	<a href="http://www.chromtech.net.au/pdf2/vici09_171.pdf">http://www.chromtech.net.au/pdf2/vici09_171.pdf</a>
C75-16R8	1/16inx250um 8 Position,Replace Rotor	Cheminert Valves	\$122.00	<a href="http://www.chromtech.net.au/pdf2/vici09_171.pdf">http://www.chromtech.net.au/pdf2/vici09_171.pdf</a>
C75-16R0	1/16inx250um 10 Position,Replace Rotor	Cheminert Valves	\$122.00	<a href="http://www.chromtech.net.au/pdf2/vici09_171.pdf">http://www.chromtech.net.au/pdf2/vici09_171.pdf</a>
C75-1C96	1/16inx250um 6 Position,Replace Stator	Cheminert Valves	\$912.00	<a href="http://www.chromtech.net.au/pdf2/vici09_171.pdf">http://www.chromtech.net.au/pdf2/vici09_171.pdf</a>
C75-1C98	1/16inx250um 8 Position,Replace Stator	Cheminert Valves	\$975.00	<a href="http://www.chromtech.net.au/pdf2/vici09_171.pdf">http://www.chromtech.net.au/pdf2/vici09_171.pdf</a>
C75-1C90	1/16inx250um 10 Position,Replace Stator	Cheminert Valves	\$1,040.00	<a href="http://www.chromtech.net.au/pdf2/vici09_171.pdf">http://www.chromtech.net.au/pdf2/vici09_171.pdf</a>
C75H-1690E	1/16inx250um 10 Position,w/Std Electric actuator	Cheminert Valves	\$2,264.00	<a href="http://www.chromtech.net.au/pdf2/vici09_171.pdf">http://www.chromtech.net.au/pdf2/vici09_171.pdf</a>
C75H-1696E	1/16inx250um 6 Position,w/Std Electric actuator	Cheminert Valves	\$2,136.00	<a href="http://www.chromtech.net.au/pdf2/vici09_171.pdf">http://www.chromtech.net.au/pdf2/vici09_171.pdf</a>
C75-16R6	1/16inx250um 6 Position,Replace Rotor	Cheminert Valves	\$122.00	<a href="http://www.chromtech.net.au/pdf2/vici09_171.pdf">http://www.chromtech.net.au/pdf2/vici09_171.pdf</a>
C75H-1696EMH	1/16inx250um 6 Position,w/Microelectric actuator	Cheminert Valves	\$2,340.00	<a href="http://www.chromtech.net.au/pdf2/vici09_171.pdf">http://www.chromtech.net.au/pdf2/vici09_171.pdf</a>
C75H-1698E	1/16inx250um 8 Position,w/Std Electric actuator	Cheminert Valves	\$2,199.00	<a href="http://www.chromtech.net.au/pdf2/vici09_171.pdf">http://www.chromtech.net.au/pdf2/vici09_171.pdf</a>
C75H-1690A	1/16inx250um 10 Position,w/Pneumatic actuator	Cheminert Valves	\$1,746.00	<a href="http://www.chromtech.net.au/pdf2/vici09_171.pdf">http://www.chromtech.net.au/pdf2/vici09_171.pdf</a>
C75H-1698A	1/16inx250um 8 Position,w/Pneumatic actuator	Cheminert Valves	\$1,683.00	<a href="http://www.chromtech.net.au/pdf2/vici09_171.pdf">http://www.chromtech.net.au/pdf2/vici09_171.pdf</a>
C75H-1696A	1/16inx250um 6 Position,w/Pneumatic actuator	Cheminert Valves	\$1,620.00	<a href="http://www.chromtech.net.au/pdf2/vici09_171.pdf">http://www.chromtech.net.au/pdf2/vici09_171.pdf</a>
C75N-6C90	1/32inx150um 10 Position,Replace Stator	Cheminert Valves	\$1,294.00	<a href="http://www.chromtech.net.au/pdf2/vici09_171.pdf">http://www.chromtech.net.au/pdf2/vici09_171.pdf</a>
C75N-66R0	1/32inx150um 10 Position,Replace Rotor	Cheminert Valves	\$122.00	<a href="http://www.chromtech.net.au/pdf2/vici09_170.pdf">http://www.chromtech.net.au/pdf2/vici09_170.pdf</a>
C75NX-6698D	1/32inx150um 8 Position,Replace Valve	Cheminert Valves	\$1,612.00	<a href="http://www.chromtech.net.au/pdf2/vici09_170.pdf">http://www.chromtech.net.au/pdf2/vici09_170.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C75NX-6690D	1/32inx150um 10 Position,Replace Valve	Cheminert Valves	\$1,677.00	<a href="http://www.chromtech.net.au/pdf2/vici09_170.pdf">http://www.chromtech.net.au/pdf2/vici09_170.pdf</a>
C75NX-6696A	1/32inx150um 6 Position,w/Pneumatic actuator	Cheminert Valves	\$2,066.00	<a href="http://www.chromtech.net.au/pdf2/vici09_170.pdf">http://www.chromtech.net.au/pdf2/vici09_170.pdf</a>
C75NX-6690EMT	1/32inx150um 10 Position,w/Microelectric actuator	Cheminert Valves	\$2,990.00	<a href="http://www.chromtech.net.au/pdf2/vici09_170.pdf">http://www.chromtech.net.au/pdf2/vici09_170.pdf</a>
C75NX-6698EMT	1/32inx150um 8 Position,w/Microelectric actuator	Cheminert Valves	\$2,926.00	<a href="http://www.chromtech.net.au/pdf2/vici09_170.pdf">http://www.chromtech.net.au/pdf2/vici09_170.pdf</a>
C75NX-6696EMH	1/32inx150um 6 Position,w/Microelectric actuator	Cheminert Valves	\$2,786.00	<a href="http://www.chromtech.net.au/pdf2/vici09_170.pdf">http://www.chromtech.net.au/pdf2/vici09_170.pdf</a>
C75NX-6696D	1/32inx150um 6 Position,Replace Valve	Cheminert Valves	\$1,550.00	<a href="http://www.chromtech.net.au/pdf2/vici09_170.pdf">http://www.chromtech.net.au/pdf2/vici09_170.pdf</a>
C75N-66R8	1/32inx150um 8 Position,Replace Rotor	Cheminert Valves	\$122.00	<a href="http://www.chromtech.net.au/pdf2/vici09_170.pdf">http://www.chromtech.net.au/pdf2/vici09_170.pdf</a>
C75N-6C96	1/32inx150um 6 Position,Replace Stator	Cheminert Valves	\$1,167.00	<a href="http://www.chromtech.net.au/pdf2/vici09_170.pdf">http://www.chromtech.net.au/pdf2/vici09_170.pdf</a>
C75N-6C98	1/32inx150um 8 Position,Replace Stator	Cheminert Valves	\$1,230.00	<a href="http://www.chromtech.net.au/pdf2/vici09_170.pdf">http://www.chromtech.net.au/pdf2/vici09_170.pdf</a>
C75NX-6690E	1/32inx150um 10 Position,w/Std Electric actuator	Cheminert Valves	\$2,709.00	<a href="http://www.chromtech.net.au/pdf2/vici09_170.pdf">http://www.chromtech.net.au/pdf2/vici09_170.pdf</a>
C75NX-6698E	1/32inx150um 8 Position,w/Std Electric actuator	Cheminert Valves	\$2,646.00	<a href="http://www.chromtech.net.au/pdf2/vici09_170.pdf">http://www.chromtech.net.au/pdf2/vici09_170.pdf</a>
C75NX-6696E	1/32inx150um 6 Position,w/Std Electric actuator	Cheminert Valves	\$2,582.00	<a href="http://www.chromtech.net.au/pdf2/vici09_170.pdf">http://www.chromtech.net.au/pdf2/vici09_170.pdf</a>
C75NX-6690A	1/32inx150um 10 Position,w/Pneumatic actuator	Cheminert Valves	\$2,193.00	<a href="http://www.chromtech.net.au/pdf2/vici09_170.pdf">http://www.chromtech.net.au/pdf2/vici09_170.pdf</a>
C75NX-6698A	1/32inx150um 8 Position,w/Pneumatic actuator	Cheminert Valves	\$2,130.00	<a href="http://www.chromtech.net.au/pdf2/vici09_170.pdf">http://www.chromtech.net.au/pdf2/vici09_170.pdf</a>
C75N-66R6	1/32inx150um 6 Position,Replace Rotor	Cheminert Valves	\$122.00	<a href="http://www.chromtech.net.au/pdf2/vici09_170.pdf">http://www.chromtech.net.au/pdf2/vici09_170.pdf</a>
FSR.5-5	Valcon Polyimide 1/32in Removable adaptors tube OD 0.4mm to 0.5mm	Valco Fittings	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_17.pdf">http://www.chromtech.net.au/pdf2/vici09_17.pdf</a>
FS1L.2-5	Valcon Polyimide 1/16in Replacement liners tube OD less than 0.2mm	Valco Fittings	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_17.pdf">http://www.chromtech.net.au/pdf2/vici09_17.pdf</a>
FS1R.4-5	Valcon Polyimide 1/16in Removable adaptors tube OD 0.25mm to 0.4mm	Valco Fittings	\$108.00	<a href="http://www.chromtech.net.au/pdf2/vici09_17.pdf">http://www.chromtech.net.au/pdf2/vici09_17.pdf</a>
LZCN1-10		Valco Fittings	\$57.00	<a href="http://www.chromtech.net.au/pdf2/vici09_17.pdf">http://www.chromtech.net.au/pdf2/vici09_17.pdf</a>
FS1L.15PK-5	PEEK 1/16in Replacement liners tube OD less than 0.15mm	Valco Fittings	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_17.pdf">http://www.chromtech.net.au/pdf2/vici09_17.pdf</a>
ZN.5-10	Stainless steel 1/32in nut	Valco Fittings	\$38.00	<a href="http://www.chromtech.net.au/pdf2/vici09_17.pdf">http://www.chromtech.net.au/pdf2/vici09_17.pdf</a>
ZF1V-5	Valcon Polyimide 1/16in Adapters tube OD 0.5mm to 0.8mm	Valco Fittings	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_17.pdf">http://www.chromtech.net.au/pdf2/vici09_17.pdf</a>
ZF.5V-5	Valcon Polyimide 1/32in Adapters tube OD 0.5mm to 0.8mm	Valco Fittings	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_17.pdf">http://www.chromtech.net.au/pdf2/vici09_17.pdf</a>
FS1R.2PK-5	PEEK 1/16in Removable adaptors tube OD 0.15mm to 0.2mm	Valco Fittings	\$108.00	<a href="http://www.chromtech.net.au/pdf2/vici09_17.pdf">http://www.chromtech.net.au/pdf2/vici09_17.pdf</a>
FS1R.4PK-5	PEEK 1/16in Removable adaptors tube OD 0.25mm to 0.4mm	Valco Fittings	\$108.00	<a href="http://www.chromtech.net.au/pdf2/vici09_17.pdf">http://www.chromtech.net.au/pdf2/vici09_17.pdf</a>
FS1R.15PK-5	PEEK 1/16in Removable adaptors tube OD less than 0.15mm	Valco Fittings	\$108.00	<a href="http://www.chromtech.net.au/pdf2/vici09_17.pdf">http://www.chromtech.net.au/pdf2/vici09_17.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

FS1R.8PK-5	PEEK 1/16in Removable adaptors tube OD 0.5mm to 0.8mm	Valco Fittings	\$108.00	<a href="http://www.chromtech.net.au/pdf2/vici09_17.pdf">http://www.chromtech.net.au/pdf2/vici09_17.pdf</a>
FSR.4-5	Valcon Polyimide 1/32in Removable adaptors tube OD 0.25mm to 0.4mm	Valco Fittings	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_17.pdf">http://www.chromtech.net.au/pdf2/vici09_17.pdf</a>
FSR.3-5	Valcon Polyimide 1/32in Removable adaptors tube OD 0.3mm to 0.35mm	Valco Fittings	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_17.pdf">http://www.chromtech.net.au/pdf2/vici09_17.pdf</a>
FS1R.2-5	Valcon Polyimide 1/16in Removable adaptors tube OD 0.15mm to 0.2mm	Valco Fittings	\$108.00	<a href="http://www.chromtech.net.au/pdf2/vici09_17.pdf">http://www.chromtech.net.au/pdf2/vici09_17.pdf</a>
FS1R.5-5	Valcon Polyimide 1/16in Removable adaptors tube OD 0.4mm to 0.5mm	Valco Fittings	\$108.00	<a href="http://www.chromtech.net.au/pdf2/vici09_17.pdf">http://www.chromtech.net.au/pdf2/vici09_17.pdf</a>
FS1R.5PK-5	PEEK 1/16in Removable adaptors tube OD 0.4mm to 0.5mm	Valco Fittings	\$108.00	<a href="http://www.chromtech.net.au/pdf2/vici09_17.pdf">http://www.chromtech.net.au/pdf2/vici09_17.pdf</a>
FS1R.8-5	Valcon Polyimide 1/16in Removable adaptors tube OD 0.5mm to 0.8mm	Valco Fittings	\$108.00	<a href="http://www.chromtech.net.au/pdf2/vici09_17.pdf">http://www.chromtech.net.au/pdf2/vici09_17.pdf</a>
FSR.25-5	Valcon Polyimide 1/32in Removable adaptors tube OD less than 0.25mm	Valco Fittings	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_17.pdf">http://www.chromtech.net.au/pdf2/vici09_17.pdf</a>
FS1R1.0-5	Removable Fused Silica Adaptor 1/16" One Piece Valcon 0.9-1.0mm ID	Valco Fittings	\$108.00	<a href="http://www.chromtech.net.au/pdf2/vici09_17.pdf">http://www.chromtech.net.au/pdf2/vici09_17.pdf</a>
FS1L.4-5	Valcon Polyimide 1/16in Replacement liners tube OD 0.2mm to 0.4mm	Valco Fittings	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_17.pdf">http://www.chromtech.net.au/pdf2/vici09_17.pdf</a>
FS1L.4PK-5	PEEK 1/16in Replacement liners tube OD 0.2mm to 0.4mm	Valco Fittings	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_17.pdf">http://www.chromtech.net.au/pdf2/vici09_17.pdf</a>
FS1L.5-5	Valcon Polyimide 1/16in liners tube OD 0.4mm to 0.5mm	Valco Fittings	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_17.pdf">http://www.chromtech.net.au/pdf2/vici09_17.pdf</a>
FS1L.5PK-5	PEEK 1/16in liners tube OD 0.4mm to 0.5mm	Valco Fittings	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_17.pdf">http://www.chromtech.net.au/pdf2/vici09_17.pdf</a>
FS1L.8-5	Valcon Polyimide 1/16in liners tube OD 0.5mm to 0.8mm	Valco Fittings	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_17.pdf">http://www.chromtech.net.au/pdf2/vici09_17.pdf</a>
FS1L.8PK-5	PEEK 1/16in liners tube OD 0.5mm to 0.8mm	Valco Fittings	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_17.pdf">http://www.chromtech.net.au/pdf2/vici09_17.pdf</a>
FS1L1.0-5	Valcon Polyimide 1/16in liners tube OD 0.9mm to 1.0mm	Valco Fittings	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_17.pdf">http://www.chromtech.net.au/pdf2/vici09_17.pdf</a>
FS1L1.0PK-5	PEEK 1/16in liners tube OD 0.9mm to 1.0mm	Valco Fittings	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_17.pdf">http://www.chromtech.net.au/pdf2/vici09_17.pdf</a>
ZF1PK-10	PEEK Polymeric Ferrules 1/16in	Valco Fittings	\$42.00	<a href="http://www.chromtech.net.au/pdf2/vici09_17.pdf">http://www.chromtech.net.au/pdf2/vici09_17.pdf</a>
FS1L.2PK-5	PEEK 1/16in Replacement liners tube OD less than 0.2mm	Valco Fittings	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_17.pdf">http://www.chromtech.net.au/pdf2/vici09_17.pdf</a>
ZCN1-10		Valco Fittings	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_17.pdf">http://www.chromtech.net.au/pdf2/vici09_17.pdf</a>
FS1R1.0PK-5	PEEK 1/16in Removable adaptors tube OD 0.9mm to 1.0mm	Valco Fittings	\$108.00	<a href="http://www.chromtech.net.au/pdf2/vici09_17.pdf">http://www.chromtech.net.au/pdf2/vici09_17.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C24-2184-.5E	0.5ul Electric actuator 1/4-28 1/16in tubing 0.50mm ports	Valco Valves	\$1,059.00	<a href="http://www.chromtech.net.au/pdf2/vici09_167.pdf">http://www.chromtech.net.au/pdf2/vici09_167.pdf</a>
C24-2184-2A	2ul Pneumatic actuator 1/4-28 1/16in tubing 0.50mm ports	Valco Valves	\$651.00	<a href="http://www.chromtech.net.au/pdf2/vici09_167.pdf">http://www.chromtech.net.au/pdf2/vici09_167.pdf</a>
C24-2184-2	2ul Manual 1/4-28 1/16in tubing 0.50mm ports	Valco Valves	\$332.00	<a href="http://www.chromtech.net.au/pdf2/vici09_167.pdf">http://www.chromtech.net.au/pdf2/vici09_167.pdf</a>
C24-2184-1E	1ul Electric actuator 1/4-28 1/16in tubing 0.50mm ports	Valco Valves	\$106.00	<a href="http://www.chromtech.net.au/pdf2/vici09_167.pdf">http://www.chromtech.net.au/pdf2/vici09_167.pdf</a>
C24-2184-2E	2ul Electric actuator 1/4-28 1/16in tubing 0.50mm ports	Valco Valves	\$1,059.00	<a href="http://www.chromtech.net.au/pdf2/vici09_167.pdf">http://www.chromtech.net.au/pdf2/vici09_167.pdf</a>
C24-2184-1EH	1ul Microelectric actuator 1/4-28 1/16in tubing 0.50mm ports	Valco Valves	\$1,275.00	<a href="http://www.chromtech.net.au/pdf2/vici09_167.pdf">http://www.chromtech.net.au/pdf2/vici09_167.pdf</a>
C24-2184-2EH	2ul Microelectric actuator 1/4-28 1/16in tubing 0.50mm ports	Valco Valves	\$1,275.00	<a href="http://www.chromtech.net.au/pdf2/vici09_167.pdf">http://www.chromtech.net.au/pdf2/vici09_167.pdf</a>
C24-2184-.5D	0.5ul Replacement valve 1/4-28 1/16in tubing 0.50mm ports	Valco Valves	\$332.00	<a href="http://www.chromtech.net.au/pdf2/vici09_167.pdf">http://www.chromtech.net.au/pdf2/vici09_167.pdf</a>
C24-2184-1D	1ul Replacement valve 1/4-28 1/16in tubing 0.50mm ports	Valco Valves	\$332.00	<a href="http://www.chromtech.net.au/pdf2/vici09_167.pdf">http://www.chromtech.net.au/pdf2/vici09_167.pdf</a>
C24-2184-2D	2ul Replacement valve 1/4-28 1/16in tubing 0.50mm ports	Valco Valves	\$332.00	<a href="http://www.chromtech.net.au/pdf2/vici09_167.pdf">http://www.chromtech.net.au/pdf2/vici09_167.pdf</a>
C24-2184-1A	1ul Pneumatic actuator 1/4-28 1/16in tubing 0.50mm ports	Valco Valves	\$651.00	<a href="http://www.chromtech.net.au/pdf2/vici09_167.pdf">http://www.chromtech.net.au/pdf2/vici09_167.pdf</a>
C24-2184-.5A	0.5ul Pneumatic actuator 1/4-28 1/16in tubing 0.50mm ports	Valco Valves	\$651.00	<a href="http://www.chromtech.net.au/pdf2/vici09_167.pdf">http://www.chromtech.net.au/pdf2/vici09_167.pdf</a>
C24-2184-.5EH	0.5ul Microelectric actuator 1/4-28 1/16in tubing 0.50mm ports	Valco Valves	\$1,275.00	<a href="http://www.chromtech.net.au/pdf2/vici09_167.pdf">http://www.chromtech.net.au/pdf2/vici09_167.pdf</a>
C24-2184-.5	0.5ul Manual 1/4-28 1/16in tubing 0.50mm ports	Valco Valves	\$332.00	<a href="http://www.chromtech.net.au/pdf2/vici09_167.pdf">http://www.chromtech.net.au/pdf2/vici09_167.pdf</a>
C24-10R-.5	0.5ul Replacement rotor 1/16in ZDV fittings 0.40mm ports	Valco Valves	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_167.pdf">http://www.chromtech.net.au/pdf2/vici09_167.pdf</a>
C24-10R-1	1ul Replacement rotor 1/4-28 1/16in tubing 0.50mm ports	Valco Valves	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_167.pdf">http://www.chromtech.net.au/pdf2/vici09_167.pdf</a>
C24-10R-2	2ul Replacement rotor 1/4-28 1/16in tubing 0.50mm ports	Valco Valves	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_167.pdf">http://www.chromtech.net.au/pdf2/vici09_167.pdf</a>
C24-1C8	1ul Replacement stator 1/4-28 1/16in tubing 0.50mm ports	Valco Valves	\$192.00	<a href="http://www.chromtech.net.au/pdf2/vici09_167.pdf">http://www.chromtech.net.au/pdf2/vici09_167.pdf</a>
C24-1C8	1ul Replacement stator 1/4-28 1/16in tubing 0.50mm ports	Valco Valves	\$192.00	<a href="http://www.chromtech.net.au/pdf2/vici09_167.pdf">http://www.chromtech.net.au/pdf2/vici09_167.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C24-1C8	1ul Replacement stator 1/4-28 1/16in tubing 0.50mm ports	Valco Valves	\$192.00	<a href="http://www.chromtech.net.au/pdf2/vici09_167.pdf">http://www.chromtech.net.au/pdf2/vici09_167.pdf</a>
C24-2184-1	1ul Manual 1/4-28 1/16in tubing 0.50mm ports	Valco Valves	\$332.00	<a href="http://www.chromtech.net.au/pdf2/vici09_167.pdf">http://www.chromtech.net.au/pdf2/vici09_167.pdf</a>
C24-10R-.5	0.5ul Replacement rotor 1/16in ZDV fittings 0.40mm ports	Low Pressure Valves	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_166.pdf">http://www.chromtech.net.au/pdf2/vici09_166.pdf</a>
C24Z-2184-.2D	0.2ul Replacement valve 1/16in ZDV fittings 0.40mm ports	Low Pressure Valves	\$332.00	<a href="http://www.chromtech.net.au/pdf2/vici09_166.pdf">http://www.chromtech.net.au/pdf2/vici09_166.pdf</a>
C24Z-2184-.2E	0.2ul Electric actuator 1/16in ZDV fittings 0.40mm ports	Low Pressure Valves	\$1,059.00	<a href="http://www.chromtech.net.au/pdf2/vici09_166.pdf">http://www.chromtech.net.au/pdf2/vici09_166.pdf</a>
C24Z-2184-.5D	0.5ul Replacement valve 1/16in ZDV fittings 0.40mm ports	Low Pressure Valves	\$332.00	<a href="http://www.chromtech.net.au/pdf2/vici09_166.pdf">http://www.chromtech.net.au/pdf2/vici09_166.pdf</a>
C24Z-2184-.5A	0.5ul Pneumatic actuator 1/16in ZDV fittings 0.40mm ports	Low Pressure Valves	\$651.00	<a href="http://www.chromtech.net.au/pdf2/vici09_166.pdf">http://www.chromtech.net.au/pdf2/vici09_166.pdf</a>
C24Z-2184-.2A	0.2ul Pneumatic actuator 1/16in ZDV fittings 0.40mm ports	Low Pressure Valves	\$651.00	<a href="http://www.chromtech.net.au/pdf2/vici09_166.pdf">http://www.chromtech.net.au/pdf2/vici09_166.pdf</a>
C24Z-2184-.5	0.5ul Manual 1/16in ZDV fittings 0.40mm ports	Low Pressure Valves	\$332.00	<a href="http://www.chromtech.net.au/pdf2/vici09_166.pdf">http://www.chromtech.net.au/pdf2/vici09_166.pdf</a>
C24-10R-.2	0.2ul Replacement rotor 1/16in ZDV fittings 0.40mm ports	Low Pressure Valves	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_166.pdf">http://www.chromtech.net.au/pdf2/vici09_166.pdf</a>
C24Z-2184-.2EH	0.2ul Microelectricactuator 1/16in ZDV fittings 0.40mm ports	Low Pressure Valves	\$1,275.00	<a href="http://www.chromtech.net.au/pdf2/vici09_166.pdf">http://www.chromtech.net.au/pdf2/vici09_166.pdf</a>
C24-10R-1	1ul Replacement rotor 1/4-28 1/16in tubing 0.50mm ports	Low Pressure Valves	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_166.pdf">http://www.chromtech.net.au/pdf2/vici09_166.pdf</a>
C24Z-1C8	0.5ul Replacement stator 1/16in ZDV fittings 0.40mm ports	Low Pressure Valves	\$192.00	<a href="http://www.chromtech.net.au/pdf2/vici09_166.pdf">http://www.chromtech.net.au/pdf2/vici09_166.pdf</a>
C24Z-1C8	0.5ul Replacement stator 1/16in ZDV fittings 0.40mm ports	Low Pressure Valves	\$192.00	<a href="http://www.chromtech.net.au/pdf2/vici09_166.pdf">http://www.chromtech.net.au/pdf2/vici09_166.pdf</a>
C24Z-1C8	0.5ul Replacement stator 1/16in ZDV fittings 0.40mm ports	Low Pressure Valves	\$192.00	<a href="http://www.chromtech.net.au/pdf2/vici09_166.pdf">http://www.chromtech.net.au/pdf2/vici09_166.pdf</a>
C24Z-2184-.2	0.2ul Manual 1/16in ZDV fittings 0.40mm ports	Low Pressure Valves	\$332.00	<a href="http://www.chromtech.net.au/pdf2/vici09_166.pdf">http://www.chromtech.net.au/pdf2/vici09_166.pdf</a>
C24Z-2184-1D	1ul Replacement valve 1/16in ZDV fittings 0.40mm ports	Valco Valves	\$332.00	<a href="http://www.chromtech.net.au/pdf2/vici09_166.pdf">http://www.chromtech.net.au/pdf2/vici09_166.pdf</a>
C24Z-2184-1A	1ul Pneumatic actuator 1/16in ZDV fittings 0.40mm ports	Valco Valves	\$651.00	<a href="http://www.chromtech.net.au/pdf2/vici09_166.pdf">http://www.chromtech.net.au/pdf2/vici09_166.pdf</a>
C24Z-2184-.5EH	0.5ul Microelectricactuator 1/16in ZDV fittings 0.40mm ports	Valco Valves	\$1,275.00	<a href="http://www.chromtech.net.au/pdf2/vici09_166.pdf">http://www.chromtech.net.au/pdf2/vici09_166.pdf</a>
C24Z-2184-1	1ul Manual 1/16in ZDV fittings 0.40mm ports	Valco Valves	\$332.00	<a href="http://www.chromtech.net.au/pdf2/vici09_166.pdf">http://www.chromtech.net.au/pdf2/vici09_166.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C24Z-2184-1EH	1ul Microelectricactuator 1/16in ZDV fittings 0.40mm ports	Valco Valves	\$1,275.00	<a href="http://www.chromtech.net.au/pdf2/vici09_166.pdf">http://www.chromtech.net.au/pdf2/vici09_166.pdf</a>
C24Z-2184-.5E	0.5ul Electric actuator 1/16in ZDV fittings 0.40mm ports	Valco Valves	\$1,059.00	<a href="http://www.chromtech.net.au/pdf2/vici09_166.pdf">http://www.chromtech.net.au/pdf2/vici09_166.pdf</a>
C24Z-2184-1E	1ul Electric actuator 1/16in ZDV fittings 0.40mm ports	Valco Valves	\$1,059.00	<a href="http://www.chromtech.net.au/pdf2/vici09_166.pdf">http://www.chromtech.net.au/pdf2/vici09_166.pdf</a>
C22-6184EH	Electric actuator 4 Port 1/8in 1/4-28 fittings 1.5mm ports	Low Pressure Valves	\$1,275.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
C22-618	Replacement rotor 8 Port 1/8in 1/4-28 fittings 1.5mm ports	Low Pressure Valves	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
C22-616	Replacement rotor 6 Port 1/8in 1/4-28 fittings 1.5mm ports	Low Pressure Valves	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
C22-614	Replacement rotor 4 Port 1/8in 1/4-28 fittings 1.5mm ports	Low Pressure Valves	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
C22-6188D	Replacement valve 8 Port 1/8in 1/4-28 fittings 1.5mm ports	Low Pressure Valves	\$447.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
C22-6186D	Replacement valve 6 Port 1/8in 1/4-28 fittings 1.5mm ports	Low Pressure Valves	\$332.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
C22-6184D	Replacement valve 4 Port 1/8in 1/4-28 fittings 1.5mm ports	Low Pressure Valves	\$332.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
CFSL1KFEP	FEP Sample Loops-1ml	Low Pressure Valves	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
C22-6186EH	Electric actuator 6 Port 1/8in 1/4-28 fittings 1.5mm ports	Low Pressure Valves	\$1,275.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
C22-680	Replacement stator 10 Port 1/8in 1/4-28 fittings 1.5mm ports	Low Pressure Valves	\$338.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
C22-6188E	Electric actuator 8 Port 1/8in 1/4-28 fittings 1.5mm ports	Low Pressure Valves	\$1,173.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
C22-6186E	Electric actuator 6 Port 1/8in 1/4-28 fittings 1.5mm ports	Low Pressure Valves	\$1,059.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
C22-6184E	Electric actuator 4 Port 1/8in 1/4-28 fittings 1.5mm ports	Low Pressure Valves	\$1,059.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
C22-6188A	Pneumatic actuator 8 Port 1/8in 1/4-28 fittings 1.5mm ports	Low Pressure Valves	\$765.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
C22-6186A	Pneumatic actuator 6 Port 1/8in 1/4-28 fittings 1.5mm ports	Low Pressure Valves	\$651.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
C22-6184A	Pneumatic actuator 4 Port 1/8in 1/4-28 fittings 1.5mm ports	Low Pressure Valves	\$651.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>



## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C22-6188EH	Electric actuator 8 Port 1/8in 1/4-28 fittings 1.5mm ports	Low Pressure Valves	\$1,389.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
CFSL50PK	PEEK Sample Loops-50ul	Low Pressure Valves	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
CFSL500PK	PEEK Sample Loops-500ul	Low Pressure Valves	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
CFSL500TF	PTFE Sample Loops-500ul	Low Pressure Valves	\$26.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
CFSL500FEP	FEP Sample Loops-250ul	Low Pressure Valves	\$26.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
CFSL250PK	PEEK Sample Loops-250ul	Low Pressure Valves	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
CFSL250TF	PTFE Sample Loops-250ul	Low Pressure Valves	\$22.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
CFSL250FEP	FEP Sample Loops-250ul	Low Pressure Valves	\$22.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
C22-610	Replacement rotor 10 Port 1/8in 1/4-28 fittings 1.5mm ports	Low Pressure Valves	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
CFSL100TF	PTFE Sample Loops-100ul	Low Pressure Valves	\$22.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
C22-684	Replacement stator 4 Port 1/8in 1/4-28 fittings 1.5mm ports	Low Pressure Valves	\$192.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
CFSL50TF	PTFE Sample Loops-50ul	Low Pressure Valves	\$22.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
CFSL50FEP	FEP Sample Loops-50ul	Low Pressure Valves	\$22.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
CFSL20PK	PEEK Sample Loops-20ul	Low Pressure Valves	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
CFSL20TF	PTFE Sample Loops-20ul	Low Pressure Valves	\$22.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
CFSL20FEP	FEP Sample Loops-20ul	Low Pressure Valves	\$22.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
C22-6186	6 Port 1/8in 1/4-28 fittings 1.5mm ports,Manual	Low Pressure Valves	\$332.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
CFSL100PK	PEEK Sample Loops-100ul	Low Pressure Valves	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
C22-3186	Manual 6 Port 1/8in 1/4-28 fittings 1.5mm ports	Low Pressure Valves	\$332.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
C22-3186E	Electric actuator 6 Port 1/16in 1/4-28 fittings 0.75mm ports	Low Pressure Valves	\$1,059.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
C22-3184E	Electric actuator 4 Port 1/16in 1/4-28 fittings 0.75mm ports	Low Pressure Valves	\$1,059.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
C22-3180A	Pneumatic actuator 10 Port 1/16in 1/4-28 fittings 0.75mm ports	Low Pressure Valves	\$796.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
C22-3188A	Pneumatic actuator 8 Port 1/16in 1/4-28 fittings 0.75mm ports	Low Pressure Valves	\$765.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
C22-3186A	Pneumatic actuator 6 Port 1/16in 1/4-28 fittings 0.75mm ports	Low Pressure Valves	\$651.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
C22-3184A	Pneumatic actuator 4 Port 1/16in 1/4-28 fittings 0.75mm ports	Low Pressure Valves	\$651.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
C22-6180	10 Port 1/8in 1/4-28 fittings 1.5mm ports,Manual	Low Pressure Valves	\$478.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
C22-3188	Manual 8 Port 1/16in 1/4-28 fittings 0.75mm ports	Low Pressure Valves	\$447.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C22-3184EH	Electric actuator 4 Port 1/16in 1/4-28 fittings 0.75mm ports	Low Pressure Valves	\$1,275.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
C22-3184	Manual 4 Port 1/16in 1/4-28 fittings 0.75mm ports	Low Pressure Valves	\$332.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
CFSL2KPK	PEEK Sample Loops-2ml	Low Pressure Valves	\$80.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
CFSL2KTF	FEP Sample Loops-2ml	Low Pressure Valves	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
CFSL2KFEP	FEP Sample Loops-2ml	Low Pressure Valves	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
CFSL1KPK	PEEK Sample Loops-1ml	Low Pressure Valves	\$51.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
CFSL1KTF	PTFE Sample Loops-1ml	Low Pressure Valves	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
C22-3180	Manual 10 Port 1/16in 1/4-28 fittings 0.75mm ports	Low Pressure Valves	\$478.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
C22-3180D	Replacement valve 10 Port 1/16in 1/4-28 fittings 0.75mm ports	Low Pressure Valves	\$478.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
C22-386	Replacement stator 6 Port 1/16in 1/4-28 fittings 0.75mm ports	Low Pressure Valves	\$192.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
C22-6184	4 Port 1/8in 1/4-28 fittings 1.5mm ports,Manual	Low Pressure Valves	\$332.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
C22-380	Replacement stator 10 Port 1/16in 1/4-28 fittings 0.75mm ports	Low Pressure Valves	\$338.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
C22-388	Replacement stator 8 Port 1/16in 1/4-28 fittings 0.75mm ports	Low Pressure Valves	\$306.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
C22-384	Replacement stator 4 Port 1/16in 1/4-28 fittings 0.75mm ports	Low Pressure Valves	\$192.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
C22-318	Replacement rotor 8 Port 1/16in 1/4-28 fittings 0.75mm ports	Low Pressure Valves	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
C22-3188E	Electric actuator 8 Port 1/16in 1/4-28 fittings 0.75mm ports	Low Pressure Valves	\$1,173.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
C22-314	Replacement rotor 4 Port1/16in 1/4-28 fittings 0.75mm ports	Low Pressure Valves	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
C22-3180E	Electric actuator 10 Port 1/16in 1/4-28 fittings 0.75mm ports	Low Pressure Valves	\$1,204.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
C22-3188D	Replacement valve 8 Port 1/16in 1/4-28 fittings 0.75mm ports	Low Pressure Valves	\$447.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
C22-3186D	Replacement valve 6 Port 1/16in 1/4-28 fittings 0.75mm ports	Low Pressure Valves	\$332.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
C22-3184D	Replacement valve 4 Port 1/16in 1/4-28 fittings 0.75mm ports	Low Pressure Valves	\$332.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
C22-3180EH	Electric actuator 10 Port 1/16in 1/4-28 fittings 0.75mm ports	Low Pressure Valves	\$1,422.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C22-3188EH	Electric actuator 8 Port 1/16in 1/4-28 fittings 0.75mm ports	Low Pressure Valves	\$1,389.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
C22-3186EH	Electric actuator 6 Port 1/16in 1/4-28 fittings 0.75mm ports	Low Pressure Valves	\$1,275.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
C22-6188	8 Port 1/8in 1/4-28 fittings 1.5mm ports,Manual	Low Pressure Valves	\$447.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
C22-316	Replacement rotor 6 Port 1/16in 1/4-28 fittings 0.75mm ports	Low Pressure Valves	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
C22-310	Replacement rotor 10 Port 1/16in 1/4-28 fittings 0.75mm ports	Low Pressure Valves	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
C22-6180A	Pneumatic actuator 10 Port 1/8in 1/4-28 fittings 1.5mm ports	Valco Valves	\$796.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
C22-6180E	Electric actuator 10 Port 1/8in 1/4-28 fittings 1.5mm ports	Valco Valves	\$1,204.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
C22-6180EH	Electric actuator 10 Port 1/8in 1/4-28 fittings 1.5mm ports	Valco Valves	\$1,422.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
C22-688	Replacement stator 8 Port 1/8in 1/4-28 fittings 1.5mm ports	Valco Valves	\$306.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
CFSL100FEP	FEP Sample Loops-100ul	Valco Valves	\$22.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
C22-6180D	Replacement valve 10 Port 1/8in 1/4-28 fittings 1.5mm ports	Valco Valves	\$478.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
C22-686	Replacement stator 6 Port 1/8in 1/4-28 fittings 1.5mm ports	Valco Valves	\$192.00	<a href="http://www.chromtech.net.au/pdf2/vici09_165.pdf">http://www.chromtech.net.au/pdf2/vici09_165.pdf</a>
CZSL20TF	20ul PTFE Sample loop	Low Pressure Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
CZSL50PK	50ul PEEK Sample loop	Low Pressure Valves	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
CZSL5TF	5ul PTFE Sample loop	Low Pressure Valves	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
CZSL50TF	50ul PTFE Sample loop	Low Pressure Valves	\$24.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
CZSL5PK	5ul PEEK Sample loops for C1 C2 and C3 inc fits	Low Pressure Valves	\$40.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
CZSL10PK	10ul PEEK Sample loops for C1 C2 and C3 inc fits	Low Pressure Valves	\$40.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
CZSL20PK	20ul PEEK Sample loops for C1 C2 and C3 inc fits	Low Pressure Valves	\$34.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
CZSL10FEP	100ul TEF Sample loop	Low Pressure Valves	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
CZSL20FEP	10ul FEP Sample loop	Low Pressure Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
CZSL100FEP	100ul FEP Sample loop	Low Pressure Valves	\$24.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
CZSL500FEP	500ul FEP Sample loop	Low Pressure Valves	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
CZSL10TF	10ul TEF Sample loop	Low Pressure Valves	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
CZSL50FEP	50ul FEP Sample loop	Low Pressure Valves	\$24.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
CZSL100TF	100ul PTFE Sample loop	Low Pressure Valves	\$24.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

CZSL100PK	100ul PEEK Sample loop	Low Pressure Valves	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
CZSL250FEP	250ul FEP Sample loop	Low Pressure Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
CZSL250PK	250ul PEEK Sample loops for C1 C2 and C3 inc fits	Low Pressure Valves	\$40.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
CZSL500TF	500ul PTFE Sample loop	Low Pressure Valves	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
CZSL500PK	500ul PEEK Sample loop	Low Pressure Valves	\$48.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
CZSL1KFEP	1ml FEP Sample loop	Low Pressure Valves	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
CZSL1KTF	1ml PTFE Sample loop	Low Pressure Valves	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
CZSL1KPK	1ml PEEK Sample loop	Low Pressure Valves	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
CZSL5FEP	5ul FEP Sample loop	Low Pressure Valves	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
C22Z-3184E	Electric actuator 4 Port 1/16in ZDV fittings 0.75mm ports	Low Pressure Valves	\$1,059.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
CZSL2KTF	2ml PTFE Sample loop	Low Pressure Valves	\$48.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
CZSL2KFEP	2ml FEP Sample loop	Low Pressure Valves	\$48.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
CZSL250TF	250ul PTFE Sample loop	Low Pressure Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
C22Z-3188EH	Electric actuator 8 Port 1/16in ZDV fittings 0.75mm ports	Low Pressure Valves	\$1,314.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
C22Z-3186EH	Electric actuator 6 Port 1/16in ZDV fittings 0.75mm ports	Low Pressure Valves	\$1,275.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
CZSL2KPK	2ml PEEK Sample loop	Low Pressure Valves	\$87.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
C22Z-3184	Manual 4 Port 1/16in ZDV fittings 0.75mm ports	Low Pressure Valves	\$332.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
C22Z-3186	Manual 6 Port 1/16in ZDV fittings 0.75mm ports	Low Pressure Valves	\$332.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
C22Z-3188	Manual 8 Port 1/16in ZDV fittings 0.75mm ports	Low Pressure Valves	\$369.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
C22Z-3180	Manual 10 Port 1/16in ZDV fittings 0.75mm ports	Low Pressure Valves	\$402.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
C22Z-3184A	Pneumatic actuator 4 Port 1/16in ZDV fittings 0.75mm ports	Low Pressure Valves	\$651.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
C22Z-3186A	Pneumatic actuator 6 Port 1/16in ZDV fittings 0.75mm ports	Low Pressure Valves	\$651.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
C22Z-3188A	Pneumatic actuator 8 Port 1/16in ZDV fittings 0.75mm ports	Low Pressure Valves	\$688.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
C22Z-3180A	Pneumatic actuator 10 Port 1/16in ZDV fittings 0.75mm ports	Low Pressure Valves	\$720.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
C22Z-3186E	Electric actuator 6 Port 1/16in ZDV fittings 0.75mm ports	Low Pressure Valves	\$1,059.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
C22Z-3188E	Electric actuator 8 Port 1/16in ZDV fittings 0.75mm ports	Low Pressure Valves	\$1,096.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C22Z-3184EH	Electric actuator 4 Port 1/16in ZDV fittings 0.75mm ports	Low Pressure Valves	\$1,275.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
C22Z-380	Replacement stator 10 Port 1/16in ZDV fittings 0.75mm ports	Low Pressure Valves	\$261.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
C22Z-3184D	Replacement valve 4 Port 1/16in ZDV fittings 0.75mm ports	Low Pressure Valves	\$332.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
C22Z-3186D	Replacement valve 6 Port 1/16in ZDV fittings 0.75mm ports	Low Pressure Valves	\$332.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
C22Z-3188D	Replacement valve 8 Port 1/16in ZDV fittings 0.75mm ports	Low Pressure Valves	\$369.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
C22Z-3180D	Replacement valve 10 Port 1/16in ZDV fittings 0.75mm ports	Low Pressure Valves	\$402.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
C12-314	Replacement rotor 4 Port 1/16in ZDV fittings 0.75mm ports	Low Pressure Valves	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
C12-316	Replacement rotor 6 Port 1/16in ZDV fittings 0.75mm ports	Low Pressure Valves	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
C12-318	Replacement rotor 8 Port 1/16in ZDV fittings 0.75mm ports	Low Pressure Valves	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
C12-310	Replacement rotor 10 Port 1/16in ZDV fittings 0.75mm ports	Low Pressure Valves	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
C22Z-384	Replacement stator 4 Port 1/16in ZDV fittings 0.75mm ports	Low Pressure Valves	\$192.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
C22Z-386	Replacement stator 6 Port 1/16in ZDV fittings 0.75mm ports	Low Pressure Valves	\$192.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
C22Z-388	Replacement stator 8 Port 1/16in ZDV fittings 0.75mm ports	Low Pressure Valves	\$230.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
C22Z-3180E	Electric actuator 10 Port 1/16in ZDV fittings 0.75mm ports	Low Pressure Valves	\$1,128.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
C22Z-3180EH	Electric actuator 10 Port 1/16in ZDV fittings 0.75mm ports	Low Pressure Valves	\$1,346.00	<a href="http://www.chromtech.net.au/pdf2/vici09_164.pdf">http://www.chromtech.net.au/pdf2/vici09_164.pdf</a>
C6-2C06	1/16inx0.4mm Continuous Flow 6 Port N60 SS Stator,w/Microelectric actuator	Analytical HPLC Injectors	\$345.00	<a href="http://www.chromtech.net.au/pdf2/vici09_163.pdf">http://www.chromtech.net.au/pdf2/vici09_163.pdf</a>
C6-2346A	1/16inx0.4mm Continuous Flow 6 Port PAEK stator,Manual	Analytical HPLC Injectors	\$1,020.00	<a href="http://www.chromtech.net.au/pdf2/vici09_163.pdf">http://www.chromtech.net.au/pdf2/vici09_163.pdf</a>
C6-2006E	1/16inx0.4mm Continuous Flow 6 Port N60 SS Stator,w/Std Electric actuator	Analytical HPLC Injectors	\$1,287.00	<a href="http://www.chromtech.net.au/pdf2/vici09_163.pdf">http://www.chromtech.net.au/pdf2/vici09_163.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C6-2346E	1/16inx0.4mm Continuous Flow 6 Port PAEK stator,Replace Valve	Analytical HPLC Injectors	\$1,428.00	<a href="http://www.chromtech.net.au/pdf2/vici09_163.pdf">http://www.chromtech.net.au/pdf2/vici09_163.pdf</a>
C6-2006EH	1/16inx0.4mm Continuous Flow 6 Port N60 SS Stator,w/Microelectric actuator	Analytical HPLC Injectors	\$1,504.00	<a href="http://www.chromtech.net.au/pdf2/vici09_163.pdf">http://www.chromtech.net.au/pdf2/vici09_163.pdf</a>
C6-2346EH	1/16inx0.4mm Continuous Flow 6 Port PAEK stator,w/Std Electric actuator	Analytical HPLC Injectors	\$1,644.00	<a href="http://www.chromtech.net.au/pdf2/vici09_163.pdf">http://www.chromtech.net.au/pdf2/vici09_163.pdf</a>
C6-2C46	1/16inx0.4mm Continuous Flow 6 Port PAEK stator,w/Microelectric actuator	Analytical HPLC Injectors	\$484.00	<a href="http://www.chromtech.net.au/pdf2/vici09_163.pdf">http://www.chromtech.net.au/pdf2/vici09_163.pdf</a>
C2-23R6	Analytical 20ul 0.40mm ports PAEK stator Replacement rotor 1/16in fittings	Analytical HPLC Injectors	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_163.pdf">http://www.chromtech.net.au/pdf2/vici09_163.pdf</a>
C6-2006A	1/16inx0.4mm Continuous Flow 6 Port N60 SS Stator,Pneumatic actuator	Analytical HPLC Injectors	\$879.00	<a href="http://www.chromtech.net.au/pdf2/vici09_163.pdf">http://www.chromtech.net.au/pdf2/vici09_163.pdf</a>
C6-2006D	1/16inx0.4mm Continuous Flow 6 Port N60 SS Stator,Replace Valve	Analytical HPLC Injectors	\$561.00	<a href="http://www.chromtech.net.au/pdf2/vici09_163.pdf">http://www.chromtech.net.au/pdf2/vici09_163.pdf</a>
C6-2346D	1/16inx0.4mm Continuous Flow 6 Port PAEK stator,Pneumatic actuator	Analytical HPLC Injectors	\$702.00	<a href="http://www.chromtech.net.au/pdf2/vici09_163.pdf">http://www.chromtech.net.au/pdf2/vici09_163.pdf</a>
C-2C46	PAEK Stator Replacement stator 6 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$420.00	<a href="http://www.chromtech.net.au/pdf2/vici09_163.pdf">http://www.chromtech.net.au/pdf2/vici09_163.pdf</a>
C1-2006	1/16in fittings -N60 stainless stator Analytical 0.40mm ports	Analytical HPLC Injectors	\$688.00	<a href="http://www.chromtech.net.au/pdf2/vici09_163.pdf">http://www.chromtech.net.au/pdf2/vici09_163.pdf</a>
C1-2346	1/16in fittings -PAEK stator Analytical 0.40mm ports	Analytical HPLC Injectors	\$828.00	<a href="http://www.chromtech.net.au/pdf2/vici09_163.pdf">http://www.chromtech.net.au/pdf2/vici09_163.pdf</a>
C2-20R6	Titanium stator Replacement rotor 6 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_163.pdf">http://www.chromtech.net.au/pdf2/vici09_163.pdf</a>
C6-2346	1/16inx0.4mm Continuous Flow 6 Port N60 SS Stator,Manual	Analytical HPLC Injectors	\$702.00	<a href="http://www.chromtech.net.au/pdf2/vici09_163.pdf">http://www.chromtech.net.au/pdf2/vici09_163.pdf</a>
C-2C06	1/16in fittings -Replacement stator Analytical 0.40mm ports	Analytical HPLC Injectors	\$280.00	<a href="http://www.chromtech.net.au/pdf2/vici09_163.pdf">http://www.chromtech.net.au/pdf2/vici09_163.pdf</a>
C1CF-2006	1/16inx0.4mm Continuous Thru-Handle 6 Port N60 SS Stator	Analytical HPLC Injectors	\$75.00	<a href="http://www.chromtech.net.au/pdf2/vici09_163.pdf">http://www.chromtech.net.au/pdf2/vici09_163.pdf</a>
C1CF-2346	1/16inx0.4mm Continuous Thru-Handle 6 Port PAEK Stator	Analytical HPLC Injectors	\$898.00	<a href="http://www.chromtech.net.au/pdf2/vici09_163.pdf">http://www.chromtech.net.au/pdf2/vici09_163.pdf</a>
C6-2006	1/16inx0.4mm Continuous Flow 6 Port N60 SS Stator,Manual	Analytical HPLC Injectors	\$561.00	<a href="http://www.chromtech.net.au/pdf2/vici09_163.pdf">http://www.chromtech.net.au/pdf2/vici09_163.pdf</a>
C1-20R6	1/16in fittings -Replacement rotor Analytical 0.40mm ports	Analytical HPLC Injectors	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_163.pdf">http://www.chromtech.net.au/pdf2/vici09_163.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C1-23R6	1/16in fittings -Replacement rotor Analytical 0.40mm ports	Analytical HPLC Injectors	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_163.pdf">http://www.chromtech.net.au/pdf2/vici09_163.pdf</a>
C-261	Injector Fitting for Thru the Handle Manual HPLC Injector valve	Analytical HPLC Injectors	\$50.00	<a href="http://www.chromtech.net.au/pdf2/vici09_163.pdf">http://www.chromtech.net.au/pdf2/vici09_163.pdf</a>
C4-1004-.1E	0.1ul N60 stainless stator Electric actuator 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$1,428.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-1344-.2	0.1ul PAEK stator Manual 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$842.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-1344-.1	0.1ul PAEK stator Manual 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$842.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-1C0	0.5ul N60 stainless stator Replacement stator 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$490.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-1C0	0.5ul N60 stainless stator Replacement stator 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$490.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-10R-.1	0.1ul N60 stainless stator Replacement rotor 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-1004-.2D	0.2ul N60 stainless stator Replacement valve 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$702.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-1004-.1D	0.1ul N60 stainless stator Replacement valve 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$702.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-1004-.5EH	0.5ul N60 stainless stator Microelectricactuator 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$1,644.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-1004-.2EH	0.2ul N60 stainless stator Microelectricactuator 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$1,644.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-1004-.1EH	0.1ul N60 stainless stator Microelectricactuator 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$1,644.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-1034-.1EH	0.1ul Titanium stator Microelectric actuator 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$2,014.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-1004-.2E	0.2ul N60 stainless stator Electric actuator 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$1,428.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-1344-.2A	0.1ul PAEK stator Pneumatic actuator 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$1,161.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-1004-.5A	0.5ul N60 stainless stator Pneumatic actuator 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$1,020.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-1004-.2A	0.2ul N60 stainless stator Pneumatic actuator 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$1,020.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C4-1004-.1A	0.1ul N60 stainless stator Pneumatic actuator 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$1,020.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-1004-.5	0.5ul N60 stainless stator Manual 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$702.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-1004-.2	0.2ul N60 stainless stator Manual 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$702.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-1004-.1	0.1ul N60 stainless stator Manual 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$702.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-1034-.1E	0.1ul Titanium stator Electric actuator 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$1,926.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-1034-.2E	0.2ul Titanium stator Electric actuator 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$1,926.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-1034-.5E	0.5ul Titanium stator Electric actuator 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$1,926.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-1004-.5E	0.5ul N60 stainless stator Electric actuator 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$1,428.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-1344-.5D	0.5ul PAEK stator Replacement valve 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$842.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-1034-.5A	0.5ul Titanium stator Pneumatic actuator 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$1,389.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-1034-.2A	0.2ul Titanium stator Pneumatic actuator 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$1,389.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-1034-.1A	0.1ul Titanium stator Pneumatic actuator 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$1,389.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-1034-.5	0.5ul Titanium stator Manual N60 stainless stator 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$1,071.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-1034-.2	0.2ul Titanium stator Manual 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$1,071.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-1034-.1	0.1ul Titanium stator Manual 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$1,071.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-1C4	0.1ul PAEK stator Replacement stator 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$632.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-1C4	0.1ul PAEK stator Replacement stator 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$632.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-1C4	0.1ul PAEK stator Replacement stator 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$632.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>



## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C4-13R-.5	0.5ul PAEK stator Replacement rotor 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-1344-.5	0.1ul PAEK stator Pneumatic actuator 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$842.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-13R-.1	0.1ul PAEK stator Replacement rotor 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-1344-.1A	0.1ul PAEK stator Pneumatic actuator 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$1,161.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-1344-.2D	0.1ul PAEK stator Replacement valve 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$842.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-1344-.1D	0.1ul PAEK stator Replacement valve 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$842.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-1344-.5EH	0.5ul PAEK stator Microelectricactuator 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$1,912.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-1344-.2EH	0.1ul PAEK stator Microelectric actuator 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$1,912.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-1344-.1EH	0.1ul PAEK stator Microelectric actuator 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$1,912.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-1344-.5E	0.5ul PAEK stator Electric actuator 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$1,569.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-1344-.2E	0.1ul PAEK stator Electric actuator 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$1,569.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-1344-.1E	0.1ul PAEK stator Electric actuator 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$1,569.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-1344-.5A	0.5ul PAEK stator Pneumatic actuator 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$1,161.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-1004-.5D	0.5ul N60 stainless stator Replacement valve 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$702.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-13R-.2	0.2ul PAEK stator Replacement rotor 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-1C0	0.5ul N60 stainless stator Replacement stator 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$490.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-1034-.2EH	0.2ul Titanium stator Microelectricactuator 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$2,014.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-1C3	0.1ul Titanium stator Replacement stator 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$861.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C4-1C3	0.1ul Titanium stator Replacement stator 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$861.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-1C3	0.1ul Titanium stator Replacement stator 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$861.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-10R-.5	0.5ul N60 stainless stator Replacement rotor 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-1034-.5D	0.5ul Titanium stator Replacement valve 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$1,071.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-1034-.2D	0.2ul Titanium stator Replacement valve 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$1,071.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-1034-.1D	0.1ul Titanium stator Replacement valve 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$1,071.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-1034-.5EH	0.5ul Titanium stator Microelectric actuator 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$2,014.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C4-10R-.2	0.1ul N60 stainless stator Replacement rotor 1/16in fittings 0.25mm ports	Analytical HPLC Injectors	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_162.pdf">http://www.chromtech.net.au/pdf2/vici09_162.pdf</a>
C-2C44	PAEK Stator Replacement stator 4 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$420.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2-2034A	Titanium stator Pneumatic actuator 4 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$1,179.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2H-2030	Titanium Stator Manual 10 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$1,000.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2H-2038	Titanium Stator Manual 8 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$930.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2-2036	Titanium stator Manual 4 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$861.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2-2034	Titanium stator Manual 4 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$861.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C-2C40H	PAEK Stator Replacement stator 10 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$561.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2H-2038A	Titanium Stator Pneumatic actuator 8 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$1,250.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C-2C46	PAEK Stator Replacement stator 6 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$420.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2H-2030A	Titanium Stator Pneumatic actuator 10 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$1,320.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C2-23R0H	PAEK Stator 10 Port 1/16in fittings 0.40mm ports, Replace Rotor	Analytical HPLC Injectors	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2-23R8H	PAEK Stator Replacement rotor 8 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2-23R6	PAEK Stator Replacement rotor 6 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2-23R4	PAEK Stator Replacement rotor 4 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2H-2340D	PAEK Stator Replacement valve 10 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$771.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2H-2348D	PAEK Stator Replacement valve 8 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$702.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2-2346D	PAEK Stator Replacement valve 6 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$632.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2-2344D	PAEK Stator Replacement valve 4 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$632.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C-2C48H	PAEK Stator Replacement stator 8 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$490.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2-2034D	Titanium stator Replacement valve 4 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$861.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C-2C38H	Titanium stator Replacement stator 8 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$720.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C-2C36	Titanium stator Replacement stator 6 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$651.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C-2C34	Titanium stator Replacement stator 4 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$651.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2-20R0H	Titanium stator 6 Port 1/16in fittings 0.40mm ports, Replace Rotor	Analytical HPLC Injectors	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2-20R8H	Titanium stator Replacement rotor 8 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2-20R4	Titanium stator Replacement rotor 4 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2H-2030D	Titanium Stator Replacement valve 10 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$1,000.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2-2036A	Titanium stator Pneumatic actuator 6 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$1,179.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C2-2036D	Titanium stator Replacement valve 6 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$861.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2-2346EH	PAEK Stator Microelectric actuator 6 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$1,575.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2H-2030EH	Titanium Stator Microelectric actuator 10 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$1,944.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2H-2038EH	Titanium Stator Microelectric actuator 8 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$2,001.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2-2036EH	Titanium stator Microelectric actuator 6 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$1,932.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2-2034EH	Titanium stator Microelectric actuator 4 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$1,932.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2H-2030E	Titanium Stator Electric actuator 10 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$1,728.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2H-2038E	Titanium Stator Electric actuator 8 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$1,658.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2-2036E	Titanium stator Electric actuator 6 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$1,587.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2-2034E	Titanium stator Electric actuator 4 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$1,587.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2H-2038D	Titanium Stator Replacement valve 8 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$930.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C-2C30H	N60 SS Stator Replacement stator 10 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$790.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C-2C04	N60 SS Stator Replacement stator 4 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$280.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2H-2348EH	PAEK Stator Microelectric actuator 8 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$1,644.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2-20R4	Titanium stator Replacement rotor 4 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2H-2340EH	PAEK Stator Microelectric actuator 10 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$1,714.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2H-2008D	N60 SS Stator Stator Replacement valve 8 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$561.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C-2C06	N60 SS Stator Replacement stator 6 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$280.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C2H-2008EH	N60 SS Stator Stator Microelectric actuator 8 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$1,504.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2-20R6	Titanium stator Replacement rotor 6 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2H-2008A	N60 SS Stator Stator Pneumatic actuator 8 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$879.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2H-2000E	N60 SS Stator Stator Electric actuator 10 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$1,358.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2H-2008E	N60 SS Stator Stator Electric actuator 8 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$1,287.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2H-2000A	N60 SS Stator Stator Pneumatic actuator 10 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$950.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2H-2000	N60 SS Stator Stator Manual 10 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$632.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2H-2008	N60 SS Stator Stator Manual 8 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$561.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2H-2000EH	N60 SS Stator Stator Microelectric actuator 10 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$1,575.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2-2346A	PAEK Stator Pneumatic actuator 6 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$950.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2H-2348E	PAEK Stator Electric actuator 8 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$1,428.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2-2346E	PAEK Stator Electric actuator 6 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$1,358.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2-2344E	PAEK Stator Electric actuator 4 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$1,358.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2H-2000D	N60 SS Stator Stator Replacement valve 10 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$632.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2-2344EH	PAEK Stator Microelectric actuator 4 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$1,575.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C-2C08H	N60 SS Stator Replacement stator 8 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$351.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2H-2340E	PAEK Stator Electric actuator 10 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$1,498.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2H-2348A	PAEK Stator Pneumatic actuator 8 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$1,020.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C2H-2340A	PAEK Stator Pneumatic actuator 10 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$1,090.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2-2344A	PAEK Stator Pneumatic actuator 4 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$950.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2H-2340	PAEK Stator Manual 10 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$771.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2H-2348	PAEK Stator Manual 8 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$702.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2-2346	PAEK Stator Manual 4 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$632.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2-2344	PAEK Stator Manual 4 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$632.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C-2C00H	N60 SS Stator Replacement rotor 10 Port 1/16in fittings 0.40mm ports	Analytical HPLC Injectors	\$420.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2-2006	N60 stainless stator Manual 4 Port 1/16in fittings 0.40mm ports	Microbore HPLC 1/16	\$490.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2-2004A	N60 stainless stator Pneumatic actuator 4 Port 1/16in fittings 0.40mm ports	Microbore HPLC 1/16	\$810.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2-2004	N60 stainless stator Manual 4 Port 1/16in fittings 0.40mm ports	Microbore HPLC 1/16	\$490.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2-2004E	N60 stainless stator Electric actuator 4 Port 1/16in fittings 0.40mm ports	Microbore HPLC 1/16	\$1,218.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2-2006E	N60 stainless stator Electric actuator 6 Port 1/16in fittings 0.40mm ports	Microbore HPLC 1/16	\$1,218.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2-2006EH	N60 stainless stator Microelectric actuator 6 Port 1/16in fittings 0.40mm ports	Microbore HPLC 1/16	\$1,434.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2-2004D	N60 stainless stator Replacement valve 4 Port 1/16in fittings 0.40mm ports	Microbore HPLC 1/16	\$490.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2-2006D	N60 stainless stator Replacement valve 6 Port 1/16in fittings 0.40mm ports	Microbore HPLC 1/16	\$490.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2-2004EH	N60 stainless stator Microelectric actuator 4 Port 1/16in fittings 0.40mm ports	Microbore HPLC 1/16	\$1,434.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C2-2006A	N60 stainless stator Pneumatic actuator 4 Port 1/16in fittings 0.40mm ports	Microbore HPLC 1/16	\$810.00	<a href="http://www.chromtech.net.au/pdf2/vici09_161.pdf">http://www.chromtech.net.au/pdf2/vici09_161.pdf</a>
C1-1006	6 Port 1/16inx0.25mm Thru-Handle Injector,N60 SS Stator	Microbore HPLC 1/16	\$898.00	<a href="http://www.chromtech.net.au/pdf2/vici09_160.pdf">http://www.chromtech.net.au/pdf2/vici09_160.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C6-1006A	1/16x0.25mm Cont Flow,N6 OSS Stator,w/Pneumatic actuator	Microbore HPLC 1/16	\$1,090.00	<a href="http://www.chromtech.net.au/pdf2/vici09_160.pdf">http://www.chromtech.net.au/pdf2/vici09_160.pdf</a>
C2-10R6	Microbore 5ul 0.25mm ports Titanium stator Replacement rotor 1/16in fittings	Microbore HPLC 1/16	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_160.pdf">http://www.chromtech.net.au/pdf2/vici09_160.pdf</a>
C2-13R6	Microbore 5ul 0.25mm ports PEEK stator Replacement rotor 1/16in fittings	Microbore HPLC 1/16	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_160.pdf">http://www.chromtech.net.au/pdf2/vici09_160.pdf</a>
C6-1006D	1/16x0.25mm Cont Flow,N6 OSS Stator,Replace Valve	Microbore HPLC 1/16	\$771.00	<a href="http://www.chromtech.net.au/pdf2/vici09_160.pdf">http://www.chromtech.net.au/pdf2/vici09_160.pdf</a>
C6-1346EH	1/16x0.25mm Cont Flow,PAEK Stator,w/MicroElectric actuator	Microbore HPLC 1/16	\$1,983.00	<a href="http://www.chromtech.net.au/pdf2/vici09_160.pdf">http://www.chromtech.net.au/pdf2/vici09_160.pdf</a>
C6-1006EH	1/16x0.25mm Cont Flow,N6 OSS Stator,w/Microelectric actuator	Microbore HPLC 1/16	\$1,714.00	<a href="http://www.chromtech.net.au/pdf2/vici09_160.pdf">http://www.chromtech.net.au/pdf2/vici09_160.pdf</a>
C6-1346E	1/16x0.25mm Cont Flow,PAEK Stator,w/Std Electric actuator	Microbore HPLC 1/16	\$1,638.00	<a href="http://www.chromtech.net.au/pdf2/vici09_160.pdf">http://www.chromtech.net.au/pdf2/vici09_160.pdf</a>
C6-1006E	1/16x0.25mm Cont Flow,N6 OSS Stator,w/Std Electric actuator	Microbore HPLC 1/16	\$1,498.00	<a href="http://www.chromtech.net.au/pdf2/vici09_160.pdf">http://www.chromtech.net.au/pdf2/vici09_160.pdf</a>
C6-1346D	1/16x0.25mm Cont Flow,PAEK Stator,Replace Valve	Microbore HPLC 1/16	\$912.00	<a href="http://www.chromtech.net.au/pdf2/vici09_160.pdf">http://www.chromtech.net.au/pdf2/vici09_160.pdf</a>
C1-10R6	6 Port 1/16inx0.25mm Thru-Handle Injector,N60 SS Stator;Replace Rotor	Microbore HPLC 1/16	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_160.pdf">http://www.chromtech.net.au/pdf2/vici09_160.pdf</a>
C-1C06	6 Port 1/16inx0.25mm Thru-Handle Injector,Replace N60SS Stator	Microbore HPLC 1/16	\$490.00	<a href="http://www.chromtech.net.au/pdf2/vici09_160.pdf">http://www.chromtech.net.au/pdf2/vici09_160.pdf</a>
C-1C46	1/16in fittings -Replacement stator Microbore 0.25mm ports	Microbore HPLC 1/16	\$632.00	<a href="http://www.chromtech.net.au/pdf2/vici09_160.pdf">http://www.chromtech.net.au/pdf2/vici09_160.pdf</a>
C-261	Injector Fittin for Thru the Handle Manual HPLC Injector valve	Microbore HPLC 1/16	\$50.00	<a href="http://www.chromtech.net.au/pdf2/vici09_160.pdf">http://www.chromtech.net.au/pdf2/vici09_160.pdf</a>
C1CF-1006	6 Port 1/16inx0.25mm Thru-Handle Injector,Cont Flow,N60 SS Stator	Microbore HPLC 1/16	\$969.00	<a href="http://www.chromtech.net.au/pdf2/vici09_160.pdf">http://www.chromtech.net.au/pdf2/vici09_160.pdf</a>
C1CF-1346	6 Port 1/16inx0.25mm Thru-Handle Injector,Cont Flow,PAEL Stator Stator	Microbore HPLC 1/16	\$1,110.00	<a href="http://www.chromtech.net.au/pdf2/vici09_160.pdf">http://www.chromtech.net.au/pdf2/vici09_160.pdf</a>
C1-13R6	6 Port 1/16inx0.25mm Thru-Handle Injector,PAEK Stator,Replace Rotor,	Microbore HPLC 1/16	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_160.pdf">http://www.chromtech.net.au/pdf2/vici09_160.pdf</a>
C6-1C46	1/16x0.25mm Cont Flow,Replace PEEK Stator	Microbore HPLC 1/16	\$702.00	<a href="http://www.chromtech.net.au/pdf2/vici09_160.pdf">http://www.chromtech.net.au/pdf2/vici09_160.pdf</a>
C6-1346	1/16x0.25mm Cont Flow,PAEK Stator,Manual	Microbore HPLC 1/16	\$912.00	<a href="http://www.chromtech.net.au/pdf2/vici09_160.pdf">http://www.chromtech.net.au/pdf2/vici09_160.pdf</a>
C6-1C06	1/16x0.25mm Cont Flow,Replace N60 SS Stator	Microbore HPLC 1/16	\$561.00	<a href="http://www.chromtech.net.au/pdf2/vici09_160.pdf">http://www.chromtech.net.au/pdf2/vici09_160.pdf</a>
C6-1346A	1/16x0.25mm Cont Flow,PAEK Stator,Pneumatic actuator	Microbore HPLC 1/16	\$1,230.00	<a href="http://www.chromtech.net.au/pdf2/vici09_160.pdf">http://www.chromtech.net.au/pdf2/vici09_160.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C-261	Injector Fittin for Thru the Handle Manual HPLC Injector valve	Microbore HPLC 1/16	\$50.00	<a href="http://www.chromtech.net.au/pdf2/vici09_160.pdf">http://www.chromtech.net.au/pdf2/vici09_160.pdf</a>
C1-1346	6 Port 1/16inx0.25mm Thru-Handle Injector, PEEK Stator	Microbore HPLC 1/16	\$1,040.00	<a href="http://www.chromtech.net.au/pdf2/vici09_160.pdf">http://www.chromtech.net.au/pdf2/vici09_160.pdf</a>
FS1.8PK-5	PEEK 1/16in Adapters tube OD 0.5mm to 0.8mm	Valco Fittings	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_16.pdf">http://www.chromtech.net.au/pdf2/vici09_16.pdf</a>
FS1.3PK-5	1/16",0.30to0.35mm,PEEK,5pk	Valco Fittings	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_16.pdf">http://www.chromtech.net.au/pdf2/vici09_16.pdf</a>
FS1.3V1-5	1/16",0.30to0.35mm,PolyI,5pk	Valco Fittings	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_16.pdf">http://www.chromtech.net.au/pdf2/vici09_16.pdf</a>
FS1.4-5	Valcon Polyimide 1/16in Adapters tube OD 0.25mm to 0.4mm	Valco Fittings	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_16.pdf">http://www.chromtech.net.au/pdf2/vici09_16.pdf</a>
FS1.4PK-5	PEEK 1/16in Adapters tube OD 0.25mm to 0.4mm	Valco Fittings	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_16.pdf">http://www.chromtech.net.au/pdf2/vici09_16.pdf</a>
FS1.4V1-5	Virgin Polyimide 1/16in Adapters tube OD 0.25mm to 0.4mm	Valco Fittings	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_16.pdf">http://www.chromtech.net.au/pdf2/vici09_16.pdf</a>
FS1.5-5	Fused Silica Adaptor 1/16" One Piece 0.4-0.5mm ID	Valco Fittings	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_16.pdf">http://www.chromtech.net.au/pdf2/vici09_16.pdf</a>
FS1.3-5	1/16",0.30to0.35mm,PolyI/Valcon,5pk	Valco Fittings	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_16.pdf">http://www.chromtech.net.au/pdf2/vici09_16.pdf</a>
FS1.5V1-5	Virgin Polyimide 1/16in Adapters tube OD 0.4mm to 0.5mm	Valco Fittings	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_16.pdf">http://www.chromtech.net.au/pdf2/vici09_16.pdf</a>
FS1.8-5	1/16",0.50to0.80mm,PolyI/Valcon,5pk	Valco Fittings	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_16.pdf">http://www.chromtech.net.au/pdf2/vici09_16.pdf</a>
FS1.9-5	1/16",0.80to0.90mm,PolyI/Valcon,5pk	Valco Fittings	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_16.pdf">http://www.chromtech.net.au/pdf2/vici09_16.pdf</a>
FS1.9PK-5	1/16",0.80to0.90mm,PEEK,5pk	Valco Fittings	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_16.pdf">http://www.chromtech.net.au/pdf2/vici09_16.pdf</a>
FS1.9V1-5	1/16"0.80to0.90mm,PolyI,5pk	Valco Fittings	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_16.pdf">http://www.chromtech.net.au/pdf2/vici09_16.pdf</a>
FS11.0-5	Fused Silica Adaptor 1/16" One Piece 0.9-1.0mm ID	Valco Fittings	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_16.pdf">http://www.chromtech.net.au/pdf2/vici09_16.pdf</a>
FS11.0PK-5	1/16",0.90to1.0mm,PEEK,5pk	Valco Fittings	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_16.pdf">http://www.chromtech.net.au/pdf2/vici09_16.pdf</a>
FS11.0V1-5	1/16"0.90to1.00mm,PolyI,5pk	Valco Fittings	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_16.pdf">http://www.chromtech.net.au/pdf2/vici09_16.pdf</a>
FS1.5PK-5	PEEK 1/16in Adapters tube OD 0.4mm to 0.5mm	Valco Fittings	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_16.pdf">http://www.chromtech.net.au/pdf2/vici09_16.pdf</a>
FS.5-5	Valcon Polyimide 1/32in Adapters tube OD 0.25mm to 0.5mm 5-pk	Valco Fittings	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_16.pdf">http://www.chromtech.net.au/pdf2/vici09_16.pdf</a>
FS.25-5	Valcon Polyimide One-piece 1/32in Adapters tube OD to 0.25mm	Valco Fittings	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_16.pdf">http://www.chromtech.net.au/pdf2/vici09_16.pdf</a>
FS1.8V1-5	Virgin Polyimide 1/16in Adapters tube OD 0.5mm to 0.8mm	Valco Fittings	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_16.pdf">http://www.chromtech.net.au/pdf2/vici09_16.pdf</a>
FS1.25V1-5	Virgin Polyimide 1/16in Adapters tube OD to 0.25mm	Valco Fittings	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_16.pdf">http://www.chromtech.net.au/pdf2/vici09_16.pdf</a>
FS.25PK-5	PEEK 1/32in Adapters tube OD to 0.25mm	Valco Fittings	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_16.pdf">http://www.chromtech.net.au/pdf2/vici09_16.pdf</a>



## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

FS.25V1-5	Virgin Polyimide 1/32in Adapters tube OD to 0.25mm	Valco Fittings	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_16.pdf">http://www.chromtech.net.au/pdf2/vici09_16.pdf</a>
FS.4-5	Valcon Polyimide 1/32in Adapters tube OD 0.25mm to 0.4mm	Valco Fittings	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_16.pdf">http://www.chromtech.net.au/pdf2/vici09_16.pdf</a>
FS.4V1-5	Virgin Polyimide 1/32in Adapters tube OD 0.25mm to 0.4mm	Valco Fittings	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_16.pdf">http://www.chromtech.net.au/pdf2/vici09_16.pdf</a>
FS.5PK-5	PEEK 1/32in Adapters tube OD 0.4mm to 0.5mm	Valco Fittings	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_16.pdf">http://www.chromtech.net.au/pdf2/vici09_16.pdf</a>
FS1.2PK-5	1/16", <0.25mm, PEEK, 5pk	Valco Fittings	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_16.pdf">http://www.chromtech.net.au/pdf2/vici09_16.pdf</a>
FS1.25PK-5	PEEK 1/16in Adapters tube OD to 0.25mm	Valco Fittings	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_16.pdf">http://www.chromtech.net.au/pdf2/vici09_16.pdf</a>
FS1.25-5	Valcon Polyimide 1/16in Adapters tube OD to 0.25mm	Valco Fittings	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_16.pdf">http://www.chromtech.net.au/pdf2/vici09_16.pdf</a>
FS.4PK-5	PEEK 1/32in Adapters tube OD 0.25mm to 0.4mm	Valco Fittings	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_16.pdf">http://www.chromtech.net.au/pdf2/vici09_16.pdf</a>
FS1.2V1-5	1/16", <0.25mm, PolyI, 5pk	Valco Fittings	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_16.pdf">http://www.chromtech.net.au/pdf2/vici09_16.pdf</a>
FS.5V1-5	Virgin Polyimide 1/32in Adapters tube OD 0.4mm to 0.5mm	Valco Fittings	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_16.pdf">http://www.chromtech.net.au/pdf2/vici09_16.pdf</a>
FS.25-5	Valcon Polyimide One-piece 1/32in Adapters tube OD to 0.25mm	Valco Fittings	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_16.pdf">http://www.chromtech.net.au/pdf2/vici09_16.pdf</a>
ZF.5V1-5	Virgin Polyimide 1/32in Adapters tube OD 0.5mm to 0.8mm	Valco Fittings	\$74.00	<a href="http://www.chromtech.net.au/pdf2/vici09_16.pdf">http://www.chromtech.net.au/pdf2/vici09_16.pdf</a>
ZF.5PK-5	PEEK 1/32in Adapters tube OD 0.5mm to 0.8mm	Valco Fittings	\$27.00	<a href="http://www.chromtech.net.au/pdf2/vici09_16.pdf">http://www.chromtech.net.au/pdf2/vici09_16.pdf</a>
ZF.5V-5	Valcon Polyimide 1/32in Adapters tube OD 0.5mm to 0.8mm	Valco Fittings	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_16.pdf">http://www.chromtech.net.au/pdf2/vici09_16.pdf</a>
CSL50	50ul Stainless Steel Sample loops for C1 C2 and C3 inc fits	Microbore HPLC 1/16	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
CZSL20PK	20ul PEEK Sample loops for C1 C2 and C3 inc fits	Microbore HPLC 1/16	\$34.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
CSL20	20ul 1/16in Stainless steel loops inc fits	Microbore HPLC 1/16	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
CSL10TI	10ul Titanium Sample loops for C1 C2 and C3 inc fits	Microbore HPLC 1/16	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
CZSL50PK	50ul PEEK Sample loop	Microbore HPLC 1/16	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
CZSL10PK	10ul PEEK Sample loops for C1 C2 and C3 inc fits	Microbore HPLC 1/16	\$40.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
CSL20TI	20ul Titanium Sample loops for C1 C2 and C3 inc fits	Microbore HPLC 1/16	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
C4-0C4	0.02ul PAEK stator Replacement stator 1/16in fittings 0.15mm ports	Microbore HPLC 1/16	\$842.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
CZSL5PK	5ul PEEK Sample loops for C1 C2 and C3 inc fits	Microbore HPLC 1/16	\$40.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

CSL5	5ul Stainless Steel Sample loops for C1 C2 and C3 inc fits	Microbore HPLC 1/16	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
CZSL2PK	2ul PEEK Sample loops for C1 C2 and C3 inc fits	Microbore HPLC 1/16	\$40.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
CSL2	2ul Stainless Steel Sample loops for C1 C2 and C3 inc fits	Microbore HPLC 1/16	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
CSL50TI	50ul Sample Loops, Tit for C1,C2,C2V,C3,C6 valves	Microbore HPLC 1/16	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
CSL1K	1ml Stainless Steel Sample loops for C1 C2 and C3 inc fits	Microbore HPLC 1/16	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
C4-0C4	0.02ul PAEK stator Replacement stator 1/16in fittings 0.15mm ports	Microbore HPLC 1/16	\$842.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
C4-03R-.05	0.05ul PAEK stator Replacement rotor 1/16in fittings 0.15mm ports	Microbore HPLC 1/16	\$105.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
C4-0C4	0.02ul PAEK stator Replacement stator 1/16in fittings 0.15mm ports	Microbore HPLC 1/16	\$842.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
CSL10	10ul Stainless Steel Sample loops for C1 C2 and C3 inc fits	Microbore HPLC 1/16	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
CZSL1KPK	1ml PEEK Sample loop	Microbore HPLC 1/16	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
CSL10K	10ml 1/16in Stainless steel loops inc fits	Microbore HPLC 1/16	\$87.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
CZSL5KPK	5ml PEEK Sample loops for C1 C2 and C3 inc fits	Microbore HPLC 1/16	\$135.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
C4-0004-.05D	0.05ul Replacement valve N60 stainless stator 1/16in fittings 0.15mm ports	Microbore HPLC 1/16	\$912.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
C4-03R-.02	0.02ul PAEK stator Replacement rotor 1/16in fittings 0.15mm ports	Microbore HPLC 1/16	\$105.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
CSL5K	5ml 1/16in Stainless steel loops inc fits	Microbore HPLC 1/16	\$68.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
CZSL2KPK	2ml PEEK Sample loop	Microbore HPLC 1/16	\$87.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
CZSL500PK	500ul PEEK Sample loop	Microbore HPLC 1/16	\$48.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
CSL1KTI	1ml Titanium Sample loops for C1 C2 and C3 inc fits	Microbore HPLC 1/16	\$302.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
CSL100	100ul 1/16in Stainless steel loops inc fits	Microbore HPLC 1/16	\$40.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
CSL500TI	500ul 1/16in Titanium loops inc fits	Microbore HPLC 1/16	\$192.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
CSL500	500ul 1/16in Stainless steel loops inc fits	Microbore HPLC 1/16	\$40.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
CSL250TI	250ul Titanium Sample loops for C1 C2 and C3 inc fits	Microbore HPLC 1/16	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
CZSL250PK	250ul PEEK Sample loops for C1 C2 and C3 inc fits	Microbore HPLC 1/16	\$40.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
CSL250	250ul Stainless Steel Sample loops for C1 C2 and C3 inc fits	Microbore HPLC 1/16	\$40.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

CSL100TI	100ul Titanium Sample loops for C1 C2 and C3 inc fits	Microbore HPLC 1/16	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
CZSL100PK	100ul PEEK Sample loop	Microbore HPLC 1/16	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
CSL2K	2ml 1/16in Stainless steel loops inc fits	Microbore HPLC 1/16	\$57.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
C4-0004-.01E	0.01ul Electric actuator N60 stainless stator 1/16in fittings 0.15mm ports	Microbore HPLC 1/16	\$1,638.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
C4-00R-.02	0.02ul Replacement rotor N60 stainless stator 1/16in fittings 0.15mm ports	Microbore HPLC 1/16	\$105.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
C4-00R-.01	0.01ul Replacement rotor N60 stainless stator 1/16in fittings 0.15mm ports	Microbore HPLC 1/16	\$105.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
C4-0004-.02D	0.02ul Replacement valve N60 stainless stator 1/16in fittings 0.15mm ports	Microbore HPLC 1/16	\$912.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
C4-0004-.01D	0.01ul Replacement valve N60 stainless stator 1/16in fittings 0.15mm ports	Microbore HPLC 1/16	\$912.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
C4-0004-.05EH	0.05 icroelectricactuator N60 stainless stator 1/16in fittings 0.15mm ports	Microbore HPLC 1/16	\$1,983.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
C4-0004-.02EH	0.02 Microelectricactuator N60 stainless stator 1/16in fittings 0.15mm ports	Microbore HPLC 1/16	\$1,983.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
C4-0C0	0.02ul Replacement stator N60 stainless stator 1/16in fittings 0.15mm ports	Microbore HPLC 1/16	\$702.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
C4-0004-.05E	0.05ul Electric actuator N60 stainless stator 1/16in fittings 0.15mm ports	Microbore HPLC 1/16	\$1,638.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
C4-0C0	0.02ul Replacement stator N60 stainless stator 1/16in fittings 0.15mm ports	Microbore HPLC 1/16	\$702.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
C4-0004-.02A	0.02ul Pneumatic actuator N60 stainless stator 1/16in fittings 0.15mm ports	Microbore HPLC 1/16	\$1,230.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
C4-0004-.01A	0.01ul Pneumatic actuator N60 stainless stator 1/16in fittings 0.15mm ports	Microbore HPLC 1/16	\$1,230.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
C4-0004-.05	0.05 Manual N60 stainless stator 1/16in fittings 0.15mm ports	Microbore HPLC 1/16	\$912.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
C4-0004-.02	0.02 Manual N60 stainless stator 1/16in fittings 0.15mm ports	Microbore HPLC 1/16	\$912.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
C4-0004-.01	0.01 Manual N60 stainless stator 1/16in fittings 0.15mm ports	Microbore HPLC 1/16	\$912.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
C4-0004-.02E	0.02ul Electric actuator N60 stainless stator 1/16in fittings 0.15mm ports	Microbore HPLC 1/16	\$1,638.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C1-13R6	1/16in fittings -Replacement rotor Microbore 0.25mm ports	Microbore HPLC 1/16	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
C4-0004-.01EH	0.01 Microelectricactuator N60 stainless stator 1/16in fittings 0.15mm ports	Microbore HPLC 1/16	\$1,983.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
C4-0344-.01E	0.01ul PAEK stator Electric actuator 1/16in fittings 0.15mm ports	Microbore HPLC 1/16	\$1,779.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
C4-0344-.05D	0.05ul PAEK stator Replacement valve 1/16in fittings 0.15mm ports	Microbore HPLC 1/16	\$1,052.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
C4-0344-.02D	0.02ul PAEK stator Replacement valve 1/16in fittings 0.15mm ports	Microbore HPLC 1/16	\$1,052.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
C4-0344-.01D	0.01ul PAEK stator Replacement valve 1/16in fittings 0.15mm ports	Microbore HPLC 1/16	\$1,052.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
C4-0344-.05EH	0.05ul PAEK stator Microelectricactuator 1/16in fittings 0.15mm ports	Microbore HPLC 1/16	\$1,995.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
C4-0344-.02EH	0.02ul PAEK stator Microelectricactuator 1/16in fittings 0.15mm ports	Microbore HPLC 1/16	\$1,995.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
C4-0344-.01EH	0.01ul PAEK stator Microelectricactuator 1/16in fittings 0.15mm ports	Microbore HPLC 1/16	\$1,995.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
C4-00R-.05	0.05ul Replacement rotor N60 stainless stator 1/16in fittings 0.15mm ports	Microbore HPLC 1/16	\$105.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
C4-0344-.02E	0.02ul PAEK stator Electric actuator 1/16in fittings 0.15mm ports	Microbore HPLC 1/16	\$1,779.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
C4-03R-.01	0.01ul PAEK stator Replacement rotor 1/16in fittings 0.15mm ports	Microbore HPLC 1/16	\$105.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
C4-0344-.05A	0.05ul PAEK stator Pneumatic actuator 1/16in fittings 0.15mm ports	Microbore HPLC 1/16	\$1,371.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
C4-0344-.02A	0.02ul PAEK stator Pneumatic actuator 1/16in fittings 0.15mm ports	Microbore HPLC 1/16	\$1,371.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
C4-0344-.01A	0.01ul PAEK stator Pneumatic actuator 1/16in fittings 0.15mm ports	Microbore HPLC 1/16	\$1,371.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
C4-0344-.05	0.05ul PAEK stator Manual 1/16in fittings 0.15mm ports	Microbore HPLC 1/16	\$1,052.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
C4-0344-.02	0.02ul PAEK stator Manual 1/16in fittings 0.15mm ports	Microbore HPLC 1/16	\$1,052.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
C4-0344-.01	0.01ul PAEK stator Manual 1/16in fittings 0.15mm ports	Microbore HPLC 1/16	\$1,052.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C4-0C0	0.02ul Replacement stator N60 stainless stator 1/16in fittings 0.15mm ports	Microbore HPLC 1/16	\$702.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
C4-0344-.05E	0.05ul PAEK stator Electric actuator 1/16in fittings 0.15mm ports	Microbore HPLC 1/16	\$1,779.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
C4-0004-.05A	0.05ul Pneumatic actuator N60 stainless stator 1/16in fittings 0.15mm ports	Microbore HPLC 1/16	\$1,230.00	<a href="http://www.chromtech.net.au/pdf2/vici09_159.pdf">http://www.chromtech.net.au/pdf2/vici09_159.pdf</a>
C-1C48H	8 Port 1/16inx0.25mm,Replace PAEK Stator	Microbore HPLC 1/16	\$702.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2-1034A	4 Port 1/16inx0.25mm,w/Pneumatic actuator Tit Stator	Microbore HPLC 1/16	\$1,389.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2H-1340E	10 Port 1/16inx0.25mm,w/Std Electric actuator,PAEK Stator	Microbore HPLC 1/16	\$1,708.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2H-1038	8 Port 1/16inx0.25mm,Manual Tit Stator	Microbore HPLC 1/16	\$1,142.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2-1036	6 Port 1/16inx0.25mm,Manual Tit SS Stator	Microbore HPLC 1/16	\$1,071.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2-1034	4 Port 1/16inx0.25mm,Manual Tit Stator	Microbore HPLC 1/16	\$1,071.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2-1004A	4 Port 1/16inx0.25mm,w/Pneumatic actuator N60 SS Stator	Microbore HPLC 1/16	\$1,020.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2H-1000	10 Port 1/16inx0.25mm,Manual, N60 SS Stator	Microbore HPLC 1/16	\$842.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2H-1008	8 Port 1/16inx0.25mm,Manual, N60 SS Stator	Microbore HPLC 1/16	\$771.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C-1C44	4 Port 1/16inx0.25mm,Replace PAEK Stator	Microbore HPLC 1/16	\$632.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C-1C40H	10 Port 1/16inx0.25mm,Replace PAEK Stator	Microbore HPLC 1/16	\$771.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2-13R0H	10 Port 1/16inx0.25mm,w/Replace Rotor, PAEK Stator	Microbore HPLC 1/16	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C-1C46	6 Port 1/16inx0.25mm,Replace PAEK Stator	Microbore HPLC 1/16	\$632.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2-1006	6 Port 1/16inx0.25mm,Manual N60 SS Stator	Microbore HPLC 1/16	\$702.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2H-1030	10 Port 1/16inx0.25mm,Manual, Tit Stator	Microbore HPLC 1/16	\$1,212.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2H-1348E	8 Port 1/16inx0.25mm,w/Std Electric actuator,PAEK Stator	Microbore HPLC 1/16	\$1,638.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C-1C30H	10 Port 1/16inx0.25mm,Replace Tit Stator	Microbore HPLC 1/16	\$1,000.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C-1C38H	8 Port 1/16inx0.25mm,Replace Tit Stator	Microbore HPLC 1/16	\$930.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C-1C36	6 Port 1/16inx0.25mm,Replace Tit Stator	Microbore HPLC 1/16	\$861.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C-1C34	4 Port 1/16inx0.25mm,Replace Tit Stator	Microbore HPLC 1/16	\$861.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2-10R0H	10 Port 1/16inx0.25mm,Replace N60 SS Stator	Microbore HPLC 1/16	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2-10R6	6 Port 1/16inx0.25mm,Replacement rotor, Tit Stator	Microbore HPLC 1/16	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2H-1030D	10 Port 1/16inx0.25mm,Replace Valve,Tit Stator	Microbore HPLC 1/16	\$1,212.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2H-1038D	8 Port 1/16inx0.25mm,Replace Valve,Tit Stator	Microbore HPLC 1/16	\$1,142.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C2H-1038EH	8 Port 1/16inx0.25mm,w/Microelectric actuator,Tit Stator	Microbore HPLC 1/16	\$2,085.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2-1034D	4 Port 1/16inx0.25mm,w/Replace valve Tit Stator	Microbore HPLC 1/16	\$1,071.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2H-1038A	8 Port 1/16inx0.25mm,w/Pneumatic actuator,Tit Stator	Microbore HPLC 1/16	\$1,460.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2H-1030EH	10 Port 1/16inx0.25mm,w/Microelectric actuator,Tit Stator	Microbore HPLC 1/16	\$2,154.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2-1006A	6 Port 1/16inx0.25mm,w/Pneumatic actuator N60 SS Stator	Microbore HPLC 1/16	\$1,020.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2-1036EH	6 Port 1/16inx0.25mm,w/Microelectric actuator Tit Stator	Microbore HPLC 1/16	\$2,014.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2-1034EH	4 Port 1/16inx0.25mm,w/Microelectric actuator Tit Stator	Microbore HPLC 1/16	\$2,014.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2H-1030E	10 Port 1/16inx0.25mm,w/Std Electric actuator,Tit Stator	Microbore HPLC 1/16	\$1,938.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2H-1038E	8 Port 1/16inx0.25mm,w/Std Electric actuator,Tit Stator	Microbore HPLC 1/16	\$1,995.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2-1036E	6 Port 1/16inx0.25mm,w/Std Electric actuator Tit Stator	Microbore HPLC 1/16	\$1,926.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2-1034E	4 Port 1/16inx0.25mm,w/Std Electric actuator Tit Stator	Microbore HPLC 1/16	\$1,926.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2H-1030A	10 Port 1/16inx0.25mm,w/Pneumatic actuator,Tit Stator	Microbore HPLC 1/16	\$1,530.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2-1036D	6 Port 1/16inx0.25mm,w/Replace valve Tit Stator	Microbore HPLC 1/16	\$1,071.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2-13R4	4 Port 1/16inx0.25mm,w/Replace Rotor, PAEK Stator	Microbore HPLC 1/16	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2H-1348	8 Port 1/16inx0.25mm,Manual, PAEK Stator	Microbore HPLC 1/16	\$912.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2-1004	4 Port 1/16inx0.25mm,Manual N60 SS Stator	Microbore HPLC 1/16	\$702.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2H-1340	10 Port 1/16inx0.25mm,Manual, PAEK Stator	Microbore HPLC 1/16	\$981.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2-1344A	4 Port 1/16inx0.25mm,w/Pneumatic actuator PAEK Rotor	Microbore HPLC 1/16	\$1,161.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2-1346A	6 Port 1/16inx0.25mm,w/Pneumatic actuator PAEK Rotor	Microbore HPLC 1/16	\$1,161.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2H-1348A	8 Port 1/16inx0.25mm,w/Pneumatic actuator,PAEK Stator	Microbore HPLC 1/16	\$1,230.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2H-1340A	10 Port 1/16inx0.25mm,w/Pneumatic actuator,PAEK Stator	Microbore HPLC 1/16	\$1,300.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C2-1346E	6 Port 1/16inx0.25mm,w/Std Electric actuator PAEK Rotor	Microbore HPLC 1/16	\$1,569.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2-1346	6 Port 1/16inx0.25mm,Manual PAEK Rotor	Microbore HPLC 1/16	\$842.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2-13R6	6 Port 1/16inx0.25mm,w/Replace Rotor, PAEK Stator	Microbore HPLC 1/16	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2-1346D	6 Port 1/16inx0.25mm,w/Replace Valve PAEK Rotor	Microbore HPLC 1/16	\$842.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2H-1340D	10 Port 1/16inx0.25mm,Replace Valve,PAEK Stator	Microbore HPLC 1/16	\$981.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2H-1348D	8 Port 1/16inx0.25mm,Replace Valve,PAEK Stator	Microbore HPLC 1/16	\$912.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2-1344D	4 Port 1/16inx0.25mm,w/Replace Valve PAEK Rotor	Microbore HPLC 1/16	\$842.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2H-1348EH	8 Port 1/16inx0.25mm,w/Microelectric actuator,PAEK Stator	Microbore HPLC 1/16	\$1,983.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2-1346EH	6 Port 1/16inx0.25mm,w/Microelectric actuator PAEK Rotor	Microbore HPLC 1/16	\$1,912.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2-1344EH	4 Port 1/16inx0.25mm,w/Microelectric actuator PAEK Rotor	Microbore HPLC 1/16	\$1,912.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2H-1008A	8 Port 1/16inx0.25mm,w/Pneumatic actuator, N60 SS Stator	Microbore HPLC 1/16	\$1,090.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2-1036A	6 Port 1/16inx0.25mm,w/Pneumatic actuator Tit Stator	Microbore HPLC 1/16	\$1,389.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2-13R8H	8 Port 1/16inx0.25mm,w/Replace Rotor, PAEK Stator	Microbore HPLC 1/16	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2H-1000E	10 Port 1/16inx0.25mm,w/Std Electric actuator, N60 SS Stator	Microbore HPLC 1/16	\$1,569.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2H-1000A	10 Port 1/16inx0.25mm,w/Pneumatic actuator, N60 SS Stator	Microbore HPLC 1/16	\$1,161.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2-1004E	4 Port 1/16inx0.25mm,w/Std Electric actuator N60 SS Stator	Microbore HPLC 1/16	\$1,428.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2H-1340EH	10 Port 1/16inx0.25mm,w/Microelectric actuator,PAEK Stator	Microbore HPLC 1/16	\$1,926.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2-1344	4 Port 1/16inx0.25mm,Manual PAEK Stator	Microbore HPLC 1/16	\$842.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2-1006E	6 Port 1/16inx0.25mm,w/Std Electric actuator N60 SS Stator	Microbore HPLC 1/16	\$1,428.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2H-1008E	8 Port 1/16inx0.25mm,Replace Valve, N60 SS Stator	Microbore HPLC 1/16	\$1,498.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C2-1004EH	4 Port 1/16inx0.25mm,w/Microelectric actuator N60 SS Stator	Microbore HPLC 1/16	\$1,644.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2-1006EH	6 Port 1/16inx0.25mm,w/Microelectric actuator N60 SS Stator	Microbore HPLC 1/16	\$1,644.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2H-1008EH	8 Port 1/16inx0.25mm,w/Microelectric actuator, N60 SS Stator	Microbore HPLC 1/16	\$1,714.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2H-1000EH	10 Port 1/16inx0.25mm,w/Microelectric actuator, N60 SS Stator	Microbore HPLC 1/16	\$1,912.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2-1004D	4 Port 1/16inx0.25mm,w/Replace valve N60 SS Stator	Microbore HPLC 1/16	\$702.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C-1C08H	8 Port 1/16inx0.25mm,Replace N60SS Stator	Microbore HPLC 1/16	\$561.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C-1C06	6 Port 1/16inx0.25mm,Replace N60SS Stator	Microbore HPLC 1/16	\$490.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2-1006D	6 Port 1/16inx0.25mm,w/Replace valve N60 SS Stator	Microbore HPLC 1/16	\$702.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C-1C04	4 Port 1/16inx0.25mm,Replace N60SS Stator	Microbore HPLC 1/16	\$490.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2-10R4	4 Port 1/16inx0.25mm,Replacement rotor, Tit Stator	Microbore HPLC 1/16	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2-10R8H	8 Port 1/16inx0.25mm,Replacement rotor, Tit Stator	Microbore HPLC 1/16	\$92.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2H-1000D	10 Port 1/16inx0.25mm,Replace Valve, N60 SS Stator	Microbore HPLC 1/16	\$842.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2-1344E	4 Port 1/16inx0.25mm,w/Std Electric actuator PAEK Rotor	Microbore HPLC 1/16	\$1,569.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C2H-1008D	8 Port 1/16inx0.25mm,Replace Valve, N60 SS Stator	Microbore HPLC 1/16	\$771.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C-1C00H	10 Port 1/16inx0.25mm,Replace N60SS Stator	Microbore HPLC 1/16	\$632.00	<a href="http://www.chromtech.net.au/pdf2/vici09_158.pdf">http://www.chromtech.net.au/pdf2/vici09_158.pdf</a>
C74X-1694-.05EH	1/16inx0.25m Internal Sample 20nl,w/Microelectric actuator,10,000psi	Microbore UHPLC 1/16	\$2,142.00	<a href="http://www.chromtech.net.au/pdf2/vici09_157.pdf">http://www.chromtech.net.au/pdf2/vici09_157.pdf</a>
C74H-1694-.02A	1/16inx0.25m Internal Sample 20nl,w/Pneumatic actuator,10,000psi	Microbore UHPLC 1/16	\$1,326.00	<a href="http://www.chromtech.net.au/pdf2/vici09_157.pdf">http://www.chromtech.net.au/pdf2/vici09_157.pdf</a>
C74H-1694-.01A	1/16inx0.25m Internal Sample 10nl,w/Pneumatic actuator,10,000psi	Microbore UHPLC 1/16	\$1,326.00	<a href="http://www.chromtech.net.au/pdf2/vici09_157.pdf">http://www.chromtech.net.au/pdf2/vici09_157.pdf</a>
C74X-1694-.01E	1/16inx0.25m Internal Sample 10nl,w/Std Electric actuator,10,000psi	Microbore UHPLC 1/16	\$1,926.00	<a href="http://www.chromtech.net.au/pdf2/vici09_157.pdf">http://www.chromtech.net.au/pdf2/vici09_157.pdf</a>
C74H-1694-.02	1/16inx0.25m Internal Sample 20nl,w/Pneumatic actuator,10,000psi	Microbore UHPLC 1/16	\$1,008.00	<a href="http://www.chromtech.net.au/pdf2/vici09_157.pdf">http://www.chromtech.net.au/pdf2/vici09_157.pdf</a>



## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C74-16R-.01	1/16inx0.25m Internal Sample 10nl,Replace Rotor,15,000psi	Microbore UHPLC 1/16	\$122.00	<a href="http://www.chromtech.net.au/pdf2/vici09_157.pdf">http://www.chromtech.net.au/pdf2/vici09_157.pdf</a>
C74X-1694-.05D	1/16inx0.25m Internal Sample 50nl,w/Replace Valve,10,000psi	Microbore UHPLC 1/16	\$1,198.00	<a href="http://www.chromtech.net.au/pdf2/vici09_157.pdf">http://www.chromtech.net.au/pdf2/vici09_157.pdf</a>
C74H-1694-.05A	1/16inx0.25m Internal Sample 50nl,w/Pneumatic actuator,10,000psi	Microbore UHPLC 1/16	\$1,326.00	<a href="http://www.chromtech.net.au/pdf2/vici09_157.pdf">http://www.chromtech.net.au/pdf2/vici09_157.pdf</a>
C74X-1694-.01D	1/16inx0.25m Internal Sample 10nl,w/Replace Valve,10,000psi	Microbore UHPLC 1/16	\$1,198.00	<a href="http://www.chromtech.net.au/pdf2/vici09_157.pdf">http://www.chromtech.net.au/pdf2/vici09_157.pdf</a>
C74H-1694-.05	1/16inx0.25m Internal Sample 50nl,w/Pneumatic actuator,10,000psi	Microbore UHPLC 1/16	\$1,008.00	<a href="http://www.chromtech.net.au/pdf2/vici09_157.pdf">http://www.chromtech.net.au/pdf2/vici09_157.pdf</a>
C74X-1694-.02EH	1/16inx0.25m Internal Sample 20nl,w/Microelectric actuator,10,000psi	Microbore UHPLC 1/16	\$2,142.00	<a href="http://www.chromtech.net.au/pdf2/vici09_157.pdf">http://www.chromtech.net.au/pdf2/vici09_157.pdf</a>
C74X-1694-.01EH	1/16inx0.25m Internal Sample 10nl,w/Microelectric actuator,10,000psi	Microbore UHPLC 1/16	\$2,142.00	<a href="http://www.chromtech.net.au/pdf2/vici09_157.pdf">http://www.chromtech.net.au/pdf2/vici09_157.pdf</a>
C74X-1694-.05E	1/16inx0.25m Internal Sample 50nl,w/Std Electric actuator,10,000psi	Microbore UHPLC 1/16	\$1,926.00	<a href="http://www.chromtech.net.au/pdf2/vici09_157.pdf">http://www.chromtech.net.au/pdf2/vici09_157.pdf</a>
C74X-1694-.02E	1/16inx0.25m Internal Sample 20nl,w/Std Electric actuator,10,000psi	Microbore UHPLC 1/16	\$1,926.00	<a href="http://www.chromtech.net.au/pdf2/vici09_157.pdf">http://www.chromtech.net.au/pdf2/vici09_157.pdf</a>
C74H-1694-.05D	1/16inx0.25m Internal Sample 50nl,w/Replace Valve,10,000psi	Microbore UHPLC 1/16	\$1,008.00	<a href="http://www.chromtech.net.au/pdf2/vici09_157.pdf">http://www.chromtech.net.au/pdf2/vici09_157.pdf</a>
C74-16R-.02	1/16inx0.25m Internal Sample 20nl,Replace Rotor,15,000psi	Microbore UHPLC 1/16	\$122.00	<a href="http://www.chromtech.net.au/pdf2/vici09_157.pdf">http://www.chromtech.net.au/pdf2/vici09_157.pdf</a>
C74-16R-.05	1/16inx0.25m Internal Sample 50nl,Replace Rotor,15,000psi	Microbore UHPLC 1/16	\$122.00	<a href="http://www.chromtech.net.au/pdf2/vici09_157.pdf">http://www.chromtech.net.au/pdf2/vici09_157.pdf</a>
C74-1C9	1/16inx0.25m Internal Sample 10nl,Replace Stator,15,000psi	Microbore UHPLC 1/16	\$816.00	<a href="http://www.chromtech.net.au/pdf2/vici09_157.pdf">http://www.chromtech.net.au/pdf2/vici09_157.pdf</a>
C74X-1694-.02D	1/16inx0.25m Internal Sample 20nl,w/Replace Valve,10,000psi	Microbore UHPLC 1/16	\$1,198.00	<a href="http://www.chromtech.net.au/pdf2/vici09_157.pdf">http://www.chromtech.net.au/pdf2/vici09_157.pdf</a>
C74X-1694-.02	1/16inx0.25m Internal Sample 20nl,w/Pneumatic actuator,10,000psi	Microbore UHPLC 1/16	\$1,198.00	<a href="http://www.chromtech.net.au/pdf2/vici09_157.pdf">http://www.chromtech.net.au/pdf2/vici09_157.pdf</a>
C74H-1694-.01	1/16inx0.25m Internal Sample 10nl,Manual,10,000psi	Microbore UHPLC 1/16	\$1,008.00	<a href="http://www.chromtech.net.au/pdf2/vici09_157.pdf">http://www.chromtech.net.au/pdf2/vici09_157.pdf</a>
C74X-1694-.01	1/16inx0.25m Internal Sample 10nl,w/Pneumatic actuator,10,000psi	Microbore UHPLC 1/16	\$1,198.00	<a href="http://www.chromtech.net.au/pdf2/vici09_157.pdf">http://www.chromtech.net.au/pdf2/vici09_157.pdf</a>
C74H-1694-.01E	1/16inx0.25m Internal Sample 10nl,w/Std Electric actuator,10,000psi	Microbore UHPLC 1/16	\$1,734.00	<a href="http://www.chromtech.net.au/pdf2/vici09_157.pdf">http://www.chromtech.net.au/pdf2/vici09_157.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C74X-1694-.05	1/16inx0.25m Internal Sample 50nl,w/Pneumatic actuator,10,000psi	Microbore UHPLC 1/16	\$1,198.00	<a href="http://www.chromtech.net.au/pdf2/vici09_157.pdf">http://www.chromtech.net.au/pdf2/vici09_157.pdf</a>
C74X-1694-.01A	1/16inx0.25m Internal Sample 10nl,w/Pneumatic actuator,10,000psi	Microbore UHPLC 1/16	\$1,518.00	<a href="http://www.chromtech.net.au/pdf2/vici09_157.pdf">http://www.chromtech.net.au/pdf2/vici09_157.pdf</a>
C74X-1694-.02A	1/16inx0.25m Internal Sample 20nl,w/Pneumatic actuator,10,000psi	Microbore UHPLC 1/16	\$1,518.00	<a href="http://www.chromtech.net.au/pdf2/vici09_157.pdf">http://www.chromtech.net.au/pdf2/vici09_157.pdf</a>
C74X-1694-.05A	1/16inx0.25m Internal Sample 50nl,w/Pneumatic actuator,10,000psi	Microbore UHPLC 1/16	\$1,518.00	<a href="http://www.chromtech.net.au/pdf2/vici09_157.pdf">http://www.chromtech.net.au/pdf2/vici09_157.pdf</a>
C74H-1694-.02D	1/16inx0.25m Internal Sample 20nl,w/Replace Valve,10,000psi	Microbore UHPLC 1/16	\$1,008.00	<a href="http://www.chromtech.net.au/pdf2/vici09_157.pdf">http://www.chromtech.net.au/pdf2/vici09_157.pdf</a>
C74H-1694-.01D	1/16inx0.25m Internal Sample 10nl,w/Replace Valve,10,000psi	Microbore UHPLC 1/16	\$1,008.00	<a href="http://www.chromtech.net.au/pdf2/vici09_157.pdf">http://www.chromtech.net.au/pdf2/vici09_157.pdf</a>
C74H-1694-.05EH	1/16inx0.25m Internal Sample 50nl,w/Microelectric actuator,10,000psi	Microbore UHPLC 1/16	\$1,950.00	<a href="http://www.chromtech.net.au/pdf2/vici09_157.pdf">http://www.chromtech.net.au/pdf2/vici09_157.pdf</a>
C74H-1694-.02EH	1/16inx0.25m Internal Sample 20nl,w/Microelectric actuator,10,000psi	Microbore UHPLC 1/16	\$1,950.00	<a href="http://www.chromtech.net.au/pdf2/vici09_157.pdf">http://www.chromtech.net.au/pdf2/vici09_157.pdf</a>
C74H-1694-.01EH	1/16inx0.25m Internal Sample 10nl,w/Microelectric actuator,10,000psi	Microbore UHPLC 1/16	\$1,950.00	<a href="http://www.chromtech.net.au/pdf2/vici09_157.pdf">http://www.chromtech.net.au/pdf2/vici09_157.pdf</a>
C74H-1694-.05E	1/16inx0.25m Internal Sample 50nl,w/Std Electric actuator,10,000psi	Microbore UHPLC 1/16	\$1,734.00	<a href="http://www.chromtech.net.au/pdf2/vici09_157.pdf">http://www.chromtech.net.au/pdf2/vici09_157.pdf</a>
C74H-1694-.02E	1/16inx0.25m Internal Sample 20nl,w/Std Electric actuator,10,000psi	Microbore UHPLC 1/16	\$1,734.00	<a href="http://www.chromtech.net.au/pdf2/vici09_157.pdf">http://www.chromtech.net.au/pdf2/vici09_157.pdf</a>
C72X-1690	10 Port 1/16mmx 0.25mm Manual,15,000psi	Microbore UHPLC 1/16	\$1,326.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
C72X-1698ED	8 Port 1/16mmx 0.25mm w/Microelectric actuator,10,000psi	Microbore UHPLC 1/16	\$2,283.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
C72X-1696EH	6 Port 1/16mmx 0.25mm w/Microelectric actuator,05,000psi	Microbore UHPLC 1/16	\$2,142.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
C72X-1694EH	4 Port 1/16mmx 0.25mm w/Microelectric actuator,15,000psi	Microbore UHPLC 1/16	\$2,142.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
C72X-1690E	10 Port 1/16mmx 0.25mm w/Std Electric actator,15,000psi	Microbore UHPLC 1/16	\$2,052.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
C72X-1698E	8 Port 1/16mmx 0.25mm w/Std Electric actator,10,000psi	Microbore UHPLC 1/16	\$1,989.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
C72X-1696E	6 Port 1/16mmx 0.25mm w/Std Electric actator,10,000psi	Microbore UHPLC 1/16	\$1,926.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
C72X-1694E	6 Port 1/16mmx 0.25mm w/Std Electric actator,15,000psi	Microbore UHPLC 1/16	\$1,926.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C72X-1690A	10 Port 1/16mmx 0.25mm Replace Valve,15,000psi	Microbore UHPLC 1/16	\$1,644.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
C72X-1698A	8 Port 1/16mmx 0.25mm w/Pneumatic actuator,10,000psi	Microbore UHPLC 1/16	\$1,581.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
C72X-1694A	6 Port 1/16mmx 0.25mm w/Pneumatic actuator,15,000psi	Microbore UHPLC 1/16	\$1,518.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
C72X-1696D	6 Port 1/16mmx 0.25mm Replace Valve,10,000psi	Microbore UHPLC 1/16	\$1,198.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
C72X-1698	8 Port 1/16mmx 0.25mm Manual,10,000psi	Microbore UHPLC 1/16	\$1,263.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
C72X-1696	6 Port 1/16mmx 0.25mm Manual,10,000psi	Microbore UHPLC 1/16	\$1,198.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
C72X-1694	6 Port 1/16mmx 0.25mm Manual,10,000psi	Microbore UHPLC 1/16	\$1,198.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
C72X-1696A	6 Port 1/16mmx 0.25mm w/Pneumatic actuator,10,000psi	Microbore UHPLC 1/16	\$1,518.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
CSL1K	1ml Stainless Steel Sample loops for C1 C2 and C3 inc fits	Microbore UHPLC 1/16	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
C72X-1690ED	10 Port 1/16mmx 0.25mm w/Microelectric actuator,15,000psi	Microbore UHPLC 1/16	\$2,346.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
C72H-1694D	4 Port 1/16mmx 0.25mm Replace Valve,15,000psi	Microbore UHPLC 1/16	\$1,008.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
C72H-1696D	6 Port 1/16mmx 0.25mm Replace Valve,15,000psi	Microbore UHPLC 1/16	\$1,008.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
C72H-1698D	8 Port 1/16mmx 0.25mm Replace Valve,15,000psi	Microbore UHPLC 1/16	\$1,071.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
C72H-1690D	10 Port 1/16mmx 0.25mm Replace Valve,15,000psi	Microbore UHPLC 1/16	\$1,134.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
C72-16R0	10 Port 1/16mmx 0.25mm Replace Rotor,15,000psi	Microbore UHPLC 1/16	\$122.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
C72-1C94	4 Port 1/16mmx 0.25mm Replace Stator,15,000psi	Microbore UHPLC 1/16	\$816.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
C72-1C98	8 Port 1/16mmx 0.25mm Replace Rotor,15,000psi	Microbore UHPLC 1/16	\$879.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
C72H-1698ED	8 Port 1/16mmx 0.25mm w/Microelectric actuator,15,000psi	Microbore UHPLC 1/16	\$2,091.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
CSL50	50ul Stainless Steel Sample loops for C1 C2 and C3 inc fits	Microbore UHPLC 1/16	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
C72H-1696EH	6 Port 1/16mmx 0.25mm w/Microelectric actuator,15,000psi	Microbore UHPLC 1/16	\$1,950.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
CSL5	5ul Stainless Steel Sample loops for C1 C2 and C3 inc fits	Microbore UHPLC 1/16	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
CSL100	100ul 1/16in Stainless steel loops inc fits	Microbore UHPLC 1/16	\$40.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
CSL2K	2ml 1/16in Stainless steel loops inc fits	Microbore UHPLC 1/16	\$57.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
CSL10	10ul Stainless Steel Sample loops for C1 C2 and C3 inc fits	Microbore UHPLC 1/16	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
CSL250	250ul Stainless Steel Sample loops for C1 C2 and C3 inc fits	Microbore UHPLC 1/16	\$40.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

CSL5K	5ml 1/16in Stainless steel loops inc fits	Microbore UHPLC 1/16	\$68.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
CSL20	20ul 1/16in Stainless steel loops inc fits	Microbore UHPLC 1/16	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
CSL500	500ul 1/16in Stainless steel loops inc fits	Microbore UHPLC 1/16	\$40.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
CSL10K	10ml 1/16in Stainless steel loops inc fits	Microbore UHPLC 1/16	\$87.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
CSL2	2ul Stainless Steel Sample loops for C1 C2 and C3 inc fits	Microbore UHPLC 1/16	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
C72H-1696	6 Port 1/16mmx 0.25mm Manual,15,000psi	Microbore UHPLC 1/16	\$1,008.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
C72H-1696A	6 Port 1/16mmx 0.25mm w/Pneumatic actuator,15,000psi	Microbore UHPLC 1/16	\$1,326.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
C72X-1698D	8 Port 1/16mmx 0.25mm Replace Valve,10,000psi	Microbore UHPLC 1/16	\$1,263.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
C72X-1690D	10 Port 1/16mmx 0.25mm Replace Valve,15,000psi	Microbore UHPLC 1/16	\$1,326.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
C72-16R4	4 Port 1/16mmx 0.25mm Replace Rotor,15,000psi	Microbore UHPLC 1/16	\$122.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
C72-16R6	6 Port 1/16mmx 0.25mm Replace Rotor,15,000psi	Microbore UHPLC 1/16	\$122.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
C72-16R8	8 Port 1/16mmx 0.25mm Replace Rotor,15,000psi	Microbore UHPLC 1/16	\$122.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
	10 Port 1/16mmx 0.25mm Replace Rotor,15,000psi	Microbore UHPLC 1/16	\$122.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
C72-1C96	6 Port 1/16mmx 0.25mm Replace Rotor,15,000psi	Microbore UHPLC 1/16	\$816.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
C72H-1690ED	10 Port 1/16mmx 0.25mm w/Microelectric actuator,15,000psi	Microbore UHPLC 1/16	\$2,154.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
C72H-1694	4 Port 1/16mmx 0.25mm Manual,15,000psi	Microbore UHPLC 1/16	\$1,008.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
C72X-1694D	6 Port 1/16mmx 0.25mm Replace Valve,15,000psi	Microbore UHPLC 1/16	\$1,198.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
C72H-1698	8 Port 1/16mmx 0.25mm Manual,15,000psi	Microbore UHPLC 1/16	\$1,071.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
C72H-1690	10 Port 1/16mmx 0.25mm Manual,15,000psi	Microbore UHPLC 1/16	\$1,134.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
C72H-1694A	4 Port 1/16mmx 0.25mm w/Pneumatic actuator,15,000psi	Microbore UHPLC 1/16	\$1,326.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
C72H-1698A	8 Port 1/16mmx 0.25mm w/Pneumatic actuator,15,000psi	Microbore UHPLC 1/16	\$1,389.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
C72H-1694E	4 Port 1/16mmx 0.25mm w/Std Electric actator,15,000psi	Microbore UHPLC 1/16	\$1,734.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
C72H-1696E	6 Port 1/16mmx 0.25mm w/Std Electric actator,15,000psi	Microbore UHPLC 1/16	\$1,734.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
C72H-1698E	8 Port 1/16mmx 0.25mm w/Std Electric actator,15,000psi	Microbore UHPLC 1/16	\$1,926.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
C72H-1690E	10 Port 1/16mmx 0.25mm w/Std Electric actator,15,000psi	Microbore UHPLC 1/16	\$1,989.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
C72H-1694EH	4 Port 1/16mmx 0.25mm w/Microelectric actuator,15,000psi	Microbore UHPLC 1/16	\$1,950.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C72-1C90	10 Port 1/16mmx 0.25mm Replace Stator,15,000psi	Microbore UHPLC 1/16	\$944.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
C72H-1690A	10 Port 1/16mmx 0.25mm w/Pneumatic actuator,15,000psi	Microbore UHPLC 1/16	\$1,454.00	<a href="http://www.chromtech.net.au/pdf2/vici09_156.pdf">http://www.chromtech.net.au/pdf2/vici09_156.pdf</a>
CN4-4344-.004D	1/32 100um ports Internal Sample 4nl, Replace Valve, 5,000psi	Nanovolume Injector 1/32	\$1,346.00	<a href="http://www.chromtech.net.au/pdf2/vici09_155.pdf">http://www.chromtech.net.au/pdf2/vici09_155.pdf</a>
CN4-4344-.01D	1/32 100um ports Internal Sample 10nl, Replace Valve, 5,000psi	Nanovolume Injector 1/32	\$1,346.00	<a href="http://www.chromtech.net.au/pdf2/vici09_155.pdf">http://www.chromtech.net.au/pdf2/vici09_155.pdf</a>
CN4-4C4I	1/32 100um ports Internal Sample 4nl, Replace Stator, 5,000psi	Nanovolume Injector 1/32	\$759.00	<a href="http://www.chromtech.net.au/pdf2/vici09_155.pdf">http://www.chromtech.net.au/pdf2/vici09_155.pdf</a>
CN4-43R-.02	1/32 100um ports Internal Sample 2nl, Replace Rotor, 5,000psi	Nanovolume Injector 1/32	\$117.00	<a href="http://www.chromtech.net.au/pdf2/vici09_155.pdf">http://www.chromtech.net.au/pdf2/vici09_155.pdf</a>
CN4-43R-.01	1/32 100um ports Internal Sample 10nl, Replace Rotor, 5,000psi	Nanovolume Injector 1/32	\$117.00	<a href="http://www.chromtech.net.au/pdf2/vici09_155.pdf">http://www.chromtech.net.au/pdf2/vici09_155.pdf</a>
CN4-43R-.004	1/32 100um ports Internal Sample 4nl, Replace Rotor, 5,000psi	Nanovolume Injector 1/32	\$117.00	<a href="http://www.chromtech.net.au/pdf2/vici09_155.pdf">http://www.chromtech.net.au/pdf2/vici09_155.pdf</a>
CN4-4344-.02EH	1/32 100um ports Internal Sample 20nl, w/Microelectric actuator, 5,000psi	Nanovolume Injector 1/32	\$2,289.00	<a href="http://www.chromtech.net.au/pdf2/vici09_155.pdf">http://www.chromtech.net.au/pdf2/vici09_155.pdf</a>
CN4-4344-.01EH	1/32 100um ports Internal Sample 10nl, w/Microelectric actuator, 5,000psi	Nanovolume Injector 1/32	\$2,289.00	<a href="http://www.chromtech.net.au/pdf2/vici09_155.pdf">http://www.chromtech.net.au/pdf2/vici09_155.pdf</a>
CN4-4344-.004EH	1/32 100um ports Internal Sample 4nl, w/Microelectric actuator, 5,000psi	Nanovolume Injector 1/32	\$2,289.00	<a href="http://www.chromtech.net.au/pdf2/vici09_155.pdf">http://www.chromtech.net.au/pdf2/vici09_155.pdf</a>
CN4-4344-.02	1/32 100um ports Internal Sample 20nl, Manual,5,000psi	Nanovolume Injector 1/32	\$1,346.00	<a href="http://www.chromtech.net.au/pdf2/vici09_155.pdf">http://www.chromtech.net.au/pdf2/vici09_155.pdf</a>
CN4-4344-.01	1/32 100um ports Internal Sample 10nl, Manual,5,000psi	Nanovolume Injector 1/32	\$1,346.00	<a href="http://www.chromtech.net.au/pdf2/vici09_155.pdf">http://www.chromtech.net.au/pdf2/vici09_155.pdf</a>
CN4-4344-.004	1/32 100um ports Internal Sample 4nl, Manual,5,000psi	Nanovolume Injector 1/32	\$1,346.00	<a href="http://www.chromtech.net.au/pdf2/vici09_155.pdf">http://www.chromtech.net.au/pdf2/vici09_155.pdf</a>
CN4-4344-.02D	1/32 100um ports Internal Sample 20nl, Replace Valve, 5,000psi	Nanovolume Injector 1/32	\$1,346.00	<a href="http://www.chromtech.net.au/pdf2/vici09_155.pdf">http://www.chromtech.net.au/pdf2/vici09_155.pdf</a>
CN2-4346	10 Port 1/32in 100um bore Manual	Nanovolume Injector 1/32	\$1,346.00	<a href="http://www.chromtech.net.au/pdf2/vici09_154.pdf">http://www.chromtech.net.au/pdf2/vici09_154.pdf</a>
CNSL250PK	250nl Sample Loop PEEK, for CN2 valve	Nanovolume Injector 1/32	\$48.00	<a href="http://www.chromtech.net.au/pdf2/vici09_154.pdf">http://www.chromtech.net.au/pdf2/vici09_154.pdf</a>
CNSL5KPK	5ul Sample Loop PEEK, for CN2 valve	Nanovolume Injector 1/32	\$48.00	<a href="http://www.chromtech.net.au/pdf2/vici09_154.pdf">http://www.chromtech.net.au/pdf2/vici09_154.pdf</a>
CNSL2KPK	2ul Sample Loop PEEK, for CN2 valve	Nanovolume Injector 1/32	\$48.00	<a href="http://www.chromtech.net.au/pdf2/vici09_154.pdf">http://www.chromtech.net.au/pdf2/vici09_154.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

CNSL500PK	500nl Sample Loop PEEK, for CN2 valve	Nanovolume Injector 1/32	\$48.00	<a href="http://www.chromtech.net.au/pdf2/vici09_154.pdf">http://www.chromtech.net.au/pdf2/vici09_154.pdf</a>
CN2-4346EH	10 Port 1/32in 100um bore Replace Valve	Nanovolume Injector 1/32	\$2,289.00	<a href="http://www.chromtech.net.au/pdf2/vici09_154.pdf">http://www.chromtech.net.au/pdf2/vici09_154.pdf</a>
CN2-4340EH	6 Port 1/32in 100um borew/Microelectric actuator	Nanovolume Injector 1/32	\$2,428.00	<a href="http://www.chromtech.net.au/pdf2/vici09_154.pdf">http://www.chromtech.net.au/pdf2/vici09_154.pdf</a>
CN2-4346D	10 Port 1/32in 100um borew/Microelectric actuator	Nanovolume Injector 1/32	\$1,346.00	<a href="http://www.chromtech.net.au/pdf2/vici09_154.pdf">http://www.chromtech.net.au/pdf2/vici09_154.pdf</a>
CN2-4340D	6 Port 1/32in 100um bore Replace Valve	Nanovolume Injector 1/32	\$1,485.00	<a href="http://www.chromtech.net.au/pdf2/vici09_154.pdf">http://www.chromtech.net.au/pdf2/vici09_154.pdf</a>
CN2-43R6	6 Port 1/32in 100um bore Replace Rotor	Nanovolume Injector 1/32	\$117.00	<a href="http://www.chromtech.net.au/pdf2/vici09_154.pdf">http://www.chromtech.net.au/pdf2/vici09_154.pdf</a>
CN2-43R0	10 Port 1/32in 100um bore Replace Rotor	Nanovolume Injector 1/32	\$117.00	<a href="http://www.chromtech.net.au/pdf2/vici09_154.pdf">http://www.chromtech.net.au/pdf2/vici09_154.pdf</a>
CN2-4C46I	6 Port 1/32in 100um bore Replace Stator	Nanovolume Injector 1/32	\$759.00	<a href="http://www.chromtech.net.au/pdf2/vici09_154.pdf">http://www.chromtech.net.au/pdf2/vici09_154.pdf</a>
CN2-4C40I	10 Port 1/32in 100um bore Replace Stator	Nanovolume Injector 1/32	\$867.00	<a href="http://www.chromtech.net.au/pdf2/vici09_154.pdf">http://www.chromtech.net.au/pdf2/vici09_154.pdf</a>
CNSL1KPK	1ul Sample Loop PEEK, for CN2 valve	Nanovolume Injector 1/32	\$48.00	<a href="http://www.chromtech.net.au/pdf2/vici09_154.pdf">http://www.chromtech.net.au/pdf2/vici09_154.pdf</a>
CN2-4340	6 Port 1/32in 100um bore Manual	UHPLC NanoInjectors 1/32	\$1,485.00	<a href="http://www.chromtech.net.au/pdf2/vici09_154.pdf">http://www.chromtech.net.au/pdf2/vici09_154.pdf</a>
C74NX-6694-.004E	1/32in 150um bore Internal Sample 4nl w/Std Electric actuator, 15,000psi	UHPLC NanoInjectors 1/32	\$2,181.00	<a href="http://www.chromtech.net.au/pdf2/vici09_153.pdf">http://www.chromtech.net.au/pdf2/vici09_153.pdf</a>
C74NX-6694-.02EH	1/32in 150um bore Internal Sample 20nl w/Microelectric actuator, 15,000psi	UHPLC NanoInjectors 1/32	\$2,397.00	<a href="http://www.chromtech.net.au/pdf2/vici09_153.pdf">http://www.chromtech.net.au/pdf2/vici09_153.pdf</a>
C74NX-6694-.004D	1/32in 150um bore Internal Sample 4nl Replace Valve, 15,000psi	UHPLC NanoInjectors 1/32	\$1,454.00	<a href="http://www.chromtech.net.au/pdf2/vici09_153.pdf">http://www.chromtech.net.au/pdf2/vici09_153.pdf</a>
C74NX-6694-.01D	1/32in 150um bore Internal Sample 10nl Replace Valve, 15,000psi	UHPLC NanoInjectors 1/32	\$1,454.00	<a href="http://www.chromtech.net.au/pdf2/vici09_153.pdf">http://www.chromtech.net.au/pdf2/vici09_153.pdf</a>
C74NX-6694-.02D	1/32in 150um bore Internal Sample 20nl Replace Valve, 15,000psi	UHPLC NanoInjectors 1/32	\$1,454.00	<a href="http://www.chromtech.net.au/pdf2/vici09_153.pdf">http://www.chromtech.net.au/pdf2/vici09_153.pdf</a>
C74N-66R-.004	6 Port 1/32in 150um bore Internal Sample 4nl Manual, 15,000psi	UHPLC NanoInjectors 1/32	\$122.00	<a href="http://www.chromtech.net.au/pdf2/vici09_153.pdf">http://www.chromtech.net.au/pdf2/vici09_153.pdf</a>
C74N-66R-.01	6 Port 1/32in 150um bore Internal Sample 10nl Manual, 15,000psi	UHPLC NanoInjectors 1/32	\$122.00	<a href="http://www.chromtech.net.au/pdf2/vici09_153.pdf">http://www.chromtech.net.au/pdf2/vici09_153.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C74N-6C9	6 Port 1/32in 150um bore Internal Sample Replace Stator, 15,000psi	UHPLC NanoInjectors 1/32	\$1,071.00	<a href="http://www.chromtech.net.au/pdf2/vici09_153.pdf">http://www.chromtech.net.au/pdf2/vici09_153.pdf</a>
C72NX-6690D	10 Port 1/32in 150um bore Replace Valve, 15,000psi	UHPLC NanoInjectors 1/32	\$1,581.00	<a href="http://www.chromtech.net.au/pdf2/vici09_153.pdf">http://www.chromtech.net.au/pdf2/vici09_153.pdf</a>
C74NX-6694-.01EH	1/32in 150um bore Internal Sample 10nl w/Microelectric actuator, 15,000psi	UHPLC NanoInjectors 1/32	\$2,397.00	<a href="http://www.chromtech.net.au/pdf2/vici09_153.pdf">http://www.chromtech.net.au/pdf2/vici09_153.pdf</a>
C74N-66R-.02	6 Port 1/32in 150um bore Internal Sample 20nl Manual, 15,000psi	UHPLC NanoInjectors 1/32	\$122.00	<a href="http://www.chromtech.net.au/pdf2/vici09_153.pdf">http://www.chromtech.net.au/pdf2/vici09_153.pdf</a>
CSLN1K	1ul Sample Loop for C72NX Valve	UHPLC NanoInjectors 1/32	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_153.pdf">http://www.chromtech.net.au/pdf2/vici09_153.pdf</a>
C72NX-6696EH	6 Port 1/32in 150um bore w/Microelectric actuator, 15,000psi	UHPLC NanoInjectors 1/32	\$2,397.00	<a href="http://www.chromtech.net.au/pdf2/vici09_153.pdf">http://www.chromtech.net.au/pdf2/vici09_153.pdf</a>
C72NX-6696	6 Port 1/32in 150um bore Replace Rotor, 15,000psi	UHPLC NanoInjectors 1/32	\$1,454.00	<a href="http://www.chromtech.net.au/pdf2/vici09_153.pdf">http://www.chromtech.net.au/pdf2/vici09_153.pdf</a>
C72NX-6690	10 Port 1/32in 150um bore Replace Rotor, 15,000psi	UHPLC NanoInjectors 1/32	\$1,581.00	<a href="http://www.chromtech.net.au/pdf2/vici09_153.pdf">http://www.chromtech.net.au/pdf2/vici09_153.pdf</a>
C72NX-6696A	6 Port 1/32in 150um bore w/Pneumatic actuator, 15,000psi	UHPLC NanoInjectors 1/32	\$1,773.00	<a href="http://www.chromtech.net.au/pdf2/vici09_153.pdf">http://www.chromtech.net.au/pdf2/vici09_153.pdf</a>
C72NX-6690A	10 Port 1/32in 150um bore w/Pneumatic actuator, 15,000psi	UHPLC NanoInjectors 1/32	\$2,028.00	<a href="http://www.chromtech.net.au/pdf2/vici09_153.pdf">http://www.chromtech.net.au/pdf2/vici09_153.pdf</a>
C72NX-6696E	6 Port 1/32in 150um bore w/Std Electric actuator, 15,000psi	UHPLC NanoInjectors 1/32	\$2,181.00	<a href="http://www.chromtech.net.au/pdf2/vici09_153.pdf">http://www.chromtech.net.au/pdf2/vici09_153.pdf</a>
C72NX-6690E	10 Port 1/32in 150um bore w/Std Electric actuator, 15,000psi	UHPLC NanoInjectors 1/32	\$2,307.00	<a href="http://www.chromtech.net.au/pdf2/vici09_153.pdf">http://www.chromtech.net.au/pdf2/vici09_153.pdf</a>
C72NX-6690ED	10 Port 1/32in 150um bore w/Microelectric actuator, 15,000psi	UHPLC NanoInjectors 1/32	\$2,601.00	<a href="http://www.chromtech.net.au/pdf2/vici09_153.pdf">http://www.chromtech.net.au/pdf2/vici09_153.pdf</a>
C72N-66R6	6 Port 1/32in 150um bore Replace Rotor, 15,000psi	UHPLC NanoInjectors 1/32	\$122.00	<a href="http://www.chromtech.net.au/pdf2/vici09_153.pdf">http://www.chromtech.net.au/pdf2/vici09_153.pdf</a>
C72N-66R0	10 Port 1/32in 150um bore Replace Rotor, 15,000psi	UHPLC NanoInjectors 1/32	\$122.00	<a href="http://www.chromtech.net.au/pdf2/vici09_153.pdf">http://www.chromtech.net.au/pdf2/vici09_153.pdf</a>
C72N-6C90	10 Port 1/32in 150um bore Replace Rotor, 15,000psi	UHPLC NanoInjectors 1/32	\$1,198.00	<a href="http://www.chromtech.net.au/pdf2/vici09_153.pdf">http://www.chromtech.net.au/pdf2/vici09_153.pdf</a>
C74NX-6694-.004EH	1/32in 150um bore Internal Sample 4nl w/Microelectric actuator, 15,000psi	UHPLC NanoInjectors 1/32	\$2,397.00	<a href="http://www.chromtech.net.au/pdf2/vici09_153.pdf">http://www.chromtech.net.au/pdf2/vici09_153.pdf</a>
CSLN2K	2ul Sample Loop for C72NX Valve	UHPLC NanoInjectors 1/32	\$57.00	<a href="http://www.chromtech.net.au/pdf2/vici09_153.pdf">http://www.chromtech.net.au/pdf2/vici09_153.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

CSLN5K	5ul Sample Loop for C72NX Valve	UHPLC NanoInjectors 1/32	\$68.00	<a href="http://www.chromtech.net.au/pdf2/vici09_153.pdf">http://www.chromtech.net.au/pdf2/vici09_153.pdf</a>
C74NX-6694-.004	1/32in 150um bore Internal Sample 4nl Replace Rotor, 15,000psi	UHPLC NanoInjectors 1/32	\$1,454.00	<a href="http://www.chromtech.net.au/pdf2/vici09_153.pdf">http://www.chromtech.net.au/pdf2/vici09_153.pdf</a>
C74NX-6694-.01	1/32in 150um bore Internal Sample 10nl Manual, 15,000psi	UHPLC NanoInjectors 1/32	\$1,454.00	<a href="http://www.chromtech.net.au/pdf2/vici09_153.pdf">http://www.chromtech.net.au/pdf2/vici09_153.pdf</a>
C74NX-6694-.02	1/32in 150um bore Internal Sample 20nl w/Microelectric actuator, 15,000psi	UHPLC NanoInjectors 1/32	\$1,454.00	<a href="http://www.chromtech.net.au/pdf2/vici09_153.pdf">http://www.chromtech.net.au/pdf2/vici09_153.pdf</a>
C74NX-6694-.004A	1/32in 150um bore Internal Sample 4nl w/Pneumatic actuator, 15,000psi	UHPLC NanoInjectors 1/32	\$1,773.00	<a href="http://www.chromtech.net.au/pdf2/vici09_153.pdf">http://www.chromtech.net.au/pdf2/vici09_153.pdf</a>
C74NX-6694-.01A	1/32in 150um bore Internal Sample 10l w/Pneumatic actuator, 15,000psi	UHPLC NanoInjectors 1/32	\$1,773.00	<a href="http://www.chromtech.net.au/pdf2/vici09_153.pdf">http://www.chromtech.net.au/pdf2/vici09_153.pdf</a>
C74NX-6694-.02A	1/32in 150um bore Internal Sample 20l w/Pneumatic actuator, 15,000psi	UHPLC NanoInjectors 1/32	\$1,773.00	<a href="http://www.chromtech.net.au/pdf2/vici09_153.pdf">http://www.chromtech.net.au/pdf2/vici09_153.pdf</a>
C74NX-6694-.01E	1/32in 150um bore Internal Sample 10nl w/Std Electric actuator, 15,000psi	UHPLC NanoInjectors 1/32	\$2,181.00	<a href="http://www.chromtech.net.au/pdf2/vici09_153.pdf">http://www.chromtech.net.au/pdf2/vici09_153.pdf</a>
C74NX-6694-.02E	1/32in 150um bore Internal Sample 20nl w/Std Electric actuator, 15,000psi	UHPLC NanoInjectors 1/32	\$2,181.00	<a href="http://www.chromtech.net.au/pdf2/vici09_153.pdf">http://www.chromtech.net.au/pdf2/vici09_153.pdf</a>
C72N-6C96	6 Port 1/32in 150um bore Replace Rotor, 15,000psi	UHPLC NanoInjectors 1/32	\$1,071.00	<a href="http://www.chromtech.net.au/pdf2/vici09_153.pdf">http://www.chromtech.net.au/pdf2/vici09_153.pdf</a>
C72NX-6696D	6 Port 1/32in 150um bore Replace Valve, 15,000psi	UHPLC NanoInjectors 1/32	\$1,454.00	<a href="http://www.chromtech.net.au/pdf2/vici09_153.pdf">http://www.chromtech.net.au/pdf2/vici09_153.pdf</a>
C72M-66R6	6 Port 360um 150um bore Replace Rotor,15.000psi	UHPLC NanoInjectors 360um	\$122.00	<a href="http://www.chromtech.net.au/pdf2/vici09_152.pdf">http://www.chromtech.net.au/pdf2/vici09_152.pdf</a>
C72M-66R0	10 Port 360um 150um bore Replace Rotor,15.000psi	UHPLC NanoInjectors 360um	\$122.00	<a href="http://www.chromtech.net.au/pdf2/vici09_152.pdf">http://www.chromtech.net.au/pdf2/vici09_152.pdf</a>
C72MU-4694A	4 Port 360um 100um bore w/Pneumatic actuator, 20,000psi	UHPLC NanoInjectors 360um	\$2,283.00	<a href="http://www.chromtech.net.au/pdf2/vici09_152.pdf">http://www.chromtech.net.au/pdf2/vici09_152.pdf</a>
C72M-6C96	6 Port 360um 150um bore Replace Stator,15.000psi	UHPLC NanoInjectors 360um	\$1,198.00	<a href="http://www.chromtech.net.au/pdf2/vici09_152.pdf">http://www.chromtech.net.au/pdf2/vici09_152.pdf</a>
C72MU-4696A	6 Port 360um 100um bore w/Pneumatic actuator, 20,000psi	UHPLC NanoInjectors 360um	\$2,283.00	<a href="http://www.chromtech.net.au/pdf2/vici09_152.pdf">http://www.chromtech.net.au/pdf2/vici09_152.pdf</a>
C72M-6C90	10 Port 360um 150um bore Replace Stator,15.000psi	UHPLC NanoInjectors 360um	\$1,326.00	<a href="http://www.chromtech.net.au/pdf2/vici09_152.pdf">http://www.chromtech.net.au/pdf2/vici09_152.pdf</a>
C72MU-4694E	4 Port 360um 100um bore w/Std Electric actuator, 20,000psi	UHPLC NanoInjectors 360um	\$2,691.00	<a href="http://www.chromtech.net.au/pdf2/vici09_152.pdf">http://www.chromtech.net.au/pdf2/vici09_152.pdf</a>



## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C72MX-6690D	10 Port 360um 150um bore Replace Valve, 15,000psi	UHPLC NanoInjectors 360um	\$1,708.00	<a href="http://www.chromtech.net.au/pdf2/vici09_152.pdf">http://www.chromtech.net.au/pdf2/vici09_152.pdf</a>
C72MX-6696E	6 Port 360um 150um bore w/Std Electric actuator, 15,000psi	UHPLC NanoInjectors 360um	\$2,307.00	<a href="http://www.chromtech.net.au/pdf2/vici09_152.pdf">http://www.chromtech.net.au/pdf2/vici09_152.pdf</a>
C72MU-4696E	6 Port 360um 100um bore Replace Valve, 20,000psi	UHPLC NanoInjectors 360um	\$2,691.00	<a href="http://www.chromtech.net.au/pdf2/vici09_152.pdf">http://www.chromtech.net.au/pdf2/vici09_152.pdf</a>
C72MU-4694EH	4 Port 360um 100um bore w/Microelectric actuator, 20,000psi	UHPLC NanoInjectors 360um	\$2,907.00	<a href="http://www.chromtech.net.au/pdf2/vici09_152.pdf">http://www.chromtech.net.au/pdf2/vici09_152.pdf</a>
C72MU-4696EH	6 Port 360um 100um bore w/Microelectric actuator, 20,000psi	UHPLC NanoInjectors 360um	\$2,907.00	<a href="http://www.chromtech.net.au/pdf2/vici09_152.pdf">http://www.chromtech.net.au/pdf2/vici09_152.pdf</a>
C72MU-4696D	6 Port 360um 100um bore Replace Valve, 20,000psi	UHPLC NanoInjectors 360um	\$1,964.00	<a href="http://www.chromtech.net.au/pdf2/vici09_152.pdf">http://www.chromtech.net.au/pdf2/vici09_152.pdf</a>
C72MX-6696D	6 Port 360um 150um bore Replace Valve, 15,000psi	UHPLC NanoInjectors 360um	\$1,581.00	<a href="http://www.chromtech.net.au/pdf2/vici09_152.pdf">http://www.chromtech.net.au/pdf2/vici09_152.pdf</a>
C72MU-4696	6 Port 360um 100um bore Manual, 20,000psi	UHPLC NanoInjectors 360um	\$1,964.00	<a href="http://www.chromtech.net.au/pdf2/vici09_152.pdf">http://www.chromtech.net.au/pdf2/vici09_152.pdf</a>
C72MX-6696EH	6 Port 360um 150um bore w/Microelectric actuator, 15,000psi	UHPLC NanoInjectors 360um	\$2,524.00	<a href="http://www.chromtech.net.au/pdf2/vici09_152.pdf">http://www.chromtech.net.au/pdf2/vici09_152.pdf</a>
C72MU-4694D	4 Port 360um 100um bore Replace Valve, 20,000psi	UHPLC NanoInjectors 360um	\$1,964.00	<a href="http://www.chromtech.net.au/pdf2/vici09_152.pdf">http://www.chromtech.net.au/pdf2/vici09_152.pdf</a>
C72MX-6690A	10 Port 360um 150um bore w/Pneumatic actuator, 15,000psi	UHPLC NanoInjectors 360um	\$2,028.00	<a href="http://www.chromtech.net.au/pdf2/vici09_152.pdf">http://www.chromtech.net.au/pdf2/vici09_152.pdf</a>
C72MX-6696A	6 Port 360um 150um bore w/Pneumatic actuator, 15,000psi	UHPLC NanoInjectors 360um	\$2,028.00	<a href="http://www.chromtech.net.au/pdf2/vici09_152.pdf">http://www.chromtech.net.au/pdf2/vici09_152.pdf</a>
C72MX-6690	10 Port 360um 150um bore Manual, 15,000psi	UHPLC NanoInjectors 360um	\$1,708.00	<a href="http://www.chromtech.net.au/pdf2/vici09_152.pdf">http://www.chromtech.net.au/pdf2/vici09_152.pdf</a>
C72MX-6696	6 Port 360um 150um bore Manual, 15,000psi	UHPLC NanoInjectors 360um	\$1,581.00	<a href="http://www.chromtech.net.au/pdf2/vici09_152.pdf">http://www.chromtech.net.au/pdf2/vici09_152.pdf</a>
C72M-4C94	4 Port 360um 100um bore Replace Stator, 20,000psi	UHPLC NanoInjectors 360um	\$1,389.00	<a href="http://www.chromtech.net.au/pdf2/vici09_152.pdf">http://www.chromtech.net.au/pdf2/vici09_152.pdf</a>
C72M-46R6	6 Port 360um 100um bore Replace Rotor,20,000psi	UHPLC NanoInjectors 360um	\$122.00	<a href="http://www.chromtech.net.au/pdf2/vici09_152.pdf">http://www.chromtech.net.au/pdf2/vici09_152.pdf</a>
C72MU-4694	4 Port 360um 100um bore Manual, 20,000psi	UHPLC NanoInjectors 360um	\$1,964.00	<a href="http://www.chromtech.net.au/pdf2/vici09_152.pdf">http://www.chromtech.net.au/pdf2/vici09_152.pdf</a>
C72M-46R4	4 Port 360um 100um bore Replace Rotor, 20,000psi	UHPLC NanoInjectors 360um	\$122.00	<a href="http://www.chromtech.net.au/pdf2/vici09_152.pdf">http://www.chromtech.net.au/pdf2/vici09_152.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

C72MX-6690E	10 Port 360um 150um bore w/Std Electric actuator, 15,000psi	UHPLC NanoInjectors 360um	\$2,436.00	<a href="http://www.chromtech.net.au/pdf2/vici09_152.pdf">http://www.chromtech.net.au/pdf2/vici09_152.pdf</a>
C72M-4C96	6 Port 360um 100um bore Replace Stator, 20,000psi	UHPLC NanoInjectors 360um	\$1,389.00	<a href="http://www.chromtech.net.au/pdf2/vici09_152.pdf">http://www.chromtech.net.au/pdf2/vici09_152.pdf</a>
C72MX-6690ED	10 Port 360um 150um bore w/Microelectric actuator, 15,000psi	UHPLC NanoInjectors 360um	\$2,728.00	<a href="http://www.chromtech.net.au/pdf2/vici09_152.pdf">http://www.chromtech.net.au/pdf2/vici09_152.pdf</a>
EZRF42PK-5	1/4"to1/8" PEEK/Glass,5pk	Valco Fittings	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
EZRF42V-5	1/4"to1/8" PolyI/Valcon,5pk	Valco Fittings	\$66.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
EZRF1.5KF-5	1/16"to1/32" CTFE,5pk	Valco Fittings	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
EZRF2.5KF-5	1/8"to1/32" CTFE,5pk	Valco Fittings	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
EZRF21KF-5	1/8"to1/16" CTFE,5pk	Valco Fittings	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
EZRF41KF-5	1/4"to1/16" CTFE,5pk	Valco Fittings	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
EZRF42KF-5	1/4"to1/8" CTFE,5pk	Valco Fittings	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
RF1.5TFG-5	1/16"to1/32" PTFE/Glass,5pk	Valco Fittings	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
RF1.5PK-5	1/16"to1/32" PEEK/Glass,5pk	Valco Fittings	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
RF1.5V-5	1/16"to1/32" PolyI/Valcon,5pk	Valco Fittings	\$38.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
RF21KF-5	1/8"to1/16" CTFE,5pk	Valco Fittings	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
RF2.5V-5	1/8"to1/32" PolyI/Valcon,5pk	Valco Fittings	\$38.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
ZRF2.5V1-5	Virgin Polyimide 1/8in to 1/32in Polymeric Internal	Valco Fittings	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
RF21PK-5	1/8"to1/16" PEEK/Glass,5pk	Valco Fittings	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
RF21V-5	1/8"to1/16" PolyI/Valcon,5pk	Valco Fittings	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
RF41TFG-5	1/4"to1/16" PTFE/Glass,5pk	Valco Fittings	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
RF41PK-5	1/4"to1/16" PEEK/Glass,5pk	Valco Fittings	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
RF41V-5	1/4"to1/16" PolyI/Valcon,5pk	Valco Fittings	\$51.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
RF42TFG-5	1/4"to1/8" PTFE/Glass,5pk	Valco Fittings	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
RF42PK-5	1/4"to1/8" PEEK/Glass,5pk	Valco Fittings	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
RF42V-5	1/4"to1/8" PolyI/Valcon,5pk	Valco Fittings	\$51.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
RF1.5KF-5	1/16"to1/32" CTFE,5pk	Valco Fittings	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
RF2.5KF-5	1/8"to1/32" CTFE,5pk	Valco Fittings	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
RF2.5TFG-5	1/8"to1/32" PTFE/Glass,5pk	Valco Fittings	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
ZRF42V1-5	Virgin Polyimide 1/4in to 1/8in Polymeric Internal	Valco Fittings	\$114.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
ZRF1.5V-5	Valcon Polyimide 1/16in to 1/32in Polymeric Internal	Valco Fittings	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
ZRF2.5TFG-5	1/8"to1/32" PTFEGlass,5pk	Valco Fittings	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
ZRF2.5PK-5	1/8"to1/32" PEEK,5pk	Valco Fittings	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

ZRF2.5V-5	Valcon Polyimide 1/8in to 1/32in Polymeric Internal	Valco Fittings	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
ZRF21TFG-5	Glass-filled PTFE 1/8in to 1/16in Polymeric Internal	Valco Fittings	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
ZRF21PK-5	PEEK 1/8in to 1/16in Polymeric Internal	Valco Fittings	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
ZRF21V-5	Valcon Polyimide 1/8in to 1/16in Polymeric Internal	Valco Fittings	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
ZRF21KF-5	CTFE 1/8in to 1/16in Polymeric Internal	Valco Fittings	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
ZRF21V1-5	Virgin Polyimide 1/8in to 1/16in Polymeric Internal	Valco Fittings	\$102.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
ZRF41KF-5	CTFE 1/4in to 1/16in Polymeric Internal	Valco Fittings	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
EZRF42TFG-5	1/4"to1/8" PTFE/Glass,5pk	Valco Fittings	\$38.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
ZRF42KF-5	CTFE 1/4in to 1/8in Polymeric Internal	Valco Fittings	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
EZRF41V-5	1/4"to1/16" Poly/Valcon,5pk	Valco Fittings	\$66.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
EZRF1.5TFG-5	1/16"to1/32" PTFE/Glass,5pk	Valco Fittings	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
EZRF1.5V-5	1/16"to1/32" Poly/Valcon,5pk	Valco Fittings	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
EZRF2.5TFG-5	1/8"to1/32" PTFE/Glass,5pk	Valco Fittings	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
EZRF2.5PK-5	1/8"to1/32" PEEK/Glass,5pk	Valco Fittings	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
EZRF2.5V-5	1/8"to1/32" Poly/Valcon,5pk	Valco Fittings	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
EZRF21TFG-5	1/8"to1/16" PTFE/Glass,5pk	Valco Fittings	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
EZRF21PK-5	1/8"to1/16" PEEK/Glass,5pk	Valco Fittings	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
EZRF21V-5	1/8"to1/16" Poly/Valcon,5pk	Valco Fittings	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
EZRF41TFG-5	1/4"to1/16" PTFE/Glass,5pk	Valco Fittings	\$38.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
EZRF41PK-5	1/4"to1/16"PEEK,5pk	Valco Fittings	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
RF21TFG-5	1/8"to1/16" PTFE/Glass,5pk	Valco Fittings	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
ZRF41V1-5	Virgin Polyimide 1/4in to 1/16in Polymeric Internal	Valco Fittings	\$114.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
ZRF1.5V1-5	Virgin Polyimide 1/16in to 1/32in Polymeric Internal	Valco Fittings	\$90.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
ZRF1.5TFG-5	Glass-filled PTFE 1/16in to 1/32in Polymeric Internal	Valco Fittings	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
RF41KF-5	1/4"to1/16" CTFE,5pk	Valco Fittings	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
ZRF1.5PK-5	PEEK 1/16in to 1/32in Polymeric Internal	Valco Fittings	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
RF2.5PK-5	1/8"to1/32" PEEK/Glass,5pk	Valco Fittings	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
ZRF1.5KF-5	CTFE 1/16in to 1/32in Polymeric Internal	Valco Fittings	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

ZRF42V-5	Valcon Polyimide 1/4in to 1/8in Polymeric Internal	Valco Fittings	\$66.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
ZRF41TFG-5	Glass-filled PTFE 1/4in to 1/16in Polymeric Internal	Valco Fittings	\$38.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
ZRF42TFG-5	Glass-filled PTFE 1/4in to 1/8in Polymeric Internal	Valco Fittings	\$38.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
ZRF41V-5	Valcon Polyimide 1/4in to 1/16in Polymeric Internal	Valco Fittings	\$66.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
ZRF41PK-5	PEEK 1/4in to 1/16in Polymeric Internal	Valco Fittings	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
ZRF42PK-5	PEEK 1/4in to 1/8in Polymeric Internal	Valco Fittings	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
RF42KF-5	1/4"to1/8" CTFE,5pk	Valco Fittings	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
ZRF2.5KF-5	CTFE 1/8in to 1/32in Polymeric Internal	Valco Fittings	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>
DV22-2116	6 port sampling/switching 1/16in fittings 0.40mm ports	Valco Valves	\$1,320.00	<a href="http://www.chromtech.net.au/pdf2/vici09_142.pdf">http://www.chromtech.net.au/pdf2/vici09_142.pdf</a>
DV22-31D	Diaphragm valves 1/16in 0.75mm ports, replace Diaphragm for 0.030in	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_142.pdf">http://www.chromtech.net.au/pdf2/vici09_142.pdf</a>
DV22-21D	Diaphragm valves 1/16in 0.75mm ports, replace Diaphragm for 0.010,0.016in	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_142.pdf">http://www.chromtech.net.au/pdf2/vici09_142.pdf</a>
CSL10K	10ml 1/16in Stainless steel loops inc fits	Valco Valves	\$87.00	<a href="http://www.chromtech.net.au/pdf2/vici09_142.pdf">http://www.chromtech.net.au/pdf2/vici09_142.pdf</a>
DV12-1116	6 port sampling/switching 1/32in 0.25mm ports	Valco Valves	\$1,530.00	<a href="http://www.chromtech.net.au/pdf2/vici09_142.pdf">http://www.chromtech.net.au/pdf2/vici09_142.pdf</a>
DV12-1110	10 port multifunctional 1/32in 0.25mm ports	Valco Valves	\$1,912.00	<a href="http://www.chromtech.net.au/pdf2/vici09_142.pdf">http://www.chromtech.net.au/pdf2/vici09_142.pdf</a>
DV22-3114-1	1ul internal sample 4 port 1/16in fittings 0.75mm ports	Valco Valves	\$1,320.00	<a href="http://www.chromtech.net.au/pdf2/vici09_142.pdf">http://www.chromtech.net.au/pdf2/vici09_142.pdf</a>
DV22-22D	Diaphragm valves 1/16in 0.75mm ports, replace Diaphragm for 0.030in	Valco Valves	\$74.00	<a href="http://www.chromtech.net.au/pdf2/vici09_142.pdf">http://www.chromtech.net.au/pdf2/vici09_142.pdf</a>
DV22-2114-1	1ul internal sample 4 port 1/16in fittings 0.40mm ports	Valco Valves	\$1,320.00	<a href="http://www.chromtech.net.au/pdf2/vici09_142.pdf">http://www.chromtech.net.au/pdf2/vici09_142.pdf</a>
CSL5	5ul Stainless Steel Sample loops for C1 C2 and C3 inc fits	Valco Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_142.pdf">http://www.chromtech.net.au/pdf2/vici09_142.pdf</a>
DV22-2110	10 port multifunctional 1/16in fittings 0.40mm ports	Valco Valves	\$1,530.00	<a href="http://www.chromtech.net.au/pdf2/vici09_142.pdf">http://www.chromtech.net.au/pdf2/vici09_142.pdf</a>
DV22-2114-.5	.5ul internal sample 4 port 1/16in fittings 0.40mm ports	Valco Valves	\$1,320.00	<a href="http://www.chromtech.net.au/pdf2/vici09_142.pdf">http://www.chromtech.net.au/pdf2/vici09_142.pdf</a>
DV22-3116	6 port sampling/switching 1/16in fittings 0.75mm ports	Valco Valves	\$1,320.00	<a href="http://www.chromtech.net.au/pdf2/vici09_142.pdf">http://www.chromtech.net.au/pdf2/vici09_142.pdf</a>
DV22-3110	10 port multifunctional 1/16in fittings 0.75mm ports	Valco Valves	\$1,530.00	<a href="http://www.chromtech.net.au/pdf2/vici09_142.pdf">http://www.chromtech.net.au/pdf2/vici09_142.pdf</a>
CSL100	100ul 1/16in Stainless steel loops inc fits	Valco Valves	\$40.00	<a href="http://www.chromtech.net.au/pdf2/vici09_142.pdf">http://www.chromtech.net.au/pdf2/vici09_142.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

CSL250	250ul Stainless Steel Sample loops for C1 C2 and C3 inc fits	Valco Valves	\$40.00	<a href="http://www.chromtech.net.au/pdf2/vici09_142.pdf">http://www.chromtech.net.au/pdf2/vici09_142.pdf</a>
CSL500	500ul 1/16in Stainless steel loops inc fits	Valco Valves	\$40.00	<a href="http://www.chromtech.net.au/pdf2/vici09_142.pdf">http://www.chromtech.net.au/pdf2/vici09_142.pdf</a>
CSL1K	1ml Stainless Steel Sample loops for C1 C2 and C3 inc fits	Valco Valves	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_142.pdf">http://www.chromtech.net.au/pdf2/vici09_142.pdf</a>
CSL2K	2ml 1/16in Stainless steel loops inc fits	Valco Valves	\$57.00	<a href="http://www.chromtech.net.au/pdf2/vici09_142.pdf">http://www.chromtech.net.au/pdf2/vici09_142.pdf</a>
CSL50	50ul Stainless Steel Sample loops for C1 C2 and C3 inc fits	Valco Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_142.pdf">http://www.chromtech.net.au/pdf2/vici09_142.pdf</a>
CSL5K	5ml 1/16in Stainless steel loops inc fits	Valco Valves	\$68.00	<a href="http://www.chromtech.net.au/pdf2/vici09_142.pdf">http://www.chromtech.net.au/pdf2/vici09_142.pdf</a>
DV12-1114-1	1ul internal sample 4 port 1/32in 0.25mm ports	Valco Valves	\$1,530.00	<a href="http://www.chromtech.net.au/pdf2/vici09_142.pdf">http://www.chromtech.net.au/pdf2/vici09_142.pdf</a>
DV12-1114-.5	.5ul internal sample 4 port 1/32in 0.25mm ports	Valco Valves	\$1,530.00	<a href="http://www.chromtech.net.au/pdf2/vici09_142.pdf">http://www.chromtech.net.au/pdf2/vici09_142.pdf</a>
CSL10	10ul Stainless Steel Sample loops for C1 C2 and C3 inc fits	Valco Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_142.pdf">http://www.chromtech.net.au/pdf2/vici09_142.pdf</a>
CSL2	2ul Stainless Steel Sample loops for C1 C2 and C3 inc fits	Valco Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_142.pdf">http://www.chromtech.net.au/pdf2/vici09_142.pdf</a>
CSL20	20ul 1/16in Stainless steel loops inc fits	Valco Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_142.pdf">http://www.chromtech.net.au/pdf2/vici09_142.pdf</a>
DV22-3114-.5	.5ul internal sample 4 port 1/16in fittings 0.75mm ports	Valco Valves	\$1,320.00	<a href="http://www.chromtech.net.au/pdf2/vici09_142.pdf">http://www.chromtech.net.au/pdf2/vici09_142.pdf</a>
SL500CSTUW	500ul 1/16in Stainless steel loops	Valco Valves	\$42.00	<a href="http://www.chromtech.net.au/pdf2/vici09_133.pdf">http://www.chromtech.net.au/pdf2/vici09_133.pdf</a>
SL15CSTUW	15ul 1/16in Stainless steel loops	Valco Valves	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_133.pdf">http://www.chromtech.net.au/pdf2/vici09_133.pdf</a>
SL250CSTUW	250ul 1/16in Stainless steel loops	Valco Valves	\$36.00	<a href="http://www.chromtech.net.au/pdf2/vici09_133.pdf">http://www.chromtech.net.au/pdf2/vici09_133.pdf</a>
SL10CSTUW	10ul 1/16in Stainless steel loops	Valco Valves	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_133.pdf">http://www.chromtech.net.au/pdf2/vici09_133.pdf</a>
SSACST6UW	Replacement rotor 6 Columns or Loops 1/16in fittings 0.4mm ports	Valco Valves	\$134.00	<a href="http://www.chromtech.net.au/pdf2/vici09_133.pdf">http://www.chromtech.net.au/pdf2/vici09_133.pdf</a>
SL20CSTUW	20ul 1/16in Stainless steel loops	Valco Valves	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_133.pdf">http://www.chromtech.net.au/pdf2/vici09_133.pdf</a>
EMTCST6UW	Microelectric actuator 6 Columns or Loops 1/16in fittings 0.4mm ports	Valco Valves	\$2,301.00	<a href="http://www.chromtech.net.au/pdf2/vici09_133.pdf">http://www.chromtech.net.au/pdf2/vici09_133.pdf</a>
SSACST4UW	Replacement rotor 4 Columns or Loops 1/16in fittings 0.4mm ports	Valco Valves	\$134.00	<a href="http://www.chromtech.net.au/pdf2/vici09_133.pdf">http://www.chromtech.net.au/pdf2/vici09_133.pdf</a>
DCST4UW	Replacement valve 4 Columns or Loops 1/16in fittings 0.4mm ports	Valco Valves	\$918.00	<a href="http://www.chromtech.net.au/pdf2/vici09_133.pdf">http://www.chromtech.net.au/pdf2/vici09_133.pdf</a>
ACST4UW	Air actuator 4 Columns or Loops 1/16in fittings 0.4mm ports	Valco Valves	\$1,434.00	<a href="http://www.chromtech.net.au/pdf2/vici09_133.pdf">http://www.chromtech.net.au/pdf2/vici09_133.pdf</a>
ACST6UW	Air actuator 6 Columns or Loops 1/16in fittings 0.4mm ports	Valco Valves	\$1,504.00	<a href="http://www.chromtech.net.au/pdf2/vici09_133.pdf">http://www.chromtech.net.au/pdf2/vici09_133.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

ECST4UW	Electric actuator 4 Columns or Loops 1/16in fittings 0.4mm ports	Valco Valves	\$1,950.00	<a href="http://www.chromtech.net.au/pdf2/vici09_133.pdf">http://www.chromtech.net.au/pdf2/vici09_133.pdf</a>
ECST6UW	Electric actuator 6 Columns or Loops 1/16in fittings 0.4mm ports	Valco Valves	\$2,020.00	<a href="http://www.chromtech.net.au/pdf2/vici09_133.pdf">http://www.chromtech.net.au/pdf2/vici09_133.pdf</a>
EMTCST4UW	Microelectric actuator 4 Columns or Loops 1/16in fittings 0.4mm ports	Valco Valves	\$2,232.00	<a href="http://www.chromtech.net.au/pdf2/vici09_133.pdf">http://www.chromtech.net.au/pdf2/vici09_133.pdf</a>
DCST6UW	Replacement valve 6 Columns or Loops 1/16in fittings 0.4mm ports	Valco Valves	\$988.00	<a href="http://www.chromtech.net.au/pdf2/vici09_133.pdf">http://www.chromtech.net.au/pdf2/vici09_133.pdf</a>
SL25CSTUW	25ul 1/16in Stainless steel loops	Valco Valves	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_133.pdf">http://www.chromtech.net.au/pdf2/vici09_133.pdf</a>
SL2KCSTUW	2ml 1/16in Stainless steel loops	Valco Valves	\$58.00	<a href="http://www.chromtech.net.au/pdf2/vici09_133.pdf">http://www.chromtech.net.au/pdf2/vici09_133.pdf</a>
SL50CSTUW	50ul 1/16in Stainless steel loops	Valco Valves	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_133.pdf">http://www.chromtech.net.au/pdf2/vici09_133.pdf</a>
SL5KCSTUW	5ml 1/16in Stainless steel loops	Valco Valves	\$69.00	<a href="http://www.chromtech.net.au/pdf2/vici09_133.pdf">http://www.chromtech.net.au/pdf2/vici09_133.pdf</a>
SL100CSTUW	100ul 1/16in Stainless steel loops	Valco Valves	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_133.pdf">http://www.chromtech.net.au/pdf2/vici09_133.pdf</a>
SL10KCSTUW	10ml 1/16in Stainless steel loops	Valco Valves	\$90.00	<a href="http://www.chromtech.net.au/pdf2/vici09_133.pdf">http://www.chromtech.net.au/pdf2/vici09_133.pdf</a>
SL1KCSTUW	1ml 1/16in Stainless steel loops	Valco Valves	\$46.00	<a href="http://www.chromtech.net.au/pdf2/vici09_133.pdf">http://www.chromtech.net.au/pdf2/vici09_133.pdf</a>
DSD8UW	Replacement valve 8 Position 1/8in fittings 0.75mm ports	Valco Valves	\$981.00	<a href="http://www.chromtech.net.au/pdf2/vici09_132.pdf">http://www.chromtech.net.au/pdf2/vici09_132.pdf</a>
SD6UW	Manual 6 Position 1/8in fittings 0.75mm ports	Valco Valves	\$1,077.00	<a href="http://www.chromtech.net.au/pdf2/vici09_132.pdf">http://www.chromtech.net.au/pdf2/vici09_132.pdf</a>
SD8UW	Manual 8 Position 1/8in fittings 0.75mm ports	Valco Valves	\$1,185.00	<a href="http://www.chromtech.net.au/pdf2/vici09_132.pdf">http://www.chromtech.net.au/pdf2/vici09_132.pdf</a>
ASD6UW	Air actuator 6 Position 1/8in fittings 0.75mm ports	Valco Valves	\$1,389.00	<a href="http://www.chromtech.net.au/pdf2/vici09_132.pdf">http://www.chromtech.net.au/pdf2/vici09_132.pdf</a>
ASD8UW	Air actuator 8 Position 1/8in fittings 0.75mm ports	Valco Valves	\$1,498.00	<a href="http://www.chromtech.net.au/pdf2/vici09_132.pdf">http://www.chromtech.net.au/pdf2/vici09_132.pdf</a>
ESD4UW	Electric actuator 4 Position 1/8in fittings 0.75mm ports	Valco Valves	\$1,964.00	<a href="http://www.chromtech.net.au/pdf2/vici09_132.pdf">http://www.chromtech.net.au/pdf2/vici09_132.pdf</a>
SSASD8UW	Replacement rotor 8 Position 1/8in fittings 0.75mm ports	Valco Valves	\$134.00	<a href="http://www.chromtech.net.au/pdf2/vici09_132.pdf">http://www.chromtech.net.au/pdf2/vici09_132.pdf</a>
SSASD4UW	Replacement rotor 4 Position 1/8in fittings 0.75mm ports	Valco Valves	\$134.00	<a href="http://www.chromtech.net.au/pdf2/vici09_132.pdf">http://www.chromtech.net.au/pdf2/vici09_132.pdf</a>
SSACSD6UW	Replacement rotor 6 Position 1/16in fittings 0.4mm ports	Valco Valves	\$134.00	<a href="http://www.chromtech.net.au/pdf2/vici09_132.pdf">http://www.chromtech.net.au/pdf2/vici09_132.pdf</a>
DSD6UW	Replacement valve 6 Position 1/8in fittings 0.75mm ports	Valco Valves	\$873.00	<a href="http://www.chromtech.net.au/pdf2/vici09_132.pdf">http://www.chromtech.net.au/pdf2/vici09_132.pdf</a>
DSD4UW	Replacement valve 4 Position 1/8in fittings 0.75mm ports	Valco Valves	\$804.00	<a href="http://www.chromtech.net.au/pdf2/vici09_132.pdf">http://www.chromtech.net.au/pdf2/vici09_132.pdf</a>
EMTSD8UW	Microelectric actuator 8 Position 1/8in fittings 0.75mm ports	Valco Valves	\$2,295.00	<a href="http://www.chromtech.net.au/pdf2/vici09_132.pdf">http://www.chromtech.net.au/pdf2/vici09_132.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

EMTSD6UW	Microelectric actuator 6 Position 1/8in fittings 0.75mm ports	Valco Valves	\$2,187.00	<a href="http://www.chromtech.net.au/pdf2/vici09_132.pdf">http://www.chromtech.net.au/pdf2/vici09_132.pdf</a>
EMTSD4UW	Microelectric actuator 4 Position 1/8in fittings 0.75mm ports	Valco Valves	\$2,116.00	<a href="http://www.chromtech.net.au/pdf2/vici09_132.pdf">http://www.chromtech.net.au/pdf2/vici09_132.pdf</a>
ESD8UW	Electric actuator 8 Position 1/8in fittings 0.75mm ports	Valco Valves	\$2,014.00	<a href="http://www.chromtech.net.au/pdf2/vici09_132.pdf">http://www.chromtech.net.au/pdf2/vici09_132.pdf</a>
SSASD6UW	Replacement rotor 6 Position 1/8in fittings 0.75mm ports	Valco Valves	\$134.00	<a href="http://www.chromtech.net.au/pdf2/vici09_132.pdf">http://www.chromtech.net.au/pdf2/vici09_132.pdf</a>
ECSD6UW	Electric actuator 6 Position 1/16in fittings 0.4mm ports	Valco Valves	\$2,034.00	<a href="http://www.chromtech.net.au/pdf2/vici09_132.pdf">http://www.chromtech.net.au/pdf2/vici09_132.pdf</a>
CSD6UW	Manual 6 Position 1/16in fittings 0.4mm ports	Valco Valves	\$1,077.00	<a href="http://www.chromtech.net.au/pdf2/vici09_132.pdf">http://www.chromtech.net.au/pdf2/vici09_132.pdf</a>
CSD10UW	Manual 10 Position 1/16in fittings 0.4mm ports	Valco Valves	\$1,185.00	<a href="http://www.chromtech.net.au/pdf2/vici09_132.pdf">http://www.chromtech.net.au/pdf2/vici09_132.pdf</a>
ACSD4UW	Air actuator 4 Position 1/16in fittings 0.4mm ports	Valco Valves	\$1,320.00	<a href="http://www.chromtech.net.au/pdf2/vici09_132.pdf">http://www.chromtech.net.au/pdf2/vici09_132.pdf</a>
ECSD10UW	Electric actuator 10 Position 1/16in fittings 0.4mm ports	Valco Valves	\$2,014.00	<a href="http://www.chromtech.net.au/pdf2/vici09_132.pdf">http://www.chromtech.net.au/pdf2/vici09_132.pdf</a>
ACSD6UW	Air actuator 6 Position 1/16in fittings 0.4mm ports	Valco Valves	\$1,389.00	<a href="http://www.chromtech.net.au/pdf2/vici09_132.pdf">http://www.chromtech.net.au/pdf2/vici09_132.pdf</a>
ESD6UW	Electric actuator 6 Position 1/8in fittings 0.75mm ports	Valco Valves	\$2,034.00	<a href="http://www.chromtech.net.au/pdf2/vici09_132.pdf">http://www.chromtech.net.au/pdf2/vici09_132.pdf</a>
SD4UW	Manual 4 Position 1/8in fittings 0.75mm ports	Valco Valves	\$1,008.00	<a href="http://www.chromtech.net.au/pdf2/vici09_132.pdf">http://www.chromtech.net.au/pdf2/vici09_132.pdf</a>
ECSD4UW	Electric actuator 4 Position 1/16in fittings 0.4mm ports	Valco Valves	\$1,964.00	<a href="http://www.chromtech.net.au/pdf2/vici09_132.pdf">http://www.chromtech.net.au/pdf2/vici09_132.pdf</a>
SSACSD10UW	Replacement rotor 10 Position 1/16in fittings 0.4mm ports	Valco Valves	\$134.00	<a href="http://www.chromtech.net.au/pdf2/vici09_132.pdf">http://www.chromtech.net.au/pdf2/vici09_132.pdf</a>
EMTCSD4UW	Microelectric actuator 4 Position 1/16in fittings 0.4mm ports	Valco Valves	\$2,116.00	<a href="http://www.chromtech.net.au/pdf2/vici09_132.pdf">http://www.chromtech.net.au/pdf2/vici09_132.pdf</a>
EMTCSD10UW	Microelectric actuator 10 Position 1/16in fittings 0.4mm ports	Valco Valves	\$2,295.00	<a href="http://www.chromtech.net.au/pdf2/vici09_132.pdf">http://www.chromtech.net.au/pdf2/vici09_132.pdf</a>
DCSD4UW	Replacement valve 4 Position 1/16in fittings 0.4mm ports	Valco Valves	\$804.00	<a href="http://www.chromtech.net.au/pdf2/vici09_132.pdf">http://www.chromtech.net.au/pdf2/vici09_132.pdf</a>
DCSD6UW	Replacement valve 6 Position 1/16in fittings 0.4mm ports	Valco Valves	\$873.00	<a href="http://www.chromtech.net.au/pdf2/vici09_132.pdf">http://www.chromtech.net.au/pdf2/vici09_132.pdf</a>
DCSD10UW	Replacement valve 10 Position 1/16in fittings 0.4mm ports	Valco Valves	\$981.00	<a href="http://www.chromtech.net.au/pdf2/vici09_132.pdf">http://www.chromtech.net.au/pdf2/vici09_132.pdf</a>
SSACSD4UW	Replacement rotor 4 Position 1/16in fittings 0.4mm ports	Valco Valves	\$134.00	<a href="http://www.chromtech.net.au/pdf2/vici09_132.pdf">http://www.chromtech.net.au/pdf2/vici09_132.pdf</a>
CSD4UW	Manual 4 Position 1/16in fittings 0.4mm ports	Valco Valves	\$1,008.00	<a href="http://www.chromtech.net.au/pdf2/vici09_132.pdf">http://www.chromtech.net.au/pdf2/vici09_132.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

ACSD10UW	Air actuator 10 Position 1/16in fittings 0.4mm ports	Valco Valves	\$1,498.00	<a href="http://www.chromtech.net.au/pdf2/vici09_132.pdf">http://www.chromtech.net.au/pdf2/vici09_132.pdf</a>
EMTCSD6UW	Microelectric actuator 6 Position 1/16in fittings 0.4mm ports	Valco Valves	\$2,187.00	<a href="http://www.chromtech.net.au/pdf2/vici09_132.pdf">http://www.chromtech.net.au/pdf2/vici09_132.pdf</a>
E2STF12MWE	12 Port Selector Trapping FloThru STF w/Std Electric actuator 1/8x1.0mm ports	Valco Valves	\$2,301.00	<a href="http://www.chromtech.net.au/pdf2/vici09_131.pdf">http://www.chromtech.net.au/pdf2/vici09_131.pdf</a>
SSASTF12MWE	Replacement rotor 12 Position 1/8in fittings 1.0mm ports 4in standoff	Valco Valves	\$147.00	<a href="http://www.chromtech.net.au/pdf2/vici09_131.pdf">http://www.chromtech.net.au/pdf2/vici09_131.pdf</a>
SSASTF10MWE	Replacement rotor 10 Position 1/8in fittings 1.0mm ports 4in standoff	Valco Valves	\$147.00	<a href="http://www.chromtech.net.au/pdf2/vici09_131.pdf">http://www.chromtech.net.au/pdf2/vici09_131.pdf</a>
SSASTF6MWE	Replacement rotor 6 Position 1/8in fittings 1.0mm ports 4in standoff	Valco Valves	\$134.00	<a href="http://www.chromtech.net.au/pdf2/vici09_131.pdf">http://www.chromtech.net.au/pdf2/vici09_131.pdf</a>
DSTF12MWE	Replacement valve 12 Position 1/8in fittings 1.0mm ports 4in standoff	Valco Valves	\$1,269.00	<a href="http://www.chromtech.net.au/pdf2/vici09_131.pdf">http://www.chromtech.net.au/pdf2/vici09_131.pdf</a>
DSTF6MWE	Replacement valve 6 Position 1/8in fittings 1.0mm ports 4in standoff	Valco Valves	\$1,059.00	<a href="http://www.chromtech.net.au/pdf2/vici09_131.pdf">http://www.chromtech.net.au/pdf2/vici09_131.pdf</a>
EMT2STF16MWE	10 Port Selector Trapping FloThru STF w/Microelectric actuator 1/8x1.0mm ports	Valco Valves	\$2,582.00	<a href="http://www.chromtech.net.au/pdf2/vici09_131.pdf">http://www.chromtech.net.au/pdf2/vici09_131.pdf</a>
EMT2STF12MWE	10 Port Selector Trapping FloThru STF w/Microelectric actuator 1/8x1.0mm ports	Valco Valves	\$2,582.00	<a href="http://www.chromtech.net.au/pdf2/vici09_131.pdf">http://www.chromtech.net.au/pdf2/vici09_131.pdf</a>
EMT2STF10MWE	10 Port Selector Trapping FloThru STF w/Microelectric actuator 1/8x1.0mm ports	Valco Valves	\$2,480.00	<a href="http://www.chromtech.net.au/pdf2/vici09_131.pdf">http://www.chromtech.net.au/pdf2/vici09_131.pdf</a>
SSASTF16MWE	Replacement rotor 16 Position 1/8in fittings 1.0mm ports 4in standoff	Valco Valves	\$147.00	<a href="http://www.chromtech.net.au/pdf2/vici09_131.pdf">http://www.chromtech.net.au/pdf2/vici09_131.pdf</a>
E2STF16MWE	16 Port Selector Trapping FloThru STF w/Std Electric actuator 1/8x1.0mm ports	Valco Valves	\$2,301.00	<a href="http://www.chromtech.net.au/pdf2/vici09_131.pdf">http://www.chromtech.net.au/pdf2/vici09_131.pdf</a>
DSTF16MWE	Replacement valve 16 Position 1/8in fittings 1.0mm 4in standoff ports	Valco Valves	\$1,269.00	<a href="http://www.chromtech.net.au/pdf2/vici09_131.pdf">http://www.chromtech.net.au/pdf2/vici09_131.pdf</a>
E2STF10MWE	10 Port Selector Trapping FloThru STF w/Std Electric actuator 1/8x1.0mm ports	Valco Valves	\$2,199.00	<a href="http://www.chromtech.net.au/pdf2/vici09_131.pdf">http://www.chromtech.net.au/pdf2/vici09_131.pdf</a>
E2STF6MWE	6 Port Selector Trapping FloThru STF w/Std Electric actuator 1/8x1.0mm ports	Valco Valves	\$2,091.00	<a href="http://www.chromtech.net.au/pdf2/vici09_131.pdf">http://www.chromtech.net.au/pdf2/vici09_131.pdf</a>
A2STF16MWE	10 Port Selector Trapping FloThru STF w/Air actuator 1/8x1.0mm ports	Valco Valves	\$1,912.00	<a href="http://www.chromtech.net.au/pdf2/vici09_131.pdf">http://www.chromtech.net.au/pdf2/vici09_131.pdf</a>
A2STF12MWE	10 Port Selector Trapping FloThru STF w/Air actuator 1/8x1.0mm ports	Valco Valves	\$1,912.00	<a href="http://www.chromtech.net.au/pdf2/vici09_131.pdf">http://www.chromtech.net.au/pdf2/vici09_131.pdf</a>



## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

A2STF10MWE	10 Port Selector Trapping FloThru STF w/Air actuator 1/8x1.0mm ports	Valco Valves	\$1,683.00	<a href="http://www.chromtech.net.au/pdf2/vici09_131.pdf">http://www.chromtech.net.au/pdf2/vici09_131.pdf</a>
A2STF6MWE	10 Port Selector Trapping FloThru STF w/Air actuator 1/8x1.0mm ports	Valco Valves	\$1,575.00	<a href="http://www.chromtech.net.au/pdf2/vici09_131.pdf">http://www.chromtech.net.au/pdf2/vici09_131.pdf</a>
2STF16MWE	16 Port Selector Trapping FloThru STF Manual 1/8x1.0mm ports	Valco Valves	\$1,473.00	<a href="http://www.chromtech.net.au/pdf2/vici09_131.pdf">http://www.chromtech.net.au/pdf2/vici09_131.pdf</a>
2STF12MWE	12 Port Selector Trapping FloThru STF Manual 1/8x1.0mm ports	Valco Valves	\$1,473.00	<a href="http://www.chromtech.net.au/pdf2/vici09_131.pdf">http://www.chromtech.net.au/pdf2/vici09_131.pdf</a>
2STF10MWE	10 Port Selector Trapping FloThru STF Manual 1/8x1.0mm ports	Valco Valves	\$1,371.00	<a href="http://www.chromtech.net.au/pdf2/vici09_131.pdf">http://www.chromtech.net.au/pdf2/vici09_131.pdf</a>
2STF6MWE	6 Port Selector Trapping FloThru STF Manual 1/8x1.0mm ports	Valco Valves	\$1,263.00	<a href="http://www.chromtech.net.au/pdf2/vici09_131.pdf">http://www.chromtech.net.au/pdf2/vici09_131.pdf</a>
EMT2STF6MWE	10 Port Selector Trapping FloThru STF w/Microelectric actuator 1/8x1.0mm ports	Valco Valves	\$2,372.00	<a href="http://www.chromtech.net.au/pdf2/vici09_131.pdf">http://www.chromtech.net.au/pdf2/vici09_131.pdf</a>
DSTF10MWE	Replacement valve 10 Position 1/8in fittings 1.0mm ports 4in standoff	Valco Valves	\$1,167.00	<a href="http://www.chromtech.net.au/pdf2/vici09_131.pdf">http://www.chromtech.net.au/pdf2/vici09_131.pdf</a>
E2CSTF16MWE	16 Port Selector Trapping FloThru STF w/Std Electric 1/16x0.75mm ports	Valco Valves	\$2,301.00	<a href="http://www.chromtech.net.au/pdf2/vici09_130.pdf">http://www.chromtech.net.au/pdf2/vici09_130.pdf</a>
EMT2CSTF6MWE	6 Port Selector Trapping FloThru STF w/Microelectric 1/16x0.75mm ports	Valco Valves	\$2,372.00	<a href="http://www.chromtech.net.au/pdf2/vici09_130.pdf">http://www.chromtech.net.au/pdf2/vici09_130.pdf</a>
EMT2CSTF10MWE	10 Port Selector Trapping FloThru STF w/Microelectric 1/16x0.75mm ports	Valco Valves	\$2,480.00	<a href="http://www.chromtech.net.au/pdf2/vici09_130.pdf">http://www.chromtech.net.au/pdf2/vici09_130.pdf</a>
EMT2CSTF12MWE	12 Port Selector Trapping FloThru STF w/Microelectric 1/16x0.75mm ports	Valco Valves	\$2,582.00	<a href="http://www.chromtech.net.au/pdf2/vici09_130.pdf">http://www.chromtech.net.au/pdf2/vici09_130.pdf</a>
SSACSTF12MWE	Replacement rotor 12 Position 1/16in fittings 0.75mm ports 4in standoff	Valco Valves	\$147.00	<a href="http://www.chromtech.net.au/pdf2/vici09_130.pdf">http://www.chromtech.net.au/pdf2/vici09_130.pdf</a>
SSACSTF16MWE	Replacement rotor 16 Position 1/16in fittings 0.75mm ports 4in standoff	Valco Valves	\$147.00	<a href="http://www.chromtech.net.au/pdf2/vici09_130.pdf">http://www.chromtech.net.au/pdf2/vici09_130.pdf</a>
E2CSTF12MWE	12 Port Selector Trapping FloThru STF w/Std Electric 1/16x0.75mm ports	Valco Valves	\$2,301.00	<a href="http://www.chromtech.net.au/pdf2/vici09_130.pdf">http://www.chromtech.net.au/pdf2/vici09_130.pdf</a>
2CSTF10MWE	10 Port Selector Trapping FloThru STF manual 1/16x0.75mm ports	Valco Valves	\$1,371.00	<a href="http://www.chromtech.net.au/pdf2/vici09_130.pdf">http://www.chromtech.net.au/pdf2/vici09_130.pdf</a>
SSACSTF6MWE	Replacement rotor 6 Position 1/16in fittings 0.75mm ports 4in standoff	Valco Valves	\$134.00	<a href="http://www.chromtech.net.au/pdf2/vici09_130.pdf">http://www.chromtech.net.au/pdf2/vici09_130.pdf</a>
E2CSTF10MWE	10 Port Selector Trapping FloThru STF w/Air actuator 1/16x0.75mm ports	Valco Valves	\$2,199.00	<a href="http://www.chromtech.net.au/pdf2/vici09_130.pdf">http://www.chromtech.net.au/pdf2/vici09_130.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

E2CSTF6MWE	6 Port Selector Trapping FloThru STF w/Std Electric 1/16x0.75mm ports	Valco Valves	\$2,091.00	<a href="http://www.chromtech.net.au/pdf2/vici09_130.pdf">http://www.chromtech.net.au/pdf2/vici09_130.pdf</a>
A2CSTF16MWE	16 Port Selector Trapping FloThru STF w/Air actuator 1/16x0.75mm ports	Valco Valves	\$1,912.00	<a href="http://www.chromtech.net.au/pdf2/vici09_130.pdf">http://www.chromtech.net.au/pdf2/vici09_130.pdf</a>
A2CSTF12MWE	12 Port Selector Trapping FloThru STF w/Air actuator 1/16x0.75mm ports	Valco Valves	\$1,912.00	<a href="http://www.chromtech.net.au/pdf2/vici09_130.pdf">http://www.chromtech.net.au/pdf2/vici09_130.pdf</a>
A2CSTF10MWE	10 Port Selector Trapping FloThru STF manual 1/16x0.75mm ports	Valco Valves	\$1,683.00	<a href="http://www.chromtech.net.au/pdf2/vici09_130.pdf">http://www.chromtech.net.au/pdf2/vici09_130.pdf</a>
A2CSTF6MWE	6 Port Selector Trapping FloThru STF w/Air actuator 1/16x0.75mm ports	Valco Valves	\$1,575.00	<a href="http://www.chromtech.net.au/pdf2/vici09_130.pdf">http://www.chromtech.net.au/pdf2/vici09_130.pdf</a>
2CSTF12MWE	12 Port Selector Trapping FloThru STF manual 1/16x0.75mm ports	Valco Valves	\$1,473.00	<a href="http://www.chromtech.net.au/pdf2/vici09_130.pdf">http://www.chromtech.net.au/pdf2/vici09_130.pdf</a>
2CSTF6MWE	6 Port Selector Trapping FloThru STF manual 1/16x0.75mm ports	Valco Valves	\$1,263.00	<a href="http://www.chromtech.net.au/pdf2/vici09_130.pdf">http://www.chromtech.net.au/pdf2/vici09_130.pdf</a>
EMT2CSTF16MWE	16 Port Selector Trapping FloThru STF w/Microelectric 1/16x0.75mm ports	Valco Valves	\$2,582.00	<a href="http://www.chromtech.net.au/pdf2/vici09_130.pdf">http://www.chromtech.net.au/pdf2/vici09_130.pdf</a>
2CSTF16MWE	16 Port Selector Trapping FloThru STF manual 1/16x0.75mm ports	Valco Valves	\$1,473.00	<a href="http://www.chromtech.net.au/pdf2/vici09_130.pdf">http://www.chromtech.net.au/pdf2/vici09_130.pdf</a>
SSACSTF10MWE	Replacement rotor 10 Position 1/16in fittings 0.75mm ports 4in standoff	Valco Valves	\$147.00	<a href="http://www.chromtech.net.au/pdf2/vici09_130.pdf">http://www.chromtech.net.au/pdf2/vici09_130.pdf</a>
ZF6PK-10	PEEK 3/8in Polymeric ferrules	High Pressure PEEK Fittings	\$66.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>
ZF4PK-10	PEEK 1/4in Polymeric ferrules	High Pressure PEEK Fittings	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>
ZF8PK-10	PEEK 1/2in Polymeric ferrules	High Pressure PEEK Fittings	\$90.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>
ZF1PK-10	PEEK Polymeric Ferrules 1/16in	High Pressure PEEK Fittings	\$42.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>
ZF2PK-10	PEEK Polymeric Ferrules 1/8in	High Pressure PEEK Fittings	\$42.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>
ZF.5PK-10	PEEK 1/32in Polymeric ferrules	High Pressure PEEK Fittings	\$42.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>
ZF2PFA-10	PFA Polymeric Ferrules 1/8in	Valco Fittings	\$44.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>
ZF8TF-10	PTFE Virgin 1/2in Polymeric ferrules	Valco Fittings	\$57.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>
ZF8TFG-10	PTFE Glass-filled 1/2in Polymeric ferrules	Valco Fittings	\$57.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>
ZF8FEP-10	FEP 1/2in Polymeric ferrules	Valco Fittings	\$57.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>
ZF.5PFA-10	PFA Polymeric Ferrules 1/32in	Valco Fittings	\$51.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

ZF.5KF-10	CTFE 1/32in Polymeric ferrules	Valco Fittings	\$51.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>
ZF8PFA-10	PFA 1/2in Polymeric ferrules	Valco Fittings	\$90.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>
ZF1KF-10	CTFE 1/16in Polymeric ferrules	Valco Fittings	\$44.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>
ZF2KF-10	CTFE 1/8in Polymeric ferrules	Valco Fittings	\$44.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>
ZF4PFA-10	PFA 1/4in Polymeric ferrules	Valco Fittings	\$44.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>
ZF4KF-10	CTFE 1/4in Polymeric ferrules	Valco Fittings	\$44.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>
ZF6PFA-10	PFA 3/8in Polymeric ferrules	Valco Fittings	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>
ZF6KF-10	CTFE 3/8in Polymeric ferrules	Valco Fittings	\$63.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>
ZF1PFA-10	PFA 1/16in Polymeric ferrules	Valco Fittings	\$44.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>
ZF2TF-10	PTFE Virgin 1/8in Polymeric ferrules	Valco Fittings	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>
ZF8V1-5	1/2" PolyI,5pk	Valco Fittings	\$204.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>
ZF8V-5	1/2" PolyI/Valcon,5pk	Valco Fittings	\$141.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>
ZF8GV-5	1/2" PolyI/Graph,5pk	Valco Fittings	\$141.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>
ZF.5TF-10	PTFE Virgin 1/32in Polymeric ferrules	Valco Fittings	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>
ZF.5TFG-10	PTFE Glass-filled 1/32in Polymeric ferrules	Valco Fittings	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>
ZF.5FEP-10	FEP 1/32in Polymeric ferrules	Valco Fittings	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>
ZF1TF-10	PTFE Virgin 1/16in Polymeric ferrules	Valco Fittings	\$22.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>
ZF6FEP-10	FEP 3/8in Polymeric ferrules	Valco Fittings	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>
ZF1FEP-10	FEP 1/16in Polymeric ferrules	Valco Fittings	\$22.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>
ZF6TFG-10	PTFE Glass-filled 3/8in Polymeric ferrules	Valco Fittings	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>
ZF2TFG-10	PTFE Glass-filled 1/8in Polymeric ferrules	Valco Fittings	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>
ZF2FEP-10	FEP 1/8in Polymeric ferrules	Valco Fittings	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>
ZF4TF-10	PTFE Virgin 1/4in Polymeric ferrules	Valco Fittings	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>
ZF4TFG-10	PTFE Glass-filled 1/4in Polymeric ferrules	Valco Fittings	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>
ZF4FEP-10	FEP 1/4in Polymeric ferrules	Valco Fittings	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>
ZF6TF-10	PTFE Virgin 3/8in Polymeric ferrules	Valco Fittings	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>
ZF1GV-5	1/16" PolyI/Graph,5pk	Valco Fittings	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>
ZF1TFG-10	PTFE Glass-filled 1/16in Polymeric ferrules	Valco Fittings	\$22.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>
ZF4V-5	Polyimide Valcon 1/4in Polymeric Ferrules (5pk)	Valco Fittings	\$57.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>
ZF.5V-5	Valcon Polyimide 1/32in Adapters tube OD 0.5mm to 0.8mm	Valco Fittings	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>
ZF8KF-10	CTFE 1/2in Polymeric ferrules	Valco Fittings	\$90.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>
ZF6V-5	3/8" PolyI/Valcon,5pk	Valco Fittings	\$82.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>
ZF4V1-5	1/4" PolyI,5pk	Valco Fittings	\$153.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>
ZF4GV-5	1/4" PolyI/Graph,5pk	Valco Fittings	\$57.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

ZF2V1-5	1/8" PolyI,5pk	Valco Fittings	\$114.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>
ZF6V1-5	3/8" PolyI,5pk	Valco Fittings	\$192.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>
ZF.5V1-5	Virgin Polyimide 1/32in Adapters tube OD 0.5mm to 0.8mm	Valco Fittings	\$87.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>
ZF2GV-5	1/8" PolyI/Graph,5pk	Valco Fittings	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>
ZF1V1-5	1/16" PolyI,5pk	Valco Fittings	\$82.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>
ZF1V-5	1/16" PolyI/Valcon,5pk	Valco Fittings	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>
ZF2V-5	1/8" PolyI/Valcon,5pk	Valco Fittings	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>
ZF6GV-5	3/8" PolyI/Graph,5pk	Valco Fittings	\$82.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>
ZF.5GV-5	1/32" PolyI/Graph,5pk	Valco Fittings	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_13.pdf">http://www.chromtech.net.au/pdf2/vici09_13.pdf</a>
SL100STP	100ul 1/8in Stainless steel	Valco Valves	\$36.00	<a href="http://www.chromtech.net.au/pdf2/vici09_129.pdf">http://www.chromtech.net.au/pdf2/vici09_129.pdf</a>
DST12MWE	Replacement valve 12 Position 1/8in fittings 0.75mm ports 4in standoff	Valco Valves	\$1,236.00	<a href="http://www.chromtech.net.au/pdf2/vici09_129.pdf">http://www.chromtech.net.au/pdf2/vici09_129.pdf</a>
DST16MWE	Replacement valve 16 Position 1/8in fittings 0.75mm 4in standoff ports	Valco Valves	\$1,236.00	<a href="http://www.chromtech.net.au/pdf2/vici09_129.pdf">http://www.chromtech.net.au/pdf2/vici09_129.pdf</a>
SSAST6MWE	Replacement rotor 6 Position 1/8in fittings 0.75mm ports 4in standoff	Valco Valves	\$134.00	<a href="http://www.chromtech.net.au/pdf2/vici09_129.pdf">http://www.chromtech.net.au/pdf2/vici09_129.pdf</a>
SSAST10MWE	Replacement rotor 10 Position 1/8in fittings 0.75mm ports 4in standoff	Valco Valves	\$147.00	<a href="http://www.chromtech.net.au/pdf2/vici09_129.pdf">http://www.chromtech.net.au/pdf2/vici09_129.pdf</a>
2ST6MWE	6 Port Selector Trapping ST Manual 1/8x1.0mm ports	Valco Valves	\$1,230.00	<a href="http://www.chromtech.net.au/pdf2/vici09_129.pdf">http://www.chromtech.net.au/pdf2/vici09_129.pdf</a>
SSACST16MWE	10 Port Selector Trapping ST Replace Rotor 1/8x1.0mm ports	Valco Valves	\$147.00	<a href="http://www.chromtech.net.au/pdf2/vici09_129.pdf">http://www.chromtech.net.au/pdf2/vici09_129.pdf</a>
EMT2ST10MWE	10 Port Selector Trapping ST w/Microelectric actuator 1/8x1.0mm ports	Valco Valves	\$2,448.00	<a href="http://www.chromtech.net.au/pdf2/vici09_129.pdf">http://www.chromtech.net.au/pdf2/vici09_129.pdf</a>
SL1KSTP	1ml 1/8in Stainless steel	Valco Valves	\$42.00	<a href="http://www.chromtech.net.au/pdf2/vici09_129.pdf">http://www.chromtech.net.au/pdf2/vici09_129.pdf</a>
SL250STP	250ul 1/8in Stainless steel loops inc fittings	Valco Valves	\$36.00	<a href="http://www.chromtech.net.au/pdf2/vici09_129.pdf">http://www.chromtech.net.au/pdf2/vici09_129.pdf</a>
SL2KSTP	2ml 1/8in Stainless steel loops inc fittings	Valco Valves	\$44.00	<a href="http://www.chromtech.net.au/pdf2/vici09_129.pdf">http://www.chromtech.net.au/pdf2/vici09_129.pdf</a>
SL500STP	500ul 1/8in Stainless steel	Valco Valves	\$40.00	<a href="http://www.chromtech.net.au/pdf2/vici09_129.pdf">http://www.chromtech.net.au/pdf2/vici09_129.pdf</a>
SL5KSTP	5ml 1/8in Stainless steel	Valco Valves	\$52.00	<a href="http://www.chromtech.net.au/pdf2/vici09_129.pdf">http://www.chromtech.net.au/pdf2/vici09_129.pdf</a>
SL10KSTP	10ml 1/8in Stainless steel	Valco Valves	\$74.00	<a href="http://www.chromtech.net.au/pdf2/vici09_129.pdf">http://www.chromtech.net.au/pdf2/vici09_129.pdf</a>
SSAST12MWE	Replacement rotor 12 Position 1/8in fittings 0.75mm ports 4in standoff	Valco Valves	\$147.00	<a href="http://www.chromtech.net.au/pdf2/vici09_129.pdf">http://www.chromtech.net.au/pdf2/vici09_129.pdf</a>
EMT2ST16MWE	16 Port Selector Trapping ST w/Microelectric actuator 1/8x1.0mm ports	Valco Valves	\$2,550.00	<a href="http://www.chromtech.net.au/pdf2/vici09_129.pdf">http://www.chromtech.net.au/pdf2/vici09_129.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

EMT2ST6MWE	6 Port Selector Trapping ST w/Microelectric actuator 1/8x1.0mm ports	Valco Valves	\$2,340.00	<a href="http://www.chromtech.net.au/pdf2/vici09_129.pdf">http://www.chromtech.net.au/pdf2/vici09_129.pdf</a>
E2ST16MWE	16 Port Selector Trapping ST w/Std Electric actuator 1/8x1.0mm ports	Valco Valves	\$2,270.00	<a href="http://www.chromtech.net.au/pdf2/vici09_129.pdf">http://www.chromtech.net.au/pdf2/vici09_129.pdf</a>
E2ST12MWE	12 Port Selector Trapping ST w/Std Electric actuator 1/8x1.0mm ports	Valco Valves	\$2,270.00	<a href="http://www.chromtech.net.au/pdf2/vici09_129.pdf">http://www.chromtech.net.au/pdf2/vici09_129.pdf</a>
E2ST10MWE	10 Port Selector Trapping ST w/Std Electric actuator 1/8x1.0mm ports	Valco Valves	\$2,168.00	<a href="http://www.chromtech.net.au/pdf2/vici09_129.pdf">http://www.chromtech.net.au/pdf2/vici09_129.pdf</a>
E2ST6MWE	6 Port Selector Trapping ST w/Std Electric actuator 1/8x1.0mm ports	Valco Valves	\$2,060.00	<a href="http://www.chromtech.net.au/pdf2/vici09_129.pdf">http://www.chromtech.net.au/pdf2/vici09_129.pdf</a>
A2ST16MWE	6Port Selector Trapping ST w/Air actuator 1/8x1.0mm ports	Valco Valves	\$1,754.00	<a href="http://www.chromtech.net.au/pdf2/vici09_129.pdf">http://www.chromtech.net.au/pdf2/vici09_129.pdf</a>
A2ST12MWE	12 Port Selector Trapping ST w/Air actuator 1/8x1.0mm ports	Valco Valves	\$1,754.00	<a href="http://www.chromtech.net.au/pdf2/vici09_129.pdf">http://www.chromtech.net.au/pdf2/vici09_129.pdf</a>
A2ST10MWE	10 Port Selector Trapping ST w/Air actuator 1/8x1.0mm ports	Valco Valves	\$1,652.00	<a href="http://www.chromtech.net.au/pdf2/vici09_129.pdf">http://www.chromtech.net.au/pdf2/vici09_129.pdf</a>
A2ST6MWE	10 Port Selector Trapping ST w/Air actuator 1/8x1.0mm ports	Valco Valves	\$1,542.00	<a href="http://www.chromtech.net.au/pdf2/vici09_129.pdf">http://www.chromtech.net.au/pdf2/vici09_129.pdf</a>
2ST16MWE	16 Port Selector Trapping ST Manual 1/8x1.0mm ports	Valco Valves	\$1,440.00	<a href="http://www.chromtech.net.au/pdf2/vici09_129.pdf">http://www.chromtech.net.au/pdf2/vici09_129.pdf</a>
2ST12MWE	12 Port Selector Trapping ST Manual 1/8x1.0mm ports	Valco Valves	\$1,440.00	<a href="http://www.chromtech.net.au/pdf2/vici09_129.pdf">http://www.chromtech.net.au/pdf2/vici09_129.pdf</a>
2ST10MWE	10 Port Selector Trapping ST Manual 1/8x1.0mm ports	Valco Valves	\$1,338.00	<a href="http://www.chromtech.net.au/pdf2/vici09_129.pdf">http://www.chromtech.net.au/pdf2/vici09_129.pdf</a>
EMT2ST12MWE	12 Port Selector Trapping ST w/Microelectric actuator 1/8x1.0mm ports	Valco Valves	\$2,550.00	<a href="http://www.chromtech.net.au/pdf2/vici09_129.pdf">http://www.chromtech.net.au/pdf2/vici09_129.pdf</a>
DST10MWE	Replacement valve 10 Position 1/8in fittings 0.75mm ports 4in standoff	Valco Valves	\$1,134.00	<a href="http://www.chromtech.net.au/pdf2/vici09_129.pdf">http://www.chromtech.net.au/pdf2/vici09_129.pdf</a>
DST6MWE	Replacement valve 6 Position 1/8in fittings 0.75mm ports 4in standoff	Valco Valves	\$1,026.00	<a href="http://www.chromtech.net.au/pdf2/vici09_129.pdf">http://www.chromtech.net.au/pdf2/vici09_129.pdf</a>
SSACST16MWE	Replacement rotor 16 Position 1/16in fittings 0.75mm ports 4in standoff	Valco Valves	\$147.00	<a href="http://www.chromtech.net.au/pdf2/vici09_128.pdf">http://www.chromtech.net.au/pdf2/vici09_128.pdf</a>
2CST6MWE	6 Port Selector Trapping ST Manual 1/16x0.75mm ports	Valco Valves	\$1,230.00	<a href="http://www.chromtech.net.au/pdf2/vici09_128.pdf">http://www.chromtech.net.au/pdf2/vici09_128.pdf</a>
DCST10MWE	Replacement valve 10 Position 1/16in fittings 0.75mm ports 4in standoff	Valco Valves	\$1,134.00	<a href="http://www.chromtech.net.au/pdf2/vici09_128.pdf">http://www.chromtech.net.au/pdf2/vici09_128.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

DCST12MWE	Replacement valve 12 Position 1/16in fittings 0.75mm ports 4in standoff	Valco Valves	\$1,236.00	<a href="http://www.chromtech.net.au/pdf2/vici09_128.pdf">http://www.chromtech.net.au/pdf2/vici09_128.pdf</a>
DCST16MWE	Replacement valve 16 Position 1/16in fittings 0.75mm ports	Valco Valves	\$1,236.00	<a href="http://www.chromtech.net.au/pdf2/vici09_128.pdf">http://www.chromtech.net.au/pdf2/vici09_128.pdf</a>
SSACST6MWE	Replacement rotor 6 Position 1/16in fittings 0.75mm ports 4in standoff	Valco Valves	\$134.00	<a href="http://www.chromtech.net.au/pdf2/vici09_128.pdf">http://www.chromtech.net.au/pdf2/vici09_128.pdf</a>
EMT2CST16MWE	16 Port Selector Trapping ST w/Microelectric actuator 1/16x0.75mm ports	Valco Valves	\$2,550.00	<a href="http://www.chromtech.net.au/pdf2/vici09_128.pdf">http://www.chromtech.net.au/pdf2/vici09_128.pdf</a>
SSACST12MWE	Replacement rotor 12 Position 1/16in fittings 0.75mm ports 4in standoff	Valco Valves	\$147.00	<a href="http://www.chromtech.net.au/pdf2/vici09_128.pdf">http://www.chromtech.net.au/pdf2/vici09_128.pdf</a>
SL50CSTP	50ul 1/16in Stainless steel loops inc fittings	Valco Valves	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_128.pdf">http://www.chromtech.net.au/pdf2/vici09_128.pdf</a>
SL1KCSTP	1ml 1/16in Stainless steel loops inc fittings	Valco Valves	\$46.00	<a href="http://www.chromtech.net.au/pdf2/vici09_128.pdf">http://www.chromtech.net.au/pdf2/vici09_128.pdf</a>
SL100CSTP	100ul 1/16in Stainless steel loops inc fittings	Valco Valves	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_128.pdf">http://www.chromtech.net.au/pdf2/vici09_128.pdf</a>
SL2KCSTP	2ml 1/16in Stainless steel loops inc fittings	Valco Valves	\$58.00	<a href="http://www.chromtech.net.au/pdf2/vici09_128.pdf">http://www.chromtech.net.au/pdf2/vici09_128.pdf</a>
SL250CSTP	250ul 1/16in Stainless steel loops inc fittings	Valco Valves	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_128.pdf">http://www.chromtech.net.au/pdf2/vici09_128.pdf</a>
SL5KCSTP	5ml 1/16in Stainless steel loops inc fittings	Valco Valves	\$69.00	<a href="http://www.chromtech.net.au/pdf2/vici09_128.pdf">http://www.chromtech.net.au/pdf2/vici09_128.pdf</a>
SSACST10MWE	Replacement rotor 10 Position 1/16in fittings 0.75mm ports 4in standoff	Valco Valves	\$147.00	<a href="http://www.chromtech.net.au/pdf2/vici09_128.pdf">http://www.chromtech.net.au/pdf2/vici09_128.pdf</a>
E2CST6MWE	6 Port Selector Trapping ST w/Std Electric actuator 1/16x0.75mm ports	Valco Valves	\$2,060.00	<a href="http://www.chromtech.net.au/pdf2/vici09_128.pdf">http://www.chromtech.net.au/pdf2/vici09_128.pdf</a>
2CST10MWE	10 Port Selector Trapping ST Manual 1/16x0.75mm ports	Valco Valves	\$1,338.00	<a href="http://www.chromtech.net.au/pdf2/vici09_128.pdf">http://www.chromtech.net.au/pdf2/vici09_128.pdf</a>
2CST12MWE	12 Port Selector Trapping ST Manual 1/16x0.75mm ports	Valco Valves	\$1,440.00	<a href="http://www.chromtech.net.au/pdf2/vici09_128.pdf">http://www.chromtech.net.au/pdf2/vici09_128.pdf</a>
2CST16MWE	16 Port Selector Trapping ST Manual 1/16x0.75mm ports	Valco Valves	\$1,440.00	<a href="http://www.chromtech.net.au/pdf2/vici09_128.pdf">http://www.chromtech.net.au/pdf2/vici09_128.pdf</a>
A2CST6MWE	6 Port Selector Trapping ST w/Air actuator 1/16x0.75mm ports	Valco Valves	\$1,542.00	<a href="http://www.chromtech.net.au/pdf2/vici09_128.pdf">http://www.chromtech.net.au/pdf2/vici09_128.pdf</a>
A2CST10MWE	10 Port Selector Trapping ST w/Air actuator 1/16x0.75mm ports	Valco Valves	\$1,652.00	<a href="http://www.chromtech.net.au/pdf2/vici09_128.pdf">http://www.chromtech.net.au/pdf2/vici09_128.pdf</a>
DCST6MWE	Replacement valve 6 Position 1/16in fittings 0.75mm ports 4in standoff	Valco Valves	\$1,026.00	<a href="http://www.chromtech.net.au/pdf2/vici09_128.pdf">http://www.chromtech.net.au/pdf2/vici09_128.pdf</a>
A2CST16MWE	16 Port Selector Trapping ST w/Air actuator 1/16x0.75mm ports	Valco Valves	\$1,754.00	<a href="http://www.chromtech.net.au/pdf2/vici09_128.pdf">http://www.chromtech.net.au/pdf2/vici09_128.pdf</a>
SL10KCSTP	10ml 1/16in Stainless steel loops inc fittings	Valco Valves	\$90.00	<a href="http://www.chromtech.net.au/pdf2/vici09_128.pdf">http://www.chromtech.net.au/pdf2/vici09_128.pdf</a>
E2CST10MWE	10 Port Selector Trapping ST w/Std Electric actuator 1/16x0.75mm ports	Valco Valves	\$2,168.00	<a href="http://www.chromtech.net.au/pdf2/vici09_128.pdf">http://www.chromtech.net.au/pdf2/vici09_128.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

E2CST12MWE	12 Port Selector Trapping ST w/Std Electric actuator 1/16x0.75mm ports	Valco Valves	\$2,270.00	<a href="http://www.chromtech.net.au/pdf2/vici09_128.pdf">http://www.chromtech.net.au/pdf2/vici09_128.pdf</a>
E2CST16MWE	16 Port Selector Trapping ST w/Std Electric actuator 1/16x0.75mm ports	Valco Valves	\$2,270.00	<a href="http://www.chromtech.net.au/pdf2/vici09_128.pdf">http://www.chromtech.net.au/pdf2/vici09_128.pdf</a>
EMT2CST6MWE	6 Port Selector Trapping ST w/Microelectric actuator 1/16x0.75mm ports	Valco Valves	\$2,340.00	<a href="http://www.chromtech.net.au/pdf2/vici09_128.pdf">http://www.chromtech.net.au/pdf2/vici09_128.pdf</a>
EMT2CST10MWE	10 Port Selector Trapping ST w/Microelectric actuator 1/16x0.75mm ports	Valco Valves	\$2,448.00	<a href="http://www.chromtech.net.au/pdf2/vici09_128.pdf">http://www.chromtech.net.au/pdf2/vici09_128.pdf</a>
EMT2CST12MWE	12 Port Selector Trapping ST w/Microelectric actuator 1/16x0.75mm ports	Valco Valves	\$2,550.00	<a href="http://www.chromtech.net.au/pdf2/vici09_128.pdf">http://www.chromtech.net.au/pdf2/vici09_128.pdf</a>
A2CST12MWE	12 Port Selector Trapping ST w/Air actuator 1/16x0.75mm ports	Valco Valves	\$1,754.00	<a href="http://www.chromtech.net.au/pdf2/vici09_128.pdf">http://www.chromtech.net.au/pdf2/vici09_128.pdf</a>
SL500CSTP	500ul 1/16in Stainless steel loops inc fittings	Valco Valves	\$42.00	<a href="http://www.chromtech.net.au/pdf2/vici09_128.pdf">http://www.chromtech.net.au/pdf2/vici09_128.pdf</a>
DVLSF4MWE2	Replacement valve 4 Position 1/4in fittings 4.0mm ports 4in standoff	Valco Valves	\$1,504.00	<a href="http://www.chromtech.net.au/pdf2/vici09_127.pdf">http://www.chromtech.net.au/pdf2/vici09_127.pdf</a>
SSASF12MWE	Replacement rotor12 Position 1/8in fittings 1.0mm ports	Valco Valves	\$147.00	<a href="http://www.chromtech.net.au/pdf2/vici09_127.pdf">http://www.chromtech.net.au/pdf2/vici09_127.pdf</a>
SSASF16MWE	Replacement rotor16 Position 1/8in fittings 1.0mm ports	Valco Valves	\$147.00	<a href="http://www.chromtech.net.au/pdf2/vici09_127.pdf">http://www.chromtech.net.au/pdf2/vici09_127.pdf</a>
AH2VLSF4MWE2	4 Port Selector Flow-Thru SF w/Air actuator 1/4inx4.0mm ports	Valco Valves	\$2,307.00	<a href="http://www.chromtech.net.au/pdf2/vici09_127.pdf">http://www.chromtech.net.au/pdf2/vici09_127.pdf</a>
AH2VLSF6MWE2	6 Port Selector Flow-Thru SF w/Air actuator 1/4inx4.0mm ports	Valco Valves	\$2,378.00	<a href="http://www.chromtech.net.au/pdf2/vici09_127.pdf">http://www.chromtech.net.au/pdf2/vici09_127.pdf</a>
AH2VLSF8MWE2	8 Port Selector Flow-Thru SF w/Air actuator 1/4inx4.0mm ports	Valco Valves	\$2,487.00	<a href="http://www.chromtech.net.au/pdf2/vici09_127.pdf">http://www.chromtech.net.au/pdf2/vici09_127.pdf</a>
E2VLSF4MWE2	4 Port Selector Flow-Thru SF w/Std Electric actuator 1/4inx4.0mm ports	Valco Valves	\$2,538.00	<a href="http://www.chromtech.net.au/pdf2/vici09_127.pdf">http://www.chromtech.net.au/pdf2/vici09_127.pdf</a>
E2VLSF6MWE2	6 Port Selector Flow-Thru SF w/Std Electric actuator 1/4inx4.0mm ports	Valco Valves	\$2,607.00	<a href="http://www.chromtech.net.au/pdf2/vici09_127.pdf">http://www.chromtech.net.au/pdf2/vici09_127.pdf</a>
E2VLSF8MWE2	8 Port Selector Flow-Thru SF w/Std Electric actuator 1/4inx4.0mm ports	Valco Valves	\$2,715.00	<a href="http://www.chromtech.net.au/pdf2/vici09_127.pdf">http://www.chromtech.net.au/pdf2/vici09_127.pdf</a>
EMT2VLSF4MWE2	4 Port Selector Flow-Thru SF w/Microelectric actuator 1/4inx4.0mm ports	Valco Valves	\$2,817.00	<a href="http://www.chromtech.net.au/pdf2/vici09_127.pdf">http://www.chromtech.net.au/pdf2/vici09_127.pdf</a>
SSASF10MWE	Replacement rotor10 Position 1/8in fittings 1.0mm ports	Valco Valves	\$147.00	<a href="http://www.chromtech.net.au/pdf2/vici09_127.pdf">http://www.chromtech.net.au/pdf2/vici09_127.pdf</a>
EMT2VLSF8MWE2	8 Port Selector Flow-Thru SF w/Microelectric actuator 1/4inx4.0mm ports	Valco Valves	\$2,997.00	<a href="http://www.chromtech.net.au/pdf2/vici09_127.pdf">http://www.chromtech.net.au/pdf2/vici09_127.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

EMT2SF10MWE	10 Port Selector Flow-Thru SF w/Microelectric actuator 1/16inx1.0mm ports	Valco Valves	\$2,448.00	<a href="http://www.chromtech.net.au/pdf2/vici09_127.pdf">http://www.chromtech.net.au/pdf2/vici09_127.pdf</a>
DVLSF6MWE2	Replacement valve 6 Position 1/4in fittings 4.0mm ports 4in standoff	Valco Valves	\$1,575.00	<a href="http://www.chromtech.net.au/pdf2/vici09_127.pdf">http://www.chromtech.net.au/pdf2/vici09_127.pdf</a>
SSAVLSF4MWE2	Replacement rotor 4 Position 1/4in fittings 4.0mm ports 4in standoff	Valco Valves	\$243.00	<a href="http://www.chromtech.net.au/pdf2/vici09_127.pdf">http://www.chromtech.net.au/pdf2/vici09_127.pdf</a>
SSAVLSF8MWE2	Replacement rotor 8 Position 1/4in fittings 4.0mm ports 4in standoff	Valco Valves	\$243.00	<a href="http://www.chromtech.net.au/pdf2/vici09_127.pdf">http://www.chromtech.net.au/pdf2/vici09_127.pdf</a>
DSF12MWE	Replacement valve 12 Position 1/8in fittings 1.0mm ports	Valco Valves	\$1,236.00	<a href="http://www.chromtech.net.au/pdf2/vici09_127.pdf">http://www.chromtech.net.au/pdf2/vici09_127.pdf</a>
DSF16MWE	Replacement valve 16 Position 1/8in fittings 1.0mm ports	Valco Valves	\$1,236.00	<a href="http://www.chromtech.net.au/pdf2/vici09_127.pdf">http://www.chromtech.net.au/pdf2/vici09_127.pdf</a>
EMT2VLSF6MWE2	6 Port Selector Flow-Thru SF w/Microelectric actuator 1/4inx4.0mm ports	Valco Valves	\$2,888.00	<a href="http://www.chromtech.net.au/pdf2/vici09_127.pdf">http://www.chromtech.net.au/pdf2/vici09_127.pdf</a>
A2SF16MWE	16 Port Selector Flow-Thru SF w/Air actuator 1/8inx1.0mm ports	Valco Valves	\$1,754.00	<a href="http://www.chromtech.net.au/pdf2/vici09_127.pdf">http://www.chromtech.net.au/pdf2/vici09_127.pdf</a>
SSASF6MWE	Replacement rotor 6 Position 1/8in fittings 1.0mm ports	Valco Valves	\$134.00	<a href="http://www.chromtech.net.au/pdf2/vici09_127.pdf">http://www.chromtech.net.au/pdf2/vici09_127.pdf</a>
DVLSF8MWE2	Replacement valve 8 Position 1/4in fittings 4.0mm ports 4in standoff	Valco Valves	\$1,683.00	<a href="http://www.chromtech.net.au/pdf2/vici09_127.pdf">http://www.chromtech.net.au/pdf2/vici09_127.pdf</a>
2SF10MWE	10 Port Selector Flow-Thru SF w/Manual 1/8inx1.0mm ports	Valco Valves	\$1,338.00	<a href="http://www.chromtech.net.au/pdf2/vici09_127.pdf">http://www.chromtech.net.au/pdf2/vici09_127.pdf</a>
2SF12MWE	12 Port Selector Flow-Thru SF w/Manual 1/8inx1.0mm ports	Valco Valves	\$1,440.00	<a href="http://www.chromtech.net.au/pdf2/vici09_127.pdf">http://www.chromtech.net.au/pdf2/vici09_127.pdf</a>
2SF16MWE	16 Port Selector Flow-Thru SF w/Manual 1/8inx1.0mm ports	Valco Valves	\$1,440.00	<a href="http://www.chromtech.net.au/pdf2/vici09_127.pdf">http://www.chromtech.net.au/pdf2/vici09_127.pdf</a>
DSF10MWE	Replacement valve 10 Position 1/8in fittings 1.0mm ports	Valco Valves	\$1,134.00	<a href="http://www.chromtech.net.au/pdf2/vici09_127.pdf">http://www.chromtech.net.au/pdf2/vici09_127.pdf</a>
A2SF6MWE	6 Port Selector Flow-Thru SF w/Air actuator 1/8inx1.0mm ports	Valco Valves	\$1,542.00	<a href="http://www.chromtech.net.au/pdf2/vici09_127.pdf">http://www.chromtech.net.au/pdf2/vici09_127.pdf</a>
EMT2SF16MWE	16 Port Selector Flow-Thru SF w/Microelectric actuator 1/16inx1.0mm ports	Valco Valves	\$2,550.00	<a href="http://www.chromtech.net.au/pdf2/vici09_127.pdf">http://www.chromtech.net.au/pdf2/vici09_127.pdf</a>
A2SF12MWE	12 Port Selector Flow-Thru SF w/Air actuator 1/8inx1.0mm ports	Valco Valves	\$1,754.00	<a href="http://www.chromtech.net.au/pdf2/vici09_127.pdf">http://www.chromtech.net.au/pdf2/vici09_127.pdf</a>
DSF6MWE	Replacement valve 6 Position 1/8in fittings 1.0mm ports	Valco Valves	\$1,026.00	<a href="http://www.chromtech.net.au/pdf2/vici09_127.pdf">http://www.chromtech.net.au/pdf2/vici09_127.pdf</a>



## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

E2SF6MWE	6 Port Selector Flow-Thru SF w/Std Electric actuator 1/16inx1.0mm ports	Valco Valves	\$2,060.00	<a href="http://www.chromtech.net.au/pdf2/vici09_127.pdf">http://www.chromtech.net.au/pdf2/vici09_127.pdf</a>
E2SF10MWE	10 Port Selector Flow-Thru SF w/Std Electric actuator 1/16inx1.0mm ports	Valco Valves	\$2,168.00	<a href="http://www.chromtech.net.au/pdf2/vici09_127.pdf">http://www.chromtech.net.au/pdf2/vici09_127.pdf</a>
E2SF12MWE	12 Port Selector Flow-Thru SF w/Std Electric actuator 1/16inx1.0mm ports	Valco Valves	\$2,270.00	<a href="http://www.chromtech.net.au/pdf2/vici09_127.pdf">http://www.chromtech.net.au/pdf2/vici09_127.pdf</a>
E2SF16MWE	16 Port Selector Flow-Thru SF w/Std Electric actuator 1/16inx1.0mm ports	Valco Valves	\$2,270.00	<a href="http://www.chromtech.net.au/pdf2/vici09_127.pdf">http://www.chromtech.net.au/pdf2/vici09_127.pdf</a>
EMT2SF6MWE	6 Port Selector Flow-Thru SF w/Microelectric actuator 1/16inx1.0mm ports	Valco Valves	\$2,340.00	<a href="http://www.chromtech.net.au/pdf2/vici09_127.pdf">http://www.chromtech.net.au/pdf2/vici09_127.pdf</a>
EMT2SF12MWE	12 Port Selector Flow-Thru SF w/Microelectric actuator 1/16inx1.0mm ports	Valco Valves	\$2,550.00	<a href="http://www.chromtech.net.au/pdf2/vici09_127.pdf">http://www.chromtech.net.au/pdf2/vici09_127.pdf</a>
2SF6MWE	6 Port Selector Flow-Thru SF w/Manual 1/8inx1.0mm ports	Valco Valves	\$1,230.00	<a href="http://www.chromtech.net.au/pdf2/vici09_127.pdf">http://www.chromtech.net.au/pdf2/vici09_127.pdf</a>
A2SF10MWE	10 Port Selector Flow-Thru SF w/Air actuator 1/8inx1.0mm ports	Valco Valves	\$1,652.00	<a href="http://www.chromtech.net.au/pdf2/vici09_127.pdf">http://www.chromtech.net.au/pdf2/vici09_127.pdf</a>
SSAVLSF6MWE2	Replacement rotor 6 Position 1/4in fittings 4.0mm ports 4in standoff	Valco Valves	\$243.00	<a href="http://www.chromtech.net.au/pdf2/vici09_127.pdf">http://www.chromtech.net.au/pdf2/vici09_127.pdf</a>
2CSF12MWE	12 Port Selector Flow-Thru SF Manual 1/16inx1.0mm ports	Valco Valves	\$1,440.00	<a href="http://www.chromtech.net.au/pdf2/vici09_126.pdf">http://www.chromtech.net.au/pdf2/vici09_126.pdf</a>
EMT2CSF12MWE	12 Port Selector Flow-Thru SFw/Microelectric actuator 1/16inx1.0mm ports	Valco Valves	\$2,550.00	<a href="http://www.chromtech.net.au/pdf2/vici09_126.pdf">http://www.chromtech.net.au/pdf2/vici09_126.pdf</a>
E2CSF6MWE	6 Port Selector Flow-Thru SFw/Std Electric actuator 1/16inx1.0mm ports	Valco Valves	\$2,060.00	<a href="http://www.chromtech.net.au/pdf2/vici09_126.pdf">http://www.chromtech.net.au/pdf2/vici09_126.pdf</a>
EMT2CSF6MWE	6 Port Selector Flow-Thru SFw/Microelectric actuator 1/16inx1.0mm ports	Valco Valves	\$2,340.00	<a href="http://www.chromtech.net.au/pdf2/vici09_126.pdf">http://www.chromtech.net.au/pdf2/vici09_126.pdf</a>
E2CSF16MWE	16 Port Selector Flow-Thru SFw/Std Electric actuator 1/16inx1.0mm ports	Valco Valves	\$2,270.00	<a href="http://www.chromtech.net.au/pdf2/vici09_126.pdf">http://www.chromtech.net.au/pdf2/vici09_126.pdf</a>
EMT2CSF16MWE	16 Port Selector Flow-Thru SFw/Microelectric actuator 1/16inx1.0mm ports	Valco Valves	\$2,550.00	<a href="http://www.chromtech.net.au/pdf2/vici09_126.pdf">http://www.chromtech.net.au/pdf2/vici09_126.pdf</a>
E2CSF12MWE	12 Port Selector Flow-Thru SF w/Std Electric actuator 1/16inx1.0mm ports	Valco Valves	\$2,270.00	<a href="http://www.chromtech.net.au/pdf2/vici09_126.pdf">http://www.chromtech.net.au/pdf2/vici09_126.pdf</a>
E2CSF10MWE	10 Port Selector Flow-Thru SF w/Std Electric actuator 1/16inx1.0mm ports	Valco Valves	\$2,168.00	<a href="http://www.chromtech.net.au/pdf2/vici09_126.pdf">http://www.chromtech.net.au/pdf2/vici09_126.pdf</a>
EMT2CSF10MWE	10 Port Selector Flow-Thru SFw/Microelectric actuator 1/16inx1.0mm ports	Valco Valves	\$2,448.00	<a href="http://www.chromtech.net.au/pdf2/vici09_126.pdf">http://www.chromtech.net.au/pdf2/vici09_126.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

A2CSF16MWE	16 Port Selector Flow-Thru SFw/Air actuator 1/16inx1.0mm ports	Valco Valves	\$1,754.00	<a href="http://www.chromtech.net.au/pdf2/vici09_126.pdf">http://www.chromtech.net.au/pdf2/vici09_126.pdf</a>
A2CSF12MWE	12 Port Selector Flow-Thru SFw/Air actuator 1/16inx1.0mm ports	Valco Valves	\$1,754.00	<a href="http://www.chromtech.net.au/pdf2/vici09_126.pdf">http://www.chromtech.net.au/pdf2/vici09_126.pdf</a>
A2CSF10MWE	10 Port Selector Flow-Thru SFw/Air actuator 1/16inx1.0mm ports	Valco Valves	\$1,652.00	<a href="http://www.chromtech.net.au/pdf2/vici09_126.pdf">http://www.chromtech.net.au/pdf2/vici09_126.pdf</a>
2CSF16MWE	16 Port Selector Flow-Thru SF Manual 1/16inx1.0mm ports	Valco Valves	\$1,440.00	<a href="http://www.chromtech.net.au/pdf2/vici09_126.pdf">http://www.chromtech.net.au/pdf2/vici09_126.pdf</a>
2CSF10MWE	10 Port Selector Flow-Thru SF Manual 1/16inx1.0mm ports	Valco Valves	\$1,338.00	<a href="http://www.chromtech.net.au/pdf2/vici09_126.pdf">http://www.chromtech.net.au/pdf2/vici09_126.pdf</a>
2CSF6MWE	6 Port Selector Flow-Thru SF Manual 1/16inx1.0mm ports	Valco Valves	\$1,230.00	<a href="http://www.chromtech.net.au/pdf2/vici09_126.pdf">http://www.chromtech.net.au/pdf2/vici09_126.pdf</a>
SSACSF16MWE	Replacement rotor 16 Position 1/16in fittings 1.0mm ports 4in standoff	Valco Valves	\$147.00	<a href="http://www.chromtech.net.au/pdf2/vici09_126.pdf">http://www.chromtech.net.au/pdf2/vici09_126.pdf</a>
DCSF6MWE	Replacement valve 6 Position 1/16in fittings 1.0mm ports 4in standoff	Valco Valves	\$1,026.00	<a href="http://www.chromtech.net.au/pdf2/vici09_126.pdf">http://www.chromtech.net.au/pdf2/vici09_126.pdf</a>
A2CSF6MWE	6 Port Selector Flow-Thru SFw/Air actuator 1/16inx1.0mm ports	Valco Valves	\$1,542.00	<a href="http://www.chromtech.net.au/pdf2/vici09_126.pdf">http://www.chromtech.net.au/pdf2/vici09_126.pdf</a>
DCSF12MWE	Replacement valve 12 Position 1/16in fittings 1.0mm ports 4in standoff	Valco Valves	\$1,236.00	<a href="http://www.chromtech.net.au/pdf2/vici09_126.pdf">http://www.chromtech.net.au/pdf2/vici09_126.pdf</a>
DCSF16MWE	Replacement valve 16 Position 1/16in fittings 1.0mm ports 4in standoff	Valco Valves	\$1,236.00	<a href="http://www.chromtech.net.au/pdf2/vici09_126.pdf">http://www.chromtech.net.au/pdf2/vici09_126.pdf</a>
SSACSF6MWE	Replacement rotor 6 Position 1/16in fittings 1.0mm ports 4in standoff	Valco Valves	\$147.00	<a href="http://www.chromtech.net.au/pdf2/vici09_126.pdf">http://www.chromtech.net.au/pdf2/vici09_126.pdf</a>
SSACSF10MWE	Replacement rotor 10 Position 1/16in fittings 1.0mm ports 4in standoff	Valco Valves	\$147.00	<a href="http://www.chromtech.net.au/pdf2/vici09_126.pdf">http://www.chromtech.net.au/pdf2/vici09_126.pdf</a>
SSACSF12MWE	Replacement rotor 12 Position 1/16in fittings 1.0mm ports 4in standoff	Valco Valves	\$147.00	<a href="http://www.chromtech.net.au/pdf2/vici09_126.pdf">http://www.chromtech.net.au/pdf2/vici09_126.pdf</a>
DCSF10MWE	Replacement valve 10 Position 1/16in fittings 1.0mm ports 4in standoff	Valco Valves	\$1,134.00	<a href="http://www.chromtech.net.au/pdf2/vici09_126.pdf">http://www.chromtech.net.au/pdf2/vici09_126.pdf</a>
EMT2VLSC6MWE2	6 Port Selector LowP Common Outlet w/Microelectric actuator 1/4inx4.0mm ports	Valco Valves	\$2,742.00	<a href="http://www.chromtech.net.au/pdf2/vici09_125.pdf">http://www.chromtech.net.au/pdf2/vici09_125.pdf</a>
E2SC12MWE	12 Port Selector LowP Common Outlet w/Std Electric actuator 1/8inx1.0mm ports	Valco Valves	\$2,122.00	<a href="http://www.chromtech.net.au/pdf2/vici09_125.pdf">http://www.chromtech.net.au/pdf2/vici09_125.pdf</a>
E2SC10MWE	4 Port Selector LowP Common Outlet w/Std Electric actuator 1/8inx1.0mm ports	Valco Valves	\$2,020.00	<a href="http://www.chromtech.net.au/pdf2/vici09_125.pdf">http://www.chromtech.net.au/pdf2/vici09_125.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

E2SC6MWE	6 Port Selector LowP Common Outlet w/Std Electric actuator 1/8inx1.0mm ports	Valco Valves	\$1,912.00	<a href="http://www.chromtech.net.au/pdf2/vici09_125.pdf">http://www.chromtech.net.au/pdf2/vici09_125.pdf</a>
A2SC16MWE	16 Port Selector LowP Common Outlet w/Air actuator 1/16inx1.0mm ports	Valco Valves	\$1,606.00	<a href="http://www.chromtech.net.au/pdf2/vici09_125.pdf">http://www.chromtech.net.au/pdf2/vici09_125.pdf</a>
A2SC12MWE	12 Port Selector LowP Common Outlet w/Air actuator 1/16inx1.0mm ports	Valco Valves	\$1,606.00	<a href="http://www.chromtech.net.au/pdf2/vici09_125.pdf">http://www.chromtech.net.au/pdf2/vici09_125.pdf</a>
A2SC10MWE	10 Port Selector LowP Common Outlet w/Air actuator 1/16inx1.0mm ports	Valco Valves	\$1,504.00	<a href="http://www.chromtech.net.au/pdf2/vici09_125.pdf">http://www.chromtech.net.au/pdf2/vici09_125.pdf</a>
2SC16MWE	16 Port Selector LowP Common Outlet Manual 1/16inx1.0mm ports	Valco Valves	\$1,294.00	<a href="http://www.chromtech.net.au/pdf2/vici09_125.pdf">http://www.chromtech.net.au/pdf2/vici09_125.pdf</a>
E2SC16MWE	16 Port Selector LowP Common Outlet w/Std Electric actuator 1/8inx1.0mm ports	Valco Valves	\$2,122.00	<a href="http://www.chromtech.net.au/pdf2/vici09_125.pdf">http://www.chromtech.net.au/pdf2/vici09_125.pdf</a>
2SC6MWE	6 Port Selector LowP Common Outlet Manual 1/16inx1.0mm ports	Valco Valves	\$1,083.00	<a href="http://www.chromtech.net.au/pdf2/vici09_125.pdf">http://www.chromtech.net.au/pdf2/vici09_125.pdf</a>
A2SC6MWE	6 Port Selector LowP Common Outlet w/Air actuator 1/16inx1.0mm ports	Valco Valves	\$1,396.00	<a href="http://www.chromtech.net.au/pdf2/vici09_125.pdf">http://www.chromtech.net.au/pdf2/vici09_125.pdf</a>
EMT2VLSC8MWE2	8 Port Selector LowP Common Outlet w/Microelectric actuator 1/4inx4.0mm ports	Valco Valves	\$2,850.00	<a href="http://www.chromtech.net.au/pdf2/vici09_125.pdf">http://www.chromtech.net.au/pdf2/vici09_125.pdf</a>
DVLSC4MWE2	Replacement valve 4 Position 1/4in fittings 4.0mm ports 4in standoff	Valco Valves	\$1,358.00	<a href="http://www.chromtech.net.au/pdf2/vici09_125.pdf">http://www.chromtech.net.au/pdf2/vici09_125.pdf</a>
DVLSC6MWE2	Replacement valve 6 Position 1/4in fittings 4.0mm ports 4in standoff	Valco Valves	\$1,428.00	<a href="http://www.chromtech.net.au/pdf2/vici09_125.pdf">http://www.chromtech.net.au/pdf2/vici09_125.pdf</a>
DVLSC8MWE2	8 Port Selector LowP Common Outlet w/Std Electric actuator 1/4inx4.0mm ports	Valco Valves	\$1,536.00	<a href="http://www.chromtech.net.au/pdf2/vici09_125.pdf">http://www.chromtech.net.au/pdf2/vici09_125.pdf</a>
SSAVLSC4MWE2	Replacement rotor 4 Position 1/4in fittings 4.0mm ports 4in standoff	Valco Valves	\$243.00	<a href="http://www.chromtech.net.au/pdf2/vici09_125.pdf">http://www.chromtech.net.au/pdf2/vici09_125.pdf</a>
SSAVLSC6MWE2	Replacement rotor 6 Position 1/4in fittings 4.0mm ports 4in standoff	Valco Valves	\$243.00	<a href="http://www.chromtech.net.au/pdf2/vici09_125.pdf">http://www.chromtech.net.au/pdf2/vici09_125.pdf</a>
SSAVLSC8MWE2	Replacement rotor 8 Position 1/4in fittings 4.0mm ports 4in standoff	Valco Valves	\$243.00	<a href="http://www.chromtech.net.au/pdf2/vici09_125.pdf">http://www.chromtech.net.au/pdf2/vici09_125.pdf</a>
2SC10MWE	10 Port Selector LowP Common Outlet Manual 1/16inx1.0mm ports	Valco Valves	\$1,192.00	<a href="http://www.chromtech.net.au/pdf2/vici09_125.pdf">http://www.chromtech.net.au/pdf2/vici09_125.pdf</a>
DSC10MWE	Replacement valve 10 Position 1/8in fittings 1.0mm ports 4in standoff	Valco Valves	\$988.00	<a href="http://www.chromtech.net.au/pdf2/vici09_125.pdf">http://www.chromtech.net.au/pdf2/vici09_125.pdf</a>
2SC12MWE	12 Port Selector LowP Common Outlet Manual 1/16inx1.0mm ports	Valco Valves	\$1,294.00	<a href="http://www.chromtech.net.au/pdf2/vici09_125.pdf">http://www.chromtech.net.au/pdf2/vici09_125.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

EMT2SC16MWE	16 Port Selector LowP Common Outlet w/Microelectric actuator 1/8inx1.0mm ports	Valco Valves	\$2,403.00	<a href="http://www.chromtech.net.au/pdf2/vici09_125.pdf">http://www.chromtech.net.au/pdf2/vici09_125.pdf</a>
EMT2SC6MWE	6 Port Selector LowP Common Outlet w/Microelectric actuator1/8inx1.0mm ports	Valco Valves	\$2,193.00	<a href="http://www.chromtech.net.au/pdf2/vici09_125.pdf">http://www.chromtech.net.au/pdf2/vici09_125.pdf</a>
DSC12MWE	Replacement valve 12 Position 1/8in fittings 1.0mm ports 4in standoff	Valco Valves	\$1,090.00	<a href="http://www.chromtech.net.au/pdf2/vici09_125.pdf">http://www.chromtech.net.au/pdf2/vici09_125.pdf</a>
DSC16MWE	Replacement valve 16 Position 1/8in fittings 1.0mm ports 4in standoff	Valco Valves	\$1,090.00	<a href="http://www.chromtech.net.au/pdf2/vici09_125.pdf">http://www.chromtech.net.au/pdf2/vici09_125.pdf</a>
SSASC6MWE	Replacement rotor 6 Position1/8in fittings 1.0mm ports 4in standoff	Valco Valves	\$134.00	<a href="http://www.chromtech.net.au/pdf2/vici09_125.pdf">http://www.chromtech.net.au/pdf2/vici09_125.pdf</a>
SSASC10MWE	Replacement rotor 10 Position 1/8in fittings 1.0mm ports 4in standoff	Valco Valves	\$134.00	<a href="http://www.chromtech.net.au/pdf2/vici09_125.pdf">http://www.chromtech.net.au/pdf2/vici09_125.pdf</a>
SSASC12MWE	Replacement rotor 12 Position 1/8in fittings 1.0mm ports 4in standoff	Valco Valves	\$134.00	<a href="http://www.chromtech.net.au/pdf2/vici09_125.pdf">http://www.chromtech.net.au/pdf2/vici09_125.pdf</a>
SSASC16MWE	Replacement rotor 16 Position 1/8in fittings 1.0mm ports 4in standoff	Valco Valves	\$147.00	<a href="http://www.chromtech.net.au/pdf2/vici09_125.pdf">http://www.chromtech.net.au/pdf2/vici09_125.pdf</a>
EMT2VLSC4MWE2	4 Port Selector LowP Common Outlet w/Microelectric actuator 1/4inx4.0mm ports	Valco Valves	\$2,672.00	<a href="http://www.chromtech.net.au/pdf2/vici09_125.pdf">http://www.chromtech.net.au/pdf2/vici09_125.pdf</a>
EMT2SC12MWE	12 Port Selector LowP Common Outlet w/Microelectric actuator 1/8inx1.0mm ports	Valco Valves	\$2,403.00	<a href="http://www.chromtech.net.au/pdf2/vici09_125.pdf">http://www.chromtech.net.au/pdf2/vici09_125.pdf</a>
AH2VLSC4MWE2	4 Port Selector LowP Common Outlet w/Air actuator 1/4inx4.0mm ports	Valco Valves	\$2,162.00	<a href="http://www.chromtech.net.au/pdf2/vici09_125.pdf">http://www.chromtech.net.au/pdf2/vici09_125.pdf</a>
EMT2SC10MWE	10 Port Selector LowP Common Outlet w/Microelectric actuator 1/8inx1.0mm ports	Valco Valves	\$2,301.00	<a href="http://www.chromtech.net.au/pdf2/vici09_125.pdf">http://www.chromtech.net.au/pdf2/vici09_125.pdf</a>
E2VLSC8MWE2	8 Port Selector LowP Common Outlet w/Std Electric actuator 1/4inx4.0mm ports	Valco Valves	\$2,570.00	<a href="http://www.chromtech.net.au/pdf2/vici09_125.pdf">http://www.chromtech.net.au/pdf2/vici09_125.pdf</a>
E2VLSC6MWE2	6 Port Selector LowP Common Outlet w/Std Electric actuator 1/4inx4.0mm ports	Valco Valves	\$2,460.00	<a href="http://www.chromtech.net.au/pdf2/vici09_125.pdf">http://www.chromtech.net.au/pdf2/vici09_125.pdf</a>
E2VLSC4MWE2	4 Port Selector LowP Common Outlet w/Std Electric actuator 1/4inx4.0mm ports	Valco Valves	\$2,391.00	<a href="http://www.chromtech.net.au/pdf2/vici09_125.pdf">http://www.chromtech.net.au/pdf2/vici09_125.pdf</a>
AH2VLSC8MWE2	8 Port Selector LowP Common Outlet w/Air actuator 1/4inx4.0mm ports	Valco Valves	\$2,340.00	<a href="http://www.chromtech.net.au/pdf2/vici09_125.pdf">http://www.chromtech.net.au/pdf2/vici09_125.pdf</a>
AH2VLSC6MWE2	6 Port Selector LowP Common Outlet w/Air actuator 1/4inx4.0mm ports	Valco Valves	\$2,232.00	<a href="http://www.chromtech.net.au/pdf2/vici09_125.pdf">http://www.chromtech.net.au/pdf2/vici09_125.pdf</a>
DCSC12MWE	Replacement valve 12 Position 1/16in fittings 1.0mm ports	Valco Valves	\$1,090.00	<a href="http://www.chromtech.net.au/pdf2/vici09_124.pdf">http://www.chromtech.net.au/pdf2/vici09_124.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

EMT2CSC6MWE	1/16inx0.75mm Common Outlet SC, 6 position,w/Microelecric actuator	Valco Valves	\$2,193.00	<a href="http://www.chromtech.net.au/pdf2/vici09_124.pdf">http://www.chromtech.net.au/pdf2/vici09_124.pdf</a>
EMT2CSC10MWE	1/16inx0.75mm Common Outlet SC, 10 position, w/Microelecric actuator	Valco Valves	\$2,301.00	<a href="http://www.chromtech.net.au/pdf2/vici09_124.pdf">http://www.chromtech.net.au/pdf2/vici09_124.pdf</a>
EMT2CSC12MWE	1/16inx0.75mm Common Outlet SC, 12 position, w/Microelecric actuator	Valco Valves	\$2,403.00	<a href="http://www.chromtech.net.au/pdf2/vici09_124.pdf">http://www.chromtech.net.au/pdf2/vici09_124.pdf</a>
SSACSC16MWE	Replacement rotor 16 Position 1/16in fittings 1.0mm ports	Valco Valves	\$147.00	<a href="http://www.chromtech.net.au/pdf2/vici09_124.pdf">http://www.chromtech.net.au/pdf2/vici09_124.pdf</a>
EMT2CSC16MWE	1/16inx0.75mm Common Outlet SC, 16 position, w/Microelecric actuator	Valco Valves	\$2,403.00	<a href="http://www.chromtech.net.au/pdf2/vici09_124.pdf">http://www.chromtech.net.au/pdf2/vici09_124.pdf</a>
SSACSC10MWE	Replacement rotor 10 Position 1/16in fittings 1.0mm ports	Valco Valves	\$134.00	<a href="http://www.chromtech.net.au/pdf2/vici09_124.pdf">http://www.chromtech.net.au/pdf2/vici09_124.pdf</a>
DCSC6MWE	Replacement valve 6 Position 1/16in fittings 1.0mm ports	Valco Valves	\$879.00	<a href="http://www.chromtech.net.au/pdf2/vici09_124.pdf">http://www.chromtech.net.au/pdf2/vici09_124.pdf</a>
SSACSC6MWE	Replacement rotor 6 Position 1/16in fittings 1.0mm ports	Valco Valves	\$134.00	<a href="http://www.chromtech.net.au/pdf2/vici09_124.pdf">http://www.chromtech.net.au/pdf2/vici09_124.pdf</a>
E2CSC16MWE	1/16inx0.75mm Common Outlet SC, 16 position, w/Std Elecric actuator	Valco Valves	\$2,122.00	<a href="http://www.chromtech.net.au/pdf2/vici09_124.pdf">http://www.chromtech.net.au/pdf2/vici09_124.pdf</a>
A2CSC6MWE	1/16inx0.75mm Common Outlet SC, 6 position,w/Air actuator	Valco Valves	\$1,396.00	<a href="http://www.chromtech.net.au/pdf2/vici09_124.pdf">http://www.chromtech.net.au/pdf2/vici09_124.pdf</a>
SSACSC12MWE	Replacement rotor 12 Position 1/16in fittings 1.0mm ports	Valco Valves	\$134.00	<a href="http://www.chromtech.net.au/pdf2/vici09_124.pdf">http://www.chromtech.net.au/pdf2/vici09_124.pdf</a>
E2CSC12MWE	1/16inx0.75mm Common Outlet SC, 12 position, w/Std Elecric actuator	Valco Valves	\$2,122.00	<a href="http://www.chromtech.net.au/pdf2/vici09_124.pdf">http://www.chromtech.net.au/pdf2/vici09_124.pdf</a>
E2CSC10MWE	1/16inx0.75mm Common Outlet SC, 10 position, w/Std Elecric actuator	Valco Valves	\$2,020.00	<a href="http://www.chromtech.net.au/pdf2/vici09_124.pdf">http://www.chromtech.net.au/pdf2/vici09_124.pdf</a>
E2CSC6MWE	1/16inx0.75mm Common Outlet SC, 6 position, w/Std Elecric actuator	Valco Valves	\$1,912.00	<a href="http://www.chromtech.net.au/pdf2/vici09_124.pdf">http://www.chromtech.net.au/pdf2/vici09_124.pdf</a>
A2CSC16MWE	1/16inx0.75mm Common Outlet SC, 16 position, w/Air actuator	Valco Valves	\$1,606.00	<a href="http://www.chromtech.net.au/pdf2/vici09_124.pdf">http://www.chromtech.net.au/pdf2/vici09_124.pdf</a>
A2CSC10MWE	1/16inx0.75mm Common Outlet SC, 10 position, w/Air actuator	Valco Valves	\$1,504.00	<a href="http://www.chromtech.net.au/pdf2/vici09_124.pdf">http://www.chromtech.net.au/pdf2/vici09_124.pdf</a>
2CSC16MWE	1/16inx0.75mm Common Outlet SC, 16 position, Manual	Valco Valves	\$1,294.00	<a href="http://www.chromtech.net.au/pdf2/vici09_124.pdf">http://www.chromtech.net.au/pdf2/vici09_124.pdf</a>
2CSC12MWE	1/16inx0.75mm Common Outlet SC, 12 position, Manual	Valco Valves	\$1,294.00	<a href="http://www.chromtech.net.au/pdf2/vici09_124.pdf">http://www.chromtech.net.au/pdf2/vici09_124.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

2CSC10MWE	1/16inx0.75mm Common Outlet SC, 10 position, Manual	Valco Valves	\$1,192.00	<a href="http://www.chromtech.net.au/pdf2/vici09_124.pdf">http://www.chromtech.net.au/pdf2/vici09_124.pdf</a>
2CSC6MWE	1/16inx0.75mm Common Outlet SC, 6 position, Manual	Valco Valves	\$1,083.00	<a href="http://www.chromtech.net.au/pdf2/vici09_124.pdf">http://www.chromtech.net.au/pdf2/vici09_124.pdf</a>
DCSC10MWE	Replacement valve 10 Position 1/16in fittings 1.0mm ports	Valco Valves	\$988.00	<a href="http://www.chromtech.net.au/pdf2/vici09_124.pdf">http://www.chromtech.net.au/pdf2/vici09_124.pdf</a>
DSC6MWE	Replacement valve 6 Position 1/8in fittings 1.0mm ports 4in standoff	Valco Valves	\$879.00	<a href="http://www.chromtech.net.au/pdf2/vici09_124.pdf">http://www.chromtech.net.au/pdf2/vici09_124.pdf</a>
A2CSC12MWE	1/16inx0.75mm Common Outlet SC, 12 position, w/Air actuator	Valco Valves	\$1,606.00	<a href="http://www.chromtech.net.au/pdf2/vici09_124.pdf">http://www.chromtech.net.au/pdf2/vici09_124.pdf</a>
DCSC16MWE	Replacement valve 16 Position 1/16in fittings 1.0mm ports	Valco Valves	\$1,090.00	<a href="http://www.chromtech.net.au/pdf2/vici09_124.pdf">http://www.chromtech.net.au/pdf2/vici09_124.pdf</a>
AH2VLS6MWE2	6 Port DeadEnd Selector SD w/Air actuator 1/4inx4.0mm ports	Valco Valves	\$2,193.00	<a href="http://www.chromtech.net.au/pdf2/vici09_123.pdf">http://www.chromtech.net.au/pdf2/vici09_123.pdf</a>
DSD10MWE	Replacement valve 10 Position 1/8in fittings 1.0mm ports	Valco Valves	\$950.00	<a href="http://www.chromtech.net.au/pdf2/vici09_123.pdf">http://www.chromtech.net.au/pdf2/vici09_123.pdf</a>
DSD12MWE	Replacement valve 12 Position 1/8in fittings 1.0mm ports	Valco Valves	\$1,052.00	<a href="http://www.chromtech.net.au/pdf2/vici09_123.pdf">http://www.chromtech.net.au/pdf2/vici09_123.pdf</a>
DSD16MWE	Replacement valve 16 Position 1/8in fittings 1.0mm ports	Valco Valves	\$1,052.00	<a href="http://www.chromtech.net.au/pdf2/vici09_123.pdf">http://www.chromtech.net.au/pdf2/vici09_123.pdf</a>
SSASD6MWE	Replacement rotor 6 Position 1/8in fittings 1.0mm ports	Valco Valves	\$134.00	<a href="http://www.chromtech.net.au/pdf2/vici09_123.pdf">http://www.chromtech.net.au/pdf2/vici09_123.pdf</a>
SSASD10MWE	Replacement rotor 10 Position 1/8in fittings 1.0mm ports	Valco Valves	\$134.00	<a href="http://www.chromtech.net.au/pdf2/vici09_123.pdf">http://www.chromtech.net.au/pdf2/vici09_123.pdf</a>
SSASD12MWE	Replacement rotor 12 Position 1/8in fittings 1.0mm ports	Valco Valves	\$134.00	<a href="http://www.chromtech.net.au/pdf2/vici09_123.pdf">http://www.chromtech.net.au/pdf2/vici09_123.pdf</a>
E2VLS10MWE2	10 Port DeadEnd Selector SD w/Std Electric actuator 1/4inx4.0mm ports	Valco Valves	\$2,530.00	<a href="http://www.chromtech.net.au/pdf2/vici09_123.pdf">http://www.chromtech.net.au/pdf2/vici09_123.pdf</a>
AH2VLS4MWE2	4 Port DeadEnd Selector SD w/Air actuator 1/4inx4.0mm ports	Valco Valves	\$2,122.00	<a href="http://www.chromtech.net.au/pdf2/vici09_123.pdf">http://www.chromtech.net.au/pdf2/vici09_123.pdf</a>
EMT2SD12MWE	12 Port DeadEnd Selector SD w/Microelectric actuator 1/8inx1.0mm ports	Valco Valves	\$2,366.00	<a href="http://www.chromtech.net.au/pdf2/vici09_123.pdf">http://www.chromtech.net.au/pdf2/vici09_123.pdf</a>
AH2VLS8MWE2	8 Port DeadEnd Selector SD w/Air actuator 1/4inx4.0mm ports	Valco Valves	\$2,301.00	<a href="http://www.chromtech.net.au/pdf2/vici09_123.pdf">http://www.chromtech.net.au/pdf2/vici09_123.pdf</a>
AH2VLS10MWE2	10 Port DeadEnd Selector SD w/Air actuator 1/4inx4.0mm ports	Valco Valves	\$2,301.00	<a href="http://www.chromtech.net.au/pdf2/vici09_123.pdf">http://www.chromtech.net.au/pdf2/vici09_123.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

E2VLS4MWE2	4 Port DeadEnd Selector SD w/Std Electric actuator 1/4inx4.0mm ports	Valco Valves	\$2,352.00	<a href="http://www.chromtech.net.au/pdf2/vici09_123.pdf">http://www.chromtech.net.au/pdf2/vici09_123.pdf</a>
E2VLS6MWE2	6 Port DeadEnd Selector SD w/Std Electric actuator 1/4inx4.0mm ports	Valco Valves	\$2,422.00	<a href="http://www.chromtech.net.au/pdf2/vici09_123.pdf">http://www.chromtech.net.au/pdf2/vici09_123.pdf</a>
E2VLS8MWE2	8 Port DeadEnd Selector SD w/Std Electric actuator 1/4inx4.0mm ports	Valco Valves	\$2,530.00	<a href="http://www.chromtech.net.au/pdf2/vici09_123.pdf">http://www.chromtech.net.au/pdf2/vici09_123.pdf</a>
SSASD16MWE	Replacement rotor 16 Position 1/8in fittings 1.0mm ports	Valco Valves	\$134.00	<a href="http://www.chromtech.net.au/pdf2/vici09_123.pdf">http://www.chromtech.net.au/pdf2/vici09_123.pdf</a>
2CSD6MWE	1/16inx0.75mm Common Outlet SC, 6 position, Manual	Valco Valves	\$2,001.00	<a href="http://www.chromtech.net.au/pdf2/vici09_123.pdf">http://www.chromtech.net.au/pdf2/vici09_123.pdf</a>
2SD6MWE	6 Port DeadEnd Selector SD Manual 1/8inx1.0mm ports	Valco Valves	\$1,046.00	<a href="http://www.chromtech.net.au/pdf2/vici09_123.pdf">http://www.chromtech.net.au/pdf2/vici09_123.pdf</a>
2SD10MWE	10 Port DeadEnd Selector SD Manual 1/8inx1.0mm ports	Valco Valves	\$1,154.00	<a href="http://www.chromtech.net.au/pdf2/vici09_123.pdf">http://www.chromtech.net.au/pdf2/vici09_123.pdf</a>
2SD12MWE	12 Port DeadEnd Selector SD Manual 1/8inx1.0mm ports	Valco Valves	\$1,256.00	<a href="http://www.chromtech.net.au/pdf2/vici09_123.pdf">http://www.chromtech.net.au/pdf2/vici09_123.pdf</a>
2SD16MWE	16 Port DeadEnd Selector SD Manual 1/8inx1.0mm ports	Valco Valves	\$1,256.00	<a href="http://www.chromtech.net.au/pdf2/vici09_123.pdf">http://www.chromtech.net.au/pdf2/vici09_123.pdf</a>
A2SD6MWE	6 Port DeadEnd Selector SD w/Air actuator 1/8inx1.0mm ports	Valco Valves	\$1,358.00	<a href="http://www.chromtech.net.au/pdf2/vici09_123.pdf">http://www.chromtech.net.au/pdf2/vici09_123.pdf</a>
A2SD10MWE	10 Port DeadEnd Selector SD w/Air actuator 1/8inx1.0mm ports	Valco Valves	\$1,467.00	<a href="http://www.chromtech.net.au/pdf2/vici09_123.pdf">http://www.chromtech.net.au/pdf2/vici09_123.pdf</a>
DSD6MWE	Replacement valve 6 Position 1/8in fittings 1.0mm ports	Valco Valves	\$842.00	<a href="http://www.chromtech.net.au/pdf2/vici09_123.pdf">http://www.chromtech.net.au/pdf2/vici09_123.pdf</a>
A2SD16MWE	16 Port DeadEnd Selector SD w/Air actuator 1/8inx1.0mm ports	Valco Valves	\$1,569.00	<a href="http://www.chromtech.net.au/pdf2/vici09_123.pdf">http://www.chromtech.net.au/pdf2/vici09_123.pdf</a>
EMT2SD16MWE	16 Port DeadEnd Selector SD w/Microelectric actuator 1/8inx1.0mm ports	Valco Valves	\$2,366.00	<a href="http://www.chromtech.net.au/pdf2/vici09_123.pdf">http://www.chromtech.net.au/pdf2/vici09_123.pdf</a>
E2SD10MWE	10 Port DeadEnd Selector SD w/Std Electric actuator 1/8inx1.0mm ports	Valco Valves	\$1,983.00	<a href="http://www.chromtech.net.au/pdf2/vici09_123.pdf">http://www.chromtech.net.au/pdf2/vici09_123.pdf</a>
E2SD12MWE	12 Port DeadEnd Selector SD w/Std Electric actuator 1/8inx1.0mm ports	Valco Valves	\$2,085.00	<a href="http://www.chromtech.net.au/pdf2/vici09_123.pdf">http://www.chromtech.net.au/pdf2/vici09_123.pdf</a>
E2SD16MWE	16 Port DeadEnd Selector SD w/Std Electric actuator 1/8inx1.0mm ports	Valco Valves	\$2,085.00	<a href="http://www.chromtech.net.au/pdf2/vici09_123.pdf">http://www.chromtech.net.au/pdf2/vici09_123.pdf</a>
EMT2SD6MWE	6 Port DeadEnd Selector SD w/Microelectric actuator 1/8inx1.0mm ports	Valco Valves	\$2,154.00	<a href="http://www.chromtech.net.au/pdf2/vici09_123.pdf">http://www.chromtech.net.au/pdf2/vici09_123.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

EMT2SD10MWE	10 Port DeadEnd Selector SD w/Microelectric actuator 1/8inx1.0mm ports	Valco Valves	\$2,264.00	<a href="http://www.chromtech.net.au/pdf2/vici09_123.pdf">http://www.chromtech.net.au/pdf2/vici09_123.pdf</a>
EMT2VLSD8MWE2	8 Port DeadEnd Selector SD w/Microelectric actuator 1/4inx4.0mm ports	Valco Valves	\$2,811.00	<a href="http://www.chromtech.net.au/pdf2/vici09_123.pdf">http://www.chromtech.net.au/pdf2/vici09_123.pdf</a>
A2SD12MWE	12 Port DeadEnd Selector SD w/Air actuator 1/8inx1.0mm ports	Valco Valves	\$1,569.00	<a href="http://www.chromtech.net.au/pdf2/vici09_123.pdf">http://www.chromtech.net.au/pdf2/vici09_123.pdf</a>
SSAVLSD8MWE2	Replacement rotor 8 Position1/4in fittings 4.0mm ports	Valco Valves	\$243.00	<a href="http://www.chromtech.net.au/pdf2/vici09_123.pdf">http://www.chromtech.net.au/pdf2/vici09_123.pdf</a>
EMT2VLSD4MWE2	4 Port DeadEnd Selector SD w/Microelectric actuator 1/4inx4.0mm ports	Valco Valves	\$2,632.00	<a href="http://www.chromtech.net.au/pdf2/vici09_123.pdf">http://www.chromtech.net.au/pdf2/vici09_123.pdf</a>
EMT2VLSD6MWE2	6 Port DeadEnd Selector SD w/Microelectric actuator 1/4inx4.0mm ports	Valco Valves	\$2,703.00	<a href="http://www.chromtech.net.au/pdf2/vici09_123.pdf">http://www.chromtech.net.au/pdf2/vici09_123.pdf</a>
SSAVLSD10MWE2	Replacement rotor 10 Position1/4in fittings 4.0mm ports	Valco Valves	\$243.00	<a href="http://www.chromtech.net.au/pdf2/vici09_123.pdf">http://www.chromtech.net.au/pdf2/vici09_123.pdf</a>
SSAVLSD6MWE2	Replacement rotor 6 Position1/4in fittings 4.0mm ports	Valco Valves	\$243.00	<a href="http://www.chromtech.net.au/pdf2/vici09_123.pdf">http://www.chromtech.net.au/pdf2/vici09_123.pdf</a>
SSAVLSD4MWE2	Replacement rotor 4 Position 1/4in fittings 4.0mm ports	Valco Valves	\$243.00	<a href="http://www.chromtech.net.au/pdf2/vici09_123.pdf">http://www.chromtech.net.au/pdf2/vici09_123.pdf</a>
DVLSD10MWE2	Replacement valve 10 Position 1/4in fittings 4.0mm ports	Valco Valves	\$1,498.00	<a href="http://www.chromtech.net.au/pdf2/vici09_123.pdf">http://www.chromtech.net.au/pdf2/vici09_123.pdf</a>
DVLSD6MWE2	Replacement valve 6 Position 1/4in fittings 4.0mm ports	Valco Valves	\$1,389.00	<a href="http://www.chromtech.net.au/pdf2/vici09_123.pdf">http://www.chromtech.net.au/pdf2/vici09_123.pdf</a>
DVLSD4MWE2	Replacement valve 4 Position 1/4in fittings 4.0mm ports	Valco Valves	\$1,320.00	<a href="http://www.chromtech.net.au/pdf2/vici09_123.pdf">http://www.chromtech.net.au/pdf2/vici09_123.pdf</a>
EMT2VLSD10MWE2	10 Port DeadEnd Selector SD w/Microelectric actuator 1/4inx4.0mm ports	Valco Valves	\$2,811.00	<a href="http://www.chromtech.net.au/pdf2/vici09_123.pdf">http://www.chromtech.net.au/pdf2/vici09_123.pdf</a>
DVLSD8MWE2	Replacement valve 8 Position 1/4in fittings 4.0mm ports	Valco Valves	\$1,498.00	<a href="http://www.chromtech.net.au/pdf2/vici09_123.pdf">http://www.chromtech.net.au/pdf2/vici09_123.pdf</a>
EMT2CSD6MWE	LP 6 position SD Configuration microelectric actuator 1/16" fittings 0.75mm ports	Valco Valves	\$2,154.00	<a href="http://www.chromtech.net.au/pdf2/vici09_122.pdf">http://www.chromtech.net.au/pdf2/vici09_122.pdf</a>
2CSD12MWE	12 Port DeadEnd Selector SD Manual 1/16inx0.75mm	Valco Valves	\$1,256.00	<a href="http://www.chromtech.net.au/pdf2/vici09_122.pdf">http://www.chromtech.net.au/pdf2/vici09_122.pdf</a>
EMT2CSD10MWE	10 position SD Configuration microelectric actuator 1/16" fittings 0.75mm ports	Valco Valves	\$2,264.00	<a href="http://www.chromtech.net.au/pdf2/vici09_122.pdf">http://www.chromtech.net.au/pdf2/vici09_122.pdf</a>
E2CSD10MWE	10 Port DeadEnd Selector SD w/Std Eleelectric actuator 1/16inx0.75mm	Valco Valves	\$1,983.00	<a href="http://www.chromtech.net.au/pdf2/vici09_122.pdf">http://www.chromtech.net.au/pdf2/vici09_122.pdf</a>



## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

2CSD6MWE	6 Port DeadEnd Selector SD Manual 1/16inx0.75mm	Valco Valves	\$1,046.00	<a href="http://www.chromtech.net.au/pdf2/vici09_122.pdf">http://www.chromtech.net.au/pdf2/vici09_122.pdf</a>
2CSD10MWE	10 Port DeadEnd Selector SD Manual 1/16inx0.75mm	Valco Valves	\$1,154.00	<a href="http://www.chromtech.net.au/pdf2/vici09_122.pdf">http://www.chromtech.net.au/pdf2/vici09_122.pdf</a>
2CSD16MWE	16 Port DeadEnd Selector SD Manual 1/16inx0.75mm	Valco Valves	\$1,256.00	<a href="http://www.chromtech.net.au/pdf2/vici09_122.pdf">http://www.chromtech.net.au/pdf2/vici09_122.pdf</a>
A2CSD10MWE	10 Port DeadEnd Selector SD w/Air actuator 1/16inx0.75mm	Valco Valves	\$1,467.00	<a href="http://www.chromtech.net.au/pdf2/vici09_122.pdf">http://www.chromtech.net.au/pdf2/vici09_122.pdf</a>
A2CSD12MWE	12 Port DeadEnd Selector SD w/Air actuator 1/16inx0.75mm	Valco Valves	\$1,569.00	<a href="http://www.chromtech.net.au/pdf2/vici09_122.pdf">http://www.chromtech.net.au/pdf2/vici09_122.pdf</a>
A2CSD16MWE	16 Port DeadEnd Selector SD w/Air actuator 1/16inx0.75mm	Valco Valves	\$1,569.00	<a href="http://www.chromtech.net.au/pdf2/vici09_122.pdf">http://www.chromtech.net.au/pdf2/vici09_122.pdf</a>
A2CSD6MWE	6 Port DeadEnd Selector SD w/Air actuator 1/16inx0.75mm	Valco Valves	\$1,358.00	<a href="http://www.chromtech.net.au/pdf2/vici09_122.pdf">http://www.chromtech.net.au/pdf2/vici09_122.pdf</a>
E2CSD12MWE	12 Port DeadEnd Selector SD w/Std Eleelectric actuator 1/16inx0.75mm	Valco Valves	\$2,085.00	<a href="http://www.chromtech.net.au/pdf2/vici09_122.pdf">http://www.chromtech.net.au/pdf2/vici09_122.pdf</a>
SSACSD12MWE	Replacement rotor 12 position 1/16in fittings 0.75mm ports	Valco Valves	\$134.00	<a href="http://www.chromtech.net.au/pdf2/vici09_122.pdf">http://www.chromtech.net.au/pdf2/vici09_122.pdf</a>
EMT2CSD12MWE	6 Port DeadEnd Selector SD w/Microelectric actuator 1/16inx0.75mm	Valco Valves	\$2,366.00	<a href="http://www.chromtech.net.au/pdf2/vici09_122.pdf">http://www.chromtech.net.au/pdf2/vici09_122.pdf</a>
EMT2CSD16MWE	6 Port DeadEnd Selector SD w/Microelectric actuator 1/16inx0.75mm	Valco Valves	\$2,366.00	<a href="http://www.chromtech.net.au/pdf2/vici09_122.pdf">http://www.chromtech.net.au/pdf2/vici09_122.pdf</a>
DCSD6MWE	Replacement valve 6 position 1/16in fittings 0.75mm ports	Valco Valves	\$842.00	<a href="http://www.chromtech.net.au/pdf2/vici09_122.pdf">http://www.chromtech.net.au/pdf2/vici09_122.pdf</a>
DCSD10MWE	Replacement valve 10 position 1/16in fittings 0.75mm ports	Valco Valves	\$950.00	<a href="http://www.chromtech.net.au/pdf2/vici09_122.pdf">http://www.chromtech.net.au/pdf2/vici09_122.pdf</a>
DCSD12MWE	Replacement valve 12 position 1/16in fittings 0.75mm ports	Valco Valves	\$1,052.00	<a href="http://www.chromtech.net.au/pdf2/vici09_122.pdf">http://www.chromtech.net.au/pdf2/vici09_122.pdf</a>
DCSD16MWE	Replacement valve 16 position 1/16in fittings 0.75mm ports	Valco Valves	\$1,052.00	<a href="http://www.chromtech.net.au/pdf2/vici09_122.pdf">http://www.chromtech.net.au/pdf2/vici09_122.pdf</a>
SSACSD6MWE	Replacement rotor 6 position 1/16in fittings 0.75mm ports	Valco Valves	\$134.00	<a href="http://www.chromtech.net.au/pdf2/vici09_122.pdf">http://www.chromtech.net.au/pdf2/vici09_122.pdf</a>
SSACSD10MWE	Replacement rotor 10 position 1/16in fittings 0.75mm ports	Valco Valves	\$134.00	<a href="http://www.chromtech.net.au/pdf2/vici09_122.pdf">http://www.chromtech.net.au/pdf2/vici09_122.pdf</a>
E2CSD6MWE	6 Port DeadEnd Selector SD w/Std Eleelectric actuator 1/16inx0.75mm	Valco Valves	\$2,001.00	<a href="http://www.chromtech.net.au/pdf2/vici09_122.pdf">http://www.chromtech.net.au/pdf2/vici09_122.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

E2CSD16MWE	16 Port DeadEnd Selector SD w/Std Elelectric actuator 1/16inx0.75mm	Valco Valves	\$2,085.00	<a href="http://www.chromtech.net.au/pdf2/vici09_122.pdf">http://www.chromtech.net.au/pdf2/vici09_122.pdf</a>
ZF.5TI	Titanium 1/32in Metal ferrules	Valco Fittings	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_12.pdf">http://www.chromtech.net.au/pdf2/vici09_12.pdf</a>
ZF.5NI	Nickel 1/32in Metal ferrules	Valco Fittings	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_12.pdf">http://www.chromtech.net.au/pdf2/vici09_12.pdf</a>
ZF.5HC	Hastelloy C 1/32in Metal ferrules	Valco Fittings	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_12.pdf">http://www.chromtech.net.au/pdf2/vici09_12.pdf</a>
ZF4GP-10	Stainless Gold-plated 1/4in Metal ferrules	Valco Fittings	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_12.pdf">http://www.chromtech.net.au/pdf2/vici09_12.pdf</a>
ZF4S6-10	Stainless Type 316 1/4in Metal ferrules	Valco Fittings	\$24.00	<a href="http://www.chromtech.net.au/pdf2/vici09_12.pdf">http://www.chromtech.net.au/pdf2/vici09_12.pdf</a>
ZF2GP-10	Stainless Gold-plated 1/8in Metal ferrules	Valco Fittings	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_12.pdf">http://www.chromtech.net.au/pdf2/vici09_12.pdf</a>
ZF.5GP-10	Stainless Gold-plated 1/32in Metal ferrules	Valco Fittings	\$54.00	<a href="http://www.chromtech.net.au/pdf2/vici09_12.pdf">http://www.chromtech.net.au/pdf2/vici09_12.pdf</a>
ZF.5-10	Stainless Type 303 1/32in Metal ferrules	Valco Fittings	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_12.pdf">http://www.chromtech.net.au/pdf2/vici09_12.pdf</a>
ZF.5S6-10	Stainless Type 316 1/32in Metal ferrules	Valco Fittings	\$51.00	<a href="http://www.chromtech.net.au/pdf2/vici09_12.pdf">http://www.chromtech.net.au/pdf2/vici09_12.pdf</a>
ZF1HC	Hastelloy C 1/16in Metal ferrules	Valco Fittings	\$10.00	<a href="http://www.chromtech.net.au/pdf2/vici09_12.pdf">http://www.chromtech.net.au/pdf2/vici09_12.pdf</a>
ZF1-10	1/16" SS Ferrule	Valco Fittings	\$26.00	<a href="http://www.chromtech.net.au/pdf2/vici09_12.pdf">http://www.chromtech.net.au/pdf2/vici09_12.pdf</a>
ZF1S6-10	1/16in Valco SS ferrules	Valco Fittings	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_12.pdf">http://www.chromtech.net.au/pdf2/vici09_12.pdf</a>
ZF1GP-10	Stainless Gold-plated 1/16in Metal ferrules	Valco Fittings	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_12.pdf">http://www.chromtech.net.au/pdf2/vici09_12.pdf</a>
ZF2S6-10	Stainless Type 316 1/8in Metal ferrules	Valco Fittings	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_12.pdf">http://www.chromtech.net.au/pdf2/vici09_12.pdf</a>
ZF1NI	Nickel 1/16in Metal ferrules	Valco Fittings	\$10.00	<a href="http://www.chromtech.net.au/pdf2/vici09_12.pdf">http://www.chromtech.net.au/pdf2/vici09_12.pdf</a>
ZF1TI	Titanium 1/16in Metal ferrules	Valco Fittings	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_12.pdf">http://www.chromtech.net.au/pdf2/vici09_12.pdf</a>
ZF2HC	Hastelloy C 1/8in Metal ferrules	Valco Fittings	\$10.00	<a href="http://www.chromtech.net.au/pdf2/vici09_12.pdf">http://www.chromtech.net.au/pdf2/vici09_12.pdf</a>
ZF2NI	Nickel 1/8in Metal ferrules	Valco Fittings	\$10.00	<a href="http://www.chromtech.net.au/pdf2/vici09_12.pdf">http://www.chromtech.net.au/pdf2/vici09_12.pdf</a>
ZF2TI	Titanium 1/8in Metal ferrules	Valco Fittings	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_12.pdf">http://www.chromtech.net.au/pdf2/vici09_12.pdf</a>
ZF4NI	Nickel 1/4in Metal ferrules	Valco Fittings	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_12.pdf">http://www.chromtech.net.au/pdf2/vici09_12.pdf</a>
ZF4TI	Titanium 1/4in Metal ferrules	Valco Fittings	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_12.pdf">http://www.chromtech.net.au/pdf2/vici09_12.pdf</a>
ZF.5b-10	Brass 1/32in Metal ferrules	Valco Fittings	\$22.00	<a href="http://www.chromtech.net.au/pdf2/vici09_12.pdf">http://www.chromtech.net.au/pdf2/vici09_12.pdf</a>
ZF1b-10	Brass 1/16in Metal ferrules	Valco Fittings	\$16.00	<a href="http://www.chromtech.net.au/pdf2/vici09_12.pdf">http://www.chromtech.net.au/pdf2/vici09_12.pdf</a>
ZF2b-10	Brass 1/8in Metal ferrules	Valco Fittings	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_12.pdf">http://www.chromtech.net.au/pdf2/vici09_12.pdf</a>
ZF4b-10	Brass 1/4in Metal ferrules	Valco Fittings	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_12.pdf">http://www.chromtech.net.au/pdf2/vici09_12.pdf</a>
ZF2-10	Stainless Type 303 1/8in Metal ferrules	Valco Fittings	\$22.00	<a href="http://www.chromtech.net.au/pdf2/vici09_12.pdf">http://www.chromtech.net.au/pdf2/vici09_12.pdf</a>
ZF4HC	Hastelloy C 1/4in Metal ferrules	Valco Fittings	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_12.pdf">http://www.chromtech.net.au/pdf2/vici09_12.pdf</a>
SL1KUW	1ml 1/8in Stainless steel loops	Valco Valves	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_116.pdf">http://www.chromtech.net.au/pdf2/vici09_116.pdf</a>
SSAL8UW	Replacement rotor 8 port 1/8in fittings large bore	Valco Valves	\$134.00	<a href="http://www.chromtech.net.au/pdf2/vici09_116.pdf">http://www.chromtech.net.au/pdf2/vici09_116.pdf</a>
DL4UW	Replacement valve 4 port 1/8in fittings large bore	Valco Valves	\$734.00	<a href="http://www.chromtech.net.au/pdf2/vici09_116.pdf">http://www.chromtech.net.au/pdf2/vici09_116.pdf</a>
DL6UW	Replacement valve 6 port 1/8in fittings large bore	Valco Valves	\$804.00	<a href="http://www.chromtech.net.au/pdf2/vici09_116.pdf">http://www.chromtech.net.au/pdf2/vici09_116.pdf</a>
DL8UW	Replacement valve 8 port 1/8in fittings large bore	Valco Valves	\$873.00	<a href="http://www.chromtech.net.au/pdf2/vici09_116.pdf">http://www.chromtech.net.au/pdf2/vici09_116.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

DL10UW	Replacement valve 10 port 1/8in fittings large bore	Valco Valves	\$873.00	<a href="http://www.chromtech.net.au/pdf2/vici09_116.pdf">http://www.chromtech.net.au/pdf2/vici09_116.pdf</a>
SSAL4UW	Replacement rotor 4 port 1/8in fittings large bore	Valco Valves	\$134.00	<a href="http://www.chromtech.net.au/pdf2/vici09_116.pdf">http://www.chromtech.net.au/pdf2/vici09_116.pdf</a>
SSAL6UW	Replacement rotor 6 port 1/8in fittings large bore	Valco Valves	\$134.00	<a href="http://www.chromtech.net.au/pdf2/vici09_116.pdf">http://www.chromtech.net.au/pdf2/vici09_116.pdf</a>
SSAL10UW	Replacement rotor 10 port 1/8in fittings large bore	Valco Valves	\$134.00	<a href="http://www.chromtech.net.au/pdf2/vici09_116.pdf">http://www.chromtech.net.au/pdf2/vici09_116.pdf</a>
SL2KUW	2ml 1/8in Stainless steel loops	Valco Valves	\$42.00	<a href="http://www.chromtech.net.au/pdf2/vici09_116.pdf">http://www.chromtech.net.au/pdf2/vici09_116.pdf</a>
SL5KUW	5ml 1/8in Stainless steel loops	Valco Valves	\$51.00	<a href="http://www.chromtech.net.au/pdf2/vici09_116.pdf">http://www.chromtech.net.au/pdf2/vici09_116.pdf</a>
EDL10UW	10 Port w/ Microelectric actuator 1/8in fittings large bore	Valco Valves	\$2,020.00	<a href="http://www.chromtech.net.au/pdf2/vici09_116.pdf">http://www.chromtech.net.au/pdf2/vici09_116.pdf</a>
SL10KUW	1ml Stainless steel loops	Valco Valves	\$72.00	<a href="http://www.chromtech.net.au/pdf2/vici09_116.pdf">http://www.chromtech.net.au/pdf2/vici09_116.pdf</a>
SL100UW	100ul Stainless steel loops	Valco Valves	\$34.00	<a href="http://www.chromtech.net.au/pdf2/vici09_116.pdf">http://www.chromtech.net.au/pdf2/vici09_116.pdf</a>
SL20KUW	20ml 1/8in Stainless steel loops	Valco Valves	\$128.00	<a href="http://www.chromtech.net.au/pdf2/vici09_116.pdf">http://www.chromtech.net.au/pdf2/vici09_116.pdf</a>
SL500UW	500ul 1/8in Stainless steel loops	Valco Valves	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_116.pdf">http://www.chromtech.net.au/pdf2/vici09_116.pdf</a>
2L10UW	10 Port Manual 1/8in fittings large bore 1.0mm	Valco Valves	\$988.00	<a href="http://www.chromtech.net.au/pdf2/vici09_116.pdf">http://www.chromtech.net.au/pdf2/vici09_116.pdf</a>
SL250UW	250ul Stainless steel loops	Valco Valves	\$34.00	<a href="http://www.chromtech.net.au/pdf2/vici09_116.pdf">http://www.chromtech.net.au/pdf2/vici09_116.pdf</a>
EDL8UW	8 Port w/ Microelectric actuator 1/8in fittings large bore	Valco Valves	\$2,020.00	<a href="http://www.chromtech.net.au/pdf2/vici09_116.pdf">http://www.chromtech.net.au/pdf2/vici09_116.pdf</a>
L6UW	Manual 6 port 1/8in fittings large bore	Valco Valves	\$861.00	<a href="http://www.chromtech.net.au/pdf2/vici09_116.pdf">http://www.chromtech.net.au/pdf2/vici09_116.pdf</a>
L8UW	Manual 8 port 1/8in fittings large bore	Valco Valves	\$930.00	<a href="http://www.chromtech.net.au/pdf2/vici09_116.pdf">http://www.chromtech.net.au/pdf2/vici09_116.pdf</a>
L4UW	Manual 4 port 1/8in fittings large bore	Valco Valves	\$790.00	<a href="http://www.chromtech.net.au/pdf2/vici09_116.pdf">http://www.chromtech.net.au/pdf2/vici09_116.pdf</a>
AL4UW	Air actuated 4 port 1/8in fittings large bore	Valco Valves	\$1,052.00	<a href="http://www.chromtech.net.au/pdf2/vici09_116.pdf">http://www.chromtech.net.au/pdf2/vici09_116.pdf</a>
AL6UW	Air actuated 6 port 1/8in fittings large bore	Valco Valves	\$1,122.00	<a href="http://www.chromtech.net.au/pdf2/vici09_116.pdf">http://www.chromtech.net.au/pdf2/vici09_116.pdf</a>
EL10UW	Electric actuated 10 port 1/8in fittings large bore	Valco Valves	\$1,600.00	<a href="http://www.chromtech.net.au/pdf2/vici09_116.pdf">http://www.chromtech.net.au/pdf2/vici09_116.pdf</a>
AL10UW	Air actuated 10 port 1/8in fittings large bore	Valco Valves	\$1,192.00	<a href="http://www.chromtech.net.au/pdf2/vici09_116.pdf">http://www.chromtech.net.au/pdf2/vici09_116.pdf</a>
EL4UW	Electric actuated 4 port 1/8in fittings large bore	Valco Valves	\$1,460.00	<a href="http://www.chromtech.net.au/pdf2/vici09_116.pdf">http://www.chromtech.net.au/pdf2/vici09_116.pdf</a>
EL6UW	Electric actuated 6 port 1/8in fittings large bore	Valco Valves	\$1,530.00	<a href="http://www.chromtech.net.au/pdf2/vici09_116.pdf">http://www.chromtech.net.au/pdf2/vici09_116.pdf</a>
EL8UW	Electric actuated 8 port 1/8in fittings large bore	Valco Valves	\$1,600.00	<a href="http://www.chromtech.net.au/pdf2/vici09_116.pdf">http://www.chromtech.net.au/pdf2/vici09_116.pdf</a>
AL8UW	Air actuated port 1/8in fittings large bore	Valco Valves	\$1,192.00	<a href="http://www.chromtech.net.au/pdf2/vici09_116.pdf">http://www.chromtech.net.au/pdf2/vici09_116.pdf</a>
EDL6UW	6 Port w/ Microelectric actuator 1/8in fittings large bore	Valco Valves	\$1,950.00	<a href="http://www.chromtech.net.au/pdf2/vici09_116.pdf">http://www.chromtech.net.au/pdf2/vici09_116.pdf</a>
EDL4UW	4 Port w/ Microelectric actuator 1/8in fittings large bore	Valco Valves	\$1,754.00	<a href="http://www.chromtech.net.au/pdf2/vici09_116.pdf">http://www.chromtech.net.au/pdf2/vici09_116.pdf</a>
D6UW	Replacement valve 6 port 1/8in fittings 0.75mm	Valco Valves	\$734.00	<a href="http://www.chromtech.net.au/pdf2/vici09_115.pdf">http://www.chromtech.net.au/pdf2/vici09_115.pdf</a>
SL250UW	250ul Stainless steel loops	Valco Valves	\$34.00	<a href="http://www.chromtech.net.au/pdf2/vici09_115.pdf">http://www.chromtech.net.au/pdf2/vici09_115.pdf</a>
ED10UW	Microelectric actuated 10 port 1/8in fittings 0.75mm	Valco Valves	\$1,950.00	<a href="http://www.chromtech.net.au/pdf2/vici09_115.pdf">http://www.chromtech.net.au/pdf2/vici09_115.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

D10UW	Replacement valve 10 port 1/8in fittings 0.75mm	Valco Valves	\$804.00	<a href="http://www.chromtech.net.au/pdf2/vici09_115.pdf">http://www.chromtech.net.au/pdf2/vici09_115.pdf</a>
SSA4UW	Replacement rotor 4 port 1/8in fittings 0.75mm	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_115.pdf">http://www.chromtech.net.au/pdf2/vici09_115.pdf</a>
SSA6UW	Replacement rotor 6 port 1/8in fittings 0.75mm	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_115.pdf">http://www.chromtech.net.au/pdf2/vici09_115.pdf</a>
D8UW	Replacement valve 8 port 1/8in fittings 0.75mm	Valco Valves	\$804.00	<a href="http://www.chromtech.net.au/pdf2/vici09_115.pdf">http://www.chromtech.net.au/pdf2/vici09_115.pdf</a>
SSA8UW	Replacement rotor 8 port 1/8in fittings 0.75mm	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_115.pdf">http://www.chromtech.net.au/pdf2/vici09_115.pdf</a>
SL15UW	15ul 1/8in Stainless steel loops	Valco Valves	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_115.pdf">http://www.chromtech.net.au/pdf2/vici09_115.pdf</a>
SL10UW	10ul Stainless steel loops	Valco Valves	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_115.pdf">http://www.chromtech.net.au/pdf2/vici09_115.pdf</a>
ED6UW	Microelectric actuated 6 port 1/8in fittings 0.75mm	Valco Valves	\$1,754.00	<a href="http://www.chromtech.net.au/pdf2/vici09_115.pdf">http://www.chromtech.net.au/pdf2/vici09_115.pdf</a>
210UW	Manual 10 port 1/8in fittings 0.75mm	Valco Valves	\$918.00	<a href="http://www.chromtech.net.au/pdf2/vici09_115.pdf">http://www.chromtech.net.au/pdf2/vici09_115.pdf</a>
SL500UW	500ul 1/8in Stainless steel loops	Valco Valves	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_115.pdf">http://www.chromtech.net.au/pdf2/vici09_115.pdf</a>
SSA10UW	Replacement rotor 10 port 1/8in fittings 0.75mm	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_115.pdf">http://www.chromtech.net.au/pdf2/vici09_115.pdf</a>
ED4UW	Microelectric actuated 4 port 1/8in fittings 0.75mm	Valco Valves	\$1,683.00	<a href="http://www.chromtech.net.au/pdf2/vici09_115.pdf">http://www.chromtech.net.au/pdf2/vici09_115.pdf</a>
E10UW	Electric actuated 10 port 1/8in fittings 0.75mm	Valco Valves	\$1,530.00	<a href="http://www.chromtech.net.au/pdf2/vici09_115.pdf">http://www.chromtech.net.au/pdf2/vici09_115.pdf</a>
E8UW	Electric actuated 8 port 1/8in fittings 0.75mm	Valco Valves	\$1,530.00	<a href="http://www.chromtech.net.au/pdf2/vici09_115.pdf">http://www.chromtech.net.au/pdf2/vici09_115.pdf</a>
E6UW	Electric actuated 6 port 1/8in fittings 0.75mm	Valco Valves	\$1,460.00	<a href="http://www.chromtech.net.au/pdf2/vici09_115.pdf">http://www.chromtech.net.au/pdf2/vici09_115.pdf</a>
E4UW	Electric actuated 4 port 1/8in fittings 0.75mm	Valco Valves	\$1,389.00	<a href="http://www.chromtech.net.au/pdf2/vici09_115.pdf">http://www.chromtech.net.au/pdf2/vici09_115.pdf</a>
A10UW	Air actuated 10 port 1/8in fittings 0.75mm	Valco Valves	\$1,122.00	<a href="http://www.chromtech.net.au/pdf2/vici09_115.pdf">http://www.chromtech.net.au/pdf2/vici09_115.pdf</a>
A8UW	Air actuated port 1/8in fittings 0.75mm	Valco Valves	\$1,122.00	<a href="http://www.chromtech.net.au/pdf2/vici09_115.pdf">http://www.chromtech.net.au/pdf2/vici09_115.pdf</a>
A4UW	Air actuated 4 port 1/8in fittings 0.75mm	Valco Valves	\$981.00	<a href="http://www.chromtech.net.au/pdf2/vici09_115.pdf">http://www.chromtech.net.au/pdf2/vici09_115.pdf</a>
8UW	Manual 8 port 1/8in fittings 0.75mm	Valco Valves	\$861.00	<a href="http://www.chromtech.net.au/pdf2/vici09_115.pdf">http://www.chromtech.net.au/pdf2/vici09_115.pdf</a>
6UW	Manual 6 port 1/8in fittings 0.75mm	Valco Valves	\$790.00	<a href="http://www.chromtech.net.au/pdf2/vici09_115.pdf">http://www.chromtech.net.au/pdf2/vici09_115.pdf</a>
4UW	Manual 4 port 1/8in fittings 0.75mm	Valco Valves	\$720.00	<a href="http://www.chromtech.net.au/pdf2/vici09_115.pdf">http://www.chromtech.net.au/pdf2/vici09_115.pdf</a>
SL20UW	20ul 1/8in Stainless steel loops	Valco Valves	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_115.pdf">http://www.chromtech.net.au/pdf2/vici09_115.pdf</a>
D4UW	Replacement valve 4 port 1/8in fittings 0.75mm	Valco Valves	\$663.00	<a href="http://www.chromtech.net.au/pdf2/vici09_115.pdf">http://www.chromtech.net.au/pdf2/vici09_115.pdf</a>
A6UW	Air actuated 6 port 1/8in fittings 0.75mm	Valco Valves	\$1,052.00	<a href="http://www.chromtech.net.au/pdf2/vici09_115.pdf">http://www.chromtech.net.au/pdf2/vici09_115.pdf</a>
SL20KUW	20ml 1/8in Stainless steel loops	Valco Valves	\$128.00	<a href="http://www.chromtech.net.au/pdf2/vici09_115.pdf">http://www.chromtech.net.au/pdf2/vici09_115.pdf</a>
ED8UW	Microelectric actuated 8 port 1/8in fittings 0.75mm	Valco Valves	\$1,950.00	<a href="http://www.chromtech.net.au/pdf2/vici09_115.pdf">http://www.chromtech.net.au/pdf2/vici09_115.pdf</a>
SL10KUW	1ml Stainless steel loops	Valco Valves	\$72.00	<a href="http://www.chromtech.net.au/pdf2/vici09_115.pdf">http://www.chromtech.net.au/pdf2/vici09_115.pdf</a>
SL100UW	100ul Stainless steel loops	Valco Valves	\$34.00	<a href="http://www.chromtech.net.au/pdf2/vici09_115.pdf">http://www.chromtech.net.au/pdf2/vici09_115.pdf</a>
SL5KUW	5ml 1/8in Stainless steel loops	Valco Valves	\$51.00	<a href="http://www.chromtech.net.au/pdf2/vici09_115.pdf">http://www.chromtech.net.au/pdf2/vici09_115.pdf</a>
SL50UW	50ul Stainless steel loops	Valco Valves	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_115.pdf">http://www.chromtech.net.au/pdf2/vici09_115.pdf</a>
SL2KUW	2ml 1/8in Stainless steel loops	Valco Valves	\$42.00	<a href="http://www.chromtech.net.au/pdf2/vici09_115.pdf">http://www.chromtech.net.au/pdf2/vici09_115.pdf</a>
SL25UW	25ul 1/8in Stainless steel loops	Valco Valves	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_115.pdf">http://www.chromtech.net.au/pdf2/vici09_115.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

SL1KUW	1ml 1/8in Stainless steel loops	Valco Valves	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_115.pdf">http://www.chromtech.net.au/pdf2/vici09_115.pdf</a>
SL500CUW	500ul Stainless steel loops	Valco Valves	\$40.00	<a href="http://www.chromtech.net.au/pdf2/vici09_114.pdf">http://www.chromtech.net.au/pdf2/vici09_114.pdf</a>
DC8UW	Replacement valve 8 Port 1/16in fittings 0.75mm ports	Valco Valves	\$804.00	<a href="http://www.chromtech.net.au/pdf2/vici09_114.pdf">http://www.chromtech.net.au/pdf2/vici09_114.pdf</a>
DC10UW	Replacement valve 10 Port 1/16in fittings 0.75mm ports	Valco Valves	\$804.00	<a href="http://www.chromtech.net.au/pdf2/vici09_114.pdf">http://www.chromtech.net.au/pdf2/vici09_114.pdf</a>
SSAC10UW	Replacement rotor 10 Port 1/16in fittings 0.75mm ports	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_114.pdf">http://www.chromtech.net.au/pdf2/vici09_114.pdf</a>
SL3CUW	3ml Stainless steel loops	Valco Valves	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_114.pdf">http://www.chromtech.net.au/pdf2/vici09_114.pdf</a>
SL100CUW	100ul Stainless steel loops	Valco Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_114.pdf">http://www.chromtech.net.au/pdf2/vici09_114.pdf</a>
SL5CUW	5ul 1/16in stainless steel loop	Valco Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_114.pdf">http://www.chromtech.net.au/pdf2/vici09_114.pdf</a>
DC6UW	Replacement valve 6 Port 1/16in fittings 0.75mm ports	Valco Valves	\$734.00	<a href="http://www.chromtech.net.au/pdf2/vici09_114.pdf">http://www.chromtech.net.au/pdf2/vici09_114.pdf</a>
SL10CUW	10ul Stainless steel loops	Valco Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_114.pdf">http://www.chromtech.net.au/pdf2/vici09_114.pdf</a>
SL20CUW	20ul Stainless steel loops	Valco Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_114.pdf">http://www.chromtech.net.au/pdf2/vici09_114.pdf</a>
SL15CUW	15ul Stainless steel loops	Valco Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_114.pdf">http://www.chromtech.net.au/pdf2/vici09_114.pdf</a>
SL1KCUW	1ml 1/16in stainless steel loop	Valco Valves	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_114.pdf">http://www.chromtech.net.au/pdf2/vici09_114.pdf</a>
SL2KCUW	2ml Stainless steel loops	Valco Valves	\$57.00	<a href="http://www.chromtech.net.au/pdf2/vici09_114.pdf">http://www.chromtech.net.au/pdf2/vici09_114.pdf</a>
SL5KCUW	5ml Stainless steel loops	Valco Valves	\$68.00	<a href="http://www.chromtech.net.au/pdf2/vici09_114.pdf">http://www.chromtech.net.au/pdf2/vici09_114.pdf</a>
SL50CUW	50ul Stainless steel loops	Valco Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_114.pdf">http://www.chromtech.net.au/pdf2/vici09_114.pdf</a>
SL10KCUW	10ml Stainless steel loops	Valco Valves	\$87.00	<a href="http://www.chromtech.net.au/pdf2/vici09_114.pdf">http://www.chromtech.net.au/pdf2/vici09_114.pdf</a>
SL250CUW	250ul Stainless steel loops	Valco Valves	\$34.00	<a href="http://www.chromtech.net.au/pdf2/vici09_114.pdf">http://www.chromtech.net.au/pdf2/vici09_114.pdf</a>
C10UW	Manual 10 Port 1/16in fittings 0.75mm ports	Valco Valves	\$861.00	<a href="http://www.chromtech.net.au/pdf2/vici09_114.pdf">http://www.chromtech.net.au/pdf2/vici09_114.pdf</a>
SL25CUW	25ul 1/16in stainless steel loop	Valco Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_114.pdf">http://www.chromtech.net.au/pdf2/vici09_114.pdf</a>
DC4UW	Replacement valve 4 Port 1/16in fittings 0.75mm ports	Valco Valves	\$663.00	<a href="http://www.chromtech.net.au/pdf2/vici09_114.pdf">http://www.chromtech.net.au/pdf2/vici09_114.pdf</a>
C6UW	Manual 6 Port 1/16in fittings 0.75mm ports	Valco Valves	\$790.00	<a href="http://www.chromtech.net.au/pdf2/vici09_114.pdf">http://www.chromtech.net.au/pdf2/vici09_114.pdf</a>
C8UW	Manual 8 Port 1/16in fittings 0.75mm ports	Valco Valves	\$861.00	<a href="http://www.chromtech.net.au/pdf2/vici09_114.pdf">http://www.chromtech.net.au/pdf2/vici09_114.pdf</a>
AC4UW	Air actuated 4 Port 1/16in fittings 0.75mm ports	Valco Valves	\$981.00	<a href="http://www.chromtech.net.au/pdf2/vici09_114.pdf">http://www.chromtech.net.au/pdf2/vici09_114.pdf</a>
AC6UW	Air actuated 6 Port 1/16in fittings 0.75mm ports	Valco Valves	\$1,052.00	<a href="http://www.chromtech.net.au/pdf2/vici09_114.pdf">http://www.chromtech.net.au/pdf2/vici09_114.pdf</a>
AC8UW	Air actuated 8 Port 1/16in fittings 0.75mm ports	Valco Valves	\$1,122.00	<a href="http://www.chromtech.net.au/pdf2/vici09_114.pdf">http://www.chromtech.net.au/pdf2/vici09_114.pdf</a>
AC10UW	Air actuated 10 Port 1/16in fittings 0.75mm ports	Valco Valves	\$1,122.00	<a href="http://www.chromtech.net.au/pdf2/vici09_114.pdf">http://www.chromtech.net.au/pdf2/vici09_114.pdf</a>
EDC4UW	Electric Actuated 4 Port 1/16in fittings 0.75mm ports	Valco Valves	\$1,683.00	<a href="http://www.chromtech.net.au/pdf2/vici09_114.pdf">http://www.chromtech.net.au/pdf2/vici09_114.pdf</a>
EDC8UW	Electric Actuated 8 Port 1/16in fittings 0.75mm ports	Valco Valves	\$1,950.00	<a href="http://www.chromtech.net.au/pdf2/vici09_114.pdf">http://www.chromtech.net.au/pdf2/vici09_114.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

EC4UW	Electric Actuated 4 Port 1/16in fittings 0.75mm ports	Valco Valves	\$1,389.00	<a href="http://www.chromtech.net.au/pdf2/vici09_114.pdf">http://www.chromtech.net.au/pdf2/vici09_114.pdf</a>
EDC10UW	Electric Actuated 10 Port 1/16in fittings 0.75mm ports	Valco Valves	\$1,950.00	<a href="http://www.chromtech.net.au/pdf2/vici09_114.pdf">http://www.chromtech.net.au/pdf2/vici09_114.pdf</a>
EDC6UW	Electric Actuated 6 Port 1/16in fittings 0.75mm ports	Valco Valves	\$1,754.00	<a href="http://www.chromtech.net.au/pdf2/vici09_114.pdf">http://www.chromtech.net.au/pdf2/vici09_114.pdf</a>
C4UW	Manual 4 Port 1/16in fittings 0.75mm ports	Valco Valves	\$720.00	<a href="http://www.chromtech.net.au/pdf2/vici09_114.pdf">http://www.chromtech.net.au/pdf2/vici09_114.pdf</a>
EC10UW	Electric Actuated 10 Port 1/16in fittings 0.75mm ports	Valco Valves	\$1,530.00	<a href="http://www.chromtech.net.au/pdf2/vici09_114.pdf">http://www.chromtech.net.au/pdf2/vici09_114.pdf</a>
EC8UW	Electric Actuated 8 Port 1/16in fittings 0.75mm ports	Valco Valves	\$1,530.00	<a href="http://www.chromtech.net.au/pdf2/vici09_114.pdf">http://www.chromtech.net.au/pdf2/vici09_114.pdf</a>
EC6UW	Electric Actuated 6 Port 1/16in fittings 0.75mm ports	Valco Valves	\$1,460.00	<a href="http://www.chromtech.net.au/pdf2/vici09_114.pdf">http://www.chromtech.net.au/pdf2/vici09_114.pdf</a>
EC8W	Electric Actuated 8 ports 1/16in fittings 0.40mm ports	Valco Valves	\$1,536.00	<a href="http://www.chromtech.net.au/pdf2/vici09_113.pdf">http://www.chromtech.net.au/pdf2/vici09_113.pdf</a>
C10W	Manual 10 ports 1/16in fittings 0.40mm ports	Valco Valves	\$867.00	<a href="http://www.chromtech.net.au/pdf2/vici09_113.pdf">http://www.chromtech.net.au/pdf2/vici09_113.pdf</a>
SL25CW	25ul Stainless steel loops	Valco Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_113.pdf">http://www.chromtech.net.au/pdf2/vici09_113.pdf</a>
SL5KCW	5ml Stainless steel loops	Valco Valves	\$68.00	<a href="http://www.chromtech.net.au/pdf2/vici09_113.pdf">http://www.chromtech.net.au/pdf2/vici09_113.pdf</a>
C4W	Manual 4 ports 1/16in fittings 0.40mm ports	Valco Valves	\$726.00	<a href="http://www.chromtech.net.au/pdf2/vici09_113.pdf">http://www.chromtech.net.au/pdf2/vici09_113.pdf</a>
C6W	Manual 6 ports 1/16in fittings 0.40mm ports	Valco Valves	\$796.00	<a href="http://www.chromtech.net.au/pdf2/vici09_113.pdf">http://www.chromtech.net.au/pdf2/vici09_113.pdf</a>
SL2KCW	2ml 1/16in Stainless sample loops	Valco Valves	\$57.00	<a href="http://www.chromtech.net.au/pdf2/vici09_113.pdf">http://www.chromtech.net.au/pdf2/vici09_113.pdf</a>
C8W	Manual 8 ports 1/16in fittings 0.40mm ports	Valco Valves	\$867.00	<a href="http://www.chromtech.net.au/pdf2/vici09_113.pdf">http://www.chromtech.net.au/pdf2/vici09_113.pdf</a>
AC4W	Air Actuated 4 ports 1/16in fittings 0.40mm ports	Valco Valves	\$988.00	<a href="http://www.chromtech.net.au/pdf2/vici09_113.pdf">http://www.chromtech.net.au/pdf2/vici09_113.pdf</a>
AC6W	Air Actuated 6 ports 1/16in fittings 0.40mm ports	Valco Valves	\$1,059.00	<a href="http://www.chromtech.net.au/pdf2/vici09_113.pdf">http://www.chromtech.net.au/pdf2/vici09_113.pdf</a>
AC8W	Air Actuated 8 ports 1/16in fittings 0.40mm ports	Valco Valves	\$1,128.00	<a href="http://www.chromtech.net.au/pdf2/vici09_113.pdf">http://www.chromtech.net.au/pdf2/vici09_113.pdf</a>
AC10W	Air Actuated 10 ports 1/16in fittings 0.40mm ports	Valco Valves	\$1,128.00	<a href="http://www.chromtech.net.au/pdf2/vici09_113.pdf">http://www.chromtech.net.au/pdf2/vici09_113.pdf</a>
SL20CW	20ul Stainless steel loops	Valco Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_113.pdf">http://www.chromtech.net.au/pdf2/vici09_113.pdf</a>
EC6W	Electric Actuated 6 ports 1/16in fittings 0.40mm ports	Valco Valves	\$1,467.00	<a href="http://www.chromtech.net.au/pdf2/vici09_113.pdf">http://www.chromtech.net.au/pdf2/vici09_113.pdf</a>
SL5CW	5ul Stainless steel loops	Valco Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_113.pdf">http://www.chromtech.net.au/pdf2/vici09_113.pdf</a>
EC10W	Electric Actuated 10 ports 1/16in fittings 0.40mm ports	Valco Valves	\$1,536.00	<a href="http://www.chromtech.net.au/pdf2/vici09_113.pdf">http://www.chromtech.net.au/pdf2/vici09_113.pdf</a>
EPC4W	Microelectric Actuated 4 ports 1/16in fittings 0.40mm ports	Valco Valves	\$1,612.00	<a href="http://www.chromtech.net.au/pdf2/vici09_113.pdf">http://www.chromtech.net.au/pdf2/vici09_113.pdf</a>
EPC6W	Microelectric Actuated 6 ports 1/16in fittings 0.40mm ports	Valco Valves	\$1,683.00	<a href="http://www.chromtech.net.au/pdf2/vici09_113.pdf">http://www.chromtech.net.au/pdf2/vici09_113.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

EC4W	Electric Actuated 4 ports 1/16in fittings 0.40mm ports	Valco Valves	\$1,396.00	<a href="http://www.chromtech.net.au/pdf2/vici09_113.pdf">http://www.chromtech.net.au/pdf2/vici09_113.pdf</a>
SSAC10W	Replacement rotor 10 ports 1/16in fittings 0.40mm ports	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_113.pdf">http://www.chromtech.net.au/pdf2/vici09_113.pdf</a>
SL10KCW	10ml Stainless steel loops	Valco Valves	\$87.00	<a href="http://www.chromtech.net.au/pdf2/vici09_113.pdf">http://www.chromtech.net.au/pdf2/vici09_113.pdf</a>
EPC10W	Microelectric Actuated 10 ports 1/16in fittings 0.40mm ports	Valco Valves	\$1,754.00	<a href="http://www.chromtech.net.au/pdf2/vici09_113.pdf">http://www.chromtech.net.au/pdf2/vici09_113.pdf</a>
DC4W	Replacement valve 4 ports 1/16in fittings 0.40mm ports	Valco Valves	\$669.00	<a href="http://www.chromtech.net.au/pdf2/vici09_113.pdf">http://www.chromtech.net.au/pdf2/vici09_113.pdf</a>
DC6W	Replacement valve 6 ports 1/16in fittings 0.40mm ports	Valco Valves	\$740.00	<a href="http://www.chromtech.net.au/pdf2/vici09_113.pdf">http://www.chromtech.net.au/pdf2/vici09_113.pdf</a>
DC8W	Replacement valve 8 ports 1/16in fittings 0.40mm ports	Valco Valves	\$810.00	<a href="http://www.chromtech.net.au/pdf2/vici09_113.pdf">http://www.chromtech.net.au/pdf2/vici09_113.pdf</a>
DC10W	Replacement valve 10 ports 1/16in fittings 0.40mm ports	Valco Valves	\$810.00	<a href="http://www.chromtech.net.au/pdf2/vici09_113.pdf">http://www.chromtech.net.au/pdf2/vici09_113.pdf</a>
SSAC4W	Replacement rotor 4 ports 1/16in fittings 0.40mm ports	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_113.pdf">http://www.chromtech.net.au/pdf2/vici09_113.pdf</a>
SL250CW	250ul 1/16in Stainless sample loops	Valco Valves	\$34.00	<a href="http://www.chromtech.net.au/pdf2/vici09_113.pdf">http://www.chromtech.net.au/pdf2/vici09_113.pdf</a>
SSAC8W	Replacement rotor 8 ports 1/16in fittings 0.40mm ports	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_113.pdf">http://www.chromtech.net.au/pdf2/vici09_113.pdf</a>
SL1KCW	1ml Stainless steel loops	Valco Valves	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_113.pdf">http://www.chromtech.net.au/pdf2/vici09_113.pdf</a>
SL2CW	2ul Stainless steel loops	Valco Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_113.pdf">http://www.chromtech.net.au/pdf2/vici09_113.pdf</a>
SL100CW	100ul Stainless steel loops	Valco Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_113.pdf">http://www.chromtech.net.au/pdf2/vici09_113.pdf</a>
SL50CW	50ul Stainless steel loops	Valco Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_113.pdf">http://www.chromtech.net.au/pdf2/vici09_113.pdf</a>
EPC8W	Microelectric Actuated 8 ports 1/16in fittings 0.40mm ports	Valco Valves	\$1,754.00	<a href="http://www.chromtech.net.au/pdf2/vici09_113.pdf">http://www.chromtech.net.au/pdf2/vici09_113.pdf</a>
SL10CW	10ul Stainless steel loops	Valco Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_113.pdf">http://www.chromtech.net.au/pdf2/vici09_113.pdf</a>
SL500CW	500ul 1/16in Stainless sample loops	Valco Valves	\$40.00	<a href="http://www.chromtech.net.au/pdf2/vici09_113.pdf">http://www.chromtech.net.au/pdf2/vici09_113.pdf</a>
SL15CW	15ul Stainless steel loops	Valco Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_113.pdf">http://www.chromtech.net.au/pdf2/vici09_113.pdf</a>
SSAC6W	Replacement rotor 6 ports 1/16in fittings 0.40mm ports	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_113.pdf">http://www.chromtech.net.au/pdf2/vici09_113.pdf</a>
EPCI4W.06	Microelectric Actuated 1/16in fittings 0.40mm ports x0.25mm column port	Valco Valves	\$1,677.00	<a href="http://www.chromtech.net.au/pdf2/vici09_112.pdf">http://www.chromtech.net.au/pdf2/vici09_112.pdf</a>
SSACI4W.5	Replacement Valve 1/16in fittings 0.40mm ports x0.25mm column port	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_112.pdf">http://www.chromtech.net.au/pdf2/vici09_112.pdf</a>
SSACI4W.2	Replacement Valve 1/16in fittings 0.40mm ports x0.25mm column port	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_112.pdf">http://www.chromtech.net.au/pdf2/vici09_112.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

SSACI4W.1	Replacement Valve 1/16in fittings 0.40mm ports x0.25mm column port	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_112.pdf">http://www.chromtech.net.au/pdf2/vici09_112.pdf</a>
SSACI4W.06	Replacement Valve 1/16in fittings 0.40mm ports x0.25mm column port	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_112.pdf">http://www.chromtech.net.au/pdf2/vici09_112.pdf</a>
DCI4W.5	Replacement Valve 1/16in fittings 0.40mm ports x0.25mm column port	Valco Valves	\$734.00	<a href="http://www.chromtech.net.au/pdf2/vici09_112.pdf">http://www.chromtech.net.au/pdf2/vici09_112.pdf</a>
DCI4W.2	Replacement Valve 1/16in fittings 0.40mm ports x0.25mm column port	Valco Valves	\$734.00	<a href="http://www.chromtech.net.au/pdf2/vici09_112.pdf">http://www.chromtech.net.au/pdf2/vici09_112.pdf</a>
DCI4W.1	Replacement Valve 1/16in fittings 0.40mm ports x0.25mm column port	Valco Valves	\$734.00	<a href="http://www.chromtech.net.au/pdf2/vici09_112.pdf">http://www.chromtech.net.au/pdf2/vici09_112.pdf</a>
DCI4W.06	Replacement Valve 1/16in fittings 0.40mm ports x0.25mm column port	Valco Valves	\$734.00	<a href="http://www.chromtech.net.au/pdf2/vici09_112.pdf">http://www.chromtech.net.au/pdf2/vici09_112.pdf</a>
EPCI4W.2	Microelectric Actuated 1/16in fittings 0.40mm ports x0.25mm column port	Valco Valves	\$1,677.00	<a href="http://www.chromtech.net.au/pdf2/vici09_112.pdf">http://www.chromtech.net.au/pdf2/vici09_112.pdf</a>
CI4UW.2	14 Ports Manual 1/16in fittings 0.75mm ports 0.2ul	Valco Valves	\$784.00	<a href="http://www.chromtech.net.au/pdf2/vici09_112.pdf">http://www.chromtech.net.au/pdf2/vici09_112.pdf</a>
EDCI4UW.2	14 Ports w/Microelectric actuator 1/16in fittings x0.75mm column port 0.2ul	Valco Valves	\$1,746.00	<a href="http://www.chromtech.net.au/pdf2/vici09_112.pdf">http://www.chromtech.net.au/pdf2/vici09_112.pdf</a>
SSACI4UW2	14 Ports Replacement rotor 1/16in fittings x0.75mm column port 2ul	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_112.pdf">http://www.chromtech.net.au/pdf2/vici09_112.pdf</a>
SSACI4UW1	14 Ports Replacement rotor 1/16in fittings x0.75mm column port 1ul	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_112.pdf">http://www.chromtech.net.au/pdf2/vici09_112.pdf</a>
SSACI4UW.5	14 Ports Replacement rotor 1/16in fittings x0.75mm column port 0.5ul	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_112.pdf">http://www.chromtech.net.au/pdf2/vici09_112.pdf</a>
SSACI4UW.2	14 Ports Replacement rotor 1/16in fittings x0.75mm column port 0.2ul	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_112.pdf">http://www.chromtech.net.au/pdf2/vici09_112.pdf</a>
DCI4UW2	Replacement Valve 1/16in fittings 0.75mm ports	Valco Valves	\$726.00	<a href="http://www.chromtech.net.au/pdf2/vici09_112.pdf">http://www.chromtech.net.au/pdf2/vici09_112.pdf</a>
DCI4UW1	14 Ports Replacement valve 1/16in fittings x0.75mm column port 1ul	Valco Valves	\$726.00	<a href="http://www.chromtech.net.au/pdf2/vici09_112.pdf">http://www.chromtech.net.au/pdf2/vici09_112.pdf</a>
DCI4UW.5	14 Ports Replacement valve 1/16in fittings x0.75mm column port 0.5ul	Valco Valves	\$726.00	<a href="http://www.chromtech.net.au/pdf2/vici09_112.pdf">http://www.chromtech.net.au/pdf2/vici09_112.pdf</a>
DCI4UW.2	14 Ports Replacement valve 1/16in fittings x0.75mm column port 0.2ul	Valco Valves	\$726.00	<a href="http://www.chromtech.net.au/pdf2/vici09_112.pdf">http://www.chromtech.net.au/pdf2/vici09_112.pdf</a>
EDCI4UW2	14 Ports w/Microelectric actuator 1/16in fittings x0.75mm column port 2ul	Valco Valves	\$1,746.00	<a href="http://www.chromtech.net.au/pdf2/vici09_112.pdf">http://www.chromtech.net.au/pdf2/vici09_112.pdf</a>
CI4UW.5	14 Ports Manual 1/16in fittings 0.75mm ports 0.5ul	Valco Valves	\$784.00	<a href="http://www.chromtech.net.au/pdf2/vici09_112.pdf">http://www.chromtech.net.au/pdf2/vici09_112.pdf</a>
EDCI4UW.5	14 Ports w/Microelectric actuator 1/16in fittings x0.75mm column port 0.5ul	Valco Valves	\$1,746.00	<a href="http://www.chromtech.net.au/pdf2/vici09_112.pdf">http://www.chromtech.net.au/pdf2/vici09_112.pdf</a>



## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

CI4UW1	14 Ports Manual 1/16in fittings 0.75mm ports 1ul	Valco Valves	\$784.00	<a href="http://www.chromtech.net.au/pdf2/vici09_112.pdf">http://www.chromtech.net.au/pdf2/vici09_112.pdf</a>
ECI4UW2	Electric Actuated 1/16in fittings 0.75mm ports	Valco Valves	\$1,454.00	<a href="http://www.chromtech.net.au/pdf2/vici09_112.pdf">http://www.chromtech.net.au/pdf2/vici09_112.pdf</a>
ECI4UW1	Electric Actuated 1/16in fittings 0.75mm ports	Valco Valves	\$1,454.00	<a href="http://www.chromtech.net.au/pdf2/vici09_112.pdf">http://www.chromtech.net.au/pdf2/vici09_112.pdf</a>
ECI4UW.5	14 Ports w/Std Electric actuator 1/16in fittings x0.75mm column port 0.5ul	Valco Valves	\$1,454.00	<a href="http://www.chromtech.net.au/pdf2/vici09_112.pdf">http://www.chromtech.net.au/pdf2/vici09_112.pdf</a>
ECI4UW.2	14 Ports w/Std Electric actuator 1/16in fittings x0.75mm column port 0.2ul	Valco Valves	\$1,454.00	<a href="http://www.chromtech.net.au/pdf2/vici09_112.pdf">http://www.chromtech.net.au/pdf2/vici09_112.pdf</a>
ACI4UW2	Air Actuated 1/16in fittings 0.75mm ports 2ul	Valco Valves	\$1,046.00	<a href="http://www.chromtech.net.au/pdf2/vici09_112.pdf">http://www.chromtech.net.au/pdf2/vici09_112.pdf</a>
ACI4UW1	Air Actuated 1/16in fittings 0.75mm ports 1ul	Valco Valves	\$1,046.00	<a href="http://www.chromtech.net.au/pdf2/vici09_112.pdf">http://www.chromtech.net.au/pdf2/vici09_112.pdf</a>
ACI4UW.5	Air Actuated 1/16in fittings 0.75mm ports 0.5ul	Valco Valves	\$1,046.00	<a href="http://www.chromtech.net.au/pdf2/vici09_112.pdf">http://www.chromtech.net.au/pdf2/vici09_112.pdf</a>
ACI4UW.2	Air Actuated 1/16in fittings 0.75mm ports 0.2ul	Valco Valves	\$1,046.00	<a href="http://www.chromtech.net.au/pdf2/vici09_112.pdf">http://www.chromtech.net.au/pdf2/vici09_112.pdf</a>
CI4UW2	14 Ports Manual 1/16in fittings 0.75mm ports 2ul	Valco Valves	\$784.00	<a href="http://www.chromtech.net.au/pdf2/vici09_112.pdf">http://www.chromtech.net.au/pdf2/vici09_112.pdf</a>
EPCI4W.1	Microelectric Actuated 1/16in fittings 0.40mm ports x0.25mm column port	Valco Valves	\$1,677.00	<a href="http://www.chromtech.net.au/pdf2/vici09_112.pdf">http://www.chromtech.net.au/pdf2/vici09_112.pdf</a>
EDCI4UW1	14 Ports w/Microelectric actuator 1/16in fittings x0.75mm column port 1ul	Valco Valves	\$1,746.00	<a href="http://www.chromtech.net.au/pdf2/vici09_112.pdf">http://www.chromtech.net.au/pdf2/vici09_112.pdf</a>
CI4W.5	14 Ports Manual 1/16in fittings 0.40mm ports x0.25mm column port 0.5ul	Valco Valves	\$790.00	<a href="http://www.chromtech.net.au/pdf2/vici09_112.pdf">http://www.chromtech.net.au/pdf2/vici09_112.pdf</a>
ECI4W.5	Electric Actuated 1/16in fittings 0.40mm ports x0.25mm column port	Valco Valves	\$1,460.00	<a href="http://www.chromtech.net.au/pdf2/vici09_112.pdf">http://www.chromtech.net.au/pdf2/vici09_112.pdf</a>
EPCI4W.5	Microelectric Actuated 1/16in fittings 0.40mm ports x0.25mm column port	Valco Valves	\$1,677.00	<a href="http://www.chromtech.net.au/pdf2/vici09_112.pdf">http://www.chromtech.net.au/pdf2/vici09_112.pdf</a>
CI4W.06	14 Ports Manual 1/16in fittings 0.40mm ports x0.25mm column port 0.06ul	Valco Valves	\$790.00	<a href="http://www.chromtech.net.au/pdf2/vici09_112.pdf">http://www.chromtech.net.au/pdf2/vici09_112.pdf</a>
CI4W.2	14 Ports Manual 1/16in fittings 0.40mm ports x0.25mm column port 0.2ul	Valco Valves	\$790.00	<a href="http://www.chromtech.net.au/pdf2/vici09_112.pdf">http://www.chromtech.net.au/pdf2/vici09_112.pdf</a>
ACI4W.06	Air Actuated 1/16in fittings 0.40mm ports x0.25mm column port 0.06ul	Valco Valves	\$1,052.00	<a href="http://www.chromtech.net.au/pdf2/vici09_112.pdf">http://www.chromtech.net.au/pdf2/vici09_112.pdf</a>
ACI4W.1	Air Actuated 1/16in fittings 0.40mm ports x0.25mm column port 0.1ul	Valco Valves	\$1,052.00	<a href="http://www.chromtech.net.au/pdf2/vici09_112.pdf">http://www.chromtech.net.au/pdf2/vici09_112.pdf</a>
ACI4W.2	Air Actuated 1/16in fittings 0.40mm ports x0.25mm column port 0.2ul	Valco Valves	\$1,052.00	<a href="http://www.chromtech.net.au/pdf2/vici09_112.pdf">http://www.chromtech.net.au/pdf2/vici09_112.pdf</a>
ACI4W.5	Air Actuated 1/16in fittings 0.40mm ports x0.25mm column port 0.5ul	Valco Valves	\$1,052.00	<a href="http://www.chromtech.net.au/pdf2/vici09_112.pdf">http://www.chromtech.net.au/pdf2/vici09_112.pdf</a>
ECI4W.06	Electric Actuated 1/16in fittings 0.40mm ports x0.25mm column port	Valco Valves	\$1,460.00	<a href="http://www.chromtech.net.au/pdf2/vici09_112.pdf">http://www.chromtech.net.au/pdf2/vici09_112.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

ECI4W.1	Electric Actuated 1/16in fittings 0.40mm ports x0.25mm column port	Valco Valves	\$1,460.00	<a href="http://www.chromtech.net.au/pdf2/vici09_112.pdf">http://www.chromtech.net.au/pdf2/vici09_112.pdf</a>
ECI4W.2	Electric Actuated 1/16in fittings 0.40mm ports x0.25mm column port	Valco Valves	\$1,460.00	<a href="http://www.chromtech.net.au/pdf2/vici09_112.pdf">http://www.chromtech.net.au/pdf2/vici09_112.pdf</a>
CI4W.1	14 Ports Manual 1/16in fittings 0.40mm ports x0.25mm column port 0.1ul	Valco Valves	\$790.00	<a href="http://www.chromtech.net.au/pdf2/vici09_112.pdf">http://www.chromtech.net.au/pdf2/vici09_112.pdf</a>
A410UWT	10 Ports Manual w/standoff 1/8infittings x0.75mm ports	Valco Valves	\$1,059.00	<a href="http://www.chromtech.net.au/pdf2/vici09_111.pdf">http://www.chromtech.net.au/pdf2/vici09_111.pdf</a>
44UWT	4 Ports Manual w/standoff 1/8infittings x0.75mm ports	Valco Valves	\$714.00	<a href="http://www.chromtech.net.au/pdf2/vici09_111.pdf">http://www.chromtech.net.au/pdf2/vici09_111.pdf</a>
46UWT	6 Ports Manual 1/8infittings x0.75mm ports	Valco Valves	\$784.00	<a href="http://www.chromtech.net.au/pdf2/vici09_111.pdf">http://www.chromtech.net.au/pdf2/vici09_111.pdf</a>
SL50UW	50ul Stainless steel loops	Valco Valves	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_111.pdf">http://www.chromtech.net.au/pdf2/vici09_111.pdf</a>
410UWT	10 Ports Manual 1/8infittings x0.75mm ports	Valco Valves	\$855.00	<a href="http://www.chromtech.net.au/pdf2/vici09_111.pdf">http://www.chromtech.net.au/pdf2/vici09_111.pdf</a>
SL15UW	15ul 1/8in Stainless steel loops	Valco Valves	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_111.pdf">http://www.chromtech.net.au/pdf2/vici09_111.pdf</a>
SL20KUW	20ml 1/8in Stainless steel loops	Valco Valves	\$128.00	<a href="http://www.chromtech.net.au/pdf2/vici09_111.pdf">http://www.chromtech.net.au/pdf2/vici09_111.pdf</a>
SL10KUW	1ml Stainless steel loops	Valco Valves	\$72.00	<a href="http://www.chromtech.net.au/pdf2/vici09_111.pdf">http://www.chromtech.net.au/pdf2/vici09_111.pdf</a>
SL100UW	100ul Stainless steel loops	Valco Valves	\$34.00	<a href="http://www.chromtech.net.au/pdf2/vici09_111.pdf">http://www.chromtech.net.au/pdf2/vici09_111.pdf</a>
SL5KUW	5ml 1/8in Stainless steel loops	Valco Valves	\$51.00	<a href="http://www.chromtech.net.au/pdf2/vici09_111.pdf">http://www.chromtech.net.au/pdf2/vici09_111.pdf</a>
48UWT	8 Ports Manual 1/8infittings x0.75mm ports	Valco Valves	\$855.00	<a href="http://www.chromtech.net.au/pdf2/vici09_111.pdf">http://www.chromtech.net.au/pdf2/vici09_111.pdf</a>
SL2KUW	2ml 1/8in Stainless steel loops	Valco Valves	\$42.00	<a href="http://www.chromtech.net.au/pdf2/vici09_111.pdf">http://www.chromtech.net.au/pdf2/vici09_111.pdf</a>
SL1KUW	1ml 1/8in Stainless steel loops	Valco Valves	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_111.pdf">http://www.chromtech.net.au/pdf2/vici09_111.pdf</a>
SL25UW	25ul 1/8in Stainless steel loops	Valco Valves	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_111.pdf">http://www.chromtech.net.au/pdf2/vici09_111.pdf</a>
SL500UW	500ul 1/8in Stainless steel loops	Valco Valves	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_111.pdf">http://www.chromtech.net.au/pdf2/vici09_111.pdf</a>
A46UWT	6 Ports Manual w/standoff 1/8infittings x0.75mm ports	Valco Valves	\$988.00	<a href="http://www.chromtech.net.au/pdf2/vici09_111.pdf">http://www.chromtech.net.au/pdf2/vici09_111.pdf</a>
SL250UW	250ul Stainless steel loops	Valco Valves	\$34.00	<a href="http://www.chromtech.net.au/pdf2/vici09_111.pdf">http://www.chromtech.net.au/pdf2/vici09_111.pdf</a>
SL10UW	10ul Stainless steel loops	Valco Valves	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_111.pdf">http://www.chromtech.net.au/pdf2/vici09_111.pdf</a>
SSA10UWT	10 Ports Replacement rotor 1/8infittings x0.75mm ports	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_111.pdf">http://www.chromtech.net.au/pdf2/vici09_111.pdf</a>
SSA8UWT	8 Ports Replacement rotor 1/8infittings x0.75mm ports	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_111.pdf">http://www.chromtech.net.au/pdf2/vici09_111.pdf</a>
SSA6UWT	8 Ports Replacement rotor 1/8infittings x0.75mm ports	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_111.pdf">http://www.chromtech.net.au/pdf2/vici09_111.pdf</a>
SSA4UWT	4 Ports Replacement rotor 1/8infittings x0.75mm ports	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_111.pdf">http://www.chromtech.net.au/pdf2/vici09_111.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

D10UWT	10 Ports Replacement valve 1/8in fittings x0.75mm ports	Valco Valves	\$740.00	<a href="http://www.chromtech.net.au/pdf2/vici09_111.pdf">http://www.chromtech.net.au/pdf2/vici09_111.pdf</a>
A48UWT	8 Ports Manual w/standoff 1/8in fittings x0.75mm ports	Valco Valves	\$1,059.00	<a href="http://www.chromtech.net.au/pdf2/vici09_111.pdf">http://www.chromtech.net.au/pdf2/vici09_111.pdf</a>
SL20UW	20ul 1/8in Stainless steel loops	Valco Valves	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_111.pdf">http://www.chromtech.net.au/pdf2/vici09_111.pdf</a>
E44UWT	4 Ports w/standard electric actuator 1/8in fittings x0.75mm ports	Valco Valves	\$1,326.00	<a href="http://www.chromtech.net.au/pdf2/vici09_111.pdf">http://www.chromtech.net.au/pdf2/vici09_111.pdf</a>
D6UWT	6Ports Replacement valve 1/8in fittings x0.75mm ports	Valco Valves	\$669.00	<a href="http://www.chromtech.net.au/pdf2/vici09_111.pdf">http://www.chromtech.net.au/pdf2/vici09_111.pdf</a>
A44UWT	4 Ports Manual w/standoff 1/8in fittings x0.75mm ports	Valco Valves	\$918.00	<a href="http://www.chromtech.net.au/pdf2/vici09_111.pdf">http://www.chromtech.net.au/pdf2/vici09_111.pdf</a>
D4UWT	4 Ports Replacement valve 1/8in fittings x0.75mm ports	Valco Valves	\$600.00	<a href="http://www.chromtech.net.au/pdf2/vici09_111.pdf">http://www.chromtech.net.au/pdf2/vici09_111.pdf</a>
ED410UWT	8 Ports w/Microelectric actuator 1/8in fittings x0.75mm ports	Valco Valves	\$1,760.00	<a href="http://www.chromtech.net.au/pdf2/vici09_111.pdf">http://www.chromtech.net.au/pdf2/vici09_111.pdf</a>
ED48UWT	8 Ports w/Microelectric actuator 1/8in fittings x0.75mm ports	Valco Valves	\$1,760.00	<a href="http://www.chromtech.net.au/pdf2/vici09_111.pdf">http://www.chromtech.net.au/pdf2/vici09_111.pdf</a>
ED46UWT	6 Ports w/Microelectric actuator 1/8in fittings x0.75mm ports	Valco Valves	\$1,689.00	<a href="http://www.chromtech.net.au/pdf2/vici09_111.pdf">http://www.chromtech.net.au/pdf2/vici09_111.pdf</a>
ED44UWT	4 Ports w/Microelectric actuator 1/8in fittings x0.75mm ports	Valco Valves	\$1,620.00	<a href="http://www.chromtech.net.au/pdf2/vici09_111.pdf">http://www.chromtech.net.au/pdf2/vici09_111.pdf</a>
E410UWT	10 Ports w/standard electric actuator 1/8in fittings x0.75mm ports	Valco Valves	\$1,467.00	<a href="http://www.chromtech.net.au/pdf2/vici09_111.pdf">http://www.chromtech.net.au/pdf2/vici09_111.pdf</a>
E48UWT	8 Ports w/standard electric actuator 1/8in fittings x0.75mm ports	Valco Valves	\$1,467.00	<a href="http://www.chromtech.net.au/pdf2/vici09_111.pdf">http://www.chromtech.net.au/pdf2/vici09_111.pdf</a>
D8UWT	8 Ports Replacement valve 1/8in fittings x0.75mm ports	Valco Valves	\$740.00	<a href="http://www.chromtech.net.au/pdf2/vici09_111.pdf">http://www.chromtech.net.au/pdf2/vici09_111.pdf</a>
E46UWT	6 Ports w/standard electric actuator 1/8in fittings x0.75mm ports	Valco Valves	\$1,396.00	<a href="http://www.chromtech.net.au/pdf2/vici09_111.pdf">http://www.chromtech.net.au/pdf2/vici09_111.pdf</a>
E4VL8MWE2	8 Ports w/electric actuator Low temp 1/4in x 4.0mm ports	Valco Valves	\$2,079.00	<a href="http://www.chromtech.net.au/pdf2/vici09_110.pdf">http://www.chromtech.net.au/pdf2/vici09_110.pdf</a>
46UWE	6 Ports Manual with standoff Medium temp 1/8in fittings x0.75mm ports	Valco Valves	\$784.00	<a href="http://www.chromtech.net.au/pdf2/vici09_110.pdf">http://www.chromtech.net.au/pdf2/vici09_110.pdf</a>
8UWE	8 Ports Manual Medium temp 1/8in fittings x0.75mm ports	Valco Valves	\$796.00	<a href="http://www.chromtech.net.au/pdf2/vici09_110.pdf">http://www.chromtech.net.au/pdf2/vici09_110.pdf</a>
6UWE	6 Ports Manual Medium temp 1/8in fittings x0.75mm ports	Valco Valves	\$726.00	<a href="http://www.chromtech.net.au/pdf2/vici09_110.pdf">http://www.chromtech.net.au/pdf2/vici09_110.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

4UWE	4 Ports Manual Medium temp 1/8in fittings x0.75mm ports	Valco Valves	\$657.00	<a href="http://www.chromtech.net.au/pdf2/vici09_110.pdf">http://www.chromtech.net.au/pdf2/vici09_110.pdf</a>
E4VL6MWE2	6 Ports w/electric actuator Low temp 1/4in x 4.0mm ports	Valco Valves	\$2,008.00	<a href="http://www.chromtech.net.au/pdf2/vici09_110.pdf">http://www.chromtech.net.au/pdf2/vici09_110.pdf</a>
E4VL4MWE2	4 Ports w/electric actuator Low temp 1/4in x 4.0mm ports	Valco Valves	\$1,938.00	<a href="http://www.chromtech.net.au/pdf2/vici09_110.pdf">http://www.chromtech.net.au/pdf2/vici09_110.pdf</a>
A4VL8MWE2	4 Ports w/air actuator Low temp 1/4in x 4.0mm ports	Valco Valves	\$1,671.00	<a href="http://www.chromtech.net.au/pdf2/vici09_110.pdf">http://www.chromtech.net.au/pdf2/vici09_110.pdf</a>
A4VL6MWE2	4 Ports w/air actuator Low temp 1/4in x 4.0mm ports	Valco Valves	\$1,600.00	<a href="http://www.chromtech.net.au/pdf2/vici09_110.pdf">http://www.chromtech.net.au/pdf2/vici09_110.pdf</a>
410UWE	10 Ports Manual with standoff Medium temp 1/8in fittings x0.75mm ports	Valco Valves	\$855.00	<a href="http://www.chromtech.net.au/pdf2/vici09_110.pdf">http://www.chromtech.net.au/pdf2/vici09_110.pdf</a>
ET4VL4MWE2	4 Ports w/Microelectric actuator Low temp 1/4in x 4.0mm ports	Valco Valves	\$2,232.00	<a href="http://www.chromtech.net.au/pdf2/vici09_110.pdf">http://www.chromtech.net.au/pdf2/vici09_110.pdf</a>
48UWE	8 Ports Manual with standoff Medium temp 1/8in fittings x0.75mm ports	Valco Valves	\$855.00	<a href="http://www.chromtech.net.au/pdf2/vici09_110.pdf">http://www.chromtech.net.au/pdf2/vici09_110.pdf</a>
ET4VL6MWE2	6 Ports w/Microelectric actuator Low temp 1/4in x 4.0mm ports	Valco Valves	\$2,301.00	<a href="http://www.chromtech.net.au/pdf2/vici09_110.pdf">http://www.chromtech.net.au/pdf2/vici09_110.pdf</a>
ET4VL8MWE2	8 Ports w/Microelectric actuator Low temp 1/4in x 4.0mm ports	Valco Valves	\$2,372.00	<a href="http://www.chromtech.net.au/pdf2/vici09_110.pdf">http://www.chromtech.net.au/pdf2/vici09_110.pdf</a>
DVL4MWE2	4 Ports Replacement valve Low temp 1/4in x 4.0mm ports	Valco Valves	\$1,212.00	<a href="http://www.chromtech.net.au/pdf2/vici09_110.pdf">http://www.chromtech.net.au/pdf2/vici09_110.pdf</a>
DVL6MWE2	6 Ports Replacement valve Low temp 1/4in x 4.0mm ports	Valco Valves	\$1,281.00	<a href="http://www.chromtech.net.au/pdf2/vici09_110.pdf">http://www.chromtech.net.au/pdf2/vici09_110.pdf</a>
DVL8MWE2	8 Ports Replacement valve Low temp 1/4in x 4.0mm ports	Valco Valves	\$1,352.00	<a href="http://www.chromtech.net.au/pdf2/vici09_110.pdf">http://www.chromtech.net.au/pdf2/vici09_110.pdf</a>
SSAVL4MWE2	4 Ports Replacement rotor Low temp 1/4in x 4.0mm ports	Valco Valves	\$224.00	<a href="http://www.chromtech.net.au/pdf2/vici09_110.pdf">http://www.chromtech.net.au/pdf2/vici09_110.pdf</a>
SSAVL6MWE2	6 Ports Replacement rotor Low temp 1/4in x 4.0mm ports	Valco Valves	\$224.00	<a href="http://www.chromtech.net.au/pdf2/vici09_110.pdf">http://www.chromtech.net.au/pdf2/vici09_110.pdf</a>
SSAVL8MWE2	8 Ports Replacement rotor Low temp 1/4in x 4.0mm ports	Valco Valves	\$224.00	<a href="http://www.chromtech.net.au/pdf2/vici09_110.pdf">http://www.chromtech.net.au/pdf2/vici09_110.pdf</a>
4VL8MWE2	6 Ports Manual Low temp 1/4in x 4.0mm ports	Valco Valves	\$1,498.00	<a href="http://www.chromtech.net.au/pdf2/vici09_110.pdf">http://www.chromtech.net.au/pdf2/vici09_110.pdf</a>
WA4VL4MWE2	10 Ports Replacement rotor Low temp 1/4in x 4.0mm ports	Valco Valves	\$1,530.00	<a href="http://www.chromtech.net.au/pdf2/vici09_110.pdf">http://www.chromtech.net.au/pdf2/vici09_110.pdf</a>
D4UWE	4 Ports Replacement valve Medium temp 1/8in fittings x0.75mm ports	Valco Valves	\$600.00	<a href="http://www.chromtech.net.au/pdf2/vici09_110.pdf">http://www.chromtech.net.au/pdf2/vici09_110.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

4VL6MWE2	4 Ports Manual Low temp 1/4in x 4.0mm ports	Valco Valves	\$1,428.00	<a href="http://www.chromtech.net.au/pdf2/vici09_110.pdf">http://www.chromtech.net.au/pdf2/vici09_110.pdf</a>
4VL4MWE2	4 Ports Manual Low temp 1/4in x 4.0mm ports	Valco Valves	\$1,358.00	<a href="http://www.chromtech.net.au/pdf2/vici09_110.pdf">http://www.chromtech.net.au/pdf2/vici09_110.pdf</a>
SSA10UWE	10 Ports Replacement rotor Medium temp 1/8in fittings x0.75mm ports	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_110.pdf">http://www.chromtech.net.au/pdf2/vici09_110.pdf</a>
SSA8UWE	8 Ports Replacement rotor Medium temp 1/8in fittings x0.75mm ports	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_110.pdf">http://www.chromtech.net.au/pdf2/vici09_110.pdf</a>
SSA6UWE	8 Ports Replacement rotor Medium temp 1/8in fittings x0.75mm ports	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_110.pdf">http://www.chromtech.net.au/pdf2/vici09_110.pdf</a>
SSA4UWE	4 Ports Replacement rotor Medium temp 1/8in fittings x0.75mm ports	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_110.pdf">http://www.chromtech.net.au/pdf2/vici09_110.pdf</a>
D10UWE	10 Ports Replacement valve Medium temp 1/8in fittings x0.75mm ports	Valco Valves	\$740.00	<a href="http://www.chromtech.net.au/pdf2/vici09_110.pdf">http://www.chromtech.net.au/pdf2/vici09_110.pdf</a>
D8UWE	8 Ports Replacement valve Medium temp 1/8in fittings x0.75mm ports	Valco Valves	\$740.00	<a href="http://www.chromtech.net.au/pdf2/vici09_110.pdf">http://www.chromtech.net.au/pdf2/vici09_110.pdf</a>
44UWE	4 Ports Manual with standoff Medium temp 1/8in fittings x0.75mm ports	Valco Valves	\$714.00	<a href="http://www.chromtech.net.au/pdf2/vici09_110.pdf">http://www.chromtech.net.au/pdf2/vici09_110.pdf</a>
D6UWE	6 Ports Replacement valve Medium temp 1/8in fittings x0.75mm ports	Valco Valves	\$669.00	<a href="http://www.chromtech.net.au/pdf2/vici09_110.pdf">http://www.chromtech.net.au/pdf2/vici09_110.pdf</a>
A44UWE	4 Ports w/air actuator Medium temp 1/8in fittings x0.75mm ports	Valco Valves	\$918.00	<a href="http://www.chromtech.net.au/pdf2/vici09_110.pdf">http://www.chromtech.net.au/pdf2/vici09_110.pdf</a>
E44UWE	4 Ports w/electric actuator Medium temp 1/8in fittings x0.75mm ports	Valco Valves	\$1,326.00	<a href="http://www.chromtech.net.au/pdf2/vici09_110.pdf">http://www.chromtech.net.au/pdf2/vici09_110.pdf</a>
A46UWE	6 Ports w/air actuator Medium temp 1/8in fittings x0.75mm ports	Valco Valves	\$988.00	<a href="http://www.chromtech.net.au/pdf2/vici09_110.pdf">http://www.chromtech.net.au/pdf2/vici09_110.pdf</a>
A410UWE	10 Ports w/air actuator Medium temp 1/8in fittings x0.75mm ports	Valco Valves	\$1,059.00	<a href="http://www.chromtech.net.au/pdf2/vici09_110.pdf">http://www.chromtech.net.au/pdf2/vici09_110.pdf</a>
ED410UWE	4 Ports w/Microelectric actuator Medium temp 1/8in x0.75mm ports	Valco Valves	\$1,760.00	<a href="http://www.chromtech.net.au/pdf2/vici09_110.pdf">http://www.chromtech.net.au/pdf2/vici09_110.pdf</a>
E46UWE	6 Ports w/electric actuator Medium temp 1/8in fittings x0.75mm ports	Valco Valves	\$1,396.00	<a href="http://www.chromtech.net.au/pdf2/vici09_110.pdf">http://www.chromtech.net.au/pdf2/vici09_110.pdf</a>
E48UWE	8 Ports w/electric actuator Medium temp 1/8in fittings x0.75mm ports	Valco Valves	\$1,467.00	<a href="http://www.chromtech.net.au/pdf2/vici09_110.pdf">http://www.chromtech.net.au/pdf2/vici09_110.pdf</a>
E410UWE	10 Ports w/electric actuator Medium temp 1/8in fittings x0.75mm ports	Valco Valves	\$1,467.00	<a href="http://www.chromtech.net.au/pdf2/vici09_110.pdf">http://www.chromtech.net.au/pdf2/vici09_110.pdf</a>
ED44UWE	4 Ports w/Microelectric actuator Medium temp 1/8in x0.75mm ports	Valco Valves	\$1,620.00	<a href="http://www.chromtech.net.au/pdf2/vici09_110.pdf">http://www.chromtech.net.au/pdf2/vici09_110.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

ED46UWE	6 Ports w/Microelecric actuator Medium temp 1/8in x0.75mm ports	Valco Valves	\$1,689.00	<a href="http://www.chromtech.net.au/pdf2/vici09_110.pdf">http://www.chromtech.net.au/pdf2/vici09_110.pdf</a>
ED48UWE	4 Ports w/Microelecric actuator Medium temp 1/8in x0.75mm ports	Valco Valves	\$1,760.00	<a href="http://www.chromtech.net.au/pdf2/vici09_110.pdf">http://www.chromtech.net.au/pdf2/vici09_110.pdf</a>
A48UWE	8 Ports w/air actuator Medium temp 1/8infittings x0.75mm ports	Valco Valves	\$1,059.00	<a href="http://www.chromtech.net.au/pdf2/vici09_110.pdf">http://www.chromtech.net.au/pdf2/vici09_110.pdf</a>
1EN1K	1" external nut	Valco Valves	\$20.00	<a href="http://www.chromtech.net.au/pdf2/vici09_11.pdf">http://www.chromtech.net.au/pdf2/vici09_11.pdf</a>
EN.5	Stainless steel 1/32in external nut	Valco Fittings	\$6.00	<a href="http://www.chromtech.net.au/pdf2/vici09_11.pdf">http://www.chromtech.net.au/pdf2/vici09_11.pdf</a>
EN1	Stainless steel 1/16in external nut	Valco Fittings	\$3.00	<a href="http://www.chromtech.net.au/pdf2/vici09_11.pdf">http://www.chromtech.net.au/pdf2/vici09_11.pdf</a>
EN2	Stainless steel 1/8in external nut	Valco Fittings	\$3.00	<a href="http://www.chromtech.net.au/pdf2/vici09_11.pdf">http://www.chromtech.net.au/pdf2/vici09_11.pdf</a>
EN4	Stainless steel 1/4in external nut	Valco Fittings	\$3.00	<a href="http://www.chromtech.net.au/pdf2/vici09_11.pdf">http://www.chromtech.net.au/pdf2/vici09_11.pdf</a>
EN6	Stainless steel 3/8in external nut	Valco Fittings	\$4.00	<a href="http://www.chromtech.net.au/pdf2/vici09_11.pdf">http://www.chromtech.net.au/pdf2/vici09_11.pdf</a>
EN8	Stainless steel 1/2in external nut	Valco Fittings	\$8.00	<a href="http://www.chromtech.net.au/pdf2/vici09_11.pdf">http://www.chromtech.net.au/pdf2/vici09_11.pdf</a>
EN.5KN	Stainless steel 1/32in external nut knurled	Valco Fittings	\$6.00	<a href="http://www.chromtech.net.au/pdf2/vici09_11.pdf">http://www.chromtech.net.au/pdf2/vici09_11.pdf</a>
LZP1H	1/16" .75"	Valco Fittings	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_11.pdf">http://www.chromtech.net.au/pdf2/vici09_11.pdf</a>
ZP.5	Stainless steel 1/32in Plugs	Valco Fittings	\$10.00	<a href="http://www.chromtech.net.au/pdf2/vici09_11.pdf">http://www.chromtech.net.au/pdf2/vici09_11.pdf</a>
LZP1	Stainless steel 1/16in x.750in long Plugs	Valco Fittings	\$9.00	<a href="http://www.chromtech.net.au/pdf2/vici09_11.pdf">http://www.chromtech.net.au/pdf2/vici09_11.pdf</a>
ZP1H	Stainless steel 1/16in High pressure Plugs	Valco Fittings	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_11.pdf">http://www.chromtech.net.au/pdf2/vici09_11.pdf</a>
ZP2	Stainless steel 1/8in Plugs	Valco Fittings	\$8.00	<a href="http://www.chromtech.net.au/pdf2/vici09_11.pdf">http://www.chromtech.net.au/pdf2/vici09_11.pdf</a>
LZP2	Stainless steel 1/8in x.825in long Plugs	Valco Fittings	\$9.00	<a href="http://www.chromtech.net.au/pdf2/vici09_11.pdf">http://www.chromtech.net.au/pdf2/vici09_11.pdf</a>
ZP4	Stainless steel 1/4in Plugs	Valco Fittings	\$10.00	<a href="http://www.chromtech.net.au/pdf2/vici09_11.pdf">http://www.chromtech.net.au/pdf2/vici09_11.pdf</a>
ZP2H	Stainless steel 1/16in High pressure Plugs	Valco Fittings	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_11.pdf">http://www.chromtech.net.au/pdf2/vici09_11.pdf</a>
ZP2H	Stainless steel 1/16in High pressure Plugs	Valco Fittings	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_11.pdf">http://www.chromtech.net.au/pdf2/vici09_11.pdf</a>
ZC.5	Stainless steel 1/16in Caps	Valco Fittings	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_11.pdf">http://www.chromtech.net.au/pdf2/vici09_11.pdf</a>
ZC1	Stainless steel 1/16in Caps	Valco Fittings	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_11.pdf">http://www.chromtech.net.au/pdf2/vici09_11.pdf</a>
ZC2	Stainless steel 1/8in Caps	Valco Fittings	\$12.00	<a href="http://www.chromtech.net.au/pdf2/vici09_11.pdf">http://www.chromtech.net.au/pdf2/vici09_11.pdf</a>
ZC4	Stainless steel 1/4in Caps	Valco Fittings	\$15.00	<a href="http://www.chromtech.net.au/pdf2/vici09_11.pdf">http://www.chromtech.net.au/pdf2/vici09_11.pdf</a>
ZP.5H	Stainless steel 1/32in High pressure Plugs	Valco Fittings	\$15.00	<a href="http://www.chromtech.net.au/pdf2/vici09_11.pdf">http://www.chromtech.net.au/pdf2/vici09_11.pdf</a>
ZP1	Stainless steel 1/16in Plugs	Valco Fittings	\$8.00	<a href="http://www.chromtech.net.au/pdf2/vici09_11.pdf">http://www.chromtech.net.au/pdf2/vici09_11.pdf</a>
SL500CUW	500ul Stainless steel loops	Valco Valves	\$40.00	<a href="http://www.chromtech.net.au/pdf2/vici09_109.pdf">http://www.chromtech.net.au/pdf2/vici09_109.pdf</a>
SSAC4UWT	4 Ports Replacement rotor 1/16infittings x0.75mm ports	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_109.pdf">http://www.chromtech.net.au/pdf2/vici09_109.pdf</a>
SSAC6UWT	8 Ports Replacement rotor 1/16infittings x0.75mm ports	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_109.pdf">http://www.chromtech.net.au/pdf2/vici09_109.pdf</a>
SSAC8UWT	8 Ports Replacement rotor 1/16infittings x0.75mm ports	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_109.pdf">http://www.chromtech.net.au/pdf2/vici09_109.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

SSAC10UWT	10 Ports Replacement rotor 1/16in fittings x0.75mm ports	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_109.pdf">http://www.chromtech.net.au/pdf2/vici09_109.pdf</a>
SL5CUW	5ul 1/16in stainless steel loop	Valco Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_109.pdf">http://www.chromtech.net.au/pdf2/vici09_109.pdf</a>
SL100CUW	100ul Stainless steel loops	Valco Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_109.pdf">http://www.chromtech.net.au/pdf2/vici09_109.pdf</a>
SL10CUW	10ul Stainless steel loops	Valco Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_109.pdf">http://www.chromtech.net.au/pdf2/vici09_109.pdf</a>
SL10KCUW	10ml Stainless steel loops	Valco Valves	\$87.00	<a href="http://www.chromtech.net.au/pdf2/vici09_109.pdf">http://www.chromtech.net.au/pdf2/vici09_109.pdf</a>
SL15CUW	15ul Stainless steel loops	Valco Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_109.pdf">http://www.chromtech.net.au/pdf2/vici09_109.pdf</a>
DC6UWT	6 Ports Replacement valve 1/16in fittings x0.75mm ports	Valco Valves	\$669.00	<a href="http://www.chromtech.net.au/pdf2/vici09_109.pdf">http://www.chromtech.net.au/pdf2/vici09_109.pdf</a>
SL20CUW	20ul Stainless steel loops	Valco Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_109.pdf">http://www.chromtech.net.au/pdf2/vici09_109.pdf</a>
SL1KCUW	1ml 1/16in stainless steel loop	Valco Valves	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_109.pdf">http://www.chromtech.net.au/pdf2/vici09_109.pdf</a>
SL25CUW	25ul 1/16in stainless steel loop	Valco Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_109.pdf">http://www.chromtech.net.au/pdf2/vici09_109.pdf</a>
SL2KCUW	2ml Stainless steel loops	Valco Valves	\$57.00	<a href="http://www.chromtech.net.au/pdf2/vici09_109.pdf">http://www.chromtech.net.au/pdf2/vici09_109.pdf</a>
SL50CUW	50ul Stainless steel loops	Valco Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_109.pdf">http://www.chromtech.net.au/pdf2/vici09_109.pdf</a>
SL5KCUW	5ml Stainless steel loops	Valco Valves	\$68.00	<a href="http://www.chromtech.net.au/pdf2/vici09_109.pdf">http://www.chromtech.net.au/pdf2/vici09_109.pdf</a>
SL250CUW	250ul Stainless steel loops	Valco Valves	\$34.00	<a href="http://www.chromtech.net.au/pdf2/vici09_109.pdf">http://www.chromtech.net.au/pdf2/vici09_109.pdf</a>
E4C4UWT	4 Ports w/standard electric actuator 1/16in fittings x0.75mm ports	Valco Valves	\$1,326.00	<a href="http://www.chromtech.net.au/pdf2/vici09_109.pdf">http://www.chromtech.net.au/pdf2/vici09_109.pdf</a>
4C4UWT	1/16in fittings x0.75mm ports 4 Ports Manual w/standoff	Valco Valves	\$714.00	<a href="http://www.chromtech.net.au/pdf2/vici09_109.pdf">http://www.chromtech.net.au/pdf2/vici09_109.pdf</a>
4C6UWT	6 Ports Manual 1/16in fittings x0.75mm ports	Valco Valves	\$784.00	<a href="http://www.chromtech.net.au/pdf2/vici09_109.pdf">http://www.chromtech.net.au/pdf2/vici09_109.pdf</a>
4C8UWT	8 Ports Manual 1/16in fittings x0.75mm ports	Valco Valves	\$855.00	<a href="http://www.chromtech.net.au/pdf2/vici09_109.pdf">http://www.chromtech.net.au/pdf2/vici09_109.pdf</a>
4C10UWT	10 Ports Manual 1/16in fittings x0.75mm ports	Valco Valves	\$855.00	<a href="http://www.chromtech.net.au/pdf2/vici09_109.pdf">http://www.chromtech.net.au/pdf2/vici09_109.pdf</a>
A4C4UWT	4 Ports Manual w/standoff 1/16in fittings x0.75mm ports	Valco Valves	\$918.00	<a href="http://www.chromtech.net.au/pdf2/vici09_109.pdf">http://www.chromtech.net.au/pdf2/vici09_109.pdf</a>
A4C6UWT	6 Ports Manual w/standoff 1/16in fittings x0.75mm ports	Valco Valves	\$988.00	<a href="http://www.chromtech.net.au/pdf2/vici09_109.pdf">http://www.chromtech.net.au/pdf2/vici09_109.pdf</a>
DC10UWT	10 Ports Replacement valve 1/16in fittings x0.75mm ports	Valco Valves	\$740.00	<a href="http://www.chromtech.net.au/pdf2/vici09_109.pdf">http://www.chromtech.net.au/pdf2/vici09_109.pdf</a>
A4C10UWT	10 Ports Manual w/standoff 1/16in fittings x0.75mm ports	Valco Valves	\$1,059.00	<a href="http://www.chromtech.net.au/pdf2/vici09_109.pdf">http://www.chromtech.net.au/pdf2/vici09_109.pdf</a>
DC8UWT	8 Ports Replacement valve 1/16in fittings x0.75mm ports	Valco Valves	\$740.00	<a href="http://www.chromtech.net.au/pdf2/vici09_109.pdf">http://www.chromtech.net.au/pdf2/vici09_109.pdf</a>
E4C6UWT	6 Ports w/standard electric actuator 1/16in fittings x0.75mm ports	Valco Valves	\$1,396.00	<a href="http://www.chromtech.net.au/pdf2/vici09_109.pdf">http://www.chromtech.net.au/pdf2/vici09_109.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

E4C8UWT	8 Ports w/standard electric actuator 1/16in fittings x0.75mm ports	Valco Valves	\$1,467.00	<a href="http://www.chromtech.net.au/pdf2/vici09_109.pdf">http://www.chromtech.net.au/pdf2/vici09_109.pdf</a>
ED4C4UWT	4 Ports w/Micro electric actuator 1/16in fittings x0.75mm ports	Valco Valves	\$1,620.00	<a href="http://www.chromtech.net.au/pdf2/vici09_109.pdf">http://www.chromtech.net.au/pdf2/vici09_109.pdf</a>
ED4C8UWT	8 Ports w/Micro electric actuator 1/16in fittings x0.75mm ports	Valco Valves	\$1,760.00	<a href="http://www.chromtech.net.au/pdf2/vici09_109.pdf">http://www.chromtech.net.au/pdf2/vici09_109.pdf</a>
ED4C10UWT	10 Ports w/Micro electric actuator 1/16in fittings x0.75mm ports	Valco Valves	\$1,760.00	<a href="http://www.chromtech.net.au/pdf2/vici09_109.pdf">http://www.chromtech.net.au/pdf2/vici09_109.pdf</a>
DC4UWT	4 Ports Replacement valve 1/16in fittings x0.75mm ports	Valco Valves	\$600.00	<a href="http://www.chromtech.net.au/pdf2/vici09_109.pdf">http://www.chromtech.net.au/pdf2/vici09_109.pdf</a>
ED4C6UWT	6 Ports w/Micro electric actuator 1/16in fittings x0.75mm ports	Valco Valves	\$1,689.00	<a href="http://www.chromtech.net.au/pdf2/vici09_109.pdf">http://www.chromtech.net.au/pdf2/vici09_109.pdf</a>
A4C8UWT	8 Ports Manual w/standoff 1/16in fittings x0.75mm ports	Valco Valves	\$1,059.00	<a href="http://www.chromtech.net.au/pdf2/vici09_109.pdf">http://www.chromtech.net.au/pdf2/vici09_109.pdf</a>
E4C10UWT	10 Ports w/standard electric actuator 1/16in fittings x0.75mm ports	Valco Valves	\$1,467.00	<a href="http://www.chromtech.net.au/pdf2/vici09_109.pdf">http://www.chromtech.net.au/pdf2/vici09_109.pdf</a>
SL500CW	500ul 1/16in Stainless sample loops	Valco Valves	\$40.00	<a href="http://www.chromtech.net.au/pdf2/vici09_108.pdf">http://www.chromtech.net.au/pdf2/vici09_108.pdf</a>
SL10CW	10ul Stainless steel loops	Valco Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_108.pdf">http://www.chromtech.net.au/pdf2/vici09_108.pdf</a>
SL250CW	250ul 1/16in Stainless sample loops	Valco Valves	\$34.00	<a href="http://www.chromtech.net.au/pdf2/vici09_108.pdf">http://www.chromtech.net.au/pdf2/vici09_108.pdf</a>
SL5CW	5ul Stainless steel loops	Valco Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_108.pdf">http://www.chromtech.net.au/pdf2/vici09_108.pdf</a>
SSAC8WT	8 Ports Replacement rotor High temp 1/16in fittings x0.75mm ports	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_108.pdf">http://www.chromtech.net.au/pdf2/vici09_108.pdf</a>
SL25CW	25ul Stainless steel loops	Valco Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_108.pdf">http://www.chromtech.net.au/pdf2/vici09_108.pdf</a>
SSAC6WT	6 Ports Replacement rotor High temp 1/16in fittings x0.40mm ports	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_108.pdf">http://www.chromtech.net.au/pdf2/vici09_108.pdf</a>
SSAC10WT	10 Ports Replacement rotor High temp 1/16in fittings x0.75mm ports	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_108.pdf">http://www.chromtech.net.au/pdf2/vici09_108.pdf</a>
SL15CW	15ul Stainless steel loops	Valco Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_108.pdf">http://www.chromtech.net.au/pdf2/vici09_108.pdf</a>
SL1KCW	1ml Stainless steel loops	Valco Valves	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_108.pdf">http://www.chromtech.net.au/pdf2/vici09_108.pdf</a>
SL10KCW	10ml Stainless steel loops	Valco Valves	\$87.00	<a href="http://www.chromtech.net.au/pdf2/vici09_108.pdf">http://www.chromtech.net.au/pdf2/vici09_108.pdf</a>
SL2KCW	2ml 1/16in Stainless sample loops	Valco Valves	\$57.00	<a href="http://www.chromtech.net.au/pdf2/vici09_108.pdf">http://www.chromtech.net.au/pdf2/vici09_108.pdf</a>
SL5KCW	5ml Stainless steel loops	Valco Valves	\$68.00	<a href="http://www.chromtech.net.au/pdf2/vici09_108.pdf">http://www.chromtech.net.au/pdf2/vici09_108.pdf</a>
SL50CW	50ul Stainless steel loops	Valco Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_108.pdf">http://www.chromtech.net.au/pdf2/vici09_108.pdf</a>
SSAC4WT	4 Ports Replacement rotor High temp 1/16in fittings x0.75mm ports	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_108.pdf">http://www.chromtech.net.au/pdf2/vici09_108.pdf</a>



## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

E4C4WT	4 Ports w/standard electric actuator High temp 1/16in fittings x0.75mm ports	Valco Valves	\$1,396.00	<a href="http://www.chromtech.net.au/pdf2/vici09_108.pdf">http://www.chromtech.net.au/pdf2/vici09_108.pdf</a>
SL20CW	20ul Stainless steel loops	Valco Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_108.pdf">http://www.chromtech.net.au/pdf2/vici09_108.pdf</a>
E4C10WT	10 Ports w/standard electric actuator High temp 1/16in fittings x0.75mm ports	Valco Valves	\$1,536.00	<a href="http://www.chromtech.net.au/pdf2/vici09_108.pdf">http://www.chromtech.net.au/pdf2/vici09_108.pdf</a>
A4C6WT	6 Ports Manual w/standoff High temp 1/16in fittings x0.75mm ports	Valco Valves	\$1,059.00	<a href="http://www.chromtech.net.au/pdf2/vici09_108.pdf">http://www.chromtech.net.au/pdf2/vici09_108.pdf</a>
4C4WT	4 Ports Manual w/standoff High temp 1/16in fittings x0.75mm ports	Valco Valves	\$784.00	<a href="http://www.chromtech.net.au/pdf2/vici09_108.pdf">http://www.chromtech.net.au/pdf2/vici09_108.pdf</a>
4C6WT	6 Ports Manual High temp 1/16in fittings x0.75mm ports	Valco Valves	\$855.00	<a href="http://www.chromtech.net.au/pdf2/vici09_108.pdf">http://www.chromtech.net.au/pdf2/vici09_108.pdf</a>
4C8WT	8 Ports Manual High temp 1/16in fittings x0.75mm ports	Valco Valves	\$924.00	<a href="http://www.chromtech.net.au/pdf2/vici09_108.pdf">http://www.chromtech.net.au/pdf2/vici09_108.pdf</a>
4C10WT	10 Ports Manual High temp 1/16in fittings x0.75mm ports	Valco Valves	\$924.00	<a href="http://www.chromtech.net.au/pdf2/vici09_108.pdf">http://www.chromtech.net.au/pdf2/vici09_108.pdf</a>
A4C4WT	4 Ports Manual w/standoff High temp 1/16in fittings x0.75mm ports	Valco Valves	\$988.00	<a href="http://www.chromtech.net.au/pdf2/vici09_108.pdf">http://www.chromtech.net.au/pdf2/vici09_108.pdf</a>
E4C8WT	8 Ports w/standard electric actuator High temp 1/16in fittings x0.75mm ports	Valco Valves	\$1,536.00	<a href="http://www.chromtech.net.au/pdf2/vici09_108.pdf">http://www.chromtech.net.au/pdf2/vici09_108.pdf</a>
E4C6WT	6 Ports w/standard electric actuator High temp 1/16in fittings x0.75mm ports	Valco Valves	\$1,467.00	<a href="http://www.chromtech.net.au/pdf2/vici09_108.pdf">http://www.chromtech.net.au/pdf2/vici09_108.pdf</a>
DC10WT	10 Ports Replacement valve High temp 1/16in fittings x0.75mm ports	Valco Valves	\$810.00	<a href="http://www.chromtech.net.au/pdf2/vici09_108.pdf">http://www.chromtech.net.au/pdf2/vici09_108.pdf</a>
EH4C4WT	4 Ports w/microelectric actuator High temp 1/16in fittings x0.75mm ports	Valco Valves	\$1,612.00	<a href="http://www.chromtech.net.au/pdf2/vici09_108.pdf">http://www.chromtech.net.au/pdf2/vici09_108.pdf</a>
EH4C6WT	6 Ports w/microelectric actuator High temp 1/16in fittings x0.75mm ports	Valco Valves	\$1,683.00	<a href="http://www.chromtech.net.au/pdf2/vici09_108.pdf">http://www.chromtech.net.au/pdf2/vici09_108.pdf</a>
EH4C8WT	8 Ports w/microelectric actuator High temp 1/16in fittings x0.75mm ports	Valco Valves	\$1,754.00	<a href="http://www.chromtech.net.au/pdf2/vici09_108.pdf">http://www.chromtech.net.au/pdf2/vici09_108.pdf</a>
EH4C10WT	10 Ports w/microelectric actuator High temp 1/16in fittings x0.75mm ports	Valco Valves	\$1,754.00	<a href="http://www.chromtech.net.au/pdf2/vici09_108.pdf">http://www.chromtech.net.au/pdf2/vici09_108.pdf</a>
DC4WT	4 Ports Replacement valve High temp 1/16in fittings x0.75mm ports	Valco Valves	\$669.00	<a href="http://www.chromtech.net.au/pdf2/vici09_108.pdf">http://www.chromtech.net.au/pdf2/vici09_108.pdf</a>
DC6WT	6Ports Replacement valve High temp 1/16in fittings x0.75mm ports	Valco Valves	\$740.00	<a href="http://www.chromtech.net.au/pdf2/vici09_108.pdf">http://www.chromtech.net.au/pdf2/vici09_108.pdf</a>
DC8WT	8 Ports Replacement valve High temp 1/16in fittings x0.75mm ports	Valco Valves	\$810.00	<a href="http://www.chromtech.net.au/pdf2/vici09_108.pdf">http://www.chromtech.net.au/pdf2/vici09_108.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

A4C8WT	8 Ports Manual w/standoff High temp 1/16in fittings x0.75mm ports	Valco Valves	\$1,128.00	<a href="http://www.chromtech.net.au/pdf2/vici09_108.pdf">http://www.chromtech.net.au/pdf2/vici09_108.pdf</a>
A4C10WT	10 Ports Manual w/standoff High temp 1/16in fittings x0.75mm ports	Valco Valves	\$1,128.00	<a href="http://www.chromtech.net.au/pdf2/vici09_108.pdf">http://www.chromtech.net.au/pdf2/vici09_108.pdf</a>
SL100CUW	100ul Stainless steel loops	Valco Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_107.pdf">http://www.chromtech.net.au/pdf2/vici09_107.pdf</a>
SL10KCUW	10ml Stainless steel loops	Valco Valves	\$87.00	<a href="http://www.chromtech.net.au/pdf2/vici09_107.pdf">http://www.chromtech.net.au/pdf2/vici09_107.pdf</a>
DC8UWE	8 Ports Replacement valve 1/16in fittings x0.75mm ports	Valco Valves	\$740.00	<a href="http://www.chromtech.net.au/pdf2/vici09_107.pdf">http://www.chromtech.net.au/pdf2/vici09_107.pdf</a>
SSAC4UWE	4 Ports Replacement rotor 1/16in fittings x0.75mm ports	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_107.pdf">http://www.chromtech.net.au/pdf2/vici09_107.pdf</a>
SSAC6UWE	8 Ports Replacement rotor 1/16in fittings x0.75mm ports	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_107.pdf">http://www.chromtech.net.au/pdf2/vici09_107.pdf</a>
SSAC8UWE	8 Ports Replacement rotor 1/16in fittings x0.75mm ports	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_107.pdf">http://www.chromtech.net.au/pdf2/vici09_107.pdf</a>
SSAC10UWE	10 Ports Replacement rotor 1/16in fittings x0.75mm ports	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_107.pdf">http://www.chromtech.net.au/pdf2/vici09_107.pdf</a>
SL5CUW	5ul 1/16in stainless steel loop	Valco Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_107.pdf">http://www.chromtech.net.au/pdf2/vici09_107.pdf</a>
DC10UWE	10 Ports Replacement valve 1/16in fittings x0.75mm ports	Valco Valves	\$740.00	<a href="http://www.chromtech.net.au/pdf2/vici09_107.pdf">http://www.chromtech.net.au/pdf2/vici09_107.pdf</a>
SL10CUW	10ul Stainless steel loops	Valco Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_107.pdf">http://www.chromtech.net.au/pdf2/vici09_107.pdf</a>
SL250CUW	250ul Stainless steel loops	Valco Valves	\$34.00	<a href="http://www.chromtech.net.au/pdf2/vici09_107.pdf">http://www.chromtech.net.au/pdf2/vici09_107.pdf</a>
SL15CUW	15ul Stainless steel loops	Valco Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_107.pdf">http://www.chromtech.net.au/pdf2/vici09_107.pdf</a>
SL500CUW	500ul Stainless steel loops	Valco Valves	\$40.00	<a href="http://www.chromtech.net.au/pdf2/vici09_107.pdf">http://www.chromtech.net.au/pdf2/vici09_107.pdf</a>
SL20CUW	20ul Stainless steel loops	Valco Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_107.pdf">http://www.chromtech.net.au/pdf2/vici09_107.pdf</a>
SL1KCUW	1ml 1/16in stainless steel loop	Valco Valves	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_107.pdf">http://www.chromtech.net.au/pdf2/vici09_107.pdf</a>
SL25CUW	25ul 1/16in stainless steel loop	Valco Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_107.pdf">http://www.chromtech.net.au/pdf2/vici09_107.pdf</a>
SL2KCUW	2ml Stainless steel loops	Valco Valves	\$57.00	<a href="http://www.chromtech.net.au/pdf2/vici09_107.pdf">http://www.chromtech.net.au/pdf2/vici09_107.pdf</a>
SL5KCUW	5ml Stainless steel loops	Valco Valves	\$68.00	<a href="http://www.chromtech.net.au/pdf2/vici09_107.pdf">http://www.chromtech.net.au/pdf2/vici09_107.pdf</a>
DC6UWE	6Ports Replacement valve 1/16in fittings x0.75mm ports	Valco Valves	\$669.00	<a href="http://www.chromtech.net.au/pdf2/vici09_107.pdf">http://www.chromtech.net.au/pdf2/vici09_107.pdf</a>
ED4C4UWE	4 Ports -w/Micro electric actuator 1/16in fittings x0.75mm ports	Valco Valves	\$1,620.00	<a href="http://www.chromtech.net.au/pdf2/vici09_107.pdf">http://www.chromtech.net.au/pdf2/vici09_107.pdf</a>
SL50CUW	50ul Stainless steel loops	Valco Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_107.pdf">http://www.chromtech.net.au/pdf2/vici09_107.pdf</a>
C6UWE	6 Ports Manual Medium temp 1/16in fittings x0.75mm ports	Valco Valves	\$726.00	<a href="http://www.chromtech.net.au/pdf2/vici09_107.pdf">http://www.chromtech.net.au/pdf2/vici09_107.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

ED4C8UWE	8 Ports -w/Micro electric actuator 1/16infittings x0.75mm ports	Valco Valves	\$1,760.00	<a href="http://www.chromtech.net.au/pdf2/vici09_107.pdf">http://www.chromtech.net.au/pdf2/vici09_107.pdf</a>
DC4UWE	4 Ports Replacement valve 1/16infittings x0.75mm ports	Valco Valves	\$600.00	<a href="http://www.chromtech.net.au/pdf2/vici09_107.pdf">http://www.chromtech.net.au/pdf2/vici09_107.pdf</a>
C4UWE	4 Ports Manual Medium temp 1/16in fittings x0.75mm ports	Valco Valves	\$657.00	<a href="http://www.chromtech.net.au/pdf2/vici09_107.pdf">http://www.chromtech.net.au/pdf2/vici09_107.pdf</a>
C8UWE	8 Ports Manual Medium temp 1/16in fittings x0.75mm ports	Valco Valves	\$796.00	<a href="http://www.chromtech.net.au/pdf2/vici09_107.pdf">http://www.chromtech.net.au/pdf2/vici09_107.pdf</a>
C10UWE	10 Ports 1/16inx0.75mm ports Manual	Valco Valves	\$796.00	<a href="http://www.chromtech.net.au/pdf2/vici09_107.pdf">http://www.chromtech.net.au/pdf2/vici09_107.pdf</a>
4C4UWE	4 Ports 1/16inx0.75mm ports Manual w/Standoff	Valco Valves	\$714.00	<a href="http://www.chromtech.net.au/pdf2/vici09_107.pdf">http://www.chromtech.net.au/pdf2/vici09_107.pdf</a>
4C6UWE	6 Ports 1/16inx0.75mm ports Manual w/Standoff	Valco Valves	\$784.00	<a href="http://www.chromtech.net.au/pdf2/vici09_107.pdf">http://www.chromtech.net.au/pdf2/vici09_107.pdf</a>
4C8UWE	8 Ports 1/16inx0.75mm ports Manual w/Standoff	Valco Valves	\$855.00	<a href="http://www.chromtech.net.au/pdf2/vici09_107.pdf">http://www.chromtech.net.au/pdf2/vici09_107.pdf</a>
4C10UWE	10 Ports 1/16inx0.75mm ports Manual w/Standoff	Valco Valves	\$855.00	<a href="http://www.chromtech.net.au/pdf2/vici09_107.pdf">http://www.chromtech.net.au/pdf2/vici09_107.pdf</a>
E4C10UWE	10 Ports -w/standard electric actuator 1/16infittings x0.75mm ports	Valco Valves	\$1,467.00	<a href="http://www.chromtech.net.au/pdf2/vici09_107.pdf">http://www.chromtech.net.au/pdf2/vici09_107.pdf</a>
A4C6UWE	6 Ports -w/air actuator 1/16infittings x0.75mm ports	Valco Valves	\$988.00	<a href="http://www.chromtech.net.au/pdf2/vici09_107.pdf">http://www.chromtech.net.au/pdf2/vici09_107.pdf</a>
A4C8UWE	8 Ports -w/air actuator 1/16infittings x0.75mm ports	Valco Valves	\$1,059.00	<a href="http://www.chromtech.net.au/pdf2/vici09_107.pdf">http://www.chromtech.net.au/pdf2/vici09_107.pdf</a>
A4C10UWE	10 Ports -w/air actuator 1/16infittings x0.75mm ports	Valco Valves	\$1,059.00	<a href="http://www.chromtech.net.au/pdf2/vici09_107.pdf">http://www.chromtech.net.au/pdf2/vici09_107.pdf</a>
E4C4UWE	4 Ports -w/standard electric actuator 1/16infittings x0.75mm ports	Valco Valves	\$1,326.00	<a href="http://www.chromtech.net.au/pdf2/vici09_107.pdf">http://www.chromtech.net.au/pdf2/vici09_107.pdf</a>
E4C6UWE	6 Ports -w/standard electric actuator 1/16infittings x0.75mm ports	Valco Valves	\$1,396.00	<a href="http://www.chromtech.net.au/pdf2/vici09_107.pdf">http://www.chromtech.net.au/pdf2/vici09_107.pdf</a>
E4C8UWE	8 Ports -w/standard electric actuator 1/16infittings x0.75mm ports	Valco Valves	\$1,467.00	<a href="http://www.chromtech.net.au/pdf2/vici09_107.pdf">http://www.chromtech.net.au/pdf2/vici09_107.pdf</a>
ED4C6UWE	6 Ports -w/Micro electric actuator 1/16infittings x0.75mm ports	Valco Valves	\$1,689.00	<a href="http://www.chromtech.net.au/pdf2/vici09_107.pdf">http://www.chromtech.net.au/pdf2/vici09_107.pdf</a>
ED4C10UWE	10 Ports -w/Micro electric actuator 1/16infittings x0.75mm ports	Valco Valves	\$1,760.00	<a href="http://www.chromtech.net.au/pdf2/vici09_107.pdf">http://www.chromtech.net.au/pdf2/vici09_107.pdf</a>
A4C4UWE	4 Ports -w/air actuator 1/16infittings x0.75mm ports	Valco Valves	\$918.00	<a href="http://www.chromtech.net.au/pdf2/vici09_107.pdf">http://www.chromtech.net.au/pdf2/vici09_107.pdf</a>
SL1KCW	1ml Stainless steel loops	Valco Valves	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_106.pdf">http://www.chromtech.net.au/pdf2/vici09_106.pdf</a>
DC10WE	10 Ports Replacement valve 1/16infittings x0.40mm ports	Valco Valves	\$810.00	<a href="http://www.chromtech.net.au/pdf2/vici09_106.pdf">http://www.chromtech.net.au/pdf2/vici09_106.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

SSAC4WE	4 Ports Replacement rotor 1/16in fittings x0.40mm ports	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_106.pdf">http://www.chromtech.net.au/pdf2/vici09_106.pdf</a>
SSAC6WE	6 Ports Replacement rotor 1/16in fittings x0.40mm ports	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_106.pdf">http://www.chromtech.net.au/pdf2/vici09_106.pdf</a>
SSAC8WE	8 Ports Replacement rotor 1/16in fittings x0.40mm ports	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_106.pdf">http://www.chromtech.net.au/pdf2/vici09_106.pdf</a>
SSAC10WE	10 Ports Replacement rotor 1/16in fittings x0.40mm ports	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_106.pdf">http://www.chromtech.net.au/pdf2/vici09_106.pdf</a>
SL2CW	2ul Stainless steel loops	Valco Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_106.pdf">http://www.chromtech.net.au/pdf2/vici09_106.pdf</a>
SL100CW	100ul Stainless steel loops	Valco Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_106.pdf">http://www.chromtech.net.au/pdf2/vici09_106.pdf</a>
SL250CW	250ul 1/16in Stainless sample loops	Valco Valves	\$34.00	<a href="http://www.chromtech.net.au/pdf2/vici09_106.pdf">http://www.chromtech.net.au/pdf2/vici09_106.pdf</a>
DC8WE	8 Ports Replacement valve 1/16in fittings x0.40mm ports	Valco Valves	\$810.00	<a href="http://www.chromtech.net.au/pdf2/vici09_106.pdf">http://www.chromtech.net.au/pdf2/vici09_106.pdf</a>
SL15CW	15ul Stainless steel loops	Valco Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_106.pdf">http://www.chromtech.net.au/pdf2/vici09_106.pdf</a>
SL5CW	5ul Stainless steel loops	Valco Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_106.pdf">http://www.chromtech.net.au/pdf2/vici09_106.pdf</a>
SL20CW	20ul Stainless steel loops	Valco Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_106.pdf">http://www.chromtech.net.au/pdf2/vici09_106.pdf</a>
SL2KCW	2ml 1/16in Stainless sample loops	Valco Valves	\$57.00	<a href="http://www.chromtech.net.au/pdf2/vici09_106.pdf">http://www.chromtech.net.au/pdf2/vici09_106.pdf</a>
SL25CW	25ul Stainless steel loops	Valco Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_106.pdf">http://www.chromtech.net.au/pdf2/vici09_106.pdf</a>
SL5KCW	5ml Stainless steel loops	Valco Valves	\$68.00	<a href="http://www.chromtech.net.au/pdf2/vici09_106.pdf">http://www.chromtech.net.au/pdf2/vici09_106.pdf</a>
SL50CW	50ul Stainless steel loops	Valco Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_106.pdf">http://www.chromtech.net.au/pdf2/vici09_106.pdf</a>
SL10KCW	10ml Stainless steel loops	Valco Valves	\$87.00	<a href="http://www.chromtech.net.au/pdf2/vici09_106.pdf">http://www.chromtech.net.au/pdf2/vici09_106.pdf</a>
C6WE	6 Ports Manual 1/16in fittings x0.40mm ports	Valco Valves	\$796.00	<a href="http://www.chromtech.net.au/pdf2/vici09_106.pdf">http://www.chromtech.net.au/pdf2/vici09_106.pdf</a>
C4WE	4 Ports Manual with standoff 1/16in fittings x0.40mm ports	Valco Valves	\$726.00	<a href="http://www.chromtech.net.au/pdf2/vici09_106.pdf">http://www.chromtech.net.au/pdf2/vici09_106.pdf</a>
SL500CW	500ul 1/16in Stainless sample loops	Valco Valves	\$40.00	<a href="http://www.chromtech.net.au/pdf2/vici09_106.pdf">http://www.chromtech.net.au/pdf2/vici09_106.pdf</a>
A4C6WE	6 Ports With air actuator 1/16in fittings x0.40mm ports	Valco Valves	\$1,059.00	<a href="http://www.chromtech.net.au/pdf2/vici09_106.pdf">http://www.chromtech.net.au/pdf2/vici09_106.pdf</a>
C8WE	8 Ports Manual 1/16in fittings x0.40mm ports	Valco Valves	\$867.00	<a href="http://www.chromtech.net.au/pdf2/vici09_106.pdf">http://www.chromtech.net.au/pdf2/vici09_106.pdf</a>
SL10CW	10ul Stainless steel loops	Valco Valves	\$28.00	<a href="http://www.chromtech.net.au/pdf2/vici09_106.pdf">http://www.chromtech.net.au/pdf2/vici09_106.pdf</a>
DC6WE	6Ports Replacement valve 1/16in fittings x0.40mm ports	Valco Valves	\$740.00	<a href="http://www.chromtech.net.au/pdf2/vici09_106.pdf">http://www.chromtech.net.au/pdf2/vici09_106.pdf</a>
4C4WE	4 Ports anual with standoff 1/16in fittings x0.40mm ports	Valco Valves	\$784.00	<a href="http://www.chromtech.net.au/pdf2/vici09_106.pdf">http://www.chromtech.net.au/pdf2/vici09_106.pdf</a>
4C6WE	6 Ports Manual with standoff 1/16in fittings x0.40mm ports	Valco Valves	\$855.00	<a href="http://www.chromtech.net.au/pdf2/vici09_106.pdf">http://www.chromtech.net.au/pdf2/vici09_106.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

4C8WE	8 Ports Manual with standoff 1/16in fittings x0.40mm ports	Valco Valves	\$924.00	<a href="http://www.chromtech.net.au/pdf2/vici09_106.pdf">http://www.chromtech.net.au/pdf2/vici09_106.pdf</a>
A4C4WE	4 Ports With air actuator 1/16in fittings x0.40mm ports	Valco Valves	\$988.00	<a href="http://www.chromtech.net.au/pdf2/vici09_106.pdf">http://www.chromtech.net.au/pdf2/vici09_106.pdf</a>
C10WE	10 Ports Manual 1/16in fittings x0.40mm ports	Valco Valves	\$867.00	<a href="http://www.chromtech.net.au/pdf2/vici09_106.pdf">http://www.chromtech.net.au/pdf2/vici09_106.pdf</a>
A4C8WE	8 Ports With air actuator 1/16in fittings x0.40mm ports	Valco Valves	\$1,128.00	<a href="http://www.chromtech.net.au/pdf2/vici09_106.pdf">http://www.chromtech.net.au/pdf2/vici09_106.pdf</a>
A4C10WE	10 Ports With air actuator 1/16in fittings x0.40mm ports	Valco Valves	\$1,128.00	<a href="http://www.chromtech.net.au/pdf2/vici09_106.pdf">http://www.chromtech.net.au/pdf2/vici09_106.pdf</a>
EH4C6WE	6 Ports With microelectric actuator 1/16in fittings x0.40mm ports	Valco Valves	\$1,683.00	<a href="http://www.chromtech.net.au/pdf2/vici09_106.pdf">http://www.chromtech.net.au/pdf2/vici09_106.pdf</a>
DC4WE	4 Ports Replacement valve 1/16in fittings x0.40mm ports	Valco Valves	\$669.00	<a href="http://www.chromtech.net.au/pdf2/vici09_106.pdf">http://www.chromtech.net.au/pdf2/vici09_106.pdf</a>
EH4C10WE	10 Ports With microelectric actuator 1/16in fittings x0.40mm ports	Valco Valves	\$1,754.00	<a href="http://www.chromtech.net.au/pdf2/vici09_106.pdf">http://www.chromtech.net.au/pdf2/vici09_106.pdf</a>
4C10WE	10 Ports Manual with standoff 1/16in fittings x0.40mm ports	Valco Valves	\$924.00	<a href="http://www.chromtech.net.au/pdf2/vici09_106.pdf">http://www.chromtech.net.au/pdf2/vici09_106.pdf</a>
EH4C8WE	8 Ports With microelectric actuator 1/16in fittings x0.40mm ports	Valco Valves	\$1,754.00	<a href="http://www.chromtech.net.au/pdf2/vici09_106.pdf">http://www.chromtech.net.au/pdf2/vici09_106.pdf</a>
EH4C4WE	4 Ports With microelectric actuator 1/16in fittings x0.40mm ports	Valco Valves	\$1,612.00	<a href="http://www.chromtech.net.au/pdf2/vici09_106.pdf">http://www.chromtech.net.au/pdf2/vici09_106.pdf</a>
E4C10WE	10 Ports With standard electric actuator 1/16in fittings x0.40mm ports	Valco Valves	\$1,536.00	<a href="http://www.chromtech.net.au/pdf2/vici09_106.pdf">http://www.chromtech.net.au/pdf2/vici09_106.pdf</a>
E4C8WE	8 Ports With standard electric actuator 1/16in fittings x0.40mm ports	Valco Valves	\$1,536.00	<a href="http://www.chromtech.net.au/pdf2/vici09_106.pdf">http://www.chromtech.net.au/pdf2/vici09_106.pdf</a>
E4C6WE	6 Ports With standard electric actuator 1/16in fittings x0.40mm ports	Valco Valves	\$1,467.00	<a href="http://www.chromtech.net.au/pdf2/vici09_106.pdf">http://www.chromtech.net.au/pdf2/vici09_106.pdf</a>
E4C4WE	4 Ports With standard electric actuator 1/16in fittings x0.40mm ports	Valco Valves	\$1,396.00	<a href="http://www.chromtech.net.au/pdf2/vici09_106.pdf">http://www.chromtech.net.au/pdf2/vici09_106.pdf</a>
E4N6WT	6 Ports With standard electric actuator 1/32in fittings x0.25mm ports	Valco Valves	\$1,606.00	<a href="http://www.chromtech.net.au/pdf2/vici09_105.pdf">http://www.chromtech.net.au/pdf2/vici09_105.pdf</a>
E4N8WT	8 Ports With standard electric actuator 1/32in fittings x0.25mm ports	Valco Valves	\$1,677.00	<a href="http://www.chromtech.net.au/pdf2/vici09_105.pdf">http://www.chromtech.net.au/pdf2/vici09_105.pdf</a>
E4N10WT	10 Ports With standard electric actuator 1/32in fittings x0.25mm ports	Valco Valves	\$1,677.00	<a href="http://www.chromtech.net.au/pdf2/vici09_105.pdf">http://www.chromtech.net.au/pdf2/vici09_105.pdf</a>
EH4N4WT	4 Ports With microelectric actuator 1/32in fittings x0.25mm ports	Valco Valves	\$1,754.00	<a href="http://www.chromtech.net.au/pdf2/vici09_105.pdf">http://www.chromtech.net.au/pdf2/vici09_105.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

EH4N6WT	6 Ports With microelectric actuator 1/32infittings x0.25mm ports	Valco Valves	\$1,950.00	<a href="http://www.chromtech.net.au/pdf2/vici09_105.pdf">http://www.chromtech.net.au/pdf2/vici09_105.pdf</a>
EH4N10WT	10 Ports With microelectric actuator 1/32infittings x0.25mm ports	Valco Valves	\$2,020.00	<a href="http://www.chromtech.net.au/pdf2/vici09_105.pdf">http://www.chromtech.net.au/pdf2/vici09_105.pdf</a>
A4N4WT	4 Ports with air actuator 1/32infittings x0.25mm ports	Valco Valves	\$1,128.00	<a href="http://www.chromtech.net.au/pdf2/vici09_105.pdf">http://www.chromtech.net.au/pdf2/vici09_105.pdf</a>
EH4N8WT	8 Ports With microelectric actuator 1/32infittings x0.25mm ports	Valco Valves	\$2,020.00	<a href="http://www.chromtech.net.au/pdf2/vici09_105.pdf">http://www.chromtech.net.au/pdf2/vici09_105.pdf</a>
E4N4WT	4 Ports with standard electric actuator 1/32infittings x0.25mm ports	Valco Valves	\$1,536.00	<a href="http://www.chromtech.net.au/pdf2/vici09_105.pdf">http://www.chromtech.net.au/pdf2/vici09_105.pdf</a>
A4N10WT	10 Ports With air actuator 1/32infittings x0.25mm ports	Valco Valves	\$1,269.00	<a href="http://www.chromtech.net.au/pdf2/vici09_105.pdf">http://www.chromtech.net.au/pdf2/vici09_105.pdf</a>
A4N8WT	8 Ports With air actuator 1/32infittings x0.25mm ports	Valco Valves	\$1,269.00	<a href="http://www.chromtech.net.au/pdf2/vici09_105.pdf">http://www.chromtech.net.au/pdf2/vici09_105.pdf</a>
4N4WT	4 Ports Manual with standoff 1/32infittings x0.25mm ports	Valco Valves	\$924.00	<a href="http://www.chromtech.net.au/pdf2/vici09_105.pdf">http://www.chromtech.net.au/pdf2/vici09_105.pdf</a>
4N6WT	6 Ports Manual with standoff 1/32infittings x0.25mm ports	Valco Valves	\$994.00	<a href="http://www.chromtech.net.au/pdf2/vici09_105.pdf">http://www.chromtech.net.au/pdf2/vici09_105.pdf</a>
DN4WT	4 Ports Replacement valve 1/32infittings x0.25mm ports	Valco Valves	\$810.00	<a href="http://www.chromtech.net.au/pdf2/vici09_105.pdf">http://www.chromtech.net.au/pdf2/vici09_105.pdf</a>
4N10WT	10 Ports Manual with standoff 1/32infittings x0.25mm ports	Valco Valves	\$1,065.00	<a href="http://www.chromtech.net.au/pdf2/vici09_105.pdf">http://www.chromtech.net.au/pdf2/vici09_105.pdf</a>
SL5NW	5ul 1/32infittings Stainless steel loops	Valco Valves	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_105.pdf">http://www.chromtech.net.au/pdf2/vici09_105.pdf</a>
4N8WT	8 Ports Manual with standoff 1/32infittings x0.25mm ports	Valco Valves	\$1,065.00	<a href="http://www.chromtech.net.au/pdf2/vici09_105.pdf">http://www.chromtech.net.au/pdf2/vici09_105.pdf</a>
SL25NW	25ul 1/32infittings Stainless steel loops	Valco Valves	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_105.pdf">http://www.chromtech.net.au/pdf2/vici09_105.pdf</a>
SL500NW	500ul Stainless steel loops	Valco Valves	\$48.00	<a href="http://www.chromtech.net.au/pdf2/vici09_105.pdf">http://www.chromtech.net.au/pdf2/vici09_105.pdf</a>
A4N6WT	6 Ports With air actuator 1/32infittings x0.25mm ports	Valco Valves	\$1,198.00	<a href="http://www.chromtech.net.au/pdf2/vici09_105.pdf">http://www.chromtech.net.au/pdf2/vici09_105.pdf</a>
SL20NW	20ul 1/32infittings Stainless steel loops	Valco Valves	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_105.pdf">http://www.chromtech.net.au/pdf2/vici09_105.pdf</a>
SL250NW	250ul 1/32infittings Stainless steel loops	Valco Valves	\$40.00	<a href="http://www.chromtech.net.au/pdf2/vici09_105.pdf">http://www.chromtech.net.au/pdf2/vici09_105.pdf</a>
SL15NW	15ul 1/32infittings Stainless steel loops	Valco Valves	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_105.pdf">http://www.chromtech.net.au/pdf2/vici09_105.pdf</a>
SL100NW	100ul 1/32infittings Stainless steel loops	Valco Valves	\$34.00	<a href="http://www.chromtech.net.au/pdf2/vici09_105.pdf">http://www.chromtech.net.au/pdf2/vici09_105.pdf</a>
SL2NW	2ul 1/32infittings Stainless steel loops	Valco Valves	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_105.pdf">http://www.chromtech.net.au/pdf2/vici09_105.pdf</a>
SL50NW	50ul 1/32infittings Stainless steel loops	Valco Valves	\$34.00	<a href="http://www.chromtech.net.au/pdf2/vici09_105.pdf">http://www.chromtech.net.au/pdf2/vici09_105.pdf</a>
DN6WT	6 Ports Replacement valve 1/32infittings x0.25mm ports	Valco Valves	\$879.00	<a href="http://www.chromtech.net.au/pdf2/vici09_105.pdf">http://www.chromtech.net.au/pdf2/vici09_105.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

SSAN10WT	10 Ports Replacement rotor 1/32in fittings x0.25mm ports	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_105.pdf">http://www.chromtech.net.au/pdf2/vici09_105.pdf</a>
SSAN8WT	8 Ports Replacement rotor 1/32in fittings x0.25mm ports	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_105.pdf">http://www.chromtech.net.au/pdf2/vici09_105.pdf</a>
SSAN6WT	6 Ports Replacement rotor 1/32in fittings x0.25mm ports	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_105.pdf">http://www.chromtech.net.au/pdf2/vici09_105.pdf</a>
SSAN4WT	4 Ports Replacement rotor 1/32in fittings x0.25mm ports	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_105.pdf">http://www.chromtech.net.au/pdf2/vici09_105.pdf</a>
DN10WT	10 Ports Replacement valve 1/32in fittings x0.25mm ports	Valco Valves	\$950.00	<a href="http://www.chromtech.net.au/pdf2/vici09_105.pdf">http://www.chromtech.net.au/pdf2/vici09_105.pdf</a>
DN8WT	8 Ports Replacement valve 1/32in fittings x0.25mm ports	Valco Valves	\$950.00	<a href="http://www.chromtech.net.au/pdf2/vici09_105.pdf">http://www.chromtech.net.au/pdf2/vici09_105.pdf</a>
SL10NW	10ul 1/32in fittings Stainless steel loops	Valco Valves	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_105.pdf">http://www.chromtech.net.au/pdf2/vici09_105.pdf</a>
SL50NW	50ul 1/32in fittings Stainless steel loops	Valco Valves	\$34.00	<a href="http://www.chromtech.net.au/pdf2/vici09_104.pdf">http://www.chromtech.net.au/pdf2/vici09_104.pdf</a>
DN8WE	8 Ports Replacement valve 1/32in fittings x0.25mm ports	Valco Valves	\$950.00	<a href="http://www.chromtech.net.au/pdf2/vici09_104.pdf">http://www.chromtech.net.au/pdf2/vici09_104.pdf</a>
DN10WE	10 Ports Replacement valve 1/32in fittings x0.25mm ports	Valco Valves	\$950.00	<a href="http://www.chromtech.net.au/pdf2/vici09_104.pdf">http://www.chromtech.net.au/pdf2/vici09_104.pdf</a>
SSAN4WE	4 Ports Replacement rotor 1/32in fittings x0.25mm ports	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_104.pdf">http://www.chromtech.net.au/pdf2/vici09_104.pdf</a>
SSAN6WE	6 Ports Replacement rotor 1/32in fittings x0.25mm ports	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_104.pdf">http://www.chromtech.net.au/pdf2/vici09_104.pdf</a>
SSAN8WE	8 Ports Replacement rotor 1/32in fittings x0.25mm ports	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_104.pdf">http://www.chromtech.net.au/pdf2/vici09_104.pdf</a>
SSAN10WE	10 Ports Replacement rotor 1/32in fittings x0.25mm ports	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_104.pdf">http://www.chromtech.net.au/pdf2/vici09_104.pdf</a>
SL2NW	2ul 1/32in fittings Stainless steel loops	Valco Valves	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_104.pdf">http://www.chromtech.net.au/pdf2/vici09_104.pdf</a>
SL25NW	25ul 1/32in fittings Stainless steel loops	Valco Valves	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_104.pdf">http://www.chromtech.net.au/pdf2/vici09_104.pdf</a>
SL5NW	5ul 1/32in fittings Stainless steel loops	Valco Valves	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_104.pdf">http://www.chromtech.net.au/pdf2/vici09_104.pdf</a>
SL15NW	15ul 1/32in fittings Stainless steel loops	Valco Valves	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_104.pdf">http://www.chromtech.net.au/pdf2/vici09_104.pdf</a>
SL500NW	500ul Stainless steel loops	Valco Valves	\$48.00	<a href="http://www.chromtech.net.au/pdf2/vici09_104.pdf">http://www.chromtech.net.au/pdf2/vici09_104.pdf</a>
SL100NW	100ul 1/32in fittings Stainless steel loops	Valco Valves	\$34.00	<a href="http://www.chromtech.net.au/pdf2/vici09_104.pdf">http://www.chromtech.net.au/pdf2/vici09_104.pdf</a>
SL250NW	250ul 1/32in fittings Stainless steel loops	Valco Valves	\$40.00	<a href="http://www.chromtech.net.au/pdf2/vici09_104.pdf">http://www.chromtech.net.au/pdf2/vici09_104.pdf</a>
SL20NW	20ul 1/32in fittings Stainless steel loops	Valco Valves	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_104.pdf">http://www.chromtech.net.au/pdf2/vici09_104.pdf</a>
DN6WE	6 Ports Replacement valve 1/32in fittings x0.25mm ports	Valco Valves	\$879.00	<a href="http://www.chromtech.net.au/pdf2/vici09_104.pdf">http://www.chromtech.net.au/pdf2/vici09_104.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

SL10NW	10ul 1/32infittings Stainless steel loops	Valco Valves	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_104.pdf">http://www.chromtech.net.au/pdf2/vici09_104.pdf</a>
A4N4WE	4 Ports With air actuator 1/32infittings x0.25mm ports	Valco Valves	\$1,128.00	<a href="http://www.chromtech.net.au/pdf2/vici09_104.pdf">http://www.chromtech.net.au/pdf2/vici09_104.pdf</a>
E4N4WE	4 Ports With air actuator 1/32infittings x0.25mm ports	Valco Valves	\$1,536.00	<a href="http://www.chromtech.net.au/pdf2/vici09_104.pdf">http://www.chromtech.net.au/pdf2/vici09_104.pdf</a>
4N4WE	4 Ports Manual with standoff 1/32infittings x0.25mm ports	Valco Valves	\$924.00	<a href="http://www.chromtech.net.au/pdf2/vici09_104.pdf">http://www.chromtech.net.au/pdf2/vici09_104.pdf</a>
4N6WE	6 Ports Manual with standoff 1/32infittings x0.25mm ports	Valco Valves	\$994.00	<a href="http://www.chromtech.net.au/pdf2/vici09_104.pdf">http://www.chromtech.net.au/pdf2/vici09_104.pdf</a>
DN4WE	4 Ports Replacement valve 1/32infittings x0.25mm ports	Valco Valves	\$810.00	<a href="http://www.chromtech.net.au/pdf2/vici09_104.pdf">http://www.chromtech.net.au/pdf2/vici09_104.pdf</a>
4N10WE	10 Ports Manual with standoff 1/32infittings x0.25mm ports	Valco Valves	\$1,065.00	<a href="http://www.chromtech.net.au/pdf2/vici09_104.pdf">http://www.chromtech.net.au/pdf2/vici09_104.pdf</a>
A4N6WE	6 Ports With air actuator 1/32infittings x0.25mm ports	Valco Valves	\$1,198.00	<a href="http://www.chromtech.net.au/pdf2/vici09_104.pdf">http://www.chromtech.net.au/pdf2/vici09_104.pdf</a>
A4N8WE	8 Ports With air actuator 1/32infittings x0.25mm ports	Valco Valves	\$1,269.00	<a href="http://www.chromtech.net.au/pdf2/vici09_104.pdf">http://www.chromtech.net.au/pdf2/vici09_104.pdf</a>
EH4N8WE	8 Ports With microelectric actuator 1/32infittings x0.25mm ports	Valco Valves	\$2,020.00	<a href="http://www.chromtech.net.au/pdf2/vici09_104.pdf">http://www.chromtech.net.au/pdf2/vici09_104.pdf</a>
E4N6WE	6 Ports With standard electric actuator 1/32infittings x0.25mm ports	Valco Valves	\$1,606.00	<a href="http://www.chromtech.net.au/pdf2/vici09_104.pdf">http://www.chromtech.net.au/pdf2/vici09_104.pdf</a>
E4N8WE	8 Ports With standard electric actuator 1/32infittings x0.25mm ports	Valco Valves	\$1,677.00	<a href="http://www.chromtech.net.au/pdf2/vici09_104.pdf">http://www.chromtech.net.au/pdf2/vici09_104.pdf</a>
E4N10WE	10 Ports With standard electric actuator 1/32infittings x0.25mm ports	Valco Valves	\$1,677.00	<a href="http://www.chromtech.net.au/pdf2/vici09_104.pdf">http://www.chromtech.net.au/pdf2/vici09_104.pdf</a>
EH4N4WE	4 Ports With microelectric actuator 1/32infittings x0.25mm ports	Valco Valves	\$1,754.00	<a href="http://www.chromtech.net.au/pdf2/vici09_104.pdf">http://www.chromtech.net.au/pdf2/vici09_104.pdf</a>
EH4N6WE	6 Ports With microelectric actuator 1/32infittings x0.25mm ports	Valco Valves	\$1,950.00	<a href="http://www.chromtech.net.au/pdf2/vici09_104.pdf">http://www.chromtech.net.au/pdf2/vici09_104.pdf</a>
A4N10WE	10 Ports With air actuator 1/32infittings x0.25mm ports	Valco Valves	\$1,269.00	<a href="http://www.chromtech.net.au/pdf2/vici09_104.pdf">http://www.chromtech.net.au/pdf2/vici09_104.pdf</a>
EH4N10WE	10 Ports With microelectric actuator 1/32infittings x0.25mm ports	Valco Valves	\$2,020.00	<a href="http://www.chromtech.net.au/pdf2/vici09_104.pdf">http://www.chromtech.net.au/pdf2/vici09_104.pdf</a>
4N8WE	8 Ports Manual with standoff 1/32infittings x0.25mm ports	Valco Valves	\$1,065.00	<a href="http://www.chromtech.net.au/pdf2/vici09_104.pdf">http://www.chromtech.net.au/pdf2/vici09_104.pdf</a>
E2CI4UWE1	1/16inx0.75mm - 1/16inx0.75mm -0.5ul w/std electric actuator	Valco Valves	\$1,389.00	<a href="http://www.chromtech.net.au/pdf2/vici09_103.pdf">http://www.chromtech.net.au/pdf2/vici09_103.pdf</a>



## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

ED2I4UWE.5	1/8inx0.75mm -0. 5ul w/microelectric actuator	Valco Valves	\$1,683.00	<a href="http://www.chromtech.net.au/pdf2/vici09_103.pdf">http://www.chromtech.net.au/pdf2/vici09_103.pdf</a>
ED2I4UWE1	1/8inx0.75mm -1ul w/microelectric actuator	Valco Valves	\$1,683.00	<a href="http://www.chromtech.net.au/pdf2/vici09_103.pdf">http://www.chromtech.net.au/pdf2/vici09_103.pdf</a>
ED2I4UWE2	1/8inx0.75mm -2ul w/microelectric actuator	Valco Valves	\$1,683.00	<a href="http://www.chromtech.net.au/pdf2/vici09_103.pdf">http://www.chromtech.net.au/pdf2/vici09_103.pdf</a>
DI4UWE.2	1/8inx0.75mm - 0.2ul Replacement valve	Valco Valves	\$663.00	<a href="http://www.chromtech.net.au/pdf2/vici09_103.pdf">http://www.chromtech.net.au/pdf2/vici09_103.pdf</a>
DI4UWE.5	1/8inx0.75mm - 0.5ul Replacement valve	Valco Valves	\$663.00	<a href="http://www.chromtech.net.au/pdf2/vici09_103.pdf">http://www.chromtech.net.au/pdf2/vici09_103.pdf</a>
DI4UWE1	1/8inx0.75mm - 1ul Replacement valve	Valco Valves	\$663.00	<a href="http://www.chromtech.net.au/pdf2/vici09_103.pdf">http://www.chromtech.net.au/pdf2/vici09_103.pdf</a>
2CI4UWE.2	1/16inx0.75mm - 1/16inx0.75mm -0.1ul Manual w/standoff	Valco Valves	\$777.00	<a href="http://www.chromtech.net.au/pdf2/vici09_103.pdf">http://www.chromtech.net.au/pdf2/vici09_103.pdf</a>
SSAI4UWE.2	1/8inx0.75mm - 0.2ul Replacement rotor	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_103.pdf">http://www.chromtech.net.au/pdf2/vici09_103.pdf</a>
SSAI4UWE1	1/8inx0.75mm - 1ul Replacement rotor	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_103.pdf">http://www.chromtech.net.au/pdf2/vici09_103.pdf</a>
SSAI4UWE2	1/8inx0.75mm - 2ul Replacement rotor	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_103.pdf">http://www.chromtech.net.au/pdf2/vici09_103.pdf</a>
ED2CI4UWE.5	1/16inx0.75mm - 0.5ul w/microelectric actuator	Valco Valves	\$1,683.00	<a href="http://www.chromtech.net.au/pdf2/vici09_103.pdf">http://www.chromtech.net.au/pdf2/vici09_103.pdf</a>
ED2I4UWE.2	1/8inx0.75mm -0. 2ul w/microelectric actuator	Valco Valves	\$1,683.00	<a href="http://www.chromtech.net.au/pdf2/vici09_103.pdf">http://www.chromtech.net.au/pdf2/vici09_103.pdf</a>
E2CI4UWE2	1/16inx0.75mm - 1/16inx0.75mm -1ul w/std electric actuator	Valco Valves	\$1,389.00	<a href="http://www.chromtech.net.au/pdf2/vici09_103.pdf">http://www.chromtech.net.au/pdf2/vici09_103.pdf</a>
DI4UWE2	1/8inx0.75mm - 2ul Replacement valve	Valco Valves	\$663.00	<a href="http://www.chromtech.net.au/pdf2/vici09_103.pdf">http://www.chromtech.net.au/pdf2/vici09_103.pdf</a>
E2CI4UWE.5	1/16inx0.75mm - 1/16inx0.75mm -0.2ul w/std electric actuator	Valco Valves	\$1,389.00	<a href="http://www.chromtech.net.au/pdf2/vici09_103.pdf">http://www.chromtech.net.au/pdf2/vici09_103.pdf</a>
E2CI4UWE.2	1/16inx0.75mm - 1/16inx0.75mm -0.1ul w/std electric actuator	Valco Valves	\$1,389.00	<a href="http://www.chromtech.net.au/pdf2/vici09_103.pdf">http://www.chromtech.net.au/pdf2/vici09_103.pdf</a>
A2CI4UWE2	1/16inx0.75mm - 1/16inx0.75mm -1ul w/air actuator	Valco Valves	\$981.00	<a href="http://www.chromtech.net.au/pdf2/vici09_103.pdf">http://www.chromtech.net.au/pdf2/vici09_103.pdf</a>
A2CI4UWE1	1/16inx0.75mm - 1/16inx0.75mm -0.5ul w/air actuator	Valco Valves	\$981.00	<a href="http://www.chromtech.net.au/pdf2/vici09_103.pdf">http://www.chromtech.net.au/pdf2/vici09_103.pdf</a>
A2CI4UWE.5	1/16inx0.75mm - 1/16inx0.75mm -0.2ul w/air actuator	Valco Valves	\$981.00	<a href="http://www.chromtech.net.au/pdf2/vici09_103.pdf">http://www.chromtech.net.au/pdf2/vici09_103.pdf</a>
A2CI4UWE.2	1/16inx0.75mm - 1/16inx0.75mm -0.1ul w/air actuator	Valco Valves	\$981.00	<a href="http://www.chromtech.net.au/pdf2/vici09_103.pdf">http://www.chromtech.net.au/pdf2/vici09_103.pdf</a>
2CI4UWE2	1/16inx0.75mm - 1/16inx0.75mm -1ul Manual w/standoff	Valco Valves	\$777.00	<a href="http://www.chromtech.net.au/pdf2/vici09_103.pdf">http://www.chromtech.net.au/pdf2/vici09_103.pdf</a>
2CI4UWE1	1/16inx0.75mm - 1/16inx0.75mm -0.5ul Manual w/standoff	Valco Valves	\$777.00	<a href="http://www.chromtech.net.au/pdf2/vici09_103.pdf">http://www.chromtech.net.au/pdf2/vici09_103.pdf</a>
2CI4UWE.5	1/16inx0.75mm - 1/16inx0.75mm -0.2ul Manual w/standoff	Valco Valves	\$777.00	<a href="http://www.chromtech.net.au/pdf2/vici09_103.pdf">http://www.chromtech.net.au/pdf2/vici09_103.pdf</a>
E2I4UWE1	1/8inx0.75mm - 1ul w/std electric actuator	Valco Valves	\$1,389.00	<a href="http://www.chromtech.net.au/pdf2/vici09_103.pdf">http://www.chromtech.net.au/pdf2/vici09_103.pdf</a>
ED2CI4UWE.2	1/16inx0.75mm - 0.2ul w/microelectric actuator	Valco Valves	\$1,683.00	<a href="http://www.chromtech.net.au/pdf2/vici09_103.pdf">http://www.chromtech.net.au/pdf2/vici09_103.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

I4UWE.2	1/8inx0.75mm - 0.2ul Manual	Valco Valves	\$720.00	<a href="http://www.chromtech.net.au/pdf2/vici09_103.pdf">http://www.chromtech.net.au/pdf2/vici09_103.pdf</a>
E2I4UWE.5	1/8inx0.75mm - 0.5ul w/std electric actuator	Valco Valves	\$1,389.00	<a href="http://www.chromtech.net.au/pdf2/vici09_103.pdf">http://www.chromtech.net.au/pdf2/vici09_103.pdf</a>
E2I4UWE.2	1/8inx0.75mm - 0.2ul w/std electric actuator	Valco Valves	\$1,389.00	<a href="http://www.chromtech.net.au/pdf2/vici09_103.pdf">http://www.chromtech.net.au/pdf2/vici09_103.pdf</a>
A2I4UWE2	1/8inx0.75mm - 2ul w/air actuator	Valco Valves	\$981.00	<a href="http://www.chromtech.net.au/pdf2/vici09_103.pdf">http://www.chromtech.net.au/pdf2/vici09_103.pdf</a>
A2I4UWE1	1/8inx0.75mm - 1ul w/air actuator	Valco Valves	\$981.00	<a href="http://www.chromtech.net.au/pdf2/vici09_103.pdf">http://www.chromtech.net.au/pdf2/vici09_103.pdf</a>
A2I4UWE.5	1/8inx0.75mm - 0.5ul w/air actuator	Valco Valves	\$981.00	<a href="http://www.chromtech.net.au/pdf2/vici09_103.pdf">http://www.chromtech.net.au/pdf2/vici09_103.pdf</a>
A2I4UWE.2	1/8inx0.75mm - 0.2ul w/air actuator	Valco Valves	\$981.00	<a href="http://www.chromtech.net.au/pdf2/vici09_103.pdf">http://www.chromtech.net.au/pdf2/vici09_103.pdf</a>
2I4UWE2	1/8inx0.75mm - 2ul Manual w/standoff	Valco Valves	\$777.00	<a href="http://www.chromtech.net.au/pdf2/vici09_103.pdf">http://www.chromtech.net.au/pdf2/vici09_103.pdf</a>
2I4UWE1	1/8inx0.75mm - 1ul Manual w/standoff	Valco Valves	\$777.00	<a href="http://www.chromtech.net.au/pdf2/vici09_103.pdf">http://www.chromtech.net.au/pdf2/vici09_103.pdf</a>
2I4UWE.5	1/8inx0.75mm - 0.5ul Manual w/standoff	Valco Valves	\$777.00	<a href="http://www.chromtech.net.au/pdf2/vici09_103.pdf">http://www.chromtech.net.au/pdf2/vici09_103.pdf</a>
2I4UWE.2	1/8inx0.75mm - 0.2ul Manual w/standoff	Valco Valves	\$777.00	<a href="http://www.chromtech.net.au/pdf2/vici09_103.pdf">http://www.chromtech.net.au/pdf2/vici09_103.pdf</a>
I4UWE2	1/8inx0.75mm - 2ul Manual	Valco Valves	\$720.00	<a href="http://www.chromtech.net.au/pdf2/vici09_103.pdf">http://www.chromtech.net.au/pdf2/vici09_103.pdf</a>
SSAI4UWE.5	1/8inx0.75mm - 0.5ul Replacement rotor	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_103.pdf">http://www.chromtech.net.au/pdf2/vici09_103.pdf</a>
I4UWE.5	1/8inx0.75mm - 0.5ul Manual	Valco Valves	\$720.00	<a href="http://www.chromtech.net.au/pdf2/vici09_103.pdf">http://www.chromtech.net.au/pdf2/vici09_103.pdf</a>
SSACI4UWE2	1/16inx0.75mm - 1/16inx0.75mm -1ul Replacement rotor	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_103.pdf">http://www.chromtech.net.au/pdf2/vici09_103.pdf</a>
DCI4UWE1	1/16inx0.75mm - 1/16inx0.75mm -1ul Replacement valve	Valco Valves	\$663.00	<a href="http://www.chromtech.net.au/pdf2/vici09_103.pdf">http://www.chromtech.net.au/pdf2/vici09_103.pdf</a>
E2I4UWE2	1/8inx0.75mm - 2ul w/std electric actuator	Valco Valves	\$1,389.00	<a href="http://www.chromtech.net.au/pdf2/vici09_103.pdf">http://www.chromtech.net.au/pdf2/vici09_103.pdf</a>
ED2CI4UWE1	1/16inx0.75mm - 1ul w/microelectric actuator	Valco Valves	\$1,683.00	<a href="http://www.chromtech.net.au/pdf2/vici09_103.pdf">http://www.chromtech.net.au/pdf2/vici09_103.pdf</a>
ED2CI4UWE2	1/16inx0.75mm - 2ul w/microelectric actuator	Valco Valves	\$1,683.00	<a href="http://www.chromtech.net.au/pdf2/vici09_103.pdf">http://www.chromtech.net.au/pdf2/vici09_103.pdf</a>
I4UWE1	1/8inx0.75mm - 1ul Manual	Valco Valves	\$720.00	<a href="http://www.chromtech.net.au/pdf2/vici09_103.pdf">http://www.chromtech.net.au/pdf2/vici09_103.pdf</a>
DCI4UWE.5	1/16inx0.75mm - 1/16inx0.75mm -0.5ul Replacement valve	Valco Valves	\$663.00	<a href="http://www.chromtech.net.au/pdf2/vici09_103.pdf">http://www.chromtech.net.au/pdf2/vici09_103.pdf</a>
DCI4UWE2	1/16inx0.75mm - 1/16inx0.75mm -2ul Replacement valve	Valco Valves	\$663.00	<a href="http://www.chromtech.net.au/pdf2/vici09_103.pdf">http://www.chromtech.net.au/pdf2/vici09_103.pdf</a>
SSACI4UWE.2	1/16inx0.75mm - 1/16inx0.75mm -0.1ul Replacement rotor	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_103.pdf">http://www.chromtech.net.au/pdf2/vici09_103.pdf</a>
SSACI4UWE.5	1/16inx0.75mm - 1/16inx0.75mm -0.2ul Replacement rotor	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_103.pdf">http://www.chromtech.net.au/pdf2/vici09_103.pdf</a>
SSACI4UWE1	1/16inx0.75mm - 1/16inx0.75mm -0.5ul Replacement rotor	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_103.pdf">http://www.chromtech.net.au/pdf2/vici09_103.pdf</a>
DCI4UWE.2	1/16inx0.75mm - 1/16inx0.75mm -0.2ul Replacement valve	Valco Valves	\$663.00	<a href="http://www.chromtech.net.au/pdf2/vici09_103.pdf">http://www.chromtech.net.au/pdf2/vici09_103.pdf</a>
CI4WE.2	.2ul Manual 1/32in fittings x0.25mm ports	Valco Valves	\$790.00	<a href="http://www.chromtech.net.au/pdf2/vici09_102.pdf">http://www.chromtech.net.au/pdf2/vici09_102.pdf</a>
CI4WE.5	.5ul Manual 1/32in fittings x0.25mm ports	Valco Valves	\$790.00	<a href="http://www.chromtech.net.au/pdf2/vici09_102.pdf">http://www.chromtech.net.au/pdf2/vici09_102.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

M2CI4WE.06		Valco Valves	\$848.00	<a href="http://www.chromtech.net.au/pdf2/vici09_102.pdf">http://www.chromtech.net.au/pdf2/vici09_102.pdf</a>
2CI4WE.1	.1ul Manual w/standoff 1/32in fittings x0.25mm ports	Valco Valves	\$848.00	<a href="http://www.chromtech.net.au/pdf2/vici09_102.pdf">http://www.chromtech.net.au/pdf2/vici09_102.pdf</a>
2CI4WE.2	.2ul Manual w/standoff 1/32in fittings x0.25mm ports	Valco Valves	\$848.00	<a href="http://www.chromtech.net.au/pdf2/vici09_102.pdf">http://www.chromtech.net.au/pdf2/vici09_102.pdf</a>
2CI4WE.5	.5ul Manual w/standoff 1/32in fittings x0.25mm ports	Valco Valves	\$848.00	<a href="http://www.chromtech.net.au/pdf2/vici09_102.pdf">http://www.chromtech.net.au/pdf2/vici09_102.pdf</a>
SSACI4WE.5	.5ul Replacement rotor 1/32in fittings x0.25mm ports	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_102.pdf">http://www.chromtech.net.au/pdf2/vici09_102.pdf</a>
SSACI4WE.1	.1ul Replacement rotor 1/32in fittings x0.25mm ports	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_102.pdf">http://www.chromtech.net.au/pdf2/vici09_102.pdf</a>
SSACI4WE.06	.06ul Replacement rotor 1/32in fittings x0.25mm ports	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_102.pdf">http://www.chromtech.net.au/pdf2/vici09_102.pdf</a>
CI4WE.1	.1ul Manual 1/32in fittings x0.25mm ports	Valco Valves	\$790.00	<a href="http://www.chromtech.net.au/pdf2/vici09_102.pdf">http://www.chromtech.net.au/pdf2/vici09_102.pdf</a>
SSACI4WE.2	.2ul Replacement rotor 1/32in fittings x0.25mm ports	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_102.pdf">http://www.chromtech.net.au/pdf2/vici09_102.pdf</a>
CI4WE.06	.06ul Manual 1/32in fittings x0.25mm ports	Valco Valves	\$790.00	<a href="http://www.chromtech.net.au/pdf2/vici09_102.pdf">http://www.chromtech.net.au/pdf2/vici09_102.pdf</a>
SSANI4WE.5	.5ul Replacement rotor 1/32in fittings x0.25mm ports	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_102.pdf">http://www.chromtech.net.au/pdf2/vici09_102.pdf</a>
SSANI4WE.2	.2ul Replacement rotor 1/32in fittings x0.25mm ports	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_102.pdf">http://www.chromtech.net.au/pdf2/vici09_102.pdf</a>
SSANI4WE.1	.1ul Replacement rotor 1/32in fittings x0.25mm ports	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_102.pdf">http://www.chromtech.net.au/pdf2/vici09_102.pdf</a>
SSANI4WE.06	.06ul Replacement rotor 1/32in fittings x0.25mm ports	Valco Valves	\$96.00	<a href="http://www.chromtech.net.au/pdf2/vici09_102.pdf">http://www.chromtech.net.au/pdf2/vici09_102.pdf</a>
DNI4WE.5	.5ul Replacement valve 1/32in fittings x0.25mm ports	Valco Valves	\$873.00	<a href="http://www.chromtech.net.au/pdf2/vici09_102.pdf">http://www.chromtech.net.au/pdf2/vici09_102.pdf</a>
DCI4WE.5	.5ul Replacement valve 1/32in fittings x0.25mm ports	Valco Valves	\$734.00	<a href="http://www.chromtech.net.au/pdf2/vici09_102.pdf">http://www.chromtech.net.au/pdf2/vici09_102.pdf</a>
DNI4WE.1	.1ul Replacement valve 1/32in fittings x0.25mm ports	Valco Valves	\$873.00	<a href="http://www.chromtech.net.au/pdf2/vici09_102.pdf">http://www.chromtech.net.au/pdf2/vici09_102.pdf</a>
E2CI4WE.2	.2ul w/std electric actuator 1/32in fittings x0.25mm ports	Valco Valves	\$1,460.00	<a href="http://www.chromtech.net.au/pdf2/vici09_102.pdf">http://www.chromtech.net.au/pdf2/vici09_102.pdf</a>
DNI4WE.06	.06ul Replacement valve 1/32in fittings x0.25mm ports	Valco Valves	\$873.00	<a href="http://www.chromtech.net.au/pdf2/vici09_102.pdf">http://www.chromtech.net.au/pdf2/vici09_102.pdf</a>
CatNo	Description	Category	UD Price 1.0	WebLink
A2NI4WE.06	.06ul w/air actuator 1/32in fittings x0.25mm ports	Valco Valves	\$1,192.00	<a href="http://www.chromtech.net.au/pdf2/vici09_102.pdf">http://www.chromtech.net.au/pdf2/vici09_102.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

EP2NI4WE.1	.1ul w/std electric actuator 1/32in fittings x0.25mm ports	Valco Valves	\$1,944.00	<a href="http://www.chromtech.net.au/pdf2/vici09_102.pdf">http://www.chromtech.net.au/pdf2/vici09_102.pdf</a>
DNI4WE.2	.2ul Replacement valve 1/32in fittings x0.25mm ports	Valco Valves	\$873.00	<a href="http://www.chromtech.net.au/pdf2/vici09_102.pdf">http://www.chromtech.net.au/pdf2/vici09_102.pdf</a>
A2CI4WE.2	.2ul w/air actuator 1/32in fittings x0.25mm ports	Valco Valves	\$1,052.00	<a href="http://www.chromtech.net.au/pdf2/vici09_102.pdf">http://www.chromtech.net.au/pdf2/vici09_102.pdf</a>
2NI4WE.5	.5ul Manual w/standoff 1/32in fittings x0.25mm ports	Valco Valves	\$988.00	<a href="http://www.chromtech.net.au/pdf2/vici09_102.pdf">http://www.chromtech.net.au/pdf2/vici09_102.pdf</a>
2NI4WE.2	.2ul Manual w/standoff 1/32in fittings x0.25mm ports	Valco Valves	\$988.00	<a href="http://www.chromtech.net.au/pdf2/vici09_102.pdf">http://www.chromtech.net.au/pdf2/vici09_102.pdf</a>
2NI4WE.1	.1ul Manual w/standoff 1/32in fittings x0.25mm ports	Valco Valves	\$988.00	<a href="http://www.chromtech.net.au/pdf2/vici09_102.pdf">http://www.chromtech.net.au/pdf2/vici09_102.pdf</a>
2NI4WE.06	.06ul Manual w/standoff 1/32in fittings x0.25mm ports	Valco Valves	\$988.00	<a href="http://www.chromtech.net.au/pdf2/vici09_102.pdf">http://www.chromtech.net.au/pdf2/vici09_102.pdf</a>
A2NI4WE.1	.1ul w/air actuator 1/32in fittings x0.25mm ports	Valco Valves	\$1,192.00	<a href="http://www.chromtech.net.au/pdf2/vici09_102.pdf">http://www.chromtech.net.au/pdf2/vici09_102.pdf</a>
A2NI4WE.2	.2ul w/air actuator 1/32in fittings x0.25mm ports	Valco Valves	\$1,192.00	<a href="http://www.chromtech.net.au/pdf2/vici09_102.pdf">http://www.chromtech.net.au/pdf2/vici09_102.pdf</a>
A2NI4WE.5	.5ulw/air actuator 1/32in fittings x0.25mm ports	Valco Valves	\$1,192.00	<a href="http://www.chromtech.net.au/pdf2/vici09_102.pdf">http://www.chromtech.net.au/pdf2/vici09_102.pdf</a>
E2NI4WE.06	.06ul w/std electric actuator 1/32in fittings x0.25mm ports	Valco Valves	\$1,600.00	<a href="http://www.chromtech.net.au/pdf2/vici09_102.pdf">http://www.chromtech.net.au/pdf2/vici09_102.pdf</a>
E2NI4WE.1	.1ul w/std electric actuator 1/32in fittings x0.25mm ports	Valco Valves	\$1,600.00	<a href="http://www.chromtech.net.au/pdf2/vici09_102.pdf">http://www.chromtech.net.au/pdf2/vici09_102.pdf</a>
E2NI4WE.2	.2ul w/std electric actuator 1/32in fittings x0.25mm ports	Valco Valves	\$1,600.00	<a href="http://www.chromtech.net.au/pdf2/vici09_102.pdf">http://www.chromtech.net.au/pdf2/vici09_102.pdf</a>
EP2CI4WE.06	.06ul w/std electric actuator 1/32in fittings x0.25mm ports	Valco Valves	\$1,677.00	<a href="http://www.chromtech.net.au/pdf2/vici09_102.pdf">http://www.chromtech.net.au/pdf2/vici09_102.pdf</a>
A2CI4WE.1	.1ul w/air actuator 1/32in fittings x0.25mm ports	Valco Valves	\$1,052.00	<a href="http://www.chromtech.net.au/pdf2/vici09_102.pdf">http://www.chromtech.net.au/pdf2/vici09_102.pdf</a>
DCI4WE.2	.2ul Replacement valve 1/32in fittings x0.25mm ports	Valco Valves	\$734.00	<a href="http://www.chromtech.net.au/pdf2/vici09_102.pdf">http://www.chromtech.net.au/pdf2/vici09_102.pdf</a>
A2CI4WE.5	.5ul w/air actuator 1/32in fittings x0.25mm ports	Valco Valves	\$1,052.00	<a href="http://www.chromtech.net.au/pdf2/vici09_102.pdf">http://www.chromtech.net.au/pdf2/vici09_102.pdf</a>
E2CI4WE.06	.06ul w/std electric actuator 1/32in fittings x0.25mm ports	Valco Valves	\$1,460.00	<a href="http://www.chromtech.net.au/pdf2/vici09_102.pdf">http://www.chromtech.net.au/pdf2/vici09_102.pdf</a>
E2CI4WE.1	.1ul w/std electric actuator 1/32in fittings x0.25mm ports	Valco Valves	\$1,460.00	<a href="http://www.chromtech.net.au/pdf2/vici09_102.pdf">http://www.chromtech.net.au/pdf2/vici09_102.pdf</a>
EP2NI4WE.06	.06ul w/std electric actuator 1/32in fittings x0.25mm ports	Valco Valves	\$1,944.00	<a href="http://www.chromtech.net.au/pdf2/vici09_102.pdf">http://www.chromtech.net.au/pdf2/vici09_102.pdf</a>
E2CI4WE.5	.5ul w/std electric actuator 1/32in fittings x0.25mm ports	Valco Valves	\$1,460.00	<a href="http://www.chromtech.net.au/pdf2/vici09_102.pdf">http://www.chromtech.net.au/pdf2/vici09_102.pdf</a>

## VICI Valco Cheminert - Catalog 60 - AUD Prices at USD/AUD=1.00 - Ask for QUOTATION

E2NI4WE.5	.5ul w/std electric actuator 1/32in fittings x0.25mm ports	Valco Valves	\$1,600.00	<a href="http://www.chromtech.net.au/pdf2/vici09_102.pdf">http://www.chromtech.net.au/pdf2/vici09_102.pdf</a>
EP2CI4WE.1	.1ul w/std electric actuator 1/32in fittings x0.25mm ports	Valco Valves	\$1,677.00	<a href="http://www.chromtech.net.au/pdf2/vici09_102.pdf">http://www.chromtech.net.au/pdf2/vici09_102.pdf</a>
EP2CI4WE.2	.2ul w/std electric actuator 1/32in fittings x0.25mm ports	Valco Valves	\$1,677.00	<a href="http://www.chromtech.net.au/pdf2/vici09_102.pdf">http://www.chromtech.net.au/pdf2/vici09_102.pdf</a>
EP2CI4WE.5	.5ul w/std electric actuator 1/32in fittings x0.25mm ports	Valco Valves	\$1,677.00	<a href="http://www.chromtech.net.au/pdf2/vici09_102.pdf">http://www.chromtech.net.au/pdf2/vici09_102.pdf</a>
DCI4WE.06	.06ul Replacement valve 1/32in fittings x0.25mm ports	Valco Valves	\$734.00	<a href="http://www.chromtech.net.au/pdf2/vici09_102.pdf">http://www.chromtech.net.au/pdf2/vici09_102.pdf</a>
DCI4WE.1	.1ul Replacement valve 1/32in fittings x0.25mm ports	Valco Valves	\$734.00	<a href="http://www.chromtech.net.au/pdf2/vici09_102.pdf">http://www.chromtech.net.au/pdf2/vici09_102.pdf</a>
A2CI4WE.06	.06ul w/air actuator 1/32in fittings x0.25mm ports	Valco Valves	\$1,052.00	<a href="http://www.chromtech.net.au/pdf2/vici09_102.pdf">http://www.chromtech.net.au/pdf2/vici09_102.pdf</a>
EP2NI4WE.2	.2ul w/std electric actuator 1/32in fittings x0.25mm ports	Valco Valves	\$1,944.00	<a href="http://www.chromtech.net.au/pdf2/vici09_102.pdf">http://www.chromtech.net.au/pdf2/vici09_102.pdf</a>
ZSN1	1/16", short .30"	Valco Fittings	\$4.00	<a href="http://www.chromtech.net.au/pdf2/vici09_10.pdf">http://www.chromtech.net.au/pdf2/vici09_10.pdf</a>
ZN1Q	1/16", standard .43"	Valco Fittings	\$4.00	<a href="http://www.chromtech.net.au/pdf2/vici09_10.pdf">http://www.chromtech.net.au/pdf2/vici09_10.pdf</a>
LZN4-10	Stainless steel long nut 1/4in x1.110in	Valco Fittings	\$51.00	<a href="http://www.chromtech.net.au/pdf2/vici09_10.pdf">http://www.chromtech.net.au/pdf2/vici09_10.pdf</a>
ZN4-10	Stainless steel 1/4in nut	Valco Fittings	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_10.pdf">http://www.chromtech.net.au/pdf2/vici09_10.pdf</a>
XLZN2-10	Stainless steel 1/8in x1.075in long nut	Valco Fittings	\$42.00	<a href="http://www.chromtech.net.au/pdf2/vici09_10.pdf">http://www.chromtech.net.au/pdf2/vici09_10.pdf</a>
LZN.5-10	Stainless steel 1/32in x.450in long nut	Valco Fittings	\$45.00	<a href="http://www.chromtech.net.au/pdf2/vici09_10.pdf">http://www.chromtech.net.au/pdf2/vici09_10.pdf</a>
ZN2-10	Stainless steel 1/8in nut	Valco Fittings	\$27.00	<a href="http://www.chromtech.net.au/pdf2/vici09_10.pdf">http://www.chromtech.net.au/pdf2/vici09_10.pdf</a>
XLZN1-10	1/16" nut 1.00"	Valco Fittings	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_10.pdf">http://www.chromtech.net.au/pdf2/vici09_10.pdf</a>
LZN1-10	1/16" nut .75"	Valco Fittings	\$34.00	<a href="http://www.chromtech.net.au/pdf2/vici09_10.pdf">http://www.chromtech.net.au/pdf2/vici09_10.pdf</a>
MZN1-10	Stainless steel 1/16in x.50in long nut	Valco Fittings	\$32.00	<a href="http://www.chromtech.net.au/pdf2/vici09_10.pdf">http://www.chromtech.net.au/pdf2/vici09_10.pdf</a>
ZN1-10	Stainless steel 1/16in x.43in long nut	Valco Fittings	\$21.00	<a href="http://www.chromtech.net.au/pdf2/vici09_10.pdf">http://www.chromtech.net.au/pdf2/vici09_10.pdf</a>
LZN2-10	Stainless steel 1/8in x.825in long nut	Valco Fittings	\$34.00	<a href="http://www.chromtech.net.au/pdf2/vici09_10.pdf">http://www.chromtech.net.au/pdf2/vici09_10.pdf</a>
ZN.5-10	Stainless steel 1/32in nut	Valco Fittings	\$39.00	<a href="http://www.chromtech.net.au/pdf2/vici09_10.pdf">http://www.chromtech.net.au/pdf2/vici09_10.pdf</a>
SSACSD16MWE	Replacement rotor 16 position 1/16in fittings 0.75mm ports	Valco Valves	\$156.00	<a href="http://www.chromtech.net.au/pdf2/vici09_132.pdf">http://www.chromtech.net.au/pdf2/vici09_132.pdf</a>
EZRF1.5PK-5	1/16"to1/32" PEEK/Glass,5pk	Reducing Ferrules	\$30.00	<a href="http://www.chromtech.net.au/pdf2/vici09_15.pdf">http://www.chromtech.net.au/pdf2/vici09_15.pdf</a>