

## Nuts

### Internal nuts – stainless steel

Nuts with product numbers starting with Z are for use with all standard Valco internal fittings and most valves. They may be used with fittings from other manufacturers as well. The L (long) and XL (extra-long) types are for situations where the fitting head may be otherwise inaccessible or where interference between fittings exists, as on many Valco multiposition valves. Standard material is 300 series stainless.



Stainless nuts			
Package of 10:	Length	Prod No	Price
1/32" nut	.30"	ZN.5-10	
1/32" nut	.45"	LZN.5-10	
1/16" nut	.43"	ZN1-10	
1/16" nut	.50"	MZN1-10	
1/16" nut	.75"	LZN1-10	
1/16" nut	1.00"	XLZN1-10	
1/8" nut	.57"	ZN2-10	
1/8" nut	.82"	LZN2-10	
1/8" nut	1.07"	XLZN2-10	
1/4" nut	.70"	ZN4-10	
1/4" nut	1.11"	LZN4-10	



### NEW Specialty nuts – stainless steel

These special purpose nuts facilitate a tight bend as the tube exits the fitting, and can also help prevent kinks in very thin wall tubing. Quick bend nuts are available in standard length (.43") and in a short version (.30") for certain custom applications. Note that the short version (ZSN1) can *only* be used in certain applications. Call for more information.

Stainless nuts			
Description	Length	Prod No	Price
1/16", standard	.43"	ZN1Q	
1/16", short	.30"	ZSN1	



#### MORE INFORMATION

PEEK nuts . . . . .	page 63
HPLC column end fittings . . . . .	43-46
Reducing unions	
Internal . . . . .	29
External . . . . .	30
External/internal . . . . .	31
Internal/external . . . . .	31
Unions	
Internal . . . . .	26
External . . . . .	27
External/internal . . . . .	27

#### TECH TIP

Fittings for **360 micron** tubing are available on pages 57-58.

0.25 mm =	.010"
0.50 mm =	.020"
0.75 mm =	.030"
1.0 mm =	.040"
1.5 mm =	.060"
2.0 mm =	.080"
4.6 mm =	.180"
6.0 mm =	.236"
6.4 mm =	.253"
7.0 mm =	.275"
10.0 mm =	.400"
27.0 mm =	1.08"
1/32" =	0.8 mm
1/16" =	1.6 mm
1/8" =	3.2 mm
1/4" =	6.4 mm
3/8" =	9.5 mm
1/2" =	12.7 mm