

**Capillary GC**

**Sampling and switching valves, 1/32" fittings, 0.25 mm ports (.010")**

*W Type*

**Med temp**  
**1/32" 0.25 mm**

Includes 4" standoff.  
Manual version not available without standoff.

Standard electric actuator:  
110 VAC for USA  
110/230 VAC to 24 VDC power supply for international  
Microelectric actuator:  
24 VDC, with 110/230 VAC to 24 VDC power supply

Sample loops are not included with valves. Order separately.

**SPECS**  
**400 psi gas**  
**225°C max**  
Nitronic 60 valve body  
Valcon E rotor  
  
For 300 psi, 350°C max, see facing page.



**4 Ports**

Prod No Price



**6 Ports**

Prod No Price



**8 Ports**

Prod No Price



**10 Ports**

Prod No Price

Manual with standoff  
With air actuator

4N4WE  
A4N4WE

4N6WE  
A4N6WE

4N8WE  
A4N8WE

4N10WE  
A4N10WE

With standard electric actuator  
With microelectric actuator

E4N4WE  
EH4N4WE

E4N6WE  
EH4N6WE

E4N8WE  
EH4N8WE

E4N10WE  
EH4N10WE

Replacement valve  
Replacement rotor

DN4WE  
SSAN4WE


DN6WE  
SSAN6WE

DN8WE  
SSAN8WE

DN10WE  
SSAN10WE

**OPTIONS**

- 3 and 12 port valves available
- 2", 3", and 6" standoffs
- Materials: Hastelloy C, Inconel 600, Monel 400, Nickel 200, Nitronic 50, Titanium, Zirconium (see pages 254-255)



**W Type**  
**1/32" fittings**

**1/32" Stainless steel loops for W Type valves**

Each stainless steel loop includes two stainless nuts and two stainless ferrules. Order special fittings separately.



Volume	Prod No	Price	Volume	Prod No	Price
2 µl	SL2NW		25 µl	SL25NW	
5 µl	SL5NW		50 µl	SL50NW	
10 µl	SL10NW		100 µl	SL100NW	
15 µl	SL15NW		250 µl	SL250NW	
20 µl	SL20NW		500 µl	SL500NW	

**MORE INFORMATION**

- Actuators
- Air ..... page 195
  - Manual ..... 204
  - Microelectric .. 188-189
  - Standard elec. .... 193
- Materials
- Metals ..... 254-255
  - Polymers ..... 256
  - Valve rotors ..... 257
- Standoff assemblies ..... 205

**ABOUT LOOPS**

- Other materials available in many sizes: Electroformed Nickel, Nickel 200, PEEK, and PTFE