

Sampling and switching valves, 1/32" fittings, 0.25 mm ports (.010")

W Type

SPECS

300 psi gas
350°C max

Nitronic 60 valve body
Valcon T rotor

For 400 psi, 225°C max,
see facing page

Includes 4" standoff.
Manual version not
available without
standoff.

Standard electric actuator:
110 VAC for USA
110/230 VAC to 24 VDC power supply for
international
Microelectric actuator:
24 VDC, with 110/230 VAC to 24 VDC
power supply

Sample loops are not
included with valves.
Order separately.

High temp

1/32" 0.25 mm

OPTIONS

- 3 and 12 port valves available
- 2", 3", and 6" standoffs
- Materials: Hastelloy C, Inconel 600, Monel 400, Nickel 200, Nitronic 50, Titanium, Zirconium (see pages 254-255)



4 Ports

Prod No Price



6 Ports

Prod No Price



8 Ports

Prod No Price



10 Ports

Prod No Price

Manual with standoff
With air actuator

4N4WT \$725
A4N4WT 885

4N6WT \$780
A4N6WT 940

4N8WT \$835
A4N8WT 995

4N10WT \$835
A4N10WT 995

With standard electric actuator
With microelectric actuator

E4N4WT 1205
EH4N4WT 1375

E4N6WT 1260
EH4N6WT 1430

E4N8WT 1315
EH4N8WT 1485

E4N10WT 1315
EH4N10WT 1485


Replacement valve
Replacement rotor

DN4WT 635
SSAN4WT 75

DN6WT 690
SSAN6WT 75

DN8WT 745
SSAN8WT 75

DN10WT 745
SSAN10WT 75



W Type
1/32" fittings

1/32" Stainless steel loops

for W Type valves

Each stainless steel loop includes two stainless
nuts and two stainless ferrules. Order special
fittings separately.



Volume	Prod No	Price	Volume	Prod No	Price
2 µl	SL2NW	\$25.00	25 µl	SL25NW	\$25.00
5 µl	SL5NW	25.00	50 µl	SL50NW	27.50
10 µl	SL10NW	25.00	100 µl	SL100NW	27.50
15 µl	SL15NW	25.00	250 µl	SL250NW	31.25
20 µl	SL20NW	25.00	500 µl	SL500NW	37.50

ABOUT LOOPS

- Other materials available in many sizes: Electroformed Nickel, Nickel 200, PEEK, and PTFE