

Sampling and switching valves, 1/16" fittings, 0.75 mm ports (.030")

UW Type

SPECS
300 psi gas
330°C max

Nitronic 60 valve body
Valcon T rotor

For 400 psi, 225°C max,
see page 107.

Includes 4" standoff

Standard electric actuator:
110 VAC for USA
110/230 VAC to 24 VDC power supply for
international
Microelectric actuator:
24 VDC, with 110/230 VAC to 24 VDC
power supply

Sample loops are not
included with valves.
Order separately.

High temp

1/16" 0.75 mm

OPTIONS

- 3, 12 and 14 port valves available
- 2", 3", and 6" standoffs
- Materials: Hastelloy C, Inconel 600, Monel 400, Nickel 200, Nitronic 50, Titanium, Zirconium (see pages 254-255)
- Larger bore available



4 Ports

Prod No Price



6 Ports

Prod No Price



8 Ports

Prod No Price



10 Ports

Prod No Price

Manual with standoff
With air actuator

4C4UWT
A4C4UWT

4C6UWT
A4C6UWT

4C8UWT
A4C8UWT

4C10UWT \$670
A4C10UWT

With standard electric actuator
With microelectric actuator

E4C4UWT
ED4C4UWT

E4C6UWT
ED4C6UWT

E4C8UWT
ED4C8UWT

E4C10UWT
ED4C10UWT


Replacement valve
Replacement rotor

DC4UWT
SSAC4UWT

DC6UWT
SSAC6UWT

DC8UWT
SSAC8UWT

DC10UWT
SSAC10UWT



UW Type
1/16" fittings

1/16" Stainless steel loops

for UW Type valves

Each stainless steel loop includes two stainless nuts and two stainless ferrules. Order special fittings separately.



Volume	Prod No	Price	Volume	Prod No	Price
5 µl	SL5CUW		100 µl	SL100CUW	
10 µl	SL10CUW		250 µl	SL250CUW	
15 µl	SL15CUW		500 µl	SL500CUW	
20 µl	SL20CUW		1 ml	SL1KCUW	
25 µl	SL25CUW		2 ml	SL2KCUW	
50 µl	SL50CUW		5 ml	SL5KCUW	
			10 ml	SL10KCUW	

ABOUT LOOPS

- Other materials available in many sizes: Electroformed Nickel, Hastelloy C, Nickel 200, PEEK, PTFE, and Titanium
- Loops > 2 ml are made from 1/8" OD tubing with brazed or welded 1/16" tube ends or reducing unions.