

## Semi-Preparative HPLC

### Injectors and switching valves, 1/16" fittings, 0.75 mm ports (.030")

*UW Type*

**5,000 psi**

**Semi-prep**

**1/16" 0.75 mm**

Standard electric actuator:  
110 VAC for USA  
110/230 VAC to 24 VDC power supply for international

Microelectric actuator:  
24 VDC, with 110/230 VAC to 24 VDC power supply

Sample loops are not included with valves. Order separately.


#### SPECS

**5000 psi liq**  
**75°C max**  
Nitronic 60 valve body  
Valcon H rotor

	4 Ports		6 Ports		8 Ports		10 Ports	
	Prod No	Price	Prod No	Price	Prod No	Price	Prod No	Price
Manual	C4UW		C6UW		C8UW		C10UW	
With air actuator	AC4UW		AC6UW		AC8UW		AC10UW	
With standard electric actuator	EC4UW		EC6UW		EC8UW		EC10UW	
With microelectric actuator	EDC4UW		EDC6UW		EDC8UW		EDC10UW	
Replacement valve	DC4UW		DC6UW		DC8UW		DC10UW	
Replacement rotor	SSAC4UW		SSAC6UW		SSAC8UW		SSAC10UW	

#### OPTIONS

- 3, 12, and 14 port valves available
- 2", 3", 4", and 6" standoffs
- 1/32" and 1/16" versions available with 0.25 mm (.010") bore
- Materials: Hastelloy C, Inconel 600, Monel 400, Nickel 200, Nitronic 50, Titanium, Zirconium (see pages 254-255)
- Larger bore available.



**UW Type**  
**1/16" fittings**

### 1/16" Stainless steel loops *for UW Type valves*

Each stainless steel loop includes two stainless nuts and two stainless ferrules. Order special fittings separately.



Volume	Prod No	Price	Volume	Prod No	Price
3 µl	SL3CUW		100 µl	SL100CUW	
5 µl	SL5CUW		250 µl	SL250CUW	
10 µl	SL10CUW		500 µl	SL500CUW	
15 µl	SL15CUW		1 ml	SL1KCUW	
20 µl	SL20CUW		2 ml	SL2KCUW	
25 µl	SL25CUW		5 ml	SL5KCUW	
50 µl	SL50CUW		10 ml	SL10KCUW	

#### ABOUT LOOPS

- Other materials available in many sizes: Electroformed Nickel, Hastelloy C, Nickel 200, PEEK, PTFE, and Titanium
- Loops > 2 ml are made from 1/8" OD tubing with brazed or welded 1/16" tube ends or reducing unions.

**CHROMalytic TECHNOlogy Pty Ltd AUSTRALIAN Distributors e-mail: sales@chromtech.net.au Tel: 03 9762 2034**