

## Preparative HPLC

### Injectors and switching valves, 1/8" fittings, large bore

*UW Type*

**5,000 psi**

**Prep**

**1/8" Large bore**

Manual 10 port includes 2" standoff.

Standard electric actuator:  
110 VAC for USA  
110/230 VAC to 24 VDC power supply for international.

Microelectric actuator:  
24 VDC, with 110/230 VAC to 24 VDC power supply.

Sample loops are not included with valves. Order separately.

**SPECS**  
**5000 psi liq**  
**75°C max**  
Nitronic 60 valve body  
Valcon H rotor



**4 Ports**  
**1.7 mm (.067")**  
Prod No Price



**6 Ports**  
**1.7 mm (.067")**  
Prod No Price




**8 Ports**  
**1.3 mm (.050")**  
Prod No Price



**10 Ports**  
**1.0 mm (.040")**  
Prod No Price

Manual	L4UW	L6UW	L8UW	2L10UW
With air actuator	AL4UW	AL6UW	AL8UW	AL10UW
With standard electric actuator	EL4UW	EL6UW	EL8UW	EL10UW
With microelectric actuator	EDL4UW	EDL6UW	EDL8UW	EDL10UW
Replacement valve	DL4UW	DL6UW	DL8UW	DL10UW
Replacement rotor	SSAL4UW	SSAL6UW	SSAL8UW	SSAL10UW

- OPTIONS**
- 3 port valve available
  - 2", 3", 4", and 6" standoffs
  - Materials: Hastelloy C, Inconel 600, Monel 400, Nickel 200, Nitronic 50, Titanium, Zirconium (see pages 254-255)
  - Smaller bore available. (see page 115)



**UW Type**  
**1/8" fittings**

### MORE INFORMATION

Actuators  
Air ..... page 195  
Manual ..... 204  
Microelectric .. 188-189  
Standard elec. .... 193

Materials  
Metals ..... 254-255  
Polymers ..... 256  
Valve rotors. .... 257

Standoff assemblies ..... 205

### 1/8" Stainless steel loops *for UW Type valves*

Each stainless steel loop includes two stainless nuts and two stainless ferrules. Order special fittings separately.

Volume	Prod No	Price	Volume	Prod No	Price
100 µl	SL100UW		2 ml	SL2KUW	
250 µl	SL250UW		5 ml	SL5KUW	
500 µl	SL500UW		10 ml	SL10KUW	
1 ml	SL1KUW		20 ml	SL20KUW	



### ABOUT LOOPS

- Other materials available in many sizes: Electroformed Nickel, Hastelloy C, Nickel 200, PEEK, PTFE, and Titanium
- Loops < 100 µl are made from 1/16" OD tubing with brazed or welded 1/8" tube ends..

**CHROMalytic TECHNOlogy Pty Ltd AUSTRALIAN Distributors e-mail: sales@chromtech.net.au Tel: 03 9762 2034**