

1/8" fittings, 1.0 mm ports (.040")

MWType

**SPECS**

4-8 Positions:  
**400 psi gas**  
**200°C max**  
 10-16 Positions:  
**200 psi gas**  
**200°C max**  
 Nitronic 60 body  
 Valcon E rotor

Includes 2" standoff.  
 Ask about closemount  
 assembly if valve will  
 not be heated.

Standard electric actuators:  
 110 VAC for USA  
 110/230 VAC to 24 VDC power supply for international  
 Microelectric actuators:  
 24 VDC (includes a 110/230 VAC to 24 VDC power supply)

**OPTIONS**

- 4 and 8 positions available
- 3", 4", and 6" standoffs
- Materials: Hastelloy C, Inconel 600, Monel 400, Nickel 200, Nitronic 50, Titanium (see pages 254-255)
- Larger bore available

Low pressure

SD  
 Dead-end

1/8"

1.0 mm

	6 Position		10 Position		12 Position		16 Position	
	Prod No	Price	Prod No	Price	Prod No	Price	Prod No	Price
Manual (not recommended)	2SD6MWE		2SD10MWE		2SD12MWE		2SD16MWE	
With air actuator	A2SD6MWE		A2SD10MWE		A2SD12MWE		A2SD16MWE	
With standard electric actuator	E2SD6MWE		E2SD10MWE		E2SD12MWE		E2SD16MWE	
With microelectric actuator	EMT2SD6MWE		EMT2SD10MWE		EMT2SD12MWE		EMT2SD16MWE	
Replacement valve	DSD6MWE		DSD10MWE		DSD12MWE		DSD16MWE	
Replacement rotor	SSASD6MWE		SSASD10MWE		SSASD12MWE		SSASD16MWE	

1/4" fittings, 4.0 mm ports (.156")

MWType

**SPECS**

**100 psi gas**  
**75°C max**  
 Nitronic 60 body  
 Valcon E2 rotor

Includes 2" standoff.  
 Ask about closemount  
 assembly if valve  
 will not be heated.

Manual version not  
 available.

Standard electric actuators:  
 110 VAC for USA  
 110/230 VAC to 24 VDC power supply for  
 international  
 Microelectric actuators:  
 24 VDC (includes a 110/230 VAC to 24 VDC  
 power supply)

**OPTIONS**

- 3", 4", and 6" standoffs
- Materials: Hastelloy C, Inconel 600, Monel 400, Nickel 200, Nitronic 50, Titanium (see pages 254-255)

Low pressure

SD  
 Dead-end

1/4"

4.0 mm

	4 Position		6 Position		8 Position		10 Position	
	Prod No	Price	Prod No	Price	Prod No	Price	Prod No	Price
With air actuator	AH2VLS4MWE2		AH2VLS6MWE2		AH2VLS8MWE2		AH2VLS10MWE2	
With std electric actuator	E2VLS4MWE2		E2VLS6MWE2		E2VLS8MWE2		E2VLS10MWE2	
With microelectric actuator	EMT2VLS4MWE2		EMT2VLS6MWE2		EMT2VLS8MWE2		EMT2VLS10MWE2	
Replacement valve	DVLS4MWE2		DVLS6MWE2		DVLS8MWE2		DVLS10MWE2	
Replacement rotor	SSAVLS4MWE2		SSAVLS6MWE2		SSAVLS8MWE2		SSAVLS10MWE2	

