


## Selectors – Low Pressure

### Trapping flowpath – ST configuration

ST selectors are used for multi-column, multi-sample, or multi-trap operations, and are available for use with 4 to 16 loops, or positions. For an application suggestion, see page 137.



### 1/16" fittings, 0.75 mm ports (.030")


MW Type

<b>Low pressure</b>	Includes 2" standoff. Ask about closemount assembly if valve will not be heated.	Standard electric actuators: 110 VAC for USA 110/230 VAC to 24 VDC power supply for international	<b>SPECS</b> 200 psi gas 200°C max Nitronic 60 body Valcon E rotor
<b>ST Trapping</b>		Microelectric actuators: 24 VDC (includes a 110/230 VAC to 24 VDC power supply)	
<b>1/16" 0.75 mm</b>			

#### OPTIONS

- 4 and 8 positions available
- 3", 4", and 6" standoffs
- Materials: Hastelloy C, Inconel 600, Monel 400, Nickel 200, Nitronic 50, Titanium (see pages 254-255)

	6 Position		10 Position		12 Position		16 Position	
	Prod No	Price	Prod No	Price	Prod No	Price	Prod No	Price
Manual (not recommended)	2CST6MWE		2CST10MWE		2CST12MWE		2CST16MWE	
With air actuator	A2CST6MWE		A2CST10MWE		A2CST12MWE		A2CST16MWE	
With standard electric actuator	E2CST6MWE		E2CST10MWE		E2CST12MWE		E2CST16MWE	
With microelectric actuator	EMT2CST6MWE		EMT2CST10MWE		EMT2CST12MWE		EMT2CST16MWE	
Replacement valve	DCST6MWE		DCST10MWE		DCST12MWE		DCST16MWE	
Replacement rotor	SSACST6MWE		SSACST10MWE		SSACST12MWE		SSACST16MWE	



### 1/16" Stainless steel loops for MW Type valves

Each stainless steel loop includes two stainless nuts and two stainless ferrules. Order special fittings separately. When a set of loops is ordered, loops will be supplied from the same lot.

Volume	Prod No	Price	Volume	Prod No	Price
50 µl	SL50CSTP		1 ml	SL1KCSTP	
100 µl	SL100CSTP		2 ml	SL2KCSTP	
250 µl	SL250CSTP		5 ml	SL5KCSTP	
500 µl	SL500CSTP		10 ml	SL10KCSTP	



#### MORE INFORMATION

Application.....	page 137
Actuators	
Air .....	194
Microelectric ..	190-191
Standard elec.....	193
Materials	
Metals.....	254-255
Polymers .....	256
Valve rotors.....	257
Mounting hardware	
Closemount .....	208
Standoff.....	205

CHROMalytic TECHNOlogy Pty Ltd AUSTRALIAN Distributors e-mail: sales@chromtech.net.au Tel: 03 9762 2034