



**Polymeric ferrules**

**MORE INFORMATION**  
 PEEK ferrules . . . . . page 63  
 Grooved PEEK  
 ferrules . . . . . 63

For more detailed  
 information on polymers,  
 refer to the discussion  
 on page 256.

**POLYMERS  
 AT A GLANCE**

CTFE . . . . . KF  
*Resists all inorganic  
 corrosives.  
 Produced as Kel-F®*

FEP . . . . . FEP  
*Chemical resistance  
 equals PTFE, but lower  
 creep and higher  
 friction*

PTFE, Glass-filled . . . . . TFG  
*Inert, mechanically  
 stable*

PTFE, Virgin . . . . . TF  
*Inert; very soft, easily  
 cold flows.  
 Produced as Teflon®*

Polyimide, Graphite . . . . . GV  
*Soft, easy to form fer-  
 rules*

Polyimide, Valcon . . . . . V  
*High temp, graphite  
 reinforced*

Polyimide, Virgin . . . . . V1  
*High temp, electrical  
 insulator*

	Prod No	Price	Prod No	Price	Prod No	Price
<i>Package of 10:</i>	<b>PTFE, Virgin</b>		<b>PTFE, Glass-filled</b>		<b>FEP</b>	
1/32"	ZF.5TF-10		ZF.5TFG-10		ZF.5FEP-10	
1/16"	ZF1TF-10		ZF1TFG-10		ZF1FEP-10	
1/8"	ZF2TF-10		ZF2TFG-10		ZF2FEP-10	
1/4"	ZF4TF-10		ZF4TFG-10		ZF4FEP-10	
3/8"	ZF6TF-10		ZF6TFG-10		ZF6FEP-10	
1/2"	ZF8TF-10		ZF8TFG-10		ZF8FEP-10	

	Prod No	Price	Prod No	Price
<i>Package of 10:</i>	<b>PFA</b>		<b>CTFE</b>	
1/32"	ZF.5PFA-10		ZF.5KF-10	
1/16"	ZF1PFA-10		ZF1KF-10	
1/8"	ZF2PFA-10		ZF2KF-10	
1/4"	ZF4PFA-10		ZF4KF-10	
3/8"	ZF6PFA-10		ZF6KF-10	
1/2"	ZF8PFA-10		ZF8KF-10	

	Prod No	Price	Prod No	Price	Prod No	Price
<i>Package of 5:</i>	<b>Polyimide, Graphite</b>		<b>Polyimide, Valcon</b>		<b>Polyimide, Virgin</b>	
1/32"	ZF.5GV-5		ZF.5V-5		ZF.5V1-5	
1/16"	ZF1GV-5		ZF1V-5		ZF1V1-5	
1/8"	ZF2GV-5		ZF2V-5		ZF2V1-5	
1/4"	ZF4GV-5		ZF4V-5		ZF4V1-5	
3/8"	ZF6GV-5		ZF6V-5		ZF6V1-5	
1/2"	ZF8GV-5		ZF8V-5		ZF8V1-5	

**FERRULE  
 IDENTIFICATION**

To differentiate among  
 the most commonly  
 ordered metal ferrules,  
 ring(s) are engraved on  
 the non-sealing surfaces.

